# LLNL Nanosecond Gated Camera 2.1.2

Generated by Doxygen 1.10.0

1	Namespace Index	1
	1.1 Package List	1
2	P Hierarchical Index	3
	2.1 Class Hierarchy	3
3	3 Class Index	5
	3.1 Class List	5
4	File Index	7
	4.1 File List	7
5	i Namespace Documentation	9
	5.1 nsCamera Namespace Reference	9
	5.1.1 Detailed Description	9
	5.1.2 Variable Documentation	10
	5.1.2.1 <u>all</u>	10
	5.2 nsCamera.boards Namespace Reference	10
	5.2.1 Detailed Description	
	5.2.2 Variable Documentation	10
	5.2.2.1 <u>all</u>	10
	5.3 nsCamera.boards.LLNL_v1 Namespace Reference	
	5.3.1 Detailed Description	
	5.4 nsCamera.boards.LLNL_v4 Namespace Reference	
	5.4.1 Detailed Description	
	5.5 nsCamera.CameraAssembler Namespace Reference	11
	5.5.1 Detailed Description	
	5.6 nsCamera.comms Namespace Reference	
	5.6.1 Detailed Description	
	5.6.2 Variable Documentation	
	5.6.2.1 <u>all</u>	
	5.7 nsCamera.comms.GigE Namespace Reference	
	5.7.1 Detailed Description	
	5.8 nsCamera.comms.RS422 Namespace Reference	
	5.8.1 Detailed Description	
	5.9 nsCamera.sensors Namespace Reference	
	5.9.1 Detailed Description	
	5.9.2 Variable Documentation	
	5.9.2.1all	
	5.10 nsCamera.sensors.daedalus Namespace Reference	
	5.10.1 Detailed Description	

5.11 nsCamera.sensors.icarus Namespace Reference
5.11.1 Detailed Description
5.12 nsCamera.sensors.icarus2 Namespace Reference
5.12.1 Detailed Description
5.13 nsCamera.sensors.sensorBase Namespace Reference
5.13.1 Detailed Description
5.14 nsCamera.utils Namespace Reference
5.14.1 Detailed Description
5.14.2 Variable Documentation
5.14.2.1all
5.15 nsCamera.utils.crc16pure Namespace Reference
5.15.1 Detailed Description
5.15.2 Function Documentation
5.15.2.1 _crc16()
5.15.2.2 crc16xmodem()
5.15.3 Variable Documentation
5.15.3.1 CRC16_XMODEM_TABLE
5.16 nsCamera.utils.FlatField Namespace Reference
5.16.1 Detailed Description
5.16.2 Function Documentation
5.16.2.1 generateFF()
5.16.2.2 getFilenames()
5.16.2.3 getROlvector()
5.16.2.4 removeFF()
5.16.2.5 removeFFall()
5.16.2.6 tslopes()
5.16.3 Variable Documentation
5.16.3.1 action
5.16.3.2 args
5.16.3.3 default
5.16.3.4 dest
5.16.3.5 directory
5.16.3.6 framelist
5.16.3.7 help
5.16.3.8 nargs
5.16.3.9 parser
5.17 nsCamera.utils.misc Namespace Reference
5.17.1 Detailed Description
5.17.2 Function Documentation

5.17.2.1 bytes2str()	24
5.17.2.2 checkCRC()	25
5.17.2.3 flattenlist()	25
5.17.2.4 generateFrames()	26
5.17.2.5 getEnter()	27
5.17.2.6 loadDumpedData()	27
5.17.2.7 makeLogLabels()	30
5.17.2.8 partition()	31
5.17.2.9 plotFrames()	32
5.17.2.10 saveTiffs()	33
5.17.2.11 str2bytes()	35
5.17.2.12 str2nparray()	36
5.18 nsCamera.utils.Packet Namespace Reference	36
5.18.1 Detailed Description	36
5.19 nsCamera.utils.Subregister Namespace Reference	37
5.19.1 Detailed Description	37
Class Documentation	39
·	
·	
v	
v	
6.1.3.12 configADCs()	51
6.1.3.13 disarm()	
6.1.3.15 dumpRegisters()	
6.1.3.16 dumpStatus()	
	5.17.2.2 checkCRC() 5.17.2.3 flattenlist() 5.17.2.4 generateFrames() 5.17.2.5 getEnter() 5.17.2.6 loadDumpedData() 5.17.2.6 loadDumpedData() 5.17.2.7 makeLogLabels() 5.17.2.8 partition() 5.17.2.9 plotFrames() 5.17.2.10 saveTiffs() 5.17.2.11 str2bytes() 5.17.2.12 str2nparray() 5.18 nsCamera.utils.Packet Namespace Reference 5.18.1 Detailed Description 5.19 nsCamera.utils.Subregister Namespace Reference 5.19.1 Detailed Description 6.1 nsCamera.CameraAssembler.CameraAssembler Class Reference 6.1.1 Detailed Description 6.1.2 Constructor & Destructor Documentation 6.1.2.1init() 6.1.3 Member Function Documentation 6.1.3.1 abortReadoff() 6.1.3.2 arm() 6.1.3.3 batchAcquire() 6.1.3.4 bytes2str() 6.1.3.5 checkCRC() 6.1.3.5 checkCRG() 6.1.3.7 checkSensorVoltStat() 6.1.3.9 checkStatus() 6.1.3.1 closeDevice() 6.1.3.11 closeDevice() 6.1.3.11 closeDevice() 6.1.3.11 closeDevice() 6.1.3.13 disarm() 6.1.3.13 disarm() 6.1.3.13 disarm() 6.1.3.14 dumpNumpy()

6.1.3.17 dumpSubregisters()	53
6.1.3.18 enableLED()	54
6.1.3.19 flattenlist()	54
6.1.3.20 getBoardInfo()	54
6.1.3.21 getEnter()	55
6.1.3.22 getManualTiming()	55
6.1.3.23 getMonV()	55
6.1.3.24 getPot()	56
6.1.3.25 getPotV()	58
6.1.3.26 getPressure()	58
6.1.3.27 getRegister()	59
6.1.3.28 getSensTemp()	60
6.1.3.29 getSubregister()	60
6.1.3.30 getTemp()	61
6.1.3.31 getTimer()	61
6.1.3.32 getTiming()	61
6.1.3.33 initBoard()	61
6.1.3.34 initialize()	62
6.1.3.35 initPots()	
6.1.3.36 initPowerCheck()	63
6.1.3.37 initSensor()	64
6.1.3.38 latchPots()	64
6.1.3.39 mmReadoff()	64
6.1.3.40 parseReadoff()	65
6.1.3.41 partition()	65
6.1.3.42 plotFrames()	65
6.1.3.43 powerCheck()	66
6.1.3.44 printBoardInfo()	66
6.1.3.45 readFrames()	67
6.1.3.46 readImgs()	67
6.1.3.47 readoff()	68
6.1.3.48 readSerial()	68
6.1.3.49 readSRAM()	68
6.1.3.50 reboot()	68
6.1.3.51 reinitialize()	69
6.1.3.52 reportEdgeDetects()	69
6.1.3.53 reportStatus()	69
6.1.3.54 resetTimer()	69
6.1.3.55 resolveSubreg()	70

	6.1.3.56 saveFrames()	70
	6.1.3.57 saveHDF()	71
	6.1.3.58 saveNumpys()	72
	6.1.3.59 saveTiffs()	74
	6.1.3.60 selectOscillator()	74
	6.1.3.61 sendCMD()	74
	6.1.3.62 sensorSpecific()	74
	6.1.3.63 setArbTiming()	75
	6.1.3.64 setExtClk()	75
	6.1.3.65 setFrames()	75
	6.1.3.66 setHighFullWell()	76
	6.1.3.67 setInterlacing()	76
	6.1.3.68 setLED()	77
	6.1.3.69 setManualShutters()	77
	6.1.3.70 setManualTiming()	77
	6.1.3.71 setPhiDelay()	77
	6.1.3.72 setPot()	78
	6.1.3.73 setPotV()	79
	6.1.3.74 setPowerSave()	82
	6.1.3.75 setPPER()	82
	6.1.3.76 setRegister()	83
	6.1.3.77 setRows()	84
	6.1.3.78 setSubregister()	85
	6.1.3.79 setTiming()	86
	6.1.3.80 setTriggerDelay()	87
	6.1.3.81 setZeroDeadTime()	87
	6.1.3.82 startCapture()	87
	6.1.3.83 str2bytes()	87
	6.1.3.84 str2nparray()	87
	6.1.3.85 submitMessages()	88
	6.1.3.86 waitForSRAM()	88
	6.1.3.87 writeSerial()	89
6.1.4 N	Member Data Documentation	89
	6.1.4.1 abort	89
	6.1.4.2 arch	89
	6.1.4.3 armed	89
	6.1.4.4 board	89
	6.1.4.5 boardname	89
	6.1.4.6 commname	90

6.1.4.7 comms	 90
6.1.4.8 currtime	
6.1.4.9 cycle	 90
6.1.4.10 FPGAboardtype	 90
6.1.4.11 FPGAinterfaces	 90
6.1.4.12 FPGAinvalid	 90
6.1.4.13 FPGANum	 91
6.1.4.14 FPGArad	 91
6.1.4.15 FPGAsensor	 91
6.1.4.16 FPGAVersion	 91
6.1.4.17 inittime	 91
6.1.4.18 iplist	 91
6.1.4.19 logcrit	 91
6.1.4.20 logcritbase	 92
6.1.4.21 logdebug	 92
6.1.4.22 logdebugbase	 92
6.1.4.23 logerr	 92
6.1.4.24 logerrbase	 92
6.1.4.25 loginfo	 92
6.1.4.26 loginfobase	 92
6.1.4.27 logtag	 93
6.1.4.28 logwarn	 93
6.1.4.29 logwarnbase	 93
6.1.4.30 oldtime	 93
6.1.4.31 packageroot	 93
6.1.4.32 padToFull	 93
6.1.4.33 parsedtime	 93
6.1.4.34 payloaderror	 94
6.1.4.35 platform	 94
6.1.4.36 port	 94
6.1.4.37 PY3	 94
6.1.4.38 pyth1	 94
6.1.4.39 pyth2	 94
6.1.4.40 python	 94
6.1.4.41 read	 95
6.1.4.42 savetime	 95
6.1.4.43 sensmanual	 95
6.1.4.44 sensor	 95
6.1.4.45 sensorname	 95

6.1.4.46 senstiming
6.1.4.47 timeout
6.1.4.48 trigtime
6.1.4.49 unstringed
6.1.4.50 verblevel
6.1.4.51 verbmap
6.1.4.52 verbose
6.1.4.53 version
6.1.4.54 waited
6.2 nsCamera.sensors.daedalus.daedalus Class Reference
6.2.1 Detailed Description
6.2.2 Constructor & Destructor Documentation
6.2.2.1init()
6.2.3 Member Function Documentation
6.2.3.1 getManualTiming()
6.2.3.2 getSensTemp()
6.2.3.3 parseReadoff()
6.2.3.4 reportStatusSensor()
6.2.3.5 selectOscillator()
6.2.3.6 sensorSpecific()
6.2.3.7 setExtClk()
6.2.3.8 setHighFullWell()
6.2.3.9 setInterlacing()
6.2.3.10 setManualShutters()
6.2.3.11 setManualTiming()
6.2.3.12 setPhiDelay()
6.2.3.13 setTriggerDelay()
6.2.3.14 setZeroDeadTime()
6.2.4 Member Data Documentation
6.2.4.1 bytesperpixel
6.2.4.2 ca
6.2.4.3 columns [1/2]
6.2.4.4 columns [2/2]
6.2.4.5 detect
6.2.4.6 firstframe
6.2.4.7 firstrow
6.2.4.8 fpganumID
6.2.4.9 height
6.2.4.10 HFW [1/2]

<b>6.2.4.11 HFW</b> [2/2]	115
<b>6.2.4.12 interlacing</b> [1/2]	115
<b>6.2.4.13 interlacing</b> [2/2]	115
6.2.4.14 lastframe	115
6.2.4.15 lastrow	115
6.2.4.16 loglabel	116
6.2.4.17 maxframe	116
6.2.4.18 maxheight	116
6.2.4.19 maxwidth	116
6.2.4.20 minframe	116
6.2.4.21 nframes	116
6.2.4.22 padToFull	116
6.2.4.23 sens_registers	117
6.2.4.24 sens_subregisters	117
6.2.4.25 sensfam	117
6.2.4.26 specwarn	117
6.2.4.27 toffset	117
6.2.4.28 tslope	117
6.2.4.29 width	117
6.2.4.30 ZDT [1/2]	118
6.2.4.31 ZDT [2/2]	118
6.3 nsCamera.utils.misc.fakeCA Class Reference	118
6.3.1 Detailed Description	119
6.3.2 Constructor & Destructor Documentation	119
6.3.2.1init()	119
6.3.3 Member Function Documentation	120
6.3.3.1 partition()	120
6.3.4 Member Data Documentation	120
6.3.4.1 boardname	120
6.3.4.2 logcrit	120
6.3.4.3 logcritbase	120
6.3.4.4 logdebug	120
6.3.4.5 logdebugbase	120
6.3.4.6 logerr	121
6.3.4.7 logerrbase	121
6.3.4.8 loginfo	121
6.3.4.9 loginfobase	121
6.3.4.10 logtag	121
6.3.4.11 logwarn	121

6.3.4.12 logwarnbase
6.3.4.13 padToFull
6.3.4.14 sensor
6.3.4.15 sensorname
6.4 nsCamera.comms.GigE.GigE Class Reference
6.4.1 Detailed Description
6.4.2 Constructor & Destructor Documentation
6.4.2.1init()
6.4.3 Member Function Documentation
6.4.3.1 arm()
6.4.3.2 closeDevice()
6.4.3.3 getCardInfo()
6.4.3.4 getCardIP()
6.4.3.5 openDevice()
6.4.3.6 readFrames()
6.4.3.7 readoff()
6.4.3.8 readSerial()
6.4.3.9 sendCMD()
6.4.3.10 writeSerial()
6.4.4 Member Data Documentation
6.4.4.1 _zest
6.4.4.2 ca
6.4.4.3 CardInfo
6.4.4.4 CardInfoP
6.4.4.5 closecard
6.4.4.6 Connection
6.4.4.7 dport
6.4.4.8 logcrit
6.4.4.9 logdebug
6.4.4.10 logerr
6.4.4.11 loginfo
6.4.4.12 logwarn
6.4.4.13 mode
6.4.4.14 payloadsize
6.4.4.15 readTimeout
6.4.4.16 skipError
6.4.4.17 writeTimeout
6.4.4.18 ZCountCards
6.4.4.19 ZErrorDict

6.4.4.20 ZESTETM1_CARD_INFO
6.4.4.21 ZOpenConnection
6.4.4.22 ZReadData
6.4.4.23 ZWriteData
6.5 nsCamera.sensors.icarus.icarus Class Reference
6.5.1 Detailed Description
6.5.2 Constructor & Destructor Documentation
6.5.2.1init()
6.5.3 Member Function Documentation
6.5.3.1 checkSensorVoltStat()
6.5.3.2 sensorSpecific()
6.5.4 Member Data Documentation
6.5.4.1 bytesperpixel
6.5.4.2 ca
6.5.4.3 columns
6.5.4.4 detect
6.5.4.5 firstframe
6.5.4.6 firstrow
6.5.4.7 fpganumID
6.5.4.8 height
6.5.4.9 icarustype
6.5.4.10 interlacing
6.5.4.11 lastframe
6.5.4.12 lastrow
6.5.4.13 loglabel
6.5.4.14 maxframe
6.5.4.15 maxheight
6.5.4.16 maxwidth
6.5.4.17 minframe
6.5.4.18 nframes
6.5.4.19 padToFull
6.5.4.20 sens_registers
6.5.4.21 sens_subregisters
6.5.4.22 sensfam
6.5.4.23 specwarn
6.5.4.24 width
6.6 nsCamera.sensors.icarus2.icarus2 Class Reference
6.6.1 Detailed Description
6.6.2 Constructor & Destructor Documentation

	6.6.2.1init()	16
	6.6.3 Member Function Documentation	17
	6.6.3.1 sensorSpecific()	١7
	6.6.4 Member Data Documentation	8
	6.6.4.1 bytesperpixel	8
	6.6.4.2 ca	8
	6.6.4.3 columns	8
	6.6.4.4 detect	8
	6.6.4.5 firstframe	8
	6.6.4.6 firstrow	١9
	6.6.4.7 fpganumID	19
	6.6.4.8 height	١9
	6.6.4.9 icarustype	19
	6.6.4.10 interlacing	١9
	6.6.4.11 lastframe	١9
	6.6.4.12 lastrow	١9
	6.6.4.13 loglabel	50
	6.6.4.14 maxframe	50
	6.6.4.15 maxheight	50
	6.6.4.16 maxwidth	50
	6.6.4.17 minframe	50
	6.6.4.18 nframes	50
	6.6.4.19 padToFull	50
	6.6.4.20 sens_registers	51
	6.6.4.21 sens_subregisters	51
	6.6.4.22 sensfam	51
	6.6.4.23 specwarn	51
	6.6.4.24 width	51
6.7 n	Camera.boards.LLNL_v1.llnl_v1 Class Reference	52
	6.7.1 Detailed Description	53
	6.7.2 Constructor & Destructor Documentation	53
	6.7.2.1init()	53
	5.7.3 Member Function Documentation	55
	6.7.3.1 checkStatus()	55
	6.7.3.2 checkStatus2()	55
	6.7.3.3 clearStatus()	6
	6.7.3.4 configADCs()	6
	6.7.3.5 disarm()	57
	6.7.3.6 dumpStatus()	57

6.8

6.7	3.7 enableLED()
6.7	3.8 getPressure()
6.7	3.9 getTemp()
6.7	3.10 getTimer()
6.7	3.11 initBoard()
6.7	3.12 initPots()
6.7	3.13 initSensor()
6.7	3.14 latchPots()
6.7	3.15 readSRAM()
6.7	3.16 reportEdgeDetects()
6.7	3.17 reportStatus()
6.7	3.18 resetTimer()
6.7	3.19 setLED()
6.7	3.20 setPowerSave()
6.7	3.21 setPPER()
6.7	3.22 softReboot()
6.7	3.23 startCapture()
6.7	3.24 waitForSRAM()
6.7.4 Memb	er Data Documentation
6.7	4.1 ADC5_bipolar
6.7	4.2 ADC5_mult
6.7	4.3 ca
6.7	4.4 daedalus_monitor_controls
6.7	4.5 daedalus_subreg_aliases
6.7	4.6 icarus_monitor_controls
6.7	4.7 icarus_subreg_aliases
6.7	4.8 logcrit
6.7	4.9 logdebug
6.7	4.10 logerr
6.7	4.11 loginfo
6.7	4.12 logwarn
6.7	4.13 registers
6.7	4.14 rs422_baud
6.7	4.15 rs422_cmd_wait
6.7	4.16 subregisters
6.7	4.17 subreglist
6.7	4.18 VREF
nsCamera.bo	ards.LLNL_v4.llnl_v4 Class Reference
6.8.1 Detail	ed Description

6.8.2 Constructor & Destructor Documentation
6.8.2.1init()
6.8.3 Member Function Documentation
6.8.3.1 checkStatus()
6.8.3.2 checkStatus2()
6.8.3.3 clearStatus()
6.8.3.4 configADCs()
6.8.3.5 disarm()
6.8.3.6 dumpStatus()
6.8.3.7 enableLED()
6.8.3.8 getPressure()
6.8.3.9 getTemp()
6.8.3.10 getTimer()
6.8.3.11 initBoard()
6.8.3.12 initPots()
6.8.3.13 initSensor()
6.8.3.14 latchPots()
6.8.3.15 readSRAM()
6.8.3.16 reportEdgeDetects()
6.8.3.17 reportStatus()
6.8.3.18 resetTimer()
6.8.3.19 setLED()
6.8.3.20 setPowerSave()
6.8.3.21 setPPER()
6.8.3.22 softReboot()
6.8.3.23 startCapture()
6.8.3.24 waitForSRAM()
6.8.4 Member Data Documentation
6.8.4.1 ADC5_bipolar
6.8.4.2 ADC5_mult
6.8.4.3 ca
6.8.4.4 daedalus_monitor_controls
6.8.4.5 daedalus_subreg_aliases
6.8.4.6 defoff
6.8.4.7 defsens
6.8.4.8 icarus_monitor_controls
6.8.4.9 icarus_subreg_aliases
6.8.4.10 logcrit
6.8.4.11 logdebug

6.8.4.12 logerr	190
6.8.4.13 loginfo	190
6.8.4.14 logwarn	190
6.8.4.15 registers	191
6.8.4.16 rs422_baud	191
6.8.4.17 rs422_cmd_wait	191
6.8.4.18 subregisters	191
6.8.4.19 subreglist	191
6.8.4.20 VREF	191
6.9 nsCamera.utils.Packet.Packet Class Reference	192
6.9.1 Detailed Description	192
6.9.2 Constructor & Destructor Documentation	193
6.9.2.1init()	193
6.9.3 Member Function Documentation	193
6.9.3.1 calculateCRC()	193
6.9.3.2 checkCRC()	194
6.9.3.3 checkReadPacket()	194
6.9.3.4 checkResponsePacket()	195
6.9.3.5 checkResponseString()	195
6.9.3.6 pktStr()	196
6.9.4 Member Data Documentation	196
6.9.4.1 addr	196
6.9.4.2 cmd	197
6.9.4.3 crc	197
6.9.4.4 data	197
6.9.4.5 payload	197
6.9.4.6 payload_length	197
6.9.4.7 preamble	197
6.9.4.8 PY3	197
6.9.4.9 seqID	198
6.9.4.10 type	198
6.10 nsCamera.comms.RS422.RS422 Class Reference	198
6.10.1 Detailed Description	199
6.10.2 Constructor & Destructor Documentation	199
6.10.2.1init()	199
6.10.3 Member Function Documentation	201
6.10.3.1 arm()	201
6.10.3.2 closeDevice()	202
6.10.3.3 readFrames()	202

6.10.3.4 readoff()	203
6.10.3.5 readSerial()	204
6.10.3.6 sendCMD()	205
6.10.3.7 serialClose()	207
6.10.3.8 writeSerial()	207
6.10.4 Member Data Documentation	208
6.10.4.1 _ser	208
6.10.4.2 baud	208
6.10.4.3 ca	208
6.10.4.4 datatimeout	209
6.10.4.5 logcrit	209
6.10.4.6 logdebug	209
6.10.4.7 logerr	209
6.10.4.8 loginfo	209
6.10.4.9 logwarn	209
6.10.4.10 mode	209
6.10.4.11 par	210
6.10.4.12 payloadsize	210
6.10.4.13 port	210
6.10.4.14 PY3	210
6.10.4.15 read_timeout	210
6.10.4.16 skipError	210
6.10.4.17 stop	210
6.10.4.18 write_timeout	211
6.11 nsCamera.sensors.sensorBase.sensorBase Class Reference	211
6.11.1 Detailed Description	212
6.11.2 Constructor & Destructor Documentation	212
6.11.2.1init()	212
6.11.3 Member Function Documentation	213
6.11.3.1 checkSensorVoltStat()	213
6.11.3.2 getManualTiming()	213
6.11.3.3 getSensorStatus()	214
6.11.3.4 getSensTemp()	214
6.11.3.5 getTiming()	215
6.11.3.6 init_board_specific()	216
6.11.3.7 parseReadoff()	217
6.11.3.8 reportStatusSensor()	217
6.11.3.9 selectOscillator()	218
6.11.3.10 setArbTiming()	218

6.11.3.11 setExtClk()
6.11.3.12 setHighFullWell()
6.11.3.13 setInterlacing()
6.11.3.14 setManualShutters()
6.11.3.15 setManualTiming()
6.11.3.16 setPhiDelay()
6.11.3.17 setTiming()
6.11.3.18 setTriggerDelay()
6.11.3.19 setZeroDeadTime()
6.11.4 Member Data Documentation
6.11.4.1 ca
6.11.4.2 fpganumID
6.11.4.3 logcrit
6.11.4.4 logdebug
6.11.4.5 logerr
6.11.4.6 loginfo
6.11.4.7 logwarn
6.11.4.8 sensfam
6.12 nsCamera.utils.Subregister.SubRegister Class Reference
6.12.1 Detailed Description
6.12.2 Constructor & Destructor Documentation
6.12.2.1init()
6.12.3 Member Data Documentation
6.12.3.1 addr
6.12.3.2 max
6.12.3.3 max_value
6.12.3.4 maxV
6.12.3.5 min
6.12.3.6 minV
6.12.3.7 name
6.12.3.8 register
6.12.3.9 resolution
6.12.3.10 start_bit
6.12.3.11 value
6.12.3.12 width
6.12.3.13 writable
6.13 nsCamera.comms.GigE.GigE.ZESTETM1_CARD_INFO Class Reference
6.13.1 Detailed Description
6.13.2 Member Data Documentation

	6.13.2.1 _fields	. 232
	6.13.2.2 ubyte4	. 232
	6.13.2.3 ubyte6	. 232
7	File Documentation	233
	7.1 C:/Users/hill35/git/camera_python/nsCamera/initpy File Reference	
	7.2initpy	
	7.3 C:/Users/hill35/git/camera_python/nsCamera/boards/initpy File Reference	
	7.4initpy	
	7.5 C:/Users/hill35/git/camera_python/nsCamera/comms/initpy File Reference	
	7.6 <u>init</u> .py	
	7.7 C:/Users/hill35/git/camera_python/nsCamera/sensors/initpy File Reference	
	7.8 <u>init</u> .py	
	7.9 C:/Users/hill35/git/camera_python/nsCamera/utils/initpy File Reference	
	7.10initpy	
	7.11 C:/Users/hill35/git/camera_python/nsCamera/boards/LLNL_v1.py File Reference	
	7.12 LLNL_v1.py	. 238
	7.13 C:/Users/hill35/git/camera_python/nsCamera/boards/LLNL_v4.py File Reference	. 249
	7.14 LLNL_v4.py	. 249
	7.15 C:/Users/hill35/git/camera_python/nsCamera/CameraAssembler.py File Reference	. 260
	7.16 CameraAssembler.py	. 261
	7.17 C:/Users/hill35/git/camera_python/nsCamera/comms/GigE.py File Reference	. 287
	7.18 GigE.py	. 287
	7.19 C:/Users/hill35/git/camera_python/nsCamera/comms/RS422.py File Reference	. 294
	7.20 RS422.py	. 294
	7.21 C:/Users/hill35/git/camera_python/nsCamera/comms/ZestETM1/Data.c File Reference	. 301
	7.21.1 Class Documentation	. 302
	7.21.1.1 struct ZESTETM1_WRITE_REG_CMD	. 302
	7.21.1.2 struct ZESTETM1_WRITE_REG_RESPONSE	. 302
	7.21.1.3 struct ZESTETM1_READ_REG_CMD	. 303
	7.21.1.4 struct ZESTETM1_READ_REG_RESPONSE	. 303
	7.21.1.5 struct ZESTETM1_MAILBOX_INT_CMD	. 303
	7.21.1.6 struct ZESTETM1_MAILBOX_INT_RESPONSE	. 303
	7.21.2 Macro Definition Documentation	. 304
	7.21.2.1 _CRT_SECURE_NO_WARNINGS	. 304
	7.21.2.2 ZESTETM1_COMMAND_MAILBOX_INT	. 304
	7.21.2.3 ZESTETM1_COMMAND_READ_REG	. 304
	7.21.2.4 ZESTETM1_COMMAND_SPI	. 304
	7.21.2.5 ZESTETM1_COMMAND_WRITE_REG	. 304

7.21.3 Function Documentation
7.21.3.1 ZestETM1_CloseConnection()
7.21.3.2 ZestETM1_OpenConnection()
7.21.3.3 ZestETM1_ReadData()
7.21.3.4 ZestETM1_SendCommand()
7.21.3.5 ZestETM1_SPIReadWrite()
7.21.3.6 ZestETM1_WriteData()
7.21.3.7 ZestETM1CloseConnection()
7.21.3.8 ZestETM1OpenConnection()
7.21.3.9 ZestETM1ReadData()
7.21.3.10 ZestETM1ReadRegister()
7.21.3.11 ZestETM1SetInterrupt()
7.21.3.12 ZestETM1SPIReadWrite()
7.21.3.13 ZestETM1WriteData()
7.21.3.14 ZestETM1WriteRegister()
7.22 Data.c
7.23 C:/Users/hill35/git/camera_python/nsCamera/comms/ZestETM1/Error.c File Reference
7.23.1 Function Documentation
7.23.1.1 ZestETM1GetErrorMessage()
7.23.1.2 ZestETM1RegisterErrorHandler()
7.23.2 Variable Documentation
7.23.2.1 ZestETM1_ErrorHandler
7.23.2.2 ZestETM1_ErrorStrings
7.24 Error.c
7.25 C:/Users/hill35/git/camera_python/nsCamera/comms/ZestETM1/Main.c File Reference
7.25.1 Function Documentation
7.25.1.1 ZestETM1Close()
7.25.1.2 ZestETM1Init()
7.26 Main.c
7.27 C:/Users/hill35/git/camera_python/nsCamera/comms/ZestETM1/Private.h File Reference
7.27.1 Class Documentation
7.27.1.1 struct ZESTETM1_CONNECTION_STRUCT
7.27.2 Macro Definition Documentation
7.27.2.1 closesocket
7.27.2.2 SD_BOTH
7.27.2.3 ZESTETM1_CONNECTION_HANDLE_MAGIC
7.27.2.4 ZESTETM1_DEFAULT_TIMEOUT
7.27.2.5 ZESTETM1_ERROR
7.27.2.6 ZESTETM1_ERROR_CONN

7.27.2.7 ZESTETM1_ERROR_GENERAL
7.27.2.8 ZESTETM1_ERROR_STRING
7.27.2.9 ZESTETM1_RATE_10MHz
7.27.2.10 ZESTETM1_RATE_20MHz
7.27.2.11 ZESTETM1_RATE_40MHz
7.27.2.12 ZESTETM1_REVERSE
7.27.2.13 ZESTETM1_USER_DEVICE_ID
7.27.3 Typedef Documentation
7.27.3.1 SOCKET
7.27.4 Function Documentation
7.27.4.1 ZestETM1_CloseConnection()
7.27.4.2 ZestETM1_EraseFlashSector()
7.27.4.3 ZestETM1_OpenConnection()
7.27.4.4 ZestETM1_SendCommand()
7.27.4.5 ZestETM1_SPIReadWrite()
7.27.4.6 ZestETM1_WriteFlash()
7.27.5 Variable Documentation
7.27.5.1 ZestETM1_ErrorHandler
7.27.5.2 ZestETM1_ErrorStrings
7.28 Private.h
7.29 C:/Users/hill35/git/camera_python/nsCamera/comms/ZestETM1/UPnP.c File Reference
7.29.1 Class Documentation
7.29.1.1 struct ZESTETM1_GET_SETTINGS_CMD
7.29.1.2 struct ZESTETM1_GET_SETTINGS_RESPONSE
7.29.2 Macro Definition Documentation
7.29.2.1 _CRT_SECURE_NO_WARNINGS
7.29.2.2 _strnicmp
7.29.2.3 ZESTETM1_COMMAND_GET_SETTINGS
7.29.2.4 ZESTETM1_UPNP_ADDR
7.29.2.5 ZESTETM1_UPNP_PORT
7.29.3 Function Documentation
7.29.3.1 ZestETM1_GetAllAdapters()
7.29.3.2 ZestETM1_GetCardInfo()
7.29.3.3 ZestETM1_HTTPGet()
7.29.3.4 ZestETM1_JoinGroup()
7.29.3.5 ZestETM1_LeaveGroup()
7.29.3.6 ZestETM1_ReadSettings()
7.29.3.7 ZestETM1CountCards()
7.29.3.8 ZestETM1FreeCards()

7.29.3.9 ZestETM1GetCardInfo()
7.29.4 Variable Documentation
7.29.4.1 ZestETM1_SearchReq
7.30 UPnP.c
7.31 C:/Users/hill35/git/camera_python/nsCamera/comms/ZestETM1/ZestETM1.h File Reference 362
7.31.1 Class Documentation
7.31.1.1 struct ZESTETM1_CARD_INFO
7.31.2 Macro Definition Documentation
7.31.2.1 ZESTETM1_ERROR_BASE
7.31.2.2 ZESTETM1_INFO_BASE
7.31.2.3 ZESTETM1_VERSION_FALLBACK
7.31.2.4 ZESTETM1_WARNING_BASE
7.31.3 Typedef Documentation
7.31.3.1 ZESTETM1_CONNECTION
7.31.3.2 ZESTETM1_ERROR_FUNC
7.31.3.3 ZESTETM1_HANDLE
7.31.4 Enumeration Type Documentation
7.31.4.1 ZESTETM1_CONNECTION_TYPE
7.31.4.2 ZESTETM1_SPI_RATE
7.31.4.3 ZESTETM1_STATUS
7.31.5 Function Documentation
7.31.5.1 ZestETM1Close()
7.31.5.2 ZestETM1CloseConnection()
7.31.5.3 ZestETM1CountCards()
7.31.5.4 ZestETM1FreeCards()
7.31.5.5 ZestETM1GetCardInfo()
7.31.5.6 ZestETM1GetErrorMessage()
7.31.5.7 ZestETM1Init()
7.31.5.8 ZestETM1OpenConnection()
7.31.5.9 ZestETM1ReadData()
7.31.5.10 ZestETM1ReadRegister()
7.31.5.11 ZestETM1RegisterErrorHandler()
7.31.5.12 ZestETM1SetInterrupt()
7.31.5.13 ZestETM1SPIReadWrite()
7.31.5.14 ZestETM1WriteData()
7.31.5.15 ZestETM1WriteRegister()
7.32 ZestETM1.h
7.33 C:/Users/hill35/git/camera_python/nsCamera/sensors/daedalus.py File Reference
7 34 daedalus ny

	7.35 C:/Users/hill35/git/camera_python/nsCamera/sensors/icarus.py File Reference	. 387
	7.36 icarus.py	. 387
	7.37 C:/Users/hill35/git/camera_python/nsCamera/sensors/icarus2.py File Reference	. 390
	7.38 icarus2.py	. 390
	7.39 C:/Users/hill35/git/camera_python/nsCamera/sensors/sensorBase.py File Reference	. 392
	7.40 sensorBase.py	. 392
	7.41 C:/Users/hill35/git/camera_python/nsCamera/utils/crc16pure.py File Reference	. 401
	7.42 crc16pure.py	. 402
	7.43 C:/Users/hill35/git/camera_python/nsCamera/utils/FlatField.py File Reference	. 405
	7.44 FlatField.py	. 406
	7.45 C:/Users/hill35/git/camera_python/nsCamera/utils/misc.py File Reference	
	7.46 misc.py	. 409
	7.47 C:/Users/hill35/git/camera_python/nsCamera/utils/Packet.py File Reference	. 418
	7.48 Packet.py	. 419
	7.49 C:/Users/hill35/git/camera_python/nsCamera/utils/Subregister.py File Reference	. 421
	7.50 Subregister.py	. 422
In	dex	423

## **Namespace Index**

## 1.1 Package List

Here are the packages with brief descriptions (if available):

nsCamera
nsCamera.boards
nsCamera.boards.LLNL_v1
nsCamera.boards.LLNL_v4
nsCamera.CameraAssembler
nsCamera.comms
nsCamera.comms.GigE
nsCamera.comms.RS422
nsCamera.sensors
nsCamera.sensors.daedalus
nsCamera.sensors.icarus
nsCamera.sensors.icarus2
nsCamera.sensors.sensorBase
nsCamera.utils
nsCamera.utils.crc16pure
nsCamera.utils.FlatField
nsCamera.utils.misc
nsCamera.utils.Packet
nsCamera.utils.Subregister

2 Namespace Index

## **Hierarchical Index**

## 2.1 Class Hierarchy

This inheritance list is sorted roughly, but not completely, alphabetically:

nsCamera.CameraAssembler.CameraAssembler
nsCamera.utils.misc.fakeCA
nsCamera.comms.GigE.GigE
nsCamera.boards.LLNL_v1.llnl_v1
nsCamera.boards.LLNL_v4.llnl_v4
object
nsCamera.sensors.sensorBase.sensorBase
nsCamera.sensors.daedalus.daedalus
nsCamera.sensors.icarus2.icarus2
nsCamera.sensors.icarus.icarus
nsCamera.utils.Packet.Packet
nsCamera.comms.RS422.RS422
C.Structure
nsCamera.comms.GigE.GigE.ZESTETM1_CARD_INFO
nsCamera.comms.GigE.GigE.ZESTETM1_CARD_INFO
nsCamera.utils.Subregister.SubRegister
nsCamera.utils.Subregister.SubRegister
nsCamera.utils.Subregister.SubRegister       .228         ZESTETM1_CARD_INFO       .362         ZESTETM1_CONNECTION_STRUCT       .331
nsCamera.utils.Subregister.SubRegister       228         ZESTETM1_CARD_INFO       362         ZESTETM1_CONNECTION_STRUCT       331         ZESTETM1_GET_SETTINGS_CMD       342
nsCamera.utils.Subregister.SubRegister       228         ZESTETM1_CARD_INFO       362         ZESTETM1_CONNECTION_STRUCT       331         ZESTETM1_GET_SETTINGS_CMD       342         ZESTETM1_GET_SETTINGS_RESPONSE       342
nsCamera.utils.Subregister.SubRegister       228         ZESTETM1_CARD_INFO       362         ZESTETM1_CONNECTION_STRUCT       331         ZESTETM1_GET_SETTINGS_CMD       342         ZESTETM1_GET_SETTINGS_RESPONSE       342         ZESTETM1_MAILBOX_INT_CMD       301
nsCamera.utils.Subregister.SubRegister       228         ZESTETM1_CARD_INFO       362         ZESTETM1_CONNECTION_STRUCT       331         ZESTETM1_GET_SETTINGS_CMD       342         ZESTETM1_GET_SETTINGS_RESPONSE       342         ZESTETM1_MAILBOX_INT_CMD       301         ZESTETM1_MAILBOX_INT_RESPONSE       301
nsCamera.utils.Subregister.SubRegister       228         ZESTETM1_CARD_INFO       362         ZESTETM1_CONNECTION_STRUCT       331         ZESTETM1_GET_SETTINGS_CMD       342         ZESTETM1_GET_SETTINGS_RESPONSE       342         ZESTETM1_MAILBOX_INT_CMD       301         ZESTETM1_MAILBOX_INT_RESPONSE       301         ZESTETM1_READ_REG_CMD       301

4 Hierarchical Index

## **Class Index**

## 3.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

sCamera.CameraAssembler.CameraAssembler	39
sCamera.sensors.daedalus.daedalus	7
sCamera.utils.misc.fakeCA	8
sCamera.comms.GigE.GigE	22
sCamera.sensors.icarus.icarus	37
sCamera.sensors.icarus2.icarus2	4
sCamera.boards.LLNL_v1.llnl_v1	52
sCamera.boards.LLNL_v4.llnl_v4	1
sCamera.utils.Packet.Packet	)2
sCamera.comms.RS422.RS422	8
sCamera.sensors.sensorBase.sensorBase	1
sCamera.utils.Subregister.SubRegister	28
sCamera.comms.GigE.GigE.ZESTETM1	31

6 Class Index

## **File Index**

## 4.1 File List

Here is a list of all files with brief descriptions:

C:/Users/hill35/git/camera_python/nsCamera/initpy
C:/Users/hill35/git/camera_python/nsCamera/CameraAssembler.py
C:/Users/hill35/git/camera_python/nsCamera/boards/initpy
C:/Users/hill35/git/camera_python/nsCamera/boards/LLNL_v1.py
C:/Users/hill35/git/camera_python/nsCamera/boards/LLNL_v4.py
C:/Users/hill35/git/camera_python/nsCamera/comms/initpy
C:/Users/hill35/git/camera_python/nsCamera/comms/GigE.py
C:/Users/hill35/git/camera_python/nsCamera/comms/RS422.py
C:/Users/hill35/git/camera_python/nsCamera/comms/ZestETM1/Data.c
C:/Users/hill35/git/camera_python/nsCamera/comms/ZestETM1/Error.c
C:/Users/hill35/git/camera_python/nsCamera/comms/ZestETM1/Main.c
C:/Users/hill35/git/camera_python/nsCamera/comms/ZestETM1/Private.h
C:/Users/hill35/git/camera_python/nsCamera/comms/ZestETM1/UPnP.c
C:/Users/hill35/git/camera_python/nsCamera/comms/ZestETM1/ZestETM1.h
C:/Users/hill35/git/camera_python/nsCamera/sensors/initpy
C:/Users/hill35/git/camera_python/nsCamera/sensors/daedalus.py
C:/Users/hill35/git/camera_python/nsCamera/sensors/icarus.py
C:/Users/hill35/git/camera_python/nsCamera/sensors/icarus2.py
C:/Users/hill35/git/camera_python/nsCamera/sensors/sensorBase.py
C:/Users/hill35/git/camera_python/nsCamera/utils/initpy
C:/Users/hill35/git/camera_python/nsCamera/utils/crc16pure.py
C:/Users/hill35/git/camera_python/nsCamera/utils/FlatField.py
C:/Users/hill35/git/camera_python/nsCamera/utils/misc.py
C:/Users/hill35/git/camera_python/nsCamera/utils/Packet.py
C:/Users/hill35/git/camera_python/nsCamera/utils/Subregister.py 421

8 File Index

## **Namespace Documentation**

## 5.1 nsCamera Namespace Reference

#### **Namespaces**

- namespace boards
- namespace CameraAssembler
- · namespace comms
- namespace sensors
- · namespace utils

### **Variables**

• list \_\_all\_\_ = ["CameraAssembler"]

### 5.1.1 Detailed Description

```
Created on Tue Mar 22 15:47:43 2016
```

The Package includes a Camera object and an assembler.

The camera object will be the workhorse of the API. The assembler is used to create the Camera object.

Author: Matthew Dayton (dayton5@llnl.gov)

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved. LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL. 'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)

#### 5.1.2 Variable Documentation

```
5.1.2.1 __all__
list nsCamera.__all__ = ["CameraAssembler"] [private]
Definition at line 26 of file init .py.
```

## 5.2 nsCamera.boards Namespace Reference

### **Namespaces**

- namespace LLNL v1
- namespace LLNL v4

#### **Variables**

```
• list __all__ = ["LLNL_v1", "LLNL_v4"]
```

### 5.2.1 Detailed Description

This package is a collection of modules that represent the camera boards Each board has its own number of ADCs, POTs, and sensors. More devices can be added in the future. The list of imports will grow as we make more types of boards.

```
Author: Matthew Dayton (dayton5@llnl.gov)
```

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved. LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL. 'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)

### 5.2.2 Variable Documentation

```
5.2.2.1 __all__
list nsCamera.boards.__all__ = ["LLNL_v1", "LLNL_v4"] [private]

Definition at line 23 of file __init__.py.
```

### 5.3 nsCamera.boards.LLNL v1 Namespace Reference

#### Classes

class IInl v1

### 5.3.1 Detailed Description

LLNLv1 board definition, including monitors, pots, and other board-specific settings

Author: Jeremy Martin Hill (jerhill@llnl.gov)
Author: Matthew Dayton (dayton5@llnl.gov)

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved. LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL. 'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)

## 5.4 nsCamera.boards.LLNL\_v4 Namespace Reference

#### **Classes**

· class IInl v4

### 5.4.1 Detailed Description

LLNLv4 board definition, including monitors, DACS, and other board-specific settings

Author: Jeremy Martin Hill (jerhill@llnl.gov)
Author: Matthew Dayton (dayton5@llnl.gov)

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved. LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL. 'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)

## 5.5 nsCamera.CameraAssembler Namespace Reference

#### **Classes**

class CameraAssembler

### 5.5.1 Detailed Description

CameraAssembler assembles the separate camera parts into a camera object. This object controls a combination of three components:

```
    board: FPGA board -- LLNL_V1, LLNL_V4
    comms: communication interface -- GigE, RS422
    sensor: sensor type -- icarus, icarus2, daedalus
    Author: Jeremy Martin Hill (jerhill@llnl.gov)
    Author: Matthew Dayton (dayton5@llnl.gov)
```

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved. LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL. 'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)

### 5.6 nsCamera.comms Namespace Reference

#### **Namespaces**

- namespace GigE
- namespace RS422

#### **Variables**

• list \_\_all\_\_ = ["RS422", "GigE"]

#### 5.6.1 Detailed Description

```
Created on Tue Mar 22 15:47:43 2016
```

This package is a collection of modules for uniform handling of the nsCamera's communication systems

Author: Matthew Dayton (dayton5@llnl.gov)

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved. LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL. 'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)

### 5.6.2 Variable Documentation

```
5.6.2.1 __all__
list nsCamera.comms.__all__ = ["RS422", "GigE"] [private]
Definition at line 25 of file init .py.
```

# 5.7 nsCamera.comms.GigE Namespace Reference

#### Classes

class GigE

# 5.7.1 Detailed Description

```
Gigabit Ethernet interface for nsCamera.

Author: Jeremy Martin Hill (jerhill@llnl.gov)

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.

LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.

'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)
```

# 5.8 nsCamera.comms.RS422 Namespace Reference

### Classes

• class RS422

# 5.8.1 Detailed Description

```
RS422 driver for nsCamera

Author: Brad Funsten (funsten1@llnl.gov)

Author: Jeremy Martin Hill (jerhill@llnl.gov)

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.

LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.

'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)
```

# 5.9 nsCamera.sensors Namespace Reference

#### **Namespaces**

- · namespace daedalus
- namespace icarus
- namespace icarus2
- namespace sensorBase

#### **Variables**

```
• list __all__ = ["icarus", "icarus2", "daedalus"]
```

## 5.9.1 Detailed Description

```
This package is a collection of modules for uniform handling of the various sensors

Author: Jeremy Martin Hill (jerhill@llnl.gov)

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.

LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.

'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)
```

#### 5.9.2 Variable Documentation

```
5.9.2.1 __all__
list nsCamera.sensors.__all__ = ["icarus", "icarus2", "daedalus"] [private]

Definition at line 23 of file __init__.py.
```

# 5.10 nsCamera.sensors.daedalus Namespace Reference

#### Classes

· class daedalus

# 5.10.1 Detailed Description

Parameters and functions specific to the daedalus three-frame sensor

Author: Jeremy Martin Hill (jerhill@llnl.gov)

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved. LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL. 'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)

# 5.11 nsCamera.sensors.icarus Namespace Reference

#### Classes

· class icarus

# 5.11.1 Detailed Description

Parameters and functions specific to the icarus two-frame sensor

\*\*\*Do not use this file as a template for new code development\*\*\*

Author: Jeremy Martin Hill (jerhill@llnl.gov)

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved. LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL. 'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)

# 5.12 nsCamera.sensors.icarus2 Namespace Reference

#### **Classes**

· class icarus2

# 5.12.1 Detailed Description

Parameters and functions specific to the four-frame icarus2 sensor

Author: Jeremy Martin Hill (jerhill@llnl.gov)

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved. LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL. 'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)

# 5.13 nsCamera.sensors.sensorBase Namespace Reference

#### **Classes**

class sensorBase

# 5.13.1 Detailed Description

Superclass for nsCamera sensors

Author: Jeremy Martin Hill (jerhill@llnl.gov)

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved. LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL. 'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)

# 5.14 nsCamera.utils Namespace Reference

### **Namespaces**

- namespace crc16pure
- namespace FlatField
- · namespace misc
- namespace Packet
- namespace Subregister

#### **Variables**

• list \_\_all\_\_ = ["SubRegister", "Packet", "FlatField", "misc"]

# 5.14.1 Detailed Description

This package is a collection of utility classes for the CameraAssembler

Author: Jeremy Martin Hill (jerhill@llnl.gov)

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved. LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL. 'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)

#### 5.14.2 Variable Documentation

```
5.14.2.1 all
```

list nsCamera.utils.\_\_all\_\_ = ["SubRegister", "Packet", "FlatField", "misc"] [private]

Definition at line 30 of file \_\_init\_\_.py.

# 5.15 nsCamera.utils.crc16pure Namespace Reference

#### **Functions**

- crc16 (data, crc, table)
- crc16xmodem (data, crc=0)

#### **Variables**

list CRC16 XMODEM TABLE

# 5.15.1 Detailed Description

Pure python library for calculating CRC16 NOTE: modified slightly to combine Python 2 and Python 3 versions in single file

#### 5.15.2 Function Documentation

### 5.15.2.1 \_crc16()

```
nsCamera.utils.crc16pure._crc16 (
                data,
                crc.
                table ) [protected]
Calculate CRC16 using the given table.
              - data for calculating CRC, must be a string
'crc'
              - initial value
`table`
              - table for caclulating CRC (list of 256 integers)
Return calculated value of CRC
Definition at line 299 of file crc16pure.py.
00299 def _crc16(data, crc, table):
           """Calculate CRC16 using the given table.
          `data` - data for calculating CRC, must be a string
`crc` - initial value

`table` - table for caclulating CRC (list of 256 integers)
00302
00303
00304
          Return calculated value of CRC
00305
00306
          for byte in data:
              if sys.version_info > (3,):
                  crc = ((crc « 8) & 0xFF00) ^ table[((crc » 8) & 0xFF) ^ byte]
00308
00309
00310
                   crc = ((crc « 8) & 0xFF00) ^ table[((crc » 8) & 0xFF) ^ ord(byte)]
00311
00312
         return crc & 0xFFFF
00313
00314
```

### 5.15.2.2 crc16xmodem()

```
nsCamera.utils.crc16pure.crc16xmodem (

data,

crc = 0 )

Calculate CRC-CCITT (XModem) variant of CRC16.

'data' - data for calculating CRC, must be a string
'crc' - initial value
Return calculated value of CRC
```

#### Definition at line 315 of file crc16pure.py.

```
00315 def crc16xmodem(data, crc=0):

00316 """Calculate CRC-CCITT (XModem) variant of CRC16.

00317 'data' - data for calculating CRC, must be a string

00318 'crc' - initial value

00319 Return calculated value of CRC

00320 """

00321 return _crc16(data, crc, CRC16_XMODEM_TABLE)
```

### 5.15.3 Variable Documentation

### 5.15.3.1 CRC16\_XMODEM\_TABLE

list nsCamera.utils.crc16pure.CRC16\_XMODEM\_TABLE

Definition at line 39 of file crc16pure.py.

# 5.16 nsCamera.utils.FlatField Namespace Reference

#### **Functions**

- getFilenames (frame="Frame 1")
- getROlvector (imgfilename, roi)
- tslopes (x, y)
- generateFF (FRAMES=["Frame\_0", "Frame\_1", "Frame\_2", "Frame\_3"], roi=[0, 0, 512, 1024], directory="", ncores=-1)
- removeFF (filename, directory="", roi=[0, 0, 512, 1024])
- removeFFall (directory="", FRAMES=["Frame\_0", "Frame\_1", "Frame\_2", "Frame\_3"], roi=[0, 0, 512, 1024])

#### **Variables**

- parser = argparse.ArgumentParser()
- action
- dest
- default
- help
- nargs
- args = parser.parse\_args()
- list framelist = ["Frame\_" + str(frame) for frame in args.frames]
- directory

### 5.16.1 Detailed Description

```
Functions for batch flat-field image corrections

***Do not use this file as a template for new code development***

Author: Jeremy Martin Hill (jerhill@llnl.gov)

Author: Matthew Dayton (dayton5@llnl.gov)

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.

LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL. 'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)
```

00034 00035

00036

00037 00038 00039 onlyfiles = next(os.walk("./"))[2]

return [k for k in onlyfiles if frame in k and "tif" in k]

#### 5.16.2 Function Documentation

#### 5.16.2.1 generateFF()

```
nsCamera.utils.FlatField.generateFF (
                FRAMES = ["Frame_0", "Frame_1", "Frame_2", "Frame_3"],
                 roi = [0, 0, 512, 1024],
                 directory = "",
                ncores = -1 )
Definition at line 58 of file FlatField.py.
00063):
00064
          # TODO: documentation
00065
          # use of ROI here not compatible with use of ROI in removeFF
00066
00067
          if directory:
00068
              cwd = os.getcwd()
              newpath = os.path.join(cwd, directory)
00069
              os.chdir(newpath)
00070
00071
          if not FRAMES:
              print("No framelist provided, defaulting to four frames")
FRAMES = ["Frame_0", "Frame_1", "Frame_2", "Frame_3"]
00072
00073
00074
          for f in FRAMES:
00075
              files = getFilenames(frame=f)
00076
              imgslist = [getROIvector(fn, roi) for fn in files] # a list of flattened images
              imgsarray = np.vstack(imgslist) # turn the list into an array
npix = np.shape(imgsarray)[1] # total number of pixels
00077
00078
00079
              x = np.median(imgsarray, axis=1) # median of each image used for flat fielding
              y = []
00080
00081
              for i in range(npix):
00082
                  \# each member of y represents a pixel, as a list of magnitudes over all the
00083
                      images
00084
                  y.append(imgsarray[:, i])
00085
              # get pixel gain and offset for flatfield ff using Thiel-Sen slopes
00086
              ff = []
00087
              \label{eq:ff}  \mbox{ff = parallel.Parallel(n_jobs=ncores, verbose=5, pre\_dispatch="2 * n_jobs")() } 
00088
                  delayed(tslopes)(x, pixel) for pixel in y
00089
00090
              \# x is the dependent variable; here uses median of image as characteristic of
00091
                  noise level
00092
              m, c = zip(*ff)
                                # separate into gain and offset
00093
              m = np.array(m)
00094
              m[m < 0.1] = 0.1 # handle outliers
              m[m > 1000] = 1000 # handle outliers
00095
00096
              m = 1.0 / m
00097
              m = m.reshape(roi[3] - roi[1], roi[2] - roi[0]) # turn into matrix
00098
              c = np.array(c).reshape(roi[3] - roi[1], roi[2] - roi[0]) # turn into matrix
00099
              with open("px_gain_%s.txt" % f.replace("Frame_", "f"), "w+") as file:
00100
                  np.savetxt(file, m)
00101
00102
              with open("px_off_%s.txt" % f.replace("Frame_", "f"), "w+") as file:
00103
                  np.savetxt(file, c)
00104
00105
5.16.2.2 getFilenames()
nsCamera.utils.FlatField.getFilenames (
                frame = "Frame 1" )
get a list of tiff filenames in current working director for frame
Definition at line 32 of file FlatField.py.
00032 def getFilenames(frame="Frame 1"):
00033
          get a list of tiff filenames in current working director for frame
```

### 5.16.2.3 getROlvector()

```
nsCamera.utils.FlatField.getROIvector (
                imgfilename,
               roi )
return a numpy row vector of version of the image
Definition at line 40 of file FlatField.py.
00040 def getRoIvector(imgfilename, roi):
00041
00042
          return a numpy row vector of version of the image
00043
00044
         img = imread(imgfilename)
00045
          vroi = img[(roi[1]) : (roi[3]), (roi[0]) : (roi[2])].flattenlist()
00046
          return vroi
00047
00048
```

#### 5.16.2.4 removeFF()

```
nsCamera.utils.FlatField.removeFF (
              filename,
              directory = "",
              roi = [0, 0, 512, 1024] )
```

#### Definition at line 106 of file FlatField.py.

```
00106 def removeFF(filename, directory="", roi=[0, 0, 512, 1024]):
           if directory:
00108
               cwd = os.getcwd()
00109
               newpath = os.path.join(cwd, directory)
00110
               os.chdir(newpath)
         framenum = re.search("Frame_(\d)", filename).group(1)
gainFilename = "px_gain_f" + framenum + ".txt"
00111
00112
00113
          gainall = np.loadtxt(gainFilename)
          gain = gainall[(roi[1]) : (roi[3]), (roi[0]) : (roi[2])]
offFilename = "px_off_f" + framenum + ".txt"
00114
00115
           offsetall = np.loadtxt(offFilename, dtype="uint32")
00116
00117
           offset = offsetall[(roi[1]) : (roi[3]), (roi[0]) : (roi[2])]
00118
00119
           beforeImageall = imread(filename)
           beforeImage = beforeImageall[(roi[1]) : (roi[3]), (roi[0]) : (roi[2])]
00120
00121
          imageMed = np.median(beforeImage)
00122
00123
           flat = imageMed * gain + offset
00124
           flat = flat.clip(0)
           fix = beforeImage - flat
00125
           clipped = fix.clip(0)
fixinit = clipped.astype("uint16")
00126
00127
00128
           fiximg = Image.fromarray(fixinit)
00129
           fixFilename = filename[:-4] + "ff" + filename[-4:]
00130
           fiximg.save(fixFilename)
00131
00132
```

# 5.16.2.5 removeFFall()

```
nsCamera.utils.FlatField.removeFFall (
              directory = "",
              FRAMES = ["Frame_0", "Frame_1", "Frame_2", "Frame_3"],
              roi = [0, 0, 512, 1024] )
```

Definition at line 133 of file FlatField.py.

```
00138
          cwd = os.getcwd()
00139
          if directory:
00140
              newpath = os.path.join(cwd, directory)
00141
00142
              newpath = cwd
00143
       os.chdir(newpath)
00144
         files = next(os.walk("./"))[2]
        files = next(os.walk
filelist = []
for frame in FRAMES:
00145
00147
              filelist.extend([k for k in files if frame in k and "tif" in k])
         for fname in filelist:
00148
              removeFF(fname, directory, roi)
00150
00151
```

#### 5.16.2.6 tslopes()

```
nsCamera.utils.FlatField.tslopes ( x, \\ y \ ) theilslopes implements a method for robust linear regression. It computes the slope as the median of all slopes between paired values.
```

#### Definition at line 49 of file FlatField.py.

```
00049 def tslopes(x, y):
00050 """

00051 theilslopes implements a method for robust linear regression.
00052 It computes the slope as the median of all slopes between paired values.
00053 """

00054 val = theilslopes(x, y)
00055 return [val[0], val[1]]
00056
00057
```

# 5.16.3 Variable Documentation

#### 5.16.3.1 action

```
nsCamera.utils.FlatField.action
```

Definition at line 157 of file FlatField.py.

## 5.16.3.2 args

```
nsCamera.utils.FlatField.args = parser.parse_args()
```

Definition at line 167 of file FlatField.py.

#### 5.16.3.3 default

```
nsCamera.utils.FlatField.default
```

Definition at line 157 of file FlatField.py.

### 5.16.3.4 dest

nsCamera.utils.FlatField.dest

Definition at line 157 of file FlatField.py.

#### 5.16.3.5 directory

nsCamera.utils.FlatField.directory

Definition at line 169 of file FlatField.py.

#### 5.16.3.6 framelist

```
nsCamera.utils.FlatField.framelist = ["Frame_" + str(frame) for frame in args.frames]
```

Definition at line 168 of file FlatField.py.

### 5.16.3.7 help

nsCamera.utils.FlatField.help

Definition at line 157 of file FlatField.py.

#### 5.16.3.8 nargs

nsCamera.utils.FlatField.nargs

Definition at line 161 of file FlatField.py.

# 5.16.3.9 parser

nsCamera.utils.FlatField.parser = argparse.ArgumentParser()

Definition at line 155 of file FlatField.py.

# 5.17 nsCamera.utils.misc Namespace Reference

#### Classes

class fakeCA

#### **Functions**

- makeLogLabels (logtag, label)
- getEnter (text)
- checkCRC (rval)
- str2bytes (astring)
- bytes2str (bytesequence)
- str2nparray (valstring)
- flattenlist (x)
- generateFrames (camassem, data, columns=1)
- loadDumpedData (filename="frames.txt", path=None, filetype="txt", sensor="daedalus", firstframe=None, last-frame=None, width=None, height=None, padToFull=None, firstrow=None, lastrow=None, maxwidth=None, maxheight=None, bytesperpixel=None, interlacing=None, columns=1)
- saveTiffs (self, frames, path=None, filename="Frame", prefix=None, index=None)
- plotFrames (self, frames, index=None)
- · partition (self, frames, columns)

## 5.17.1 Detailed Description

```
Miscellaneous utilities, including batch processing of images acquired using the nsCamera. These are functions that don't require a cameraAssembler object to be instantiated before use.

Author: Jeremy Martin Hill (jerhill@llnl.gov)

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved. LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL. 'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)
```

#### 5.17.2 Function Documentation

#### 5.17.2.1 bytes2str()

```
nsCamera.utils.misc.bytes2str (

bytesequence )

Python-version-agnostic converter of bytes to hexadecimal strings

Args:

bytesequence: sequence of bytes as string (Py2) or bytes (Py3)

Returns:

hexadecimal string representation of 'bytes' without '0x'
```

```
Definition at line 123 of file misc.py.
```

```
00123 def bytes2str(bytesequence):
00125
          Python-version-agnostic converter of bytes to hexadecimal strings
00126
00127
00128
             bytesequence: sequence of bytes as string (Py2) or bytes (Py3)
00129
00130
          Returns:
          hexadecimal string representation of 'bytes' without '0x' """
00131
00132
00133
         try:
00134
              estring = binascii.b2a_hex(bytesequence)
00135
          except TypeError:
00136
              logging.error(
00137
                  "ERROR: bytes2str: Invalid byte sequence: '{bytesequence}'; returning an"
00138
                  " empty string".format(bytesequence=bytesequence)
00139
              )
00140
              return ""
         python, _, _, _, _ = sys.version_info
if python >= 3:
00141
00142
00143
             estring = str(estring)[2:-1]
00144
         return estring
00145
00146
```

#### 5.17.2.2 checkCRC()

### Definition at line 72 of file misc.py.

```
00072 def checkCRC(rval):
00073
00074
          Calculate CRC for rval[:-4] and compare with expected CRC in rval[-4:]
00075
00076
         Args:
00077
             rval: hexadecimal string
00078
00079
         Returns:
         boolean, True if CRCs match, False if they don't match or the input is invalid
00080
00081
00082
         if not isinstance(rval, str) or len(rval) < 5:</pre>
             logging.error("ERROR: checkCRC: Invalid input: {rval}".format(rval=rval))
00083
00084
              return False
00085
         data_crc = int(rval[-4:], base=16)
00086
         CRC_calc = crc16pure.crc16xmodem(str2bytes(rval[:-4]))
00087
         return CRC_calc == data_crc
00088
00089
```

#### 5.17.2.3 flattenlist()

```
nsCamera.utils.misc.flattenlist ( x\ ) Flatten list of lists recursively into single list
```

### Definition at line 179 of file misc.py.

```
00179 def flattenlist(x):
00180
00181
          Flatten list of lists recursively into single list
00182
00183
          python, _, _, _ = sys.version_info
00184
00185
              if python >= 3:
00186
                  if isinstance(x, collections.abc.Iterable):
00187
                      return [a for i in x for a in flattenlist(i)]
00188
00189
                     return [x]
00190
              else:
                  if isinstance(x, collections.Iterable):
00191
00192
                      return [a for i in x for a in flattenlist(i)]
00193
                  else:
00194
                      return [x]
00195
          except RecursionError:
00196
              logging.error(
                  "ERROR: flattenlist: input '{x}' is pathological and cannot be flattened."
00197
00198
                  " Attempting to return the input unchanged"
00199
00200
              return x
00201
00202
```

#### 5.17.2.4 generateFrames()

```
nsCamera.utils.misc.generateFrames (
              camassem,
              data,
              columns = 1)
Processes data stream from board into frames and applies sensor-specific parsing.
  Generates padded data for full-size option of setRows.
If used for offline processing, replace the 'self' object with the parameter object
  returned by loadDumpedData().
If the data stream is incomplete (e.g., from an interrupted download), the data is
  padded with zeros to the correct length.
Args:
    self: during normal operation, is the CameraAssembler object. During offline
      operation, is a parameters object as returned by loadDumpedData()
    data: text stream from board or loaded file, or numpy from loaded file
    columns: 1 for single image, 2 for separate hemisphere images
Returns: list of parsed frames
```

#### Definition at line 261 of file misc.py.

```
00261 def generateFrames(camassem, data, columns=1):
00263
          Processes data stream from board into frames and applies sensor-specific parsing.
00264
            Generates padded data for full-size option of setRows.
00265
          If used for offline processing, replace the 'self' object with the parameter object
00266
            returned by loadDumpedData().
00267
          If the data stream is incomplete (e.g., from an interrupted download), the data is
00268
            padded with zeros to the correct length.
00269
00270
00271
              self: during normal operation, is the CameraAssembler object. During offline
00272
               operation, is a parameters object as returned by loadDumpedData()
00273
              data: text stream from board or loaded file, or numpy from loaded file
00274
              columns: 1 for single image, 2 for separate hemisphere images
00275
00276
          Returns: list of parsed frames
00277
          logging.debug("DEBUG" + camassem.logtag + ": generateFrames")
00278
00279
          if isinstance(data[0], str):
```

```
00280
              allframes = str2nparray(data)
00281
         else:
00282
              allframes = data
00283
          nframes = camassem.sensor.lastframe - camassem.sensor.firstframe + 1
00284
          frames = [0] * nframes
00285
         framesize = camassem.sensor.width * (
00286
              camassem.sensor.lastrow - camassem.sensor.firstrow + 1
00287
00288
         if hasattr(camassem, "ca"):
             padIt = camassem.ca.padToFull
00290
         else:
00291
             padIt = camassem.padToFull
00292
          if padIt:
00293
              toprows = camassem.sensor.firstrow
00294
              botrows = (camassem.sensor.maxheight - 1) - camassem.sensor.lastrow
00295
              padtop = np.zeros(toprows * camassem.sensor.maxwidth, dtype=int)
              padbot = np.zeros(botrows * camassem.sensor.maxwidth, dtype=int)
00296
00297
              for n in range(nframes):
00298
                 thisframe = np.concatenate(
00299
                      (padtop, allframes[n * framesize : ((n + 1) * framesize)], padbot)
00300
00301
                  frames[n] = thisframe
00302
         else:
00303
             for n in range(nframes):
00304
                  frames[n] = allframes[n * framesize : (n + 1) * framesize]
         # self.clearStatus()
00305
00306
          parsed = camassem.sensor.parseReadoff(frames, columns)
00307
          return parsed
00308
00309
```

#### 5.17.2.5 getEnter()

#### Definition at line 58 of file misc.py.

```
00058 def getEnter(text):
00060
          Wait for enter key to be pressed.
00061
00062
          Args:
          text: message asking for keypress
00063
00064
          python, _, _, _, _ = sys.version_info
if python >= 3:
00065
00066
00067
              input(text)
00068
          else:
              raw_input(text)
00070
00071
```

#### 5.17.2.6 loadDumpedData()

```
nsCamera.utils.misc.loadDumpedData (
    filename = "frames.txt",
    path = None,
    filetype = "txt",
    sensor = "daedalus",
```

```
firstframe = None,
               lastframe = None,
               width = None,
               height = None,
               padToFull = None,
               firstrow = None,
               lastrow = None,
               maxwidth = None,
               maxheight = None,
               bytesperpixel = None,
               interlacing = None,
               columns = 1)
_summary_
    Output can be piped directly to saveTiffs:
         pars, frames=loadDumpedData(filename="Dump.npy")
         saveTiffs(pars, frames)
Aras:
    filename (str, optional): _description_. Defaults to "frames.txt".
    path (_type_, optional): _description_. Defaults to None.
    filetype (str, optional): _description_. Defaults to "txt".
    sensor (str, optional): _description_. Defaults to "daedalus".
    firstframe (_type_, optional): _description_. Defaults to None.
    lastframe (_type_, optional): _description_. Defaults to None.
    width (_type_, optional): _description_. Defaults to None.
    height (_type_, optional): _description_. Defaults to None.
    padToFull (_type_, optional): _description_. Defaults to None.
    firstrow (_type_, optional): _description_. Defaults to None.
    lastrow (_type_, optional): _description_. Defaults to None.
    maxwidth (_type_, optional): _description_. Defaults to None.
    maxheight (_type_, optional): _description_. Defaults to None.
    interlacing (_type_, optional): _description_. Defaults to None.
    columns (int, optional): _description_. Defaults to 1.
Returns:
    Tuple (parameter object, list of data frames)
Definition at line 310 of file misc.py.
00327 ):
00328
          """_summary_
00329
00330
              Output can be piped directly to saveTiffs:
00331
                 pars, frames=loadDumpedData(filename="Dump.npy")
00332
                 saveTiffs(pars, frames)
00333
          Args:
00334
             filename (str, optional): _description_. Defaults to "frames.txt".
00335
             path (_type_, optional): _description_. Defaults to None.
              filetype (str, optional): _description_. Defaults to "txt".
00336
00337
              sensor (str, optional): _description_. Defaults to "daedalus".
00338
              firstframe (_type_, optional): _description_. Defaults to None.
00339
              lastframe (_type_, optional): _description_. Defaults to None.
             width (_type_, optional): _description_. Defaults to None.
00340
             height (_type_, optional): _description_. Defaults to None.
00341
00342
             padToFull (_type_, optional): _description_. Defaults to None.
              firstrow (_type_, optional): _description_. Defaults to None.
00343
             lastrow (_type_, optional): _description_. Defaults to None.
00344
             maxwidth (_type_, optional): _description_. Defaults to None. maxheight (_type_, optional): _description_. Defaults to None.
00345
00346
00347
             interlacing (_type_, optional): _description_. Defaults to None.
00348
             columns (int, optional): _description_. Defaults to 1.
00349
00350
00351
          Returns:
             Tuple (parameter object, list of data frames)
00352
00353
```

```
00354
          logging.debug("DEBUG: loadDumpedData")
00355
          if sensor.lower() == "daedalus":
00356
              import nsCamera.sensors.daedalus as snsr
00357
          elif sensor.lower() == "icarus":
00358
              import nsCamera.sensors.icarus as snsr
00359
          elif sensor.lower() == "icarus2":
00360
              import nsCamera.sensors.icarus2 as snsr
00361
00362
              logging.error(
00363
                   "ERROR loadDumpedData: invalid sensor type provided, defaulting to icarus2"
00364
00365
              import nsCamera.sensors.icarus2 as snsr
00366
00367
          def buildEmptyFrames():
00368
              cols = [0] * 512
00369
              frame = np.array([cols] * (lastrow - firstrow + 1))
00370
              return [frame] * (lastframe - firstframe + 1)
00371
00372
          # get defaults from class declarations if not specified as parameter
00373
          if firstframe is None:
00374
              firstframe = snsr.firstframe
00375
          if lastframe is None:
00376
              lastframe = snsr.lastframe
00377
          # TODO: first frame number
00378
          if width is None:
00379
              width = snsr.width
00380
          if height is None:
              height = snsr.height
00381
          if firstrow is None:
00382
00383
              firstrow = snsr.firstrow
00384
          if lastrow is None:
00385
              lastrow = snsr.lastrow
          if maxwidth is None:
00386
              maxwidth = snsr.maxwidth
00387
          if maxheight is None:
00388
              maxheight = snsr.maxheight
00389
00390
          if bytesperpixel is None:
   bytesperpixel = snsr.bytesperpixel
00391
          if interlacing is None:
00392
00393
              interlacing = snsr.interlacing
00394
00395
          parameters = fakeCA(
00396
              sensor,
00397
              firstframe,
00398
              lastframe,
00399
              width.
00400
              height,
00401
              padToFull,
00402
              firstrow,
00403
              lastrow,
00404
              maxwidth,
00405
              maxheight,
00406
              bytesperpixel,
00407
              interlacing,
00408
00409
00410
         if path is None:
00411
              path = os.path.join(os.getcwd())
          datafile = os.path.join(path, filename)
00412
00413
         if filename[-3:].lower() == "txt":
00414
              filetype = "txt"
00415
          elif filename[-3:].lower() == "npy":
00416
              filetype = "npy"
00417
          # TODO: return empty frames if error
00418
00419
          if filetype == "npy":
00420
              expectedlength = (lastframe - firstframe + 1) * (lastrow - firstrow + 1) * width
00421
                  f = np.load(datafile)
00422
00423
                  padding = expectedlength - len(f)
00424
                  if padding:
00425
                      logging.warning(
00426
                           "{logwarn}loadDumpedData: Payload is shorter than expected."
00427
                           " Padding with '0's".format(logwarn=parameters.logwarn)
00428
                      f = np.pad(f, (0, padding), "constant", constant_values=(0))
00429
                  frames = generateFrames(parameters, f, columns)
00430
00431
                  return parameters, frames
00432
00433
              except OSError as err:
00434
                  logging.error(
```

```
00435
                        "{logerr}loadDumpedData: OS error: {err}. Returning empty"
00436
                        " frames.".format(logerr=parameters.logerr, err=err)
00437
00438
                   return parameters, buildEmptyFrames()
00439
               except:
00440
                    logging.error(
00441
                        "{logerr}loadDumpedData: Unexpected error: {err}. Returning empty"
00442
                        " frames.".format(logerr=parameters.logerr, err=str(sys.exc_info()[0]))
00443
00444
                    return parameters, buildEmptyFrames()
00445
           # if filetype is not explicitly npy, try loading as text
00446
           else:
00447
               # Payload size as string implied by provided parameters
00448
               expectedlength = (
00449
                   4 * (lastframe - firstframe + 1) * (lastrow - firstrow + 1) * width
00450
00451
00452
               try:
    f = open(datafile, "r")
00453
                   s = f.read()
00454
00455
00456
                   padding = expectedlength - len(s)
00457
                    if padding:
00458
                        logging.warning(
00459
                             "{logwarn}loadDumpedData: Payload is shorter than expected."
                             " Padding with '0's".format(logwarn=parameters.logwarn)
00460
00461
00462
                        s = s.ljust(expectedlength, "0")
00463
00464
                   frames = generateFrames(parameters, s)
00465
                   return parameters, frames
00466
00467
               except OSError as err:
00468
                   logging.error(
                        "{logerr}loadDumpedData: OS error: {err}. Returning empty"
00469
                        " frames.".format(logerr=parameters.logerr, err=err)
00470
00471
00472
                   return parameters, buildEmptyFrames()
00473
               except ValueError:
00474
                   logging.error(
00475
                        "{logerr}loadDumpedData: Could not convert data to an integer."
"Returning empty frames.".format(logerr=parameters.logerr)
00476
00477
00478
                   return parameters, buildEmptyFrames()
00479
               except:
00480
                   logging.error(
                        "{logerr}loadDumpedData: Unexpected error: {err}. Returning empty"
" frames.".format(logerr=parameters.logerr, err=str(sys.exc_info()[0]))
00481
00482
00483
00484
                    return parameters, buildEmptyFrames()
00485
00486
```

#### 5.17.2.7 makeLogLabels()

nsCamera.utils.misc.makeLogLabels (

```
logtag,
                       label )
Definition at line 39 of file misc.py.
00039 def makeLogLabels(logtag, label):
00040
              if logtag is None:
                    logtag = ""
00041
00042
              logcritbase = "CRITICAL{logtag}: ".format(logtag=logtag)
logerrbase = "ERROR{logtag}: ".format(logtag=logtag)
00043
00044
               logwarnbase = "WARNING{logtag}: ".format(logtag=logtag)
loginfobase = "INFO{logtag}: ".format(logtag=logtag)
00045
00046
               logdebugbase = "DEBUG{logtag}: ".format(logtag=logtag)
00047
00048
               logcrit = "{base}{label}".format(base=logcritbase, label=label)
logerr = "{base}{label}".format(base=logerrbase, label=label)
00049
00050
               logwarn = "{base}{label}".format(base=logwarnbase, label=label)
loginfo = "{base}{label}".format(base=loginfobase, label=label)
00051
00052
```

```
00053    logdebug = "{base}{label}".format(base=logdebugbase, label=label)
00055
00056
00057
logdebug = "{base}{label}".format(base=logdebugbase, label=label)
return logcrit, logerr, logwarn, loginfo, logdebug
00057
```

```
5.17.2.8 partition()
nsCamera.utils.misc.partition (
               self,
               frames,
               columns )
Extracts interlaced frames and divides images by hemispheres. If interlacing does
  not evenly divide the height, remainder lines will be dropped
Args:
    self: during normal operation, is sensor object. During offline
      operation, is the parameter.sensor object returned by loadDumpedData()
    frames: list of full-sized frames
    columns: 1 for single image, 2 for separate hemisphere images
Returns: list of deinterlaced frames
Definition at line 634 of file misc.py.
00634 def partition(self, frames, columns):
00635
00636
         Extracts interlaced frames and divides images by hemispheres. If interlacing does
00637
           not evenly divide the height, remainder lines will be dropped
00638
00639
00640
             self: during normal operation, is sensor object. During offline
00641
               operation, is the parameter.sensor object returned by loadDumpedData()
00642
              frames: list of full-sized frames
00643
             columns: 1 for single image, 2 for separate hemisphere images
00644
00645
         Returns: list of deinterlaced frames
00646
00647
         logging.debug(
00648
              "{logdebug}partition: columns = {columns}, interlacing = {interlacing}".format(
00649
                 logdebug=self.logdebug, columns=columns, interlacing=self.sensor.interlacing
00650
00651
00652
00653
         def unshuffle(frames, ifactor):
00654
             warntrimmed = False
00655
              if self.padToFull:
00656
                 newheight = self.sensor.maxheight // (ifactor + 1)
00657
                 if newheight != (self.sensor.maxheight / (ifactor + 1)):
                     warntrimmed = True
00658
00659
                 newheight = self.sensor.height // (ifactor + 1)
00660
00661
                 if newheight != (self.sensor.height / (ifactor + 1)):
                     warntrimmed = True
00662
00663
00664
             if warntrimmed:
00665
                 logging.warning(
00666
                      "{logwarn} partition: interlacing setting requires dropping of lines to"
                      " maintain consistent frame sizes ".format(logwarn=self.logwarn)
00667
00668
00669
             delaced = []
00670
             for frame in frames:
00671
                 for sub in range(ifactor + 1):
00672
                     current = np.zeros((newheight, self.sensor.width // columns), dtype=int)
```

for line in range (newheight):

delaced.append(current)

resorted = [None] \* len(delaced)

current[line] = frame[(ifactor + 1) \* line + sub]

nframes = self.sensor.lastframe - self.sensor.firstframe + 1

00673

00674

00675

00676

00677

```
00678
              for sub in range(ifactor + 1):
00679
                   for idx, frame in enumerate(frames):
00680
                      resorted[sub * nframes + idx] = delaced[idx * (ifactor + 1) + sub]
00681
00682
00683
          if self.sensor.interlacing[0] != self.sensor.interlacing[1]:
00684
              columns = 2 # true even if not explicitly requested by readoff
00685
          if columns == 1:
00686
              if self.sensor.interlacing == [0, 0]: # don't do anything
                  return frames
00688
              else:
00689
                  return unshuffle(frames, self.sensor.interlacing[0])
00690
          else:
00691
              # reshape frame into the proper shape, then split horizontally
00692
              if self.padToFull:
00693
                  framesab = [
00694
                      np.hsplit(frame.reshape(self.sensor.maxheight, -1), 2)
00695
                       for frame in frames
00696
                  ]
00697
              else:
00698
                  framesab = [
00699
                       np.hsplit(
00700
                           frame.reshape((self.sensor.lastrow - self.sensor.firstrow + 1), -1),
00701
00702
                       )
00703
                       for frame in frames
00704
                  1
00705
              framesa = [hemis[0] for hemis in framesab]
          framesb = [hemis[1] for hemis in framesab]
if self.sensor.interlacing == [0, 0]:
00706
00707
00708
              return framesa + framesb
00709
          else:
00710
              return unshuffle(framesa, self.sensor.interlacing[0]) + unshuffle(
00711
                  framesb, self.sensor.interlacing[1]
00712
00713
00714
00715 """
00716 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00717 LLNL-CODE-838080
00718
00719 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00720 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00721 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00722 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00723 be made under this license. 00724 """
```

#### 5.17.2.9 plotFrames()

00568

00567 def plotFrames(self, frames, index=None):

```
00569
          Plot frame or list of frames as individual graphs.
00570
00571
00572
              self: during normal operation, is cameraAssembler object. During offline
00573
                operation, is the parameter object returned by loadDumpedData()
00574
              frames: numpy array or list of numpy arrays
00575
              index: number to start frame numbering
00576
00577
          Returns:
          Error string
00578
00579
00580
          # logging.info(self.loginfo + "plotFrames: index = " + str(index))
00581
          logging.info(
00582
              "{loginfo}: plotFrames: index = {index}".format(
00583
                  loginfo=self.loginfo, index=index
00584
00585
          err = ""
00586
00587
          if index is None:
00588
              nframe = self.sensor.firstframe
00589
          else:
00590
              nframe = index
00591
00592
          if not isinstance(frames, list):
00593
              frames = [frames]
00594
00595
          # if this is a text string from fast readoff, do the numpy conversion now
00596
          if isinstance(frames[0], str):
00597
              frames = generateFrames(frames)
00598
00599
          framestemp = np.copy(frames)
00600
          for idx, frame in enumerate(framestemp):
              if idx < len(framestemp) / 2:</pre>
00601
00602
                  interlacing = self.sensor.interlacing[0]
00603
              else:
                  interlacing = self.sensor.interlacing[1]
00604
00605
              try:
00606
                  if self.padToFull:
00607
                      frame = np.reshape(
00608
                          frame, (self.sensor.maxheight // (interlacing + 1), -1)
00609
00610
                  else:
00611
                      frame = np.reshape(
00612
                          frame,
00613
00614
                               (self.sensor.lastrow - self.sensor.firstrow + 1)
00615
                               // (interlacing + 1),
00616
                               -1,
00617
00618
00619
              except:
00620
                  err = "{logerr}plotFrames: unable to plot frame".format(logerr=self.logerr)
00621
                  logging.error(err)
00622
              plt.imshow(frame, cmap="gray")
name = "Frame %d" % nframe
00623
00624
              plt.title(name)
00625
              plt.show()
00626
00627
              nframe += 1
00628
         return err
00629
00630
00631 # TODO: separate images for hemispheres with different timing
00632
00633
```

# 5.17.2.10 saveTiffs()

```
Save numpy array or list of numpy arrays or single array to disk as individual
  tiffs, with frame number appended to filename. If used for standalone, use the
  parameter object returned by loadDumpedData() as the first argument
Args:
    self: during normal operation, is cameraAssembler object. During offline
      operation, is the parameter object returned by loadDumpedData()
    frames: numpy array or list of numpy arrays
    path: save path, defaults to './output'
    filename: defaults to 'Frame' followed by frame number
    prefix: prepended to 'filename', defaults to time/date
      (e.g. '160830-124704_')
    index: number to start frame numbering
Returns:
   Error string
Definition at line 487 of file misc.py.
00494 ):
00495
00496
         Save numpy array or list of numpy arrays or single array to disk as individual
00497
           tiffs, with frame number appended to filename. If used for standalone, use the
00498
           parameter object returned by loadDumpedData() as the first argument
00499
00500
             self: during normal operation, is cameraAssembler object. During offline
00501
00502
               operation, is the parameter object returned by loadDumpedData()
00503
              frames: numpy array or list of numpy arrays
             path: save path, defaults to './output'
00504
00505
              filename: defaults to 'Frame' followed by frame number
00506
             prefix: prepended to 'filename', defaults to time/date
00507
               (e.g. '160830-124704_')
00508
             index: number to start frame numbering
00509
00510
         Returns:
         Error string
00511
00512
         # logging.info("INFO" + self.logtag + ": saveTiffs")
00513
00514
         logging.info("{loginfo}: saveTiffs".format(loginfo=self.loginfo))
00515
00516
         if path is None:
00517
             path = os.path.join(os.getcwd(), "output")
00518
         if prefix is None:
00519
             prefix = datetime.now().strftime("%y%m%d-%H%M%S%f")[:-5] + "_"
00520
         if not os.path.exists(path):
00521
             os.makedirs(path)
00522
         if index is None:
00523
             firstnum = self.sensor.firstframe
00524
         else:
00525
             firstnum = index
00526
00527
         # if this is a text string from fast readoff, do the numpy conversion now
00528
         if isinstance(frames[0], str):
00529
              frames = generateFrames(frames)
00530
         framestemp = np.copy(frames)
00531
00532
         if np.issubdtype(type(framestemp[0]), np.number):
00533
              # if type(framestemp[0]) is np.uint16:
00534
              # single frame needs to be a list containing one frame
00535
              framestemp = [framestemp]
         for idx, frame in enumerate(framestemp):
00536
00537
              if idx < len(framestemp) / 2:</pre>
00538
                 interlacing = self.sensor.interlacing[0]
00539
00540
                 interlacing = self.sensor.interlacing[1]
00541
             try:
00542
                  if self.padToFull:
00543
                     frame = np.reshape(
00544
                        frame, (self.sensor.maxheight // (interlacing + 1), -1)
00545
00546
                  else:
00547
                     frame = np.reshape(
00548
                         frame,
00549
00550
                              (self.sensor.lastrow - self.sensor.firstrow + 1)
00551
                             // (interlacing + 1).
```

```
00552
00553
                          ),
00554
00555
                  frameimg = Image.fromarray(frame, "I;16")
00556
                  namenum = filename + "_%d" % firstnum
                  tifpath = os.path.join(path, prefix + namenum + ".tif")
00557
00558
                  frameimg.save(tifpath)
00559
                 firstnum += 1
00560
             except Exception:
                 err = "saveTiffs: unable to save images"
                  # logging.error("ERROR" + self.logtag + ": " + err)
00562
00563
                  logging.error("{logerr}: {err}".format(logerr=self.logerr, err=err))
00564
         return err
00565
00566
```

#### 5.17.2.11 str2bytes()

# Definition at line 90 of file misc.py.

```
00090 def str2bytes(astring):
00091
00092
          Python-version-agnostic converter of hexadecimal strings to bytes
00093
00094
00095
              astring: hexadecimal string without '0x'
00096
00097
          byte string equivalent to input string
00098
00099
00100
          python, _, _, _, _ = sys.version_info
if python >= 3:
00101
00102
00103
00104
                   dbytes = binascii.a2b_hex(astring)
00105
               except:
00106
                   logging.error(
                        "ERROR: str2bytes: invalid input: '{astring}'; returning zero"
00107
                       " byte".format(astring=astring)
00108
00109
00110
                   dbytes = b" \x00"
00111
          else:
00112
              try:
00113
                   dbytes = astring.decode("hex")
00114
               except:
00115
                   logging.error(
                       "ERROR: str2bytes: invalid input: '{astring}'; returning zero "
" byte".format(astring=astring)
00116
00117
00118
                   dbytes = b" \x00"
00119
          return dbytes
00120
00121
00122
```

### 5.17.2.12 str2nparray()

```
nsCamera.utils.misc.str2nparray (
                valstring )
Convert string into array of uint16s
    valstring: string of hexadecimal characters
Returns:
    numpy array of uint16
Definition at line 147 of file misc.py.
00147 def str2nparray(valstring):
00149
          Convert string into array of uint16s
00150
00151
00152
              valstring: string of hexadecimal characters
00153
00154
          Returns:
          numpy array of uint16
00155
00156
         if not isinstance(valstring, str):
00158
              logging.error(
                   "ERROR: str2nparray: Invalid input: {valstring} is not a string. Returning"
00159
                  " an empty array".format(valstring=valstring)
00160
00161
00162
              return np.array([])
00163
         stringlen = len(valstring)
         arraylen = int(string)
arraylen = int(stringlen / 4)
outarray = np.empty(int(arraylen), dtype="uint16")
00164
00165
00166
00167
          for i in range(0, arraylen):
00168
                 outarray[i] = int(valstring[4 * i : 4 * i + 4], 16)
00169
00170
              except ValueError:
00171
                  logging.error(
00172
                       "ERROR: str2nparray: input string does not represent a hexadecimal"
00173
                       " integer. Returning an empty array"
00174
00175
                  return np.array([])
00176
         return outarray
00177
```

# 5.18 nsCamera.utils.Packet Namespace Reference

### Classes

00178

class Packet

# 5.18.1 Detailed Description

```
Packet object for communication with boards

Author: Brad Funsten (funstenl@llnl.gov)

Author: Jeremy Hill (hill35@llnl.gov)

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
```

LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL. 'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)

# 5.19 nsCamera.utils.Subregister Namespace Reference

#### Classes

class SubRegister

# 5.19.1 Detailed Description

Subregister object represents a subset of a full register

Author: Matthew Dayton (dayton5@llnl.gov)
Author: Jeremy Martin Hill (jerhill@llnl.gov)

Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved. LLNL-CODE-838080

This work was produced at the Lawrence Livermore National Laboratory (LLNL) under contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL. 'nsCamera' is distributed under the terms of the MIT license. All new contributions must be made under this license.

Version: 2.1.2 (February 2025)

# **Chapter 6**

# **Class Documentation**

# 6.1 nsCamera.CameraAssembler.CameraAssembler Class Reference

#### **Public Member Functions**

- \_\_init\_\_ (self, boardname="llnl\_v4", commname="GigE", sensorname="icarus2", verbose=4, port=None, ip=None, logfile=None, logtag=None, timeout=30)
- initBoard (self)

Aliases to other objects' methods TODO: properly delegate these methods.

- initPots (self)
- · latchPots (self)
- initSensor (self)
- configADCs (self)
- disarm (self)
- startCapture (self, mode)
- readSRAM (self)
- waitForSRAM (self, timeout=None)
- getTimer (self)
- resetTimer (self)
- enableLED (self, status=1)
- setLED (self, LED=1, status=1)
- setPowerSave (self, status=1)
- setPPER (self, pollperiod=None)
- getTemp (self, scale=None)
- getPressure (self, offset=None, sensitivity=None, units=None)
- clearStatus (self)
- checkStatus (self)
- checkStatus2 (self)
- reportStatus (self)
- reportEdgeDetects (self)
- dumpStatus (self)
- checkSensorVoltStat (self)
- setTiming (self, side="AB", sequence=None, delay=0)
- setArbTiming (self, side="AB", sequence=None)

40 Class Documentation

- getTiming (self, side=None, actual=None)
- setManualShutters (self, timing=None)
- setManualTiming (self, timing=None)
- getManualTiming (self)
- getSensTemp (self, scale=None, offset=None, slope=None, dec=1)
- sensorSpecific (self)
- selectOscillator (self, osc=None)
- setInterlacing (self, ifactor=None, side=None)
- setHighFullWell (self, flag=True)
- setZeroDeadTime (self, flag=True, side=None)
- setTriggerDelay (self, delay=0)
- setPhiDelay (self, side=None, delay=0)
- setExtClk (self, dilation=None, frequency=None)
- parseReadoff (self, frames, columns=1)
- sendCMD (self, pkt)
- arm (self, mode=None)
- readFrames (self, waitOnSRAM=None, timeout=0, fast=False, columns=1)
- readoff (self, waitOnSRAM=None, timeout=0, fast=None, columns=1)
- writeSerial (self, cmd, timeout=None)
- readSerial (self, size, timeout=None)
- closeDevice (self)
- saveTiffs (self, frames, path=None, filename="Frame", prefix=None, index=None)
- plotFrames (self, frames, index=None)
- getEnter (self, text)
- checkCRC (self, rval)
- str2bytes (self, astring)
- bytes2str (self, bytesequence)
- str2nparray (self, valstring)
- flattenlist (self, mylist)
- · partition (self, frames, columns)
- · initialize (self)

#### End aliases.

- · reinitialize (self)
- reboot (self)
- getBoardInfo (self)
- getRegister (self, regname)
- setRegister (self, regname, regval)
- resolveSubreg (self, srname)
- getSubregister (self, subregname)
- setSubregister (self, subregname, valstring)
- submitMessages (self, messages, errorstring="Error")
- getPot (self, potname, errflag=False)
- setPot (self, potname, value=1.0, errflag=False)
- getPotV (self, potname, errflag=False)
- setPotV (self, potname, voltage, tune=False, accuracy=0.01, iterations=20, approach=0.75, errflag=False)
- getMonV (self, monname, errflag=False)
- readImgs (self, waitOnSRAM=True, mode="Hardware")
- saveFrames (self, frames, path=None, filename="frames", prefix=None)
- saveNumpys (self, frames, path=None, filename="Frame", prefix=None, index=None)
- dumpNumpy (self, datastream, path=None, filename="Dump", prefix=None)

- checkRegSet (self, regname, teststring)
- initPowerCheck (self)
- powerCheck (self, delta=10)
- printBoardInfo (self)
- dumpRegisters (self)
- dumpSubregisters (self)
- mmReadoff (self, waitOnSRAM, variation=None)
- setFrames (self, minframe=None, maxframe=None)
- setRows (self, minrow=0, maxrow=None, padToFull=False)
- abortReadoff (self, flag=True)
- batchAcquire (self, sets=1, trig="Hardware", path=None, filename="Frame", prefix=None, showProgress=0)
- saveHDF (self, frames, path=None, filename="Acquisition", prefix=None)

#### **Public Attributes**

- version
- currtime
- oldtime
- · trigtime
- waited
- read
- unstringed
- · parsedtime
- · savetime
- cycle
- boardname
- timeout
- commname
- sensorname

For regular version.

- verbose
- port
- python
- pyth1
- pyth2
- PY3
- · platform
- arch
- FPGAVersion
- FPGANum
- FPGAboardtype
- FPGArad
- FPGAsensor
- FPGAinterfaces
- FPGAinvalid
- · iplist
- packageroot
- armed
- senstiming
- · sensmanual

42 Class Documentation

- inittime
- padToFull
- · abort
- verbmap
- logtag
- logcritbase
- logerrbase
- logwarnbase
- · loginfobase
- logdebugbase
- logcrit
- logerr
- logwarn
- · loginfo
- logdebug
- verblevel
- payloaderror
- sensor
- · comms
- board

# 6.1.1 Detailed Description

```
Code to assemble correct code to manage FPGA, frame grabber, and sensor
Exposed methods:
    initialize() - initializes board registers and pots, sets up sensor
    reinitialize() - initialize board and sensors, restore last known timer settings
    reboot() - perform software reset of board and reinitialize
    getBoardInfo() - parses FPGA_NUM register to retrieve board description
   getRegister(regname) - retrieves contents of named register
    setRegister(regname, string) - sets named register to given value
    getSubregister(subregname) - return substring of register identified in board
      attribute 'subregname'
    setSubregister(subregname, valstring) - replace substring of register identified
     in board attribute 'subregname' with 'valstring'
    submitMessages (messages) - set registers or subregisters based on list of
     destination/payload tuples
    getPot(potname) - returns float (0 < value < 1) corresponding to integer stored
     in pot or monitor 'potname'
    setPot(potname, value) - 0 < value < 1; sets named pot to fixed-point number =
      'value' * (maximum pot value)
    getPotV(potname) - returns voltage setting of 'potname'
    setPotV(potname, voltage) - sets named pot to voltage
    getMonV(monname) - returns voltage read by monitor 'monname' (or monitor
      associated with given potname)
    readImgs() - calls arm() and readoff() functions
    saveFrames(frames) - save image object as one file
    saveNumpys(frames) - save individual frames as numpy data files
    dumpNumpy(datastream) - save datastream string to numpy file
    powerCheck(delta) - check that board power has not failed
    printBoardInfo() - print board information derived from FPGA_NUM register
    dumpRegisters() - return contents of all board registers
    dumpSubregisters() - return contents of all named board subregisters
    setFrames(min, max) - select subset of frames for readoff
    setRows(min, max, padToFull) - select subset of rows for readoff
    abortReadoff() - cancel readoff in wait-for-SRAM loop
    batchAquire() - fast acquire a finite series of images
```

```
loadTextFrames() - load data sets previously saved as text and convert to frames
Includes aliases to board- and sensor- specific methods:
   Board methods
       disarm() - take camera out of waiting-for-trigger state
       clearStatus() - clear contents of status registers
       checkStatus() - print contents of status register as reversed bit string
       checkStatus2() - print contents of status register 2 as reversed bit string
       reportStatus() - print report on contents of status registers
       resetTimer() - reset on-board timer
       getTimer() - read on-board timer
       enableLED(status) - enable (default) or disable (status = 0) on-board LEDs
       setLED(LED#, status) - turn LED on (default) or off (status = 0)
       setPowerSave(status) - turn powersave functionality on (default) or off
          (status = 0)
       getTemp() - read on-board temperature sensor
       getPressure() - read on-board pressure sensor
        dumpStatus() - generate dictionary of status, register, and subregister
          contents
   Sensor methods
       checkSensorVoltStat() - checks that jumper settings match sensor selection
       setTiming(side, sequencetuple, delay) - configure high-speed timing
       setArbTiming(side, sequencelist) - configure arbitrary high-speed timing
         sequence
        getTiming(side) - returns high speed timing settings from registers
       setManualTiming() - configures manual shutter timing
       getManualTiming() - returns manual shutter settings from registers
       selectOscillator(osc) - select timing oscillator
       setInterlacing(ifactor) - sets interlacing factor
       setHighFullWell(flag) - controls High Full Well mode
       setZeroDeadTime(flag, side) - controls Zero Dead Time mode
       setTriggerDelay(delayblocks) - sets trigger delay
   Comms methods
       sendCMD(pkt) - sends packet object via serial port
       arm() - configures software buffers & arms camera
       readFrames() - waits for data ready flag, then downloads image data
        readoff() - waits for data ready flag, then downloads image data
       \verb|closeDevice()| - \verb|disconnect| interface | \verb|and release | resources|
   Miscellaneous functions (bare functions that can be called as methods)
       saveTiffs(frames) - save individual frames as tiffs
       plotFrames(frames) - plot individual frames as tiffs
Informational class variables:
   version - nsCamera software version
   FPGAVersion - firmware version (date)
   FPGANum - firmware implementation identifier
   FPGAboardtype - FPGA self-identified board type (should match 'boardname')
   FPGArad = Flag indicating radiation-tolerant FPGA build
   FPGAsensor = FPGA self-identified sensor family (should correspond to
     'sensorname')
    FPGAinterfaces = FPGA self-identified interfaces (list should include
     'commname')
   FPGAinvalid = flag indicating invalid FPGA information in register
```

Definition at line 59 of file CameraAssembler.py.

# 6.1.2 Constructor & Destructor Documentation

```
6.1.2.1 __init__()

nsCamera.CameraAssembler.CameraAssembler.__init__ (

self,
```

44 Class Documentation

```
boardname = "llnl_v4",
              commname = "GigE",
              sensorname = "icarus2",
              verbose = 4,
              port = None,
              ip = None,
              logfile = None,
              logtag = None,
              timeout = 30)
Args:
   boardname: name of FPGA board: llnl_v1, llnl_v4
    commname: name of communication interface: rs422, gige
    sensorname: name of sensor: icarus, icarus2, daedalus
    verbose: optional, sets logging level
        0: print no logging messages
        1: print CRITICAL logging messages (camera will not operate, e.g.,
          unable to connect to board)
        2: print ERROR logging messages (camera will not operate as directed,
          e.g., an attempt to set the timing mode has failed, but the camera
          is still operational)
        3: print WARNING logging messages (camera will operate as directed, but
          perhaps not as expected, e.g., ca.setTiming('A', (9, 8), 1) may be
          programmed correctly, but the actual timing generated by the board
          will be {1} [9, 8, 9, 14, 9, 8, 9])
        4: print INFO logging messages (operational messages from ordinary
          camera operation)
        5. print DEBUG logging messages (detailed messages describing specific
          operations and messages)
    port: optional integer
        When using RS422, this preselects the comport for RS422 and bypasses
          port search
        When using GigE, this preselects the OrangeTree control port for GigE
          (ignored if an ip parameter is not also provided)
    ip: optional string (e.g., '192.168.1.100')
        GigE: bypasses network search and selects particular OrangeTree board -
          required for some operating systems
    logfile: optional string, name of file to divert console output
    timeout: timeout in seconds for connecting using Gigabit Ethernet
```

### Definition at line 151 of file CameraAssembler.py.

```
00162
00163
00164
00165
                  boardname: name of FPGA board: llnl_v1, llnl_v4
00166
                  commname: name of communication interface: rs422, gige
                  sensorname: name of sensor: icarus, icarus2, daedalus
                  verbose: optional, sets logging level
00169
                      0: print no logging messages
00170
                      1: print CRITICAL logging messages (camera will not operate, e.g.,
00171
                        unable to connect to board)
00172
                      2: print ERROR logging messages (camera will not operate as directed,
00173
                        e.g., an attempt to set the timing mode has failed, but the camera
00174
                        is still operational)
00175
                      3: print WARNING logging messages (camera will operate as directed, but
00176
                        perhaps not as expected, e.g., ca.setTiming('A', (9, 8), 1) may be
                        programmed correctly, but the actual timing generated by the board
00177
00178
                        will be {1} [9, 8, 9, 14, 9, 8, 9])
00179
                      4: print INFO logging messages (operational messages from ordinary
00180
                        camera operation)
00181
                      5. print DEBUG logging messages (detailed messages describing specific
00182
                        operations and messages)
                  port: optional integer
00183
00184
                      When using RS422, this preselects the comport for RS422 and bypasses
00185
                        port search
00186
                      When using \mathrm{GigE}, this preselects the OrangeTree control port for \mathrm{GigE}
00187
                        (ignored if an ip parameter is not also provided)
00188
                  ip: optional string (e.g., '192.168.1.100')
```

```
00189
                       GigE: bypasses network search and selects particular OrangeTree board -
00190
                          required for some operating systems
00191
                   logfile: optional string, name of file to divert console output
00192
                   timeout: timeout in seconds for connecting using Gigabit Ethernet
00193
               self.version = "2.1.2"
00194
00195
               self.currtime = 0
00196
               self.oldtime = 0
00197
               self.trigtime = []
00198
               self.waited = []
               self.read = []
00199
00200
               self.unstringed = []
00201
               self.parsedtime = []
               self.savetime = []
00203
               self.cycle = []
00204
               self.boardname = boardname.lower()
00205
               self.timeout = timeout
               # TODO: parse boardname, etc. in separate method
00206
              if self.boardname in ["llnlv1", "v1", "1", 1]:
self.boardname = "llnl_v1"
00207
00208
               if self.boardname in ["llnlv4", "v4", "4", 4]:
    self.boardname = "llnl_v4"
00209
00210
00211
               self.commname = commname.lower()
00212
               if self.commname[0] == "q" or self.commname[0] == "e":
00213
                   self.commname = "gige"
               if self.commname[0] == "r":
00214
                   self.commname = "rs422"
00215
               self.sensorname = sensorname.lower()
00216
               if self.sensorname in ["i1", "ic1", "icarus1"]:
    self.sensorname = "icarus"
00217
00218
               if self.sensorname in ["i2", "ic2"]:
    self.sensorname = "icarus2"
00219
00220
               if self.sensorname == "d":
00221
                   self.sensorname = "daedalus"
00222
               self.verbose = int(verbose)
00223
00224
               self.port = port
               self.python, self.pyth1, self.pyth2, _, _ = sys.version_info
00225
00226
               self.PY3 = self.python >= 3
00227
               self.platform = platform.system()
00228
               self.arch, _ = platform.architecture()
00229
               self.FPGAVersion = ""
00230
               self.FPGANum = ""
00231
00232
               # FPGA information here and below populated during initialization using
00233
                   getBoardInfo
               self.FPGAboardtype = ""
00234
00235
               self.FPGArad = False
00236
               self.FPGAsensor = ""
00237
               self.FPGAinterfaces = []
00238
00239
               # indicates invalid FPGA information in register# (0x80000001 accepted as valid)
00240
               self.FPGAinvalid = False
00241
00242
               self.iplist = None
00243
               self.packageroot = os.path.dirname(inspect.getfile(CameraAssembler))
00244
               self.armed = False
00245
00246
               \# only one of these collections (senstiming, sensmanual) should be nonempty at
00247
00248
               self.senstiming = {} # preserve HST setting against possible power failure
               self.sensmanual = [] # preserve manual timing
00249
00250
               self.inittime = 0
00251
               self.padToFull = False
00252
               self.abort = False
00253
00254
               self.verbmap = {
00255
                   0: 99,
00256
                   1: logging.CRITICAL,
00257
                   2: logging.ERROR,
00258
                   3: logging.WARNING,
00259
                   4: logging.INFO,
00260
                   5: logging.DEBUG,
00261
00262
               if logtag is None:
00263
                   logtag = ""
00264
               self.logtag = logtag
00265
               self.logcritbase = "CRITICAL {logtag}: ".format(logtag=logtag)
00266
               self.logerrbase = "ERROR {logtag}: ".format(logtag=logtag)
00267
               self.logwarnbase = "WARNING {logtag}: ".format(logtag-logtag)
self.loginfobase = "INFO {logtag}: ".format(logtag=logtag)
00268
00269
```

46 Class Documentation

```
00270
               self.logdebugbase = "DEBUG {logtag}: ".format(logtag=logtag)
00271
00272
               self.logcrit = "{lb}[CA]".format(lb=self.logcritbase)
               self.logerr = "{lb}[CA]".format(lb=self.logerrbase)
00273
              self.logwarn = "{1b}[CA]".format(lb=self.logwarnbase)
self.loginfo = "{1b}[CA]".format(lb=self.loginfobase)
00274
00275
00276
              self.logdebug = "{lb}[CA]".format(lb=self.logdebugbase)
00277
00278
              self.verblevel = self.verbmap.get(verbose, 5) # defaults to 5 for invalid entry
00279
00280
               if logfile:
00281
                   logging.basicConfig(format="%(message)s", filename=logfile)
00282
              else:
                   logging.basicConfig(format="%(message)s")
00284
               logging.getLogger().setLevel(self.verblevel)
00285
               logging.getLogger("matplotlib.font_manager").disabled = True
00286
              logging.debug(
00287
                   "{logdebug}CameraAssembler: boardname = {boardname}; commname = {commname};"
                   " sensorname = {sensorname}; verbose = {verbose}; port = {port}; ip = {ip};"
00288
                   " logfile = {logfile}; logtag = {logtag}".format(
00289
00290
                       logdebug=self.logdebug,
00291
                       boardname=boardname,
00292
                       commname=commname,
00293
                       sensorname=sensorname,
00294
                       verbose=verbose,
                       port=port,
00295
00296
                       ip=ip,
00297
                       logfile=logfile,
00298
                       logtag=logtag,
00299
00300
              )
00301
00302
              if ip:
00303
00304
                       iphex = socket.inet_aton(ip)
                   except socket.error:
00305
00306
                       logging.critical(
00307
                            "{logcrit}CameraAssembler: invalid IP provided".format(
00308
                               logcrit=self.logcrit
00309
00310
00311
                       sys.exit(1)
                   ipnum = [0, 0, 0, 0]
00312
00313
                   for i in range(4):
00314
                       if self.PY3:
00315
                           ipnum[i] = iphex[i]
00316
00317
                           ipnum[i] = int(iphex[i].encode("hex"), 16)
00318
                   self.iplist = ipnum
00319
00320
              self.payloaderror = False
00321
00322
               # code pulled out of __init__ to facilitate reinitialization of the board
00323
                  without needing to instantiate a new CameraAssembler object
00324
              self.initialize()
00325
```

#### 6.1.3 Member Function Documentation

#### 6.1.3.1 abortReadoff()

#### Definition at line 1989 of file CameraAssembler.py.

```
def abortReadoff(self, flag=True):
01990
01991
              Simple abort command for readoff in waiting mode--does not interrupt download in
01992
                 progress. Requires external threading to function. *WARNING* if not
01993
                 intercepted by active readoff command, will terminate next readoff command
01994
                 immediately at inception.
01995
              Args:
01996
                  flag: Sets passive abort flag read by readoff command
01997
              Returns:
              boolean: updated setting of flag
01998
01999
02000
              logging.info(self.loginfo + "abortReadoff")
              self.abort = flag
02002
              return flag
02003
```

#### 6.1.3.2 arm()

# Definition at line 452 of file CameraAssembler.py.

```
00452 def arm(self, mode=None):
00453 return self.comms.arm(mode)
00454
```

#### 6.1.3.3 batchAcquire()

```
nsCamera.CameraAssembler.CameraAssembler.batchAcquire (
              self,
              sets = 1,
              trig = "Hardware",
              path = None,
              filename = "Frame",
              prefix = None,
              showProgress = 0 )
Acquire a series of images as fast as possible, then process and save to disk.
*WARNING* This method stores images in RAM, so the number of sets that can be
  acquired in a single call is limited by available memory.
Aras:
    sets: Number of acquisitions to perform
   trig: trigger type; 'hardware', 'software', or 'dual'
   path: save path, defaults to './output'
    filename: defaults to 'frames.bin'
   prefix: prepended to filename, defaults to time/date (e.g. '160830-124704_')
      DO NOT USE unless providing a varying value (a fixed prefix will cause
      overwriting)
    showProgress: if non-zero, show notice every 'showProgress' acquisitions and
      print total acquisition time
Returns:
   Time taken for acquisition (seconds)
```

48 Class Documentation

### Definition at line 2004 of file CameraAssembler.py.

```
02012
02013
02014
              Acquire a series of images as fast as possible, then process and save to disk.
02015
              *WARNING* This method stores images in RAM, so the number of sets that can be
02016
                acquired in a single call is limited by available memory.
02017
02018
02019
                  sets: Number of acquisitions to perform
                  trig: trigger type; 'hardware', 'software', or 'dual'
02020
02021
                  path: save path, defaults to './output'
                  filename: defaults to 'frames.bin'
02022
02023
                  prefix: prepended to filename, defaults to time/date (e.g. '160830-124704_')
02024
                    DO NOT USE unless providing a varying value (a fixed prefix will cause
02025
                    overwriting)
02026
                  showProgress: if non-zero, show notice every 'showProgress' acquisitions and
02027
                    print total acquisition time
02028
02029
              Returns:
              Time taken for acquisition (seconds)
02030
02031
02032
              logging.debug(
02033
                  self.logdebug
02034
                  + "batchAcquire: sets = "
02035
                  + str(sets)
                  + "; trig =
02036
02037
                  + str(trig)
                  + "; path = "
02038
02039
                  + str(path)
                  + "; filename = "
02040
02041
                  + str(filename)
                  + "; prefix =
02042
02043
                  + str(prefix)
                  + "; showProgress = "
02044
02045
                  + str(showProgress)
02046
              datalist = ["0"] * sets
timelist = [datetime.now()] * sets
02047
02048
02049
              logging.info(
02050
                  self.loginfo
02051
                  + "batchAcquire: temporarily disabling warning and information logging "
02052
02053
              logging.getLogger().setLevel(self.verbmap.get(2))
02054
              beforeread = time.time()
02055
              for i in range(sets):
02056
                  if showProgress and not (i + 1) % showProgress:
02057
                      print(self.loginfo + "batchAcquire: Acquiring set " + str(i + 1))
02058
                  self.arm(trig)
02059
                  data, datalen, data_err = self.readoff(fast=True)
                  datalist[i] = data
timelist[i] = datetime.now()
02060
02061
02062
              afterread = time.time()
02063
              if showProgress:
02064
                  print (
02065
                      self.loginfo
02066
                       + "batchAcquire: "
02067
                       + str(afterread - beforeread)
02068
                       + " seconds for "
02069
                       + str(sets)
02070
                       + " sets"
02071
02072
              setnum = 0
02073
              if path is None:
02074
                  path = os.path.join(os.getcwd(), "output")
02075
              for imset, imtime in zip(datalist, timelist):
02076
                  setnum = setnum + 1
02077
                  if showProgress and not setnum % showProgress:
                      print(self.loginfo + "batchAcquire: Saving set " + str(setnum))
02078
02079
                  parsed = generateFrames(self, imset)
02080
                  if prefix is None:
02081
                      setprefix = imtime.strftime("%y%m%d-%H%M%S%f")[:-2] + "_"
02082
                  else:
02083
                      setprefix = prefix
02084
                  self.saveTiffs(parsed, path, filename, prefix=setprefix)
              logging.getLogger().setLevel(self.verblevel)
02085
02086
              logging.info(self.loginfo + "batchAcquire: re-enabling logging")
              return afterread - beforeread
02087
02088
```

# 6.1.3.4 bytes2str()

# 6.1.3.5 checkCRC()

```
nsCamera.CameraAssembler.CameraAssembler.checkCRC ( self, \\ rval \ )
```

## Definition at line 480 of file CameraAssembler.py.

```
00480 def checkCRC(self, rval):
00481 return checkCRC(rval)
00482
```

#### 6.1.3.6 checkRegSet()

## Definition at line 1662 of file CameraAssembler.py.

```
def checkRegSet(self, regname, teststring):
01663
01664
               Quick check to confirm that data read from register matches data write
01665
01666
01667
                   regname: register to test
01668
                   teststring: value to assign to register, as integer or hexadecimal string
01669
                    with or without '0x'
01670
01671
               Returns:
              boolean, True if read and write values match \tt muu
01672
01673
01674
              self.setRegister(regname, teststring)
              # tell board to send data; wait to clear before interrogating register contents
if regname == "SRAM_CTL":
01675
01676
                   time.sleep(2)
01677
01678
                   if self.commname == "rs422":
01679
                       logging.info(
                           self.loginfo + "skipping 'SRAM_CTL' register check for RS422"
01680
```

```
01681
01682
                       return True
01683
               else:
01684
                  time.sleep(0.1)
               temp = self.getRegister(regname)
resp = temp[1].upper()
01685
01686
01687
               if resp != teststring.upper():
01688
                  logging.error(
01689
                       self.logerr
01690
                        + "checkRegSet failure: "
01691
                       + regname
01692
                        + " ; set: "
01693
                        + teststring
01694
                        + " ; read:
01695
                        + resp
01696
01697
                   return False
01698
              return True
01699
```

## 6.1.3.7 checkSensorVoltStat()

```
{\tt nsCameraAssembler.CameraAssembler.checkSensorVoltStat} \ \ ( {\tt self} \ )
```

#### Definition at line 398 of file CameraAssembler.py.

```
00398 def checkSensorVoltStat(self):
00399 return self.sensor.checkSensorVoltStat()
00400
```

#### 6.1.3.8 checkStatus()

```
{\tt nsCameraAssembler.CameraAssembler.checkStatus} \ \ ( {\tt self} \ )
```

#### Definition at line 383 of file CameraAssembler.py.

```
00383 def checkStatus(self):
00384 return self.board.checkStatus()
00385
```

# 6.1.3.9 checkStatus2()

```
{\tt nsCameraAssembler.CameraAssembler.checkStatus2} \ \ ( {\tt self} \ )
```

#### Definition at line 386 of file CameraAssembler.py.

```
00386 def checkStatus2(self):
00387 return self.board.checkStatus2()
00388
```

#### 6.1.3.10 clearStatus()

```
{\tt nsCameraAssembler.CameraAssembler.clearStatus} \ \ ( {\tt self} \ )
```

## Definition at line 380 of file CameraAssembler.py.

```
00380 def clearStatus(self):
00381 return self.board.clearStatus()
00382
```

# 6.1.3.11 closeDevice()

```
nsCamera.CameraAssembler.CameraAssembler.closeDevice (
              self )
Definition at line 468 of file CameraAssembler.py.
     def closeDevice(self):
            return self.comms.closeDevice()
00470
6.1.3.12 configADCs()
nsCamera.CameraAssembler.CameraAssembler.configADCs (
              self )
Definition at line 341 of file CameraAssembler.py.
       def configADCs(self):
00342
          return self.board.configADCs()
00343
6.1.3.13 disarm()
nsCamera.CameraAssembler.CameraAssembler.disarm (
              self )
Definition at line 344 of file CameraAssembler.py.
        def disarm(self):
00344
00345
            return self.board.disarm()
00346
6.1.3.14 dumpNumpy()
nsCamera.CameraAssembler.CameraAssembler.dumpNumpy (
              self,
              datastream,
              path = None,
              filename = "Dump",
              prefix = None )
Datastream is converted directly to numpy array and saved to disk. No attempt to
  parse headers or separate into individual frames is made. The packet header is
  removed before saving
Args:
    datastream: string to be saved
    path: save path, defaults to './output'
    filename: defaults to 'Dump'
```

prefix: prepended to 'filename', defaults to time/date

Error string

Returns:

(e.g. '160830-124704\_')

Definition at line 1610 of file CameraAssembler.py.

```
01617
01618
              Datastream is converted directly to numpy array and saved to disk. No attempt to
01619
                parse headers or separate into individual frames is made. The packet header is
01620
                removed before saving
01621
01622
              Args:
01623
                  datastream: string to be saved
01624
                  path: save path, defaults to './output'
                  filename: defaults to 'Dump' prefix: prepended to 'filename', defaults to time/date
01625
01626
01627
                    (e.g. '160830-124704_')
01628
01629
              Returns:
              Error string
01630
01631
01632
              logging.info(self.loginfo + "dumpNumpy")
01633
              logging.debug(
01634
                  self.logdebug
                   + "dumpNumpy: path = "
01635
01636
                  + str(path)
01637
                  + "; filename =
01638
                  + str(filename)
01639
                  + "; prefix =
01640
                  + str(prefix)
01641
              )
              err = ""
01642
              if path is None:
    path = os.path.join(os.getcwd(), "output")
01643
01644
01645
              if prefix is None:
                  prefix = time.strftime("%y%m%d-%H%M%S_", time.localtime())
01646
              if not os.path.exists(path):
01647
01648
                  os.makedirs(path)
              npdata = str2nparray(datastream[36:])
01649
01650
              try:
01651
                  nppath = os.path.join(path, prefix + filename + ".npy")
01652
                  np.save(nppath, npdata)
01653
              except SystemExit:
01654
01655
              except KeyboardInterrupt:
01656
                  raise
01657
              except Exception:
01658
                  err = self.logerr + "dumpNumpy: unable to save data stream"
01659
                  logging.error(err)
01660
              return err
01661
```

#### 6.1.3.15 dumpRegisters()

01771

Returns:

```
nsCamera.CameraAssembler.CameraAssembler.dumpRegisters (
              self )
*DEPRECATED* use dumpStatus() instead
List contents of all registers in board.registers. *WARNING* some status flags
  will reset when read.
Returns:
    Sorted list: [register name (register address) : register contents as
      hexadecimal string without '0x']
Definition at line 1764 of file CameraAssembler.py.
01764
         def dumpRegisters(self):
01765
01766
             *DEPRECATED* use dumpStatus() instead
01767
01768
             List contents of all registers in board.registers. *WARNING* some status flags
01769
              will reset when read.
01770
```

```
01772
                  Sorted list: [register name (register address) : register contents as
01773
                    hexadecimal string without '0x']
01774
              dump = \{\}
01775
01776
              for key in self.board.registers.keys():
01777
                  err, rval = self.getRegister(key)
01778
                  dump[key] = rval
01779
              reglistmax = int(max(self.board.registers.values()), 16)
01780
              dumplist = [0] * (reglistmax + 1)
01781
              for k, v in dump.items():
                  regnum = self.board.registers[k]
01782
01783
                  dumplist[int(regnum, 16)] = (
01784
                      "(" + regnum + ") {0:<24} {1}".format(k, v.upper())
01785
01786
              reglist = [a for a in dumplist if a]
01787
              return reglist
01788
```

## 6.1.3.16 dumpStatus()

```
{\tt nsCameraAssembler.CameraAssembler.dumpStatus} \end{\footnote} \begin{tabular}{ll} {\tt self} \end{\footnote} \begin{tabular}{ll} \begin{tabular}{ll} {\tt self} \end{\footnote} \begin{tabular}{ll} {\tt self} \begin{tabular}{ll} {\tt self} \begin{tabular}{ll} {\tt self} \begin{tabular}{ll} \begin{tabular}{ll} {\tt self} \begin{tabular}{ll} {\tt self} \begin{tabular}{ll} {\tt self} \begin{tabular}{ll} \begin{tabular}{ll} {\tt self} \begin{tabular}{ll} \begin{tabular}{
```

## Definition at line 395 of file CameraAssembler.py.

```
00395 def dumpStatus(self):
00396 return self.board.dumpStatus()
00397
```

## 6.1.3.17 dumpSubregisters()

#### Definition at line 1789 of file CameraAssembler.py.

```
01789
         def dumpSubregisters(self):
01790
01791
              *DEPRECATED* use dumpStatus() instead
01792
01793
              List contents of all subregisters in board.channel_lookups and
01794
                board.monitor_lookups.
01795
              *WARNING* some registers will reset when read; only the first subregister from
01796
               such a register will return the correct value, the remainder will return zeros
01798
              Returns:
01799
                 dictionary {subregister name : subregister contents as binary string
01800
                   without initial '0b'}
01801
01802
              dump = { } { }
01803
              for sub in self.board.subreglist:
01804
                  key = sub.name
01805
                  err, resp = self.getSubregister(key)
01806
                  if err:
01807
                     logging.warning(
                          self.logwarn + "dumpSubregisters: unable to read subregister " + key
01808
01809
                  val = hex(int(resp, 2))
01810
                  dump[key] = val
01811
01812
              return dump
01813
```

# 6.1.3.18 enableLED()

00662

00663

00664 00665

00666

00667 00668 00669 else:

```
nsCamera.CameraAssembler.CameraAssembler.enableLED (
                status = 1)
Definition at line 362 of file CameraAssembler.pv.
          def enableLED(self, status=1):
00362
00363
              return self.board.enableLED(status)
00364
6.1.3.19 flattenlist()
nsCamera.CameraAssembler.CameraAssembler.flattenlist (
                mylist )
Definition at line 492 of file CameraAssembler.py.
         def flattenlist(self, mylist):
00493
            return flattenlist(mylist)
00494
6.1.3.20 getBoardInfo()
\verb|nsCamera.CameraAssembler.CameraAssembler.getBoardInfo| (
                self )
Get board info from FPGA_NUM register. Returns error flag if register contents
  are invalid and tuple (board version number, rad tolerance flag, sensor name)
Returns:
    tuple (errorFlag, (board version, rad tolerance flag, sensor name))
Definition at line 643 of file CameraAssembler.py.
00643
          def getBoardInfo(self):
00644
00645
              Get board info from FPGA\_NUM register. Returns error flag if register contents
00646
               are invalid and tuple (board version number, rad tolerance flag, sensor name)
00647
00648
              Returns:
              tuple (errorFlag, (board version, rad tolerance flag, sensor name))
00649
00650
00651
              invalidFPGANum = False
00652
              interfaces = []
00653
00654
              # TODO: move to new method (combine with parsing from initialize)
00655
              if int(self.FPGANum[0], 16) & 8:
                  if self.FPGANum[1] == "1":
    boardtype = "LLNLv1"
00656
00657
00658
                  elif self.FPGANum[1] == "4":
                     boardtype = "LLNLv4"
00659
00660
                  else:
00661
                     boardtype = "LLNLv?"
```

invalidFPGANum = True

" supported by this software "

self.logwarn + "FPGA self-identifies as SNLrevC, which is not"

boardtype = "SNLrevC"

invalidFPGANum = True

logging.warning(

```
00670
              self.FPGAboardtype = boardtype
00671
00672
              if int(self.FPGANum[6], 16) & 1:
                 rad = True
00673
00674
              else:
00675
                 rad = False
00676
             self.FPGArad = rad
00677
00678
             if self.FPGANum[7] == "1":
00679
                 sensor = "Icarus"
00680
             elif self.FPGANum[7] == "2":
00681
                 sensor = "Daedalus"
00682
             else:
00683
                 sensor = "Undefined"
00684
                  invalidFPGANum = True
00685
             self.FPGAsensor = sensor
00686
             if int(self.FPGANum[5], 16) & 1:
00687
                  interfaces.append("RS422")
00688
00689
              if int(self.FPGANum[5], 16) & 2:
00690
                 interfaces.append("GigE")
00691
             self.FPGAinterfaces = interfaces
00692
00693
             if invalidFPGANum:
00694
                 if self.FPGANum == "80000001":
                     invalidFPGANum = False
00695
00696
                  else:
                     logging.warning(self.logwarn + "FPGA self-identification is invalid")
00697
00698
              self.FPGAinvalid = invalidFPGANum
00699
00700
             return invalidFPGANum, (boardtype, rad, sensor)
00701
```

## 6.1.3.21 getEnter()

```
nsCamera.CameraAssembler.CameraAssembler.getEnter ( self, \\ text \ )
```

## Definition at line 477 of file CameraAssembler.py.

```
00477 def getEnter(self, text):
00478 return getEnter(text)
00479
```

# 6.1.3.22 getManualTiming()

```
nsCamera.CameraAssembler.CameraAssembler.getManualTiming ( self )
```

# Definition at line 416 of file CameraAssembler.py.

```
00416 def getManualTiming(self):
00417 return self.sensor.getManualTiming()
00418
```

# 6.1.3.23 getMonV()

```
Reads voltage from monitor named or that associated with the pot named 'monname'
Aras:
    monname: name of pot or monitor, e.g., VRST or MON_CH2 found in
      board.subreg_aliases or defined in board.subregisters
    errflag: if True, return tuple with error string
Returns:
    if errflag:
         tuple: (error string, float value of voltage measured by monitor)
    else:
         float value of voltage measured by monitor
Definition at line 1396 of file CameraAssembler.pv.
01396
          def getMonV(self, monname, errflag=False):
01397
01398
              Reads voltage from monitor named or that associated with the pot named 'monname'
01399
01400
                  monname: name of pot or monitor, e.g., VRST or MON_CH2 found in
01401
                    board.subreg_aliases or defined in board.subregisters
01402
01403
                  errflag: if True, return tuple with error string
01404
01405
              Returns:
01406
                 if errflag:
01407
                      tuple: (error string, float value of voltage measured by monitor)
01408
                  else:
01409
                      float value of voltage measured by monitor
01410
01411
              logging.debug(
01412
                  self.logdebug
01413
                  + "getMonV: monname = "
01414
                  + str(monname)
01415
                  + "; errflag = "
01416
                  + str(errflag)
01417
01418
              monname = monname.upper()
01419
             if monname in self.board.subreg_aliases:
01420
                  monname = self.board.subreg_aliases[monname].upper()
01421
01422
              for key, value in self.board.monitor_controls.items():
01423
                 if value == monname:
01424
                      monname = key
              if monname not in self.board.monitor_controls:
01425
01426
                  if monname in self.board.subreglist:
01427
                      pass # no change necessary
01428
                  else:
01429
                      err = (
                          self.logerr + "getMonV: invalid lookup " + monname + ", returning 0"
01430
01431
01432
                      logging.error(err)
01433
                      if errflag:
01434
                          return err, 0
                      return 0
01435
              err, monval = self.getPot(monname, errflag=True)
01436
              logging.debug(self.logdebug + "getMonV: monval = " + str(monval))
01437
01438
              if err:
01439
                 logging.error(
01440
                      self.logerr + "getMonV: unable to read monitor value for " + monname
01441
01442
              if self.board.ADC5 bipolar:
                  if monval >= 0.5:
01443
01444
                      monval -= 1 # handle negative measurements (two's complement)
01445
                  if errflag:
                  return err, 2 * self.board.ADC5_mult * monval * self.board.VREF
return 2 * self.board.ADC5_mult * monval * self.board.VREF
01446
01447
01448
              else:
01449
                  if errflag:
                  return err, self.board.ADC5_mult * monval * self.board.VREF
return self.board.ADC5_mult * monval * self.board.VREF
01450
01451
01452
```

## 6.1.3.24 getPot()

nsCamera.CameraAssembler.CameraAssembler.getPot (

```
self,
                potname,
                errflag = False )
Retrieves value of pot or ADC monitor subregister, scaled to [0,1).
Aras:
    potname: name of pot or monitor, e.g., VRST or MON_CH2 found in
       board.subreg_aliases or defined in board.subregisters
     errflag: if True, return tuple with error string
Returns:
    if errflag:
         tuple: (error string, float value of subregister, scaled to [0,1) )
     else:
         float value of subregister, scaled to [0,1)
Definition at line 1013 of file CameraAssembler.py.
01013
          def getPot(self, potname, errflag=False):
01014
01015
              Retrieves value of pot or ADC monitor subregister, scaled to [0,1).
01016
01017
01018
                  potname: name of pot or monitor, e.g., VRST or MON_CH2 found in
01019
                    board.subreg_aliases or defined in board.subregisters
01020
                  errflag: if True, return tuple with error string
01021
01022
01023
                  if errflag:
01024
                      tuple: (error string, float value of subregister, scaled to [0,1) )
01025
01026
                       float value of subregister, scaled to [0,1)
              11 11 11
01027
01028
              logging.debug(
                   "{logdebug}getPot: potname = {potname}; errflag = {errflag}".format(
01029
01030
                       logdebug=self.logdebug, potname=potname, errflag=errflag
01031
01032
              )
01033
01034
              potname, potobj, _ = self.resolveSubreg(potname)
01035
              if not potobj:
01036
                  err = "{logerr}getPot: invalid lookup: {potname}; returning 0".format(
01037
                     logerr=self.logerr, potname=potname
01038
01039
01040
                  logging.error(err)
01041
                  if errflag:
01042
                      return err, 0
01043
                  return O
01044
              err, b_pot_value = self.getSubregister(potname)
01045
              if err:
                  err = "{logerr}getPot: unable to read subregister: {potname}".format(
01046
01047
                       logerr=self.logerr, potname=potname
01048
01049
01050
              \ensuremath{\text{\#}} convert binary string back to decimal
              f_reg_value = 1.0 * int(b_pot_value, 2)
value = (f_reg_value - potobj.min) / (potobj.max - potobj.min)
# logging.debug(self.logdebug + "getpot: value = " + str(value))
01051
01052
01053
01054
01055
              logging.debug(
01056
                   "{logdebug}getpot: value = {value}".format(
01057
                       logdebug=self.logdebug, value=value
01058
01059
              )
01060
01061
              if errflag:
01062
                  return err, value
01063
              return value
```

01064

# 6.1.3.25 getPotV()

```
nsCamera.CameraAssembler.CameraAssembler.getPotV (
                self,
                potname,
                errflag = False )
Reads voltage _setting_ (not actual voltage) of specified pot
Args:
    potname: name of pot or monitor, e.g., VRST or MON_CH2 found in
      board.subreg_aliases or defined in board.subregisters
    errflag: if True, return tuple with error string
Returns:
    if errflag:
         tuple: (error string, float value of pot voltage)
         float value of pot voltage
Definition at line 1155 of file CameraAssembler.py.
          def getPotV(self, potname, errflag=False):
01156
01157
              Reads voltage _setting_ (not actual voltage) of specified pot
01158
01159
01160
                 potname: name of pot or monitor, e.g., VRST or MON_CH2 found in
01161
                   board.subreg_aliases or defined in board.subregisters
01162
                  errflag: if True, return tuple with error string
01163
01164
              Returns:
01165
                 if errflag:
01166
                     tuple: (error string, float value of pot voltage)
01167
                  else:
01168
                      float value of pot voltage
              ....
01169
01170
              logging.debug(
01171
                  self.logdebug
                  + "getPotV: potname = "
01172
01173
                  + str(potname)
                 + "; errflag = "
01174
01175
                  + str(errflag)
01176
01177
              potname, potobj, _ = self.resolveSubreg(potname)
01178
              if not potobj:
01179
                  err = (
                     self.logerr + "getPotV: invalid lookup: " + potname + " , returning 0 "
01180
01181
01182
                  logging.error(err)
01183
                  if errflag:
01184
                      return err, 0
01185
                 return 0
              err, val = self.getPot(potname, errflag=True)
logging.debug(self.logdebug + "getPotV: val = " + str(val))
01186
01187
              if err:
01188
                  logging.error(self.logerr + "getPotV: unable to read pot " + potname)
01189
             minV = potobj.minV
maxV = potobj.maxV
01190
01191
01192
             if errflag:
01193
                  return err, val * (maxV - minV)
01194
             return val * (maxV - minV)
01195
```

# 6.1.3.26 getPressure()

```
{\tt nsCameraAssembler.CameraAssembler.getPressure} \ \ (
```

```
offset = None,
                sensitivity = None,
                units = None)
Definition at line 377 of file CameraAssembler.py.
         def getPressure(self, offset=None, sensitivity=None, units=None):
00378
              return self.board.getPressure(offset, sensitivity, units)
00379
6.1.3.27 getRegister()
\verb|nsCameraAssembler.CameraAssembler.getRegister| (
               self,
               regname )
Retrieves contents of named register as hexadecimal string without '0x'
    regname: name of register as given in ICD
Returns:
    tuple: (error string, register contents as hexadecimal string without '0x')
Definition at line 702 of file CameraAssembler.py.
00702
          def getRegister(self, regname):
00703
00704
              Retrieves contents of named register as hexadecimal string without '0x'
00705
00706
00707
                 regname: name of register as given in ICD
00708
00709
              Returns:
              tuple: (error string, register contents as hexadecimal string without '0x') _{\pi\pi\pi}
00710
00711
00712
              # logging.debug(self.logdebug + "getRegister: regname = " + str(regname))
00713
              logging.debug(
00714
                  "{logdebug}getRegister: regname = {regname}".format(
00715
                     logdebug=self.logdebug, regname=regname
00716
                 )
00717
00718
00719
              regname = regname.upper()
00720
              if regname not in self.board.registers:
00721
                  err = "{logerr}getRegister: Invalid register name: {regname}; returning"
00722
                  " zeros".format(logerr=self.logerr, regname=regname)
00723
                 logging.error(err)
00724
                  return err, "00000000"
00725
              sendpkt = Packet(cmd="1", addr=self.board.registers[regname])
00726
              err, rval = self.comms.sendCMD(sendpkt)
00727
              if err:
00728
                  logging.error(
00729
                      "{logerr}getRegister: {regname}; {err}".format(
00730
                          logerr=self.logerr, regname=regname, err=err
00731
00732
                 )
00733
00734
              retval = rval[8:16]
00735
             logging.debug(
00736
                  "{logdebug}getRegister: retval = {retval}".format(
00737
                     logdebug=self.logdebug, retval=retval
00738
```

)

return err, retval

00739

00740 00741

00742

## 6.1.3.28 getSensTemp()

```
nsCamera.CameraAssembler.CameraAssembler.getSensTemp (
                scale = None,
                offset = None,
                slope = None,
                dec = 1)
Definition at line 419 of file CameraAssembler.py.
          def getSensTemp(self, scale=None, offset=None, slope=None, dec=1):
00420
              return self.sensor.getSensTemp(scale, offset, slope, dec)
00421
6.1.3.29 getSubregister()
nsCamera.CameraAssembler.CameraAssembler.getSubregister (
                self,
                subregname )
Returns substring of register identified in board attribute 'subregname'
Aras:
    subregname: listed in board.subreg_aliases or defined in board.subregisters
Returns:
    tuple: (error string, contents of subregister as binary string without '0b')
Definition at line 836 of file CameraAssembler.py.
00836
          def getSubregister(self, subregname):
00837
00838
              Returns substring of register identified in board attribute 'subregname'
00839
00840
00841
                  subregname: listed in board.subreg_aliases or defined in board.subregisters
00842
00843
              Returns:
              tuple: (error string, contents of subregister as binary string without '0b') """
00844
00845
00846
              logging.debug(
00847
                  "{logdebug}getSubegister: subregname = {subregname}".format(
00848
                      logdebug=self.logdebug,
00849
                      subregname=subregname,
00850
00851
00852
00853
              subregname, subregobj, _ = self.resolveSubreg(subregname)
00854
              if not subregobj:
00855
                  err = "{logerr}getSubregister: invalid lookup: {subregname}; returning"
00856
                  " string of zeroes".format(logerr=self.logerr, subregname=subregname)
00857
00858
                  logging.error(err)
                  return err, "".zfill(8)
00859
00860
              err, resp = self.getRegister(subregobj.register)
00861
              if err:
                  logging.error(
00862
00863
                       "{logerr}getSubregister: unable to retrieve register setting: \setminus
00864
                      {subregname}; returning '0' string".format(
00865
                          logerr=self.logerr, subregname=subregname
00866
                      )
00867
                  )
00868
              return err, "".zfill(8)
hex_str = "0x" + resp # this should be a hexadecimalstring
00869
00870
              b_reg_value = "{0:0=32b}".format(int(hex_str, 16)) # convert to binary string
00871
```

```
00872
              # list indexing is reversed from bit string; the last bit of the string is at
00873
              # index 0 in the list (thus bit 0 is at index 0)
00874
              startindex = 31 - subregobj.start_bit
              retval = b_reg_value[startindex : startindex + subregobj.width]
00875
00876
              logging.debug(
00877
                  "{logdebug}getSubregister: retval = {retval}".format(
00878
                      logdebug=self.logdebug, retval=retval
00879
00880
00881
              return "", retval
00882
```

# 6.1.3.30 getTemp()

#### Definition at line 374 of file CameraAssembler.py.

```
00374 def getTemp(self, scale=None):
00375 return self.board.getTemp(scale)
00376
```

#### 6.1.3.31 getTimer()

```
{\tt nsCameraAssembler.CameraAssembler.getTimer} \ \ ( {\tt self} \ )
```

## Definition at line 356 of file CameraAssembler.py.

```
00356 def getTimer(self):
00357 return self.board.getTimer()
00358
```

## 6.1.3.32 getTiming()

#### Definition at line 407 of file CameraAssembler.py.

```
00407 def getTiming(self, side=None, actual=None):
00408 return self.sensor.getTiming(side, actual)
00409
```

#### 6.1.3.33 initBoard()

```
{\tt nsCameraAssembler.CameraAssembler.initBoard} \ \ ( {\tt self} \ )
```

Aliases to other objects' methods TODO: properly delegate these methods.

#### Definition at line 329 of file CameraAssembler.py.

```
00329 def initBoard(self):
00330 return self.board.initBoard()
00331
```

# 6.1.3.34 initialize()

```
{\tt nsCameraAssembler.CameraAssembler.initialize} \ \ ( {\tt self} \ )
```

#### End aliases.

Initialize board registers and set pots

## Definition at line 500 of file CameraAssembler.py.

```
def initialize(self):
00501
00502
              Initialize board registers and set pots
00503
00504
              # TODO: automate sensor and board selection from firmware info
00505
00507
00508
              # get_sensor
00509
              # TODO: pull sensor, board, comm id out to separate methods
00510
              if self.sensorname == "icarus":
0.0511
                  import nsCamera.sensors.icarus as snsr
              elif self.sensorname == "icarus2":
00512
00513
                  import nsCamera.sensors.icarus2 as snsr
              elif self.sensorname == "daedalus":
00514
00515
                 import nsCamera.sensors.daedalus as snsr
00516
              else: # catch-all for added sensors to attempt object encapsulation
                  sensormodname = ".sensors." + self.sensorname
00517
00518
                  try:
00519
                      sensormod = importlib.import_module(sensormodname, "nsCamera")
00520
                  except ImportError:
                      logging.critical(self.logcrit + "invalid sensor name")
00521
00522
                      sys.exit(1)
00523
                  snsr = getattr(sensormod, self.sensorname)
00524
              self.sensor = snsr(self)
00525
00526
              # kill existing connections (for reinitialize)
00527
              if hasattr(self, "comms"):
00528
                  self.closeDevice()
00529
00530
              # get communications interface
00531
              if self.commname == "rs422":
00532
                  import nsCamera.comms.RS422 as comms
00533
              elif self.commname == "gige":
00534
                  import nsCamera.comms.GigE as comms
00535
              else:
00536
                  commsmodname = ".comms." + self.commname
00537
00538
                      commsmod = importlib.import_module(commsmodname, "nsCamera")
                  except ImportError:
00539
00540
                      logging.critical(self.logcrit + "invalid comms name")
00541
                      sys.exit(1)
00542
                  comms = getattr(commsmod, self.commname)
00543
              self.comms = comms(self)
00544
00545
00546
              if self.boardname == "llnl_v1":
00547
                  import nsCamera.boards.LLNL_v1 as brd
00548
00549
                  self.board = brd.llnl_v1(self)
00550
              elif self.boardname == "llnl_v4":
                  import nsCamera.boards.LLNL_v4 as brd
00551
00552
00553
                  self.board = brd.llnl_v4(self)
00554
              else:
00555
                  boardmodname = ".board." + self.boardname
00556
00557
                      boardmod = importlib.import_module(boardmodname, "nsCamera")
00558
                  except ImportError:
                      logging.critical(self.logcrit + "invalid board name")
00559
00560
                      svs.exit(1)
00561
                  boardobj = getattr(boardmod, self.boardname)
                  self.board = boardobj(self)
00562
00563
00564
              # Now that board exists, initialize board-specific aliases for sensors
```

```
00565
               self.sensor.init_board_specific()
00566
00567
00568
00569
               # TODO: make cython the standard version
00570
               # ##############
00571
               # # For cython version
00572
00573
              # # get sensor
00574
             # if self.sensorname == "icarus":
00575
                     import nsCamera.sensors.icarus as snsr
00576
                     self.sensor = snsr.icarus(self)
             # elif self.sensorname == "icarus2":
# import nsCamera.sensors.icarus2 as snsr
# self.sensor = snsr.icarus2(self)
00577
00578
00579
00580
             # elif self.sensorname == "daedalus":
00581
                   import nsCamera.sensors.daedalus as snsr
00582
                    self.sensor = snsr.daedalus(self)
00583
              # # kill existing connections (for reinitialize)
00584
              # if hasattr(self, "comms"):
00585
00586
                     self.closeDevice()
00587
00588
              # # get communications interface
00589
              # if self.commname == "rs422":
00590
                    import nsCamera.comms.RS422 as comms
00591
                     self.comms = comms.RS422(self)
              # elif self.commname == "gige":
00592
                   import nsCamera.comms.GigE as comms
00593
00594
                     self.comms = comms.GigE(self)
00595
             # # get board
00596
             # if self.boardname == "llnl_v1":
# import nsCamera.boards.LLNL
00597
00598
                    import nsCamera.boards.LLNL v1 as brd
00599
                     self.board = brd.llnl_v1(self)
              # elif self.boardname == "llnl_v4":
00600
00601
                     import nsCamera.boards.LLNL_v4 as brd
00602
                     self.board = brd.llnl_v4(self)
             # ################
00603
00604
              err, rval = self.getRegister("FPGA_NUM")
if err or rval == "":
00605
00606
                   err, rval = self.getRegister("FPGA_NUM")
if err or rval == "":
00607
00608
                       logging.critical(
00609
                            self.logcrit + "Initialization failed: unable to communicate with"
00610
00611
                            " board. "
00612
                      )
00613
                   sys.exit(1)
00614
            self.initBoard()
00615
00616
               self.initPots()
00617
              self.initSensor()
00618
               self.initPowerCheck()
00619
              self.getBoardInfo()
00620
              self.printBoardInfo()
00621
```

# 6.1.3.35 initPots()

#### 6.1.3.36 initPowerCheck()

```
nsCamera.CameraAssembler.CameraAssembler.initPowerCheck ( self )
```

00334

Reset software and board timers for monitoring power status

```
Definition at line 1700 of file CameraAssembler.py.
```

```
01700 def initPowerCheck(self):
01701 """
01702 Reset software and board timers for monitoring power status
01703 """
01704 self.inittime = time.time()
01705 logging.info(self.loginfo + "resetting timer for power check function")
01706 self.resetTimer()
01707
```

#### 6.1.3.37 initSensor()

```
{\tt nsCameraAssembler.CameraAssembler.initSensor} \ \ ( {\tt self} \ )
```

#### Definition at line 338 of file CameraAssembler.py.

```
00338 def initSensor(self):
00339 return self.board.initSensor()
00340
```

## 6.1.3.38 latchPots()

```
{\tt nsCameraAssembler.CameraAssembler.latchPots} \ \ ( {\tt self} \ )
```

## Definition at line 335 of file CameraAssembler.py.

```
00335 def latchPots(self):
00336 return self.board.latchPots()
00337
```

# 6.1.3.39 mmReadoff()

```
Definition at line 1814 of file CameraAssembler.py.
```

```
def mmReadoff(self, waitOnSRAM, variation=None):
01815
01816
               Convenience function for parsing frames for use by MicroManager plugin
01817
               Args:
01818
                    waitOnSRAM: readoff wait flag
01819
                    variation: format of frames generated from readoff
01820
                        default - return first frame only
01821
                         "LastFrame" - return last frame only
                         "Average" - provide average of frames as single frame
"Landscape" - stitch frames together horizontally into single wide frame
01822
01823
01824
01825
               Returns:
               ndarray - single image frame
01826
01827
               frames, datalen, data_err = self.readoff(waitOnSRAM)
if variation == "LastFrame":
01828
01829
               return frames[self.sensor.nframes - 1]
elif variation == "Average":
01830
01831
01832
                    return np.sum(frames, axis=0) // self.sensor.nframes
01833
               elif variation == "Landscape":
01834
                    shaped = [
01835
                        np.reshape(frame, (self.sensor.maxheight, self.sensor.maxwidth))
01836
                        for frame in frames
01837
01838
                    return np.concatenate(shaped, axis=1)
01839
               else:
01840
                    return frames[0]
01841
```

#### 6.1.3.40 parseReadoff()

# Definition at line 446 of file CameraAssembler.py.

```
00446 def parseReadoff(self, frames, columns=1):
00447 return self.sensor.parseReadoff(frames, columns)
00448
```

## 6.1.3.41 partition()

# Definition at line 495 of file CameraAssembler.py.

```
00495 def partition(self, frames, columns):
00496 return partition(self, frames, columns)
00497
```

#### 6.1.3.42 plotFrames()

#### Definition at line 474 of file CameraAssembler.py.

```
00474 def plotFrames(self, frames, index=None):
00475 return plotFrames(self, frames, index)
00476
```

## 6.1.3.43 powerCheck()

```
\verb|nsCameraAssembler.CameraAssembler.powerCheck| (
               self,
               delta = 10)
Check to see if board power has persisted since powerCheck was last initialized.
  Compares time elapsed since initialization against board's timer. If the
  difference is greater than 'delta,' flag as False (power has likely failed)
    delta: difference in seconds permitted between software and board timers
Returns:
   boolean, 'True' means timer difference is less than 'delta' parameter;
              'False' indicates power failure
Definition at line 1708 of file CameraAssembler.py.
         def powerCheck(self, delta=10):
01709
01710
             Check to see if board power has persisted since powerCheck was last initialized.
01711
               Compares time elapsed since initialization against board's timer. If the
01712
               difference is greater than 'delta,' flag as False (power has likely failed)
01714
01715
                 delta: difference in seconds permitted between software and board timers
01716
01717
             Returns:
                 boolean, 'True' means timer difference is less than 'delta' parameter;
01718
01719
                          'False' indicates power failure
             ....
01720
01721
             elapsed = time.time() - self.inittime
01722
             logging.debug(self.logdebug + "powerCheck: elapsed time = " + str(elapsed))
01723
             difference = abs(elapsed - self.getTimer())
01724
             if difference > delta:
01725
                 logging.warning(
                     self.logwarn + "powerCheck function has failed; may indicate current "
01726
01727
                     "or recent power failure
01728
01729
             return difference < delta
01730
```

## 6.1.3.44 printBoardInfo()

```
nsCameraAssembler.CameraAssembler.printBoardInfo ( self \ )
```

## Definition at line 1731 of file CameraAssembler.py.

```
def printBoardInfo(self):
01732
              # TODO: add override option if logging level is above info
01733
              logging.info(
01734
                  self.loginfo
01735
                  + "Python version: "
01736
                  + str(self.python)
01737
01738
                  + str(self.pyth1)
01739
01740
                  + str(self.pvth2)
01741
01742
              logging.info(self.loginfo + "nsCamera software version: " + self.version)
              logging.info(self.loginfo + "FPGA firmware version: " + self.FPGAVersion)
01743
              logging.info(self.loginfo + "FPGA implementation: " + self.FPGANum)
01744
01745
              if self.FPGAinvalid:
                  logging.info(self.loginfo + "FPGA information unavailable")
01746
01747
              else:
01748
                  logging.info(self.loginfo + "Board type: " + self.FPGAboardtype)
```

```
logging.info(self.loginfo + "Rad-Tolerant: " + str(self.FPGArad))
logging.info(self.loginfo + "Sensor family: " + self.FPGAsensor)
logging.info(self.loginfo + "Sensor label: " + self.sensor.loglabel)
01749
01750
01751
01752
01753
                           self.loginfo + "Available interfaces: " + ", ".join(self.FPGAinterfaces)
01754
01755
                 if self.commname == "gige":
01756
                      ci = self.comms.CardInfoP.contents
01757
                      ip = ".".join(str(e) for e in [b for b in ci.IPAddr])
01758
                      logging.info(
01759
                           self.loginfo + "GigE connected to " + ip + ":" + str(self.port)
01760
01761
                 elif self.commname == "rs422":
01762
                      logging.info(self.loginfo + "RS422 connected to " + self.comms.port)
01763
```

#### 6.1.3.45 readFrames()

## Definition at line 455 of file CameraAssembler.py.

```
00455 def readFrames(self, waitOnSRAM=None, timeout=0, fast=False, columns=1):
00456 frames, _, _ = self.comms.readoff(waitOnSRAM, timeout, fast, columns)
00457 return frames
00458
```

#### 6.1.3.46 readImgs()

#### Definition at line 1453 of file CameraAssembler.py.

```
def readImgs(self, waitOnSRAM=True, mode="Hardware"):
01455
              Combines arm() and readoff() functions
01456
01457
01458
                 tuple (list of numpy arrays, length of downloaded payload, payload error
                   flag) returned by readoff
01459
01460
01461
             logging.info(self.loginfo + "readImgs")
01462
              self.arm(mode)
01463
              return self.readoff(waitOnSRAM)
01464
```

# 6.1.3.47 readoff()

```
nsCamera.CameraAssembler.CameraAssembler.readoff (
                self,
                waitOnSRAM = None,
                timeout = 0,
                fast = None,
                columns = 1)
Definition at line 459 of file CameraAssembler.py.
         def readoff(self, waitOnSRAM=None, timeout=0, fast=None, columns=1):
              return self.comms.readoff(waitOnSRAM, timeout, fast, columns)
00461
6.1.3.48 readSerial()
nsCamera.CameraAssembler.CameraAssembler.readSerial (
                self,
                size,
                timeout = None)
Definition at line 465 of file CameraAssembler.py.
          def readSerial(self, size, timeout=None):
    return self.comms.readSerial(size, timeout)
00465
00466
00467
6.1.3.49 readSRAM()
nsCamera.CameraAssembler.CameraAssembler.readSRAM (
                self )
Definition at line 350 of file CameraAssembler.py.
          def readSRAM(self):
00350
00351
              return self.board.readSRAM()
00352
6.1.3.50 reboot()
nsCamera.CameraAssembler.CameraAssembler.reboot (
                self )
Perform soft reboot on board and reinitialize
Definition at line 636 of file CameraAssembler.py.
          def reboot(self):
00636
00637
00638
              Perform soft reboot on board and reinitialize
00639
00640
              self.board.softReboot()
00641
              self.reinitialize()
00642
```

# 6.1.3.51 reinitialize()

```
nsCamera.CameraAssembler.CameraAssembler.reinitialize (
               self )
Reinitialize board registers and pots, reinitialize sensor timing (if
  previously set)
Definition at line 622 of file CameraAssembler.py.
00622
         def reinitialize(self): """
00623
             Reinitialize board registers and pots, reinitialize sensor timing (if
00624
00625
             previously set)
00626
00627
             logging.info(self.loginfo + "reinitializing")
00628
             self.initialize()
00629
00630
             for side in self.senstiming:
                 self.setTiming(side, self.senstiming[side][0], self.senstiming[side][1])
00631
00632
00633
             if self.sensmanual: # should be mutually exclusive with anything in senstiming
00634
                 self.setManualShutters(self.sensmanual)
00635
6.1.3.52 reportEdgeDetects()
```

```
nsCamera.CameraAssembler.CameraAssembler.reportEdgeDetects (
              self )
```

## Definition at line 392 of file CameraAssembler.py.

```
def reportEdgeDetects(self):
00393
             return self.board.reportEdgeDetects()
00394
```

#### 6.1.3.53 reportStatus()

```
nsCamera.CameraAssembler.CameraAssembler.reportStatus (
              self )
```

#### Definition at line 389 of file CameraAssembler.py.

```
def reportStatus(self):
            return self.board.reportStatus()
00391
```

## 6.1.3.54 resetTimer()

```
nsCamera.CameraAssembler.CameraAssembler.resetTimer (
              self )
```

#### Definition at line 359 of file CameraAssembler.pv.

```
00359
         def resetTimer(self):
00360
              return self.board.resetTimer()
00361
```

# 6.1.3.55 resolveSubreg()

```
nsCamera.CameraAssembler.CameraAssembler.resolveSubreg (
                self,
                srname )
Resolves subregister name or alias, returns object associated with subregister
  and flag indicating writability
    srname: name or alias of subregister
Returns:
    tuple (subregister name string, associated object, writable flag)
Definition at line 801 of file CameraAssembler.py.
00801
          def resolveSubreg(self, srname):
00802
00803
              Resolves subregister name or alias, returns object associated with subregister
00804
               and flag indicating writability
00805
00806
              Args:
00807
                 srname: name or alias of subregister
00808
00809
              Returns:
              tuple(subregister name string, associated object, writable flag)
00810
00811
00812
              logging.debug(
00813
                  "{logdebug}resolveSubreg: srname = {srname}".format(
00814
                     logdebug=self.logdebug,
00815
                      srname=srname,
00816
00817
00818
              writable = False
00819
              srname = srname.upper()
00820
              if srname in self.board.subreg_aliases:
00821
                 srname = self.board.subreg_aliases[srname].upper()
00822
              if srname in self.board.subreglist:
00823
                 srobj = getattr(self.board, srname)
00824
                  writable = getattr(self.board, srname).writable
00825
00826
                  # No-object error is handled by calling function
00827
                  srobj = None
00828
              logging.debug(
                  "{logdebug}resolveSubreg: srobj = {srobj}, writable={writable}".format(
00829
00830
                     logdebug=self.logdebug, srobj=srobj, writable=writable
00831
00832
00833
00834
              return srname, srobj, writable
00835
```

#### 6.1.3.56 saveFrames()

```
Save list of numpy arrays to disk. If passed an unprocessed text string, saves
  it directly to disk for postprocessing. Use 'prefix=""' for no prefix
Args:
    frames: numpy array or list of numpy arrays OR text string
    path: save path, defaults to './output'
    filename: defaults to 'frames.bin'
    prefix: prepended to filename, defaults to time/date (e.g. '160830-124704_')
Returns:
    Error string
Definition at line 1465 of file CameraAssembler.py.
         def saveFrames(self, frames, path=None, filename="frames", prefix=None):
01465
01466
01467
             Save list of numpy arrays to disk. If passed an unprocessed text string, saves
               it directly to disk for postprocessing. Use 'prefix=""' for no prefix
01468
01469
01470
                 frames: numpy array or list of numpy arrays OR text string
01471
01472
                 path: save path, defaults to './output'
01473
                 filename: defaults to 'frames.bin'
01474
                 prefix: prepended to filename, defaults to time/date (e.g. '160830-124704_')
01475
01476
             Returns:
             Error string
01477
01478
01479
             logging.debug(
01480
                self.logdebug
01481
                 + "saveFrames: path = "
01482
                 + str(path)
01483
                 + "; filename = "
01484
                 + str(filename)
01485
                 + "; prefix =
01486
                 + str(prefix)
01487
01488
             logging.info(self.loginfo + "saveFrames")
01489
             err = "
             if path is None:
01490
01491
                 path = os.path.join(os.getcwd(), "output")
01492
             if prefix is None:
01493
                 prefix = datetime.now().strftime("%y%m%d-%H%M%S%f")[:-5] + "_"
01494
             if not os.path.exists(path):
01495
                 os.makedirs(path)
01496
01497
             # TODO catch save file exceptions
01498
             if isinstance(frames[0], str):
01499
                 logging.debug(self.logdebug + "saveFrames: saving text frames")
01500
                 filename = filename + ".txt"
                 savefile = open(os.path.join(path, prefix + filename), "w+")
01501
01502
                 savefile.write(frames)
01503
             else:
01504
                 logging.debug(self.logdebug + "saveFrames: saving numerical frames")
01505
                 filename = filename + ".bin"
01506
                 stacked = np.stack(frames)
01507
                 try:
01508
                     stacked = stacked.reshape(
01509
                         (
01510
                             self.sensor.nframes,
                             self.sensor.height // (self.sensor.interlacing + 1),
01511
01512
                             self.sensor.width,
01513
                          )
01514
                     )
01515
                 except Exception as e:
                     err = self.logerr + "saveFrames: unable to save frames: " + str(e)
01516
                     logging.error(err)
01517
01518
                 stacked.tofile(os.path.join(path, prefix + filename))
01519
01520
             return err
01521
6.1.3.57 saveHDF()
```

nsCamera.CameraAssembler.CameraAssembler.saveHDF (

self.

```
frames,
path = None,
filename = "Acquisition",
prefix = None )
```

#### Definition at line 2092 of file CameraAssembler.py.

```
02098
02099
02100
              logging.info(self.loginfo + ": saveHDF")
02101
              err =
              if path is None:
02102
02103
                  path = os.path.join(os.getcwd(), "output")
              if prefix is None:
02104
                  prefix = datetime.now().strftime("%y%m%d-%H%M%S%f")[:-5] + "_"
02105
02106
              if not os.path.exists(path):
02107
                  os.makedirs(path)
02108
02109
              h5file = os.path.join(path, prefix + filename + ".hdf5")
              with h5py.File(h5file, "w") as f:
02110
                  # shotgrp = f.create_group("DATA/SHOT")
02111
02112
                  frame index = 0
02113
                  for frame in frames:
                      grp = f.create_group("DATA/SHOT/FRAME_0" + str(frame_index))
02114
02115
                      data = grp.create_dataset(
                          "DATA", (self.sensor.height, self.sensor.width), data=frame
02116
02117
02118
                      frame_index += 1
02119
02120
02121 """
02122 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
02123 LLNL-CODE-838080
02124
02125 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
02126 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
02127 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
02128 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
02129 be made under this license.
02130 """
```

#### 6.1.3.58 saveNumpys()

```
nsCamera.CameraAssembler.CameraAssembler.saveNumpys (
              self.
              frames,
              path = None,
              filename = "Frame",
              prefix = None,
              index = None)
Save numpy array or list of numpy arrays to disk as individual numpy data files,
  with frame number appended to filename.
Args:
    frames: numpy array or list of numpy arrays or single numpy array
   path: save path, defaults to './output'
    filename: defaults to 'Frame' followed by frame number
    prefix: prepended to 'filename', defaults to time/date
      (e.g. '160830-124704_')
    index: number to start frame numbering
Returns:
   Error string
```

# Definition at line 1522 of file CameraAssembler.py.

```
01529
01530
              Save numpy array or list of numpy arrays to disk as individual numpy data files,
01531
01532
               with frame number appended to filename.
01533
01534
              Args:
01535
                 frames: numpy array or list of numpy arrays or single numpy array
01536
                  path: save path, defaults to './output'
                  filename: defaults to 'Frame' followed by frame number
                  prefix: prepended to 'filename', defaults to time/date
01539
                    (e.g. '160830-124704_')
01540
                  index: number to start frame numbering
01541
01542
              Returns:
              Error string
01543
01544
01545
              logging.info(self.loginfo + "saveNumpys")
01546
              logging.debug(
01547
                 self.logdebug
                  + "saveNumpys: path = "
01548
01549
                 + str(path)
01550
                 + "; filename =
01551
                 + str(filename)
                 + "; prefix =
01552
01553
                  + str(prefix)
                  + "; index =
01554
01555
                  + str(index)
             )
01556
              err = ""
01557
             if path is None:
01558
                  path = os.path.join(os.getcwd(), "output")
01559
01560
              if prefix is None:
                  prefix = datetime.now().strftime("%y%m%d-%H%M%S%f")[:-5] + "_"
01561
              if not os.path.exists(path):
01562
01563
                 os.makedirs(path)
01564
              if index is None:
01565
                  firstnum = self.sensor.firstframe
01566
              else:
                  firstnum = index
01567
01568
              if not isinstance(frames, list):
01569
                  frames = [frames]
01570
01571
              # if this is a text string from fast readoff, do the numpy conversion now
01572
              if isinstance(frames[0], str):
01573
                  frames = generateFrames(frames)
01574
01575
              framestemp = np.copy(frames)
01576
01577
              for idx, frame in enumerate(framestemp):
01578
                  if idx < len(framestemp) / 2:</pre>
01579
                      interlacing = self.sensor.interlacing[0]
01580
01581
                      interlacing = self.sensor.interlacing[1]
01582
01583
                      if self.padToFull:
01584
                         frame = np.reshape(
01585
                              frame, (self.sensor.maxheight // (interlacing + 1), -1)
01586
01587
                      else:
                          frame = np.reshape(
01588
01589
                              frame,
01590
                              (
01591
                                  (self.sensor.lastrow - self.sensor.firstrow + 1)
01592
                                  // (interlacing + 1),
01593
                                  -1,
01594
                              ),
01595
                          )
                      namenum = filename + "_%d" % firstnum
01596
01597
                      nppath = os.path.join(path, prefix + namenum + ".npy")
01598
                      np.save(nppath, frame)
01599
                      firstnum += 1
                  except SystemExit:
01600
01601
01602
                  except KeyboardInterrupt:
01603
                      raise
01604
                  except Exception:
01605
                      err = self.logerr + "saveNumpys: unable to save arrays"
01606
                      logging.error(err)
01607
01608
              return err
```

01609

# 6.1.3.59 saveTiffs()

```
nsCamera.CameraAssembler.CameraAssembler.saveTiffs (
               self,
               frames,
               path = None,
               filename = "Frame",
               prefix = None,
               index = None)
Definition at line 471 of file CameraAssembler.py.
         def saveTiffs(self, frames, path=None, filename="Frame", prefix=None, index=None):
00471
00472
             return saveTiffs(self, frames, path, filename, prefix, index)
00473
6.1.3.60 selectOscillator()
nsCamera.CameraAssembler.CameraAssembler.selectOscillator (
               self.
               osc = None )
```

# Definition at line 425 of file CameraAssembler.py.

```
00425 def selectOscillator(self, osc=None):
00426 return self.sensor.selectOscillator(osc)
00427
```

# 6.1.3.61 sendCMD()

```
nsCamera.CameraAssembler.SendCMD ( self, \\ pkt \ )
```

#### Definition at line 449 of file CameraAssembler.py.

```
00449 def sendCMD(self, pkt):
00450 return self.comms.sendCMD(pkt)
```

# 6.1.3.62 sensorSpecific()

```
{\tt nsCameraAssembler.CameraAssembler.sensorSpecific \ (} \\ self \ )
```

#### Definition at line 422 of file CameraAssembler.py.

```
00422 def sensorSpecific(self):
00423 return self.sensor.sensorSpecific()
00424
```

# 6.1.3.63 setArbTiming()

```
nsCamera.CameraAssembler.CameraAssembler.setArbTiming (
               side = "AB",
               sequence = None )
Definition at line 404 of file CameraAssembler.py.
        def setArbTiming(self, side="AB", sequence=None):
             return self.sensor.setArbTiming(side, sequence)
00406
6.1.3.64 setExtClk()
\verb|nsCameraAssembler.CameraAssembler.setExtClk| (
               self,
               dilation = None,
               frequency = None )
Definition at line 443 of file CameraAssembler.py.
         def setExtClk(self, dilation=None, frequency=None):
             return self.sensor.setExtClk(dilation, frequency)
00444
00445
6.1.3.65 setFrames()
nsCamera.CameraAssembler.CameraAssembler.setFrames (
               self,
               minframe = None,
               maxframe = None)
Sets bounds on frames returned by board, inclusive (e.g., 0,3 returns four
frames). If called without parameters, resets to full set of frames.
Aras:
    minframe: first frame to read from board
    maxframe: last frame to read from board
Returns:
    Error string
Definition at line 1842 of file CameraAssembler.py.
         def setFrames(self, minframe=None, maxframe=None):
01843
             Sets bounds on frames returned by board, inclusive (e.g., 0,3 returns four
01844
01845
             frames). If called without parameters, resets to full set of frames.
01846
01847
             Args:
01848
                minframe: first frame to read from board
01849
                 maxframe: last frame to read from board
01850
01851
             Returns:
             Error string
01852
01853
```

logging.debug(

self.logdebug

+ str(minframe)

+ "setFrames: minframe = "

01854

01855

01856 01857

```
01858
                  + "; maxframe = "
01859
                   + str(maxframe)
01860
              if minframe is None:
01861
01862
                  minframe = self.sensor.minframe
01863
              if maxframe is None:
01864
                  maxframe = self.sensor.maxframe
01865
01866
                  not isinstance(minframe, int)
01867
                  or minframe < self.sensor.minframe</pre>
01868
                  or minframe > maxframe
01869
                  or not isinstance(maxframe, int)
01870
                  or maxframe > self.sensor.maxframe
01871
01872
                  err = (
01873
                       self.logerr + "setFrames: invalid frame limits submitted. Frame "
01874
                       "selection remains unchanged. "
01875
01876
                  logging.error(err)
01877
                  return err
01878
01879
              initframe = hex(minframe)[2:].zfill(8)
              finframe = hex(maxframe)[2:].zfill(8)
err1, _ = self.setRegister("FPA_FRAME_INITIAL", initframe)
01880
01881
01882
              err2, _ = self.setRegister("FPA_FRAME_FINAL", finframe)
01883
              self.sensor.firstframe = minframe
01884
              self.sensor.lastframe = maxframe
01885
              self.sensor.nframes = maxframe - minframe + 1
01886
              self.comms.payloadsize = (
01887
                  self.sensor.width
01888
                  * self.sensor.height
01889
                  * self.sensor.nframes
01890
                  * self.sensor.bytesperpixel
01891
              plural = ""
01892
              if self.sensor.nframes > 1:
01893
01894
                  plural = "s"
01895
              logging.info(
01896
                  self.loginfo
                  + "Readoff set to "
01897
01898
                  + str(self.sensor.nframes)
                  + " frame"
01899
                  + plural
+ " ("
01900
01901
01902
                  + str(minframe)
01903
                  + ", "
01904
                  + str(maxframe)
01905
                  + ")"
01906
01907
              err = err1 + err2
01908
              if err:
01909
                  logging.error(
01910
                      self.logerr + "setFrames may not have functioned properly: " + err
01911
01912
              return err
01913
```

# 6.1.3.66 setHighFullWell()

# 6.1.3.67 setInterlacing()

```
{\tt nsCameraAssembler.CameraAssembler.setInterlacing} \ (
```

```
ifactor = None,
side = None )
```

#### Definition at line 428 of file CameraAssembler.py.

```
00428 def setInterlacing(self, ifactor=None, side=None):
00429 return self.sensor.setInterlacing(ifactor, side)
00430
```

# 6.1.3.68 setLED()

```
nsCamera.CameraAssembler.CameraAssembler.setLED ( self, \\ LED = 1, \\ status = 1 \; )
```

# Definition at line 365 of file CameraAssembler.py.

```
00365 def setLED(self, LED=1, status=1):
00366 return self.board.setLED(LED, status)
00367
```

## 6.1.3.69 setManualShutters()

```
nsCamera.CameraAssembler.CameraAssembler.setManualShutters ( self, \\ timing = None \ )
```

#### Definition at line 410 of file CameraAssembler.py.

```
00410 def setManualShutters(self, timing=None):
00411 return self.sensor.setManualTiming(timing)
00412
```

# 6.1.3.70 setManualTiming()

## Definition at line 413 of file CameraAssembler.py.

```
00413 def setManualTiming(self, timing=None):
00414 return self.sensor.setManualTiming(timing)
00415
```

# 6.1.3.71 setPhiDelay()

#### Definition at line 440 of file CameraAssembler.py.

```
00440 def setPhiDelay(self, side=None, delay=0):
00441 return self.sensor.setPhiDelay(side, delay)
00442
```

# 6.1.3.72 setPot()

```
nsCamera.CameraAssembler.CameraAssembler.setPot (
               potname,
               value = 1.0,
                errflag = False )
Sets value of pot to value, normalized so that '1.0' corresponds with the
  fixed point maximum value of pot.
Args:
    potname: common name of pot, e.g., VRST found in board.subreg_aliases or
      defined in board.subregisters
    value: float between 0 and 1
    errflag: if True, return tuple with error string
Returns:
    if errflag:
        tuple: (error string, response packet as string)
    else:
         response packet as string
Definition at line 1065 of file CameraAssembler.pv.
          def setPot(self, potname, value=1.0, errflag=False):
01065
01066
              Sets value of pot to value, normalized so that ^{\prime}1.0^{\prime} corresponds with the
01067
01068
               fixed point maximum value of pot.
01069
01070
01071
                 potname: common name of pot, e.g., VRST found in board.subreg_aliases or
01072
                   defined in board.subregisters
01073
                 value: float between 0 and 1
01074
                 errflag: if True, return tuple with error string
01075
01076
              Returns:
                 if errflag:
01077
01078
                     tuple: (error string, response packet as string)
01079
                  else:
01080
                     response packet as string
01081
01082
              logging.debug(
01083
                  "{logdebug}setPot: potname = {potname}; value={value} errflag = {errflag}"
01084
01085
                      logdebug=self.logdebug, potname=potname, value=value, errflag=errflag
01086
01087
             )
01088
01089
              if value < 0:</pre>
01090
                 value = 0.0
01091
              if value > 1:
                  value = 1.0
01092
01093
01094
              potname, potobj, writable = self.resolveSubreg(potname)
01095
              if not potobj:
01096
                 err = "{logerr}setPot: invalid lookup: {potname}; returning '0'".format(
01097
                     logerr=self.logerr, potname=potname
01098
01099
01100
                 logging.error(err)
01101
                  if errflag:
01102
                     return err. 0
01103
                 return 0
01104
              if not writable:
                  err = "{logerr}setPot: not a writable subregister: {potname}; returning '0'"
01105
                  "".format(logerr=self.logerr, potname=potname)
01106
01107
                  logging.error(err)
01108
                 if errflag:
                 return err, "0"
return "0"
01109
01110
              setpoint = int(round(value * potobj.max_value))
01111
```

```
01112
              setpointpadded = "{num:{fill}{width}b}".format(
                  num=setpoint, fill="0", width=potobj.width
01113
01114
              logging.debug(
01115
                   "{logdebug}setpot: setpointpadded = {setpointpadded}".format(
01116
01117
                       logdebug=self.logdebug, setpointpadded=setpointpadded
01118
01119
01120
01121
              err, rval = self.setSubregister(potname, setpointpadded)
01123
                  logging.error(
01124
                       err="{logerr}setPot: unable to confirm setting of subregister:"
                         {potname} ".format(logerr=self.logerr, potname=potname)
01125
01126
              ident = potname[3:]
01128
              if ident[0].isdigit(): # numbered pot scheme
01129
                  potnumlatch = int(ident) * 2 + 1
                  potnumlatchstring = "{num:{fill}{width}x}".format(
01130
                      num=potnumlatch, fill="0", width=8
01131
01132
01133
                  err1, resp = self.setRegister("POT_CTL", potnumlatchstring)
01134
              else: # alphabetical DAC scheme
                 ident = ident.upper() # expects single character, e.g. 'A' from 'DACA'
01135
                  identnum = ord(ident) - ord("A") # DACA -> 0
potnumlatch = int(identnum) * 2 + 1
01136
01137
                  potnumlatchstring = "{num:{fill}{width}x}".format(
    num=potnumlatch, fill="0", width=8
01138
01139
01140
                  err1, resp = self.setRegister("DAC_CTL", potnumlatchstring)
01141
              if errl:
01142
                  # logging.error(self.logerr + "setPot: unable to latch register")
01143
01144
01145
                  logging.error(
01146
                       err="{logerr}setPot: unable to latch register".format(
01147
                          logerr=self.logerr
01148
01149
01150
01151
              if errflag:
01152
                  return err + err1, rval
01153
              return rval
01154
```

#### 6.1.3.73 setPotV()

```
nsCamera.CameraAssembler.CameraAssembler.setPotV (
                                                                                               self.
                                                                                               potname,
                                                                                               voltage,
                                                                                                tune = False,
                                                                                               accuracy = 0.01,
                                                                                               iterations = 20,
                                                                                                approach = 0.75,
                                                                                                errflag = False )
Sets pot to specified voltage. If tune=True, uses monitor to adjust pot to
              correct voltage. Tuning will attempt to tune to closest LSB on pot; if
               'accuracy' > LSB resolution, will only complain if tuning is unable to get
              the voltage within 'accuracy'
Args:
                           potname: common name of pot, e.g., VRST found in board.subreg_aliases or
                                        defined in board.subregisters
                            voltage: voltage bound by pot max and min (set in board object)
                           tune: if True, iterate with monitor to correct voltage % \left( 1\right) =\left( 1\right) \left( 1\right
                            accuracy: acceptable error in volts (if None, attempts to find the closest
                                         possible pot setting and warns if last iteration does not reduce error
```

```
below the resolution of the pot)
iterations: number of iteration attempts
approach: approximation parameter (>1 may cause overshoot)
errflag: if True, return tuple with error string

Returns:
   if errflag:
        tuple: (error string, response string)
   else:
        response string
```

#### Definition at line 1197 of file CameraAssembler.py.

```
01206
01207
              Sets pot to specified voltage. If tune=True, uses monitor to adjust pot to
01208
                correct voltage. Tuning will attempt to tune to closest LSB on pot; if
01209
                'accuracy' > LSB resolution, will only complain if tuning is unable to get the voltage within 'accuracy'
01210
01211
01212
01213
              Args:
                  potname: common name of pot, e.g., VRST found in board.subreg_aliases or
01214
01215
                    defined in board.subregisters
01216
                  voltage: voltage bound by pot max and min (set in board object)
01217
                  tune: if True, iterate with monitor to correct voltage
01218
                  accuracy: acceptable error in volts (if None, attempts to find the closest
01219
                    possible pot setting and warns if last iteration does not reduce error
01220
                    below the resolution of the pot)
01221
                  iterations: number of iteration attempts
01222
                  approach: approximation parameter (>1 may cause overshoot)
01223
                  errflag: if True, return tuple with error string
01224
01225
              Returns:
01226
                  if errflag:
01227
                      tuple: (error string, response string)
01228
                  else:
01229
                      response string
01230
01231
              logging.debug(
01232
                  self.logdebug
01233
                  + "setPotV: potname = "
01234
                  + str(potname)
01235
                  + "; voltage = "
01236
                  + str(voltage)
01237
                  + "; tune = '
01238
                  + str(tune)
01239
                  + "; accuracy = "
01240
                  + str(accuracy)
01241
                  + "; iterations =
01242
                  + str(iterations)
01243
                  + "; approach =
01244
                  + str(approach)
01245
                  + "; errflag = "
01246
                  + str(errflag)
01247
01248
              potname, potobj, writable = self.resolveSubreg(potname)
01249
              if not potobj:
01250
                  err = (
01251
                      self.logerr
01252
                       + "setPotV: invalid lookup: "
                       + potname
+ " , returning zero "
01253
01254
01255
01256
                  logging.error(err)
01257
                  if errflag:
01258
                      return err, 0
01259
                  return 0
01260
              if not writable:
01261
                  err = (
                      self.logerr
01262
                      + "setPotV: not a writable subregister: "
01263
01264
                      + potname
01265
                      + "; returning zero"
01266
01267
                  logging.error(err)
01268
                  if errflag:
01269
                      return err, 0
01270
                  return 0
```

```
01271
               if voltage < potobj.minV:</pre>
                    voltage = potobj.minV
01272
01273
               if voltage > potobj.maxV:
                   voltage = potobj.maxV
01274
               setting = (voltage - potobj.minV) / (potobj.maxV - potobj.minV)
logging.debug(self.logdebug + "setPotV: setting = " + str(setting))
01275
01276
01277
               err, rval = self.setPot(potname, setting, errflag=True)
01278
               time.sleep(0.1)
01279
               # TODO: refactor tuning to separate method
01280
               if tune:
                    logging.debug(self.logdebug + "setPotV: beginning tuning")
01281
                    if potname not in self.board.monitor_controls.values():
01282
01283
                        err = (
                            self.logerr
01285
                             + "setPotV: pot '"
01286
                             + potname
01287
                             + "' does not have a corresponding monitor"
01288
01289
                        logging.error(err)
01290
                        if errflag:
01291
                            return err, rval
01292
                        return rval
01293
                    self.setPot(potname, 0.65)
01294
                   time.sleep(0.2)
01295
                    err1, mon65 = self.getMonV(potname, errflag=True)
01296
                    self.setPot(potname, 0.35)
01297
                    time.sleep(0.2)
01298
                    err2, mon35 = self.getMonV(potname, errflag=True)
                   # theoretical voltage range assuming linearity
potrange = (mon65 - mon35) / 0.3
01299
01300
01301
                    stepsize = potrange / (potobj.max_value + 1)
                    err += err1 + err2
01302
01303
                    if err or potrange < 1:</pre>
                        err += self.logerr + "setPotV: unable to tune pot " + potname
01304
                        if potrange < 1: # potrange should be on the order of 3.3 or 5 volts
    err += "; monitor shows insufficient change with pot variation; "</pre>
01305
01306
                             "retrying setPotV with tune=False"
01307
01308
                        logging.warning(err)
                        err, rval = self.setPotV(
01309
01310
                            potname=potname, voltage=voltage, tune=False, errflag=True
01311
01312
                        if errflag:
01313
                            return err, rval
01314
                        return rval
                   potzero = 0.35 - (mon35 / potrange)
potone = 1.65 - (mon65 / potrange)
01315
01316
01317
                   if potzero < 0:</pre>
                        potzero = 0
01318
01319
                    if potone > 1:
01320
                        potone = 1
01321
01322
                    if accuracy > stepsize:
01323
                        mindiff = accuracy
01324
01325
                        mindiff = stepsize
01326
                    setting = potzero + (voltage / potone)
01327
                    self.setPot(potname, setting)
01328
                    lastdiff = 0
01329
                    smalladjust = 0
01330
                    err3 = ""
                    for _ in range(iterations):
01331
01332
                        err3i, measured = self.getMonV(potname, errflag=True)
                        if err3i:
01333
01334
                            err3 = err3 + err3i + " "
                        diff = voltage - measured
01335
                        if abs(diff - lastdiff) < stepsize / 2:</pre>
01336
01337
                            if (
01338
                                 smalladjust > 12
01339
                                # magic number for now; if it doesn't converge after several
                                     tries, it never will, usually because the setting is pinned
01340
01341
                                     to 0 or 1 and adjust can't change it
01342
                                 logging.warning(
01343
                                     self.logwarn
                                      + "setPotV: Tuning converged too slowly: pot "
01344
01345
                                      + potname
                                      + " set to "
01346
01347
                                      + str(voltage)
                                      + "V, monitor returns "
01348
01349
                                      + str(measured)
                                      + "\mbox{\em V}; if this value is incorrect, consider trying "
01350
                                      + "tune=False"
01351
```

```
01352
01353
                                logging.debug(self.logdebug + "setPotV: tuning complete")
                                if errflag:
return "", rval
01354
01355
01356
                                return rval
01357
                           smalladjust += 1
01358
                       if not int (2 * diff / stepsize):
01359
                           # TODO: is this check redundant with the first one?
01360
                           logging.debug(self.logdebug + "setPotV: tuning complete")
                           if errflag:
    return "", rval
01361
01362
01363
                           return rval
01364
                       adjust = approach * (diff / potrange)
01365
                       setting += adjust
01366
                       if setting > 1:
01367
                           setting = 1
                       elif setting < 0:</pre>
01368
01369
                          setting = 0
                       err1, rval = self.setPot(potname, setting, True)
01370
01371
                       lastdiff = diff
01372
                       time.sleep(0.2)
01373
                   err4, measured = self.getMonV(potname, errflag=True)
                   diff = voltage - measured
01374
01375
                   # code will try to get to within one stepsize, but will only complain if it
                   # doesn't get within mindiff
if int(diff / mindiff):
01376
01377
01378
                       logging.warning(
01379
                           self.logwarn
                           + "setPotV: pot "
01380
01381
                           + potname
+ " set to "
01382
01383
                           + str(voltage)
                           + "V, monitor returns "
01384
01385
                           + str(measured)
01386
01387
01388
                   err += err1 + err2 + err3 + err4
01389
               if err:
                  logging.error(self.logerr + "setPotV: errors occurred: " + err)
01390
01391
               if errflag:
01392
                   return err, rval
               logging.debug(self.logdebug + "setPotV: tuning complete")
01393
01394
               return rval
01395
```

## 6.1.3.74 setPowerSave()

```
nsCamera.CameraAssembler.CameraAssembler.setPowerSave ( self, \\ status = 1 \ )
```

## Definition at line 368 of file CameraAssembler.py.

```
00368 def setPowerSave(self, status=1):
00369 return self.board.setPowerSave(status)
00370
```

## 6.1.3.75 setPPER()

## Definition at line 371 of file CameraAssembler.py.

```
00371 def setPPER(self, pollperiod=None):
00372 return self.board.setPPER(pollperiod)
00373
```

## 6.1.3.76 setRegister()

```
nsCamera.CameraAssembler.CameraAssembler.setRegister (
               self.
                regname,
               regval )
Sets named register to given value
Args:
    regname: name of register as given in ICD
    regval: value to assign to register, as integer or hexadecimal string
      with or without '0x'
Returns:
    tuple: (error string, response string)
Definition at line 743 of file CameraAssembler.py.
00743
          def setRegister(self, regname, regval):
00744
00745
              Sets named register to given value
00746
00747
              Args:
                 regname: name of register as given in ICD
00748
00749
                 regval: value to assign to register, as integer or hexadecimal string
00750
                   with or without '0x'
00751
00752
              Returns:
              tuple: (error string, response string)
00753
00754
00755
              logging.debug(
00756
                  "{logdebug}setRegister: regname = {regname}; regval = {regval}".format(
00757
                      logdebug=self.logdebug, regname=regname, regval=regval
00758
                 )
00759
             )
00760
00761
             regname = regname.upper()
00762
              if regname not in self.board.registers:
                  err = "{logerr}setRegister: Invalid register name: {regname}".format(
00763
                     logerr=self.logerr, regname=regname
00764
00765
00766
                 logging.error(err)
00767
                  return err, "00000000"
00768
              if isinstance(regval, int):
00769
                 regval = hex(regval)
00770
              try:
00771
                 if regval[0:2] == "0x":
00772
                      regval = regval[2:]
00773
              except TypeError:
00774
                  err = "{logerr}setRegister: invalid register value parameter".format(
00775
                     logerr=self.logerr
00776
00777
                 logging.error(err)
00778
                  return err, "00000000"
00779
              pkt = Packet(addr=self.board.registers[regname], data=regval)
00780
              err, rval = self.comms.sendCMD(pkt)
00781
              if err:
00782
                  logging.error(
00783
                      "{logerr}setRegister: {regname}: {err}".format(
00784
                          logerr=self.logerr, regname=regname, err=err
00785
00786
00787
              if len(rval) < 32:</pre>
00788
                  logging.debug(
00789
                      "{logdebug}SetRegister: rval = {rval}".format(
00790
                          logdebug=self.logdebug, rval=rval
00791
00792
                 )
00793
              else:
00794
                  logging.debug(
00795
                      "{logdebug}SetRegister: rval (truncated) = {rval}".format(
00796
                          logdebug=self.logdebug, rval=rval[0:32]
00797
00798
                 )
00799
              return err, rval
00800
```

# 6.1.3.77 setRows()

```
\verb|nsCameraAssembler.CameraAssembler.setRows| (
                minrow = 0.
                maxrow = None,
                padToFull = False )
Sets bounds on rows returned by board, inclusive (e.g., 0,1023 returns all 1024
  rows). If called without parameters, resets to full image size.
Args:
    minrow: first row to return from board
    maxrow: last row to return from board
    padToFull: if True, generate full size frames, padding collected rows with
      zeroes if necessary
Definition at line 1914 of file CameraAssembler.py.
01914
          def setRows(self, minrow=0, maxrow=None, padToFull=False):
01915
              Sets bounds on rows returned by board, inclusive (e.g., 0,1023 returns all 1024
01916
01917
                rows). If called without parameters, resets to full image size.
01918
01919
              Args:
                  minrow: first row to return from board
01920
01921
                  maxrow: last row to return from board
01922
                  padToFull: if True, generate full size frames, padding collected rows with
01923
                    zeroes if necessary
01924
              logging.debug(
01925
01926
                  self.logdebug
01927
                  + "setRows: minrow = "
01928
                  + str(minrow)
                 + "; maxrow = "
01929
01930
                  + str(maxrow)
                  + "; padToFull = "
01931
01932
                  + str(padToFull)
01933
01934
             if maxrow is None:
01935
                  maxrow = self.sensor.maxheight - 1
01936
01937
                  not isinstance(minrow, int)
01938
                  or minrow < 0
                  or minrow > maxrow
01939
01940
                  or not isinstance(maxrow, int)
01941
                  or maxrow >= self.sensor.maxheight
01942
01943
01944
                      self.logerr + "setRows: invalid row arguments submitted. Frame size"
01945
                      " remains unchanged.
01946
01947
                  logging.error(err)
01948
                  return err
01949
01950
              initrow = hex(minrow)[2:].zfill(8)
01951
              finrow = hex(maxrow)[2:].zfill(8)
              err1, _ = self.setRegister("FPA_ROW_INITIAL", initrow)
err2, _ = self.setRegister("FPA_ROW_FINAL", finrow)
01952
01953
01954
              self.sensor.firstrow = minrow
01955
              self.sensor.lastrow = maxrow
01956
              self.sensor.height = maxrow - minrow + 1
01957
              self.comms.payloadsize = (
01958
                  self.sensor.width
01959
                  * self.sensor.height
01960
                  * self.sensor.nframes
01961
                  * self.sensor.bytesperpixel
01962
              )
01963
              if self.commname == "rs422":
01964
01965
                  self.comms.datatimeout = (
                      (1.0 * self.sensor.height / self.sensor.maxheight)
01966
01967
                      * 5e7
```

nsCamera.CameraAssembler.CameraAssembler.setSubregister (

```
01968
                      * self.sensor.nframes
01969
                      / self.comms.baud
01970
01971
              self.padToFull = padToFull
01972
              logging.info(
01973
                  self.loginfo
01974
                  + "Readoff set to "
01975
                  + str(self.sensor.height)
01976
                  + " rows ("
01977
                  + str(minrow)
01978
01979
                  + str(maxrow)
01980
                  + ")"
01981
01982
              err = err1 + err2
01983
             if err:
01984
                 logging.error(
01985
                      self.logerr + "setRows may not have functioned properly: " + err
01986
01987
             return err
01988
```

#### 6.1.3.78 setSubregister()

self,

```
valstring )
Sets substring of register identified in board attribute 'subregname' to
  valstring if subregister is writable
Aras:
    subregname: listed in board.subreg_aliases or defined in board.subregisters
    valstring: integer or binary string with or without '0b'
Returns:
    tuple: (error, packet response string) from setRegister
Definition at line 883 of file CameraAssembler.py.
00883
         def setSubregister(self, subregname, valstring):
00884
00885
             Sets substring of register identified in board attribute 'subregname' to
00886
               valstring if subregister is writable
00887
00888
00889
                 subregname: listed in board.subreg_aliases or defined in board.subregisters
00890
                 valstring: integer or binary string with or without '0b'
00891
00892
             tuple: (error, packet response string) from setRegister
00893
00894
00895
             logging.debug(
00896
                  "{logdebug}setSubegister: subregname = {subregname}; valstring ="
00897
                 " {valstring}".format(
00898
                     logdebug=self.logdebug, subregname=subregname, valstring=valstring
00899
00900
             )
00901
             subregname, subregobj, writable = self.resolveSubreg(subregname)
00902
00903
             if not subregobj:
00904
                 err = "{logerr}getSubregister: invalid lookup: {subregname}".format(
00905
                     logerr=self.logerr, subregname=subregname
00906
00907
00908
                 logging.error(err)
00909
                 return err, "0'
00910
             if not writable:
00911
                 err = "{logerr}getSubregister: not a writable subregister: {subregname}"
```

```
00912
                    "".format(logerr=self.logerr, subregname=subregname)
00913
                    logging.error(err)
00914
                    return err, "0"
               if isinstance(valstring, int):
00915
00916
                   valstring = bin(valstring)[2:]
00917
00918
                   if valstring[0:2] == "0b":
00919
                        valstring = valstring[2:]
00920
               except TypeError:
00921
                   err = "{logerr}getSubregister: invalid subregister value parameter".format(
00922
                       logerr=self.logerr
00923
00924
00925
                   logging.error(err)
00926
                    return err, "0"
               if len(str(valstring)) > subregobj.width:
00928
                   err = "{logerr}getSubregister: ialue string is too long".format(
00929
                        logerr=self.logerr
00930
00931
00932
                   logging.error(err)
                   return err, "0'
00933
               # read current value of register data
00934
00935
               err, resp = self.getRegister(subregobj.register)
00936
               if err:
00937
                   logging.error(
00938
                         "{logerr}getSubregister: unable to retrieve register setting; setting"
                        " of {subregname} likely failed ".format(
00939
00940
                             {\tt logerr=self.logerr,\ subregname=subregname}
00941
00942
                   )
00943
               return err, "0"
hex_str = "0x" + resp
00944
00945
               b_reg_value = "[0:0=32b]".format(int(hex_str, 16))  # convert to binary  # list indexing is reversed from bit string; the last bit of the string is at
00946
00947
               \# index 0 in the list (thus bit 0 is at index 0)
00948
00949
               startindex = 31 - subregobj.start_bit
               valstringpadded = str(valstring).zfill(subregobj.width)
00950
00951
               fullreg = list(b_reg_value)
00952
               fullreg[startindex : startindex + subregobj.width] = valstringpadded
               # convert binary string back to hexadecimal string for writing
new_reg_value = "".join(fullreg)
h_reg_value = "{num:{fill}{width}x}".format(
00953
00954
00955
00956
                    num=int(new_reg_value, 2), fill="0", width=8
00957
00958
               err, retval = self.setRegister(subregobj.register, h_reg_value)
00959
               # logging.debug(self.logdebug + "retval = " + str(retval))
00960
                if len(retval) < 32:</pre>
00961
                    logging.debug(
00962
                        "{logdebug}setSubregister: retval = {retval}".format(
00963
                             logdebug=self.logdebug, retval=retval
00964
00965
                   )
00966
               else:
00967
                    logging.debug(
00968
                        "{logdebug}setSubregister: retval (truncated) = {retval}".format(
00969
                             logdebug=self.logdebug, retval=retval[0:32]
00970
00971
                    )
00972
00973
               return err, retval
00974
```

#### 6.1.3.79 setTiming()

```
6.1.3.80 setTriggerDelay()
nsCamera.CameraAssembler.CameraAssembler.setTriggerDelay (
               self,
               delay = 0 )
Definition at line 437 of file CameraAssembler.py.
        def setTriggerDelay(self, delay=0):
00438
             return self.sensor.setTriggerDelay(delay)
00439
6.1.3.81 setZeroDeadTime()
nsCamera.CameraAssembler.CameraAssembler.setZeroDeadTime (
               self.
               flag = True,
               side = None )
Definition at line 434 of file CameraAssembler.py.
       def setZeroDeadTime(self, flag=True, side=None):
00434
             return self.sensor.setZeroDeadTime(flag, side)
00435
00436
6.1.3.82 startCapture()
\verb|nsCameraAssembler.CameraAssembler.startCapture| (
               self,
               mode )
Definition at line 347 of file CameraAssembler.py.
00347
        def startCapture(self, mode):
00348
            return self.board.startCapture(mode)
00349
6.1.3.83 str2bytes()
```

```
nsCamera.CameraAssembler.CameraAssembler.str2bytes (
              self,
             astring )
```

# Definition at line 483 of file CameraAssembler.py.

```
def str2bytes(self, astring):
            return str2bytes(astring)
00485
```

### 6.1.3.84 str2nparray()

```
nsCamera.CameraAssembler.CameraAssembler.str2nparray (
              self,
              valstring )
```

### Definition at line 489 of file CameraAssembler.py.

```
def str2nparray(self, valstring):
00489
00490
              return str2nparray(valstring)
00491
```

### 6.1.3.85 submitMessages()

```
nsCamera.CameraAssembler.CameraAssembler.submitMessages (
               self,
               messages,
               errorstring = "Error" )
Serially set multiple register / subregister values
Args:
    messages: list of tuples (register name, integer or hexadecimal string with
      or without '0x') and/or (subregister name, integer or binary string with
      or without '0b')
    errorstring: error message to print in case of failure
Returns:
    tuple (accumulated error string, response string of final message)
Definition at line 975 of file CameraAssembler.py.
          def submitMessages(self, messages, errorstring="Error"):
00976
00977
             Serially set multiple register / subregister values
00978
00979
             Args:
00980
                 messages: list of tuples (register name, integer or hexadecimal string with
                   or without '0x') and/or (subregister name, integer or binary string with
00981
                   or without '0b')
00982
00983
                 errorstring: error message to print in case of failure
00984
00985
             Returns:
             tuple (accumulated error string, response string of final message)
00986
00987
00988
             logging.debug(
00989
                  "{logdebug}submitMessages: messages = {messages}; errorstring ="
00990
                  " {errorstring}".format(
00991
                     logdebug=self.logdebug, messages=messages, errorstring=errorstring
00992
                 )
00993
             )
00994
             errs = ""
00995
             err = ""
00996
             rval = ""
00997
00998
             for m in messages:
00999
                 if m[0].upper() in self.board.registers:
01000
                     err, rval = self.setRegister(m[0].upper(), m[1])
01001
                 elif m[0].upper() in self.board.subreglist:
01002
                     err, rval = self.setSubregister(m[0].upper(), m[1])
01003
                 else:
01004
                     err = "{logerr}submitMessages: Invalid register/subregister:"
01005
                     " {errorstring}:{m0}; ".format(
01006
                         logerr=self.logerr, errorstring=errorstring, m0=m[0]
01007
01008
01009
                     logging.error(err)
01010
                 errs = errs + err
01011
             return err, rval
01012
6.1.3.86 waitForSRAM()
nsCamera.CameraAssembler.CameraAssembler.waitForSRAM (
               self,
               timeout = None )
Definition at line 353 of file CameraAssembler.py.
         def waitForSRAM(self, timeout=None):
    return self.board.waitForSRAM(timeout)
00353
00354
00355
```

# 6.1.3.87 writeSerial()

# Definition at line 462 of file CameraAssembler.py.

```
00462 def writeSerial(self, cmd, timeout=None):
00463 return self.comms.writeSerial(cmd, timeout)
00464
```

#### 6.1.4 Member Data Documentation

#### 6.1.4.1 abort

nsCamera.CameraAssembler.CameraAssembler.abort

Definition at line 252 of file CameraAssembler.py.

#### 6.1.4.2 arch

nsCamera.CameraAssembler.CameraAssembler.arch

Definition at line 228 of file CameraAssembler.py.

### 6.1.4.3 armed

 $\verb|nsCameraAssembler.CameraAssembler.armed| \\$ 

Definition at line 244 of file CameraAssembler.py.

#### 6.1.4.4 board

 $\verb|nsCamera.CameraAssembler.CameraAssembler.board|$ 

Definition at line 549 of file CameraAssembler.py.

# 6.1.4.5 boardname

nsCamera.CameraAssembler.CameraAssembler.boardname

Definition at line 204 of file CameraAssembler.py.

# 6.1.4.6 commname

nsCamera.CameraAssembler.CameraAssembler.commname

Definition at line 211 of file CameraAssembler.py.

#### 6.1.4.7 comms

nsCamera.CameraAssembler.CameraAssembler.comms

Definition at line 543 of file CameraAssembler.py.

### 6.1.4.8 currtime

nsCamera.CameraAssembler.CameraAssembler.currtime

Definition at line 195 of file CameraAssembler.py.

# 6.1.4.9 cycle

nsCamera.CameraAssembler.CameraAssembler.cycle

Definition at line 203 of file CameraAssembler.py.

# 6.1.4.10 FPGAboardtype

 $\verb|nsCamera.CameraAssembler.CameraAssembler.FPGAboardtype|$ 

Definition at line 234 of file CameraAssembler.py.

# 6.1.4.11 FPGAinterfaces

nsCamera.CameraAssembler.CameraAssembler.FPGAinterfaces

Definition at line 237 of file CameraAssembler.py.

### 6.1.4.12 FPGAinvalid

nsCamera.CameraAssembler.CameraAssembler.FPGAinvalid

Definition at line 240 of file CameraAssembler.py.

# 6.1.4.13 FPGANum

nsCamera.CameraAssembler.CameraAssembler.FPGANum

Definition at line 231 of file CameraAssembler.py.

#### 6.1.4.14 FPGArad

nsCamera.CameraAssembler.CameraAssembler.FPGArad

Definition at line 235 of file CameraAssembler.py.

### 6.1.4.15 FPGAsensor

nsCamera.CameraAssembler.CameraAssembler.FPGAsensor

Definition at line 236 of file CameraAssembler.py.

#### 6.1.4.16 FPGAVersion

nsCamera.CameraAssembler.CameraAssembler.FPGAVersion

Definition at line 230 of file CameraAssembler.py.

### 6.1.4.17 inittime

 $\verb|nsCameraAssembler.CameraAssembler.inittime| \\$ 

Definition at line 250 of file CameraAssembler.py.

# 6.1.4.18 iplist

nsCamera.CameraAssembler.CameraAssembler.iplist

Definition at line 242 of file CameraAssembler.py.

# 6.1.4.19 logcrit

nsCamera.CameraAssembler.CameraAssembler.logcrit

Definition at line 272 of file CameraAssembler.py.

# 6.1.4.20 logcritbase

 $\verb|nsCameraAssembler.CameraAssembler.logcritbase|\\$ 

Definition at line 266 of file CameraAssembler.py.

# 6.1.4.21 logdebug

nsCamera.CameraAssembler.CameraAssembler.logdebug

Definition at line 276 of file CameraAssembler.py.

### 6.1.4.22 logdebugbase

nsCamera.CameraAssembler.CameraAssembler.logdebugbase

Definition at line 270 of file CameraAssembler.py.

# 6.1.4.23 logerr

nsCamera.CameraAssembler.CameraAssembler.logerr

Definition at line 273 of file CameraAssembler.py.

#### 6.1.4.24 logerrbase

 $\verb|nsCamera.CameraAssembler.CameraAssembler.logerrbase|$ 

Definition at line 267 of file CameraAssembler.py.

# 6.1.4.25 loginfo

 $\verb|nsCameraAssembler.CameraAssembler.loginfo|\\$ 

Definition at line 275 of file CameraAssembler.py.

# 6.1.4.26 loginfobase

nsCamera.CameraAssembler.CameraAssembler.loginfobase

Definition at line 269 of file CameraAssembler.py.

# 6.1.4.27 logtag

 $\verb|nsCamera|. Camera| Assembler. Camera| Assembler. logtag$ 

Definition at line 264 of file CameraAssembler.py.

# 6.1.4.28 logwarn

nsCamera.CameraAssembler.CameraAssembler.logwarn

Definition at line 274 of file CameraAssembler.py.

### 6.1.4.29 logwarnbase

nsCamera.CameraAssembler.CameraAssembler.logwarnbase

Definition at line 268 of file CameraAssembler.py.

#### 6.1.4.30 oldtime

nsCamera.CameraAssembler.CameraAssembler.oldtime

Definition at line 196 of file CameraAssembler.py.

#### 6.1.4.31 packageroot

 $\verb|nsCamera.CameraAssembler.CameraAssembler.packageroot|\\$ 

Definition at line 243 of file CameraAssembler.py.

# 6.1.4.32 padToFull

 $\verb|nsCameraAssembler.CameraAssembler.padToFull|$ 

Definition at line 251 of file CameraAssembler.py.

### 6.1.4.33 parsedtime

nsCamera.CameraAssembler.CameraAssembler.parsedtime

Definition at line 201 of file CameraAssembler.py.

# 6.1.4.34 payloaderror

 $\verb|nsCamera Assembler.Camera Assembler.payloaderror|\\$ 

Definition at line 320 of file CameraAssembler.py.

# 6.1.4.35 platform

nsCamera.CameraAssembler.CameraAssembler.platform

Definition at line 227 of file CameraAssembler.py.

# 6.1.4.36 port

nsCamera.CameraAssembler.CameraAssembler.port

Definition at line 224 of file CameraAssembler.py.

### 6.1.4.37 PY3

nsCamera.CameraAssembler.CameraAssembler.PY3

Definition at line 226 of file CameraAssembler.py.

#### 6.1.4.38 pyth1

 $\verb|nsCameraAssembler.CameraAssembler.pyth1|\\$ 

Definition at line 225 of file CameraAssembler.py.

# 6.1.4.39 pyth2

nsCamera.CameraAssembler.CameraAssembler.pyth2

Definition at line 225 of file CameraAssembler.py.

# 6.1.4.40 python

 $\verb|nsCamera.CameraAssembler.CameraAssembler.python|\\$ 

Definition at line 225 of file CameraAssembler.py.

# 6.1.4.41 read

nsCamera.CameraAssembler.CameraAssembler.read

Definition at line 199 of file CameraAssembler.py.

### 6.1.4.42 savetime

nsCamera.CameraAssembler.CameraAssembler.savetime

Definition at line 202 of file CameraAssembler.py.

### 6.1.4.43 sensmanual

nsCamera.CameraAssembler.CameraAssembler.sensmanual

Definition at line 249 of file CameraAssembler.py.

#### 6.1.4.44 sensor

 $\verb|nsCamera|. Camera| Assembler. Camera| Assembler. sensor$ 

Definition at line 524 of file CameraAssembler.py.

### 6.1.4.45 sensorname

nsCamera.CameraAssembler.CameraAssembler.sensorname

For regular version.

Definition at line 216 of file CameraAssembler.py.

### 6.1.4.46 senstiming

nsCamera.CameraAssembler.CameraAssembler.senstiming

Definition at line 248 of file CameraAssembler.py.

#### 6.1.4.47 timeout

nsCamera.CameraAssembler.CameraAssembler.timeout

Definition at line 205 of file CameraAssembler.py.

# 6.1.4.48 trigtime

nsCamera.CameraAssembler.CameraAssembler.trigtime

Definition at line 197 of file CameraAssembler.py.

### 6.1.4.49 unstringed

nsCamera.CameraAssembler.CameraAssembler.unstringed

Definition at line 200 of file CameraAssembler.py.

#### 6.1.4.50 verblevel

 $\verb|nsCamera.CameraAssembler.CameraAssembler.verblevel|\\$ 

Definition at line 278 of file CameraAssembler.py.

# 6.1.4.51 verbmap

nsCamera.CameraAssembler.CameraAssembler.verbmap

Definition at line 254 of file CameraAssembler.py.

# 6.1.4.52 verbose

nsCamera.CameraAssembler.CameraAssembler.verbose

Definition at line 223 of file CameraAssembler.py.

### 6.1.4.53 version

nsCamera.CameraAssembler.CameraAssembler.version

Definition at line 194 of file CameraAssembler.py.

### 6.1.4.54 waited

nsCamera.CameraAssembler.CameraAssembler.waited


Definition at line 198 of file CameraAssembler.py.

The documentation for this class was generated from the following file:

C:/Users/hill35/git/camera\_python/nsCamera/CameraAssembler.py

# 6.2 nsCamera.sensors.daedalus.daedalus Class Reference

Inheritance diagram for nsCamera.sensors.daedalus.daedalus:



#### **Public Member Functions**

- \_\_init\_\_ (self, ca)
- sensorSpecific (self)
- setInterlacing (self, ifactor=None, side=None)
- setHighFullWell (self, flag)
- setZeroDeadTime (self, flag=True, side=None)
- selectOscillator (self, osc=None)
- setTriggerDelay (self, delay=0)
- setPhiDelay (self, side=None, delay=0)
- setExtClk (self, dilation=None, frequency=None)
- setManualShutters (self, timing=None)
- setManualTiming (self, timing=None)
- getManualTiming (self)
- getSensTemp (self, scale=None, offset=None, slope=None, dec=1)
- parseReadoff (self, frames, columns)
- reportStatusSensor (self, statusbits, statusbits2)

### Public Member Functions inherited from nsCamera.sensors.sensorBase.sensorBase

- init\_board\_specific (self)
- · checkSensorVoltStat (self)
- setTiming (self, side="AB", sequence=None, delay=0)
- setArbTiming (self, side="AB", sequence=None)
- getTiming (self, side, actual)
- getSensorStatus (self)

#### **Public Attributes**

- ca
- sens\_registers
- sens\_subregisters
- interlacing
- columns
- HFW
- ZDT

# Public Attributes inherited from nsCamera.sensors.sensorBase.sensorBase

- ca
- logcrit
- logerr
- logwarn
- loginfo
- logdebug
- fpganumID
- sensfam

#### **Static Public Attributes**

- str specwarn = ""
- int minframe = 0
- int maxframe = 2
- int maxwidth = 512
- int maxheight = 1024
- int bytesperpixel = 2
- int fpganumID = 2
- str detect = "DAEDALUS\_DET"
- str sensfam = "Daedalus"
- str loglabel = "[Daedalus] "
- bool ZDT = False
- bool HFW = False
- int firstframe = 0
- int lastframe = 2
- int nframes = 3
- int width = 512
- int height = 1024
- int firstrow = 0
- int lastrow = 1023
- list interlacing = [0, 0]
- int columns = 1
- bool padToFull = True
- float toffset = -165.76
- float tslope = 81.36

# 6.2.1 Detailed Description

Definition at line 30 of file daedalus.py.

### 6.2.2 Constructor & Destructor Documentation

### 6.2.2.1 \_\_init\_\_()

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

#### Definition at line 56 of file daedalus.py.

```
def __init__(self, ca):
00057
               self.ca = ca
00058
              super(daedalus, self).__init__(ca)
00059
              self.sens_registers = OrderedDict(
00060
00061
00062
                       "HST_READBACK_A_LO": "018",
00063
                       "HST_READBACK_A_HI": "019",
                       "HST_READBACK_B_LO": "01A",
00065
                       "HST_READBACK_B_HI": "01B"
00066
                       "HSTALLWEN_WAIT_TIME": "03F",
00067
                       "VRESET_HIGH_VALUE": "04A",
                       "FRAME_ORDER_SEL": "04B",
00069
                       "EXT_PHI_CLK_SHO_ON": "050"
00070
                       "EXT_PHI_CLK_SHO_OFF": "051",
00071
                       "EXT_PHI_CLK_SH1_ON": "052",
                       "EXT_PHI_CLK_SH1_OFF": "053",
00072
00073
                       "EXT_PHI_CLK_SH2_ON": "054",
00074
                       "HST_TRIGGER_DELAY_DATA_LO": "120",
00075
                       "HST_TRIGGER_DELAY_DATA_HI": "121",
                       "HST_PHI_DELAY_DATA": "122",
00076
                       "HST_EXT_CLK_HALF_PER": "129",
00077
00078
                       "HST_COUNT_TRIG": "130",
                       "HST_DELAY_EN": "131",
00079
                       "RSL_HFW_MODE_EN": "133"
00080
                       "RSL_ZDT_MODE_B_EN": "135",
"RSL_ZDT_MODE_A_EN": "136",
00081
00082
                       "BGTRIMA": "137",
"BGTRIMB": "138",
00083
00084
00085
                       "COLUMN TEST EN": "139",
                       "RSL_CONFIG_DATA_BO": "140",
00086
                       "RSL_CONFIG_DATA_B1": "141",
00087
                       "RSL_CONFIG_DATA_B2": "142",
00088
                       "RSL_CONFIG_DATA_B3": "143",
00089
00090
                       "RSL_CONFIG_DATA_B4": "144",
                       "RSL_CONFIG_DATA_B5": "145",
00091
                       "RSL_CONFIG_DATA_B6": "146",
00092
                       "RSL_CONFIG_DATA_B7": "147"
00093
                       "RSL_CONFIG_DATA_B8": "148",
00094
                       "RSL_CONFIG_DATA_B9": "149",
00095
00096
                       "RSL_CONFIG_DATA_B10": "14A"
                       "RSL_CONFIG_DATA_B11": "14B"
00097
                       "RSL_CONFIG_DATA_B12": "14C"
00098
                       "RSL_CONFIG_DATA_B13": "14D",
00099
                       "RSL_CONFIG_DATA_B14": "14E",
00100
                       "RSL_CONFIG_DATA_B15": "14F"
00101
                       "RSL_CONFIG_DATA_B16": "150"
00102
                       "RSL_CONFIG_DATA_B17": "151",
00103
00104
                       "RSL_CONFIG_DATA_B18": "152",
                       "RSL_CONFIG_DATA_B19": "153"
00105
00106
                       "RSL_CONFIG_DATA_B20": "154",
00107
                       "RSL_CONFIG_DATA_B21": "155",
00108
                       "RSL_CONFIG_DATA_B22": "156",
00109
                       "RSL_CONFIG_DATA_B23": "157",
00110
                       "RSL_CONFIG_DATA_B24": "158",
00111
                       "RSL_CONFIG_DATA_B25": "159",
00112
                       "RSL_CONFIG_DATA_B26": "15A",
00113
                       "RSL_CONFIG_DATA_B27": "15B",
00114
                       "RSL_CONFIG_DATA_B28": "15C",
                       "RSL_CONFIG_DATA_B29": "15D",
00115
00116
                       "RSL_CONFIG_DATA_B30": "15E",
                       "RSL_CONFIG_DATA_B31": "15F",
00117
                       "RSL_CONFIG_DATA_A0": "160",
00118
```

```
00119
                             "RSL_CONFIG_DATA_A1": "161",
00120
                             "RSL_CONFIG_DATA_A2": "162",
00121
                             "RSL_CONFIG_DATA_A3": "163",
                             "RSL_CONFIG_DATA_A4": "164",
00122
00123
                             "RSL_CONFIG_DATA_A5": "165",
                            "RSL_CONFIG_DATA_A6": "166",
00124
00125
                             "RSL_CONFIG_DATA_A7": "167"
00126
                            "RSL_CONFIG_DATA_A8": "168",
00127
                             "RSL_CONFIG_DATA_A9": "169",
                            "RSL_CONFIG_DATA_A10": "16A"
00128
                             "RSL_CONFIG_DATA_A11": "16B",
00129
                            "RSL_CONFIG_DATA_A12": "16C",
00130
00131
                             "RSL_CONFIG_DATA_A13": "16D",
                            "RSL_CONFIG_DATA_A14": "16E",
00132
00133
                             "RSL_CONFIG_DATA_A15": "16F",
                            "RSL_CONFIG_DATA_A16": "170",
00135
                             "RSL_CONFIG_DATA_A17": "171",
                             "RSL_CONFIG_DATA_A18": "172",
00136
00137
                             "RSL_CONFIG_DATA_A19": "173",
00138
                             "RSL_CONFIG_DATA_A20": "174",
                             "RSL_CONFIG_DATA_A21": "175",
00139
00140
                             "RSL_CONFIG_DATA_A22": "176",
00141
                             "RSL_CONFIG_DATA_A23": "177",
                             "RSL_CONFIG_DATA_A24": "178",
00142
                             "RSL_CONFIG_DATA_A25": "179",
00143
                             "RSL CONFIG DATA A26": "17A",
00144
                             "RSL_CONFIG_DATA_A27": "17B",
00145
                             "RSL_CONFIG_DATA_A28": "17C",
00146
                             "RSL_CONFIG_DATA_A29": "17D",
00147
                             "RSL_CONFIG_DATA_A30": "17E",
00148
                             "RSL_CONFIG_DATA_A31": "17F",
00149
00150
00151
                 )
00152
                  self.sens_subregisters = [
00153
00154
                       ("HST_MODE", "HS_TIMING_CTL", 0, 1, True), ("SLOWREADOFF_O", "CTRL_REG", 4, 1, True), ("SLOWREADOFF_1", "CTRL_REG", 5, 1, True), ("MANSHUT_MODE", "CTRL_REG", 8, 1, True),
00156
00157
00158
00159
                       ("INTERLACING_EN", "CTRL_REG", 9, 1, True),
00160
                       ("HFW", "RSL_HFW_MODE_EN", 0, 1, True),
("ZDT_A", "RSL_ZDT_MODE_A_EN", 0, 1, True),
("ZDT_B", "RSL_ZDT_MODE_B_EN", 0, 1, True),
00161
00162
00163
                       ("HST_DEL_EN", "HST_DELAY_EN", 0, 1, True), ("PHI_DELAY_A", "HST_PHI_DELAY_DATA", 9, 10, True), ("PHI_DELAY_B", "HST_PHI_DELAY_DATA", 29, 10, True),
00164
00165
00166
00167
                       # Assume that daedalus is not to be used with v1 board
00168
                       ("VRESET_HIGH", "VRESET_HIGH_VALUE", 15, 16, True),
00169
                       ("STAT_SHORISEUR", "STAT_REG", 3, 1, False),
("STAT_SHOFALLUR", "STAT_REG", 4, 1, False),
("STAT_RSLNALLWENA", "STAT_REG", 12, 1, False),
("STAT_RSLNALLWENB", "STAT_REG", 15, 1, False),
00172
00173
00174
00175
00176
                        # ("STAT_CONFIGHSTDONE", "STAT_REG", 16, 1, False),
00177
00178
```

# 6.2.3 Member Function Documentation

### 6.2.3.1 getManualTiming()

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

Definition at line 659 of file daedalus.py.

```
def getManualTiming(self):
00659
00660
00661
              Read off manual shutter timing settings
00662
              list of manual timing intervals
00663
00664
00665
              timing = []
              for reg in [
00666
00667
                  "EXT_PHI_CLK_SH0_ON",
00668
                  "EXT_PHI_CLK_SHO_OFF",
00669
                  "EXT_PHI_CLK_SH1_ON",
00670
                  "EXT_PHI_CLK_SH1_OFF",
00671
                  "EXT_PHI_CLK_SH2_ON",
00672
00673
                   _, reghex = self.ca.getRegister(reg)
00674
                  timing.append(25 * int(reghex, 16))
00675
              return timing
00676
```

### 6.2.3.2 getSensTemp()

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

Definition at line 677 of file daedalus.py.

```
00677
          def getSensTemp(self, scale=None, offset=None, slope=None, dec=1):
00678
00679
              Read temperature sensor located on the Daedalus sensor
00680
              Args:
00681
                  scale: temperature scale to report (defaults to C, options are F and K)
                  offset: offset of linear fit of sensor response (defaults to self.toffset)
00682
00683
                  slope: slope of linear fit of sensor response (defaults to self.tslope)
00684
                  dec: round to 'dec' digits after the decimal point
00685
00686
              Returns:
              temperature as float on given scale, rounded to .1 degree _{\tt m,m}
00687
00688
00689
              err, rval = self.ca.getMonV("MON_TSENSE_OUT", errflag=True)
00690
              if err:
00691
                  logging.error(
00692
                      self.logerr + "unable to retrieve temperature information ("
00693
                      'getTemp), returning "0" '
00694
00695
                  return 0.0
00696
              if offset is None:
                  offset = self.toffset
00697
00698
              if slope is None:
                  slope = self.tslope
00699
```

```
00700
00701
              ctemp = offset + slope * rval
00702
              if scale == "K":
00703
                  temp = round(ctemp + 273.15, dec)
00704
              elif scale == "F":
                  temp = round(1.8 \star ctemp + 32, dec)
00705
00706
              else:
00707
                  temp = round(ctemp, dec)
00708
              return temp
00709
```

#### 6.2.3.3 parseReadoff()

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

# Definition at line 710 of file daedalus.py.

```
00710
           def parseReadoff(self, frames, columns):
00711
00712
               Parses frames from board into images
00713
                   frames: list of data arrays (frames) returned from board
00714
00715
                   columns: 1 (full width image) or 2 (hemispheres generate distinct images)
               Returns:
00716
               list of data arrays (frames) reordered and deinterlaced
00717
00718
00719
               logging.debug(self.logdebug + "parseReadoff")
00720
               w = self.width
00721
               if hasattr(self, "ca"): # TODO: this may no longer be necessary
00722
                   padIt = self.ca.padToFull
00723
                   padIt = self.padToFull
00724
               if padIt:
00725
00726
                   rows = self.maxheight
00727
00728
                   rows = self.lastrow - self.firstrow + 1
00729
               parsed = []
00730
               for frame in frames:
00731
                   current = np.zeros((rows, w), dtype=np.uint16)
00732
                   mapped = np.zeros((rows, w), dtype=np.uint16)
00733
                   frame = frame.reshape(rows, w)
00734
00735
                    for entry in range(int(w / 2)):
00736
                        col = 32 * (entry % 8) + entry // 8 # lookup from daedlookup.xls
                        for row in range(rows):
00737
                            current[row][col] = frame[row][2 * entry]
current[row][col + 256] = frame[row][2 * entry + 1]
00738
00739
00740
00741
                   for row in range (rows):
                        mapped[row][0:32] = current[row][320:352]
mapped[row][32:64] = current[row][352:384]
00742
00743
00744
                        mapped[row][64:96] = current[row][192:224]
                        mapped[row][96:128] = current[row][160:192]
00745
00746
                        mapped[row][128:160] = current[row][256:288]
                        mapped[row][160:192] = current[row][288:320]
00747
                        mapped[row][192:224] = current[row][416:448]
mapped[row][224:256] = current[row][32:64]
00748
00749
                        mapped[row][256:288] = current[row][128:160]
00750
```

```
00751
                      mapped[row][288:320] = current[row][224:256]
00752
                      mapped[row][320:352] = current[row][384:416]
00753
                      mapped[row][352:384] = current[row][448:480]
00754
                      mapped[row][384:416] = current[row][480:512]
00755
                      mapped[row][416:448] = current[row][0:32]
00756
                      mapped[row][448:480] = current[row][64:96]
00757
                      mapped[row][480:512] = current[row][96:128]
00758
                  parsed.append(mapped)
00759
00760
              images = self.ca.partition(parsed, columns)
00761
              flatimages = [flattenlist(x) for x in images]
00762
              return flatimages
00763
```

### 6.2.3.4 reportStatusSensor()

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

### Definition at line 764 of file daedalus.py.

```
00764
          def reportStatusSensor(self, statusbits, statusbits2):
00765
00766
              Print status messages from sensor-specific bits of status register or object
00767
                status flags
00768
00769
00770
                  statusbits: result of checkStatus()
00771
                  statusbits2: result of checkStatus2()
00772
00773
              if int(statusbits[3]):
00774
                  print(self.loginfo + "SHO_rise_B_edge detected")
00775
              if int(statusbits[4]):
00776
                  print(self.loginfo + "SHO_fall_B_edge detected")
00777
              if int(statusbits[12]):
00778
                  print(self.loginfo + "RSLNALLWENB detected")
00779
              if int(statusbits[15]):
00780
                  print(self.loginfo + "RSLNALLWENA detected")
00781
              if self.HFW:
00782
                  print(self.loginfo + "High Full Well mode active")
00783
              # TODO: handle two hemispheres for ZDT
00784
              elif self.ZDT:
00785
                  print(self.loginfo + "Zero Dead Time mode active")
00786
              elif self.interlacing != [0, 0]:
00787
                  print(
00788
                      "{loginfo}Interlacing active: {interlacing}".format(
00789
                          loginfo=self.loginfo, interlacing=str(self.interlacing)
00790
00791
                  )
00792
              if self.ca.sensmanual == []:
00793
                  print(
00794
                       "{loginfo}High-speed timing: A:{Atiming}, B:{Btiming}".format(
00795
                          loginfo=self.loginfo,
00796
                          Atiming=self.getTiming(side="A", actual=True),
00797
                          Btiming=self.getTiming(side="B", actual=True),
00798
00799
00800
              else:
00801
                  print(
```

```
00802
                      "{loginfo}Manual timing set to {timing}".format(
00803
                          loginfo=self.loginfo, timing=self.getManualTiming()
00804
00805
00806
00807
00808 """
00809 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00810 LLNL-CODE-838080
00812 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00813 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00814 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00815 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00816 be made under this license.
00817 """
```

#### 6.2.3.5 selectOscillator()

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

#### Definition at line 444 of file daedalus.py.

```
def selectOscillator(self, osc=None):
00445
00446
               Selects oscillator to control sensor timing
00447
               Aras:
00448
                  osc: 500|100|'ring'|external', defaults to 500 MHz
00449
00450
               Returns:
               error message as string
00451
00452
00453
               logging.info(self.loginfo + "selectOscillator; osc = " + str(osc))
00454
               if osc is None:
00455
                   osc = 500
               osc = str(osc)
00456
               if osc[:3] == "500":
00457
                   payload = "00"
00458
               elif osc[:3] == "100":
payload = "01"
00459
00460
               elif osc.upper()[:3] == "RIN":
    payload = "10"
00461
00462
               elif osc.upper()[:3] in ["EXT"]:
    payload = "11"
00463
00464
00465
               else:
00466
                   err = (
00467
                        self.logerr + "selectOscillator: invalid parameter supplied. "
00468
                        "Oscillator selection is unchanged."
00469
00470
                   logging.error(err)
00471
                   return err
               self.ca.setSubregister("OSC_SELECT", payload)
00472
00473
```

### 6.2.3.6 sensorSpecific()

```
nsCamera.sensors.daedalus.daedalus.sensorSpecific (
                      self )
Returns:
      list of tuples, (Sensor-specific register, default setting)
Definition at line 180 of file daedalus.py.
             def sensorSpecific(self):
00180
00181
00182
                   Returns:
                   list of tuples, (Sensor-specific register, default setting)
00183
00184
00185
                   return [
                         ("FPA_FRAME_INITIAL", "00000000"),
("FFA_FRAME_FINAL", "00000002"),
("FPA_ROW_INITIAL", "00000000"),
("FPA_ROW_FINAL", "000003FF"),
00186
00187
00188
00189
                         ("HS_TIMING_DATA_ALO", "00006666"), # 0db6 = 2-1; 6666 = 2-2 ("HS_TIMING_DATA_AHI", "00000000"),
00190
00191
                         ("HS_TIMING_DATA_BLO", "00006666"), ("HS_TIMING_DATA_BHI", "00000000"),
00192
00193
                         ("FRAME_ORDER_SEL", "00000000"),
("RSL_HFW_MODE_EN", "00000000"),
00194
00195
                         ("RSL_HFW_MODE_EN", "00000000"),
("RSL_ZDT_MODE_B_EN", "00000000"),
("RSL_ZDT_MODE_A_EN", "00000000"),
("RSL_CONFIG_DATA_B0", "00000000"),
("RSL_CONFIG_DATA_B1", "00000000"),
00196
00197
00198
00199
                         ("RSL_CONFIG_DATA_B2", "00000000"), ("RSL_CONFIG_DATA_B3", "00000000"),
00200
00201
                         ("RSL_CONFIG_DATA_B4", "00000000"), ("RSL_CONFIG_DATA_B5", "00000000"),
00202
00203
                         ("RSL_CONFIG_DATA_B6", "00000000"),
00204
00205
                         ("RSL_CONFIG_DATA_B7", "00000000"),
                         ("RSL_CONFIG_DATA_B8", "00000000"), ("RSL_CONFIG_DATA_B9", "00000000"),
00206
00207
                         ("RSL_CONFIG_DATA_B10", "00000000"), ("RSL_CONFIG_DATA_B11", "00000000"), ("RSL_CONFIG_DATA_B12", "00000000"), ("RSL_CONFIG_DATA_B12", "00000000"), ("RSL_CONFIG_DATA_B13", "00000000"),
00208
00209
00210
00211
                         ("RSL_CONFIG_DATA_B14", "00000000"),
00212
                         ("RSL_CONFIG_DATA_B15", "00000000"), ("RSL_CONFIG_DATA_B16", "000000000"),
00213
00214
00215
                          ("RSL_CONFIG_DATA_B17", "00000000"),
                         ("RSL_CONFIG_DATA_B18", "00000000"), ("RSL_CONFIG_DATA_B19", "000000000"),
00216
00217
                         ("RSL_CONFIG_DATA_B20", "00000000"),
("RSL_CONFIG_DATA_B21", "00000000"),
00218
00219
00220
                         ("RSL_CONFIG_DATA_B22", "00000000"),
00221
                         ("RSL_CONFIG_DATA_B23", "00000000"),
00222
                         ("RSL_CONFIG_DATA_B24", "00000000"),
00223
                         ("RSL_CONFIG_DATA_B25", "00000000"),
                         ("RSL_CONFIG_DATA_B26", "00000000"),
                         ("RSL_CONFIG_DATA_B27", "00000000"),
00225
                         ("RSL_CONFIG_DATA_B28", "00000000"),
00226
00227
                         ("RSL_CONFIG_DATA_B29", "00000000"),
                         ("RSL_CONFIG_DATA_B30", "00000000"),
                         ("RSL_CONFIG_DATA_B31", "00000000"), ("RSL_CONFIG_DATA_A0", "00000000"),
00229
00230
                         ("RSL_CONFIG_DATA_A1", "00000000"),
00231
00232
                         ("RSL_CONFIG_DATA_A2",
                                                          "00000000"),
                                                          "00000000"),
00233
                         ("RSL_CONFIG_DATA_A3",
00234
                         ("RSL_CONFIG_DATA_A4",
                                                          "00000000"),
                                                          "00000000"),
00235
                         ("RSL_CONFIG_DATA_A5",
00236
                         ("RSL_CONFIG_DATA_A6", "00000000"),
00237
                         ("RSL_CONFIG_DATA_A7",
                                                          "00000000"),
00238
                         ("RSL_CONFIG_DATA_A8", "00000000"),
                         ("RSL_CONFIG_DATA_A9", "00000000"), ("RSL_CONFIG_DATA_A10", "00000000"),
00239
00240
                         ("RSL_CONFIG_DATA_A11", "00000000"),
00241
                         ("RSL_CONFIG_DATA_A12", "00000000"),
00242
                         ("RSL_CONFIG_DATA_A13", "00000000"), ("RSL_CONFIG_DATA_A14", "000000000"),
00243
00244
                          ("RSL_CONFIG_DATA_A15", "00000000"),
00245
```

```
00246
                                   ("RSL_CONFIG_DATA_A16", "00000000"),
                                   ("RSL_CONFIG_DATA_A17", "00000000"), ("RSL_CONFIG_DATA_A18", "00000000"), ("RSL_CONFIG_DATA_A19", "00000000"),
00247
00248
00249
                                   ("RSL_CONFIG_DATA_A20", "00000000"), ("RSL_CONFIG_DATA_A21", "00000000"),
00250
00251
                                   ("RSL_CONFIG_DATA_A22", "00000000"), ("RSL_CONFIG_DATA_A23", "000000000"),
00252
00253
                                  ("RSL_CONFIG_DATA_A24", "000000000"),
("RSL_CONFIG_DATA_A25", "000000000"),
("RSL_CONFIG_DATA_A25", "00000000"),
("RSL_CONFIG_DATA_A26", "00000000"),
00254
00255
00256
00257
                                  ("RSL_CONFIG_DATA_A28", "00000000"),
("RSL_CONFIG_DATA_A29", "00000000"),
("RSL_CONFIG_DATA_A30", "00000000"),
("RSL_CONFIG_DATA_A31", "00000000"),
00258
00260
00261
                                   ("HST_TRIGGER_DELAY_DATA_LO", "00000000"), ("HST_TRIGGER_DELAY_DATA_HI", "00000000"),
00262
00263
                                   ("HST_PHI_DELAY_DATA", "00000000"),
00264
                                   ("SLOWREADOFF_0", "0"),
00265
                                   ("SLOWREADOFF_1", "0"),
00266
00267
                          ]
00268
```

### 6.2.3.7 setExtClk()

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

### Definition at line 556 of file daedalus.py.

```
00556
          def setExtClk(self, dilation=None, frequency=None):
00557
00558
              Override the standard board clock with the external clock.
00559
00560
                  dilation: ratio of base frequency (500 MHz) to desired external clock
                   frequency. Default is 25. Overridden if frequency parameter is provided
00561
00562
                  frequency: Desired frequency for phi clock.
00563
              Returns:
             error message as string
00564
00565
00566
              logging.debug(
00567
                 self.logdebug
                  + "setExtClk; dilation = "
00568
00569
                  + str(dilation)
                  + "; frequency = "
00570
00571
                  + str(frequency)
00572
00573
              if not (isinstance(frequency, int) or isinstance(frequency, float)):
00574
                  err = (
00575
                     self.logerr
00576
                      + "invalid external clock frequency submitted. Clock is not "
                      + "operating"
00577
00578
00579
                  logging.error(err)
00580
                  return err
              self.ca.selectOscillator("external")
00581
00582
              if not dilation:
```

```
00583
                  dilation = 25
00584
              if not frequency:
                  frequency = 5e7 / float(dilation)
00585
00586
              count = 2e7 / float(frequency) - 1 # base phi clock is 20 MHz?
00587
              if count < 0:
00588
00589
                  warn = (
00590
                     self.logwarn
00591
                      + "external clock frequency exceeding maximum. Frequency set to "
00592
                      + "maximum (20 MHz)"
00593
00594
                  logging.warning(warn)
00595
              if count > 0xFFFFFFFF:
                  count = 0xFFFFFFFF
00596
00597
              counthex = hex(int(count))[2:].zfill(8)
00598
              self.ca.setRegister("HST_EXT_CLK_HALF_PER", counthex)
00599
```

### 6.2.3.8 setHighFullWell()

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

```
Definition at line 354 of file daedalus.py.
```

```
00354
          def setHighFullWell(self, flag):
00355
00356
              Activates High Full Well mode. All frames are acquired simultaneously. Zero Dead
00357
                Time mode and interlacing will be automatically deactivated and column number
00358
                will be reset to 0. NOTE: after deactivating HFW, the board remains in
00359
                uninterlaced mode (interlacing = 0)
00360
00361
00362
                 flag: True to activate HFW mode, False to deactivate
00363
00364
              Returns:
              Error message
00365
00366
00367
              logging.debug(self.logdebug + "setHighFullWell; flag = " + str(flag))
00368
              err0 = ""
00369
              if flag:
00370
                  if self.ZDT:
00371
                      logging.warning(
00372
                          self.logwarn + "ZDT mode will be disengaged because of HFW "
00373
                          "setting "
00374
00375
                      err0 = self.setZeroDeadTime(False)
                  err1, _ = self.ca.setSubregister("HFW", "1")
00376
                  self.HFW = False # preclude HFW deactivation message in setInterlacing
00377
00378
                  self.setInterlacing(0)
00379
                  self.HFW = True
00380
                  logging.info(self.loginfo + "High Full Well mode active")
00381
              else:
00382
                 self.HFW = False
00383
                  errl, _ = self.ca.setSubregister("HFW", "0")
                  logging.info(self.loginfo + "High Full Well mode inactivate")
00384
00385
              err = err0 + err1
00386
              if err:
00387
                  logging.error(self.logerr + "HFW option may not be set correctly ")
00388
              return err
00389
```

### 6.2.3.9 setInterlacing()

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

### Definition at line 269 of file daedalus.py.

```
def setInterlacing(self, ifactor=None, side=None):
00269
00270
00271
              Sets interlacing factor. NOTE: if called directly when HFW or ZDT mode is active,
00272
              this will disengage those modes automatically. If hemispheres have different
00273
              factors when the image is acquired, the resulting frames are separated into
00274
              half-width images
00275
00276
              Args:
00277
                  ifactor: number of interlaced lines (generates ifactor + 1 images per frame)
00278
                    defaults to 0 (no interlacing)
00279
                  side: identify particular hemisphere (A or B) to control. If left blank,
00280
                    control both hemispheres
00281
00282
              Returns:
              integer: active interlacing factor (unchanged if error) _{\mbox{\scriptsize mum}}
00283
00284
00285
              logging.debug(self.logdebug + "setInterlacing; ifactor = " + str(ifactor))
00286
              if ifactor is None:
00287
                  ifactor = 0
00288
00289
                  not isinstance(ifactor, int)
                  or ifactor < 0</pre>
00290
00291
                  or ifactor > (self.maxheight - 1)
00292
00293
00294
                       self.logerr + "invalid interlacing factor submitted. "
00295
                       "Interlacing remains unchanged. "
00296
00297
                  logging.error(err)
00298
                  return self.interlacing
              if self.HFW:
00299
00300
                  logging.warning(
00301
                      self.logwarn + "HFW mode will be disengaged because of new "
00302
                       "interlacing setting '
00303
00304
                  self.setHighFullWell(False)
00305
              if self.ZDT:
00306
                  logging.warning(
00307
                       self.logwarn + "ZDT mode will be disengaged because of new "
00308
                       "interlacing setting '
00309
                  self.setZeroDeadTime(False)
00310
00311
              if ifactor == 0:
00312
                  bitscheme = self.maxheight * [0]
                  # deactivating one side shouldn't turn off enable for both sides
00313
00314
                  # TODO: is it a problem if sides are set separately, so interlacing is zero
00315
                       but still enabled?
```

```
00316
                    if side is None:
00317
                        self.ca.setSubregister("INTERLACING_EN", "0")
00318
               else:
                pattern = [0] + ifactor * [1]
00319
00320
                    reps = 1 + self.maxheight // (ifactor + 1)
00321
                   bitscheme = (reps * pattern)[0 : self.maxheight]
00322
                   self.ca.setSubregister("INTERLACING_EN", "1")
00323
00324
               for regnum in range(32):
                    regbits = bitscheme[32 * regnum : 32 * (regnum + 1)]
00325
00326
                    logging.debug(self.logdebug + "regbits = " + str(regbits))
                    # generated pattern is reverse order from placement in register (element 0
00327
00328
                       of the list is the LSB of the register)
                 bitsrev = regbits[::-1]
00329
00330
                   s = [str(i) for i in bitsrev]
00331
                   b = "".join(s) # assemble as binary number for processing
                   hexval = "%x" % int(b, 2)
00332
                   val = hexval.zfill(8)
00333
                   err0 = ""
00334
                   err1 = ""
00335
                   if side is None or side.lower() == "a":
    lname = "RSL_CONFIG_DATA_A" + str(regnum)
00336
00337
00338
                        err1, _ = self.ca.setRegister(lname, val)
00339
                        self.interlacing[1] = ifactor
                   if side is None or side.lower() == "b":
    rname = "RSL_CONFIG_DATA_B" + str(regnum)
00340
00341
                        err0, _ = self.ca.setRegister(rname, val)
00342
00343
                        self.interlacing[0] = ifactor
00344
                   err = err + err0 + err1
00345
             if err:
               logging.error(self.logerr + "interlacing may not be set correctly: " + err) logging.info(self.loginfo + "Interlacing set to " + str(self.interlacing))
00346
00347
               if self.interlacing[0] == self.interlacing[1]:
00348
00349
                   self.columns = 1
00350
               else:
00351
                   self.columns = 2
00352
               return self.interlacing
00353
```

# 6.2.3.10 setManualShutters()

```
nsCamera.sensors.daedalus.daedalus.setManualShutters (
              self.
              timing = None)
Legacy alias for setManualTiming()
```

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

```
Definition at line 602 of file daedalus.py.

Od602 def setManualShutters(self, timing=None):
00603
00604
                  Legacy alias for setManualTiming()
00605
00606
                  self.setManualTiming(timing)
00607
```

# 6.2.3.11 setManualTiming()

```
nsCamera.sensors.daedalus.daedalus.setManualTiming (
              self.
              timing = None )
```

```
Manual shutter timing, five intervals given in nanoseconds, e.g.,
  [100,50,100,50,100] for frame 0 open for 100 ns, an interframe pause of 50 ns,
  frame 1 open for 100 ns, etc. Timing is set for both hemispheres.
The actual timing is rounded down to the nearest multiple of 25 ns. (Each
  count = 25 ns. e.g., a request for 140 ns rounds down to a count of '5',
  which corresponds to 125 ns))
   - Minimum timing is 75 ns
   - Maximum is 25 * 2^30 ns (approximately 27 seconds)
Args:
    timing: 5-element list in nanoseconds
Returns:
    tuple (error string, response string from final message)
```

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

```
Definition at line 608 of file daedalus.py.
```

```
00608
            def setManualTiming(self, timing=None):
00609
                Manual shutter timing, five intervals given in nanoseconds, e.g.,
00610
00611
                   [100,50,100,50,100] for frame 0 open for 100 ns, an interframe pause of 50 ns,
00612
                    frame 1 open for 100 ns, etc. Timing is set for both hemispheres.
00613
                The actual timing is rounded down to the nearest multiple of 25 ns. (Each
00614
                   count = 25 ns. e.g., a request for 140 ns rounds down to a count of '5',
00615
00616
                   which corresponds to 125 ns))
                     - Minimum timing is 75 ns - Maximum is 25 \star 2^30 ns (approximately 27 seconds)
00617
00618
00619
00620
                Args:
00621
                     timing: 5-element list in nanoseconds
00622
00623
                Returns:
                tuple (error string, response string from final message)
00624
00625
00626
                if timing is None:
00627
                     logging.info(
00628
                          self.loginfo
00629
                           + "No manual timing setting provided, defaulting to (100, 150, 100, "
00630
                          " 150, 100, 150, 100) for both hemispheres"
00631
00632
                     timing = [(100, 150, 100, 150, 100)]
00633
00634
                logging.info(self.loginfo + "Manual shutter sequence: " + str(timing))
00635
                 flattened = flattenlist(timing)
00636
00637
                     len(flattened) != 5
00638
                     or not all(isinstance(x, (int, float)) for x in flattened)
00639
                     or not all(x \ge 25 for x in flattened)
00640
00641
                     err = self.logerr + "Invalid manual shutter timing list: " + str(timing)
00642
                     logging.error(err + "; timing settings unchanged")
                     return err, "00000000"
00643
00644
00645
                timecounts = [int(a // 25) for a in flattened]
                self.ca.sensmanual = timing
00646
00647
                self.ca.senstiming = {} # clear HST settings from ca object
00648
                control_messages = [
00649
00650
                     ("MANSHUT_MODE", "1"),
                     ("EXT_PHI_CLK_SHO_ON", "{0:#0{1}x}".format(timecounts[0], 10)[2:10]),
("EXT_PHI_CLK_SHO_OFF", "{0:#0{1}x}".format(timecounts[1], 10)[2:10]),
("EXT_PHI_CLK_SH1_ON", "{0:#0{1}x}".format(timecounts[2], 10)[2:10]),
("EXT_PHI_CLK_SH1_OFF", "{0:#0{1}x}".format(timecounts[3], 10)[2:10]),
("EXT_PHI_CLK_SH2_ON", "{0:#0{1}x}".format(timecounts[4], 10)[2:10]),
00651
00652
00653
00654
00655
00656
00657
                return self.ca.submitMessages(control_messages, " setManualShutters: ")
00658
```

### 6.2.3.12 setPhiDelay()

nsCamera.sensors.daedalus.daedalus.setPhiDelay (

```
self,
    side = None,
    delay = 0 )

Use phi delay timer. Actual delay is rounded down to multiple of .15 ns, up to a
    maximum delay of 1.5 ns
Args:
    side: hemisphere to delay; if None, delay both hemispheres
    delay: phi delay in ns

Returns:
    String of errors, if any
```

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

Definition at line 516 of file daedalus.py.

```
00516
            def setPhiDelay(self, side=None, delay=0):
00517
                  Use phi delay timer. Actual delay is rounded down to multiple of .15 ns, up to a maximum delay of 1.5 ns \,
00518
00519
00520
                  Args:
                       side: hemisphere to delay; if None, delay both hemispheres
00521
00522
                       delay: phi delay in ns
00523
00524
                  Returns:
                  String of errors, if any
00525
00526
                  logging.debug(self.logdebug + "setPhiDelay; delay = " + str(delay))
00527
00528
                       not (isinstance(delay, int) or isinstance(delay, float))
00529
00530
                       or delay < 0
00531
                       or delay > 1.5
00532
                 ):
00533
                       err = (
                            self.logerr + "invalid phi delay submitted. Delay remains " "unchanged. "  
00534
00535
00536
00537
                       logging.error(err)
00538
                       return err
00539
                  delayblocks = int(delay / 0.15)
00540
                  if delayblocks < 0:</pre>
00541
                       delayblocks = 0
00542
                  if delayblocks > 10:
00543
                       delayblocks = 10
                  \label{eq:delayseq} \begin{array}{lll} \texttt{delayseq} = (10 - \texttt{delayblocks}) \; \star \; \texttt{[0]} \; + \; \texttt{delayblocks} \; \star \; \texttt{[1]} \\ \texttt{seqstr} = \; \texttt{"".join}(\texttt{str}(\texttt{x}) \; \; \texttt{for} \; \; \texttt{x} \; \; \texttt{in} \; \; \texttt{delayseq}) \end{array}
00544
00545
                 err1 = ""
err2 = ""
00546
00547
00548
                  if side is None or side.upper() == "A":
00549
                       err1, _ = self.ca.setSubregister("PHI_DELAY_A", seqstr)
00550
                  if side is None or side.upper() == "B":
                  err2, _ = self.ca.setSubregister("PHI_DELAY_B", seqstr)
delayed = delayblocks * 0.15
00551
00553
                  logging.info(self.loginfo + "Actual phi delay = " + str(delayed) + " ns")
00554
                  return err1 + err2
00555
```

### 6.2.3.13 setTriggerDelay()

```
nsCamera.sensors.daedalus.daedalus.setTriggerDelay ( self, \\ delay = 0 \ ) Use trigger delay timer. Actual delay is rounded down to multiple of .15 ns, up to a maximum delay of 6 ns
```

```
Args:
delay: trigger delay in ns

Returns:
String of errors, if any
```

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

```
Definition at line 474 of file daedalus.py.
```

```
00474
          def setTriggerDelay(self, delay=0):
00475
00476
              Use trigger delay timer. Actual delay is rounded down to multiple of .15 ns, up
00477
                to a maximum delay of 6 ns
00478
00479
00480
                  delay: trigger delay in ns
00481
00482
              Returns:
              String of errors, if any
00483
00484
00485
              logging.debug(self.logdebug + "setTriggerDelay; delay = " + str(delay))
00486
00487
                  not (isinstance(delay, int) or isinstance(delay, float))
                  or delay < 0</pre>
00488
00489
                  or delay > 6
00490
00491
00492
                       self.logerr + "invalid trigger delay submitted. Delay remains "
00493
00494
00495
                  logging.error(err)
00496
                   return err
00497
              delayblocks = int(delay / 0.15)
00498
              if delayblocks < 0:
00499
                  delayblocks = 0
00500
              if delayblocks > 40:
00501
                  delayblocks = 40
00502
              delayseq = (40 - delayblocks) * [0] + delayblocks * [1]
              seqstr = "".join(str(x) for x in delayseq)
seqhex = "%x" % int(seqstr, 2)
00503
00504
              logging.debug(self.logdebug + "seqhex = " + str(seqhex))
00505
              highpart = seqhex[-10:-8].zfill(8)
00506
00507
              lowpart = seghex[-8:].zfill(8)
00508
              self.ca.setSubregister("HST_DEL_EN", "1")
              err0, _ = self.ca.setRegister("HST_TRIGGER_DELAY_DATA_LO", lowpart)
00509
              err1, _ = self.ca.setRegister("HST_TRIGGER_DELAY_DATA_HI", highpart)
00510
              err2, _ = self.ca.setSubregister("HST_MODE", "1")
00511
00512
              delayed = delayblocks * 0.15
              logging.info(self.loginfo + "Actual trigger delay = " + str(delayed) + " ns")
00513
00514
              return err0 + err1 + err2
00515
```

### 6.2.3.14 setZeroDeadTime()

```
control both hemispheres
Returns:
    Error message
```

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

Definition at line 390 of file daedalus.py.

```
def setZeroDeadTime(self, flag=True, side=None):
00391
              Activates Zero Dead Time mode. Even rows follow the assigned HST schedule; odd rows are acquired while the 'shutter' for the even rows are closed. High Full
00392
00393
00394
                Well mode and interlacing will be automatically deactivated.
00395
              *NOTE* after deactivating ZDT, the board reverts to uninterlaced mode
00396
                (interlacing = 0)
00397
00398
              Args:
00399
                  flag: True to activate ZDT mode, False to deactivate
00400
                  side: identify particular hemisphere (A or B) to control. If left blank,
                    control both hemispheres
00401
00402
00403
              Returns:
              Error message
00404
00405
00406
              logging.debug(self.logdebug + "setZeroDeadTime; flag = " + str(flag))
00407
              err0 =
              err1 = ""
00408
              err2 = ""
00409
00410
              if flag:
00411
                  if self.HFW:
00412
                      logging.warning(
                           self.logwarn + "HFW mode will be disengaged because of ZDT "
00413
00414
                           "setting "
00415
00416
                      err0 = self.setHighFullWell(False)
00417
                  if side is None or side.lower() == "a":
                      err2, _ = self.ca.setSubregister("ZDT_A", "1")
00418
00419
                      self.interlacing[0] = 1
                  if side is None or side.lower() == "b":
00420
                      errl, _ = self.ca.setSubregister("ZDT_B", "1")
00421
00422
                      self.interlacing[1] = 1
00423
                  # self.ZDT = False # preclude ZDT deactivation message in setInterlacing
00424
                  # if self.interlacing != [0, 0]:
00425
                        self.setInterlacing(0)
                  # TODO: need to handle flags when ZDT active for just one side
00426
00427
                  self.ZDT = True
00428
                  logging.info(
                      self.loginfo + "Zero Dead Time mode active; actual interlacing = 1"
00429
00430
00431
             else:
00432
                 self.ZDT = False
00433
                  if side is None or side.lower() == "a":
00434
                      err2, _ = self.ca.setSubregister("ZDT_A", "0")
00435
                  if side is None or side.lower() == "b":
00436
                      err1,
                              = self.ca.setSubregister("ZDT_B", "0")
00437
                  self.setInterlacing(0)
00438
                  logging.info(self.loginfo + "Zero Dead Time mode inactivate")
00439
              err = err0 + err1 + err2
00440
00441
                  logging.error(self.logerr + "ZDT option may not be set correctly ")
00442
              return err
00443
```

### 6.2.4 Member Data Documentation

### 6.2.4.1 bytesperpixel

```
int nsCamera.sensors.daedalus.daedalus.bytesperpixel = 2 [static]
```

Definition at line 36 of file daedalus.py.

# 6.2.4.2 ca

nsCamera.sensors.daedalus.daedalus.ca

Definition at line 57 of file daedalus.py.

### 6.2.4.3 columns [1/2]

```
int nsCamera.sensors.daedalus.daedalus.columns = 1 [static]
```

Definition at line 51 of file daedalus.py.

# 6.2.4.4 columns [2/2]

nsCamera.sensors.daedalus.daedalus.columns

Definition at line 349 of file daedalus.py.

#### 6.2.4.5 detect

```
str nsCamera.sensors.daedalus.daedalus.detect = "DAEDALUS_DET" [static]
```

Definition at line 38 of file daedalus.py.

### 6.2.4.6 firstframe

```
int nsCamera.sensors.daedalus.daedalus.firstframe = 0 [static]
```

Definition at line 43 of file daedalus.py.

# 6.2.4.7 firstrow

```
int nsCamera.sensors.daedalus.daedalus.firstrow = 0 [static]
```

Definition at line 48 of file daedalus.py.

# 6.2.4.8 fpganumID

```
int nsCamera.sensors.daedalus.daedalus.fpganumID = 2 [static]
```

Definition at line 37 of file daedalus.py.

# 6.2.4.9 height

```
int nsCamera.sensors.daedalus.daedalus.height = 1024 [static]
```

Definition at line 47 of file daedalus.py.

### 6.2.4.10 HFW [1/2]

```
bool nsCamera.sensors.daedalus.daedalus.HFW = False [static]
```

Definition at line 42 of file daedalus.py.

### 6.2.4.11 HFW [2/2]

nsCamera.sensors.daedalus.daedalus.HFW

Definition at line 377 of file daedalus.py.

# 6.2.4.12 interlacing [1/2]

```
list nsCamera.sensors.daedalus.daedalus.interlacing = [0, 0] [static]
```

Definition at line 50 of file daedalus.py.

#### 6.2.4.13 interlacing [2/2]

nsCamera.sensors.daedalus.daedalus.interlacing

Definition at line 347 of file daedalus.py.

### 6.2.4.14 lastframe

```
int nsCamera.sensors.daedalus.daedalus.lastframe = 2 [static]
```

Definition at line 44 of file daedalus.py.

### 6.2.4.15 lastrow

```
int nsCamera.sensors.daedalus.daedalus.lastrow = 1023 [static]
```

Definition at line 49 of file daedalus.py.

# 6.2.4.16 loglabel

```
str nsCamera.sensors.daedalus.daedalus.loglabel = "[Daedalus] " [static]
```

Definition at line 40 of file daedalus.py.

#### 6.2.4.17 maxframe

```
int nsCamera.sensors.daedalus.daedalus.maxframe = 2 [static]
```

Definition at line 33 of file daedalus.py.

### 6.2.4.18 maxheight

```
int nsCamera.sensors.daedalus.daedalus.maxheight = 1024 [static]
```

Definition at line 35 of file daedalus.py.

#### 6.2.4.19 maxwidth

```
int nsCamera.sensors.daedalus.daedalus.maxwidth = 512 [static]
```

Definition at line 34 of file daedalus.py.

### 6.2.4.20 minframe

```
int nsCamera.sensors.daedalus.daedalus.minframe = 0 [static]
```

Definition at line 32 of file daedalus.py.

# 6.2.4.21 nframes

```
int nsCamera.sensors.daedalus.daedalus.nframes = 3 [static]
```

Definition at line 45 of file daedalus.py.

# 6.2.4.22 padToFull

```
bool nsCamera.sensors.daedalus.daedalus.padToFull = True [static]
```

Definition at line 52 of file daedalus.py.

# 6.2.4.23 sens\_registers

nsCamera.sensors.daedalus.daedalus.sens\_registers

Definition at line 60 of file daedalus.py.

# 6.2.4.24 sens\_subregisters

nsCamera.sensors.daedalus.daedalus.sens\_subregisters

Definition at line 153 of file daedalus.py.

### 6.2.4.25 sensfam

```
str nsCamera.sensors.daedalus.daedalus.sensfam = "Daedalus" [static]
```

Definition at line 39 of file daedalus.py.

# 6.2.4.26 specwarn

```
str nsCamera.sensors.daedalus.daedalus.specwarn = "" [static]
```

Definition at line 31 of file daedalus.py.

### 6.2.4.27 toffset

```
float nsCamera.sensors.daedalus.daedalus.toffset = -165.76 [static]
```

Definition at line 53 of file daedalus.py.

# 6.2.4.28 tslope

```
float nsCamera.sensors.daedalus.daedalus.tslope = 81.36 [static]
```

Definition at line 54 of file daedalus.py.

### 6.2.4.29 width

```
int nsCamera.sensors.daedalus.daedalus.width = 512 [static]
```

Definition at line 46 of file daedalus.py.

# 6.2.4.30 ZDT [1/2]

bool nsCamera.sensors.daedalus.daedalus.ZDT = False [static]

Definition at line 41 of file daedalus.py.

### 6.2.4.31 ZDT [2/2]

nsCamera.sensors.daedalus.daedalus.ZDT

Definition at line 427 of file daedalus.py.

The documentation for this class was generated from the following file:

• C:/Users/hill35/git/camera\_python/nsCamera/sensors/daedalus.py

# 6.3 nsCamera.utils.misc.fakeCA Class Reference

#### **Public Member Functions**

- \_\_init\_\_ (self, sensorname="icarus2", firstframe=0, lastframe=3, width=512, height=1024, padToFull=True, firstrow=0, lastrow=1023, maxwidth=512, maxheight=1024, bytesperpixel=2, interlacing=None, columns=1, logtag=None)
- partition (self, parsed, columns)

### **Public Attributes**

- sensorname
- boardname
- padToFull
- logtag
- logcritbase
- logerrbase
- logwarnbase
- loginfobase
- logdebugbase
- logcrit
- logerr
- logwarn
- loginfo
- logdebug
- sensor

# 6.3.1 Detailed Description

```
Fake 'cameraAssembler' object to use as a parameter object in offline functions. Returned by generateFrames(), it contains the frame details required to save and plot images.
```

Definition at line 203 of file misc.py.

### 6.3.2 Constructor & Destructor Documentation

### 6.3.2.1 init ()

```
nsCamera.utils.misc.fakeCA.__init__ (
              self,
              sensorname = "icarus2",
              firstframe = 0,
              lastframe = 3,
              width = 512,
              height = 1024,
              padToFull = True,
              firstrow = 0,
              lastrow = 1023,
              maxwidth = 512,
              maxheight = 1024,
              bytesperpixel = 2,
              interlacing = None,
              columns = 1,
              logtag = None )
```

### Definition at line 210 of file misc.py.

```
00226
00227
              self.sensorname = sensorname
00228
              self.boardname = None
              self.padToFull = padToFull
00229
00230
              if logtag is None:
                  self.logtag = ""
00231
00232
              else:
00233
                   self.logtag = logtag
00234
              self.logcritbase = "CRITICAL" + self.logtag + ": "
self.logerrbase = "ERROR" + self.logtag + ": "
00235
00236
              self.logwarnbase = "WARNING" + self.logtag + ": "
00237
00238
              self.loginfobase = "INFO" + self.logtag + ": "
              self.logdebugbase = "DEBUG" + self.logtag + ": "
00239
00240
00241
              self.logcrit = self.logcritbase + "[FS] "
00242
              self.logerr = self.logerrbase + "[FS]
00243
              self.logwarn = self.logwarnbase + "[FS] "
00244
              self.loginfo = self.loginfobase + "[FS] "
00245
              self.logdebug = self.logdebugbase + "[FS] "
00246
              if self.sensorname == "icarus":
00247
00248
                   import nsCamera.sensors.icarus as snsr
00249
              elif self.sensorname == "icarus2":
00250
                   import nsCamera.sensors.icarus2 as snsr
              elif self.sensorname == "daedalus":
00251
00252
                   import nsCamera.sensors.daedalus as snsr
00253
00254
              self.sensor = snsr(self)
00255
```

# 6.3.3 Member Function Documentation

#### 6.3.3.1 partition()

```
nsCamera.utils.misc.fakeCA.partition (

self,

parsed,

columns)

Definition at line 256 of file misc.py.

00256 def partition(self, parsed, columns):
00257 # return nsCamera.utils.misc.partition(self, parsed, columns)

return partition(self, parsed, columns)

00258

00259

00260
```

### 6.3.4 Member Data Documentation

### 6.3.4.1 boardname

nsCamera.utils.misc.fakeCA.boardname

Definition at line 228 of file misc.py.

# 6.3.4.2 logcrit

nsCamera.utils.misc.fakeCA.logcrit

Definition at line 241 of file misc.py.

### 6.3.4.3 logcritbase

```
nsCamera.utils.misc.fakeCA.logcritbase
```

Definition at line 235 of file misc.py.

### 6.3.4.4 logdebug

```
nsCamera.utils.misc.fakeCA.logdebug
```

Definition at line 245 of file misc.py.

# 6.3.4.5 logdebugbase

nsCamera.utils.misc.fakeCA.logdebugbase

Definition at line 239 of file misc.py.

# 6.3.4.6 logerr

nsCamera.utils.misc.fakeCA.logerr

Definition at line 242 of file misc.py.

# 6.3.4.7 logerrbase

nsCamera.utils.misc.fakeCA.logerrbase

Definition at line 236 of file misc.py.

### 6.3.4.8 loginfo

nsCamera.utils.misc.fakeCA.loginfo

Definition at line 244 of file misc.py.

# 6.3.4.9 loginfobase

nsCamera.utils.misc.fakeCA.loginfobase

Definition at line 238 of file misc.py.

### 6.3.4.10 logtag

nsCamera.utils.misc.fakeCA.logtag

Definition at line 231 of file misc.py.

# 6.3.4.11 logwarn

nsCamera.utils.misc.fakeCA.logwarn

Definition at line 243 of file misc.py.

# 6.3.4.12 logwarnbase

 $\verb|nsCamera.utils.misc.fake CA.logwarn base|\\$ 

Definition at line 237 of file misc.py.

# 6.3.4.13 padToFull

```
nsCamera.utils.misc.fakeCA.padToFull
```

Definition at line 229 of file misc.py.

#### 6.3.4.14 sensor

```
nsCamera.utils.misc.fakeCA.sensor
```

Definition at line 254 of file misc.py.

#### 6.3.4.15 sensorname

```
nsCamera.utils.misc.fakeCA.sensorname
```

Definition at line 227 of file misc.py.

The documentation for this class was generated from the following file:

· C:/Users/hill35/git/camera\_python/nsCamera/utils/misc.py

# 6.4 nsCamera.comms.GigE.GigE Class Reference

### Classes

• class ZESTETM1\_CARD\_INFO

### **Public Member Functions**

- \_\_init\_\_ (self, camassem)
- sendCMD (self, pkt)
- arm (self, mode)
- readFrames (self, waitOnSRAM, timeout=0, fast=False, columns=1)
- readoff (self, waitOnSRAM, timeout=0, fast=False, columns=1)
- writeSerial (self, outstring, timeout)
- readSerial (self, size, timeout=None)
- openDevice (self, timeout=30)
- closeDevice (self)
- getCardIP (self)
- getCardInfo (self)

#### **Public Attributes**

- ca
- · logcrit
- logerr
- logwarn
- loginfo
- logdebug
- mode
- writeTimeout
- readTimeout
- · payloadsize
- skipError
- ZErrorDict
- dport
- closecard
- CardInfo
- CardInfoP
- ZCountCards
- ZESTETM1\_CARD\_INFO
- ZOpenConnection
- ZWriteData
- ZReadData
- Connection

#### **Protected Attributes**

• \_zest

### 6.4.1 Detailed Description

```
single OT card; to use multiple cards, instantiate multiple cameraAssembler
objects, each specifying the unique IPs of the corresponding OT card.

Note: Orange Tree card must be configured before use. See the README for details

Exposed methods:
    arm() - puts camera into wait state for external trigger
    readFrames() - waits for data ready register flag, then copies camera image data
    into numpy arrays
    readoff() - waits for data ready register flag, then copies camera image data
    into numpy arrays; returns payload, payload size, and error message
    sendCMD(pkt) - sends packet object via serial port
    readSerial(size, timeout) - read 'size' bytes from connection
    writeSerial(outstring) - submits string 'outstring' over connection
```

closeDevice() - close connections and free resources

getCardInfo() - prints report of details of OT card and connection

getCardIP() - returns IP address of OT card

Code to manage Gigabit Ethernet connection to board. Each GigE object manages a

Definition at line 27 of file GigE.py.

### 6.4.2 Constructor & Destructor Documentation

### 6.4.2.1 \_\_init\_\_()

```
nsCamera.comms.GigE.GigE.__init__ (
                 self,
                 camassem )
Args:
    camassem: parent cameraAssembler object
Definition at line 49 of file GigE.py.
00049
           def __init__(self, camassem):
00050
00051
               Aras:
               camassem: parent cameraAssembler object
00052
00053
00054
               self.ca = camassem
00055
               self.logcrit = self.ca.logcritbase + "[GigE] "
               self.logerr = self.ca.logerrbase + "[GigE]
00056
               self.logwarn = self.ca.logwarnbase + "[GigE] "
self.loginfo = self.ca.loginfobase + "[GigE] "
00057
00058
00059
               self.logdebug = self.ca.logdebugbase + "[GigE] "
               logging.info(self.loginfo + "Initializing GigE comms object")
00060
00061
               self.mode = 1
               self.writeTimeout = 10000
00062
               self.readTimeout = 10000
00063
00064
               self.payloadsize = (
00065
                   self.ca.sensor.width
00066
                    * self.ca.sensor.height
00067
                    * self.ca.sensor.nframes
00068
                    * self.ca.sensor.bytesperpixel
00069
00070
               logging.debug(
00071
                    self.logdebug + "Payload size: " + str(self.payloadsize) + " bytes"
00072
00073
               self.skipError = False
00074
00075
               self.ZErrorDict = {
                   0x8000: "Socket Error",
0x8001: "Internal Error",
00076
00077
                    0x8002: "Illegal Status Code",
00078
                    0x8003: "Null Parameter",
00079
08000
                    0x8004: "Out of Memory",
                   0x8005: "Invalid Connection Type",
0x8006: "Illegal Connection",
00081
00082
                    0x8007: "Socket Closed Unexpectedly",
00083
                    0x8008: "Timeout",
0x8009: "Illegal Parameter",
00084
00085
00086
00088
               if self.ca.port:
00089
                   logging.debug(
00090
                        self.logdebug + "Port supplied to GigE.py: " + str(self.ca.port)
00091
00092
                    if isinstance(self.ca.port, int) and 0 < self.ca.port < 65536:</pre>
00093
                        self.dport = self.ca.port
00094
                    else:
00095
                        logging.error(
                             self.logerr + "Invalid port number supplied, defaulting to "
00096
00097
                             "20482 "
00098
                        self.dport = 20482
00099
00100
               else:
00101
                   self.dport = 20482 # default
00102
               self.ca.port = self.dport
00103
               logging.debug(self.logdebug + "Port used by GigE.py: " + str(self.dport))
00104
00105
               logging.debug(self.logdebug + "CPU architecture: " + str(self.ca.arch))
if self.ca.arch == "64bit":
00106
00107
                   arch = "64"
00108
```

```
00109
              else:
00110
                  arch = "32"
00111
00112
              logging.debug(self.logdebug + "Operating system: " + str(self.ca.platform))
00113
              if self.ca.platform == "Windows":
00114
                  lib_name = "ZestETM1.dll"
00115
              elif self.ca.platform == "Linux" or self.ca.platform == "Darwin":
00116
                  lib_name = "libZestETM1.so"
00117
              else:
00118
                  logging.warning(
00119
                      self.logwarn + "System does not self-identify as Linux, Windows, "
00120
                      "or Mac. Assuming posix-style libraries "
00121
00122
                  lib_name = "libZestETM1.so"
00123
00124
              self.closecard = False
00125
00126
              libpath = os.path.join(self.ca.packageroot, "comms", "ZestETM1", arch, lib_name)
00127
              self._zest = C.CDLL(libpath)
00128
00129
              self.CardInfo = self.ZESTETM1_CARD_INFO()
00130
              self.CardInfoP = C.pointer(self.CardInfo)
00131
00132
              # functions
00133
              self.ZCountCards = self._zest.ZestETM1CountCards
00134
              self.ZCountCards.argtypes = [
00135
                  C.POINTER(C.c_ulong),
                  C.POINTER(C.POINTER(self.ZESTETM1_CARD_INFO)),
00136
00137
                  C.c_int,
00138
              ]
00139
00140
              self.ZOpenConnection = self._zest.ZestETM1OpenConnection
00141
              self.ZOpenConnection.argtypes = |
                  C.POINTER(self.ZESTETM1_CARD_INFO),
00142
00143
                  C.c_int,
00144
                  C.c_ushort,
00145
                  C.c ushort,
00146
                  C.POINTER(C.c_void_p),
00147
             ]
00148
00149
              self.ZWriteData = self._zest.ZestETM1WriteData
00150
              self.ZWriteData.argtypes = [
00151
                  C.c_void_p,
00152
                  C.c_void_p,
00153
                  C.c_ulong,
00154
                  C.POINTER(C.c_ulong),
00155
                  C.c_ulong,
00156
00157
00158
              self.ZReadData = self._zest.ZestETM1ReadData
00159
              self.ZReadData.argtypes = [
00160
                  C.c_void_p,
00161
                  C.c_void_p,
                  C.c_ulong,
00162
00163
                  C.POINTER(C.c_ulong),
00164
                  C.c_ulong,
00165
00166
00167
              self.Connection = C.c_void_p()
00168
              self.openDevice(self.ca.timeout)
00169
```

### 6.4.3 Member Function Documentation

### 6.4.3.1 arm()

```
nsCamera.comms.GigE.GigE.arm ( self, \\ mode \ )
```

Puts camera into wait state for trigger. Mode determines source; defaults to

```
'Hardware'
Args:
             {\tt 'Software'} \mid {\tt 'S'} \text{ activates software, disables hardware triggering}
    mode:
              'Hardware' | 'H' activates hardware, disables software triggering
                Hardware is the default
Returns:
    tuple (error, response string)
Definition at line 221 of file GigE.py.
00221
          def arm(self, mode):
00222
00223
              Puts camera into wait state for trigger. Mode determines source; defaults to
00224
                'Hardware'
00225
00226
              Args:
00227
                 mode:
                          'Software' | 'S' activates software, disables hardware triggering
00228
                          'Hardware' | 'H' activates hardware, disables software triggering
00229
                            Hardware is the default
00230
00231
              Returns:
              tuple (error, response string)
00232
00233
00234
              if not mode:
00235
                  mode = "Hardware"
00236
                  logging.info(self.loginfo + "arm")
00237
              logging.debug(self.logdebug + "arming mode: " + str(mode))
00238
              self.ca.clearStatus()
00239
              self.ca.latchPots()
00240
              err, resp = self.ca.startCapture(mode)
00241
              if err:
00242
                  logging.error(self.logerr + "unable to arm camera")
00243
              else:
00244
                  self.ca.armed = True
00245
                  self.skipError = True
00246
              return err, resp
00247
```

#### 6.4.3.2 closeDevice()

```
{\tt nsCamera.comms.GigE.GigE.closeDevice} \ \ ( {\tt self} \ )
```

Close connection to Orange Tree card and free resources

# Definition at line 467 of file GigE.py.

```
00467
          def closeDevice(self):
00468
              Close connection to Orange Tree card and free resources
00469
00470
              logging.debug(self.logdebug + "Closing connection to Orange Tree card")
00471
00472
              self._zest.ZestETM1CloseConnection(self.Connection)
00473
              if self.closecard:
00474
                  try:
00475
                      self._zest.ZestETM1FreeCards(self.CardInfoP)
00476
                  except SystemExit:
00477
                      raise
00478
                  except KeyboardInterrupt:
00479
                      rais
00480
                  except Exception:
                      logging.error(self.logerr + "Error reported in OT card closure")
00481
00482
              self._zest.ZestETM1Close()
00483
```

### 6.4.3.3 getCardInfo()

```
nsCamera.comms.GigE.GigE.getCardInfo (
                    self )
Prints status message with information returned by OT card
Definition at line 493 of file GigE.py.
            def getCardInfo(self):
00495
                 Prints status message with information returned by OT card
00496
                 ci = self.CardInfoP.contents
00497
                 print("GigE Card Status:")
00498
                 print("----")
00499
                 print("IP: " + ".".join(str(e) for e in [b for b in ci.IPAddr]))
00500
                 print("ControlPort: " + str(ci.ControlPort))
00501
                 print("Timeout: " + str(ci.Timeout))
00502
                 print("Timeout: " + str(ci.Timeout))
print("HTTPPort: " + str(ci.HTTPPort))
print("MACAddr: " + ".".join(format(e, "02X") for e in [b for b in ci.MACAddr]))
print("SubNet: " + ".".join(str(e) for e in [b for b in ci.SubNet]))
print("Gateway: " + ".".join(str(e) for e in [b for b in ci.Gateway]))
print("SerialNumber: " + str(ci.SerialNumber))
print("FirmwareVersion: " + str(ci.FirmwareVersion))
print("HardwareVersion: " + str(ci.HardwareVersion))
print("------")
00503
00504
00505
00506
00507
00508
00509
                 print("----")
00510
00511
6.4.3.4 getCardIP()
nsCamera.comms.GigE.GigE.getCardIP (
                   self )
Query IP address of OT card
Returns: address of OT card as list of bytes
Definition at line 484 of file GigE.py.
            def getCardIP(self):
00485
00486
                 Query IP address of OT card
00487
00488
                 Returns: address of OT card as list of bytes
00489
00490
                 return self.CardInfo.IPAddr
00491
6.4.3.5 openDevice()
nsCamera.comms.GigE.GigE.openDevice (
                    self,
                    timeout = 30)
Find Orange Tree card and open a connection; if IP is supplied as parameter for
   the CameraAssembler, bypass network search and connect directly to indicated
   IP address
```

timeout: timeout in seconds for attempting to connect to a card

Args:

Definition at line 402 of file GigE.py.

```
00402
          def openDevice(self, timeout=30):
00403
00404
              Find Orange Tree card and open a connection; if IP is supplied as parameter for
00405
                the CameraAssembler, bypass network search and connect directly to indicated
00406
00407
00408
              timeout: timeout in seconds for attempting to connect to a card \tt^{\tt nnn}
00409
00410
00411
              err = self._zest.ZestETM1Init()
00412
              if err:
00413
                  logging.critical(self.logcrit + "ZestETM1Init failure")
00414
                  sys.exit(1)
00415
              logging.info(self.loginfo + "searching for Orange Tree cards")
00416
              NumCards = C.c_ulong(0)
00417
00418
              if self.ca.iplist:
00419
                  ubyte4 = C.c_ubyte * 4
                  self.CardInfo.IPAddr = ubyte4(*self.ca.iplist)
00420
                  self.CardInfo.ControlPort = C.c_ushort(self.dport)
self.CardInfo.Timeout = C.c_ulong(self.writeTimeout)
00421
00422
00423
                  self.closecard = False
00424
              else:
00425
                  wait = 0
                  while True:
00426
00427
                      logging.debug(
00428
                           self.logdebug + "openDevice: connection wait = " + str(wait)
00429
                       err = self.ZCountCards(C.byref(NumCards), C.byref(self.CardInfoP), 1000)
00430
00431
                       if err:
00432
                           logging.critical(self.logcrit + "CountCards failure")
00433
                           svs.exit(1)
00434
                       if NumCards.value > 0:
00435
                           break
                       if wait == timeout:
00436
00437
                           logging.critical(
00438
                                "{}No Orange Tree cards found in {} seconds".format(
00439
                                   self.logcrit, timeout
00440
00441
                           )
00442
                           sys.exit(1)
                       elif not wait % 5:
00443
00444
                           logging.info(
00445
                               "{}Still trying to connect after {} seconds...".format(
00446
                                    self.loginfo, wait
00447
                               )
00448
                          )
00449
                       wait += 1
00450
                  logging.info(
                       self.loginfo + str(NumCards.value) + " Orange Tree card(s) found"
00451
00452
                  ) # TODO: add check for GigE bit in board description
00453
              err = self.ZOpenConnection(
00454
                  self.CardInfoP, 0, self.dport, 0, C.byref(self.Connection)
00455
00456
              if err:
00457
                  if err == 0x4000:
00458
                      logging.warning(
00459
                           self.logerr + "OT Card emitted an undefined warning message"
00460
00461
                  else:
00462
                       logging.critical(
00463
                           self.logcrit + "OpenConnection failure: " + self.ZErrorDict[err]
00464
00465
                  sys.exit(1)
00466
```

#### 6.4.3.6 readFrames()

```
Copies image data from board into numpy arrays.
Aras:
    waitOnSRAM: if True, wait until SRAM_READY flag is asserted to begin copying
    timeout: passed to waitForSRAM; after this many seconds begin copying data
      irrespective of SRAM_READY status; 'zero' means wait indefinitely
      WARNING: If acquisition fails, the SRAM will not contain a current image,
        but the code will copy the data anyway
    fast: if False, parse and convert frames to numpy arrays; if True, return
      unprocessed text stream
    columns: 1 for single image per frame, 2 for separate hemisphere images
Returns:
    list of numpy arrays OR raw text stream
Definition at line 248 of file GigE.py.
         def readFrames(self, waitOnSRAM, timeout=0, fast=False, columns=1):
00248
00249
00250
             Copies image data from board into numpy arrays.
00251
00252
             Aras:
00253
                 waitOnSRAM: if True, wait until SRAM_READY flag is asserted to begin copying
00254
                  data
                 timeout: passed to waitForSRAM; after this many seconds begin copying data irrespective of SRAM_READY status; 'zero' means wait indefinitely
00255
00256
00257
                   WARNING: If acquisition fails, the SRAM will not contain a current image,
00258
                     but the code will copy the data anyway
00259
                 fast: if False, parse and convert frames to numpy arrays; if True, return
00260
                   unprocessed text stream
                 columns: 1 for single image per frame, 2 for separate hemisphere images
00261
00262
00263
             Returns:
00264
                 list of numpy arrays OR raw text stream
00265
00266
00267
             frames, _, _ = self.readoff(waitOnSRAM, timeout, fast, columns)
00268
             return frames
00269
6.4.3.7 readoff()
nsCamera.comms.GigE.GigE.readoff (
               self.
               waitOnSRAM,
               timeout = 0,
               fast = False,
               columns = 1)
Copies image data from board into numpy arrays; returns data, length of data,
and error messages. Use 'readFrames()' unless you require this additional
information
Args:
    waitOnSRAM: if True, wait until SRAM_READY flag is asserted to begin copying
    timeout: passed to waitForSRAM; after this many seconds begin copying data
      irrespective of SRAM_READY status; 'zero' means wait indefinitely
      WARNING: If acquisition fails, the SRAM will not contain a current image,
        but the code will copy the data anyway
```

fast: if False, parse and convert frames to numpy arrays; if True, return

columns: 1 for single image per frame, 2 for separate hemisphere images

unprocessed text stream

```
Returns:

tuple (list of numpy arrays OR raw text stream, length of downloaded payload

in bytes, payload error flag) since CRC check is handled by TCP/IP,

payload error flag is always False for GigE
```

Definition at line 270 of file GigE.pv. 00270 def readoff(self, waitOnSRAM, timeout=0, fast=False, columns=1): 00271 00272 Copies image data from board into numpy arrays; returns data, length of data, 00273 and error messages. Use 'readFrames()' unless you require this additional 00274 information 00275 00276 Args: 00277 waitOnSRAM: if True, wait until SRAM\_READY flag is asserted to begin copying 00278 data timeout: passed to waitForSRAM; after this many seconds begin copying data
irrespective of SRAM\_READY status; 'zero' means wait indefinitely 00279 00280 WARNING: If acquisition fails, the SRAM will not contain a current image, 00281 00282 but the code will copy the data anyway 00283 fast: if False, parse and convert frames to numpy arrays; if True, return 00284 unprocessed text stream 00285 columns: 1 for single image per frame, 2 for separate hemisphere images 00286 00287 Returns: 00288 tuple (list of numpy arrays OR raw text stream, length of downloaded payload 00289 in bytes, payload error flag) since CRC check is handled by TCP/IP, 00290 payload error flag is always False for GigE 00291 00292 logging.info(self.loginfo + "readoff") 00293 logging.debug( 00294 self.logdebug 00295 + "readoff: waitonSRAM = " 00296 + str(waitOnSRAM) 00297 + "; timeout = 00298 + str(timeout) 00299 + "; fast = 00300 + str(fast) 00301 00302 # Wait for data to be ready on board 00303 # Skip wait only if explicitly tagged 'False' ('None' defaults to True) 00304 if waitOnSRAM is not False: 00305 self.ca.waitForSRAM(timeout) 00306 self.skipError = False 00307 err, rval = self.ca.readSRAM() 00308 if err: 00309 logging.error(self.logerr + "Error detected in readSRAM") elif self.ca.boardname == "llnl\_v4": 00310 00311 # self.ca.setSubregister('SWACK','1') 00312 00313 # extract the data. Remove header; the FPGA returns a packet without the CRC 00314 # suffix 00315 logging.debug(self.logdebug + "readoff: first 64 chars: " + str(rval[0:64])) 00316 data = rval[32:1]00317 if fast: 00318 return data, len(data) // 2, bool(err) 00319 else: parsed = generateFrames(self.ca, data, columns) 00320 00321 return parsed, len(data) // 2, bool(err)

#### 6.4.3.8 readSerial()

00322

```
size: number of bytes to read
   timeout: serial timeout in sec (defaults to self.readTimeout)
Returns:
   tuple (error string, string read from serial port)
Definition at line 363 of file GigE.py.
          def readSerial(self, size, timeout=None):
00364
00365
              Read bytes from the serial port. Does not verify packets.
00366
00367
              Args:
00368
                size: number of bytes to read
00369
                 timeout: serial timeout in sec (defaults to self.readTimeout)
00370
00371
              tuple (error string, string read from serial port)
00372
00373
00374
              logging.debug(
00375
                  self.logdebug
00376
                  + "readSerial: size = "
00377
                  + str(size)
00378
                 + "; timeout = "
00379
                  + str(timeout)
00380
00381
              if not timeout:
00382
                  timeout = self.readTimeout
00383
              inbuff = C.create_string_buffer(size + 1)
00384
              inbuffp = C.pointer(inbuff)
              readlen = C.c_ulong(0)
00385
00386
              err = self.ZReadData(self.Connection, inbuffp, size, C.byref(readlen), timeout)
00387
              if err:
00388
                  if self.skipError:
00389
                     logging.debug(
                          self.logdebug + "readSerial: skipped error: " + self.ZErrorDict[err]
00390
00391
00392
                      self.skipError = False
00393
                  elif err == 0x4000:
00394
                      logging.warning(
                          self.logerr + "OT Card emitted an undefined warning message"
00395
00396
00397
                  else:
                      logging.error(self.logerr + "readSerial error: " + self.ZErrorDict[err])
00398
00399
              return bytes2str(inbuff.raw)[:-2]
00400
```

# 6.4.3.9 sendCMD()

```
nsCamera.comms.GigE.GigE.sendCMD (
               self.
               pkt )
Submit packet and verify the response packet.
Packet communications with FPGA omit CRC suffix, so adds fake CRC bytes to
  response
Aras:
    pkt: Packet object
Returns:
    tuple (error, response string)
Definition at line 170 of file GigE.py.
         def sendCMD(self, pkt):
00170
00171
00172
             Submit packet and verify the response packet.
00173
             Packet communications with FPGA omit CRC suffix, so adds fake CRC bytes to
```

```
00174
                 response
00175
00176
               Args:
                  pkt: Packet object
00177
00178
00179
               tuple (error, response string)
00180
00181
00182
00183
               pktStr = pkt.pktStr()[0:16]
00184
               logging.debug(self.logdebug + "sendCMD packet: " + str(pktStr))
00185
00186
               self.ca.writeSerial(pktStr)
00187
               if (
                   hasattr(self.ca, "board")
and pktStr[4] == "0"
00188
00189
00190
                   and pktStr[5:8] == self.ca.board.registers["SRAM_CTL"]
00191
00192
                   bufsize = self.payloadsize + 16
                   resptext = self.readSerial(bufsize)
00193
00194
00195
                   if len(resptext) < 32:</pre>
00196
                       logging.debug(self.logdebug + "sendCMD resptext = " + str(resptext))
00197
                   else:
00198
                       logging.debug(
00199
                            self.logdebug
                            + "sendCMD resptext (truncated) = "
00200
                            + str(resptext)[0:32]
00201
00202
                       )
00203
00204
                   if len(resptext) < bufsize + 16:</pre>
00205
                       err += (
                            self.logerr + "sendCMD- packet too small, payload may be incomplete"
00206
00207
00208
                       logging.error(err)
00209
               else:
00210
                   # workaround for initial setup before board object has been initialized
00211
                   resp = self.readSerial(8)
                   logging.debug(self.logdebug + "sendCMD response: " + str(resp))
00212
00213
                   if len(resp) < 8:</pre>
                       err += self.logerr + "sendCMD- response too small, returning zeros" resptext = "000000000000000000000"
00214
00215
00216
                       logging.error(err)
00217
                   else:
00218
                       resptext = resp + "0000"
00219
               return err, resptext
00220
```

#### 6.4.3.10 writeSerial()

```
nsCamera.comms.GigE.GigE.writeSerial (
               self.
               outstring,
               timeout )
Transmit string to board
Args:
    outstring: string to write
    timeout: serial timeout in sec (defaults to self.writeTimeout)
Returns:
    integer number of bytes written
Definition at line 323 of file GigE.py.
00323
         def writeSerial(self, outstring, timeout):
00324
00325
             Transmit string to board
00326
00327
             Args:
```

```
00328
                 outstring: string to write
00329
                 timeout: serial timeout in sec (defaults to self.writeTimeout)
00330
00331
             Returns:
             integer number of bytes written
00332
00333
00334
             logging.debug(
00335
                self.logdebug
00336
                 + "writeSerial: outstring = "
00337
                 + str(outstring)
00338
                 + "; timeout =
00339
                 + str(timeout)
00340
00341
            if not timeout:
00342
                 timeout = self.writeTimeout
00343
            outstring = str2bytes(outstring)
00344
             outbuff = C.create_string_buffer(outstring)
            outbuffp = C.pointer(outbuff)
00345
             outbufflen = len(outstring)
00346
00347
             writelen = C.c_ulong(0)
             err = self.ZWriteData(
00348
00349
                 self.Connection, outbuffp, outbufflen, C.byref(writelen), timeout
00350
00351
             if err:
00352
                if err == 0x4000:
00353
                     logging.warning(
                         self.logerr + "OT Card emitted an undefined warning message"
00354
00355
00356
                 else:
00357
                     logging.error(
00358
                         self.logerr + "writeSerial error: " + self.ZErrorDict[err]
00359
             logging.debug(self.logdebug + "writeSerial: writelen = " + str(writelen))
00360
0.0361
             return writelen
00362
```

# 6.4.4 Member Data Documentation

# 6.4.4.1 \_zest

```
nsCamera.comms.GigE.GigE._zest [protected]
```

Definition at line 127 of file GigE.py.

### 6.4.4.2 ca

```
nsCamera.comms.GigE.GigE.ca
```

Definition at line 54 of file GigE.py.

#### 6.4.4.3 CardInfo

```
\verb|nsCamera.comms.GigE.GigE.CardInfo|\\
```

Definition at line 129 of file GigE.py.

# 6.4.4.4 CardInfoP

nsCamera.comms.GigE.GigE.CardInfoP

Definition at line 130 of file GigE.py.

#### 6.4.4.5 closecard

nsCamera.comms.GigE.GigE.closecard

Definition at line 124 of file GigE.py.

# 6.4.4.6 Connection

nsCamera.comms.GigE.GigE.Connection

Definition at line 167 of file GigE.py.

# 6.4.4.7 dport

nsCamera.comms.GigE.GigE.dport

Definition at line 93 of file GigE.py.

### 6.4.4.8 logcrit

nsCamera.comms.GigE.GigE.logcrit

Definition at line 55 of file GigE.py.

# 6.4.4.9 logdebug

nsCamera.comms.GigE.GigE.logdebug

Definition at line 59 of file GigE.py.

# 6.4.4.10 logerr

nsCamera.comms.GigE.GigE.logerr

Definition at line 56 of file GigE.py.

# 6.4.4.11 loginfo

nsCamera.comms.GigE.GigE.loginfo

Definition at line 58 of file GigE.py.

# 6.4.4.12 logwarn

nsCamera.comms.GigE.GigE.logwarn

Definition at line 57 of file GigE.py.

# 6.4.4.13 mode

nsCamera.comms.GigE.GigE.mode

Definition at line 61 of file GigE.py.

# 6.4.4.14 payloadsize

nsCamera.comms.GigE.GigE.payloadsize

Definition at line 64 of file GigE.py.

#### 6.4.4.15 readTimeout

 ${\tt nsCamera.comms.GigE.GigE.readTimeout}$ 

Definition at line 63 of file GigE.py.

# 6.4.4.16 skipError

nsCamera.comms.GigE.GigE.skipError

Definition at line 73 of file GigE.py.

### 6.4.4.17 writeTimeout

nsCamera.comms.GigE.GigE.writeTimeout

Definition at line 62 of file GigE.py.

# 6.4.4.18 ZCountCards

```
\verb|nsCamera.comms.GigE.GigE.ZCountCards| \\
```

Definition at line 133 of file GigE.py.

#### 6.4.4.19 ZErrorDict

```
nsCamera.comms.GigE.GigE.ZErrorDict
```

Definition at line 75 of file GigE.py.

# 6.4.4.20 ZESTETM1\_CARD\_INFO

```
nsCamera.comms.GigE.GigE.ZESTETM1_CARD_INFO
```

Definition at line 136 of file GigE.py.

# 6.4.4.21 ZOpenConnection

```
nsCamera.comms.GigE.GigE.ZOpenConnection
```

Definition at line 140 of file GigE.py.

#### 6.4.4.22 ZReadData

 ${\tt nsCamera.comms.GigE.GigE.ZReadData}$ 

Definition at line 158 of file GigE.py.

# 6.4.4.23 ZWriteData

nsCamera.comms.GigE.GigE.ZWriteData


Definition at line 149 of file GigE.py.

The documentation for this class was generated from the following file:

• C:/Users/hill35/git/camera\_python/nsCamera/comms/GigE.py

# 6.5 nsCamera.sensors.icarus.icarus Class Reference

Inheritance diagram for nsCamera.sensors.icarus.icarus:



### **Public Member Functions**

- \_\_init\_\_ (self, ca)
- · checkSensorVoltStat (self)
- sensorSpecific (self)

### Public Member Functions inherited from nsCamera.sensors.sensorBase.sensorBase

- init\_board\_specific (self)
- setInterlacing (self, ifactor)
- setHighFullWell (self, flag)
- setZeroDeadTime (self, flag)
- setTriggerDelay (self, delay)
- setPhiDelay (self, delay)
- setExtClk (self, delay)
- setTiming (self, side="AB", sequence=None, delay=0)
- setArbTiming (self, side="AB", sequence=None)
- getTiming (self, side, actual)
- setManualShutters (self, timing=None)
- setManualTiming (self, timing=None)
- getManualTiming (self)
- getSensTemp (self, scale=None, offset=None, slope=None, dec=None)
- selectOscillator (self, osc=None)
- parseReadoff (self, frames, columns)
- getSensorStatus (self)
- reportStatusSensor (self, statusbits, statusbits2)

#### **Public Attributes**

- ca
- sens\_registers
- sens\_subregisters

# Public Attributes inherited from nsCamera.sensors.sensorBase.sensorBase

- ca
- logcrit
- logerr
- logwarn
- loginfo
- logdebug
- fpganumID
- sensfam

#### **Static Public Attributes**

```
• str specwarn = " and the use of the lcarus model 1 sensor"
```

- int minframe = 1
- int maxframe = 2
- int maxwidth = 512
- int maxheight = 1024
- int bytesperpixel = 2
- int icarustype = 1
- int fpganumID = 1
- str detect = "ICARUS DET"
- str sensfam = "lcarus"
- str loglabel = "[lcarus1] "
- int firstframe = 1
- int lastframe = 2
- int nframes = 2
- int width = 512
- int height = 1024
- int firstrow = 0
- int lastrow = 1023
- list interlacing = [0, 0]
- int columns = 1
- bool padToFull = True

# 6.5.1 Detailed Description

Definition at line 27 of file icarus.py.

### 6.5.2 Constructor & Destructor Documentation

#### 6.5.2.1 init ()

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

```
Definition at line 53 of file icarus.py.
```

```
00053
              def __init__(self, ca):
00054
                    self.ca = ca
00055
                    super(icarus, self).__init__(ca)
00056
00057
                    self.sens registers = OrderedDict(
00058
00059
                                "VRESET_WAIT_TIME": "03E",
                                "ICARUS_VER_SEL": "041",
"VRESET_HIGH_VALUE": "04A"
00060
00061
                                "MISC_SENSOR_CTL": "04C"
00062
                                "MANUAL_SHUTTERS_MODE": "050",
00063
                                "W0_INTEGRATION": "051",
"W0_INTERFRAME": "052",
00064
00065
                                "W1_INTEGRATION": "053",
"W1_INTERFRAME": "054",
00066
00067
                                "W2_INTEGRATION": "055",
00068
                                "W2_INTERFRAME": "056",
00069
                                "W3_INTEGRATION": "057"
00070
00071
                                "WO_INTEGRATION_B": "058",
                                "WO_INTERFRAME_B": "059",
00072
                                "W1_INTEGRATION_B": "05A",
"W1_INTERFRAME_B": "05B",
"W2_INTEGRATION_B": "05C",
00073
00074
00075
                                "W2_INTERFRAME_B": "05D",
"W3_INTEGRATION_B": "05E",
00076
00077
00078
                                "TIME_ROW_DCD": "05F",
00079
                          }
00080
                   )
00081
00082
                    self.sens_subregisters = [
00083
00085
                          ("MANSHUT_MODE", "MANUAL_SHUTTERS_MODE", 0, 1, True),
00086
                           ("REVREAD", "CTRL_REG", 4, 1, True),
                          ("PDBIAS_LOW", "CTRL_REG", 6, 1, True), ("ROWDCD_CTL", "CTRL_REG", 7, 1, True),
00087
00088
                          ("ROWDED_CTL", "CTRL_REG", 7, 1, True),

("ACCUMULATION_CTL", "MISC_SENSOR_CTL", 0, 1, True),

("HST_TST_ANRST_EN", "MISC_SENSOR_CTL", 1, 1, True),

("HST_TST_BNRST_EN", "MISC_SENSOR_CTL", 2, 1, True),

("HST_TST_ANRST_IN", "MISC_SENSOR_CTL", 3, 1, True),

("HST_TST_BNRST_IN", "MISC_SENSOR_CTL", 4, 1, True),

("HST_PXL_RST_EN", "MISC_SENSOR_CTL", 5, 1, True),

("HST_CONT_MODE", "MISC_SENSOR_CTL", 6, 1, True),

("COLL DCD_EN", "MISC_SENSOR_CTL", 7, 1, True),
00089
00090
00091
00092
00093
00094
00095
00096
                          ("COL_DCD_EN", "MISC_SENSOR_CTL", 7, 1, True),
00097
                          ("COL_READOUT_EN", "MISC_SENSOR_CTL", 8, 1, True),
00098
                          ("READOFF_DELAY_EN", "TRIGGER_CTL", 4, 1, True),
00099
                          ("STAT_W3TOPAEDGE1", "STAT_REG", 3, 1, False), ("STAT_W3TOPBEDGE1", "STAT_REG", 4, 1, False),
00102
00103
                          ("STAT_HST_ALL_W_EN_DETECTED", "STAT_REG", 12, 1, False), ("PDBIAS_UNREADY", "STAT_REG2", 5, 1, False),
00104
00105
00106
                    ]
00107
00108
                    if self.ca.boardname == "llnl_v1":
00109
                          self.sens_subregisters.append(
                                ("VRESET_HIGH", "VRESET_HIGH_VALUE", 7, 8, True)
00110
00111
                         )
00112
                    else:
00113
                         self.sens_subregisters.extend(
00114
                                      ("VRESET_HIGH", "VRESET_HIGH_VALUE", 15, 16, True),
00115
                                      ("READOFF_DELAY_EN", "TRIGGER_CTL", 4, 1, True),
00116
00117
00118
00119
                          self.sens registers.update({"DELAY ASSERTION ROWDCD EN": "04F"})
00120
```

### 6.5.3 Member Function Documentation

#### 6.5.3.1 checkSensorVoltStat()

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

```
Definition at line 121 of file icarus.py.
```

```
00121
         def checkSensorVoltStat(self):
00123
              Checks register tied to sensor select jumpers to confirm match with sensor
00124
                object
00125
00126
              Returns:
              boolean, True if jumpers select for Icarus sensor
00127
00128
              logging.debug(self.logdebug + "checkSensorVoltStat")
00129
00130
              err, status = self.ca.getSubregister("ICARUS_DET")
00131
              if err:
00132
                  logging.error(self.logerr + "unable to confirm sensor status")
00133
                  return False
00134
              if not int(status):
                  logging.error(self.logerr + "Icarus sensor not detected")
00135
00136
                  return False
00137
              return True
00138
```

#### 6.5.3.2 sensorSpecific()

```
nsCamera.sensors.icarus.icarus.sensorSpecific ( self \ ) Returns: list of tuples, (Sensor-specific register, default setting)
```

### Definition at line 139 of file icarus.py.

```
def sensorSpecific(self):
00140
00141
                     list of tuples, (Sensor-specific register, default setting)
00142
00143
00144
                     icarussettings = [
00145
                          ("ICARUS_VER_SEL", "00000001"),
                           ("IFPA_ROW_FINAL", "00000001"),

("FPA_ROME_INITIAL", "00000002"),

("FPA_ROW_INITIAL", "00000000"),

("FPA_ROW_FINAL", "000003FF"),
00146
00147
00148
00149
                           ("VRESET_WAIT_TIME", "000003FF ),
("NS_TIMING_DATA_BHI", "00000000"),
("HS_TIMING_DATA_BLO", "00006666"),
("HS_TIMING_DATA_AHI", "00000000"),
00150
00151
                                                                                    \# 0db6 = 2-1; 6666 = 2-2
00152
00153
                            ("HS_TIMING_DATA_ALO", "00006666"),
00154
00155
00156
                     if self.ca.boardname == "llnl_v1":
```

```
00157
                 icarussettings.append(
00158
                     ("VRESET_HIGH_VALUE", "000000D5") # 3.3 V (FF = 3.96)
00159
00160
00161
                 icarussettings.append(("VRESET_HIGH_VALUE", "0000FFFF"))
00162
              return icarussettings
00163
00164
00165 """
00166 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00169 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00170 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00171 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00172 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00173 be made under this license.
00174 """
```

### 6.5.4 Member Data Documentation

### 6.5.4.1 bytesperpixel

```
int nsCamera.sensors.icarus.icarus.bytesperpixel = 2 [static]
```

Definition at line 36 of file icarus.py.

#### 6.5.4.2 ca

```
nsCamera.sensors.icarus.icarus.ca
```

Definition at line 54 of file icarus.py.

### 6.5.4.3 columns

```
int nsCamera.sensors.icarus.icarus.columns = 1 [static]
```

Definition at line 50 of file icarus.py.

# 6.5.4.4 detect

```
str nsCamera.sensors.icarus.icarus.detect = "ICARUS_DET" [static]
```

Definition at line 39 of file icarus.py.

#### 6.5.4.5 firstframe

```
int nsCamera.sensors.icarus.icarus.firstframe = 1 [static]
```

Definition at line 42 of file icarus.py.

# 6.5.4.6 firstrow

```
int nsCamera.sensors.icarus.icarus.firstrow = 0 [static]
```

Definition at line 47 of file icarus.py.

# 6.5.4.7 fpganumID

```
int nsCamera.sensors.icarus.icarus.fpganumID = 1 [static]
```

Definition at line 38 of file icarus.py.

# 6.5.4.8 height

```
int nsCamera.sensors.icarus.icarus.height = 1024 [static]
```

Definition at line 46 of file icarus.py.

# 6.5.4.9 icarustype

```
int nsCamera.sensors.icarus.icarus.icarustype = 1 [static]
```

Definition at line 37 of file icarus.py.

### 6.5.4.10 interlacing

```
list nsCamera.sensors.icarus.icarus.interlacing = [0, 0] [static]
```

Definition at line 49 of file icarus.py.

# 6.5.4.11 lastframe

```
int nsCamera.sensors.icarus.icarus.lastframe = 2 [static]
```

Definition at line 43 of file icarus.py.

### 6.5.4.12 lastrow

```
int nsCamera.sensors.icarus.icarus.lastrow = 1023 [static]
```

Definition at line 48 of file icarus.py.

# 6.5.4.13 loglabel

```
str nsCamera.sensors.icarus.icarus.loglabel = "[Icarus1] " [static]
```

Definition at line 41 of file icarus.py.

#### 6.5.4.14 maxframe

```
int nsCamera.sensors.icarus.icarus.maxframe = 2 [static]
```

Definition at line 30 of file icarus.py.

### 6.5.4.15 maxheight

```
int nsCamera.sensors.icarus.icarus.maxheight = 1024 [static]
```

Definition at line 35 of file icarus.py.

#### 6.5.4.16 maxwidth

```
int nsCamera.sensors.icarus.icarus.maxwidth = 512 [static]
```

Definition at line 34 of file icarus.py.

#### 6.5.4.17 minframe

```
int nsCamera.sensors.icarus.icarus.minframe = 1 [static]
```

Definition at line 29 of file icarus.py.

### 6.5.4.18 nframes

```
int nsCamera.sensors.icarus.icarus.nframes = 2 [static]
```

Definition at line 44 of file icarus.py.

# 6.5.4.19 padToFull

```
bool nsCamera.sensors.icarus.icarus.padToFull = True [static]
```

Definition at line 51 of file icarus.py.

### 6.5.4.20 sens\_registers

```
nsCamera.sensors.icarus.icarus.sens_registers
```

Definition at line 57 of file icarus.py.

### 6.5.4.21 sens\_subregisters

```
nsCamera.sensors.icarus.icarus.sens_subregisters
```

Definition at line 82 of file icarus.py.

#### 6.5.4.22 sensfam

```
str nsCamera.sensors.icarus.icarus.sensfam = "Icarus" [static]
```

Definition at line 40 of file icarus.py.

### 6.5.4.23 specwarn

```
str nsCamera.sensors.icarus.icarus.specwarn = " and the use of the Icarus model 1 sensor" [static]
```

Definition at line 28 of file icarus.py.

### 6.5.4.24 width

```
int nsCamera.sensors.icarus.icarus.width = 512 [static]
```


Definition at line 45 of file icarus.py.

The documentation for this class was generated from the following file:

C:/Users/hill35/git/camera\_python/nsCamera/sensors/icarus.py

# 6.6 nsCamera.sensors.icarus2.icarus2 Class Reference

Inheritance diagram for nsCamera.sensors.icarus2.icarus2:



### **Public Member Functions**

- init (self, ca)
- sensorSpecific (self)

# Public Member Functions inherited from nsCamera.sensors.sensorBase.sensorBase

- init board specific (self)
- checkSensorVoltStat (self)
- setInterlacing (self, ifactor)
- setHighFullWell (self, flag)
- setZeroDeadTime (self, flag)
- setTriggerDelay (self, delay)
- setPhiDelay (self, delay)
- setExtClk (self, delay)
- setTiming (self, side="AB", sequence=None, delay=0)
- setArbTiming (self, side="AB", sequence=None)
- getTiming (self, side, actual)
- setManualShutters (self, timing=None)
- setManualTiming (self, timing=None)
- getManualTiming (self)
- getSensTemp (self, scale=None, offset=None, slope=None, dec=None)
- selectOscillator (self, osc=None)
- parseReadoff (self, frames, columns)
- getSensorStatus (self)
- reportStatusSensor (self, statusbits, statusbits2)

### **Public Attributes**

- ca
- sens\_registers
- sens\_subregisters

### Public Attributes inherited from nsCamera.sensors.sensorBase.sensorBase

- ca
- logcrit
- logerr
- logwarn
- loginfo
- logdebug
- fpganumID
- sensfam

### **Static Public Attributes**

```
str specwarn = ""
int minframe = 0
int maxframe = 3
int maxwidth = 512
int maxheight = 1024
int bytesperpixel = 2
int icarustype = 0
int fpganumID = 1
str detect = "ICARUS_DET"
str sensfam = "Icarus"
str loglabel = "[Icarus2] "
int firstframe = 0
int lastframe = 3
int nframes = 4
int width = 512
```

- int height = 1024int firstrow = 0
- int lastrow = 1023
- list interlacing = [0, 0]
- int columns = 1
- bool padToFull = True

# 6.6.1 Detailed Description

Definition at line 24 of file icarus2.py.

### 6.6.2 Constructor & Destructor Documentation

# 6.6.2.1 init ()

Reimplemented from nsCamera.sensors.sensorBase.sensorBase.

#### Definition at line 47 of file icarus2.py.

```
def __init__(self, ca):
00047
00048
                self.ca = ca
00049
                super(icarus2, self).__init__(ca)
00050
00051
                self.sens_registers = OrderedDict(
00052
00053
                          "VRESET_WAIT_TIME": "03E",
                          "ICARUS_VER_SEL": "041",
"MISC_SENSOR_CTL": "04C",
00054
00055
                          "MANUAL_SHUTTERS_MODE": "050",
00056
                          "WO_INTEGRATION": "051",
00057
                          "WO_INTERFRAME": "052",
"W1_INTEGRATION": "053",
00058
00059
                          "W1_INTERFRAME": "054",
00060
```

```
00061
                                      "W2_INTEGRATION": "055",
00062
                                      "W2_INTERFRAME": "056",
00063
                                      "W3_INTEGRATION": "057"
                                      "W0_INTEGRATION_B": "058",
00064
00065
                                      "WO_INTERFRAME_B": "059",
                                      "W1_INTEGRATION_B": "05A",
00066
00067
                                      "W1_INTERFRAME_B": "05B",
00068
                                      "W2_INTEGRATION_B": "05C",
00069
                                      "W2_INTERFRAME_B": "05D",
00070
                                      "W3_INTEGRATION_B": "05E",
00071
                                      "TIME_ROW_DCD": "05F",
00072
00073
                       )
00074
00075
                       self.sens_subregisters = [
00076
00078
                                ("MANSHUT_MODE", "MANUAL_SHUTTERS_MODE", 0, 1, True),
                               ("REVREAD", "CTRL_REG", 4, 1, True), ("PDBIAS_LOW", "CTRL_REG", 6, 1, True), ("ROWDCD_CTL", "CTRL_REG", 7, 1, True),
00079
08000
00081
                              ("ROWDCD_CTL", "CTRL_REG", 7, 1, True),

("ACCUMULATION_CTL", "MISC_SENSOR_CTL", 0, 1, True),

("HST_TST_ANRST_EN", "MISC_SENSOR_CTL", 1, 1, True),

("HST_TST_BNRST_EN", "MISC_SENSOR_CTL", 2, 1, True),

("HST_TST_ANRST_IN", "MISC_SENSOR_CTL", 3, 1, True),

("HST_TST_BNRST_IN", "MISC_SENSOR_CTL", 4, 1, True),

("HST_PXL_RST_EN", "MISC_SENSOR_CTL", 5, 1, True),

("HST_CONT_MODE", "MISC_SENSOR_CTL", 6, 1, True),

("COL_DCD_EN", "MISC_SENSOR_CTL", 7, 1, True),
00082
00083
00084
00085
00086
00087
00088
00089
                               ("COL_READOUT_EN", "MISC_SENSOR_CTL", 8, 1, True),
00090
00091
                               ("STAT_W3TOPAEDGE1", "STAT_REG", 3, 1, False),
("STAT_W3TOPBEDGE1", "STAT_REG", 4, 1, False),
("STAT_HST_ALL_W_EN_DETECTED", "STAT_REG", 12, 1, False),
("PDBIAS_UNREADY", "STAT_REG2", 5, 1, False),
00094
00095
00096
00097
00098
                      1
00099
00100
                       if self.ca.boardname == "llnl_v4":
                               self.sens_subregisters.append(
    ("READOFF_DELAY_EN", "TRIGGER_CTL", 4, 1, True)
00101
00102
00103
00104
                               self.sens_registers.update({"DELAY_ASSERTION_ROWDCD_EN": "04F"})
00105
```

#### 6.6.3 Member Function Documentation

# 6.6.3.1 sensorSpecific()

```
nsCamera.sensors.icarus2.icarus2.sensorSpecific (
                    self )
Returns:
      list of tuples, (Sensor-specific register, default setting)
Definition at line 107 of file icarus2.py.
            def sensorSpecific(self):
00108
00109
                  Returns:
                  list of tuples, (Sensor-specific register, default setting)
00110
00111
00112
                  return [
                       ("ICARUS_VER_SEL", "00000000"),
00113
                       ("FPA_FRAME_INITIAL", "00000000"),
("FPA_FRAME_FINAL", "00000003"),
("FPA_ROW_INITIAL", "00000000"),
("FPA_ROW_FINAL", "000003FF"),
00114
00115
00116
00117
                       ("HS_TIMING_DATA_BHI", "00000000"),
("HS_TIMING_DATA_BLO", "00006666"),
00118
                                                                        \# 0db6 = 2-1; 6666 = 2-2
00119
                       ("HS_TIMING_DATA_AHI", "00000000"),
("HS_TIMING_DATA_ALO", "00006666"),
00120
00121
00122
                  1
```

```
00123
00124
00125 """
00126 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00127 LLNL-CODE-838080
00128
00129 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00130 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00131 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00132 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00133 be made under this license.
```

### 6.6.4 Member Data Documentation

### 6.6.4.1 bytesperpixel

```
int nsCamera.sensors.icarus2.icarus2.bytesperpixel = 2 [static]
```

Definition at line 30 of file icarus2.py.

#### 6.6.4.2 ca

```
nsCamera.sensors.icarus2.icarus2.ca
```

Definition at line 48 of file icarus2.py.

### 6.6.4.3 columns

```
int nsCamera.sensors.icarus2.icarus2.columns = 1 [static]
```

Definition at line 44 of file icarus2.py.

### 6.6.4.4 detect

```
str nsCamera.sensors.icarus2.icarus2.detect = "ICARUS_DET" [static]
```

Definition at line 33 of file icarus2.py.

#### 6.6.4.5 firstframe

```
int nsCamera.sensors.icarus2.icarus2.firstframe = 0 [static]
```

Definition at line 36 of file icarus2.py.

# 6.6.4.6 firstrow

```
int nsCamera.sensors.icarus2.icarus2.firstrow = 0 [static]
```

Definition at line 41 of file icarus2.py.

# 6.6.4.7 fpganumID

```
int nsCamera.sensors.icarus2.icarus2.fpganumID = 1 [static]
```

Definition at line 32 of file icarus2.py.

### 6.6.4.8 height

```
int nsCamera.sensors.icarus2.icarus2.height = 1024 [static]
```

Definition at line 40 of file icarus2.py.

# 6.6.4.9 icarustype

```
int nsCamera.sensors.icarus2.icarus2.icarustype = 0 [static]
```

Definition at line 31 of file icarus2.py.

#### 6.6.4.10 interlacing

```
list nsCamera.sensors.icarus2.icarus2.interlacing = [0, 0] [static]
```

Definition at line 43 of file icarus2.py.

### 6.6.4.11 lastframe

```
int nsCamera.sensors.icarus2.icarus2.lastframe = 3 [static]
```

Definition at line 37 of file icarus2.py.

### 6.6.4.12 lastrow

```
int nsCamera.sensors.icarus2.icarus2.lastrow = 1023 [static]
```

Definition at line 42 of file icarus2.py.

# 6.6.4.13 loglabel

```
str nsCamera.sensors.icarus2.icarus2.loglabel = "[Icarus2] " [static]
```

Definition at line 35 of file icarus2.py.

#### 6.6.4.14 maxframe

```
int nsCamera.sensors.icarus2.icarus2.maxframe = 3 [static]
```

Definition at line 27 of file icarus2.py.

### 6.6.4.15 maxheight

```
int nsCamera.sensors.icarus2.icarus2.maxheight = 1024 [static]
```

Definition at line 29 of file icarus2.py.

#### 6.6.4.16 maxwidth

```
int nsCamera.sensors.icarus2.icarus2.maxwidth = 512 [static]
```

Definition at line 28 of file icarus2.py.

#### 6.6.4.17 minframe

```
int nsCamera.sensors.icarus2.icarus2.minframe = 0 [static]
```

Definition at line 26 of file icarus2.py.

### 6.6.4.18 nframes

```
int nsCamera.sensors.icarus2.icarus2.nframes = 4 [static]
```

Definition at line 38 of file icarus2.py.

# 6.6.4.19 padToFull

```
bool nsCamera.sensors.icarus2.icarus2.padToFull = True [static]
```

Definition at line 45 of file icarus2.py.

# 6.6.4.20 sens\_registers

```
nsCamera.sensors.icarus2.icarus2.sens_registers
```

Definition at line 51 of file icarus2.py.

# 6.6.4.21 sens\_subregisters

```
nsCamera.sensors.icarus2.icarus2.sens_subregisters
```

Definition at line 75 of file icarus2.py.

### 6.6.4.22 sensfam

```
str nsCamera.sensors.icarus2.icarus2.sensfam = "Icarus" [static]
```

Definition at line 34 of file icarus2.py.

# 6.6.4.23 specwarn

```
str nsCamera.sensors.icarus2.icarus2.specwarn = "" [static]
```

Definition at line 25 of file icarus2.py.

#### 6.6.4.24 width

```
int nsCamera.sensors.icarus2.icarus2.width = 512 [static]
```

Definition at line 39 of file icarus2.py.

The documentation for this class was generated from the following file:

• C:/Users/hill35/git/camera\_python/nsCamera/sensors/icarus2.py

# 6.7 nsCamera.boards.LLNL v1.llnl v1 Class Reference

#### **Public Member Functions**

- \_\_init\_\_ (self, camassem)
- initBoard (self)
- initPots (self)
- latchPots (self)
- initSensor (self)
- configADCs (self)
- softReboot (self)
- disarm (self)
- startCapture (self, mode="Hardware")
- readSRAM (self)
- waitForSRAM (self, timeout)
- getTimer (self)
- resetTimer (self)
- enableLED (self, status)
- setLED (self, LED, status)
- setPowerSave (self, status)
- setPPER (self, pollperiod)
- getTemp (self, scale=None, offset=None, slope=None)
- getPressure (self, offset, sensitivity, units)
- · clearStatus (self)
- checkStatus (self)
- · checkStatus2 (self)
- reportStatus (self)
- reportEdgeDetects (self)
- · dumpStatus (self)

### **Public Attributes**

- ca
- logcrit
- logerr
- logwarn
- loginfo
- logdebug
- VREF
- ADC5 mult
- ADC5\_bipolar
- rs422 baud
- rs422\_cmd\_wait
- icarus\_subreg\_aliases
- icarus\_monitor\_controls
- daedalus\_subreg\_aliases
- · daedalus\_monitor\_controls
- · subreglist

#### **Static Public Attributes**

- · registers
- · list subregisters

# 6.7.1 Detailed Description

```
Livermore LLNL v1.0 board

Compatible communication protocols: RS422, GigE

Compatible sensors: icarus, icarus2, daedalus
```

Definition at line 28 of file LLNL\_v1.py.

#### 6.7.2 Constructor & Destructor Documentation

### 6.7.2.1 \_\_init\_\_()

```
nsCamera.boards.LLNL_v1.llnl_v1.__init__ (
                  self,
                  camassem )
Definition at line 192 of file LLNL_v1.py.
00192
           def __init__(self, camassem):
00193
                self.ca = camassem
00194
                self.logcrit = self.ca.logcritbase + "[LLNL_v1] "
                self.logerr = self.ca.logerrbase + "[LLNL_v1]
00195
                self.logwarn = self.ca.logwarnbase + "[LLNL_v1] "
self.loginfo = self.ca.loginfobase + "[LLNL_v1] "
00196
00197
00198
                self.logdebug = self.ca.logdebugbase + "[LLNL_v1] "
00199
                logging.info(self.loginfo + "initializing board object")
                self.VREF = 2.5  # default
self.ADC5_mult = 2  # i.e., monmax = 2 * VREF
00200
00201
00202
                \# False => monitor range runs 0 to monmax, True => +/- monmax
00203
                self.ADC5_bipolar = True
00204
                self.rs422\_baud = 921600
00205
                self.rs422\_cmd\_wait = 0.3
00206
                fpgaNum_pkt = Packet(cmd="1", addr=self.registers["FPGA_NUM"])
fpgaRev_pkt = Packet(cmd="1", addr=self.registers["FPGA_REV"])
00207
00208
00209
                _, _ = self.ca.sendCMD(fpgaNum_pkt) # dummy duplicate call err, rval = self.ca.sendCMD(fpgaNum_pkt)
00210
00211
00212
                self.ca.FPGANum = rval[8:16]
00213
00214
                err, rval = self.ca.sendCMD(fpgaRev_pkt)
00215
                self.ca.FPGAVersion = rval[8:16]
00216
00217
                # map channels to signal names for abstraction at the camera assembler level;
00218
                    each requires a corresponding entry in 'subregisters'
00219
00220
                self.icarus_subreg_aliases = OrderedDict(
00221
00222
                          "COL_BOT_IBIAS_IN": "POT1",
                          "HST_A_PDELAY": "POT2",
"HST_B_NDELAY": "POT3",
00223
00224
00225
                          "HST_RO_IBIAS": "POT4",
                          "HST OSC VREF IN": "POT5",
00226
                          "HST_B_PDELAY": "POT6",
"HST_OSC_CTL": "POT7",
00227
00228
                          "HST_A_NDELAY": "POT8",
00229
                          "COL_TOP_IBIAS_IN": "POT9",
00230
00231
                          "HST_OSC_R_BIAS": "POT10",
```

```
00232
                         "VAB": "POT11",
                         "HST_RO_NC_IBIAS": "POT12",
00233
00234
                         "VRST": "POT13",
                         "MON_HST_A_PDELAY": "MON_CH2",
00235
                         "MON_HST_B_NDELAY": "MON_CH3",
"MON_HST_RO_IBIAS": "MON_CH4",
00236
00237
00238
                         "MON_HST_OSC_VREF_IN": "MON_CH5",
00239
                         "MON_HST_B_PDELAY": "MON_CH6",
00240
                         "MON_HST_OSC_CTL": "MON_CH7"
00241
                         "MON_HST_A_NDELAY": "MON_CH8",
00242
00243
00244
               # Read-only; identifies controls corresponding to monitors
00245
               self.icarus_monitor_controls = OrderedDict(
00246
                   {
00247
                         "MON_CH2": "POT2",
                         "MON_CH3": "POT3",
00248
                         "MON_CH4": "POT4",
00249
                         "MON_CH5": "POT5",
00250
                         "MON_CH6": "POT6",
00251
                         "MON_CH7": "POT7",
"MON_CH8": "POT8",
00252
00253
                        # Note: VRST is not measured across the pot; it will read a voltage # approximately 1 Volt lower than pot13's actual output
00254
00255
00256
                         "MON_VRST": "POT13",
00257
                    }
00258
               )
00259
00260
               self.daedalus_subreg_aliases = OrderedDict(
00261
00262
                         "HST OSC CTL": "POT4",
                         "HST_RO_NC_IBIAS": "POT5",
00263
                         "HST_OSC_VREF_IN": "POT6",
"VAB": "POT11",
00264
00265
                         "MON_TSENSEOUT": "MON_CH2",
00266
                         "MON_BGREF": "MON_CH3",
00267
                         "MON_HST_OSC_CTL": "MON_CH4",
00268
                         "MON_HST_RO_NC_IBIAS": "MON_CH5",
"MON_HST_OSC_VREF_IN": "MON_CH6",
"MON_COL_TST_IN": "MON_CH7",
00269
00270
00271
                         "MON_HST_OSC_PBIAS_PAD": "MON_CH8",
00272
00273
00274
               # Read-only; identifies controls corresponding to monitors
00275
00276
               self.daedalus_monitor_controls = OrderedDict(
00277
00278
                         "MON_CH4": "POT4", "MON_CH5": "POT5",
00279
                         "MON_CH6": "POT6",
00280
00281
                         # Note: VRST is not measured across the pot; it will read a voltage
00282
                            lower than pot13's actual output
00283
                         "MON_VRST": "POT13",
00284
00285
               )
00286
00287
               self.subreglist = []
00288
               for s in self.subregisters:
00289
                    self.subreglist.append(s[0].upper())
00290
                    sr = SubRegister(
00291
                        self,
00292
                        name=s[0].upper(),
00293
                         register=s[1].upper(),
00294
                        start_bit=s[2],
00295
                         width=s[3],
00296
                        writable=s[4],
00297
00298
                    setattr(self, s[0].upper(), sr)
00299
00300
                # set voltage ranges for all pots
00301
               for n in range(1, 13):
00302
                    potname = "POT" + str(n)
00303
                    potobj = getattr(self, potname)
00304
                    potobj.minV = 0
                    potobj.maxV = 3.3
00305
00306
                    # resolution is approximately .0129 V / LSB
                    potobj.resolution = (1.0 * potobj.maxV - potobj.minV) / potobj.max_value
00307
00308
               self.POT13.minV = 0
               self.POT13.maxV = 3.96
00309
                # POT13 resolution is approximately .0155 V / LSB
00310
               self.POT13.resolution = (
00311
                    1.0 * self.POT13.maxV - self.POT13.minV
00312
```

```
00313 ) / self.POT13.max_value 00314
```

# 6.7.3 Member Function Documentation

#### 6.7.3.1 checkStatus()

```
nsCamera.boards.LLNL_v1.llnl_v1.checkStatus (
               self )
Check status register, convert to reverse-order bit stream (i.e., bit 0 is
  statusbits[0])
Returns:
    bit string (no '0b') in reversed order
Definition at line 729 of file LLNL v1.py.
00729
         def checkStatus(self):
00730
00731
              Check status register, convert to reverse-order bit stream (i.e., bit {\tt 0} is
00732
00733
00734
             Returns:
             bit string (no '0b') in reversed order
00735
00736
00737
              logging.debug(self.logdebug + "checkStatus")
00738
              err, rval = self.ca.getRegister("STAT_REG")
             if not rval:
00739
00740
                 logging.error(
00741
                     self.logerr + "Unable to check status register (zeroes returned)"
00742
00743
                 rval = "0"
00744
             rvalbits = bin(int(rval, 16))[2:].zfill(32)
             statusbits = rvalbits[::-1]
00745
00746
             return statusbits # TODO: add error handling
00747
```

#### 6.7.3.2 checkStatus2()

### Definition at line 748 of file LLNL\_v1.py.

```
00748
         def checkStatus2(self):
00749
00750
              Check second status register, convert to reverse-order bit stream (i.e., bit 0
00751
                is statusbits[0])
00752
00753
              Returns: bit string (no '0b') in reversed order
00754
              logging.debug(self.logdebug + "checkStatus2")
00755
00756
              err, rval = self.ca.getRegister("STAT_REG2")
00757
              if not rval:
00758
                 logging.error(
00759
                      self.logerr + "Unable to check status register 2 (zeroes returned)"
00760
00761
                  rval = "0"
              rvalbits = bin(int(rval, 16))[2:].zfill(5)
00762
00763
              statusbits = rvalbits[::-1]
00764
              return statusbits # TODO: add error handling
00765
```

### 6.7.3.3 clearStatus()

nsCamera.boards.LLNL\_v1.llnl\_v1.clearStatus (

```
self )
Check status registers to clear them
Returns:
     error string
Definition at line 714 of file LLNL v1.py.
00714
            def clearStatus(self):
00715
00716
                 Check status registers to clear them
00717
00718
                 Returns:
                 error string
00719
00720
00721
                 logging.debug(self.logdebug + "clearStatus")
00722
                 err1, rval = self.ca.getRegister("STAT_REG_SRC")
00723
                 err2, rval = self.ca.getRegister("STAT_REG2_SRC")
00724
                 err = err1 + err2
00725
00726
                     logging.error(self.logerr + "clearStatus failed")
00727
                 return err
00728
6.7.3.4 configADCs()
nsCamera.boards.LLNL_v1.llnl_v1.configADCs (
                   self )
Sets default ADC configuration (does not latch settings)
Returns:
     tuple (error string, response string) from final control message
Definition at line 437 of file LLNL_v1.py.
00437
            def configADCs(self):
00438
00439
                 Sets default ADC configuration (does not latch settings)
00440
00441
                 tuple (error string, response string) from final control message
00442
00443
00444
                 logging.info(self.loginfo + "configADCs")
00445
00446
                 control messages = [
00447
                     # just in case ADC_RESET was set (pull all ADCs out # of reset)
00448
                      ("ADC_RESET", "00000000"),
                      # workaround for uncertain behavior after previous readoff
00449
                      """ workaround for uncertain behavi
("ADC1_CONFIG_DATA", "FFFFFFFF"),
("ADC3_CONFIG_DATA", "FFFFFFFF"),
("ADC4_CONFIG_DATA", "FFFFFFFFF"),
00450
00451
00452
00453
00454
                      ("ADC_CTL", "FFFFFFFF"),
                      ("ADC_CIL", "FFFFFFFF"),

("ADC1_CONFIG_DATA", "81A801FF"), # ext Vref 1.25V

("ADC2_CONFIG_DATA", "81A801FF"), # ext Vref 1.25V

("ADC3_CONFIG_DATA", "81A801FF"), # ext Vref 1.25V

("ADC4_CONFIG_DATA", "81A801FF"), # ext Vref 1.25V

("ADC5_CONFIG_DATA", "81A883FF"), # int Vref 2.50V
00455
00456
00457
00458
00459
00460
00461
                 return self.ca.submitMessages(control_messages, " configADCs: ")
00462
```

#### 6.7.3.5 disarm()

```
nsCamera.boards.LLNL_v1.llnl_v1.disarm (
               self )
Takes camera out of trigger wait state. Has no effect if camera is not already
  in wait state.
Returns:
    tuple (error string, response string) from final control message
Definition at line 475 of file LLNL v1.py.
00475
         def disarm(self):
00476
00477
             Takes camera out of trigger wait state. Has no effect if camera is not already
00478
               in wait state.
00479
00480
             Returns:
             tuple (error string, response string) from final control message _{\tt nun}
00481
00482
00483
             logging.info(self.loginfo + "disarm")
00484
             self.ca.clearStatus()
00485
             self.ca.armed = False
00486
             control messages = [
                 ("HW_TRIG_EN", "0"),
00487
00488
                  ("SW_TRIG_EN", "0"),
00489
00490
             return self.ca.submitMessages(control_messages, " disarm: ")
00491
6.7.3.6 dumpStatus()
nsCamera.boards.LLNL_v1.llnl_v1.dumpStatus (
               self )
Create dictionary of status values, DAC settings, monitor values, and register
WARNING: the behavior of self-resetting subregisters may be difficult to predict
  and may generate contradictory results
    dictionary of system diagnostic values
Definition at line 824 of file LLNL v1.py.
00824
         def dumpStatus(self):
00825
00826
             Create dictionary of status values, DAC settings, monitor values, and register
00827
00828
00829
             WARNING: the behavior of self-resetting subregisters may be difficult to predict
00830
               and may generate contradictory results
00831
00832
             Returns:
             dictionary of system diagnostic values
00833
00834
00835
             statusbits = self.checkStatus()
00836
             statusbits2 = self.checkStatus2()
             temp = self.ca.getTemp()
00837
00838
             statDict = OrderedDict(
00839
00840
                      "Temperature reading": "{0:1.2f}".format(temp) + " C",
00841
```

"Sensor read complete": str(statusbits[0]),

00842

```
"Coarse trigger detected": str(statusbits[1]),
00843
00844
                       "Fine trigger detected": str(statusbits[2]),
00845
                       "Sensor readout in progress": str(statusbits[5]),
                       "Sensor readout complete": str(statusbits[6]),
00846
00847
                       "SRAM readout started": str(statusbits[7]),
00848
                       "SRAM readout complete": str(statusbits[8]),
00849
                       "High-speed timing configured": str(statusbits[9]),
00850
                       "All ADCs configured": str(statusbits[10]),
00851
                       "All pots configured": str(statusbits[11])
00852
                       "HST_All_W_En detected": str(statusbits[12]),
00853
                       "Timer has reset": str(statusbits[13]),
00854
                       "Camera is Armed": str(statusbits[14]),
00855
                       "FPA_IF_TO": str(statusbits2[0]),
                       "SRAM_RO_TO": str(statusbits2[1]),
00857
                       "PixelRd Timeout Error": str(statusbits2[2]),
00858
                       "UART_TX_TO_RST": str(statusbits2[3]),
00859
                       "UART_RX_TO_RST": str(statusbits2[4]),
00860
                  }
00861
              )
00862
              POTDict = OrderedDict()
MonDict = OrderedDict()
00863
00864
00865
              for entry in self.subreg_aliases:
00866
                   if self.subreg_aliases[entry][0] == "P":
                      val = str(round(self.ca.getPotV(entry), 3)) + " V"
POTDict["POT_" + entry] = val
00867
00868
00869
                  else:
00870
                      val = str(round(self.ca.getMonV(entry), 3)) + " V"
00871
                       MonDict[entry] = val
00872
00873
              regDict = OrderedDict()
00874
              for key in self.registers.keys():
00875
                  err, rval = self.ca.getRegister(key)
00876
                   regDict[key] = rval
00877
00878
              dumpDict = OrderedDict()
              for x in [statDict, MonDict, POTDict, regDict]:
00879
00880
                  dumpDict.update(x)
00881
              return dumpDict
00882
00883
00884 """
00885 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00886 LLNL-CODE-838080
00887
00888 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00889 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DDE)
00890 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00891 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00892 be made under this license. 00893 """
```

## 6.7.3.7 enableLED()

```
00606
00607
00608
                 status: 0 for disabled, 1 for enabled
00609
00610
             tuple: (error string, response string from setSubregister()
00611
00612
00613
             logging.info(self.loginfo + "enableLED")
00614
             if status:
                 status = 1
00615
00616
             return self.ca.setSubregister("LED_EN", str(status))
00617
```

#### 6.7.3.8 getPressure()

```
nsCamera.boards.LLNL_v1.llnl_v1.getPressure (
              self,
              offset,
              sensitivity,
              units )
Read pressure sensor
Currently unimplemented
Returns:
   0 as float
```

#### Definition at line 700 of file LLNL\_v1.py.

```
def getPressure(self, offset, sensitivity, units):
00700
00701
00702
              Read pressure sensor
00703
00704
              Currently unimplemented
00705
00706
              Returns:
              0 as float
00707
00708
00709
              logging.warning(
00710
                  "WARNING: [LLNL_v1] 'getPressure' is not implemented on the LLNLv1 board"
00711
00712
             return 0.0
00713
```

## 6.7.3.9 getTemp()

```
nsCamera.boards.LLNL_v1.llnl_v1.getTemp (
             self,
              scale = None,
              offset = None,
             slope = None )
Read temperature sensor
Args:
   scale: temperature scale to report (defaults to C, options are F and K)
   offset: unused
   slope: unused
Returns:
   temperature as float on given scale
```

```
Definition at line 672 of file LLNL_v1.py.
```

```
00672
          def getTemp(self, scale=None, offset=None, slope=None):
00673
              Read temperature sensor
00674
00675
              Args:
00676
                 scale: temperature scale to report (defaults to C, options are F and K)
00677
                  offset: unused
00678
                  slope: unused
00679
              Returns:
              temperature as float on given scale
00680
00681
00682
              logging.debug(self.logdebug + "getTemp: scale = " + str(scale))
              err, rval = self.ca.getRegister("TEMP_SENSE_DATA")
00683
00684
00685
                  logging.error(
00686
                     self.logerr + "unable to retrieve temperature information ("
00687
                      'getTemp), returning "0" '
00688
00689
                  return 0.0
00690
00691
              ctemp = int(rval[-3:], 16) / 16.0
00692
              if scale == "K":
00693
                 temp = ctemp + 273.15
00694
              elif scale == "F":
00695
                 temp = 1.8 * ctemp + 32
00696
              else:
00697
                 temp = ctemp
00698
              return temp
00699
```

#### 6.7.3.10 getTimer()

```
nsCamera.boards.LLNL_v1.llnl_v1.getTimer (
              self )
Read value of on-board timer
Returns:
   timer value as integer
```

#### Definition at line 575 of file LLNL v1.py. 00575

```
def getTimer(self):
00576
00577
              Read value of on-board timer
00578
00579
              Returns:
              timer value as integer
00580
00581
              logging.info(self.loginfo + "getTimer")
00582
              err, rval = self.ca.getRegister("TIMER_VALUE")
00583
00584
              if err:
                       self.logerr + "unable to retrieve timer information (getTimer), " 'returning "0" '
00585
                  logging.error(
00586
00587
00588
                  )
00589
                  return 0
00590
              return int(rval, 16)
00591
```

#### 6.7.3.11 initBoard()

```
nsCamera.boards.LLNL_v1.llnl_v1.initBoard (
              self )
```

```
Register and reset board, set up firmware for sensor
Returns:
    tuple (error string, response string) from final control message
Definition at line 315 of file LLNL v1.py.
           def initBoard(self):
00316
00317
                Register and reset board, set up firmware for sensor
00318
00319
                Returns:
                tuple (error string, response string) from final control message
00320
00321
                logging.info(self.loginfo + "initBoard LLNLv1")
00322
00323
                control_messages = [("LED_EN", "1")]
00324
00325
                self.clearStatus()
00326
                self.configADCs()
00327
                err, resp = self.ca.getSubregister("ADC5_VREF3")
00328
00329
                if err:
                    logging.error(self.logerr + "unable to read 'ADC5_VREF3'")
00330
                if int(resp, 2): # check to see if Vref is 3 or 2.5 volts
   vrefmax = 3.0
00331
00332
00333
                else:
00334
                   vrefmax = 2.5
00335
                err, resp = self.ca.getSubregister("ADC5_VREF")
00336
                if err:
                    logging.error(self.logerr + "unable to read 'ADC5_VREF'")
00337
                self.VREF = vrefmax * int(resp, 2) / 1024.0
00338
                err, multmask = self.ca.getSubregister("ADC5_MULT")
00339
00340
                if err:
00341
                    logging.error(self.logerr + "unable to read 'ADC5_MULT'")
00342
                 \begin{tabular}{ll} if & multmask[0] & and & multmask[1] & and & multmask[3] & and & multmask[5] : \\ \end{tabular} 
00343
                    self.ADC5 mult = 2
                 \textbf{elif not } (\texttt{multmask[0] or } \texttt{multmask[1] or } \texttt{multmask[3] or } \texttt{multmask[5]):} \\
00344
00345
                    self.ADC5_mult = 4
00346
                    logging.error(self.logerr + "inconsistent mode settings on ADC5")
00347
00348
                return self.ca.submitMessages(control_messages, " initBoard: ")
00349
6.7.3.12 initPots()
nsCamera.boards.LLNL_v1.llnl_v1.initPots (
                  self )
Configure default pot settings before image acquisition
    tuple (error string, response string) from final control message
Definition at line 350 of file LLNL v1.pv.
           def initPots(self):
00351
00352
                Configure default pot settings before image acquisition
00353
00354
                Returns:
               tuple (error string, response string) from final control message
00355
00356
00357
                logging.info(self.loginfo + "initPots")
00358
                if self.ca.sensorname == "icarus" or self.ca.sensorname == "icarus2":
                    err0, _ = self.ca.setPotV("HST_A_PDELAY", 0, errflag=True)
err1, _ = self.ca.setPotV("HST_B_NDELAY", 3.3, errflag=True)
err2, _ = self.ca.setPotV("HST_RO_IBIAS", 2.5, tune=True, errflag=True)
err3, _ = self.ca.setPotV("HST_SC_VREF_IN", 2.9, tune=True, errflag=True)
00359
00360
00361
```

err4, \_ = self.ca.setPot("HST\_B\_PDELAY", 0, errflag=True)
err5, \_ = self.ca.setPotV("HST\_OSC\_CTL", 1.45, tune=True, errflag=True)
err6, \_ = self.ca.setPotV("HST\_A\_NDELAY", 3.3, errflag=True)

00362

00363 00364 00365

```
00366
                      err7, _ = self.ca.setPotV("VAB", 0.5, errflag=True)
                      err8, _ = self.ca.setPotV("HST_RO_NC_IBIAS", 2.5, errflag=True)
00367
                      err9, _ = self.ca.setPotV("VRST", 0.3, tune=True, errflag=True)
err = err0 + err1 + err2 + err3 + err4 + err5 + err6 + err7 + err8 + err9
00368
00369
                 else: # Daedalus
00370
00371
                      err0, _ = self.ca.setPotV("HST_OSC_CTL", 1.0, tune=True, errflag=True)
                      err1, _ = self.ca.setPotV("HST_RO_NC_IBIAS", 1.0, errflag=True)
err2, _ = self.ca.setPotV("HST_OSC_VREF_IN", 1.0, tune=True, errflag=True)
00372
00373
00374
                      err3, _ = self.ca.setPotV("VAB", 0.5, errflag=True)
00375
                      err = err0 + err1 + err2 + err3
00376
                 return err, ""
00377
```

#### 6.7.3.13 initSensor()

#### Definition at line 404 of file LLNL v1.py.

```
def initSensor(self):
00405
00406
               Register sensor, set default timing settings
00407
00408
               Returns:
              tuple (error string, response string) from final control message
00409
00410
               logging.info(self.loginfo + "initSensor")
00411
00412
               if int(self.ca.FPGANum[7]) != self.ca.sensor.fpganumID:
00413
                   logging.error(
00414
                       self.logerr + "unable to confirm sensor compatibility with FPGA"
00415
00416
               self.registers.update(self.ca.sensor.sens_registers)
00417
               self.subregisters.extend(self.ca.sensor.sens_subregisters)
00418
               for s in self.ca.sensor.sens_subregisters:
00419
                  sr = SubRegister(
00420
                       self,
00421
                       name=s[0].upper(),
                       register=s[1].upper(),
00422
                       start_bit=s[2],
00423
                       width=s[3],
00424
00425
                       writable=s[4],
00426
00427
                   setattr(self, s[0].upper(), sr)
00428
                   self.subreglist.append(s[0])
              self.ca.checkSensorVoltStat()
00429
00430
               control_messages = self.ca.sensorSpecific() + [
                  # ring w/caps=01, relax=00, ring w/o caps = 02
("OSC_SELECT", "00"),
("FPA_DIVCLK_EN_ADDR", "00000001"), # TODO Make this a subregister
00431
00432
00433
00434
              return self.ca.submitMessages(control_messages, " initSensor: ")
00435
00436
```

## 6.7.3.14 latchPots()

```
nsCamera.boards.LLNL_v1.llnl_v1.latchPots ( self \ )
```

```
Latch pot settings into sensor
Returns:
     tuple (error string, response string) from final control message
Definition at line 378 of file LLNL v1.py.
            def latchPots(self):
00380
                 Latch pot settings into sensor
00381
00382
                 Returns:
                 tuple (error string, response string) from final control message
00383
00384
00385
                 logging.info(self.loginfo + "latchPots")
00386
00387
                 control_messages = [
                      ("POT_CTL", "00000003"),  # latches register settings for pot 1 ("POT_CTL", "00000005"),
00388
00389
00390
                      ("POT_CTL", "00000007"),
                      ("POT_CTL", "00000009"),
00391
                     ("POT_CTL", "00000009"),
("POT_CTL", "0000000B"),
("POT_CTL", "00000000"),
("POT_CTL", "00000001"),
("POT_CTL", "00000011"),
00392
00393
00394
00395
                      ("POT_CTL", "00000013"),
("POT_CTL", "00000015"),
00396
00397
                      ("POT_CTL", "00000017"),
("POT_CTL", "00000019"),
00398
00399
                      ("POT_CTL", "0000001B"),
00400
00401
00402
                 return self.ca.submitMessages(control_messages, " latchPots: ")
00403
```

#### 6.7.3.15 readSRAM()

```
nsCamera.boards.LLNL_v1.llnl_v1.readSRAM (
               self )
Start readoff of SRAM
Returns:
   tuple (error string, response string from register set)
Definition at line 529 of file LLNL v1.py.
         def readSRAM(self):
00530
00531
             Start readoff of SRAM
00532
00533
             Returns:
             tuple (error string, response string from register set)
00534
00535
             logging.info(self.loginfo + "readSRAM")
00537
             control_messages = [("READ_SRAM", "1")]
00538
             return self.ca.submitMessages(control_messages, " readSRAM: ")
00539
```

## 6.7.3.16 reportEdgeDetects()

```
\label{local_norm} \verb| nsCamera.boards.LLNL_v1.llnl_v1.reportEdgeDetects ( | self | )
```

Unimplemented

## Definition at line 815 of file LLNL\_v1.py.

```
00815 def reportEdgeDetects(self):
00816 """
00817 Unimplemented
00818 """
00819 logging.warning(
00820 self.logwarn + "'reportEdgeDetects' is not implemented on the LLNLv1 "
00821 "board "
00822
00823
```

#### 6.7.3.17 reportStatus()

```
\label{local_continuous_continuous} $$\operatorname{nsCamera.boards.LLNL\_v1.llnl\_v1.reportStatus} \ ($$\operatorname{\it self}$ )
```

Check contents of status register, print relevant messages

#### Definition at line 766 of file LLNL v1.py.

```
00766
          def reportStatus(self):
00767
00768
              Check contents of status register, print relevant messages
00769
00770
              statusbits = self.checkStatus()
              statusbits2 = self.checkStatus2()
00771
00772
              logging.info(self.loginfo + "Status report:")
00773
              if int(statusbits[0]):
00774
                  print(self.loginfo + "Sensor read complete")
00775
              if int(statusbits[1]):
00776
                  print(self.loginfo + "Coarse trigger detected")
00777
              if int(statusbits[2]):
00778
                  print(self.loginfo + "Fine trigger detected")
00779
              if int(statusbits[5]):
                  print(self.loginfo + "Sensor readout in progress")
00780
00781
              if int(statusbits[6]):
00782
                  print(self.loginfo + "Sensor readout complete")
00783
              if int(statusbits[7]):
00784
                  print(self.loginfo + "SRAM readout started")
00785
              if int(statusbits[8]):
00786
                  print(self.loginfo + "SRAM readout complete")
00787
              if int(statusbits[9]):
00788
                  print(self.loginfo + "High-speed timing configuration started")
00789
              if int(statusbits[10]):
                  print(self.loginfo + "All ADCs configured")
00790
00791
              if int(statusbits[11]):
00792
                  print(self.loginfo + "All pots configured")
00793
              if int(statusbits[13]):
                  print(self.loginfo + "Timer has reset")
00794
00795
              if int(statusbits[14]):
00796
                  print(self.loginfo + "Camera is Armed")
00797
              \verb|self.ca.sensor.reportStatusSensor(statusbits, statusbits2)|\\
00798
              temp = int(statusbits[27:15:-1], 2) / 16.0
00799
              logging.info(
00800
                  self.loginfo + "Temperature reading: " + "{0:1.2f}".format(temp) + " C"
00801
00802
              # press = int(statusbits[:27:-1], 2)
00803
              # logging.info(self.loginfo + "Pressure reading: " + "{0:1.2f}".format(press))
00804
              if int(statusbits2[0]):
00805
                  print(self.loginfo + "FPA_IF_TO")
00806
              if int(statusbits2[1]):
00807
                  print(self.loginfo + "SRAM_RO_TO")
00808
              if int(statusbits2[2]):
00809
                  print(self.loginfo + "PixelRd Timeout Error")
00810
              if int(statusbits2[3]):
                  print(self.loginfo + "UART_TX_TO_RST")
00811
              if int(statusbits2[4]):
00812
00813
                  print(self.loginfo + "UART_RX_TO_RST")
00814
```

## 6.7.3.18 resetTimer()

```
nsCamera.boards.LLNL_v1.llnl_v1.resetTimer (
               self )
Reset on-board timer
Returns:
    tuple (error string, response string from register set)
Definition at line 592 of file LLNL_v1.py.
         def resetTimer(self):
00592
00593
00594
              Reset on-board timer
00595
00596
              Returns:
             tuple (error string, response string from register set)
00597
00598
00599
             logging.info(self.loginfo + "resetTimer")
             control_messages = [("RESET_TIMER", "1"), ("RESET_TIMER", "0")]
00600
00601
              return self.ca.submitMessages(control_messages, " resetTimer: ")
00602
6.7.3.19 setLED()
nsCamera.boards.LLNL_v1.llnl_v1.setLED (
               self.
               LED,
               status )
Illuminate on-board LED
Args:
    LED: LED number (1-8)
    status: 0 is off, 1 is on
Returns:
    tuple: (error string, response string from setSubregister()
Definition at line 618 of file LLNL_v1.py.
00618
         def setLED(self, LED, status):
00619
00620
             Illuminate on-board LED
00621
00622
                 LED: LED number (1-8)
00623
00624
                 status: 0 is off, 1 is on
00625
00626
             Returns:
             tuple: (error string, response string from setSubregister()
00627
00628
             logging.info(self.loginfo + "setLED")
key = "LED" + str(LED)
00629
00630
             return self.ca.setSubregister(key, str(status))
00631
00632
```

## 6.7.3.20 setPowerSave()

```
nsCamera.boards.LLNL_v1.llnl_v1.setPowerSave (
               self,
               status )
Select powersave option
    status: setting for powersave option (1 is enabled)
Returns:
    tuple (error string, response string from setSubregister()
Definition at line 633 of file LLNL_v1.py.
00633
          def setPowerSave(self, status):
00634
             Select powersave option
00635
00636
00637
00638
                 status: setting for powersave option (1 is enabled)
00639
00640
              Returns:
             tuple (error string, response string from setSubregister()
00641
00642
             logging.info(self.loginfo + "setPowerSave")
00643
00644
              if status:
00645
                 status = 1
             return self.ca.setSubregister("POWERSAVE", str(status))
00646
00647
6.7.3.21 setPPER()
nsCamera.boards.LLNL_v1.llnl_v1.setPPER (
               self,
               pollperiod )
Set polling period for ADCs.
Args:
    pollperiod: milliseconds, between 1 and 255, defaults to 50
Returns:
    tuple (error string, response string from setSubregister OR invalid time
      setting string)
Definition at line 648 of file LLNL v1.py.
00648
         def setPPER(self, pollperiod):
00649
00650
              Set polling period for ADCs.
00651
00652
                 pollperiod: milliseconds, between 1 and 255, defaults to 50
00653
00654
              Returns:
00655
                 tuple (error string, response string from setSubregister OR invalid time
00656
                   setting string)
00657
00658
              logging.debug(self.logdebug + "setPPER: time = " + str(pollperiod))
              if pollperiod is None:
00659
                 pollperiod = 50
00660
              if not isinstance(pollperiod, int) or pollperiod < 1 or pollperiod > 255:
00661
00662
                 err = (
                     self.logerr + "invalid poll period submitted. Setting remains "
"unchanged."
00663
00664
00665
                 logging.error(err)
00666
00667
                 return err, str(pollperiod)
00668
              else:
00669
                 binset = bin(pollperiod)[2:].zfill(8)
                 return self.ca.setSubregister("PPER", binset)
00670
00671
```

## 6.7.3.22 softReboot()

```
nsCamera.boards.LLNL_v1.llnl_v1.softReboot (
               self )
Perform software reboot of board. WARNING: board reboot will likely prevent
  correct communication reponses and therefore will generate an error message
Returns:
    tuple (error string, response string) from final control message
Definition at line 463 of file LLNL_v1.py.
00463
         def softReboot(self):
00464
00465
             Perform software reboot of board. WARNING: board reboot will likely prevent
00466
               correct communication reponses and therefore will generate an error message
00467
00468
             tuple (error string, response string) from final control message
00469
00470
00471
             logging.info(self.loginfo + "reboot")
00472
             control_messages = [("RESET", "1")]
00473
             return self.ca.submitMessages(control_messages, " disarm: ")
00474
6.7.3.23 startCapture()
nsCamera.boards.LLNL_v1.llnl_v1.startCapture (
               self,
               mode = "Hardware" )
Selects trigger mode and enables board for image capture
    mode: trigger mode ("hardware"|"software"|"dual|"h"|"s"|"d" , is case-
      insensitive)
Returns:
    tuple (error string, response string) from final control message
Definition at line 492 of file LLNL v1.py.
00492
         def startCapture(self, mode="Hardware"):
00493
00494
             Selects trigger mode and enables board for image capture
00495
00496
00497
                 mode: trigger mode ("hardware"|"software"|"dual|"h"|"s"|"d", is case-
00498
                   insensitive)
00499
00500
             Returns:
             tuple (error string, response string) from final control message
00501
00502
00503
             logging.info(self.loginfo + "startCapture")
00504
             if self.ca.sensmanual:
                 timingReg = "MANSHUT_MODE"
00505
00506
             else:
00507
                 timingReg = "HST_MODE"
00508
00509
             if mode.upper()[0] == "S": # SOFTWARE
00510
                 trigmess = [
                     ("HW_TRIG_EN", "0"),
("SW_TRIG_EN", "1"),
00511
00512
                      ("SW_TRIG_START", "1"),
00513
```

```
00514
              else: # HARDWARE
00515
00516
                  trigmess = [
00517
                      ("SW_TRIG_EN", "0"),
00518
                      ("HW_TRIG_EN", "1"),
00519
00520
00521
              control_messages = [
00522
                  ("ADC_CTL", "0000001F"), # configure all ADCs
                  (timingReg, "1"),
00524
00525
00526
              control_messages.extend(trigmess)
              return self.ca.submitMessages(control_messages, " startCapture: ")
00528
```

## 6.7.3.24 waitForSRAM()

## Definition at line 540 of file LLNL\_v1.py.

```
00540
          def waitForSRAM(self, timeout):
00541
              Wait until subreg 'SRAM_READY' flag is true or timeout is exceeded;
00542
00543
                timeout = None or zero means wait indefinitely
00544
00545
                  timeout - time in seconds before readoff proceeds automatically without
00546
00547
                   waiting for SRAM_READY flag
00548
00549
              Returns:
              error string
00550
00551
00552
              logging.info(self.loginfo + "waitForSRAM")
00553
              waiting = True
00554
              starttime = time.time()
00555
              err = ""
00556
              while waiting:
00557
                  err, status = self.ca.getSubregister("SRAM_READY")
00558
00559
00560
                          self.logerr + "error in register read: " + err + " (waitForSRAM)"
00561
00562
                  if int(status):
00563
                      waiting = False
00564
                      logging.info(self.loginfo + "SRAM ready")
00565
                  if self.ca.abort:
00566
                      waiting = False
                      logging.info(self.loginfo + "readoff aborted by user")
00567
                      self.ca.abort = False
                  if timeout and time.time() - starttime > timeout:
00570
                      err += self.logerr + "SRAM timeout; proceeding with download attempt"
00571
                      logging.error(err)
00572
                      return err
00573
             return err
00574
```

## 6.7.4 Member Data Documentation

#### 6.7.4.1 ADC5 bipolar

```
nsCamera.boards.LLNL_v1.llnl_v1.ADC5_bipolar
```

Definition at line 203 of file LLNL\_v1.py.

## 6.7.4.2 ADC5\_mult

```
nsCamera.boards.LLNL_v1.llnl_v1.ADC5_mult
```

Definition at line 201 of file LLNL\_v1.py.

#### 6.7.4.3 ca

```
nsCamera.boards.LLNL_v1.llnl_v1.ca
```

Definition at line 193 of file LLNL v1.py.

#### 6.7.4.4 daedalus\_monitor\_controls

```
nsCamera.boards.LLNL_v1.llnl_v1.daedalus_monitor_controls
```

Definition at line 276 of file LLNL\_v1.py.

## 6.7.4.5 daedalus\_subreg\_aliases

```
nsCamera.boards.LLNL_v1.llnl_v1.daedalus_subreg_aliases
```

Definition at line 260 of file LLNL\_v1.py.

## 6.7.4.6 icarus\_monitor\_controls

```
\verb|nsCamera.boards.LLNL_v1.llnl_v1.icarus_monitor\_controls|\\
```

Definition at line 245 of file LLNL\_v1.py.

## 6.7.4.7 icarus\_subreg\_aliases

```
nsCamera.boards.LLNL_v1.llnl_v1.icarus_subreg_aliases
```

Definition at line 220 of file LLNL\_v1.py.

## 6.7.4.8 logcrit

```
nsCamera.boards.LLNL_v1.llnl_v1.logcrit
```

Definition at line 194 of file LLNL\_v1.py.

#### 6.7.4.9 logdebug

```
nsCamera.boards.LLNL_v1.llnl_v1.logdebug
```

Definition at line 198 of file LLNL\_v1.py.

#### 6.7.4.10 logerr

```
nsCamera.boards.LLNL_v1.llnl_v1.logerr
```

Definition at line 195 of file LLNL\_v1.py.

## 6.7.4.11 loginfo

```
nsCamera.boards.LLNL_v1.llnl_v1.loginfo
```

Definition at line 197 of file LLNL\_v1.py.

#### 6.7.4.12 logwarn

```
nsCamera.boards.LLNL_v1.llnl_v1.logwarn
```

Definition at line 196 of file LLNL\_v1.py.

## **6.7.4.13** registers

```
nsCamera.boards.LLNL_v1.llnl_v1.registers [static]
```

Definition at line 37 of file LLNL\_v1.py.

## 6.7.4.14 rs422\_baud

```
nsCamera.boards.LLNL_v1.llnl_v1.rs422_baud
```

Definition at line 204 of file LLNL\_v1.py.

## 6.7.4.15 rs422\_cmd\_wait

```
nsCamera.boards.LLNL_v1.llnl_v1.rs422_cmd_wait
```

Definition at line 205 of file LLNL\_v1.py.

## 6.7.4.16 subregisters

```
list nsCamera.boards.LLNL_v1.llnl_v1.subregisters [static]
```

Definition at line 98 of file LLNL v1.py.

#### 6.7.4.17 subreglist

```
nsCamera.boards.LLNL_v1.llnl_v1.subreglist
```

Definition at line 287 of file LLNL\_v1.py.

#### 6.7.4.18 VREF

```
nsCamera.boards.LLNL_v1.llnl_v1.VREF
```

Definition at line 200 of file LLNL\_v1.py.

The documentation for this class was generated from the following file:

C:/Users/hill35/git/camera\_python/nsCamera/boards/LLNL\_v1.py

# 6.8 nsCamera.boards.LLNL\_v4.llnl\_v4 Class Reference

#### **Public Member Functions**

- \_\_init\_\_ (self, camassem)
- · initBoard (self)
- initPots (self)
- · latchPots (self)
- initSensor (self)
- configADCs (self)
- softReboot (self)
- disarm (self)
- startCapture (self, mode="Hardware")
- readSRAM (self)
- waitForSRAM (self, timeout)
- getTimer (self)

- resetTimer (self)
- enableLED (self, status)
- setLED (self, LED, status)
- setPowerSave (self, status)
- setPPER (self, pollperiod)
- getTemp (self, scale=None)
- getPressure (self, offset=None, sensitivity=None, units=None)
- · clearStatus (self)
- checkStatus (self)
- checkStatus2 (self)
- reportStatus (self)
- reportEdgeDetects (self)
- dumpStatus (self)

#### **Public Attributes**

- ca
- logcrit
- logerr
- logwarn
- loginfo
- logdebug
- VREF
- ADC5\_mult
- ADC5\_bipolar
- rs422 baud
- rs422\_cmd\_wait
- defoff
- defsens
- icarus\_subreg\_aliases
- icarus\_monitor\_controls
- · daedalus\_subreg\_aliases
- daedalus\_monitor\_controls
- · subreglist

## **Static Public Attributes**

- registers
- list subregisters

## 6.8.1 Detailed Description

Livermore LLNL v4.0 board

Compatible communication protocols: RS422, GigE

Compatible sensors: icarus, icarus2, daedalus

Definition at line 29 of file LLNL\_v4.py.

## 6.8.2 Constructor & Destructor Documentation

#### 6.8.2.1 \_\_init\_\_()

```
nsCamera.boards.LLNL_v4.llnl_v4.__init__ (
                  self.
                  camassem )
Definition at line 181 of file LLNL v4.py.
00181
           def __init__(self, camassem):
00182
                self.ca = camassem
                self.logcrit = self.ca.logcritbase + "[LLNL_v4] "
00183
                self.logerr = self.ca.logerrbase + "[LLNL_v4]
00184
                self.logwarn = self.ca.logwarnbase + "[LLNL_v4] "
00185
                self.loginfo = self.ca.loginfobase + "[LLNL_v4] " self.logdebug = self.ca.logdebugbase + "[LLNL_v4] "
00186
00187
                logging.info(self.loginfo + "Initializing board object")
00188
00189
                self.VREF = 3.3 \# must be supplied externally for ADC128S102
                self.ADC5_mult = 1
00190
00191
00192
                \# ADC128S102; False => monitor range runs 0 to monmax, True => +/- monmax
00193
                self.ADC5_bipolar = False
00194
                self.rs422\_baud = 921600
00195
                self.rs422\_cmd\_wait = 0.3
00196
00197
                fpgaNum_pkt = Packet(cmd="1", addr=self.registers["FPGA_NUM"])
                fpgaRev_pkt = Packet(cmd="1", addr=self.registers["FPGA_REV"])
00198
00199
                _, _ = self.ca.sendCMD(fpgaNum_pkt) # dummy duplicate call
00200
00201
                err, rval = self.ca.sendCMD(fpgaNum_pkt)
00202
                self.ca.FPGANum = rval[8:16]
00203
00204
                err, rval = self.ca.sendCMD(fpgaRev_pkt)
00205
                self.ca.FPGAVersion = rval[8:16]
00206
00207
                self.defoff = 34.5 # default pressure sensor offset
                self.defsens = 92.5 # default pressure sensor sensitivity
00208
00209
00210
                # TODO: move to sensor scripts?
00211
                # map channels to signal names for abstraction at the camera assembler level;
00212
                   each requires a corresponding entry in 'subregisters'
00213
                self.icarus_subreg_aliases = OrderedDict(
00214
00215
                         "HST_A_PDELAY": "DACA",
                         "HST_A_NDELAY": "DACB",
00216
                         "HST_B_PDELAY": "DACC",
00217
                         "HST_B_NDELAY": "DACD",
00218
00219
                         "HST_RO_IBIAS": "DACE",
00220
                         "HST_RO_NC_IBIAS": "DACE",
00221
                          "HST_OSC_CTL": "DACF",
                         "VAB": "DACG",
00222
                         "VRST": "DACH",
00223
                         "MON_PRES_MINUS": "MON_CH1",
"MON_PRES_PLUS": "MON_CH2",
00224
00225
                          "MON_TEMP": "MON_CH3",
00226
                         "MON_COL_TOP_IBIAS_IN": "MON_CH4",
"MON_HST_OSC_R_BIAS": "MON_CH5",
00227
00228
                         "MON_UAB": "MON_CH6",
"MON_HST_RO_IBIAS": "MON_CH7",
00229
00230
                         "MON_HST_RO_NC_IBIAS": "MON_CH7",
"MON_VRST": "MON_CH8",
00231
                         "MON_COL_BOT_IBIAS_IN": "MON_CH9",
00233
                         "MON_HST_A_PDELAY": "MON_CH10",
"MON_HST_B_NDELAY": "MON_CH11",
00234
00235
                         "DOSIMETER": "MON_CH12",
00236
                         "MON_HST_OSC_CTL": "MON_CH13",
"MON_HST_OSC_CTL": "MON_CH15",
00237
00238
00239
                         "MON_HST_A_NDELAY": "MON_CH16",
"MON_CHA": "MON_CH10",
"MON_CHB": "MON_CH16",
00240
00241
00242
                         "MON_CHC": "MON_CH14",
00243
                         "MON_CHD": "MON_CH11",
00244
                         "MON_CHE": "MON_CH7",
"MON_CHF": "MON_CH15",
00245
00246
```

```
00247
                           "MON_CHG": "MON_CH6",
00248
                           "MON_CHH": "MON_CH8",
00249
00250
00251
                 # Read-only; identifies controls corresponding to monitors
00252
                 self.icarus_monitor_controls = OrderedDict(
00253
                     {
00254
                           "MON_CH10": "DACA",
00255
                           "MON_CH16": "DACB",
                           "MON_CH14": "DACC",
00256
00257
                           "MON_CH11": "DACD",
00258
                           "MON_CH7": "DACE",
                           "MON_CH15": "DACF",
"MON_CH6": "DACG",
00259
00261
                           "MON_CH8": "DACH",
00262
00263
00264
                 self.daedalus_subreg_aliases = OrderedDict(
00265
                           "HST_OSC_VREF_IN": "DACC",
00266
                           "HST_OSC_CTL": "DACE",
"COL_TST_IN": "DACF",
00267
00268
                           "VAB": "DACG",
"VRST": "DACH",
00269
00270
00271
                           "MON_PRES_MINUS": "MON_CH1",
                           "MON_PRES_PLUS": "MON_CH2",
00272
                           "MON_TEMP": "MON_CH3",
"MON_VAB": "MON_CH6",
00273
00274
                           "MON_HST_OSC_CTL": "MON_CH10",
"MON_TSENSE_OUT": "MON_CH10",
"MON_BGREF": "MON_CH11",
"DOSIMETER": "MON_CH12",
00275
00276
00277
00278
                           "MON_HST_RO_NC_IBIAS": "MON_CH13",
"MON_HST_OSC_VREF_IN": "MON_CH14",
00279
00280
                           "MON_COL_TST_IN": "MON_CH15",
"MON_HST_OSC_PBIAS_PAD": "MON_CH16",
00281
00282
                           "MON_CHC": "MON_CH14",
"MON_CHE": "MON_CH7",
00283
00284
                           "MON_CHF": "MON_CH15",
"MON_CHG": "MON_CH6",
00285
00286
                           "MON_CHH": "MON_CH8",
00287
00288
00289
00290
                 # Read-only; identifies controls corresponding to monitors
00291
                 self.daedalus\_monitor\_controls = OrderedDict(
00292
                           "MON_CH14": "DACC",
"MON_CH7": "DACE",
00293
00294
                           "MON_CH15": "DACF",
"MON_CH6": "DACG",
"MON_CH8": "DACH",
00295
00296
00297
00298
00299
00300
                 self.subreglist = []
00301
                 for s in self.subregisters:
00302
                      self.subreglist.append(s[0].upper())
00303
                      sr = SubRegister(
00304
                           self,
00305
                           name=s[0].upper(),
00306
                           register=s[1].upper(),
                           start_bit=s[2],
00307
00308
                           width=s[3],
00309
                           writable=s[4],
00310
00311
                      setattr(self, s[0].upper(), sr)
00312
00313
                 # set voltage ranges for all DACs - WARNING: actual output voltage limited to
00314
                    external supply (3.3 V)
                 # setpot('potx', n) will generate 3.3 V for all n > .66
00315
                 for n in range(0, 8):
00316
                      potname = "DAC" + string.ascii_uppercase[n]
00317
00318
                      potobj = getattr(self, potname)
                      potobj.minV = 0
00319
                      potobj.maxV = 5 #
00320
                      potobj.resolution = (
00321
                     1.0 * potobj.maxV - potobj.minV) / potobj.max_value # 76 uV / LSB
00322
00323
00324
```

## 6.8.3 Member Function Documentation

## 6.8.3.1 checkStatus()

```
nsCamera.boards.LLNL_v4.llnl_v4.checkStatus (
               self )
Check status register, convert to reverse-order bit stream (i.e., bit 0 is
  statusbits[0])
Returns:
    bit string (no '0b') in reversed order
Definition at line 706 of file LLNL v4.py.
         def checkStatus(self):
00707
00708
             Check status register, convert to reverse-order bit stream (i.e., bit 0 is
00709
               statusbits[0])
00710
00711
             Returns:
             bit string (no '0b') in reversed order
00712
00713
00714
             err, rval = self.ca.getRegister("STAT_REG")
00715
             rvalbits = bin(int(rval, 16))[2:].zfill(32)
00716
             statusbits = rvalbits[::-1]
00717
             return statusbits
00718
```

## 6.8.3.2 checkStatus2()

#### Definition at line 719 of file LLNL v4.py.

```
00719
          def checkStatus2(self):
00720
00721
              Check second status register, convert to reverse-order bit stream (i.e., bit \boldsymbol{0}
00722
                is statusbits[0])
00723
              Returns: bit string (no '0b') in reversed order
00724
00725
00726
              err, rval = self.ca.getRegister("STAT_REG2")
00727
              rvalbits = bin(int(rval, 16))[2:].zfill(6)
00728
              statusbits = rvalbits[::-1]
              return statusbits
00729
00730
```

#### 6.8.3.3 clearStatus()

nsCamera.boards.LLNL\_v4.llnl\_v4.clearStatus (

```
self )
Check status registers to clear them
Returns:
     error string
Definition at line 692 of file LLNL v4.py.
           def clearStatus(self):
00692
00693
00694
                Check status registers to clear them
00695
00696
                Returns:
                error string
00697
00698
                err1, rval = self.ca.getRegister("STAT_REG_SRC")
err2, rval = self.ca.getRegister("STAT_REG2_SRC")
00699
00700
00701
                err = err1 + err2
                if err:
00702
00703
                     logging.error(self.logerr + "clearStatus failed")
00704
                return err
00705
6.8.3.4 configADCs()
nsCamera.boards.LLNL_v4.llnl_v4.configADCs (
                  self )
Sets default ADC configuration (does not latch settings)
     tuple (error string, response string) from final control message
Definition at line 402 of file LLNL v4.py.
00402
           def configADCs(self):
00403
00404
                Sets default ADC configuration (does not latch settings)
00405
00406
                tuple (error string, response string) from final control message
00407
00408
00409
                logging.info(self.loginfo + "configADCs")
00410
00411
                control_messages = [
00412
                     # just in case ADC_RESET was set on any of the ADCs (pull all ADCs out of
00413
00414
                     ("ADC_RESET", "00000000"),
00415
                     # workaround for uncertain behavior after previous readoff
                     ("ADC2_CONFIG_DATA", "FFFFFFFF"),
("ADC3_CONFIG_DATA", "FFFFFFFF"),
("ADC3_CONFIG_DATA", "FFFFFFFF"),
("ADC4_CONFIG_DATA", "FFFFFFFFF"),
00416
00417
00418
00419
                     ("ADC_CTL", "FFFFFFFF"),
00420
                     ("ADC1_CONFIG_DATA", "81A801FF"), # ext Vref 1.25V
("ADC2_CONFIG_DATA", "81A801FF"), # ext Vref 1.25V
("ADC3_CONFIG_DATA", "81A801FF"), # ext Vref 1.25V
00421
00422
00423
                     ("ADC4_CONFIG_DATA", "81A801FF"), # ext Vref 1.25V
00424
00425
00426
                return self.ca.submitMessages(control messages, " configADCs: ")
00427
```

## 6.8.3.5 disarm()

```
nsCamera.boards.LLNL_v4.llnl_v4.disarm (
               self )
Takes camera out of trigger wait state. Has no effect if camera is not in wait
  state.
Returns:
    tuple (error string, response string) from final control message
Definition at line 440 of file LLNL v4.py.
00440
         def disarm(self):
00441
00442
              Takes camera out of trigger wait state. Has no effect if camera is not in wait
00443
00444
00445
             tuple (error string, response string) from final control message
00446
00447
00448
              logging.info(self.loginfo + "disarm")
00449
              self.ca.clearStatus()
             self.ca.armed = False
00450
00451
             control_messages = [
                 ("HW_TRIG_EN", "0"),
00452
00453
                  ("SW_TRIG_EN", "0"),
00454
00455
              self.ca.comms.skipError = False
00456
              return self.ca.submitMessages(control_messages, " disarm: ")
00457
6.8.3.6 dumpStatus()
nsCamera.boards.LLNL_v4.llnl_v4.dumpStatus (
               self )
Create dictionary of status values, DAC settings, monitor values, and register
  values.
Returns:
    dictionary of system diagnostic values
Definition at line 821 of file LLNL v4.py.
         def dumpStatus(self):
00822
00823
              Create dictionary of status values, DAC settings, monitor values, and register
00824
               values.
00825
00826
              Returns:
             dictionary of system diagnostic values
00827
00828
00829
              statusbits = self.checkStatus()
00830
              statusbits2 = self.checkStatus2()
00831
00832
              temp = int(statusbits[23:16:-1], 2) \star 3.3 \star 1000 / 4096
00833
              press = int(statusbits[:23:-1], 2) \star 3.3 \star 1000 / 4096
00834
00835
              statDict = OrderedDict(
00836
                 {
00837
                      "Temperature sensor reading (Deg C)": "{0:1.2f}".format(temp),
                      "Pressure reading (Torr)": str(round(self.ca.getPressure(), 3)),
00838
```

"Pressure sensor reading (mV)": " $\{0:1.2f\}$ ".format(press),

"Sensor read complete": str(statusbits[0]),

"Coarse trigger detected": str(statusbits[1]),

00839

00840

00841

```
00842
                       "Fine trigger detected": str(statusbits[2]),
00843
                       "Sensor readout in progress": str(statusbits[5]),
00844
                       "Sensor readout complete": str(statusbits[6]),
                       "SRAM readout started": str(statusbits[7]),
00845
00846
                       "SRAM readout complete": str(statusbits[8]),
00847
                       "High-speed timing configured": str(statusbits[9]),
00848
                       "All ADCs configured": str(statusbits[10]),
00849
                       "All DACs configured": str(statusbits[11]),
00850
                       "Timer has reset": str(statusbits[13]),
                       "Camera is Armed": str(statusbits[14]),
00851
00852
                       "FPA_IF_TO": str(statusbits2[0]),
00853
                       "SRAM_RO_TO": str(statusbits2[1]),
00854
                       "PixelRd Timeout Error": str(statusbits2[2]),
                       "UART_TX_TO_RST": str(statusbits2[3]),
00856
                       "UART_RX_TO_RST": str(statusbits2[4]),
00857
                       "PDBIAS Unready": str(statusbits2[5]),
00858
                   }
00859
00860
00861
               if self.ca.sensorname == "icarus" or self.ca.sensorname == "icarus2":
00862
                   senslabs = {
00863
                       3: "W3_Top_A_Edge1 detected",
                       4: "W3_Top_B_Edgel detected",
12: "HST_All_W_En detected",
00864
00865
00866
00867
              else:
00868
                   senslabs = {
                       3: "RSLROWOUTA",
00869
00870
                       4: "RSLROWOUTB"
00871
                       12: "RSLNALLWENB",
00872
                       15: "RSLNALLWENA",
                       16: "Config HST is done",
00873
00874
00875
              sensDict = {senslabs[x]: str(statusbits[x]) for x in senslabs.keys()}
00876
00877
              DACDict = OrderedDict()
              MonDict = OrderedDict()
00878
00879
              for entry in self.subreg_aliases:
00880
                   if self.subreg_aliases[entry][0] == "D":
00881
                       val = str(round(self.ca.getPotV(entry), 3))
                       DACDict["DAC_" + entry] = val
00882
00883
00884
                       val = str(round(self.ca.getMonV(entry), 3))
00885
                       MonDict[entry] = val
00886
00887
               regDict = OrderedDict()
00888
               for key in self.registers.keys():
                   # Load in all registers except for the read-clear status registers.
if key == "STAT_REG_SRC" or key == "STAT_REG2_SRC":
00889
00890
00891
                   else:
00892
00893
                       err, rval = self.ca.getRegister(key)
00894
                       regDict[key] = rval
00895
00896
               dumpDict = OrderedDict()
00897
               for x in [statDict, sensDict, MonDict, DACDict, regDict]:
                   dumpDict.update(x)
00898
00899
               return dumpDict
00900
00901
00902 """
00903 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00904 LLNL-CODE-838080
00906 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00907 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00908 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00909 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00910 be made under this license.
00911 """
```

#### 6.8.3.7 enableLED()

```
Dummy function; feature is not implemented on LLNL_V4 board
Returns:
    tuple: dummy of (error string, response string from setSubregister())
Definition at line 571 of file LLNL v4.py.
          def enableLED(self, status):
00573
              Dummy function; feature is not implemented on LLNL_V4 board
00574
00575
              tuple: dummy of (error string, response string from setSubregister())
00576
00577
              del status return "", "0"
00578
00579
00580
6.8.3.8 getPressure()
nsCamera.boards.LLNL_v4.llnl_v4.getPressure (
                self,
               offset = None,
                sensitivity = None,
               units = None)
Read pressure sensor. Uses default offset and sensitivity defined in init
  function unless alternatives are specified. NOTE: to reset defaults, reassign
  board.defoff and board.defsens explicitly
Aras:
    offset: non-default offset in mv/V
    sensitivity: non-default sensitivity in {\rm mV/V/span}
    units: units to report pressure (defaults to Torr, options are psi, bar,
      inHq, atm)
    Pressure as float in chosen units, defaults to torr
Definition at line 653 of file LLNL_v4.py.
00653
          def getPressure(self, offset=None, sensitivity=None, units=None):
00654
00655
              Read pressure sensor. Uses default offset and sensitivity defined in init
00656
                function unless alternatives are specified. NOTE: to reset defaults, reassign
00657
                board.defoff and board.defsens explicitly
00658
00659
00660
                 offset: non-default offset in mv/V
00661
                  sensitivity: non-default sensitivity in mV/V/span
00662
                 units: units to report pressure (defaults to Torr, options are psi, bar,
00663
                   inHg, atm)
00664
00665
              Returns:
              Pressure as float in chosen units, defaults to torm _{\tt nun}
00666
00667
00668
              if offset is None:
00669
                 offset = self.defoff
00670
              if sensitivity is None:
00671
                 sensitivity = self.defsens
00672
              if units is None:
00673
                 units = "torr"
              pplus = self.ca.getMonV("MON_PRES_PLUS")
00674
00675
              pminus = self.ca.getMonV("MON_PRES_MINUS")
00676
              delta = 1000 * (pplus - pminus)
              ratio = sensitivity / 30 # nominal is 21/30 psi = (delta - offset) / ratio
00677
00678
00679
              if units.lower() == "psi":
```

```
00680
                 press = psi
00681
             elif units.lower() == "bar":
00682
                 press = psi / 14.504
              elif units.lower() == "atm":
00683
00684
                 press = psi / 14.695
00685
              elif units.lower() == "inHg":
00686
                 press = psi * 2.036
00687
00688
                 press = 51.715 * psi # default to Torr
00690
             return press
00691
```

#### 6.8.3.9 getTemp()

```
nsCamera.boards.LLNL_v4.llnl_v4.getTemp (
                self,
                scale = None )
Read temperature sensor
Args:
    scale: temperature scale to report (defaults to C, options are F and K)
Returns:
    temperature as float on given scale
Definition at line 628 of file LLNL_v4.py.
          def getTemp(self, scale=None):
00628
00629
              Read temperature sensor
00630
00631
                  scale: temperature scale to report (defaults to C, options are F and K)
00632
00633
00634
              Returns:
              temperature as float on given scale
00635
00636
              err, rval = self.ca.getMonV("MON_TEMP", errflag=True)
00637
00638
              if err:
00639
                 logging.error(
                      self.logerr + "unable to retrieve temperature information ("
    'getTemp), returning "0" '
00640
00641
00642
00643
                  return 0.0
             ctemp = rval * 1000 - 273.15
if scale == "K":
00644
00645
00646
                  temp = ctemp + 273.15
00647
              elif scale == "F":
00648
                  temp = 1.8 * ctemp + 32
00649
00650
                  temp = ctemp
00651
             return temp
00652
```

#### 6.8.3.10 getTimer()

```
Definition at line 543 of file LLNL_v4.py.
00543
         def getTimer(self):
00544
00545
              Read value of on-board timer
00546
00547
              Returns:
              timer value as integer
00548
00549
00550
              logging.info(self.loginfo + "getTimer")
00551
              err, rval = self.ca.getRegister("TIMER_VALUE")
00552
00553
                 logging.error(
00554
                     self.logerr + "unable to retrieve timer information (getTimer), "
00555
                     returning "0" '
00556
                 return 0
00558
              return int(rval, 16)
00559
6.8.3.11 initBoard()
nsCamera.boards.LLNL_v4.llnl_v4.initBoard (
               self )
Register and reset board, set up firmware for sensor
Returns:
    tuple (error string, response string) from final control message
Definition at line 325 of file LLNL_v4.py.
         def initBoard(self):
              Register and reset board, set up firmware for sensor
00327
00328
00329
             tuple (error string, response string) from final control message
00330
00331
00332
              logging.info(self.loginfo + "initBoard LLNLv4")
00333
             control_messages = []
00334
              self.clearStatus()
00335
              self.configADCs()
00336
              return self.ca.submitMessages(control_messages, " initBoard: ")
00337
6.8.3.12 initPots()
nsCamera.boards.LLNL_v4.llnl_v4.initPots (
               self )
Dummy function; initial DAC values are set by firmware at startup
Returns:
    tuple (empty string, empty string)
Definition at line 338 of file LLNL_v4.py.
00338
         def initPots(self):
00339
00340
              Dummy function; initial DAC values are set by firmware at startup
00341
00342
              Returns:
              tuple (empty string, empty string)
```

logging.debug(self.logdebug + "InitPots")
return "", ""

00343 00344

00345 00346 00347

## 6.8.3.13 initSensor()

```
nsCamera.boards.LLNL_v4.llnl_v4.initSensor (
                self )
Register sensor, set default timing settings
Returns:
    tuple (error string, response string) from final control message
Definition at line 368 of file LLNL v4.pv.
00368
          def initSensor(self):
00369
00370
              Register sensor, set default timing settings
00371
00372
              Returns:
              tuple (error string, response string) from final control message
00373
00374
00375
              logging.info(self.loginfo + "initSensor")
00376
              if int(self.ca.FPGANum[7]) != self.ca.sensor.fpganumID:
00377
                  logging.warning(
                      self.logwarn + "unable to confirm sensor compatibility with FPGA"
00378
00379
00380
              self.registers.update(self.ca.sensor.sens_registers)
00381
              self.subregisters.extend(self.ca.sensor.sens_subregisters)
00382
              for s in self.ca.sensor.sens_subregisters:
00383
                  sr = SubRegister(
00384
                      self,
00385
                      name=s[0].upper(),
00386
                      register=s[1].upper(),
00387
                      start_bit=s[2],
00388
                      width=s[3],
00389
                      writable=s[4],
00390
00391
                  setattr(self, s[0].upper(), sr)
00392
                  self.subreglist.append(s[0])
00393
              # TODO: self.ca.checkSensorVoltStat() # SENSOR_VOLT_STAT and SENSOR_VOLT_CTL are
00394
                  deactivated for v4 icarus and daedalus firmware for now, is this permanent?
00395
              control_messages = self.ca.sensorSpecific() + [
                  # ring w/caps=01, relax=00, ring w/o caps = 02
("OSC_SELECT", "00"),
00396
00397
                  ("FPA_DIVCLK_EN_ADDR", "00000001"),
00398
00399
00400
              return self.ca.submitMessages(control_messages, " initSensor: ")
00401
```

## 6.8.3.14 latchPots()

00352

00353 00354 00355 Returns:

logging.info(self.loginfo + "latchPots")

tuple (error string, response string) from final control message

```
00356
                       control_messages = [
                          ("DAC_CTL", "00000001"), # latches register settings for DACA ("DAC_CTL", "00000003"),
00357
00358
                             ("DAC_CTL", "00000003"),

("DAC_CTL", "00000005"),

("DAC_CTL", "00000007"),

("DAC_CTL", "00000009"),

("DAC_CTL", "0000000B"),

("DAC_CTL", "0000000D"),

("DAC_CTL", "0000000F"),
00359
00360
00361
00362
00363
00364
00366
                       return self.ca.submitMessages(control_messages, " latchPots: ")
00367
6.8.3.15 readSRAM()
```

```
nsCamera.boards.LLNL_v4.llnl_v4.readSRAM (
              self )
Start readoff of SRAM
Returns:
   tuple (error string, response string from register set)
```

## Definition at line 495 of file LLNL\_v4.py.

```
def readSRAM(self):
00495
00496
00497
              Start readoff of SRAM
00498
00499
              Returns:
              tuple (error string, response string from register set)
00500
0.0501
              logging.info(self.loginfo + "readSRAM")
00502
00503
              control_messages = [("READ_SRAM", "1")]
              return self.ca.submitMessages(control_messages, " readSRAM: ")
00504
00505
```

#### 6.8.3.16 reportEdgeDetects()

```
nsCamera.boards.LLNL_v4.llnl_v4.reportEdgeDetects (
              self )
Report edge detects
```

#### Definition at line 788 of file LLNL\_v4.py.

```
def reportEdgeDetects(self):
00790
                 Report edge detects
00791
00792
                 err, rval = self.ca.getRegister("STAT_EDGE_DETECTS")
# shift to left to fake missing edge detect
00793
00794
                 edgebits = bin(int(rval, 16) « 1)[2:].zfill(32)
00795
                 # reverse to get order matching assignment
00796
                 bitsrev = edgebits[::-1]
                 detdict = {}
00797
00798
                 bitidx = 0
                 for frame in range(4):
    for vert in ("TOP", "BOT"):
00799
00800
                           for edge in range(1, 3):
    for hor in ("A", "B"):
00801
00802
00803
                                    detname = (
00804
00805
                                          + str(frame)
```

```
00806
00807
                                   + vert
00808
00809
                                   + hor
00810
                                   + "_EDGE"
00811
                                   + str(edge)
00812
00813
                               detdict[detname] = bitsrev[bitidx]
00814
              # remove faked detect
00815
              del detdict["W0_TOP_A_EDGE1"]
00816
              logging.info(self.loginfo + "Edge detect report:")
00817
00818
              for key, val in detdict.items():
                  logging.info(self.loginfo + key + ": " + val)
00820
```

#### 6.8.3.17 reportStatus()

Check contents of status register, print relevant messages

#### Definition at line 731 of file LLNL v4.py.

```
00731
          def reportStatus(self):
00732
00733
              Check contents of status register, print relevant messages
00734
00735
              statusbits = self.checkStatus()
00736
              statusbits2 = self.checkStatus2()
00737
00738
              logging.info(self.loginfo + "Status report:")
00739
              if int(statusbits[0]):
                  print(self.loginfo + "Sensor read complete")
00740
00741
              if int(statusbits[1]):
00742
                  print(self.loginfo + "Coarse trigger detected")
00743
              if int(statusbits[2]):
00744
                  print(self.loginfo + "Fine trigger detected")
00745
              if int(statusbits[5]):
00746
                  print(self.loginfo + "Sensor readout in progress")
00747
              if int(statusbits[6]):
00748
                  print(self.loginfo + "Sensor readout complete")
00749
              if int(statusbits[7]):
00750
                  print(self.loginfo + "SRAM readout started")
00751
              if int(statusbits[8]):
00752
                  print(self.loginfo + "SRAM readout complete")
              if int(statusbits[9]):
00753
00754
                  print(self.loginfo + "High-speed timing configuration started")
00755
              if int(statusbits[10]):
00756
                  print(self.loginfo + "All ADCs configured")
00757
              if int(statusbits[11]):
00758
                  print(self.loginfo + "All DACs configured")
00759
              if int(statusbits[13]):
00760
                  print(self.loginfo + "Timer has reset")
00761
              if int(statusbits[14]):
00762
                  print(self.loginfo + "Camera is Armed")
00763
              if int(statusbits[16]):
00764
                  print(self.loginfo + "High-speed timing configuration complete")
00765
              self.ca.sensor.reportStatusSensor(statusbits, statusbits2)
              temp = int(statusbits[23:16:-1], 2) * 3.3 * 1000 / 4096
00766
00767
              logging.info(
                  self.loginfo + "Temperature reading: " + "{0:1.2f}".format(temp) + " C"
00768
00769
00770
              press = int(statusbits[:23:-1], 2) * 3.3 * 1000 / 4096
00771
              logging.info(
00772
                  self.loginfo
                  + "Pressure sensor reading: '
+ "{0:1.2f}".format(press)
00773
00774
                  + " mV"
00775
00776
00777
              if int(statusbits2[0]):
                  print(self.loginfo + "FPA_IF_TO")
00778
```

```
00779
              if int(statusbits2[1]):
00780
                 print(self.loginfo + "SRAM_RO_TO")
00781
              if int(statusbits2[2]):
00782
                  print(self.loginfo + "PixelRd Timeout Error")
00783
              if int(statusbits2[3]):
00784
                  print(self.loginfo + "UART_TX_TO_RST")
00785
              if int(statusbits2[4]):
00786
                  print(self.loginfo + "UART_RX_TO_RST")
00787
```

#### 6.8.3.18 resetTimer()

#### Definition at line 560 of file LLNL v4.py.

```
00560 def resetTimer(self):
00561 """

00562 Reset on-board timer

00563 
00564 Returns:
00565 tuple (error string, response string from register set)

00566 """

00567 logging.info(self.loginfo + "resetTimer")
00568 control_messages = [("RESET_TIMER", "1"), ("RESET_TIMER", "0")]
00569 return self.ca.submitMessages(control_messages, " resetTimer: ")
```

#### 6.8.3.19 setLED()

```
00581 def setLED(self, LED, status):
00582 """
00583 Dummy function; feature is not implemented on LLNL_V4 board
00585 Returns:
00586 tuple: dummy of (error string, response string from setSubregister())
00587 """
00588 del LED, status
00589 return "", "0"
```

## 6.8.3.20 setPowerSave()

```
nsCamera.boards.LLNL_v4.llnl_v4.setPowerSave (
               self,
               status )
Select powersave option
    status: setting for powersave option (1 is enabled)
Returns:
    tuple (error string, response string from setSubregister())
Definition at line 591 of file LLNL_v4.py.
00591
         def setPowerSave(self, status):
00592
00593
              Select powersave option
00594
00595
00596
                status: setting for powersave option (1 is enabled)
00597
00598
             tuple (error string, response string from setSubregister())
00599
00600
00601
              if status:
                 status = 1
00603
              return self.ca.setSubregister("POWERSAVE", str(status))
00604
6.8.3.21 setPPER()
nsCamera.boards.LLNL_v4.llnl_v4.setPPER (
               self,
               pollperiod )
Set polling period for ADCs.
Args:
    pollperiod: milliseconds, between 1 and 255; defaults to 50
Returns:
    tuple (error string, response string from setSubregister() OR invalid time
      setting string)
Definition at line 605 of file LLNL v4.py.
00605
         def setPPER(self, pollperiod):
00606
00607
              Set polling period for ADCs.
00608
00609
                 pollperiod: milliseconds, between 1 and 255; defaults to 50
00610
00611
              Returns:
00612
                 tuple (error string, response string from setSubregister() OR invalid time
00613
                   setting string)
00614
00615
              if pollperiod is None:
00616
                 pollperiod = 50
              if not isinstance(pollperiod, int) or pollperiod < 1 or pollperiod > 255:
00617
00618
                 err = (
                     self.logerr + "invalid poll period submitted. Setting remains "
"unchanged."
00619
00620
00621
                 logging.error(err)
00622
                 return err, str(pollperiod)
00623
00624
              else:
00625
                 binset = bin(pollperiod)[2:].zfill(8)
                 return self.ca.setSubregister("PPER", binset)
00626
00627
```

## 6.8.3.22 softReboot()

```
nsCamera.boards.LLNL_v4.llnl_v4.softReboot (
                self )
Perform software reboot of board. WARNING: board reboot will likely prevent
  correct response and therefore will generate an error message
Returns:
    tuple (error string, response string) from final control message
Definition at line 428 of file LLNL_v4.py.
00428
         def softReboot(self):
00429
00430
              Perform software reboot of board. WARNING: board reboot will likely prevent
00431
               correct response and therefore will generate an error message
00432
00433
             tuple (error string, response string) from final control message
00434
00435
00436
             logging.info(self.loginfo + "reboot")
00437
              control_messages = [("RESET", "0")]
00438
              return self.ca.submitMessages(control_messages, " disarm: ")
00439
6.8.3.23 startCapture()
nsCamera.boards.LLNL_v4.llnl_v4.startCapture (
               self,
               mode = "Hardware" )
Selects trigger mode and enables board for image capture
    mode: trigger mode ("hardware"|"software"|"dual|"h"|"s"|"d" , is case-
      insensitive)
Returns:
    tuple (error string, response string) from final control message
Definition at line 458 of file LLNL v4.py.
         def startCapture(self, mode="Hardware"):
00459
00460
              Selects trigger mode and enables board for image capture
00461
00462
00463
                 mode: trigger mode ("hardware"|"software"|"dual|"h"|"s"|"d", is case-
00464
                   insensitive)
00465
00466
              Returns:
              tuple (error string, response string) from final control message _{\tt mum}
00467
00468
00469
              logging.info(self.loginfo + "startCapture")
00470
              if self.ca.sensmanual:
                 timingReg = "MANSHUT_MODE"
00471
00472
              else:
00473
                 timingReg = "HST_MODE"
00474
00475
              if mode.upper()[0] == "S": # SOFTWARE
00476
                 trigmess = [
                      ("HW_TRIG_EN", "0"),
("SW_TRIG_EN", "1"),
00477
00478
                      ("SW_TRIG_START", "1"),
00479
```

```
00480
              else: # HARDWARE
00481
00482
                  trigmess = [
00483
                      ("SW_TRIG_EN", "0"),
00484
                      ("HW_TRIG_EN", "1"),
00485
00486
00487
              control_messages = [
00488
                  ("ADC_CTL", "0000000F"), # configure all ADCs
00489
                  (timingReg, "1"),
00490
00491
00492
              control_messages.extend(trigmess)
00493
              return self.ca.submitMessages(control_messages, " startCapture: ")
00494
```

#### 6.8.3.24 waitForSRAM()

# Definition at line 506 of file LLNL\_v4.py.

```
00506
         def waitForSRAM(self, timeout):
00507
              Wait until subreg 'SRAM_READY' flag is true or timeout is exceeded;
00508
00509
                timeout = None or zero means wait indefinitely
00510
00511
00512
                 timeout - time in seconds before readoff proceeds automatically without
00513
                   waiting for SRAM_READY flag
00514
00515
              Returns:
              error string
00516
00517
00518
              logging.info(self.loginfo + "waitForSRAM, timeout = " + str(timeout))
00519
              waiting = True
00520
              starttime = time.time()
              err = ""
00521
00522
              while waiting:
00523
                 err, status = self.ca.getSubregister("SRAM_READY")
00524
00525
                      err = self.logerr + "error in register read: " + err + " (waitForSRAM)"
00526
                      logging.error(err)
00527
                  if int(status):
00528
                      waiting = False
00529
                      logging.info(self.loginfo + "SRAM ready")
00530
                  if self.ca.abort:
00531
                      waiting = False
00532
                      logging.info(self.loginfo + "readoff aborted by user")
00533
                      self.ca.abort = False
00534
                  if timeout and time.time() - starttime > timeout:
00535
                      err += self.logerr + "SRAM timeout; proceeding with download attempt"
00536
                      logging.error(err)
00537
                      return err
                  # Slow down for debugging (avoid thousands of messages)
00538
00539
                  if self.ca.verbose >= 5:
00540
                     time.sleep(0.5)
00541
              return err
00542
```

## 6.8.4 Member Data Documentation

#### 6.8.4.1 ADC5\_bipolar

```
nsCamera.boards.LLNL_v4.llnl_v4.ADC5_bipolar
```

Definition at line 193 of file LLNL\_v4.py.

## 6.8.4.2 ADC5\_mult

```
nsCamera.boards.LLNL_v4.llnl_v4.ADC5_mult
```

Definition at line 190 of file LLNL\_v4.py.

#### 6.8.4.3 ca

```
nsCamera.boards.LLNL_v4.llnl_v4.ca
```

Definition at line 182 of file LLNL\_v4.py.

#### 6.8.4.4 daedalus\_monitor\_controls

```
\verb|nsCamera.boards.LLNL_v4.llnl_v4.daedalus_monitor_controls|\\
```

Definition at line 291 of file LLNL\_v4.py.

## 6.8.4.5 daedalus\_subreg\_aliases

```
nsCamera.boards.LLNL_v4.llnl_v4.daedalus_subreg_aliases
```

Definition at line 264 of file LLNL\_v4.py.

## 6.8.4.6 defoff

```
\verb|nsCamera.boards.LLNL_v4.llnl_v4.defoff|
```

Definition at line 207 of file LLNL\_v4.py.

## 6.8.4.7 defsens

```
nsCamera.boards.LLNL_v4.llnl_v4.defsens
```

Definition at line 208 of file LLNL\_v4.py.

## 6.8.4.8 icarus\_monitor\_controls

```
nsCamera.boards.LLNL_v4.llnl_v4.icarus_monitor_controls
```

Definition at line 252 of file LLNL\_v4.py.

## 6.8.4.9 icarus\_subreg\_aliases

```
nsCamera.boards.LLNL_v4.llnl_v4.icarus_subreg_aliases
```

Definition at line 213 of file LLNL\_v4.py.

#### 6.8.4.10 logcrit

```
nsCamera.boards.LLNL_v4.llnl_v4.logcrit
```

Definition at line 183 of file LLNL\_v4.py.

## 6.8.4.11 logdebug

```
nsCamera.boards.LLNL_v4.llnl_v4.logdebug
```

Definition at line 187 of file LLNL\_v4.py.

#### 6.8.4.12 logerr

```
nsCamera.boards.LLNL_v4.llnl_v4.logerr
```

Definition at line 184 of file LLNL\_v4.py.

## 6.8.4.13 loginfo

```
nsCamera.boards.LLNL_v4.llnl_v4.loginfo
```

Definition at line 186 of file LLNL\_v4.py.

## 6.8.4.14 logwarn

nsCamera.boards.LLNL\_v4.llnl\_v4.logwarn

Definition at line 185 of file LLNL\_v4.py.

## 6.8.4.15 registers

```
nsCamera.boards.LLNL_v4.llnl_v4.registers [static]
```

Definition at line 38 of file LLNL\_v4.py.

#### 6.8.4.16 rs422\_baud

```
nsCamera.boards.LLNL_v4.llnl_v4.rs422_baud
```

Definition at line 194 of file LLNL\_v4.py.

#### 6.8.4.17 rs422 cmd wait

```
nsCamera.boards.LLNL_v4.llnl_v4.rs422_cmd_wait
```

Definition at line 195 of file LLNL\_v4.py.

## 6.8.4.18 subregisters

```
list nsCamera.boards.LLNL_v4.llnl_v4.subregisters [static]
```

Definition at line 101 of file LLNL v4.py.

#### 6.8.4.19 subreglist

```
{\tt nsCamera.boards.LLNL\_v4.llnl\_v4.subreglist}
```

Definition at line 300 of file LLNL\_v4.py.

## 6.8.4.20 VREF

```
nsCamera.boards.LLNL_v4.llnl_v4.VREF
```

Definition at line 189 of file LLNL\_v4.py.

The documentation for this class was generated from the following file:

C:/Users/hill35/git/camera\_python/nsCamera/boards/LLNL\_v4.py

## 6.9 nsCamera.utils.Packet.Packet Class Reference

#### **Public Member Functions**

- \_\_init\_\_ (self, preamble="aaaa", cmd="0", addr="", data="00000000", seqID="", payload\_length="", payload="", crc="")
- pktStr (self)
- calculateCRC (self)
- checkCRC (self)
- checkReadPacket (self, resppkt)
- checkResponsePacket (self, resppkt)
- checkResponseString (self, respstr)

#### **Public Attributes**

- PY3
- preamble
- cmd
- addr
- data
- seqID
- payload\_length
- payload
- crc
- type

## 6.9.1 Detailed Description

Definition at line 29 of file Packet.py.

### 6.9.2 Constructor & Destructor Documentation

### 6.9.2.1 \_\_init\_\_()

```
nsCamera.utils.Packet.Packet.__init__ (
                   self,
                   preamble = "aaaa",
                   cmd = "0",
                   addr = "",
                   data = "000000000",
                   seqID = "",
                   payload_length = "",
                   payload = "",
                   crc = "" )
Definition at line 52 of file Packet.py.
00063
                 self.PY3 = sys.version_info > (3,)
00065
                 self.preamble = preamble # 16 bit packet preamble
                 self.cmd = str(cmd) # 4 bit command packet
00066
                self.addr = addr.zfill(3) # 12 bit address packet
self.data = data.zfill(8) # 32 bit data packet
# 16 bit sequence ID packet (only Read Burst)
00067
00068
00069
00070
                 self.seqID = seqID
                # 16 bit payload packet (only Read Burst)
self.payload_length = payload_length
00071
00072
                # variable payload packet (only Read Burst) for now it's 16 bits
self.payload = payload
# 16 bit CRC-CCIT (XModem) packet
00073
00074
00075
00076
                self.crc = crc
                 self.type = ""
00077
                 if self.crc == "": # check if packet to be sent needs crc appended
00078
00079
                      self.crc = self.calculateCRC()
00080
```

#### 6.9.3 Member Function Documentation

nsCamera.utils.Packet.Packet.calculateCRC (

# 6.9.3.1 calculateCRC()

```
self )
Calculate CRC-CCIT (XModem) (2 bytes) from 8 byte packet for send and rcv
Returns:
   CRC as hexadecimal string without '0x'
Definition at line 108 of file Packet.py.
         def calculateCRC(self):
00108
00110
             Calculate CRC-CCIT (XModem) (2 bytes) from 8 byte packet for send and rcv
00111
00112
             CRC as hexadecimal string without '0x'
00113
00114
             preamble = self.preamble
00115
             crc = self.crc
00116
             self.crc = ""
00117
             self.preamble = ""
00118
00119
```

```
00120
              CRC_dec = crc16pure.crc16xmodem(str2bytes(self.pktStr()))
00121
              # input = int type decimal, output = hex string with 0x at the beginning
00122
              CRC_hex_0x = "0x%0.4X" % CRC_dec
              # make all hex letters lower case for comparison
00123
00124
              CRC_hex = CRC_hex_0x.lower()
00125
              # input = hex string with 0x at the beginning, output = hex str with 0x removed
00126
              CRC_hex = CRC_hex[2:]
00127
              self.preamble = preamble
00128
              self.crc = crc
00129
              return CRC_hex
00130
```

### 6.9.3.2 checkCRC()

### 6.9.3.3 checkReadPacket()

## Definition at line 137 of file Packet.py.

```
00137
           def checkReadPacket(self, resppkt):
00138
00139
               Confirm that Read Single occurred without error
00140
00141
                    resppkt: response packet
00142
00143
               Returns:
               tuple (error string, response packet as string)
00144
00145
00146
               err = ""
00147
               if int(resppkt.cmd.upper(), 16) - int(self.cmd.upper(), 16) != 0x8:
00148
                    err = "invalid command; "
               if resppkt.addr.upper() != self.addr.upper():
    err += "invalid address; "
00149
00150
               if resppkt.crc.upper() != resppkt.calculateCRC().upper():
    err += "invalid CRC; "
00151
00152
00153
               return err, resppkt.pktStr()
00154
```

### 6.9.3.4 checkResponsePacket()

```
nsCamera.utils.Packet.Packet.checkResponsePacket (
                resppkt )
Confirm that Write Single occurred without error
Aras:
    resppkt: response packet
Returns:
    tuple (error string, response packet as string)
Definition at line 155 of file Packet.py.
          def checkResponsePacket(self, resppkt):
00155
00156
00157
              Confirm that Write Single occurred without error
00158
              Args:
00159
                  resppkt: response packet
00160
00161
              Returns:
              tuple (error string, response packet as string)
00162
00163
              err = ""
00164
00165
              if int(resppkt.data, 16) & 1:
00166
                  err += "Checksum error; "
00167
              if int(resppkt.data, 16) & 2:
00168
              err += "Invalid command / command not executed; "
err1, rval = self.checkReadPacket(resppkt)
00169
00170
              err += err1
00171
              return err, rval
00172
```

#### 6.9.3.5 checkResponseString()

```
nsCamera.utils.Packet.Packet.checkResponseString (
               self,
               respstr )
Checks response string for error indicators
Args:
    respstr: packet as hexadecimal string
    tuple (error string, response packet string)
Definition at line 173 of file Packet.py.
         def checkResponseString(self, respstr):
00174
00175
             Checks response string for error indicators
00176
             Args:
00177
                respstr: packet as hexadecimal string
00178
00179
             Returns:
             tuple (error string, response packet string)
00180
00181
00182
             respstring = respstr.decode(encoding="UTF-8")
             resppkt = Packet(
00183
                 preamble=respstring[0:4],
00184
00185
                 cmd=respstring[4],
00186
                 addr=respstring[5:8],
00187
                 data=respstring[8:16],
```

)

00188

```
00189
 00190
                                           if resppkt.cmd == "8":
                                           # verify response to write command
 00191
 00192
                                                        err, rval = self.checkResponsePacket(resppkt)
 00193
                                          elif resppkt.cmd == "9":
 00194
                                                     err, rval = self.checkReadPacket(resppkt) # verify response to read command
 00195
                                           else:
 00196
                                            err = "Packet command invalid; "
 00197
                                                        rval = ""
 00198
                                           return err, rval
 00199
00200
00201 """
 00202 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
 00203 LLNL-CODE-838080
00205 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00206 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00207 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00208 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00209 be made under this license. 00210 """
6.9.3.6 pktStr()
nsCamera.utils.Packet.Packet.pktStr (
                                                 self )
Generate hexadecimal string form of packet
Returns:
           packet as hexadecimal string without '0x'
Definition at line 81 of file Packet.py.
                              def pktStr(self):
 00083
                                            Generate hexadecimal string form of packet
 00084
 00085
                                           packet as hexadecimal string without '0x' """ % \left( \frac{1}{2}\right) =\frac{1}{2}\left( \frac{1}{2}\right) =\frac{1}{2}\left
 00086
00087
 00088
                                            if self.seqID != "":
 00089
                                                        # Read burst response
00090
                                                        packetparts = [
00091
                                                                  self.preamble,
 00092
                                                                    self.cmd,
 00093
                                                                    self.seqID,
00094
                                                                    self.payload_length,
00095
                                                                    self.payload,
 00096
                                                                    self.crc,
 00097
                                                      1
 00098
                                            else:
                                                     # Single Command/Response response
00099
 00100
                                                        packetparts = [self.preamble, self.cmd, self.addr, self.data, self.crc]
 00101
                                            stringparts = [
                                                        part.decode("ascii") if isinstance(part, bytes) else part
00102
00103
                                                        for part in packetparts
00104
                                           out = "".join(stringparts)
0.0105
00106
                                           return out
```

# 6.9.4 Member Data Documentation

#### 6.9.4.1 addr

00107

nsCamera.utils.Packet.Packet.addr

Definition at line 67 of file Packet.py.

# 6.9.4.2 cmd

nsCamera.utils.Packet.Packet.cmd

Definition at line 66 of file Packet.py.

#### 6.9.4.3 crc

nsCamera.utils.Packet.Packet.crc

Definition at line 76 of file Packet.py.

# 6.9.4.4 data

nsCamera.utils.Packet.Packet.data

Definition at line 68 of file Packet.py.

# 6.9.4.5 payload

nsCamera.utils.Packet.Packet.payload

Definition at line 74 of file Packet.py.

### 6.9.4.6 payload\_length

 ${\tt nsCamera.utils.Packet.Packet.payload\_length}$ 

Definition at line 72 of file Packet.py.

# 6.9.4.7 preamble

nsCamera.utils.Packet.Packet.preamble

Definition at line 65 of file Packet.py.

### 6.9.4.8 PY3

nsCamera.utils.Packet.Packet.PY3

Definition at line 64 of file Packet.py.

## 6.9.4.9 seqID

```
nsCamera.utils.Packet.Packet.seqID
```

Definition at line 70 of file Packet.py.

## 6.9.4.10 type

```
nsCamera.utils.Packet.Packet.type
```

Definition at line 77 of file Packet.py.

The documentation for this class was generated from the following file:

C:/Users/hill35/git/camera\_python/nsCamera/utils/Packet.py

# 6.10 nsCamera.comms.RS422.RS422 Class Reference

#### **Public Member Functions**

- \_\_init\_\_ (self, camassem, baud=921600, par="O", stop=1)
- serialClose (self)
- sendCMD (self, pkt)
- arm (self, mode)
- readFrames (self, waitOnSRAM, timeout=0, fast=False, columns=1)
- readoff (self, waitOnSRAM, timeout, fast, columns=1)
- writeSerial (self, outstring, timeout)
- readSerial (self, size, timeout=None)
- closeDevice (self)

### **Public Attributes**

- ca
- logcrit
- logerr
- logwarn
- loginfo
- logdebug
- mode
- baud
- par
- stop
- read\_timeout
- write\_timeout
- datatimeout
- PY3
- skipError
- port
- · payloadsize

#### **Protected Attributes**

ser

### 6.10.1 Detailed Description

```
Code to manage RS422 connection. Will automatically query available COM interfaces until a board is found. Use the 'port=x' parameter in cameraAssembler call to specify a particular COM interface.

Exposed methods:

arm() - Puts camera into wait state for external trigger readFrames() - waits for data ready register flag, then copies camera image data into numpy arrays readoff() - waits for data ready register flag, then copies camera image data into numpy arrays; returns payload, payload size, and error message sendCMD(pkt) - sends packet object via serial port readSerial(size, timeout) - read 'size' bytes from serial port writeSerial(cmd) - submits string 'cmd' (assumes string is preformed packet) closeDevice() - close serial connections
```

Definition at line 30 of file RS422.py.

#### 6.10.2 Constructor & Destructor Documentation

#### 6.10.2.1 \_\_init\_\_()

```
nsCamera.comms.RS422.RS422.__init__ (
                self,
                camassem,
                baud = 921600,
                par = "0",
                stop = 1 )
Args:
    camassem: parent cameraAssembler object
    baud: bits per second
    par: parity type
    stop: number of stop bits
Definition at line 48 of file RS422.py.
          def __init__(self, camassem, baud=921600, par="0", stop=1):
00050
              Args:
00051
                   camassem: parent cameraAssembler object
00052
                  baud: bits per second
00053
                  par: parity type
                  stop: number of stop bits
00054
              ....
00055
00056
              self.ca = camassem
              self.logcrit = self.ca.logcritbase + "[RS422] "
00057
              self.logerr = self.ca.logerrbase + "[RS422]
00058
              self.logwarn = self.ca.logwarnbase + "[RS422] '
00059
              self.loginfo = self.ca.loginfobase + "[RS422] "
00060
              self.logdebug = self.ca.logdebugbase + "[RS422] "
logging.info(self.loginfo + "initializing RS422 comms object")
00061
00062
00063
              logging.debug(
```

```
00064
                   self.logdebug
00065
                   + "Init: baud = "
00066
                   + str(baud)
                   + "; par = "
00067
00068
                   + str(par)
00069
                   + "; stop = "
00070
                   + str(stop)
00071
00072
              self.mode = 0
00073
              self.baud = baud # Baud rate (bits/second)
00074
              self.par = par # Parity bit
00075
              self.stop = stop # Number of stop bits
00076
              self.read_timeout = 1 # default timeout for ordinary packets
00077
              self.write\_timeout = 1
00078
               # TODO: make datatimeout a cameraAssembler parameter
00079
              self.datatimeout = 60 # timeout for data read
00080
              logging.debug(
00081
                   self.logdebug + "Data timeout = " + str(self.datatimeout) + " seconds"
00082
00083
              self.PY3 = sys.version_info > (3,)
              self.skipError = False
port = ""
00084
00085
              ports = list(serial.tools.list_ports.comports())
logging.debug(self.logdebug + "Comports: " + str(ports))
00086
00087
00088
               for p, desc, add in ports:
                   if self.ca.port is None or p == "COM" + str(self.ca.port):
00089
                       logging.info(self.loginfo + "found comm port " + p)
00090
00091
                       try:
00092
                           with serial.Serial(
00093
00094
                               self.baud.
00095
                               parity=self.par,
00096
                               timeout=0.01.
00097
                               write timeout=0.01,
00098
                           ) as ser:
                               ser.write(str2bytes("aaaa1000000000001a84"))
00099
00100
                               time.sleep(1)
00101
                               s = ser.read(10)
                               resp = bytes2str(s)
00103
                               logging.debug(self.logdebug + "Init response: " + str(resp))
00104
                                    resp[0:5].lower() == "aaaa9"
00105
                               ): \mbox{\tt\#} TODO: add check for RS422 bit in board description
00106
00107
                                   boardid = resp[8:10]
00108
                                    if boardid == "00":
00109
                                        logging.critical(
                                            self.logcrit + "SNLrevC board detected - not "
00110
00111
                                            "compatible with nsCamera >= 2.0"
00112
00113
                                        sys.exit(1)
                                    elif boardid == "81":
00114
00115
                                        logging.info(self.loginfo + "LLNLv1 board detected")
00116
                                    elif boardid == "84":
                                        logging.info(self.loginfo + "LLNLv4 board detected")
00117
00118
00119
                                        logging.info(
00120
                                            self.loginfo + "unidentified board detected"
00121
00122
                                    logging.info(self.loginfo + "connected to " + p)
00123
                                    port = p
00124
                                    ser.reset_input_buffer()
00125
                                    ser.reset_output_buffer()
00126
00127
                       except Exception as e:
                          logging.error(self.logerr + "port identification: " + str(e))
00128
              if port == "":
00129
00130
                   if self.ca.port:
00131
                      logging.critical(
                           self.logcrit + "No usable board found at port " + str(self.ca.port)
00132
00133
00134
                       svs.exit(1)
00135
                   else:
00136
                      logging.critical(self.logcrit + "No usable board found")
00137
                       sys.exit(1)
              self.port = port # COM port to use for RS422 link
00138
00139
              self.ca.port = port[3:] # re-extract port number from com name
00140
00141
              self. ser = serial.Serial( # Class RS422
00142
                   port=self.port,
00143
                   baudrate=self.baud.
00144
                   parity=self.par,
```

```
00145
                  stopbits=self.stop,
00146
                  timeout=self.read_timeout, # timeout for serial read
00147
                  bytesize=serial.EIGHTBITS,
00148
00149
             self.payloadsize = (
00150
                  self.ca.sensor.width
00151
                  * self.ca.sensor.height
00152
                  * self.ca.sensor.nframes
00153
                  * self.ca.sensor.bytesperpixel
00154
00155
              logging.debug(
00156
                 self.logdebug + "Payload size: " + str(self.payloadsize) + " bytes"
00157
              self._ser.flushInput()
00159
              if not self._ser.is_open:
                  logging.critical(self.logcrit + "Unable to open serial connection")
00161
                  sys.exit(1)
00162
```

### 6.10.3 Member Function Documentation

#### 6.10.3.1 arm()

```
nsCamera.comms.RS422.RS422.arm (
               self,
               mode )
Puts camera into wait state for trigger. Mode determines source; defaults to
  'Hardware'
Args:
    mode:
             'Software' | 'S' activates software, disables hardware triggering
             'Hardware' | 'H' activates hardware, disables software triggering
               Hardware is the default
Returns:
    tuple (error, response string)
Definition at line 319 of file RS422.py.
         def arm(self, mode):
00320
00321
             Puts camera into wait state for trigger. Mode determines source; defaults to
00322
               'Hardware'
00323
00324
             Args:
                        'Software' | 'S' activates software, disables hardware triggering
00325
                 mode:
                         'Hardware' | 'H' activates hardware, disables software triggering
00326
                           Hardware is the default
00327
00328
00329
             Returns:
             tuple (error, response string)
00330
00331
00332
             if not mode:
                 mode = "Hardware"
00333
             logging.info(self.loginfo + "arm")
00334
00335
             logging.debug(self.logdebug + "arming mode: " + str(mode))
00336
             self.ca.clearStatus()
00337
             self.ca.latchPots()
00338
             err, resp = self.ca.startCapture(mode)
00339
             if err:
                 logging.error(self.logerr + "unable to arm camera")
00340
00341
             else:
00342
                 self.ca.armed = True
                 self.skipError = True
00343
00344
             return err, resp
00345
```

### 6.10.3.2 closeDevice()

nsCamera.comms.RS422.RS422.closeDevice (

```
self )
Close primary serial interface
Definition at line 496 of file RS422.pv.
00496
         def closeDevice(self):
00497
00498
             Close primary serial interface
00499
             logging.debug(self.logdebug + "Closing RS422 connection")
00500
00501
             self._ser.close()
00502
00503
00504 """
00505 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00506 LLNL-CODE-838080
00507
00508 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00509 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00510 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00511 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00512 be made under this license. 00513 """ \,
6.10.3.3 readFrames()
nsCamera.comms.RS422.RS422.readFrames (
               self.
               waitOnSRAM,
               timeout = 0,
               fast = False,
               columns = 1)
Copies image data from board into numpy arrays.
Args:
    waitOnSRAM: if True, wait until SRAM_READY flag is asserted to begin copying
    timeout: passed to waitForSRAM; after this many seconds begin copying data
      irrespective of SRAM_READY status; 'zero' means wait indefinitely
      WARNING: If acquisition fails, the SRAM will not contain a current image,
        but the code will copy the data anyway
    fast: if False, parse and convert frames to numpy arrays; if True, return
      unprocessed text stream
    columns: 1 for single image per frame, 2 for separate hemisphere images
Returns:
    list of numpy arrays OR raw text stream
Definition at line 346 of file RS422.py.
00346
         def readFrames(self, waitOnSRAM, timeout=0, fast=False, columns=1):
00347
00348
             Copies image data from board into numpy arrays.
00349
00350
             Args:
00351
                 waitOnSRAM: if True, wait until SRAM_READY flag is asserted to begin copying
00352
00353
                 timeout: passed to waitForSRAM; after this many seconds begin copying data
```

```
00354
                    irrespective of SRAM_READY status; 'zero' means wait indefinitely
00355
                    WARNING: If acquisition fails, the SRAM will not contain a current image,
00356
                     but the code will copy the data anyway
                  fast: if False, parse and convert frames to numpy arrays; if True, return
00357
00358
                   unprocessed text stream
00359
                  columns: 1 for single image per frame, 2 for separate hemisphere images
00360
00361
              Returns:
00362
                 list of numpy arrays OR raw text stream
00363
00364
00365
              frames, _, _ = self.readoff(waitOnSRAM, timeout, fast, columns)
00366
              return frames
```

#### 6.10.3.4 readoff()

```
nsCamera.comms.RS422.RS422.readoff (
              self.
              waitOnSRAM,
              timeout,
              fast,
              columns = 1)
Copies image data from board into numpy arrays; returns data, length of data,
and error messages. Use 'readFrames()' unless you require this additional
information
Args:
    waitOnSRAM: if True, wait until SRAM_READY flag is asserted to begin copying
      data
    timeout: passed to waitForSRAM; after this many seconds begin copying data
      irrespective of SRAM_READY status; 'zero' means wait indefinitely
      WARNING: If acquisition fails, the SRAM will not contain a current image,
        but the code will copy the data anyway
    fast: if False, parse and convert frames to numpy arrays; if True, return
      unprocessed text stream
    columns: 1 for single image per frame, 2 for separate hemisphere images
Returns:
   tuple (list of numpy arrays OR raw text stream, length of downloaded payload
      in bytes, payload error flag)
    NOTE: This reduces readoff by <1 second, so will have no noticeable impact
```

### Definition at line 368 of file RS422.py.

when using RS422

```
def readoff(self, waitOnSRAM, timeout, fast, columns=1):
00369
00370
              Copies image data from board into numpy arrays; returns data, length of data,
00371
              and error messages. Use 'readFrames()' unless you require this additional
00372
              information
00373
00374
              Args:
00375
                  waitOnSRAM: if True, wait until SRAM_READY flag is asserted to begin copying
00376
00377
                  timeout: passed to waitForSRAM; after this many seconds begin copying data
                    irrespective of SRAM_READY status; 'zero' means wait indefinitely
00378
                    WARNING: If acquisition fails, the SRAM will not contain a current image,
00379
00380
                     but the code will copy the data anyway
00381
                  fast: if False, parse and convert frames to numpy arrays; if True, return
00382
                   unprocessed text stream
00383
                  columns: 1 for single image per frame, 2 for separate hemisphere images
00384
00385
              Returns:
00386
                  tuple (list of numpy arrays OR raw text stream, length of downloaded payload
00387
                    in bytes, payload error flag)
```

```
00388
                   NOTE: This reduces readoff by <1 second, so will have no noticeable impact
00389
                     when using RS422
00390
               logging.info(self.loginfo + "readoff")
00391
00392
               logging.debug(
00393
                   self.logdebug
00394
                   + "readoff: waitonSRAM = "
00395
                   + str(waitOnSRAM)
00396
                   + "; timeout = "
00397
                   + str(timeout)
00398
                   + "; fast = '
00399
                   + str(fast)
00400
00401
               errortemp = False
00402
00403
               # Wait for data to be ready on board, turns off error messaging
               # Skip wait only if explicitly tagged 'False' ('None' defaults to True)
if waitOnSRAM is not False:
00404
00405
00406
                   logging.getLogger().setLevel(logging.CRITICAL)
00407
                   self.ca.waitForSRAM(timeout)
00408
                   logging.getLogger().setLevel(self.ca.verblevel)
00409
00410
               # Retrieve data
00411
               err, rval = self.ca.readSRAM()
00412
               if err:
                   logging.error(self.logerr + "Error detected in readSRAM")
00413
00414
               time.sleep(0.3)
               logging.debug(self.logdebug + "readoff: first 64 chars: " + str(rval[0:64]))
00415
               \ensuremath{\text{\#}} extract only the read burst data. Remove header & CRC footer
00416
00417
               read_burst_data = rval[36:-4]
00418
00419
               # Payload size as string implied by provided parameters
00420
               expectedlength = (
00421
                   4
                   * (self.ca.sensor.lastframe - self.ca.sensor.firstframe + 1)
* (self.ca.sensor.lastrow - self.ca.sensor.firstrow + 1)
00422
00423
00424
                   * self.ca.sensor.width
00425
00426
               padding = expectedlength - len(read_burst_data)
00427
               if padding:
00428
                   logging.warning(
00429
                        "{logwarn}readoff: Payload is shorter than expected."
                        " Padding with '0's".format(logwarn=self.logwarn)
00430
00431
00432
                   read_burst_data = read_burst_data.ljust(expectedlength, "0")
00433
00434
               if fast:
00435
                   return read_burst_data, len(read_burst_data) // 2, errortemp
00436
00437
                   parsed = generateFrames(self.ca, read_burst_data, columns)
00438
                   return parsed, len(read_burst_data) // 2, errortemp
00439
```

#### 6.10.3.5 readSerial()

Definition at line 465 of file RS422.py.

```
00465
          def readSerial(self, size, timeout=None):
00466
00467
              Read bytes from the serial port. Does not verify packets.
00468
00469
              Args:
00470
                 size: number of bytes to read
00471
                 timeout: serial timeout in sec
00472
00473
              tuple (error string, string read from serial port) _{\tt mum}
00474
00475
00476
              logging.debug(
00477
                  self.logdebug
00478
                  + "readSerial: size = "
00479
                  + str(size)
00480
                  + "; timeout = "
00481
                  + str(timeout)
00482
              err = ""
00483
00484
              if timeout:
00485
                  self._ser.timeout = timeout
00486
              else:
00487
                  self. ser.timeout = self.read timeout
00488
              resp = self._ser.read(size)
00489
              if len(resp) < 10: # bytes</pre>
00490
                  err += (
00491
                      self.logerr + "readSerial : packet too small: '" + bytes2str(resp) + "'"
00492
00493
                  logging.error(err)
00494
              return err, bytes2str(resp)
00495
```

### 6.10.3.6 sendCMD()

# Definition at line 170 of file RS422.py.

```
00170
         def sendCMD(self, pkt):
00171
00172
              Submit packet and verify response packet. Recognizes readoff packet and adjusts
00173
              read size and timeout appropriately
00174
00175
              Args:
                 pkt: Packet object
00176
00177
00178
              Returns:
              tuple (error, response string)
00179
00180
00181
              pktStr = pkt.pktStr()
00182
              logging.debug(self.logdebug + "sendCMD packet: " + str(pktStr))
00183
              self._ser.flushInput()
00184
              time.sleep(0.01) # wait 10 ms in between flushing input and output buffers
00185
              self. ser.flushOutput()
00186
              self.ca.writeSerial(pktStr)
              err0 = ""
00187
              err = ""
00188
              resp = ""
00189
              tries = 3 # TODO: make a function parameter?
00190
00191
```

```
00192
              if (
                  hasattr(self.ca, "board")
and pktStr[4] == "0"
00193
00194
00195
                  and pktStr[5:8] == self.ca.board.registers["SRAM_CTL"]
00196
                   # download data payload
00197
00198
00199
                      self.loginfo + "Payload size (bytes) = " + str(self.payloadsize)
00200
00201
                  crcresp0 = ""
00202
                  crcresp1 = ""
                  smallresp = ""
00203
00204
                  emptyResponse = False
00205
                  wrongSize = False
00206
                   # TODO: refactor payload error management to another method
00207
                  for i in range(tries):
00208
                      err, resp = self.readSerial(
00209
                          self.payloadsize + 20, timeout=self.datatimeout
00210
00211
                       if err:
00212
                          logging.error(
00213
                               self.logerr + "sendCMD: read payload failed " + pktStr + err
00214
00215
                          self.ca.payloaderror = True
00216
                       else:
00217
                           if not len(resp):
00218
                               err0 = self.logerr + "sendCMD: empty response from board"
00219
                               logging.error(err0)
00220
                               emptyResponse = True
                               self.ca.payloaderror = True
00221
00222
                           elif len(resp) != 2 * (self.payloadsize + 20):
00223
                               err0 = (
                                   self.logerr
00224
                                   + "sendCMD: incorrect response; expected "
00225
                                   + str(self.payloadsize + 20)
00226
                                   + " bytes, received "
00227
                                   + str(len(resp) // 2)
00228
00229
00230
                               logging.error(err0)
00231
                               wrongSize = True
                               smallresp = resp
00232
                               self.ca.payloaderror = True
00233
                           elif not checkCRC(resp[4:20]):
00234
00235
                               err0 = (
00236
                                   self.logerr
00237
                                   + "sendCMD: '
                                   + pktStr
+ " - payload preface CRC fail"
00238
00239
00240
00241
                               logging.error(err0)
00242
                               self.ca.payloaderror = True
00243
                               crcresp1 = resp
00244
                           elif not checkCRC(resp[24:]):
00245
                               err0 = (
00246
                                   self.logerr + "sendCMD: " + pktStr + " - payload CRC fail"
00247
00248
                               logging.error(err0)
00249
                               self.ca.payloaderror = True
00250
                               crcresp0 = resp
00251
                          err += err0
00252
                      time.sleep(5)
00253
                       if self.ca.payloaderror:
00254
                           # keep best results over multiple tries; e.g., if first try is
00255
                           # bad CRC and second try is an incomplete payload, use the
00256
                               first payload
00257
                           if i == tries - 1:
00258
                               if crcresp0:
00259
                                   logging.error(
                                       self.logerr + "sendCMD: Unable to acquire "
00260
00261
                                        "CRC-confirmed payload after
00262
                                       + str(tries)
00263
                                       + " attempts. Continuing with unconfirmed payload"
00264
                                   )
00265
                                   resp = crcresp0
00266
                               elif crcresp1:
00267
                                   logging.error(
00268
                                       self.logerr + "sendCMD: Unable to acquire "
00269
                                        "CRC-confirmed readoff after
00270
                                       + str(tries)
00271
                                       + " attempts. Continuing with unconfirmed payload"
00272
```

```
00273
                                   resp = crcresp1
00274
                               elif wrongSize:
00275
                                   logging.error(
00276
                                       self.logerr + "sendCMD: Unable to acquire complete "
00277
                                        "payload after "
00278
                                        + str(tries)
00279
                                       + " attempts. Dumping datastream to file."
00280
00281
                                   resp = smallresp
00282
                                   self.ca.dumpNumpy(resp)
00283
                               elif emptyResponse:
00284
                                   logging.error(
                                       self.logerr + "sendCMD: Unable to acquire any "
00285
                                        "payload after " + str(tries) + " attempts.
00287
00288
                           else:
00289
                               logging.info(
00290
                                   self.loginfo + "Retrying download, attempt #" + str(i + 1)
00291
00292
                               err = ""
                               err0 = ""
00293
00294
                               self.ca.payloaderror = False
00295
                               self.ca.writeSerial(pktStr)
00296
                       else:
                           logging.info(self.loginfo + "Download successful")
if self.ca.boardname == "llnl_v4":
00297
00298
00299
                               # self.ca.setSubregister('SWACK','1')
00300
00301
                           break
00302
00303
              else:
00304
                  # non-payload messages and workaround for initial setup before board object
00305
                     has been initialized
00306
                  time.sleep(0.03)
00307
                  self.\_ser.timeout = 0.02
00308
                  err, resp = self.readSerial(10)
00309
                  logging.debug(self.logdebug + "sendCMD response: " + str(resp))
00310
                  if err:
00311
                      logging.error(
                           self.logerr + "sendCMD: readSerial failed (regular packet) " + err
00312
00313
00314
                  elif not checkCRC(resp[4:20]):
                       err = self.logerr + "sendCMD- regular packet CRC fail: " + resp
00315
00316
                      logging.error(err)
00317
              return err, resp
00318
```

### 6.10.3.7 serialClose()

```
nsCamera.comms.RS422.RS422.serialClose ( self )
```

#### Definition at line 163 of file RS422.py.

```
00163 def serialClose(self):
00164 """
00165 Close serial interface
00166 """
00167 logging.debug(self.logdebug + "serialclose")
00168 self._ser.close() # close serial interface COM port
00169
```

# 6.10.3.8 writeSerial()

```
nsCamera.comms.RS422.RS422.writeSerial (
self,
outstring,
timeout )
```

```
Transmit string to board
Args:
    outstring: string to write
    timeout: serial timeout in sec
Returns:
    integer length of string written to serial port
Definition at line 440 of file RS422.py.
         def writeSerial(self, outstring, timeout):
00441
00442
              Transmit string to board
00443
00444
00445
                 outstring: string to write
00446
                 timeout: serial timeout in sec
00447
             Returns:
             . integer length of string written to serial port _{\tt mum}
00448
00449
00450
             logging.debug(
00451
                 self.logdebug
00452
                  + "writeSerial: outstring = "
00453
                 + str(outstring)
00454
                 + "; timeout = "
00455
                 + str(timeout)
00456
00457
             if timeout:
00458
                 self._ser.timeout = timeout
00459
00460
                 self._ser.timeout = self.write_timeout
00461
             lengthwritten = self._ser.write(str2bytes(outstring))
             self._ser.timeout = self.read_timeout # reset if changed above
00463
             return lengthwritten
00464
```

#### 6.10.4 Member Data Documentation

### 6.10.4.1 \_ser

```
nsCamera.comms.RS422.RS422._ser [protected]
```

Definition at line 141 of file RS422.py.

### 6.10.4.2 baud

```
nsCamera.comms.RS422.RS422.baud
```

Definition at line 73 of file RS422.py.

### 6.10.4.3 ca

```
nsCamera.comms.RS422.RS422.ca
```

Definition at line 56 of file RS422.py.

# 6.10.4.4 datatimeout

nsCamera.comms.RS422.RS422.datatimeout

Definition at line 79 of file RS422.py.

## 6.10.4.5 logcrit

nsCamera.comms.RS422.RS422.logcrit

Definition at line 57 of file RS422.py.

# 6.10.4.6 logdebug

nsCamera.comms.RS422.RS422.logdebug

Definition at line 61 of file RS422.py.

# 6.10.4.7 logerr

nsCamera.comms.RS422.RS422.logerr

Definition at line 58 of file RS422.py.

## 6.10.4.8 loginfo

nsCamera.comms.RS422.RS422.loginfo

Definition at line 60 of file RS422.py.

# 6.10.4.9 logwarn

nsCamera.comms.RS422.RS422.logwarn

Definition at line 59 of file RS422.py.

# 6.10.4.10 mode

nsCamera.comms.RS422.RS422.mode

Definition at line 72 of file RS422.py.

## 6.10.4.11 par

```
nsCamera.comms.RS422.RS422.par
```

Definition at line 74 of file RS422.py.

## 6.10.4.12 payloadsize

```
nsCamera.comms.RS422.RS422.payloadsize
```

Definition at line 149 of file RS422.py.

### 6.10.4.13 port

```
nsCamera.comms.RS422.RS422.port
```

Definition at line 138 of file RS422.py.

### 6.10.4.14 PY3

nsCamera.comms.RS422.RS422.PY3

Definition at line 83 of file RS422.py.

## 6.10.4.15 read\_timeout

```
nsCamera.comms.RS422.RS422.read_timeout
```

Definition at line 76 of file RS422.py.

# 6.10.4.16 skipError

nsCamera.comms.RS422.RS422.skipError

Definition at line 84 of file RS422.py.

# 6.10.4.17 stop

nsCamera.comms.RS422.RS422.stop

Definition at line 75 of file RS422.py.

# 6.10.4.18 write\_timeout

nsCamera.comms.RS422.RS422.write\_timeout


Definition at line 77 of file RS422.py.

The documentation for this class was generated from the following file:

C:/Users/hill35/git/camera\_python/nsCamera/comms/RS422.py

# 6.11 nsCamera.sensors.sensorBase.sensorBase Class Reference

Inheritance diagram for nsCamera.sensors.sensorBase.sensorBase:



#### **Public Member Functions**

- \_\_init\_\_ (self, camassem)
- init board specific (self)
- checkSensorVoltStat (self)
- · setInterlacing (self, ifactor)
- setHighFullWell (self, flag)
- setZeroDeadTime (self, flag)
- setTriggerDelay (self, delay)
- setPhiDelay (self, delay)
- setExtClk (self, delay)
- setTiming (self, side="AB", sequence=None, delay=0)
- setArbTiming (self, side="AB", sequence=None)
- getTiming (self, side, actual)
- setManualShutters (self, timing=None)
- setManualTiming (self, timing=None)
- getManualTiming (self)
- getSensTemp (self, scale=None, offset=None, slope=None, dec=None)
- selectOscillator (self, osc=None)
- parseReadoff (self, frames, columns)
- getSensorStatus (self)
- reportStatusSensor (self, statusbits, statusbits2)

### **Public Attributes**

- ca
- logcrit
- logerr
- logwarn
- loginfo
- logdebug
- fpganumID
- · sensfam

# 6.11.1 Detailed Description

Base class for sensors.  $^{\prime}$  Virtual  $^{\prime}$  methods below default to Icarus behavior. daedalus.py overrides these methods as necessary

Definition at line 25 of file sensorBase.py.

## 6.11.2 Constructor & Destructor Documentation

### 6.11.2.1 \_\_init\_\_()

Reimplemented in nsCamera.sensors.idaedalus.daed

### Definition at line 31 of file sensorBase.py.

```
00031
           def __init__(self, camassem):
00032
                self.ca = camassem
               # skip board settings if no board object exists
if hasattr(self.ca, "board"):
00033
00034
                    self.init_board_specific()
00035
00036
00037
                    self.logcrit,
00038
00039
                    self.logerr,
00040
                    self.logwarn,
00041
                    self.loginfo,
00042
                    self.logdebug,
00043
               ) = makeLogLabels(self.ca.logtag, self.loglabel)
00044
               # skip assignment if no comms object exists
if hasattr(self.ca, "comms"):
00045
00046
00047
                    self.ca.comms.payloadsize = (
00048
                         self.width * self.height * self.nframes * self.bytesperpixel
00049
00050
                logging.info(self.loginfo + "Initializing sensor object")
00051
00052
```

### 6.11.3 Member Function Documentation

#### 6.11.3.1 checkSensorVoltStat()

Reimplemented in nsCamera.sensors.icarus.icarus.

## Definition at line 64 of file sensorBase.py.

```
00064
          def checkSensorVoltStat(self):
00065
00066
              Checks register tied to sensor select jumpers to confirm match with sensor
00067
00068
00069
              Returns:
              boolean, True if jumpers select for Icarus sensor _{\pi\pi\pi}
00070
00071
00072
              logging.debug(self.logdebug + "checkSensorVoltStat")
00073
              err, status = self.ca.getSubregister(self.detect)
00074
00075
                  logging.error(self.logerr + "Unable to confirm sensor status")
00076
00077
              if not int(status):
00078
                  logging.error(self.logerr + self.sensfam + " sensor not detected")
00079
00080
              return True
00081
```

# 6.11.3.2 getManualTiming()

```
nsCamera.sensors.sensorBase.sensorBase.getManualTiming ( self \ ) Read off manual shutter timing settings Overridden in daedalus.py Returns:  list \ of \ 2 \ lists \ of \ timing \ from \ A \ and \ B \ sides, \ respectively
```

Reimplemented in nsCamera.sensors.daedalus.daedalus.

### Definition at line 611 of file sensorBase.py.

```
00611
          def getManualTiming(self):
00612
00613
               Read off manual shutter timing settings
00614
               Overridden in daedalus.py
00615
               Returns:
               list of 2 lists of timing from A and B sides, respectively """ \ensuremath{\text{"""}}
00616
00617
               aside = []
00618
               bside = []
00619
00620
               for reg in [
                    "WO_INTEGRATION",
00621
                    "WO_INTERFRAME",
00622
                   "W1_INTEGRATION",
00623
```

```
00624
                  "W1_INTERFRAME",
00625
                  "W2_INTEGRATION",
00626
                  "W2_INTERFRAME",
00627
                  "W3_INTEGRATION",
00628
00629
                  _, reghex = self.ca.getRegister(reg)
00630
                  aside.append(25 * int(reghex, 16))
00631
              for reg in [
00632
                  "WO_INTEGRATION_B",
00633
                  "WO_INTERFRAME_B"
00634
                  "W1_INTEGRATION_B",
00635
                  "W1_INTERFRAME_B",
00636
                  "W2_INTEGRATION_B",
00637
                  "W2_INTERFRAME_B",
00638
                  "W3_INTEGRATION_B",
00639
            ]:
00640
                  _, reghex = self.ca.getRegister(reg)
                  bside.append(25 * int(reghex, 16))
00641
00642
              return [aside, bside]
00643
```

### 6.11.3.3 getSensorStatus()

```
nsCamera.sensors.sensorBase.sensorBase.getSensorStatus ( self \ )
```

Wrapper for reportSensorStatus so that the user doesn't have to query statusbits

### Definition at line 691 of file sensorBase.py.

```
00691  def getSensorStatus(self):
00692     """
00693     Wrapper for reportSensorStatus so that the user doesn't have to query statusbits
00694     """
00695     sb1 = self.ca.board.checkstatus()
00696     sb2 = self.ca.board.checkstatus2()
00697     self.reportStatusSensor(sb1, sb2)
```

### 6.11.3.4 getSensTemp()

Virtual method (Temperature sensor is not present on Icarus sensors). Returns 0. Overridden by Daedalus method

Reimplemented in nsCamera.sensors.daedalus.daedalus.

#### Definition at line 644 of file sensorBase.py.

```
00644 def getSensTemp(self, scale=None, offset=None, slope=None, dec=None):
00645 """

00646 Virtual method (Temperature sensor is not present on Icarus sensors). Returns 0.
00647 Overridden by Daedalus method
00648 """

00649 return 0
00650
```

nsCamera.sensors.sensorBase.sensorBase.getTiming (

## 6.11.3.5 getTiming()

```
self,
               side,
               actual )
actual = True: returns actual high speed intervals that will be generated by the
             FPGA as list
          False: Returns high speed timing settings as set by setTiming. Assumes
             that timing was set via the setTiming method--it will not accurately
              report arbitrary timings set by direct register sets or manual
              shutter control
Args:
     side: Hemisphere 'A' or 'B'
    actual: False: return HST settings
             True: calculate and return actual HST behavior
Returns:
    actual= True: list of shutter intervals;
                  icarus: [delay, open0, closed0, open1, closed1, open2, closed2,
                             open3]
                  daedalus: [delay, open0, closed0, open1, closed1, open2]
             False: tuple (hemisphere label,
                               'open shutter' in ns,
                                'closed shutter' in ns,
                                initial delay in ns)
Definition at line 434 of file sensorBase.py.
00434
          def getTiming(self, side, actual):
    """
00435
00436
              actual = True: returns actual high speed intervals that will be generated by the
00437
                          FPGA as list
00438
                       False: Returns high speed timing settings as set by setTiming. Assumes
00439
                          that timing was set via the setTiming method--it will not accurately
00440
                          report arbitrary timings set by direct register sets or manual
00441
                          shutter control
00442
00443
00444
                 side: Hemisphere 'A' or 'B'
00445
                  actual: False: return HST settings
00446
                         True: calculate and return actual HST behavior
00447
00448
00449
                 actual= True: list of shutter intervals;
00450
                              icarus: [delay, open0, closed0, open1, closed1, open2, closed2,
00451
                                       open31
00452
                              daedalus: [delay, open0, closed0, open1, closed1, open2]
00453
                          False: tuple (hemisphere label,
00454
                                          'open shutter' in ns,
00455
                                          'closed shutter' in ns,
00456
                                          initial delay in ns)
00457
00458
00459
              logging.info("{}getTiming".format(self.loginfo))
00460
              if side is None:
00461
                 side = "A"
00462
              logging.info(self.loginfo + "get timing, side " + side.upper())
00463
              if side.upper() == "A":
                 lowreg = "HS_TIMING_DATA_ALO"
00464
             highreg = "HS_TIMING_DATA_AHI"
elif side.upper() == "B":
00465
00466
                  lowreg = "HS_TIMING_DATA_BLO"
00467
                 highreg = "HS_TIMING_DATA_BHI"
00468
00469
              else:
00470
                 logging.error(
00471
                      self.logerr
+ "Invalid sensor side: "
00472
00473
                      + side
```

```
00474
                      + "; timing settings unchanged"
00475
00476
                  return "", 0, 0, 0
              err, lowpart = self.ca.getRegister(lowreg)
00477
00478
              err1, highpart = self.ca.getRegister(highreg)
00479
              if err or err1:
00480
                  logging.error(
00481
                      self.logerr + "Unable to retrieve timing setting (getTiming), "
00482
                       "returning zeroes "
00483
00484
                   return side.upper(), 0, 0, 0
              full40hex = highpart[-2:] + lowpart.zfill(8)
00485
00486
              full40bin = "{0:0=40b}".format(int(full40hex, 16))
              logging.debug(self.logdebug + "full40bin = " + str(full40bin))
00488
              if actual:
00489
                  if full40bin == "0" * 40: # all-zero timing
00490
                      if self.fpganumID == 2:
00491
                          times = [0] * 6
00492
                      else:
00493
                         times = [0] * 8
00494
                  else:
00495
                      full160 = 4 * full40bin
00496
                      gblist = [[k, len(list(g))] for k, g in itertools.groupby(full160)]
00497
                       if self.fpganumID == 2:
00498
                          times = [int(x[1]) for x in gblist[:-7:-1]]
00499
                      else:
                          times = [int(x[1]) for x in gblist[:-9:-1]]
00500
                      times[0] = times[0] - 1
00501
                  if self.ca.sensorname == "icarus":
# get timing for frames 1 and 2, keep delay as offset
00502
00503
                      # TODO: should this give a 'correct' offset from frame 0?
00504
                      times12 = [times[0]] + times[3:6]
00505
00506
                      return times12
00507
                  return times
00508
              else:
                  if full40bin == "0" * 40: # all-zero timing
00509
00510
                      timeon, timeoff, delay = (0, 0, 0)
00511
                  else:
00512
                      gblist = [[k, len(list(g))] for k, g in itertools.groupby(full40bin)]
00513
                      delay = gblist[-1][1] - 1
                      timeon = gblist[-2][1]
00514
00515
00516
                      if self.ca.sensorname == "icarus":
                           if len(gblist) == 2: # 39,1 corner case
00517
00518
                              timeoff = 1
00519
                           elif len(gblist) == 3: # sequence fits only once
00520
                              timeoff = 40 - timeon
00521
                           else:
00522
                               timeoff = gblist[-3][1]
00523
00524
                           if len(gblist) < self.nframes: # sequence fits only once</pre>
00525
                               timeoff = 40 - timeon
00526
00527
                               # TODO: confirm '-3' works for daedalus
00528
                               timeoff = gblist[-3][1]
00529
                  return side.upper(), timeon, timeoff, delay
00530
```

### 6.11.3.6 init\_board\_specific()

```
nsCamera.sensors.sensorBase.sensorBase.init_board_specific ( self \ )
```

Initialize aliases and subregisters specific to the current board and sensor.

# Definition at line 53 of file sensorBase.py.

```
00053
         def init board specific(self):
00054
              """Initialize aliases and subregisters specific to the current board and sensor."""
00055
00056
              if self.ca.sensorname == "icarus" or self.ca.sensorname == "icarus2":
                  self.ca.board.subreq aliases = self.ca.board.icarus subreq aliases
00057
00058
                  self.ca.board.monitor_controls = self.ca.board.icarus_monitor_controls
00059
              else:
00060
                  self.ca.board.subreg_aliases = self.ca.board.daedalus_subreg_aliases
00061
                  self.ca.board.monitor_controls = self.ca.board.daedalus_monitor_controls
00062
```

## 6.11.3.7 parseReadoff()

Reimplemented in nsCamera.sensors.daedalus.daedalus.

#### Definition at line 683 of file sensorBase.py.

```
00683 def parseReadoff(self, frames, columns):
00684 """
00685 Virtual method (Order parsing is unnecessary for Icarus, continue to hemisphere
00686 parsing.)
00687 Overridden by Daedalus method
00688 """
00689 return self.ca.partition(frames, columns)
```

#### 6.11.3.8 reportStatusSensor()

Reimplemented in nsCamera.sensors.daedalus.daedalus.

# Definition at line 699 of file sensorBase.py.

```
00699
          def reportStatusSensor(self, statusbits, statusbits2):
00700
00701
              Print status messages from sensor-specific bits of status register, default for
00702
                Icarus family sensors
00703
00704
                  statusbits: result of checkStatus()
00705
                  statusbits2: result of checkStatus2()
00706
00707
              if int(statusbits[3]):
00708
                  print(self.loginfo + "W3_Top_A_Edge1 detected")
00709
              if int(statusbits[4]):
00710
                  print(self.loginfo + "W3_Top_B_Edge1 detected")
00711
              if int(statusbits[12]):
              print(self.loginfo + "HST_All_W_En detected")
if self.ca.boardname == "llnl_v4" and int(statusbits2[5]):
00712
00713
                  print(self.loginfo + "PDBIAS Unready")
00714
00715
00716
00717 # TODO: add function to control TIME_ROW_DCD delay
00718
00719 """
00720 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00721 LLNL-CODE-838080
00722
00723 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00724 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00725 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00726 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00727 be made under this license.
00728 """
```

### 6.11.3.9 selectOscillator()

Reimplemented in nsCamera.sensors.daedalus.daedalus.

```
Definition at line 651 of file sensorBase.py.
```

```
00651
          def selectOscillator(self, osc=None):
00652
00653
              Selects oscillator to control sensor timing
00654
              Overridden in daedalus.py
00655
00656
                  osc: 'relaxation'|'ring'|'ringnoosc'|'external', defaults to relaxation
00657
00658
              Returns:
              error message as string
00659
00660
00661
              logging.info(self.loginfo + "selectOscillator; osc = " + str(osc))
00662
              if osc is None:
                 osc = "rel"
00663
00664
              osc = str(osc)
00665
              if osc.upper()[:3] == "REL":
00666
                  payload = "00"
              elif osc.upper()[:3] == "RIN":
00667
                  if "NO" in osc.upper() or "0" in osc:
    payload = "10"
00668
00669
00670
                  else:
00671
                     payload = "01"
00672
              elif osc.lower()[:3] in ["ext", "phi"]:
                  payload = "11"
00673
00674
              else:
00675
                  err = (
00676
                       self.logerr + "selectOscillator: invalid parameter supplied. "
00677
                       "Oscillator selection is unchanged."
00678
00679
                  logging.error(err)
00680
00681
              self.ca.setSubregister("OSC_SELECT", payload)
00682
```

#### 6.11.3.10 setArbTiming()

```
for interlacing or ZDT, you should populate the entire 40-bit register,
       e.g., [0,3,2,3,2,3,2,3,2,3,2,3,2,3,2,3,2]
     \star \text{NOTE} \star \text{Icarus} sensors generally cannot use 1 ns timing, so should use at
       least 2 ns for frames 2 and 3 integration and interframe times (an initial
       delay of only 1 ns is acceptable)
     *NOTE* although the Icarus model 1 only images the middle two frames, timing
       entries must be provided for all four frames; to implement frame 1 open \,
       for X ns, shutter closed for Y ns, and frame 2 open for Z ns, use the
       sequence [0,1,1,X,Y,Z,1,1]
     *WARNING* arbitrary timings will not be restored after a board power cycle
Returns:
    list: Actual timing results
Definition at line 308 of file sensorBase.py.
00308
          def setArbTiming(self, side="AB", sequence=None):
00309
00310
              Set arbitrary high-speed timing sequence.
00311
              Args:
00312
                  side: Hemisphere 'A', 'B', 'AB'
00313
                  sequence: list of arbitrary timing intervals, beginning with initial delay.
00314
                     The conventional timing (3,2) with delay = 0 would be represented by
                     [0,3,2,3,2,3,2,3] on icarus devices, [0,3,2,3,2,3] on daedalus. If used
00315
                     for interlacing or ZDT, you should populate the entire 40-bit register,
00316
00317
                    e.g., [0,3,2,3,2,3,2,3,2,3,2,3,2,3,2,3,2]
00318
                  *NOTE* Icarus sensors generally cannot use 1 ns timing, so should use at least 2 ns for frames 2 and 3 integration and interframe times (an initial
00319
00320
00321
                    delay of only 1 ns is acceptable)
00322
00323
                  *NOTE* although the Icarus model 1 only images the middle two frames, timing
                    entries must be provided for all four frames; to implement frame 1 open
00324
                    for X ns, shutter closed for Y ns, and frame 2 open for Z ns, use the
00325
00326
                    sequence [0,1,1,X,Y,Z,1,1]
00327
00328
                  *WARNING* arbitrary timings will not be restored after a board power cycle
00329
00330
              Returns:
              list: Actual timing results
00331
00332
00333
              logging.info(
00334
                   "{}setArbTiming; side={}, sequence={}".format(self.loginfo, side, sequence)
00335
00336
              if sequence is None:
00337
                  if self.sensfam == "Daedalus":
00338
                      sequence = [0, 2, 3, 4, 5, 6]
00339
                      sequence = [0, 2, 3, 4, 5, 6, 7, 8]
00340
00341
              logging.info(
                  self.loginfo + "HST side " + side.upper() + " (arbitrary): " + str(sequence)
00342
00343
              if side.upper() == "AB":
00344
                  err1, _ = self.setArbTiming(side="A", sequence=sequence)
00345
                  err2, actual = self.setArbTiming(side="B", sequence=sequence)
00346
              return err1 + err2, actual
if side.upper() == "A":
00347
00348
                  lowreg = "HS_TIMING_DATA_ALO"
highreg = "HS_TIMING_DATA_AHI"
00349
00350
00351
              elif side.upper() == "B":
00352
                  lowreg = "HS_TIMING_DATA_BLO"
00353
                  highreg = "HS_TIMING_DATA_BHI"
00354
              else:
00355
                  err = (
00356
                       self.logerr
00357
                       + "Invalid sensor side: "
00358
00359
                       + "; timing settings unchanged"
00360
00361
                  logging.error("{}setArbTiming: {}".format(self.logerr, err))
00362
                  return err, "0000000000"
00363
00364
              full40 = [0] * 40
```

bitlist = []

00365

```
00366
              flag = 0 # similar to setTiming, but starts with delay
00367
00368
              for a in sequence:
00369
                  add = [flag] * a
00370
                  bitlist += add
00371
                  if flag:
00372
                      flag = 0
00373
00374
                      flag = 1
00375
00376
              logging.debug(self.logdebug + "bitlist = " + str(bitlist))
00377
              reversedlist = bitlist[39::-1]
00378
              full40[-(len(reversedlist) + 1) : -1] = reversedlist
00379
              full40bin = "".join(str(x) for x in full40)
00380
              logging.debug(self.logdebug + "full40bin = " + str(full40bin))
00381
              full40hex = "%x" % int(full40bin, 2)
              logging.debug(self.logdebug + "full40hex = " + str(full40hex))
00382
00383
              highpart = full40hex[-10:-8].zfill(8)
              lowpart = full40hex[-8:].zfill(8)
00384
00385
              self.ca.setRegister(lowreg, lowpart)
              # deactivates manual shutter mode if previously engaged
00386
00387
              self.ca.setSubregister("MANSHUT_MODE", "0")
self.ca.setSubregister("HST_MODE", "1")
00388
00389
00390
              actual = self.getTiming(side, actual=True)
00391
              f0delay = sequence[1] + sequence[2]
00392
00393
              if self.ca.sensorname == "icarus":
00394
                  if actual != sequence[:1] + sequence[3:6]:
00395
                       logging.warning(
00396
                           self.logwarn + "Due to sequence length and use of the Icarus model "
00397
                           "1 sensor, the actual timing sequence for side "
00398
                           + side
                          + " will be "
+ "{"
00399
00400
00401
                           + str(actual[0] + f0delay)
00402
00403
00404
                           + str(actual[1 : 2 * self.nframes])
00405
                      )
00406
                  else:
00407
                      logging.warning(
                           " timing sequence for side "
00408
00409
                          + side
+ " will be "
00410
00411
                           + "{"
00412
00413
                           + str(actual[0] + f0delay)
                           + "}"
+ ""
00414
00415
00416
                           + str(actual[1 : 2 * self.nframes])
00417
00418
              else:
                  if actual != sequence:
00419
00420
                       logging.warning(
00421
                           self.logwarn + "Due to sequence length, actual timing sequence "
00422
                           "for side "
00423
                           + side
00424
                           + " will be "
                           + "{"
00425
00426
                           + str(actual[0])
                           + "}"
00427
00428
00429
                           + str(actual[1 : 2 * self.nframes])
00430
              return "", actual
00431
00432
```

### 6.11.3.11 setExtClk()

Reimplemented in nsCamera.sensors.daedalus.daedalus.

Definition at line 136 of file sensorBase.py.

### 6.11.3.12 setHighFullWell()

Reimplemented in nsCamera.sensors.daedalus.daedalus.

```
Definition at line 96 of file sensorBase.py.
```

```
00096
         def setHighFullWell(self, flag):
00097
00098
              Virtual function; feature is not implemented on Icarus
00099
              Overridden in daedalus.py
00100
00101
              if flag:
                 logging.warning(
00102
                      self.logwarn + "HighFullWell mode is not supported by Icarus sensors. "
00103
00104
00105
```

## 6.11.3.13 setInterlacing()

Reimplemented in nsCamera.sensors.daedalus.daedalus.

#### Definition at line 82 of file sensorBase.py.

```
00082
         def setInterlacing(self, ifactor):
00083
00084
              Virtual function; feature is not implemented on Icarus
00085
              Overridden in daedalus.py
00086
00087
              Returns:
              integer 0
00088
00089
00090
              if ifactor:
00091
                  logging.warning(
00092
                     self.logwarn + "Interlacing is not supported by Icarus sensors. "
00093
                  )
00094
              return 0
00095
```

### 6.11.3.14 setManualShutters()

Reimplemented in nsCamera.sensors.daedalus.daedalus.

#### Definition at line 531 of file sensorBase.py.

```
00531 def setManualShutters(self, timing=None):
00532 """
00533 Legacy alias for setManualTiming()
00534 """
00535 self.setManualTiming(timing)
```

### 6.11.3.15 setManualTiming()

```
nsCamera.sensors.sensorBase.sensorBase.setManualTiming (
              self,
              timing = None)
Manual shutter timing, seven intervals to assign to both hemispheres, e.g.,
  [(100,150,100,150,100,150,100)] for frame 0 open for 100 ns, an interframe
  pause of 50 ns, frame 1 open for 100 ns, etc.
Provide two sets of seven intervals, e.g., [(100,150,100,150,100,150,100),
  (200,250,200,250,200,250,200)] to program the A and B hemispheres
  independently
Overridden in daedalus.py
The timing list is flattened before processing; the suggested tuple structure is
  just for clarity (first tuple is A, second is B) and is optional.
The actual timing is rounded down to the nearest multiple of 25 ns. (Each
  count = 25 ns. e.g., a request for 140 ns rounds down to a count of '5',
  which corresponds to 125 ns))
    - Minimum timing is 75 ns
    - Maximum is 25 * 2^30 ns (approximately 27 seconds)
Aras:
   timing: 7- or 14-element list (substructure optional) in nanoseconds
Returns:
    tuple (error string, response string from final message)
```

Reimplemented in nsCamera.sensors.daedalus.daedalus.

### Definition at line 537 of file sensorBase.py.

```
00537 def setManualTiming(self, timing=None):
00538 """
00539 Manual shutter timing, seven intervals to assign to both hemispheres, e.g.,
00540 [(100,150,100,150,100,150,100)] for frame 0 open for 100 ns, an interframe
00541 pause of 50 ns,frame 1 open for 100 ns, etc.
00542 Provide two sets of seven intervals, e.g., [(100,150,100,150,100,150,100),
00543 (200,250,200,250,200,250,200)] to program the A and B hemispheres
00544 independently
```

```
00545
00546
                               Overridden in daedalus.pv
00547
                               The timing list is flattened before processing; the suggested tuple structure is
00548
00549
                                   just for clarity (first tuple is A, second is B) and is optional.
00550
00551
                               The actual timing is rounded down to the nearest multiple of 25\ \mathrm{ns}. (Each
00552
                                   count = 25 ns. e.g., a request for 140 ns rounds down to a count of '5',
00553
                                   which corresponds to 125 ns))
00554
                                       - Minimum timing is 75 ns
00555
                                        - Maximum is 25 * 2^30 ns (approximately 27 seconds)
00556
00557
                               Args:
00558
                                       timing: 7- or 14-element list (substructure optional) in nanoseconds
00559
00560
                               tuple (error string, response string from final message) """
00561
00562
                               if timing is None:
00563
00564
                                       logging.info(
00565
                                                 self.loginfo
00566
                                                  + "No manual timing setting provided, defaulting to (100, 150, 100, "
00567
                                                  " 150, 100, 150, 100) for both hemispheres"
00568
00569
                                       timing = [
00570
                                                  (100, 150, 100, 150, 100, 150, 100),
00571
                                                  (100, 150, 100, 150, 100, 150, 100),
00572
00573
                               logging.info(self.loginfo + "Manual shutter sequence: " + str(timing))
00574
                               flattened = flattenlist(timing)
00575
                               if len(flattened) == 7:
00576
                                        flattened = 2 * flattened
00577
                               if (
00578
                                        len(flattened) != 14
00579
                                        or not all(isinstance(x, numbers.Real) for x in flattened)
00580
                                        or not all(x \ge 75 for x in flattened)
                                        or not all(x <= 26843545600 for x in flattened)
00581
00582
                                       logging.error(err + "; timing settings unchanged")
return err, "00000000"
                                        err = self.logerr + "Invalid manual shutter timing list: " + str(timing)
00583
00584
00585
00586
00587
                               timecounts = [int(a // 25) for a in flattened]
00588
                               self.ca.sensmanual = timing
                               self.ca.senstiming = {} # clear HST settings from ca object
00589
00590
00591
                               control_messages = [
                                         [col_messages = [ ("WO_INTEGRATION", "\{0:\#0\{1\}x\}".format(timecounts[0], 10)[2:10]), ("WO_INTERFRAME", "\{0:\#0\{1\}x\}".format(timecounts[1], 10)[2:10]), ("WI_INTEGRATION", "\{0:\#0\{1\}x\}".format(timecounts[2], 10)[2:10]), ("WI_INTERFRAME", "\{0:\#0\{1\}x\}".format(timecounts[3], 10)[2:10]), ("W2_INTEGRATION", "\{0:\#0\{1\}x\}".format(timecounts[4], 10)[2:10]), ("W2_INTEGRATION", "\{0:#0\{1\}x\}".format(timecounts[5], 10)[2:10]), ("W3_INTEGRATION", "\{0:#0\{1\}x\}".format(timecounts[6], 10)[2:10]), ("W3_INTEGRATION", "\{0:#0\{1\}x\}".format(timecounts[6], 10)[2:10]), ("W3_INTEGRATION R". "\{0:#0\{1\}x\}".format(timecounts[7], 10)[2:10]), ("W3_INTEGRATION R". "(0:#0\{1\}x]".format(timecounts[7], 10)[2:10]), ("W
00592
00593
00594
00595
00596
00597
00598
                                        ("W0_INTEGRATION", "{0:#0{1}x}".format(timecounts[6], 10)[2:10]),

("W0_INTEGRATION_B", "{0:#0{1}x}".format(timecounts[7], 10)[2:10]),

("W0_INTERFRAME_B", "{0:#0{1}x}".format(timecounts[8], 10)[2:10]),

("W1_INTEGRATION_B", "{0:#0{1}x}".format(timecounts[9], 10)[2:10]),

("W1_INTEGRATION_B", "{0:#0{1}x}".format(timecounts[10], 10)[2:10]),

("W2_INTEGRATION_B", "{0:#0{1}x}".format(timecounts[11], 10)[2:10]),

("W2_INTEGRATION_B", "{0:#0{1}x}".format(timecounts[12], 10)[2:10]),
00599
00600
00601
00602
00603
00604
                                         ("W3_INTEGRATION_B", "{0:#0{1}x}".format(timecounts[13], 10)[2:10]),
00605
                                         ("HST_MODE", "0"),
00606
                                         ("MANSHUT_MODE", "1"),
00607
00608
00609
                               return self.ca.submitMessages(control_messages, " setManualShutters: ")
00610
```

#### 6.11.3.16 setPhiDelay()

Reimplemented in nsCamera.sensors.daedalus.daedalus.

```
Definition at line 126 of file sensorBase.py.
```

### 6.11.3.17 setTiming()

```
nsCamera.sensors.sensorBase.sensorBase.setTiming (
              self,
              side = "AB",
              sequence = None,
              delay = 0 )
Sets timing registers based on 'sequence.' Requesting (0,0) timing will clear the
 timing register.
*WARNING* if the entire sequence does not fit into the 40-bit register space,
 then the actual timings generated may differ from those requested. If the
  timing sequence fits only once into the register space (i.e., for a single
  frame, open + closed > 20 ns), then the actual timing will be (n, 40-n),
 irrespective of the setting of second parameter, e.g. (35,1) will actually
  result in (35,5) timing.
*NOTE* Icarus sensors generally cannot use 1 ns timing, so all values for these
  devices (besides the delay) should be at least 2 ns
    side: Hemisphere 'A', 'B', 'AB'
    sequence: two-element tuple of timing durations in ns, e.g., '(5,2)'
    delay: initial delay in ns (1 ns delay is acceptable)
   tuple (error string, 10-character hexadecimal representation of timing
     sequence)
```

#### Definition at line 149 of file sensorBase.py.

```
def setTiming(self, side="AB", sequence=None, delay=0):
00149
00150
00151
              Sets timing registers based on 'sequence.' Requesting (0,0) timing will clear the
00152
                timing register.
00153
              *WARNING* if the entire sequence does not fit into the 40-bit register space,
00154
                then the actual timings generated may differ from those requested. If the
00155
                timing sequence fits only once into the register space (i.e., for a single
00156
                frame, open + closed > 20 ns), then the actual timing will be (n, 40-n),
00157
                irrespective of the setting of second parameter, e.g. (35,1) will actually
00158
                result in (35,5) timing.
00159
              *NOTE* Icarus sensors generally cannot use 1 ns timing, so all values for these
00160
                devices (besides the delay) should be at least 2 ns
00161
00162
00163
                  side: Hemisphere 'A', 'B', 'AB'
00164
                  sequence: two-element tuple of timing durations in ns, e.g., '(5,2)'
00165
                  delay: initial delay in ns (1 ns delay is acceptable)
00166
00167
              Returns:
00168
                 tuple (error string, 10-character hexadecimal representation of timing
00169
                   sequence)
00170
00171
              logging.info(
```

```
00172
                                    "{}setTiming; side={}, sequence={}, delay={}".format(
00173
                                           self.loginfo, side, sequence, delay
00174
00175
00176
                           if sequence is None:
00177
                                   sequence = (3, 2)
00178
                           if delay is None:
00179
                                   delay = 0
00180
                           logging.info(
00181
                                 self.loginfo
                                   + "HST side "
00182
00183
                                   + side.upper()
00184
                                  + ": "
                                   + str(sequence)
00186
                                   + "; delay =
                                   + str(delay)
00188
                           )
                           err = ""
00189
00190
                           if len(sequence) != 2:
00191
                                   err = (
00192
                                           self.logerr
                                            + "Invalid sequence setting for side: "
00193
00194
                                            + side
00195
                                            + "; timing settings are unchanged"
00196
00197
                                   logging.error(err)
                           return err, "0000000000"

if side.upper() == "AB":
00198
00199
                                                 _ = self.setTiming(side="A", sequence=sequence, delay=delay)
00200
                                   err1, _ = self.setTiming(side="A", sequence=sequence, delay=delay)
err2, full40hex = self.setTiming(side="B", sequence=sequence, delay=delay)
00201
                                    return err1 + err2, full40hex
00202
00203
                           if side.upper() == "A":
                                   lowreg = "HS_TIMING_DATA_ALO"
highreg = "HS_TIMING_DATA_AHI"
00204
00205
                           elif side.upper() == "B":
00206
                                  lowreg = "HS_TIMING_DATA_BLO"
00207
                                   highreg = "HS_TIMING_DATA_BHI"
00208
00209
                           else:
00210
                                  err = (
                                           self.logerr
00211
00212
                                            + "setTiming: Invalid sensor side: "
00213
                                           + side
00214
                                            + "; timing settings unchanged"
00215
00216
                                   logging.error(err)
00217
                                    return err, "0000000000"
00218
                           if (sequence[0] + sequence[1]) + delay > 40:
00219
                                   err = (
00220
                                           self.logerr
00221
                                            + "setTiming: Timing sequence is too long to be implemented; "
00222
                                            "timing settings unchanged "
00223
00224
                                   logging.error(err)
00225
                                    return err, "0000000000"
00226
00227
                           self.ca.senstiming[side.upper()] = (sequence, delay)
00228
                           self.ca.sensmanual = [] # clear manual settings from ca
00229
00230
                           full40 = [0] * 40
00231
                           bitlist = []
00232
                           flag = 1
00233
                           sequence = sequence[:2] # TODO: is this redundant?
00234
                           for a in sequence:
                                   add = [flag] * a
00235
                                   bitlist += add
00236
00237
                                   if flag:
00238
                                           flag = 0
00239
                                   else:
00240
                                        flag = 1
00241
                           # automatically truncates sequence to 39 characters
                           logging.debug(self.logdebug + "bitlist = " + str(bitlist))
00242
                           if bitlist: # skip this if timing is [0,0]
    reversedlist = bitlist[39::-1]
00243
00244
00245
                                   trunclist = reversedlist[:]
00246
                                   while trunclist[0] == 0:
00247
                                           trunclist.pop(0)
                                    # fullrepeat counts open/closed cycles, doesn't include final frame
00248
                                   "Introduct Control of 
00249
00250
                                   # Pattern from sequence repeated to fit inside 40 bits
repeated = trunclist + reversedlist * fullrepeats
00251
00252
```

```
00253
                   full40[-(len(repeated) + delay + 1) : -(delay + 1)] = repeated
00254
              else:
00255
                  logging.warning(self.logwarn + "setTiming: all-zero timing supplied")
                   fullrepeats = self.nframes
00256
00257
               full 40bin = "".join(str(x) for x in full 40)
00258
               logging.debug(self.logdebug + "full40bin = " + str(full40bin))
00259
               full40hex = "%x" % int(full40bin, 2)
00260
               logging.debug(self.logdebug + "full40hex = " + str(full40hex))
00261
              highpart = full40hex[-10:-8].zfill(8)
              lowpart = full40hex[-8:].zfill(8)
00262
00263
              err0, _ = self.ca.setRegister(lowreg, lowpart)
              err1, _ = self.ca.setRegister(highreg, highpart)
00264
              err2, _ = self.ca.setSubregister("MANSHUT_MODE", "0")
err3, _ = self.ca.setSubregister("HST_MODE", "1")
00265
00266
00267
              err = err0 + err1 + err2 + err3
00268
              if err:
00269
                  logging.error(
00270
                       self.logerr + "setTiming: Timing may not have been set correctly"
00271
00272
              if fullrepeats < self.nframes - 1:</pre>
00273
                  actual = self.getTiming(side, actual=True)
00274
                   if self.fpganumID == 2:
00275
                      expected = [delay] + 2 * list(sequence) + [sequence[0]]
00276
                   else:
00277
                      expected = [delay] + 3 * list(sequence) + [sequence[0]]
00278
                   if actual != expected:
00279
                       logging.warning(
00280
                           self.logwarn
00281
                           + "setTiming: Due to sequence length"
00282
                           + self.specwarn
00283
                            + ", the actual timing "
00284
                           "sequence for side '
00285
                           + side
                           + " will be "
+ "{"
00286
00287
00288
                           + str(actual[0])
                           + "}"
00289
00290
00291
                           + str(actual[1 : 2 * self.nframes])
00292
00293
              elif self.ca.sensorname == "icarus":
00294
                   f0delay = sequence[0] + sequence[1]
00295
                   logging.warning(
00296
                       self.logwarn + "setTiming: Due to use of the Icarus model 1 sensor, the"
00297
                       " initial delay for side "
00298
                       + side
00299
                       + " will actually be "
00300
                       + str(delay + f0delay)
00301
                       + " nanoseconds"
00302
00303
              return err, full40hex
00304
```

### 6.11.3.18 setTriggerDelay()

Reimplemented in nsCamera.sensors.daedalus.daedalus.

# Definition at line 116 of file sensorBase.py.

```
00116
          def setTriggerDelay(self, delay):
00117
00118
              Virtual function; feature is not implemented on Icarus
00119
              Overridden in daedalus.py
00120
              if delay:
00121
00122
                  logging.warning(
                      self.logwarn + "Trigger Delay is not supported by Icarus sensors. "
00123
00124
00125
```

## 6.11.3.19 setZeroDeadTime()

```
nsCamera.sensors.sensorBase.sensorBase.setZeroDeadTime ( self, \\ flag \ ) Virtual function; feature is not implemented on Icarus Overridden in daedalus.py
```

Reimplemented in nsCamera.sensors.daedalus.daedalus.

### Definition at line 106 of file sensorBase.py.

```
def setZeroDeadTime(self, flag):
00106
00107
00108
             Virtual function; feature is not implemented on Icarus
             Overridden in daedalus.py
00109
00110
00111
            if flag:
00112
                logging.warning(
                     self.logwarn + "ZeroDeadTime mode is not supported by Icarus sensors. "
00113
00114
00115
```

### 6.11.4 Member Data Documentation

#### 6.11.4.1 ca

nsCamera.sensors.sensorBase.sensorBase.ca

Definition at line 32 of file sensorBase.py.

## 6.11.4.2 fpganumID

nsCamera.sensors.sensorBase.sensorBase.fpganumID

Definition at line 274 of file sensorBase.py.

#### 6.11.4.3 logcrit

nsCamera.sensors.sensorBase.sensorBase.logcrit

Definition at line 38 of file sensorBase.py.

## 6.11.4.4 logdebug

nsCamera.sensors.sensorBase.sensorBase.logdebug

Definition at line 42 of file sensorBase.py.

## 6.11.4.5 logerr

nsCamera.sensors.sensorBase.sensorBase.logerr

Definition at line 39 of file sensorBase.py.

### 6.11.4.6 loginfo

nsCamera.sensors.sensorBase.sensorBase.loginfo

Definition at line 41 of file sensorBase.py.

### 6.11.4.7 logwarn

nsCamera.sensors.sensorBase.sensorBase.logwarn

Definition at line 40 of file sensorBase.py.

#### 6.11.4.8 sensfam

nsCamera.sensors.sensorBase.sensorBase.sensfam

Definition at line 337 of file sensorBase.py.

The documentation for this class was generated from the following file:

C:/Users/hill35/git/camera\_python/nsCamera/sensors/sensorBase.py

# 6.12 nsCamera.utils.Subregister.SubRegister Class Reference

#### **Public Member Functions**

• \_\_init\_\_ (self, board, name, register, start\_bit=31, width=8, writable=False, value=255, minV=0, maxV=5)

#### **Public Attributes**

- name
- register
- addr
- · start bit
- width
- value
- max\_value
- min
- max
- · writable
- minV
- maxV
- · resolution

### 6.12.1 Detailed Description

```
Represents a subset of a 32-bit register [31..0] starting at 'start_bit' consisting of 'width' bits. Consistent with the ICD usage, start_bit is MSB e.g., for [7..0], the start_bit is '7'.
```

Definition at line 21 of file Subregister.py.

#### 6.12.2 Constructor & Destructor Documentation

#### 6.12.2.1 init ()

### Definition at line 28 of file Subregister.py.

```
00039
00040
               self.name = name
00041
               self.register = register
00042
               self.addr = board.registers[register]
00043
               self.start_bit = start_bit
               self.width = width
00045
               self.value = value
00046
               self.max\_value = 2**width - 1 # used to normalize the input values to 1
00047
               self.min = 0
              self.max = self.max_value
00048
00049
               self.writable = writable
00050
              self.minV = minV
00051
               self.maxV = maxV
00052
               # resolution should be reset after init if actual min and max are different
00053
               self.resolution = (1.0 * maxV - minV) / self.max_value
00054
00055
00056 """
00057 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00058 LLNL-CODE-838080
00059
00060 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00061 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE) 00062 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00063 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00064 be made under this license. 00065 """
```

### 6.12.3 Member Data Documentation

#### 6.12.3.1 addr

nsCamera.utils.Subregister.SubRegister.addr

Definition at line 42 of file Subregister.py.

230 Class Documentation

#### 6.12.3.2 max

 $\verb|nsCamera.utils.Subregister.SubRegister.max| \\$ 

Definition at line 48 of file Subregister.py.

#### 6.12.3.3 max\_value

nsCamera.utils.Subregister.SubRegister.max\_value

Definition at line 46 of file Subregister.py.

#### 6.12.3.4 maxV

nsCamera.utils.Subregister.SubRegister.maxV

Definition at line 51 of file Subregister.py.

#### 6.12.3.5 min

nsCamera.utils.Subregister.SubRegister.min

Definition at line 47 of file Subregister.py.

#### 6.12.3.6 minV

nsCamera.utils.Subregister.SubRegister.minV

Definition at line 50 of file Subregister.py.

#### 6.12.3.7 name

nsCamera.utils.Subregister.SubRegister.name

Definition at line 40 of file Subregister.py.

### 6.12.3.8 register

nsCamera.utils.Subregister.SubRegister.register

Definition at line 41 of file Subregister.py.

#### 6.12.3.9 resolution

nsCamera.utils.Subregister.SubRegister.resolution

Definition at line 53 of file Subregister.py.

### 6.12.3.10 start\_bit

nsCamera.utils.Subregister.SubRegister.start\_bit

Definition at line 43 of file Subregister.py.

#### 6.12.3.11 value

nsCamera.utils.Subregister.SubRegister.value

Definition at line 45 of file Subregister.py.

#### 6.12.3.12 width

nsCamera.utils.Subregister.SubRegister.width

Definition at line 44 of file Subregister.py.

#### 6.12.3.13 writable

nsCamera.utils.Subregister.SubRegister.writable


Definition at line 49 of file Subregister.py.

The documentation for this class was generated from the following file:

C:/Users/hill35/git/camera\_python/nsCamera/utils/Subregister.py

# 6.13 nsCamera.comms.GigE.GigE.ZESTETM1\_CARD\_INFO Class Reference

Inheritance diagram for nsCamera.comms.GigE.GigE.ZESTETM1\_CARD\_INFO:



232 Class Documentation

#### **Static Public Attributes**

```
int ubyte4 = C.c_ubyte * 4int ubyte6 = C.c ubyte * 6
```

#### **Static Protected Attributes**

list \_fields\_

#### 6.13.1 Detailed Description

Definition at line 512 of file GigE.py.

#### 6.13.2 Member Data Documentation

#### 6.13.2.1 \_fields\_

```
list nsCamera.comms.GigE.GigE.ZESTETM1_CARD_INFO._fields_ [static], [protected]
```

#### Initial value:

```
("IPAddr", ubyte4),
    ("ControlPort", C.c_ushort),
    ("Timeout", C.c_ulong),
    ("HTTPPort", C.c_ushort),
    ("MACAddr", ubyte6),
    ("SubNet", ubyte4),
    ("Gateway", ubyte4),
    ("SerialNumber", C.c_ulong),
    ("FirmwareVersion", C.c_ulong),
    ("HardwareVersion", C.c_ulong),
```

Definition at line 515 of file GigE.py.

#### 6.13.2.2 ubyte4

```
int nsCamera.comms.GigE.GigE.ZESTETM1_CARD_INFO.ubyte4 = C.c_ubyte * 4 [static]
```

Definition at line 513 of file GigE.py.

#### 6.13.2.3 ubyte6

```
int nsCamera.comms.GigE.GigE.ZESTETM1_CARD_INFO.ubyte6 = C.c_ubyte * 6 [static]
```

Definition at line 514 of file GigE.py.

The documentation for this class was generated from the following file:

C:/Users/hill35/git/camera\_python/nsCamera/comms/GigE.py

# **Chapter 7**

# **File Documentation**

# 7.1 C:/Users/hill35/git/camera python/nsCamera/ init .py File Reference

#### **Namespaces**

namespace nsCamera

#### **Variables**

list nsCamera.\_\_all\_\_ = ["CameraAssembler"]

# 7.2 \_\_init\_\_.py

#### Go to the documentation of this file.

```
00001 # -*- coding: utf-8 -*
00002 """
00003 Created on Tue Mar 22 15:47:43 2016
00005 The Package includes a Camera object and an assembler.
00007 The camera object will be the workhorse of the API. The assembler is used
00008 to create the Camera object.
00010 Author: Matthew Dayton (dayton5@llnl.gov)
00012 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00015 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00016 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00017 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00018 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00019 be made under this license.
00020
00021 Version: 2.1.2 (February 2025)
00022 """
00023
00024 from nsCamera import CameraAssembler
00025
00026 __all__ = ["CameraAssembler"]
00027
00028 """
```

```
00029 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00030 LLNL-CODE-838080
00031
00032 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00033 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00034 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00035 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00036 be made under this license.
```

# 7.3 C:/Users/hill35/git/camera\_python/nsCamera/boards/\_\_init\_\_.py File Reference

#### **Namespaces**

- namespace nsCamera
- namespace nsCamera.boards

#### **Variables**

list nsCamera.boards.\_\_all\_\_ = ["LLNL\_v1", "LLNL\_v4"]

## 7.4 \_\_init\_\_.py

#### Go to the documentation of this file.

```
00001 # -*- coding: utf-8 -*-
00002 """
00003 This package is a collection of modules that represent the camera boards
00004 Each board has its own number of ADCs, POTs, and sensors. More devices can be added in
00005 the future. The list of imports will grow as we make more types of boards.
00006
00007 Author: Matthew Dayton (dayton5@llnl.gov)
00008
00009 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00010 LLNL-CODE-838080
00011
00012 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00013 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00014 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00015 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00016 be made under this license.
00017
00018 Version: 2.1.2 (February 2025)
00019 """
00020 from .LLNL_v1 import llnl_v1
00021 from .LLNL_v4 import llnl_v4
00022
00023 __all__ = ["LLNL_v1", "LLNL_v4"]
00024
00025 """
00026 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00029 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00030 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00031 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00032 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00033 be made under this license.
00034 """
```

# 7.5 C:/Users/hill35/git/camera\_python/nsCamera/comms/\_\_init\_\_.py File Reference

#### **Namespaces**

- · namespace nsCamera
- namespace nsCamera.comms

#### **Variables**

• list nsCamera.comms.\_\_all\_\_ = ["RS422", "GigE"]

## 7.6 init .py

```
Go to the documentation of this file.
```

```
00001 # -*- coding: utf-8 -*-
00002 """
00003 Created on Tue Mar 22 15:47:43 2016
00004
00005 This package is a collection of modules for uniform handling of the nsCamera's
00006
       communication systems
00007
00008 Author: Matthew Dayton (dayton5@llnl.gov)
00009
00010 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00011 LLNL-CODE-838080
00012
00013 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00014 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00015 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00016 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00017 be made under this license.
00018
00019 Version: 2.1.2 (February 2025)
00020 """
00021
00022 from .GigE import GigE
00023 from .RS422 import RS422
00024
00025 __all__ = ["RS422", "GigE"]
00026
00027 """
00028 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00029 LLNL-CODE-838080
00031 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00032 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00033 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00034 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00035 be made under this license.
00036 """
```

# 7.7 C:/Users/hill35/git/camera\_python/nsCamera/sensors/\_\_init\_\_.py File Reference

#### **Namespaces**

- namespace nsCamera
- namespace nsCamera.sensors

#### **Variables**

• list nsCamera.sensors. all = ["icarus", "icarus2", "daedalus"]

## 7.8 init .py

```
Go to the documentation of this file.
```

```
00001 # -*- coding: utf-8 -*- 00002 """
00003 This package is a collection of modules for uniform handling of the various sensors
00004
00005 Author: Jeremy Martin Hill (jerhill@llnl.gov)
00006
00007 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00008 LLNL-CODE-838080
00009
00010 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00011 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00012 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00013 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00014 be made under this license.
00015
00016 Version: 2.1.2 (February 2025)
00017 """
00018
00019 from .daedalus import daedalus
00020 from .icarus import icarus
00021 from .icarus2 import icarus2
00023 __all__ = ["icarus", "icarus2", "daedalus"]
00024
00026 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00027 LLNL-CODE-838080
00029 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00030 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00031 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00032 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00033 be made under this license.
00034 """
```

# 7.9 C:/Users/hill35/git/camera\_python/nsCamera/utils/\_\_init\_\_.py File Reference

#### **Namespaces**

- · namespace nsCamera
- namespace nsCamera.utils

#### **Variables**

• list nsCamera.utils.\_\_all\_\_ = ["SubRegister", "Packet", "FlatField", "misc"]

7.10 init .py 237

## 7.10 \_\_init\_\_.py

```
Go to the documentation of this file.
```

```
00001 # -*- coding: utf-8 -*-
00003 This package is a collection of utility classes for the CameraAssembler
00004
00005 Author: Jeremy Martin Hill (jerhill@llnl.gov)
00006
00007 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00008 LLNL-CODE-838080
00009
00010 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00011 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00012 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00013 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00014 be made under this license.
00015
00016 Version: 2.1.2 (February 2025) 00017 """
00018
00019 from .Packet import Packet
00020 from .Subregister import SubRegister
00021
00022 try:
00023
         from misc import *
00024
          from .crc16pure import *
00025
          from .FlatField import *
00026
00027 except:
00028
00029
00030 __all__ = ["SubRegister", "Packet", "FlatField", "misc"]
00031
00032 """
00033 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00034 LLNL-CODE-838080
00036 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00037 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00038 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00039 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00040 be made under this license.
00041 """
```

# 7.11 C:/Users/hill35/git/camera\_python/nsCamera/boards/LLNL\_v1.py File Reference

#### Classes

· class nsCamera.boards.LLNL\_v1.llnl\_v1

#### **Namespaces**

- namespace nsCamera
- namespace nsCamera.boards
- namespace nsCamera.boards.LLNL v1

## 7.12 LLNL v1.py

#### Go to the documentation of this file.

```
00001 # -*- coding: utf-8 -*
00002 """
00003 LLNLv1 board definition, including monitors, pots, and other board-specific settings
00005 Author: Jeremy Martin Hill (jerhill@llnl.gov)
00006 Author: Matthew Dayton (dayton5@llnl.gov)
00008 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00009 LLNL-CODE-838080
00011 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00012 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00013 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00014 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00015 be made under this license.
00016
00017 Version: 2.1.2 (February 2025)
00018 """
00019
00020 import logging
00021 import time
00022 from collections import OrderedDict
00023
00024 from nsCamera.utils.Packet import Packet
00025 from nsCamera.utils.Subregister import SubRegister
00026
00027
00028 class llnl_v1:
00029
00030
           Livermore LLNL v1.0 board
00031
00032
           Compatible communication protocols: RS422, GigE
00033
           Compatible sensors: icarus, icarus2, daedalus
00034
00035
00036
           # FPGA register map - use '.upper()' on keys to ensure uppercase lookup
00037
           registers = OrderedDict(
00038
                    "FPGA_NUM": "000", "FPGA_REV": "001",
00039
00040
                    "HS_TIMING_CTL": "010",
"HS_TIMING_DATA_ALO": "013",
00041
00042
                    "HS_TIMING_DATA_AHI": "014",
00043
00044
                    "HS_TIMING_DATA_BLO": "015"
                    "HS_TIMING_DATA_BHI": "016",
00045
00046
                    "SW_TRIGGER_CONTROL": "017",
00047
                    "STAT_REG": "024",
                   "CTRL_REG": "025",
"POT_CTL": "026",
"POT_REG4_TO_1": "027",
00048
00049
00050
00051
                    "POT_REG8_TO_5": "028",
00052
                    "POT_REG12_TO_9": "029",
                    "POT_REG13": "02A",
00053
00054
                    "LED_GP": "02B",
                    "SW_RESET": "02D"
                    "HST_SETTINGS": "02E",
"STAT_REG_SRC": "02F",
00056
00057
00058
                    "STAT_REG2": "030",
                    "STAT_REG2_SRC": "031",
                    "ADC_BYTECOUNTER": "032",
                    "RBP_PIXEL_CNTR": "033",
                    "DIAG_MAX_CNT_0": "034",
                    "DIAG_MAX_CNT_1": "035",
00063
                    "DIAG_CNTR_VAL_0": "036",
00064
00065
                    "DIAG_CNTR_VAL_1": "037",
                    "TRIGGER_CTL": "03A",
00066
                    "SRAM_CTL": "03B",
00067
                    "TIMER_CTL": "03C"
00068
00069
                    "TIMER_VALUE": "03D",
00070
                    "HSTALLWEN_WAIT_TIME": "03F",
00071
                    "FPA ROW INITIAL": "042",
                    "FPA_ROW_FINAL": "043",
00072
00073
                    "FPA_FRAME_INITIAL": "044",
                    "FPA_FRAME_FINAL": "045",
"FPA_DIVCLK_EN_ADDR": "046",
00074
00075
                    "FPA_OSCILLATOR_SEL_ADDR": "047",
00076
```

7.12 LLNL\_v1.py 239

```
00077
                                              "FRAME_ORDER_SEL": "04B",
00078
                                              "SENSOR_VOLT_STAT": "082",
00079
                                              "SENSOR_VOLT_CTL": "083",
00080
                                              "ADC_CTL": "090",
                                              "ADC1_CONFIG_DATA": "091",
00081
00082
                                              "ADC2_CONFIG_DATA": "092",
00083
                                              "ADC3_CONFIG_DATA": "093",
                                              "ADC4_CONFIG_DATA": "094",
00084
00085
                                              "ADC5_CONFIG_DATA": "095",
                                              "ADC5_DATA_1": "096",
00086
                                              "ADC5_DATA_2": "097",
00087
                                              "ADC5_DATA_3": "098",
00088
00089
                                              "ADC5_DATA_4": "099",
                                              "ADC5_PPER": "09A",
00090
00091
                                              "ADC_STANDBY": "09B", # board version <= AD
00092
                                              "ADC_RESET": "09B",
                                                                                                  # board version > AD
                                              "TEMP_SENSE_PPER": "OAO",
00093
                                              "TEMP_SENSE_DATA": "OA1"
00094
00095
00096
                       )
00097
00098
                         subregisters = [
00099
00103
                                    ("HST_MODE", "HS_TIMING_CTL", 0, 1, True),
                                    ("SW_TRIG_START", "SW_TRIGGER_CONTROL", 0, 1, True), ("LED_EN", "CTRL_REG", 1, 1, True),
00104
00105
                                    ("COLQUENCHEN", "CTRL_REG", 2, 1, True),
("POWERSAVE", "CTRL_REG", 3, 1, True),
00106
                                  ("COLQUENCHEN", "CTRL_REG", 2, 1, True),
("POWERSAVE", "CTRL_REG", 3, 1, True),
("POT1", "POT_REG4_TO_1", 7, 8, True),
("POT2", "POT_REG4_TO_1", 15, 8, True),
("POT3", "POT_REG4_TO_1", 23, 8, True),
("POT4", "POT_REG4_TO_1", 31, 8, True),
("POT5", "POT_REG8_TO_5", 7, 8, True),
("POT6", "POT_REG8_TO_5", 15, 8, True),
("POT7", "POT_REG8_TO_5", 31, 8, True),
("POT7", "POT_REG8_TO_5", 31, 8, True),
("POT9", "POT_REG12_TO_9", 7, 8, True),
("POT10", "POT_REG12_TO_9", 15, 8, True),
("POT11", "POT_REG12_TO_9", 15, 8, True),
("POT12", "POT_REG12_TO_9", 31, 8, True),
("POT13", "POT_REG13", 7, 8, True),
("LED1", "LED_GP", 0, 1, True),
("LED4", "LED_GP", 1, 1, True),
("LED4", "LED_GP", 3, 1, True),
("LED5", "LED_GP", 4, 1, True),
("LED6", "LED_GP", 5, 1, True),
("LED7", "LED_GP", 6, 1, True),
("LED8", "LED_GP", 6, 1, True),
("LED8", "LED_GP", 7, 1, True),
00107
00108
00109
00110
00111
00112
00113
00114
00115
00116
00117
00118
00119
00120
00121
00122
00123
00124
00125
00126
00127
                                    ("LED8", "LED_GP", 7, 1, True),
("RESET", "SW_RESET", 0, 1, True),
00128
00129
                                    ("HST_SW_CTL_EN", "HST_SETTINGS", 0, 1, True), ("SW_HSTALLWEN", "HST_SETTINGS", 1, 1, True),
00130
00131
                                    ("MAXERR FIT", "DIAG_MAX_CNT_0", 31, 16, True), ("MAXERR_SRT", "DIAG_MAX_CNT_0", 7, 8, True), ("MAXERR_UTTR", "DIAG_MAX_CNT_1", 31, 16, True), ("MAXERR_UTTR", "DIAG_MAX_CNT_1", 15, 16, True),
00132
00133
00134
00135
                                    ("HM_TRIG_EN", "TRIGGER_CTL", 0, 1, True), ("SW_TRIG_EN", "TRIGGER_CTL", 2, 1, True), ("READ_SRAM", "SRAM_CTL", 0, 1, True),
00136
00137
00138
00139
                                    ("RESET_TIMER", "TIMER_CTL", 0, 1, True),
                                    ("OSC_SELECT", "FPA_OSCILLATOR_SEL_ADDR", 1, 2, True),
("ADC5_VREF", "ADC5_CONFIG_DATA", 9, 10, True),
("ADC5_VREF3", "ADC5_CONFIG_DATA", 13, 1, True),
("ADC5_INT", "ADC5_CONFIG_DATA", 15, 1, True),
("ADC5_MULT", "ADC5_CONFIG_DATA", 24, 6, True),
00140
00141
00142
00143
00144
                                    ("PPER", "ADC5_PPER", 7, 8, True),
00145
00146
                                    ("SRAM_READY", "STAT_REG", 0, 1, False), ("STAT_COARSE", "STAT_REG", 1, 1, False), ("STAT_FINE", "STAT_REG", 2, 1, False),
00151
00152
00153
00154
                                    ("STAT_SENSREADIP", "STAT_REG", 5, 1, False),
                                    ("STAT_SENSREADIP", "STAT_REG", 5, 1, False),
("STAT_SENSREADDONE", "STAT_REG", 6, 1, False),
("STAT_SRAMREADSTART", "STAT_REG", 7, 1, False),
("STAT_SRAMREADDONE", "STAT_REG", 8, 1, False),
("STAT_HSTCONFIGSTART", "STAT_REG", 9, 1, False),
("STAT_ADCSCONFIGURED", "STAT_REG", 10, 1, False),
("STAT_POTSCONFIGURED", "STAT_REG", 11, 1, False),
00155
00156
00157
00158
00159
00160
                                    ("STAT_POISCONFIGURED", "STAT_REG", 11, 1, False),
("STAT_TIMERCOUNTERRESET", "STAT_REG", 13, 1, False),
("STAT_ARMED", "STAT_REG", 14, 1, False),
("STAT_TEMP", "STAT_REG", 27, 11, False),
("STAT_PRESS", "STAT_REG", 31, 4, False),
00161
00162
00163
00164
```

```
("FPA_IF_TO", "STAT_REG2", 0, 1, False),
("SRAM_RO_TO", "STAT_REG2", 1, 1, False),
("PIXELRD_TOUT_ERR", "STAT_REG2", 2, 1, False),
("UART_TX_TO_RST", "STAT_REG2", 3, 1, False),
("UART_RX_TO_RST", "STAT_REG2", 4, 1, False),
("SENSOR_POSN", "SENSOR_VOLT_STAT", 0, 1, False),
("SENSOR_NEGP", "SENSOR_VOLT_STAT", 1, 1, False),
("ICARUS_DET", "SENSOR_VOLT_STAT", 2, 1, False),
("DAPPALUS_DET", "SENSOR_VOLT_STAT", 3, 1, False),
00165
00166
00167
00168
00169
00170
00171
00172
00173
                   ("DAEDALUS_DET", "SENSOR_VOLT_STAT", 3, 1, False),
                   ("HORUS_DET", "SENSOR_VOLT_STAT", 4, 1, False),
00174
00175
                   ("SENSOR_POWER", "SENSOR_VOLT_STAT", 5, 1, False),
                   ("FIT_COUNT", "DIAG_CNTR_VAL_0", 31, 16, False), ("SRT_COUNT", "DIAG_CNTR_VAL_0", 7, 8, False), ("UTTR_COUNT", "DIAG_CNTR_VAL_1", 31, 16, False), ("UTTR_COUNT", "DIAG_CNTR_VAL_1", 15, 16, False),
00176
00177
00178
00179
00180
                   # monitor ADC channels defined here - the poll period will need to be set during
00181
                        camera initialization (x98)
                   ("MON_CH2", "ADCS_DATA_1", 15, 16, False), ("MON_CH3", "ADCS_DATA_1", 31, 16, False),
00182
00183
                   ("MON_CH4", "ADC5_DATA_2", 15, 16, False), ("MON_CH5", "ADC5_DATA_2", 31, 16, False), ("MON_CH6", "ADC5_DATA_3", 15, 16, False),
00184
00185
00186
                   ("MON_CH7", "ADC5_DATA_3", 31, 16, False),
("MON_CH8", "ADC5_DATA_4", 15, 16, False),
00187
00188
00189
                   ("MON_VRST", "ADC5_DATA_4", 31, 16, False),
00190
             1
00191
00192
             def __init__(self, camassem):
00193
                   self.ca = camassem
                  self.logcrit = self.ca.logcritbase + "[LLNL_v1] "
00194
00195
                   self.logerr = self.ca.logerrbase + "[LLNL_v1]
                  self.logwarn = self.ca.logwarnbase + "[LLNL_v1] "
00196
                   self.loginfo = self.ca.loginfobase + "[LLNL_v1] "
00197
                  self.logdebug = self.ca.logdebugbase + "[LLNL_v1] "
00198
                   logging.info(self.loginfo + "initializing board object")
00199
00200
                  self.VREF = 2.5 # default
                   self.ADC5_mult = 2 # i.e., monmax = 2 * VREF
00202
                   # False => monitor range runs 0 to monmax, True => +/- monmax
                  self.ADC5_bipolar = True
self.rs422_baud = 921600
00204
00205
                  self.rs422\_cmd\_wait = 0.3
00206
                  fpgaNum_pkt = Packet(cmd="1", addr=self.registers["FPGA_NUM"])
fpgaRev_pkt = Packet(cmd="1", addr=self.registers["FPGA_REV"])
00207
00208
00209
00210
                      _ = self.ca.sendCMD(fpgaNum_pkt) # dummy duplicate call
00211
                   err, rval = self.ca.sendCMD(fpgaNum_pkt)
00212
                   self.ca.FPGANum = rval[8:16]
00213
00214
                   err, rval = self.ca.sendCMD(fpgaRev_pkt)
00215
                  self.ca.FPGAVersion = rval[8:16]
00216
00217
                   # map channels to signal names for abstraction at the camera assembler level;
00218
                        each requires a corresponding entry in 'subregisters'
00219
00220
                   self.icarus_subreg_aliases = OrderedDict(
00221
                       {
00222
                              "COL_BOT_IBIAS_IN": "POT1",
                             "HST_A_PDELAY": "POT2",
"HST_B_NDELAY": "POT3",
00223
00224
                              "HST_RO_IBIAS": "POT4"
00225
00226
                              "HST_OSC_VREF_IN": "POT5",
00227
                              "HST_B_PDELAY": "POT6",
                              "HST_OSC_CTL": "POT7"
00228
                              "HST_A_NDELAY": "POT8",
                              "COL_TOP_IBIAS_IN": "POT9",
00230
                              "HST_OSC_R_BIAS": "POT10",
00231
                              "VAB": "POT11",
00232
00233
                              "HST_RO_NC_IBIAS": "POT12",
                              "VRST": "POT13",
00234
                              "MON_HST_A_PDELAY": "MON_CH2",
00235
00236
                              "MON_HST_B_NDELAY": "MON_CH3",
                              "MON_HST_RO_IBIAS": "MON_CH4",
00237
00238
                              "MON_HST_OSC_VREF_IN": "MON_CH5",
                              "MON_HST_B_PDELAY": "MON_CH6",
00239
00240
                              "MON_HST_OSC_CTL": "MON_CH7",
                              "MON_HST_A_NDELAY": "MON_CH8",
00241
00242
00243
00244
                  # Read-only; identifies controls corresponding to monitors
00245
                  self.icarus_monitor_controls = OrderedDict(
```

7.12 LLNL v1.py 241

```
00246
00247
                       "MON_CH2": "POT2",
00248
                       "MON_CH3": "POT3",
                       "MON_CH4": "POT4",
00249
00250
                       "MON_CH5": "POT5",
                       "MON_CH6": "POT6",
00251
00252
                       "MON_CH7": "POT7",
00253
                       "MON_CH8": "POT8",
00254
                       \# Note: VRST is not measured across the pot; it will read a voltage
00255
                          approximately 1 Volt lower than pot13's actual output
00256
                       "MON_VRST": "POT13",
00257
00258
              )
00260
              self.daedalus_subreg_aliases = OrderedDict(
00261
00262
                       "HST_OSC_CTL": "POT4",
                       "HST_RO_NC_IBIAS": "POT5",
00263
                       "HST_OSC_VREF_IN": "POT6",
00264
                       "VAB": "POT11",
00265
                       "MON_TSENSEOUT": "MON_CH2",
00266
00267
                       "MON_BGREF": "MON_CH3",
00268
                       "MON_HST_OSC_CTL": "MON_CH4",
00269
                       "MON_HST_RO_NC_IBIAS": "MON_CH5",
                       "MON_HST_OSC_VREF_IN": "MON_CH6",
"MON_COL_TST_IN": "MON_CH7",
00270
00271
00272
                       "MON_HST_OSC_PBIAS_PAD": "MON_CH8",
00273
00274
00275
               # Read-only; identifies controls corresponding to monitors
00276
              self.daedalus_monitor_controls = OrderedDict(
00277
                  {
                       "MON_CH4": "POT4",
00278
                       "MON_CH4: POT4,
"MON_CH5": "POT5",
"MON_CH6": "POT6",
00279
00280
                       \# Note: VRST is not measured across the pot; it will read a voltage
00281
00282
                          lower than pot13's actual output
                       "MON_VRST": "POT13",
00283
00284
                   }
00285
              )
00286
00287
              self.subreglist = []
00288
              for s in self.subregisters:
00289
                  self.subreglist.append(s[0].upper())
00290
                   sr = SubRegister(
00291
                       self.
00292
                       name=s[0].upper(),
00293
                       register=s[1].upper(),
00294
                       start_bit=s[2],
00295
                       width=s[3],
00296
                       writable=s[4],
00297
00298
                   setattr(self, s[0].upper(), sr)
00299
00300
               # set voltage ranges for all pots
00301
               for n in range(1, 13):
00302
                  potname = "POT" + str(n)
00303
                   potobj = getattr(self, potname)
00304
                   potobj.minV = 0
00305
                   potobj.maxV = 3.3
00306
                   # resolution is approximately .0129 V / LSB
00307
                   potobj.resolution = (1.0 * potobj.maxV - potobj.minV) / potobj.max_value
00308
               self.POT13.minV = 0
              self.POT13.maxV = 3.96
00309
               # POT13 resolution is approximately .0155 V / LSB
00310
00311
              self.POT13.resolution = (
00312
                  1.0 * self.POT13.maxV - self.POT13.minV
              ) / self.POT13.max_value
00313
00314
00315
          def initBoard(self):
00316
00317
              Register and reset board, set up firmware for sensor
00318
00319
              Returns:
               tuple (error string, response string) from final control message
00320
00321
              logging.info(self.loginfo + "initBoard LLNLv1")
00322
00323
              control_messages = [("LED_EN", "1")]
00324
00325
              self.clearStatus()
00326
              self.configADCs()
```

```
00327
00328
                  err, resp = self.ca.getSubregister("ADC5_VREF3")
00329
                       logging.error(self.logerr + "unable to read 'ADC5_VREF3'")
00330
00331
                  if int(resp, 2): # check to see if Vref is 3 or 2.5 volts
00332
                       vrefmax = 3.0
00333
                  else:
00334
                     vrefmax = 2.5
00335
                  err, resp = self.ca.getSubregister("ADC5_VREF")
00336
                       logging.error(self.logerr + "unable to read 'ADC5_VREF'")
00337
                  self.VREF = vrefmax * int(resp, 2) / 1024.0
00338
00339
                  err, multmask = self.ca.getSubregister("ADC5_MULT")
00340
00341
                       logging.error(self.logerr + "unable to read 'ADC5_MULT'")
00342
                  if multmask[0] and multmask[1] and multmask[3] and multmask[5]:
00343
                       self.ADC5_mult = 2
00344
                  elif not (multmask[0] or multmask[1] or multmask[3] or multmask[5]):
00345
                      self.ADC5_mult = 4
00346
00347
                       logging.error(self.logerr + "inconsistent mode settings on ADC5")
00348
                  return self.ca.submitMessages(control_messages, " initBoard: ")
00349
00350
            def initPots(self):
00351
00352
                  Configure default pot settings before image acquisition
00353
00354
                  Returns:
                  tuple (error string, response string) from final control message
00355
00356
                  logging.info(self.loginfo + "initPots")
00357
                  if self.ca.sensorname == "icarus" or self.ca.sensorname == "icarus2":
00358
                      err0, _ = self.ca.setPotV("HST_A_PDELAY", 0, errflag=True)
err1, _ = self.ca.setPotV("HST_B_NDELAY", 3.3, errflag=True)
err2, _ = self.ca.setPotV("HST_RO_IBIAS", 2.5, tune=True, errflag=True)
err3, _ = self.ca.setPotV("HST_OSC_VREF_IN", 2.9, tune=True, errflag=True)
00359
00360
00361
00362
                      err4, _ = self.ca.setPotV("HST_DELAY", 0, errflag=True)
err5, _ = self.ca.setPotV("HST_OSC_CTL", 1.45, tune=True, errflag=True)
err6, _ = self.ca.setPotV("HST_A_NDELAY", 3.3, errflag=True)
err7, _ = self.ca.setPotV("VAB", 0.5, errflag=True)
00363
00364
00365
00366
                       err8, _ = self.ca.setPotV("HST_RO_NC_IBIAS", 2.5, errflag=True)
00367
                      err9, _ = self.ca.setPotV("VRST", 0.3, tune=True, errflag=True)
err = err0 + err1 + err2 + err3 + err4 + err5 + err6 + err7 + err8 + err9
00368
00369
00370
                  else: # Daedalus
                      err0, _ = self.ca.setPotV("HST_OSC_CTL", 1.0, tune=True, errflag=True)
err1, _ = self.ca.setPotV("HST_RO_NC_IBIAS", 1.0, errflag=True)
err2, _ = self.ca.setPotV("HST_OSC_VREF_IN", 1.0, tune=True, errflag=True)
err3, _ = self.ca.setPotV("VAB", 0.5, errflag=True)
00371
00372
00373
00374
00375
                       err = err0 + err1 + err2 + err3
00376
                  return err, ""
00377
00378
             def latchPots(self):
00379
00380
                  Latch pot settings into sensor
00381
00382
                  Returns:
                  tuple (error string, response string) from final control message
00383
00384
00385
                  logging.info(self.loginfo + "latchPots")
00386
00387
                  control_messages = [
                       ("POT_CTL", "00000003"),  # latches register settings for pot 1 ("POT_CTL", "00000005"),
00388
00389
                       ("POT_CTL", "00000007"),
00390
                       ("POT_CTL", "00000009"),
00391
                       ("POT_CTL", "0000000B"),
00392
                       ("POT_CTL", "0000000D"),
00393
                       ("POT_CTL", "0000000F"),
00394
                       ("POT_CTL", "00000011"),
("POT_CTL", "00000013"),
00395
00396
00397
                       ("POT_CTL", "00000015"),
                       ("POT_CTL", "00000017"),
("POT_CTL", "00000019"),
00398
00399
                       ("POT_CTL", "0000001B"),
00400
00401
00402
                  return self.ca.submitMessages(control messages, " latchPots: ")
00403
00404
             def initSensor(self):
00405
00406
                  Register sensor, set default timing settings
00407
```

7.12 LLNL v1.py 243

```
00408
                tuple (error string, response string) from final control message
00409
00410
00411
                logging.info(self.loginfo + "initSensor")
00412
                if int(self.ca.FPGANum[7]) != self.ca.sensor.fpganumID:
00413
                    logging.error(
00414
                         self.logerr + "unable to confirm sensor compatibility with FPGA"
00415
00416
                self.registers.update(self.ca.sensor.sens_registers)
00417
                self.subregisters.extend(self.ca.sensor.sens_subregisters)
00418
                for s in self.ca.sensor.sens_subregisters:
00419
                    sr = SubRegister(
00420
                         self,
00421
                         name=s[0].upper(),
00422
                         register=s[1].upper(),
00423
                         start_bit=s[2],
00424
                         width=s[3],
00425
                         writable=s[4],
00426
00427
                    setattr(self, s[0].upper(), sr)
00428
                    self.subreglist.append(s[0])
00429
                self.ca.checkSensorVoltStat()
00430
                control messages = self.ca.sensorSpecific() + [
00431
                    # ring w/caps=01, relax=00, ring w/o caps = 02
00432
                     ("OSC_SELECT", "00"),
                     ("FPA DIVCLK EN ADDR", "00000001"), # TODO Make this a subregister
00433
00434
00435
                return self.ca.submitMessages(control_messages, " initSensor: ")
00436
00437
           def configADCs(self):
00438
00439
                Sets default ADC configuration (does not latch settings)
00440
00441
                tuple (error string, response string) from final control message
00442
00443
00444
                logging.info(self.loginfo + "configADCs")
00445
00446
                control_messages = [
                    # just in case ADC_RESET was set (pull all ADCs out # of reset)
00447
                     ("ADC_RESET", "00000000"),
00448
00449
                     # workaround for uncertain behavior after previous readoff
                     ("ADC1_CONFIG_DATA", "FFFFFFFF"), ("ADC2_CONFIG_DATA", "FFFFFFFF"),
00450
00451
                     ("ADC3_CONFIG_DATA", "FFFFFFFF"),
("ADC4_CONFIG_DATA", "FFFFFFFF"),
00452
00453
00454
                     ("ADC_CTL", "FFFFFFFF"),
                    ("ADC1_CONFIG_DATA", "81A801FF"), # ext Vref 1.25V

("ADC2_CONFIG_DATA", "81A801FF"), # ext Vref 1.25V

("ADC3_CONFIG_DATA", "81A801FF"), # ext Vref 1.25V

("ADC4_CONFIG_DATA", "81A801FF"), # ext Vref 1.25V

("ADC5_CONFIG_DATA", "81A883FF"), # int Vref 2.50V
00455
00456
00457
00458
00459
00460
00461
                return self.ca.submitMessages(control_messages, " configADCs: ")
00462
00463
           def softReboot(self):
00464
00465
                Perform software reboot of board. WARNING: board reboot will likely prevent
00466
                  correct communication reponses and therefore will generate an error message
00467
00468
                tuple (error string, response string) from final control message _{\tt mum}
00469
00470
00471
                logging.info(self.loginfo + "reboot")
00472
                control_messages = [("RESET", "1")]
00473
                return self.ca.submitMessages(control_messages, " disarm: ")
00474
00475
           def disarm(self):
00476
                ....
00477
                Takes camera out of trigger wait state. Has no effect if camera is not already
00478
                 in wait state.
00479
00480
                Returns:
                tuple (error string, response string) from final control message
00481
00482
00483
                logging.info(self.loginfo + "disarm")
00484
                self.ca.clearStatus()
00485
                self.ca.armed = False
00486
                control messages = [
                    ("HW_TRIG_EN", "0"),
("SW_TRIG_EN", "0"),
00487
00488
```

```
00489
00490
              return self.ca.submitMessages(control_messages, " disarm: ")
00491
          def startCapture(self, mode="Hardware"):
00492
00493
00494
              Selects trigger mode and enables board for image capture
00495
00496
00497
                  mode: trigger mode ("hardware"|"software"|"dual|"h"|"s"|"d" , is case-
00498
                    insensitive)
00499
00500
              tuple (error string, response string) from final control message
00501
00502
00503
              logging.info(self.loginfo + "startCapture")
00504
              if self.ca.sensmanual:
                  timingReg = "MANSHUT_MODE"
00505
00506
              else:
                  timingReg = "HST_MODE"
00507
00508
00509
              if mode.upper()[0] == "S": # SOFTWARE
00510
                   trigmess = [
                       ("HW_TRIG_EN", "0"),
("SW_TRIG_EN", "1"),
00511
00512
00513
                       ("SW_TRIG_START", "1"),
00514
                  1
              else: # HARDWARE
00515
00516
                  trigmess = [
                       ("SW_TRIG_EN", "0"),
("HW_TRIG_EN", "1"),
00517
00518
00519
                  1
00520
              control_messages = [
   ("ADC_CTL", "0000001F"), # configure all ADCs
   (timingReg, "1"),
00521
00522
00523
00524
00525
00526
              control_messages.extend(trigmess)
              return self.ca.submitMessages(control_messages, " startCapture: ")
00527
00528
          def readSRAM(self):
00529
00530
              Start readoff of SRAM
00531
00532
00533
              Returns:
               tuple (error string, response string from register set) _{\mbox{\tiny NIII}}
00534
00535
              logging.info(self.loginfo + "readSRAM")
00536
00537
              control_messages = [("READ_SRAM", "1")]
00538
               return self.ca.submitMessages(control_messages, " readSRAM: ")
00539
00540
          def waitForSRAM(self, timeout):
00541
              Wait until subreg 'SRAM_READY' flag is true or timeout is exceeded;
00542
00543
                timeout = None or zero means wait indefinitely
00544
00545
00546
                  timeout - time in seconds before readoff proceeds automatically without
00547
                     waiting for SRAM_READY flag
00548
00549
              Returns:
              error string
00550
00551
00552
              logging.info(self.loginfo + "waitForSRAM")
00553
              waiting = True
              starttime = time.time()
00554
              err = ""
00555
00556
              while waiting:
00557
                 err, status = self.ca.getSubregister("SRAM_READY")
00558
                   if err:
00559
                       logging.error(
00560
                           self.logerr + "error in register read: " + err + " (waitForSRAM)"
00561
00562
                   if int(status):
00563
                       waiting = False
00564
                       logging.info(self.loginfo + "SRAM ready")
00565
                   if self.ca.abort:
00566
                       waiting = False
00567
                       logging.info(self.loginfo + "readoff aborted by user")
00568
                       self.ca.abort = False
00569
                   if timeout and time.time() - starttime > timeout:
```

7.12 LLNL v1.py 245

```
00570
                       err += self.logerr + "SRAM timeout; proceeding with download attempt"
00571
                       logging.error(err)
00572
00573
               return err
00574
00575
          def getTimer(self):
00576
00577
              Read value of on-board timer
00578
00579
              Returns:
              timer value as integer
00580
00581
00582
               logging.info(self.loginfo + "getTimer")
              err, rval = self.ca.getRegister("TIMER_VALUE")
00583
00584
               if err:
00585
                  logging.error(
                       self.logerr + "unable to retrieve timer information (getTimer), "
00586
                       'returning "0" '
00587
00588
00589
                  return 0
00590
              return int (rval, 16)
00591
00592
          def resetTimer(self):
00593
00594
              Reset on-board timer
00595
00596
              Returns:
               tuple (error string, response string from register set)
00597
00598
00599
              logging.info(self.loginfo + "resetTimer")
              control_messages = [("RESET_TIMER", "1"), ("RESET_TIMER", "0")]
return self.ca.submitMessages(control_messages, " resetTimer: ")
00600
00601
00602
00603
          def enableLED(self, status):
00604
              Enable/disable on-board LEDs
00605
00606
00607
00608
                  status: 0 for disabled, 1 for enabled
00609
00610
               tuple: (error string, response string from setSubregister() """  
00611
00612
              logging.info(self.loginfo + "enableLED")
00613
               if status:
00614
                   status = 1
00615
00616
              return self.ca.setSubregister("LED_EN", str(status))
00617
00618
          def setLED(self, LED, status):
00619
00620
              Illuminate on-board LED
00621
00622
00623
                  LED: LED number (1-8)
00624
                   status: 0 is off, 1 is on
00625
00626
              tuple: (error string, response string from setSubregister()
00627
00628
00629
              logging.info(self.loginfo + "setLED")
00630
              key = "LED" + str(LED)
00631
               return self.ca.setSubregister(key, str(status))
00632
00633
          def setPowerSave(self, status):
00634
00635
              Select powersave option
00636
00637
              Args:
00638
                  status: setting for powersave option (1 is enabled)
00639
00640
               tuple (error string, response string from setSubregister() _{\rm m\,m\,m}
00641
00642
00643
              logging.info(self.loginfo + "setPowerSave")
00644
              if status:
00645
                  status = 1
00646
              return self.ca.setSubregister("POWERSAVE", str(status))
00647
00648
          def setPPER(self, pollperiod):
00649
00650
              Set polling period for ADCs.
```

```
00651
              Args:
                  pollperiod: milliseconds, between 1 and 255, defaults to 50
00652
00653
00654
              Returns:
00655
                  tuple (error string, response string from setSubregister OR invalid time
00656
                   setting string)
00657
00658
              logging.debug(self.logdebug + "setPPER: time = " + str(pollperiod))
00659
              if pollperiod is None:
                  pollperiod = 50
00660
00661
              if not isinstance(pollperiod, int) or pollperiod < 1 or pollperiod > 255:
00662
                  err = (
00663
                      self.logerr + "invalid poll period submitted. Setting remains "
00664
                       "unchanged.
00665
00666
                  logging.error(err)
00667
                  return err, str(pollperiod)
00668
              else:
00669
                  binset = bin(pollperiod)[2:].zfill(8)
00670
                  return self.ca.setSubregister("PPER", binset)
00671
00672
          def getTemp(self, scale=None, offset=None, slope=None):
00673
00674
              Read temperature sensor
00675
              Args:
00676
                  scale: temperature scale to report (defaults to C, options are F and K)
00677
                  offset: unused
00678
                  slope: unused
00679
              Returns:
              temperature as float on given scale
00680
00681
00682
              logging.debug(self.logdebug + "getTemp: scale = " + str(scale))
              err, rval = self.ca.getRegister("TEMP_SENSE_DATA")
00683
00684
              if err:
00685
                  logging.error(
                      self.logerr + "unable to retrieve temperature information ("
'getTemp), returning "0" '
00686
00687
00688
00689
                  return 0.0
00690
              ctemp = int(rval[-3:], 16) / 16.0
00691
00692
              if scale == "K":
                  temp = ctemp + 273.15
00693
              elif scale == "F":
00694
00695
                  temp = 1.8 * ctemp + 32
00696
              else:
00697
                  temp = ctemp
00698
              return temp
00699
00700
          def getPressure(self, offset, sensitivity, units):
00701
00702
              Read pressure sensor
00703
00704
              Currently unimplemented
00705
00706
              Returns:
              0 as float
00707
00708
00709
              logging.warning(
00710
                  "WARNING: [LLNL_v1] 'getPressure' is not implemented on the LLNLv1 board"
00711
00712
              return 0.0
00713
00714
          def clearStatus(self):
00715
00716
              Check status registers to clear them
00717
00718
              Returns:
              error string
00719
00720
00721
              logging.debug(self.logdebug + "clearStatus")
00722
              err1, rval = self.ca.getRegister("STAT_REG_SRC")
              err2, rval = self.ca.getRegister("STAT_REG2_SRC")
00723
00724
              err = err1 + err2
00725
00726
                 logging.error(self.logerr + "clearStatus failed")
00727
              return err
00728
          def checkStatus(self):
00729
00730
00731
              Check status register, convert to reverse-order bit stream (i.e., bit 0 is
```

7.12 LLNL v1.py 247

```
00732
                statusbits[0])
00733
00734
              Returns:
              bit string (no '0b') in reversed order """
00735
00736
00737
              logging.debug(self.logdebug + "checkStatus")
00738
              err, rval = self.ca.getRegister("STAT_REG")
00739
              if not rval:
00740
                  logging.error(
00741
                      self.logerr + "Unable to check status register (zeroes returned)"
00742
00743
                  rval = "0"
00744
              rvalbits = bin(int(rval, 16))[2:].zfill(32)
00745
              statusbits = rvalbits[::-1]
00746
              return statusbits # TODO: add error handling
00747
00748
          def checkStatus2(self):
00749
00750
              Check second status register, convert to reverse-order bit stream (i.e., bit 0
00751
                is statusbits[0])
00752
00753
              Returns: bit string (no '0b') in reversed order
00754
00755
              logging.debug(self.logdebug + "checkStatus2")
00756
              err, rval = self.ca.getRegister("STAT_REG2")
00757
              if not rval:
00758
                  logging.error(
00759
                      self.logerr + "Unable to check status register 2 (zeroes returned)"
00760
                  rval = "0"
00761
00762
              rvalbits = bin(int(rval, 16))[2:].zfill(5)
00763
              statusbits = rvalbits[::-1]
00764
              return statusbits # TODO: add error handling
00765
00766
          def reportStatus(self):
00767
00768
              Check contents of status register, print relevant messages
00769
00770
              statusbits = self.checkStatus()
statusbits2 = self.checkStatus2()
00771
00772
              logging.info(self.loginfo + "Status report:")
00773
              if int(statusbits[0]):
                  print(self.loginfo + "Sensor read complete")
00774
00775
              if int(statusbits[1]):
00776
                  print(self.loginfo + "Coarse trigger detected")
00777
              if int(statusbits[2]):
00778
                  print(self.loginfo + "Fine trigger detected")
00779
              if int(statusbits[5]):
00780
                  print(self.loginfo + "Sensor readout in progress")
00781
              if int(statusbits[6]):
00782
                  print(self.loginfo + "Sensor readout complete")
00783
              if int(statusbits[7]):
00784
                  print(self.loginfo + "SRAM readout started")
00785
              if int(statusbits[8]):
00786
                  print(self.loginfo + "SRAM readout complete")
00787
              if int(statusbits[9]):
00788
                  print(self.loginfo + "High-speed timing configuration started")
00789
              if int(statusbits[10]):
00790
                  print(self.loginfo + "All ADCs configured")
00791
              if int(statusbits[11]):
00792
                  print(self.loginfo + "All pots configured")
00793
              if int(statusbits[13]):
00794
                  print(self.loginfo + "Timer has reset")
00795
              if int(statusbits[14]):
00796
                  print(self.loginfo + "Camera is Armed")
00797
              self.ca.sensor.reportStatusSensor(statusbits, statusbits2)
00798
              temp = int(statusbits[27:15:-1], 2) / 16.0
00799
              logging.info(
00800
                  self.loginfo + "Temperature reading: " + "{0:1.2f}".format(temp) + " C"
00801
00802
              # press = int(statusbits[:27:-1], 2)
00803
              # logging.info(self.loginfo + "Pressure reading: " + "{0:1.2f}".format(press))
00804
              if int(statusbits2[0]):
00805
                  print(self.loginfo + "FPA_IF_TO")
00806
              if int(statusbits2[1]):
00807
                  print(self.loginfo + "SRAM_RO_TO")
00808
              if int(statusbits2[2]):
00809
                  print(self.loginfo + "PixelRd Timeout Error")
00810
              if int(statusbits2[3]):
                  print(self.loginfo + "UART_TX_TO_RST")
00811
00812
              if int(statusbits2[4]):
```

```
00813
                    print(self.loginfo + "UART_RX_TO_RST")
00814
00815
           def reportEdgeDetects(self):
00816
00817
               Unimplemented
00818
00819
               logging.warning(
00820
                    self.logwarn + "'reportEdgeDetects' is not implemented on the LLNLv1 "
00821
00822
00823
00824
           def dumpStatus(self):
00825
00826
               Create dictionary of status values, DAC settings, monitor values, and register
00827
                 values
00828
00829
               WARNING: the behavior of self-resetting subregisters may be difficult to predict
00830
                 and may generate contradictory results
00831
00832
               dictionary of system diagnostic values _{\mbox{\scriptsize mum}}
00833
00834
               statusbits = self.checkStatus()
statusbits2 = self.checkStatus2()
00835
00836
00837
               temp = self.ca.getTemp()
00838
00839
               statDict = OrderedDict(
00840
                   {
                        "Temperature reading": "{0:1.2f}".format(temp) + " C", "Sensor read complete": str(statusbits[0]),
00841
00842
00843
                         "Coarse trigger detected": str(statusbits[1]),
                         "Fine trigger detected": str(statusbits[2]),
00844
                        "Sensor readout in progress": str(statusbits[5]),
00845
                         "Sensor readout complete": str(statusbits[6]),
00846
                        "SRAM readout started": str(statusbits[7]),
"SRAM readout complete": str(statusbits[8]),
00847
00848
                        "High-speed timing configured": str(statusbits[9]),
00849
                         "All ADCs configured": str(statusbits[10]),
00850
                        "All pots configured": str(statusbits[11]),
"HST_All_W_En detected": str(statusbits[12]),
00851
00852
                        "Timer has reset": str(statusbits[13]),
00853
                        "Camera is Armed": str(statusbits[14]),
00854
                        "FPA_IF_TO": str(statusbits2[0]),
00855
                        "SRAM_RO_TO": str(statusbits2[1]),
00856
                        "PixelRd Timeout Error": str(statusbits2[2]),
"UART_TX_TO_RST": str(statusbits2[3]),
00857
00858
00859
                        "UART_RX_TO_RST": str(statusbits2[4]),
00860
                    }
00861
               )
00862
00863
               POTDict = OrderedDict()
00864
               MonDict = OrderedDict()
00865
               for entry in self.subreg_aliases:
00866
                    if self.subreg_aliases[entry][0] == "P":
00867
                        val = str(round(self.ca.getPotV(entry), 3)) + " V"
00868
                        POTDict["POT_" + entry] = val
00869
00870
                        val = str(round(self.ca.getMonV(entry), 3)) + " V"
00871
                        MonDict[entry] = val
00872
00873
               regDict = OrderedDict()
00874
               for key in self.registers.keys():
00875
                    err, rval = self.ca.getRegister(key)
00876
                    regDict[key] = rval
00877
00878
               dumpDict = OrderedDict()
00879
               for x in [statDict, MonDict, POTDict, regDict]:
00880
                   dumpDict.update(x)
00881
               return dumpDict
00882
00883
00884 """
00885 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00886 LLNL-CODE-838080
00888 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00889 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00890 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00891 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00892 be made under this license.
```

# 7.13 C:/Users/hill35/git/camera\_python/nsCamera/boards/LLNL\_v4.py File Reference

#### Classes

class nsCamera.boards.LLNL v4.llnl v4

#### **Namespaces**

- namespace nsCamera
- namespace nsCamera.boards
- · namespace nsCamera.boards.LLNL v4

## 7.14 LLNL v4.py

#### Go to the documentation of this file.

```
00001 # -*- coding: utf-8 -*-
00002 """
00003 LLNLv4 board definition, including monitors, DACS, and other board-specific settings
00004
00005 Author: Jeremy Martin Hill (jerhill@llnl.gov)
00006 Author: Matthew Dayton (dayton5@llnl.gov)
00007
00008 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00009 LLNL-CODE-838080
00010
00011 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00012 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DDE)
00013 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00014 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00015 be made under this license.
00016
00017 Version: 2.1.2 (February 2025)
00018 """
00019
00020 import logging
00021 import string
00022 import time
00023 from collections import OrderedDict
00025 from nsCamera.utils.Packet import Packet
00026 from nsCamera.utils.Subregister import SubRegister
00027
00029 class llnl_v4:
00030
00031
          Livermore LLNL v4.0 board
00033
          Compatible communication protocols: RS422, GigE
          Compatible sensors: icarus, icarus2, daedalus
00035
00036
00037
          # FPGA register map - use '.upper()' on keys to ensure uppercase lookup
00038
          registers = OrderedDict(
00039
00040
                  "FPGA_NUM": "000",
                   "FPGA_REV": "001",
00041
                   "HS_TIMING_CTL": "010",
00042
                   "HS_TIMING_DATA_ALO": "013",
00043
                   "HS TIMING DATA AHI": "014",
00044
                  "HS_TIMING_DATA_BLO": "015",
00045
                   "HS_TIMING_DATA_BHI": "016",
00046
                  "SW_TRIGGER_CONTROL": "017",
"SW_COARSE_CONTROL": "01C",
00047
00048
00049
                  "STAT_REG": "024",
```

```
00050
                               "CTRL_REG": "025",
00051
                               "DAC_CTL": "026",
                               "DAC_REG_A_AND_B": "027",
00052
00053
                               "DAC_REG_C_AND_D": "028",
00054
                               "DAC_REG_E_AND_F": "029",
00055
                               "DAC_REG_G_AND_H": "02A",
00056
                               "SW_RESET": "02D",
00057
                               "HST_SETTINGS": "02E",
00058
                               "STAT_REG_SRC": "02F",
                               "STAT_REG2": "030",
                               "STAT_REG2_SRC": "031",
00060
                               "ADC_BYTECOUNTER": "032",
00061
                               "RBP_PIXEL_CNTR": "033",
"DIAG_MAX_CNT_0": "034",
00062
00063
00064
                               "DIAG_MAX_CNT_1": "035"
                               "DIAG_CNTR_VAL_0": "036",
00065
                               "DIAG_CNTR_VAL_1": "037"
00066
                               "STAT_EDGE_DETECTS": "038",
00067
                               "TRIGGER_CTL": "03A",
00068
                               "SRAM_CTL": "03B",
00069
                               "TIMER_CTL": "03C",
00070
00071
                               "TIMER_VALUE": "03D",
00072
                               "HSTALLWEN_WAIT_TIME": "03F",
00073
                               "FPA_ROW_INITIAL": "042",
00074
                               "FPA_ROW_FINAL": "043",
00075
                               "FPA_FRAME_INITIAL": "044",
                               "FPA_FRAME_FINAL": "045",
00076
                               "FPA_DIVCLK_EN_ADDR": "046",
00077
                               "FPA_OSCILLATOR_SEL_ADDR": "047",
"SUSPEND_TIME": "04D",
00078
00079
                               "FPA_INTERFACE_STATE": "04E",
"DELAY_READOFF": "04F",
00080
00081
                                "STAT_REG_SEC": "060",
00082
                                "ADC_CTL": "090",
00083
                               "ADC1_CONFIG_DATA": "091",
00084
                               "ADC2_CONFIG_DATA": "092",
00085
                               "ADC3_CONFIG_DATA": "093",
"ADC4_CONFIG_DATA": "094",
00086
00087
                               "ADC5_DATA_1": "095",
00088
                               "ADC5_DATA_2": "096",
00089
                               "ADC5_DATA_3": "097",
00090
                               "ADC5_DATA_4": "098",
"ADC6_DATA_1": "099",
00091
00092
                               "ADC6_DATA_2": "09A",
00093
                               "ADC6_DATA_3": "09B",
00094
                               "ADC6_DATA_4": "09C",
"ADC_PPER": "09D",
"ADC_RESET": "09E",
00095
00096
00097
00098
00099
00100
00101
                 subregisters = [
00102
00104
                         ("HST_MODE", "HS_TIMING_CTL", 0, 1, True),
00105
                         ("SW_TRIG_START", "SW_TRIGGER_CONTROL", 0, 1, True),
                         ("SW_COARSE_TRIGGER", "SW_COARSE_CONTROL", 0, 1, True), ("LED_EN", "CTRL_REG", 1, 1, True),
00106
00107
                        ("COLQUENCHEN", "CTRL_REG", 2, 1, True), ("POWERSAVE", "CTRL_REG", 3, 1, True), ("PDBIAS_LOW", "CTRL_REG", 6, 1, True),
00108
00109
00110
                       ("PDBIAS_LOW", "CTRL_REG", 6, 1, True),
# ("SWACK", "CTRL_REG", 10, 1, True),
("DACA", "DAC_REG_A_AND_B", 31, 16, True),
("DACB", "DAC_REG_A_AND_B", 15, 16, True),
("DACC", "DAC_REG_C_AND_D", 31, 16, True),
("DACD", "DAC_REG_C_AND_D", 15, 16, True),
("DACE", "DAC_REG_E_AND_F", 31, 16, True),
("DACF", "DAC_REG_E_AND_F", 31, 16, True),
("DACG", "DAC_REG_G_AND_H", 15, 16, True),
("DACH", "DAC_REG_G_AND_H", 31, 16, True),
("DACH", "DAC_REG_G_AND_H", 15, 16, True),
("RESET", "SW_RESET", 0, 1, True),
("HST_SW_CTL_EN", "HST_SETINGS", 0, 1, Tr
00111
00112
00113
00114
00116
00117
00118
00119
00120
                        ("RESET", "SW_RESET", 0, 1, True),
("HST_SW_CTL_EN", "HST_SETTINGS", 0, 1, True),
("SW_HSTALLWEN", "HST_SETTINGS", 1, 1, True),
("MAXERR_FIT", "DIAG_MAX_CNT_0", 31, 16, True),
("MAXERR_SRT", "DIAG_MAX_CNT_0", 7, 8, True),
("MAXERR_UTTR", "DIAG_MAX_CNT_1", 31, 16, True),
("MAXERR_UTTR", "DIAG_MAX_CNT_1", 15, 16, True),
("HM_TRIG_EN", "TRIGGER_CTL", 0, 1, True),
("SW_TRIG_EN", "TRIGGER_CTL", 2, 1, True),
("DEPADOFF_DELAY_EN", "TRIGGER_CTL", 4, 1, True),
00121
00122
00123
00124
00125
00126
00127
00128
                         ("READOFF_DELAY_EN", "TRIGGER_CTL", 4, 1, True),
00129
                         ("READ_SRAM", "SRAM_CTL", 0, 1, True),
("RESET_TIMER", "TIMER_CTL", 0, 1, True),
00130
00131
```

7.14 LLNL v4.py 251

```
00132
                          ("OSC_SELECT", "FPA_OSCILLATOR_SEL_ADDR", 1, 2, True),
                          ("PPER", "ADC_PPER", 7, 8, True),
00133
00134
00137
                          ("SRAM_READY", "STAT_REG", 0, 1, False),
                          ("STAT_COARSE", "STAT_REG", 1, 1, False), ("STAT_FINE", "STAT_REG", 2, 1, False),
00138
00139
                         ("STAT_SENSREADIP", "STAT_REG", 5, 1, False),
("STAT_SENSREADDONE", "STAT_REG", 6, 1, False),
("STAT_SRAMREADSTART", "STAT_REG", 7, 1, False),
("STAT_SRAMREADDONE", "STAT_REG", 8, 1, False),
00140
00141
00142
00143
                          ("STAT_HSTCONFIGSTART", "STAT_REG", 9, 1, False), ("STAT_ADCSCONFIGURED", "STAT_REG", 10, 1, False),
00144
00145
00146
                          ("STAT_DACSCONFIGURED", "STAT_REG", 11, 1, False),
                          ("STAT_TIMERCOUNTERRESET", "STAT_REG", 13, 1, False),
00148
                          ("STAT_HSTCONFIGDONE", "STAT_REG", 16, 1, False),
                         ("STAT_HSICONFIGDONE", "STAT_REG", 10, 1, False), ("STAT_ARMED", "STAT_REG", 14, 1, False), ("STAT_TEMP", "STAT_REG", 23, 7, False), ("STAT_PRESS", "STAT_REG", 31, 8, False), ("FPA_IF_TO", "STAT_REG2", 0, 1, False), ("SRAM_RO_TO", "STAT_REG2", 1, 1, False),
00149
00150
00151
00152
00153
                         ("SRAM_RO_TO", "STAT_REG2", 1, 1, False),

("PIXELRD_TOUT_ERR", "STAT_REG2", 2, 1, False),

("UART_TX_TO_RST", "STAT_REG2", 3, 1, False),

("UART_RX_TO_RST", "STAT_REG2", 4, 1, False),

("FIT_COUNT", "DIAG_CNTR_VAL_0", 31, 16, False),

("SRT_COUNT", "DIAG_CNTR_VAL_0", 7, 8, False),

("UTTR_COUNT", "DIAG_CNTR_VAL_1", 31, 16, False),

("URTR_COUNT", "DIAG_CNTR_VAL_1", 15, 16, False),

("URTR_COUNT", "DIAG_CNTR_VAL_1", 15, 16, False),
00154
00155
00156
00157
00158
00159
00160
                         # monitor ADC channels defined here - the poll period will need to be set
00161
                         # monitor ADC channels defined here - the post during camera initialization (x98)

("MON_CH1", "ADC5_DATA_1", 11, 12, False),

("MON_CH2", "ADC5_DATA_1", 23, 12, False),

("MON_CH3", "ADC5_DATA_2", 11, 12, False),

("MON_CH4", "ADC5_DATA_2", 23, 12, False),

("MON_CH6", "ADC5_DATA_3", 11, 12, False),

("MON_CH6", "ADC5_DATA_3", 23, 12, False),

("MON_CH7", "ADC5_DATA_4", 11, 12, False),

("MON_CH8", "ADC5_DATA_4", 11, 12, False),

("MON_CH9", "ADC6_DATA_1", 11, 12, False),

("MON_CH10", "ADC6_DATA_1", 23, 12, False),

("MON_CH11", "ADC6_DATA_2", 11, 12, False),

("MON_CH12", "ADC6_DATA_2", 12, False),
00162
00163
00164
00165
00166
00167
00168
00169
00170
00171
00172
00173
                         ("MON_CH11", "ADC6_DATA_2", 11, 12, False), ("MON_CH12", "ADC6_DATA_2", 23, 12, False), ("MON_CH13", "ADC6_DATA_3", 11, 12, False), ("MON_CH14", "ADC6_DATA_3", 23, 12, False), ("MON_CH15", "ADC6_DATA_4", 11, 12, False), ("MON_CH16", "ADC6_DATA_4", 23, 12, False),
00174
00175
00176
00177
00178
00179
00180
00181
                  def __init__(self, camassem):
00182
                         self.ca = camassem
00183
                         self.logcrit = self.ca.logcritbase + "[LLNL_v4] "
                         self.logerr = self.ca.logerrbase + "[LLNL_v4] "
00184
00185
                         self.logwarn = self.ca.logwarnbase + "[LLNL_v4] "
00186
                         self.loginfo = self.ca.loginfobase + "[LLNL_v4] "
00187
                         self.logdebug = self.ca.logdebugbase + "[LLNL_v4] "
00188
                         logging.info(self.loginfo + "Initializing board object")
00189
                         self.VREF = 3.3 # must be supplied externally for ADC128S102
00190
                         self.ADC5_mult = 1
00191
00192
                          # ADC128S102; False => monitor range runs 0 to monmax, True => +/- monmax
                         self.ADC5_bipolar = False
00193
                         self.rs422\_baud = 921600
00194
00195
                         self.rs422 cmd wait = 0.3
00196
                         fpgaNum_pkt = Packet(cmd="1", addr=self.registers["FPGA_NUM"])
fpgaRev_pkt = Packet(cmd="1", addr=self.registers["FPGA_REV"])
00197
00198
00199
                         _, _ = self.ca.sendCMD(fpgaNum_pkt) # dummy duplicate call
00200
00201
                         err, rval = self.ca.sendCMD(fpgaNum_pkt)
00202
                         self.ca.FPGANum = rval[8:16]
00203
00204
                         err, rval = self.ca.sendCMD(fpgaRev_pkt)
00205
                         self.ca.FPGAVersion = rval[8:16]
00206
                         self.defoff = 34.5  # default pressure sensor offset
self.defsens = 92.5  # default pressure sensor sensitivity
00207
00208
00209
00210
                          # TODO: move to sensor scripts?
00211
                          # map channels to signal names for abstraction at the camera assembler level;
                          # each requires a corresponding entry in 'subregisters'
00212
00213
                         self.icarus_subreg_aliases = OrderedDict(
00214
```

```
00215
                           "HST_A_PDELAY": "DACA",
00216
                           "HST_A_NDELAY": "DACB",
                           "HST_B_PDELAY": "DACC",
"HST_B_NDELAY": "DACD",
00217
00218
00219
                           "HST_RO_IBIAS": "DACE",
00220
                           "HST_RO_NC_IBIAS": "DACE",
00221
                           "HST_OSC_CTL": "DACF",
                           "VAB": "DACG",
00222
00223
                           "VRST": "DACH",
                           "MON_PRES_MINUS": "MON_CH1",
00224
00225
                           "MON_PRES_PLUS": "MON_CH2",
                           "MON_TEMP": "MON_CH3",
00226
                           "MON_COL_TOP_IBIAS_IN": "MON_CH4",
"MON_HST_OSC_R_BIAS": "MON_CH5",
00227
                           "MON_VAB": "MON_CH6",
"MON_HST_RO_IBIAS": "MON_CH7",
00229
00230
                           "MON_HST_RO_NC_IBIAS": "MON_CH7",
00231
                           "MON_VRST": "MON_CH8",
00232
                           "MON_COL_BOT_IBIAS_IN": "MON_CH9",
00233
                           "MON_HST_A_PDELAY": "MON_CH10",
00234
                           "MON_HST_B_NDELAY": "MON_CH11",
00235
00236
                           "DOSIMETER": "MON_CH12",
                           "MON_HST_OSC_VREF_IN": "MON_CH13",
"MON_HST_B_PDELAY": "MON_CH14",
00237
00238
                           "MON_HST_OSC_CTL": "MON_CH15",
"MON_HST_A_NDELAY": "MON_CH16",
00239
00240
                           "MON_CHA": "MON_CH10",
"MON_CHB": "MON_CH16",
00241
00242
                           "MON_CHC": "MON_CH14",
00243
                           "MON_CHD": "MON_CH11",
00244
                           "MON_CHE": "MON_CH7",
00245
                           "MON_CHF": "MON_CH15",
00246
                           "MON_CHG": "MON_CH6",
00247
                           "MON_CHH": "MON_CH8",
00248
00249
00250
00251
                 # Read-only; identifies controls corresponding to monitors
00252
                 self.icarus_monitor_controls = OrderedDict(
00253
                     {
00254
                           "MON CH10": "DACA",
                           "MON_CH10": "DACA",
"MON_CH16": "DACB",
"MON_CH14": "DACC",
"MON_CH11": "DACC",
"MON_CH17": "DACE",
"MON_CH15": "DACF",
"MON_CH6": "DACG",
00255
00256
00257
00258
00259
00260
                           "MON_CH8": "DACH",
00261
00262
                      }
00263
00264
                 self.daedalus_subreg_aliases = OrderedDict(
00265
                           "HST_OSC_VREF_IN": "DACC",
00266
                           "HST_OSC_CTL": "DACE",
"COL_TST_IN": "DACF",
00267
00268
00269
                           "VAB": "DACG",
"VRST": "DACH",
00270
00271
                           "MON_PRES_MINUS": "MON_CH1",
00272
                           "MON_PRES_PLUS": "MON_CH2",
00273
                           "MON_TEMP": "MON_CH3",
00274
                           "MON_VAB": "MON_CH6",
00275
                           "MON_HST_OSC_CTL": "MON_CH7",
00276
                           "MON_TSENSE_OUT": "MON_CH10",
00277
                           "MON_BGREF": "MON_CH11",
                           "DOSIMETER": "MON_CH12",
00278
00279
                           "MON_HST_RO_NC_IBIAS": "MON_CH13",
                           "MON_HST_OSC_VREF_IN": "MON_CH14",
                           "MON_COL_TST_IN": "MON_CH15",
                           "MON_HST_OSC_PBIAS_PAD": "MON_CH16",
00282
                           "MON_CHC": "MON_CH14",
00283
                           "MON_CHE": "MON_CH7",
00284
00285
                           "MON_CHF": "MON_CH15",
                           "MON_CHG": "MON_CH6",
"MON_CHH": "MON_CH8",
00286
00287
00288
                      }
00289
00290
                 # Read-only; identifies controls corresponding to monitors
                 self.daedalus_monitor_controls = OrderedDict(
00291
00292
                      {
                           "MON_CH14": "DACC",
"MON_CH7": "DACE",
"MON_CH15": "DACF",
00293
00294
00295
```

7.14 LLNL v4.py 253

```
00296
                       "MON_CH6": "DACG",
00297
                       "MON_CH8": "DACH",
00298
00299
00300
               self.subreglist = []
00301
               for s in self.subregisters:
00302
                   self.subreglist.append(s[0].upper())
00303
                   sr = SubRegister(
00304
                       self.
                       name=s[0].upper(),
00305
00306
                       register=s[1].upper(),
00307
                       start_bit=s[2],
00308
                       width=s[3],
00309
                       writable=s[4],
00310
00311
                   setattr(self, s[0].upper(), sr)
00312
00313
               # set voltage ranges for all DACs - WARNING: actual output voltage limited to
               # external supply (3.3 V)
# setpot('potx', n) will generate 3.3 V for all n > .66
00314
00315
00316
               for n in range(0, 8):
                   potname = "DAC" + string.ascii_uppercase[n]
00317
00318
                   potobj = getattr(self, potname)
00319
                   potobj.minV = 0
00320
                   potobj.maxV = 5 #
                   potobj.resolution = (
00321
                   1.0 * potobj.maxV - potobj.minV) / potobj.max_value # 76 uV / LSB
00322
00323
00324
          def initBoard(self):
00325
00326
00327
               Register and reset board, set up firmware for sensor
00328
00329
               Returns:
               tuple (error string, response string) from final control message
00330
00331
               logging.info(self.loginfo + "initBoard LLNLv4")
00332
00333
               control_messages = []
00334
               self.clearStatus()
00335
               self.configADCs()
00336
               return self.ca.submitMessages(control_messages, " initBoard: ")
00337
00338
          def initPots(self):
00339
00340
               Dummy function; initial DAC values are set by firmware at startup
00341
00342
               Returns:
               tuple (empty string, empty string)
00343
00344
00345
               logging.debug(self.logdebug + "InitPots")
00346
00347
00348
          def latchPots(self):
00349
00350
               Latch DAC settings into sensor
00351
00352
               tuple (error string, response string) from final control message
00353
00354
00355
               logging.info(self.loginfo + "latchPots")
00356
               control_messages = [
                   ("DAC_CTL", "00000001"),  # latches register settings for DACA ("DAC_CTL", "00000003"),
00357
00358
                   ("DAC_CTL", "00000005"),
00359
                   ("DAC_CTL", "00000007"),
00360
                   ("DAC_CTL", "00000009"),
00361
                   ("DAC_CTL", "0000000B"),
("DAC_CTL", "0000000D"),
00362
00363
                   ("DAC_CTL", "0000000F"),
00364
00365
00366
               return self.ca.submitMessages(control_messages, " latchPots: ")
00367
00368
          def initSensor(self):
00369
00370
               Register sensor, set default timing settings
00371
00372
               Returns:
               tuple (error string, response string) from final control message
00373
00374
00375
               logging.info(self.loginfo + "initSensor")
00376
               if int(self.ca.FPGANum[7]) != self.ca.sensor.fpganumID:
```

```
00377
                     logging.warning(
00378
                         self.logwarn + "unable to confirm sensor compatibility with FPGA"
00379
00380
                self.registers.update(self.ca.sensor.sens_registers)
00381
                self.subregisters.extend(self.ca.sensor.sens_subregisters)
00382
                for s in self.ca.sensor.sens_subregisters:
00383
                    sr = SubRegister(
00384
                         self,
00385
                         name=s[0].upper(),
00386
                         register=s[1].upper(),
00387
                         start_bit=s[2],
00388
                         width=s[3],
00389
                         writable=s[4],
00390
00391
                     setattr(self, s[0].upper(), sr)
00392
                     self.subreglist.append(s[0])
00393
                # TODO: self.ca.checkSensorVoltStat() # SENSOR_VOLT_STAT and SENSOR_VOLT_CTL are
00394
                # deactivated for v4 icarus and daedalus firmware for now, is this permanent?
00395
                control_messages = self.ca.sensorSpecific() + [
00396
                    # ring w/caps=01, relax=00, ring w/o caps = 02
                     ("OSC_SELECT", "00"),
00397
00398
                     ("FPA_DIVCLK_EN_ADDR", "00000001"),
00399
00400
                return self.ca.submitMessages(control_messages, " initSensor: ")
00401
           def configADCs(self):
00402
00403
00404
                Sets default ADC configuration (does not latch settings)
00405
00406
                tuple (error string, response string) from final control message
00407
00408
00409
                logging.info(self.loginfo + "configADCs")
00410
00411
                control messages = [
                    # just in case ADC_RESET was set on any of the ADCs (pull all ADCs out of
00412
00413
                        reset)
                     ("ADC_RESET", "00000000"),
00414
                    ( ADC_RESEIT, "UUUUUUUU"),
# workaround for uncertain behavior after previous readoff
("ADC1_CONFIG_DATA", "FFFFFFFF"),
("ADC3_CONFIG_DATA", "FFFFFFFF"),
("ADC3_CONFIG_DATA", "FFFFFFFF"),
("ADC4_CONFIG_DATA", "FFFFFFFF"),
("ADC_CTL", "FFFFFFFF"),
("ADC_CCL", "FFFFFFFFF"),
00415
00416
00417
00418
00419
00420
                     ("ADC1_CONFIG_DATA", "81A801FF"), # ext Vref 1.25V

("ADC2_CONFIG_DATA", "81A801FF"), # ext Vref 1.25V

("ADC3_CONFIG_DATA", "81A801FF"), # ext Vref 1.25V

("ADC4_CONFIG_DATA", "81A801FF"), # ext Vref 1.25V
00421
00422
00423
00424
00425
00426
                return self.ca.submitMessages(control_messages, " configADCs: ")
00427
00428
           def softReboot(self):
00429
                ....
00430
                Perform software reboot of board. WARNING: board reboot will likely prevent
00431
                 correct response and therefore will generate an error message
00432
00433
                tuple (error string, response string) from final control message
00434
00435
00436
                logging.info(self.loginfo + "reboot")
                control_messages = [("RESET", "0")]
00437
00438
                return self.ca.submitMessages(control_messages, " disarm: ")
00439
00440
           def disarm(self):
00441
00442
                Takes camera out of trigger wait state. Has no effect if camera is not in wait
00443
                  state.
00444
00445
                Returns:
                tuple (error string, response string) from final control message
00446
00447
00448
                logging.info(self.loginfo + "disarm")
00449
                self.ca.clearStatus()
00450
                self.ca.armed = False
00451
                control_messages = [
00452
                     ("HW_TRIG_EN", "0"),
                     ("SW TRIG EN", "0"),
00453
00454
00455
                self.ca.comms.skipError = False
00456
                return self.ca.submitMessages(control_messages, " disarm: ")
00457
```

7.14 LLNL v4.py 255

```
00458
          def startCapture(self, mode="Hardware"):
00459
00460
              Selects trigger mode and enables board for image capture
00461
00462
              Args:
00463
                  mode: trigger mode ("hardware"|"software"|"dual|"h"|"s"|"d" , is case-
00464
                    insensitive)
00465
00466
              tuple (error string, response string) from final control message
00467
00468
00469
              logging.info(self.loginfo + "startCapture")
00470
              if self.ca.sensmanual:
00471
                  timingReg = "MANSHUT_MODE"
00472
              else:
                  timingReg = "HST_MODE"
00473
00474
00475
              if mode.upper()[0] == "S": # SOFTWARE
00476
                  trigmess = [
00477
                      ("HW_TRIG_EN", "0"),
                      ("SW_TRIG_EN", "1")
00478
00479
                      ("SW_TRIG_START", "1"),
00480
                  1
00481
              else: # HARDWARE
00482
                  trigmess = [
                      ("SW_TRIG_EN", "0"),
00483
00484
                      ("HW_TRIG_EN", "1"),
00485
                  1
00486
00487
              control_messages = [
                  ("ADC_CTL", "0000000F"), # configure all ADCs (timingReg, "1"),
00488
00489
00490
00491
              control_messages.extend(trigmess)
00492
00493
              return self.ca.submitMessages(control_messages, " startCapture: ")
00494
00495
          def readSRAM(self):
00496
00497
              Start readoff of SRAM
00498
00499
              Returns:
              tuple (error string, response string from register set)
00500
00501
00502
              logging.info(self.loginfo + "readSRAM")
00503
              control_messages = [("READ_SRAM", "1")]
00504
              return self.ca.submitMessages(control_messages, " readSRAM: ")
00505
00506
          def waitForSRAM(self, timeout):
00507
              Wait until subreg 'SRAM_READY' flag is true or timeout is exceeded;
00508
00509
                timeout = None or zero means wait indefinitely
00510
00511
00512
                  timeout - time in seconds before readoff proceeds automatically without
00513
                    waiting for SRAM_READY flag
00514
00515
              Returns:
              error string
00516
00517
00518
              logging.info(self.loginfo + "waitForSRAM, timeout = " + str(timeout))
              waiting = True
00519
00520
              starttime = time.time()
00521
              err = ""
              while waiting:
00522
00523
                  err, status = self.ca.getSubregister("SRAM_READY")
00524
                  if err:
00525
                      err = self.logerr + "error in register read: " + err + " (waitForSRAM)"
00526
                      logging.error(err)
00527
                  if int(status):
00528
                      waiting = False
00529
                      logging.info(self.loginfo + "SRAM ready")
                  if self.ca.abort:
00530
00531
                      waiting = False
00532
                      logging.info(self.loginfo + "readoff aborted by user")
00533
                      self.ca.abort = False
                  if timeout and time.time() - starttime > timeout:
00534
                      err += self.logerr + "SRAM timeout; proceeding with download attempt"
00535
00536
                      logging.error(err)
00537
                      return err
00538
                  # Slow down for debugging (avoid thousands of messages)
```

```
00539
                   if self.ca.verbose >= 5:
00540
                      time.sleep(0.5)
00541
              return err
00542
00543
          def getTimer(self):
00544
00545
              Read value of on-board timer
00546
00547
              timer value as integer
00548
00549
00550
              logging.info(self.loginfo + "getTimer")
00551
              err, rval = self.ca.getRegister("TIMER_VALUE")
00552
00553
                  logging.error(
00554
                      self.logerr + "unable to retrieve timer information (getTimer), "
                       'returning "0" '
00555
00556
00557
                  return 0
00558
              return int (rval, 16)
00559
00560
          def resetTimer(self):
00561
00562
              Reset on-board timer
00563
00564
              Returns:
              tuple (error string, response string from register set)
00565
00566
              logging.info(self.loginfo + "resetTimer")
control_messages = [("RESET_TIMER", "1"), ("RESET_TIMER", "0")]
00567
00568
              return self.ca.submitMessages(control_messages, " resetTimer: ")
00569
00570
00571
          def enableLED(self, status):
00572
00573
              Dummy function; feature is not implemented on LLNL_V4 board
00574
00575
              tuple: dummy of (error string, response string from setSubregister()) _{\mbox{\tiny mun}}
00576
00577
              del status
00578
              return "", "0"
00579
00580
00581
          def setLED(self, LED, status):
00582
00583
              Dummy function; feature is not implemented on LLNL_V4 board
00584
00585
              tuple: dummy of (error string, response string from setSubregister())
00586
00587
              del LED, status return "", "0"
00588
00589
00590
00591
          def setPowerSave(self, status):
00592
00593
              Select powersave option
00594
00595
00596
                  status: setting for powersave option (1 is enabled)
00597
00598
              tuple (error string, response string from setSubregister())
00599
00600
00601
              if status:
00602
                  status = 1
              return self.ca.setSubregister("POWERSAVE", str(status))
00603
00604
00605
          def setPPER(self, pollperiod):
00606
00607
              Set polling period for ADCs.
00608
                  pollperiod: milliseconds, between 1 and 255; defaults to 50
00609
00610
00611
              Returns:
00612
                  tuple (error string, response string from setSubregister() OR invalid time
00613
                    setting string)
00614
00615
              if pollperiod is None:
00616
                  pollperiod = 50
              if not isinstance(pollperiod, int) or pollperiod < 1 or pollperiod > 255:
00617
00618
                  err = (
                      self.logerr + "invalid poll period submitted. Setting remains "
00619
```

7.14 LLNL v4.py 257

```
00620
                        "unchanged. "
00621
00622
                   logging.error(err)
                   return err, str(pollperiod)
00623
00624
               else:
00625
                   binset = bin(pollperiod)[2:].zfill(8)
00626
                   return self.ca.setSubregister("PPER", binset)
00627
00628
           def getTemp(self, scale=None):
00629
00630
               Read temperature sensor
00631
               Args:
00632
                   scale: temperature scale to report (defaults to C, options are F and K)
00633
00634
               Returns:
               temperature as float on given scale
00635
00636
00637
               err, rval = self.ca.getMonV("MON_TEMP", errflag=True)
00638
               if err:
00639
                   logging.error(
                       self.logerr + "unable to retrieve temperature information ("
'getTemp), returning "0" '
00640
00641
00642
00643
                   return 0.0
00644
               ctemp = rval * 1000 - 273.15
               if scale == "K":
00645
                   temp = ctemp + 273.15
00646
00647
               elif scale == "F":
00648
                   temp = 1.8 * ctemp + 32
00649
               else:
00650
                   temp = ctemp
00651
               return temp
00652
00653
           def getPressure(self, offset=None, sensitivity=None, units=None):
00654
00655
               Read pressure sensor. Uses default offset and sensitivity defined in init
00656
                 function unless alternatives are specified. NOTE: to reset defaults, reassign
00657
                 board.defoff and board.defsens explicitly
00658
00659
               Args:
00660
                   offset: non-default offset in mv/V
                   sensitivity: non-default sensitivity in mV/V/span
00661
00662
                   units: units to report pressure (defaults to Torr, options are psi, bar,
00663
                     inHg, atm)
00664
00665
               Returns:
               Pressure as float in chosen units, defaults to torm \ensuremath{{\bf n}}\ensuremath{{\bf n}}\ensuremath{{\bf n}}\ensuremath{{\bf n}}
00666
00667
00668
               if offset is None:
00669
                   offset = self.defoff
00670
               if sensitivity is None:
00671
                   sensitivity = self.defsens
00672
               if units is None:
                   units = "torr"
00673
00674
               pplus = self.ca.getMonV("MON_PRES_PLUS")
00675
               pminus = self.ca.getMonV("MON_PRES_MINUS")
00676
               delta = 1000 * (pplus - pminus)
00677
               ratio = sensitivity / 30 # nominal is 21/30
00678
               psi = (delta - offset) / ratio
00679
               if units.lower() == "psi":
00680
                   press = psi
00681
               elif units.lower() == "bar":
00682
                   press = psi / 14.504
00683
               elif units.lower() == "atm":
                   press = psi / 14.695
00684
               elif units.lower() == "inHg":
00685
                  press = psi * 2.036
00686
               else:
00687
00688
                   press = 51.715 * psi # default to Torr
00689
00690
               return press
00691
00692
           def clearStatus(self):
00693
00694
               Check status registers to clear them
00695
00696
               Returns:
               error string
00697
00698
               err1, rval = self.ca.getRegister("STAT_REG_SRC")
err2, rval = self.ca.getRegister("STAT_REG2_SRC")
00699
00700
```

```
00701
              err = err1 + err2
00702
              if err:
00703
                  logging.error(self.logerr + "clearStatus failed")
00704
00705
00706
         def checkStatus(self):
00707
00708
              Check status register, convert to reverse-order bit stream (i.e., bit 0 is
00709
                statusbits[0])
00710
00711
              Returns:
              bit string (no '0b') in reversed order
00712
00713
00714
              err, rval = self.ca.getRegister("STAT_REG")
00715
              rvalbits = bin(int(rval, 16))[2:].zfill(32)
00716
              statusbits = rvalbits[::-1]
00717
              return statusbits
00718
00719
         def checkStatus2(self):
00720
00721
              Check second status register, convert to reverse-order bit stream (i.e., bit 0
00722
                is statusbits[0])
00723
00724
              Returns: bit string (no 'Ob') in reversed order
00725
00726
              err, rval = self.ca.getRegister("STAT_REG2")
00727
              rvalbits = bin(int(rval, 16))[2:].zfill(6)
              statusbits = rvalbits[::-1]
00728
00729
              return statusbits
00730
00731
          def reportStatus(self):
00732
              Check contents of status register, print relevant messages
00733
00734
00735
              statusbits = self.checkStatus()
              statusbits2 = self.checkStatus2()
00736
00737
00738
              logging.info(self.loginfo + "Status report:")
              if int(statusbits[0]):
00739
                  print(self.loginfo + "Sensor read complete")
00740
00741
              if int(statusbits[1]):
00742
                  print(self.loginfo + "Coarse trigger detected")
              if int(statusbits[2]):
00743
00744
                  print(self.loginfo + "Fine trigger detected")
00745
              if int(statusbits[5]):
00746
                  print(self.loginfo + "Sensor readout in progress")
00747
              if int(statusbits[6]):
00748
                  print(self.loginfo + "Sensor readout complete")
00749
              if int(statusbits[7]):
00750
                  print(self.loginfo + "SRAM readout started")
              if int(statusbits[8]):
00751
00752
                  print(self.loginfo + "SRAM readout complete")
00753
              if int(statusbits[9]):
00754
                  print(self.loginfo + "High-speed timing configuration started")
00755
              if int(statusbits[10]):
00756
                  print(self.loginfo + "All ADCs configured")
              if int(statusbits[11]):
00757
00758
                  print(self.loginfo + "All DACs configured")
00759
              if int(statusbits[13]):
00760
                  print(self.loginfo + "Timer has reset")
00761
              if int(statusbits[14]):
00762
                  print(self.loginfo + "Camera is Armed")
00763
              if int(statusbits[16]):
00764
                  print(self.loginfo + "High-speed timing configuration complete")
00765
              self.ca.sensor.reportStatusSensor(statusbits, statusbits2)
00766
              temp = int(statusbits[23:16:-1], 2) * 3.3 * 1000 / 4096
00767
              logging.info(
00768
                  self.loginfo + "Temperature reading: " + "{0:1.2f}".format(temp) + " C"
00769
00770
              press = int(statusbits[:23:-1], 2) \star 3.3 \star 1000 / 4096
00771
              logging.info(
00772
                  self.loginfo
00773
                  + "Pressure sensor reading: "
00774
                  + "{0:1.2f}".format(press)
00775
                  + " mV"
00776
00777
              if int(statusbits2[0]):
00778
                  print(self.loginfo + "FPA_IF_TO")
00779
              if int(statusbits2[1]):
00780
                  print(self.loginfo + "SRAM_RO_TO")
00781
              if int(statusbits2[2]):
```

7.14 LLNL v4.py 259

```
00782
                   print(self.loginfo + "PixelRd Timeout Error")
00783
               if int(statusbits2[3]):
00784
                   print(self.loginfo + "UART_TX_TO_RST")
               if int(statusbits2[4]):
00785
                   print(self.loginfo + "UART_RX_TO_RST")
00786
00787
00788
          def reportEdgeDetects(self):
00789
00790
               Report edge detects
00791
00792
               err, rval = self.ca.getRegister("STAT_EDGE_DETECTS")
00793
               # shift to left to fake missing edge detect
00794
               edgebits = bin(int(rval, 16) « 1)[2:].zfill(32)
00795
               # reverse to get order matching assignment
00796
               bitsrev = edgebits[::-1]
00797
               detdict = {}
00798
               bitidx = 0
00799
               for frame in range(4):
                   for vert in ("TOP", "BOT"):
00800
                       for edge in range(1, 3):
00801
                           for hor in ("A", "B"):
00802
00803
                                detname = (
00804
                                    "W"
00805
                                    + str(frame)
00806
                                    + "_"
00807
                                    + vert
00808
00809
                                    + hor
00810
                                    + " EDGE"
00811
                                    + str(edge)
00812
00813
                                detdict[detname] = bitsrev[bitidx]
00814
                                bitidx += 1
00815
               # remove faked detect
               del detdict["W0_TOP_A_EDGE1"]
00816
               logging.info(self.loginfo + "Edge detect report:")
00817
00818
               for key, val in detdict.items():
                   logging.info(self.loginfo + key + ": " + val)
00819
00820
00821
          def dumpStatus(self):
00822
00823
               Create dictionary of status values, DAC settings, monitor values, and register
00824
                values.
00825
00826
               Returns:
               dictionary of system diagnostic values
00827
00828
00829
               statusbits = self.checkStatus()
00830
               statusbits2 = self.checkStatus2()
00831
               temp = int(statusbits[23:16:-1], 2) * 3.3 * 1000 / 4096 press = int(statusbits[:23:-1], 2) * 3.3 * 1000 / 4096
00832
00833
00834
00835
               statDict = OrderedDict(
00836
                  {
00837
                       "Temperature sensor reading (Deg C)": "\{0:1.2f\}".format(temp),
00838
                       "Pressure reading (Torr)": str(round(self.ca.getPressure(), 3)),
                       "Pressure sensor reading (mV)": "{0:1.2f}".format(press),
00839
00840
                       "Sensor read complete": str(statusbits[0]),
00841
                       "Coarse trigger detected": str(statusbits[1]),
00842
                       "Fine trigger detected": str(statusbits[2]),
00843
                       "Sensor readout in progress": str(statusbits[5]),
00844
                       "Sensor readout complete": str(statusbits[6]),
00845
                       "SRAM readout started": str(statusbits[7]),
00846
                       "SRAM readout complete": str(statusbits[8]),
                       "High-speed timing configured": str(statusbits[9]),
00847
00848
                       "All ADCs configured": str(statusbits[10]),
                       "All DACs configured": str(statusbits[11]),
00849
                       "Timer has reset": str(statusbits[13]),
00850
                       "Camera is Armed": str(statusbits[14]),
00851
00852
                       "FPA_IF_TO": str(statusbits2[0]),
00853
                       "SRAM_RO_TO": str(statusbits2[1]),
00854
                       "PixelRd Timeout Error": str(statusbits2[2]),
                       "UART_TX_TO_RST": str(statusbits2[3]),
"UART_RX_TO_RST": str(statusbits2[4]),
00855
00856
                       "PDBIAS Unready": str(statusbits2[5]),
00857
00858
00859
               )
00860
00861
               if self.ca.sensorname == "icarus" or self.ca.sensorname == "icarus2":
00862
                   senslabs = {
```

```
00863
                       3: "W3_Top_A_Edge1 detected",
00864
                       4: "W3_Top_B_Edge1 detected",
00865
                       12: "HST_All_W_En detected",
00866
00867
              else:
00868
                  senslabs = {
00869
                      3: "RSLROWOUTA",
00870
                       4: "RSLROWOUTB",
00871
                       12: "RSLNALLWENB",
00872
                       15: "RSLNALLWENA",
00873
                       16: "Config HST is done",
00874
00875
              sensDict = {senslabs[x]: str(statusbits[x]) for x in senslabs.keys()}
00876
00877
              DACDict = OrderedDict()
00878
              MonDict = OrderedDict()
00879
              for entry in self.subreq_aliases:
                   if self.subreg_aliases[entry][0] == "D":
00880
00881
                       val = str(round(self.ca.getPotV(entry), 3))
                      DACDict["DAC_" + entry] = val
00882
00883
                  else:
00884
                      val = str(round(self.ca.getMonV(entry), 3))
00885
                      MonDict[entry] = val
00886
00887
              regDict = OrderedDict()
00888
              for key in self.registers.keys():
                  # Load in all registers except for the read-clear status registers.
if key == "STAT_REG_SRC" or key == "STAT_REG2_SRC":
00889
00890
00891
00892
                  else:
00893
                      err, rval = self.ca.getRegister(key)
00894
                      regDict[key] = rval
00895
              dumpDict = OrderedDict()
00896
              for x in [statDict, sensDict, MonDict, DACDict, regDict]:
00897
00898
                 dumpDict.update(x)
00899
              return dumpDict
00900
00901
00902 """
00903 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00904 LLNL-CODE-838080
00905
00906 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00907 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DDE)
00908 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00909 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00910 be made under this license.
00911 """
```

# 7.15 C:/Users/hill35/git/camera\_python/nsCamera/CameraAssembler.py File Reference

#### Classes

· class nsCamera.CameraAssembler.CameraAssembler

#### **Namespaces**

- namespace nsCamera
- namespace nsCamera.CameraAssembler

# 7.16 CameraAssembler.py

#### Go to the documentation of this file.

```
00001 # -*- coding: utf-8 00002 """
00003 CameraAssembler assembles the separate camera parts into a camera object. This object
00004 controls a combination of three components:
00006 1. board : FPGA board -- LLNL_V1, LLNL_V4
00007 2. comms: communication interface -- GigE, RS422
00008 3. sensor : sensor type -- icarus, icarus2, daedalus
00010 Author: Jeremy Martin Hill (jerhill@llnl.gov)
00011 Author: Matthew Dayton (dayton5@llnl.gov)
00013 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00014 LLNL-CODE-838080
00015
00016 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00017 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00018 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00019 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00020 be made under this license.
00021
00022 Version: 2.1.2 (February 2025)
00023 """
00024
00025 from __future__ import absolute_import
00026
00027 import importlib
00028 import inspect
00029 import logging
00030 import os
00031 import platform
00032 import socket
00033 import sys
00034 import time
00035 import h5py
00036 from datetime import datetime
00037
00038 import numpy as np
00039
00040 from nsCamera.utils.misc import (
00041
        bytes2str,
00042
          checkCRC,
00043
         flattenlist,
00044
         generateFrames,
00045
          getEnter,
00046
         partition,
00047
         plotFrames,
00048
          saveTiffs,
00049
          str2bytes,
00050
          str2nparray,
00051)
00052
00053 from nsCamera.utils.Packet import Packet
00054
00055 # TODO: move to Sphinx documentation
00056 # TODO: add pytest and tox scripts
00058
00059 class CameraAssembler:
          Code to assemble correct code to manage FPGA, frame grabber, and sensor
00063
          Exposed methods:
00064
              initialize() - initializes board registers and pots, sets up sensor
00065
              reinitialize() - initialize board and sensors, restore last known timer settings
              reboot() - perform software reset of board and reinitialize
00067
              getBoardInfo() - parses FPGA_NUM register to retrieve board description
00068
              getRegister(regname) - retrieves contents of named register
00069
              setRegister (regname, string) - sets named register to given value
00070
              getSubregister(subregname) - return substring of register identified in board
00071
               attribute 'subregname'
00072
              setSubregister(subregname, valstring) - replace substring of register identified
               in board attribute 'subregname' with 'valstring'
00073
00074
              submitMessages(messages) - set registers or subregisters based on list of
00075
               destination/payload tuples
00076
              getPot(potname) - returns float (0 < value < 1) corresponding to integer stored
```

```
00077
                in pot or monitor 'potname'
00078
               setPot(potname, value) - 0 < value < 1; sets named pot to fixed-point number =
00079
                 'value' * (maximum pot value)
08000
               getPotV(potname) - returns voltage setting of 'potname'
00081
               setPotV(potname, voltage) - sets named pot to voltage
00082
               getMonV(monname) - returns voltage read by monitor 'monname' (or monitor
00083
                associated with given potname)
00084
               readImgs() - calls arm() and readoff() functions
00085
               saveFrames(frames) - save image object as one file
               saveNumpys(frames) - save individual frames as numpy data files
00086
00087
               dumpNumpy(datastream) - save datastream string to numpy file
00088
              powerCheck(delta) - check that board power has not failed
00089
               printBoardInfo() - print board information derived from FPGA_NUM register
               dumpRegisters() - return contents of all board registers
00090
00091
               dumpSubregisters() - return contents of all named board subregisters
00092
               setFrames(min, max) - select subset of frames for readoff
00093
               setRows(min, max, padToFull) - select subset of rows for readoff
               abortReadoff() - cancel readoff in wait-for-SRAM loop
00094
               batchAquire() - fast acquire a finite series of images
00095
00096
               loadTextFrames() - load data sets previously saved as text and convert to frames
00097
00098
          Includes aliases to board- and sensor- specific methods:
00099
              Board methods
00100
                  disarm() - take camera out of waiting-for-trigger state
00101
                   clearStatus() - clear contents of status registers
                   checkStatus() - print contents of status register as reversed bit string
00102
                   checkStatus2() - print contents of status register 2 as reversed bit string reportStatus() - print report on contents of status registers
00103
00104
00105
                   resetTimer() - reset on-board timer
00106
                   getTimer() - read on-board timer
00107
                   enableLED(status) - enable (default) or disable (status = 0) on-board LEDs
00108
                   setLED(LED#, status) - turn LED on (default) or off (status = 0)
00109
                   setPowerSave(status) - turn powersave functionality on (default) or off
00110
                    (status = 0)
                   getTemp() - read on-board temperature sensor
00111
00112
                   getPressure() - read on-board pressure sensor
00113
                   dumpStatus() - generate dictionary of status, register, and subregister
00114
                    contents
00115
              Sensor methods
00116
                  \verb|checkSensorVoltStat()| - \verb|checks| that jumper settings| match sensor selection|\\
00117
                   setTiming(side, sequencetuple, delay) - configure high-speed timing
00118
                   setArbTiming(side, sequencelist) - configure arbitrary high-speed timing
00119
00120
                   getTiming(side) - returns high speed timing settings from registers
                   setManualTiming() - configures manual shutter timing
getManualTiming() - returns manual shutter settings from registers
00121
00122
00123
                   selectOscillator(osc) - select timing oscillator
00124
                   setInterlacing(ifactor) - sets interlacing factor
00125
                   setHighFullWell(flag) - controls High Full Well mode
                   setZeroDeadTime(flag, side) - controls Zero Dead Time mode
setTriggerDelay(delayblocks) - sets trigger delay
00126
00127
00128
              Comms methods
00129
                   sendCMD(pkt) - sends packet object via serial port
00130
                   arm() - configures software buffers & arms camera
00131
                   readFrames() - waits for data ready flag, then downloads image data
00132
                   readoff() - waits for data ready flag, then downloads image data
00133
                   closeDevice() - disconnect interface and release resources
00134
              Miscellaneous functions (bare functions that can be called as methods)
00135
                   saveTiffs(frames) - save individual frames as tiffs
00136
                   plotFrames(frames) - plot individual frames as tiffs
00137
00138
          Informational class variables:
00139
               version - nsCamera software version
00140
               FPGAVersion - firmware version (date)
00141
               FPGANum - firmware implementation identifier
00142
               FPGAboardtype - FPGA self-identified board type (should match 'boardname')
               FPGArad = Flag indicating radiation-tolerant FPGA build
00143
00144
              FPGAsensor = FPGA self-identified sensor family (should correspond to
00145
                 'sensorname')
00146
               FPGAinterfaces = FPGA self-identified interfaces (list should include
00147
                 'commname')
00148
              FPGAinvalid = flag indicating invalid FPGA information in register
00149
00150
00151
          def __init__(
00152
              self.
              boardname="11n1 v4",
00153
              commname="GigE",
00154
00155
              sensorname="icarus2",
00156
              verbose=4.
00157
              port=None,
```

```
00158
              ip=None,
00159
              logfile=None,
00160
              logtag=None,
00161
              timeout=30,
00162
          ):
00163
00164
              Args:
00165
                  boardname: name of FPGA board: llnl_v1, llnl_v4
00166
                  commname: name of communication interface: rs422, gige
00167
                   sensorname: name of sensor: icarus, icarus2, daedalus
00168
                   verbose: optional, sets logging level
00169
                      0: print no logging messages
00170
                       1: print CRITICAL logging messages (camera will not operate, e.g.,
00171
                        unable to connect to board)
00172
                       2: print ERROR logging messages (camera will not operate as directed,
00173
                        e.g., an attempt to set the timing mode has failed, but the camera
00174
                         is still operational)
00175
                       3: print WARNING logging messages (camera will operate as directed, but
                         perhaps not as expected, e.g., ca.setTiming('A', (9, 8), 1) may be
00176
                         programmed correctly, but the actual timing generated by the board
00177
00178
                         will be {1} [9, 8, 9, 14, 9, 8, 9])
00179
                       4: print INFO logging messages (operational messages from ordinary
00180
                         camera operation)
00181
                       5. print DEBUG logging messages (detailed messages describing specific
00182
                         operations and messages)
00183
                  port: optional integer
00184
                       When using RS422, this preselects the comport for RS422 and bypasses
00185
                         port search
00186
                       When using \mathrm{GigE}, this preselects the OrangeTree control port for \mathrm{GigE}
00187
                         (ignored if an ip parameter is not also provided)
                  ip: optional string (e.g., '192.168.1.100')

GigE: bypasses network search and selects particular OrangeTree board -
00188
00189
00190
                         required for some operating systems
                  logfile: optional string, name of file to divert console output
00191
00192
                  timeout: timeout in seconds for connecting using Gigabit Ethernet
00193
              self.version = "2.1.2"
00194
00195
              self.currtime = 0
              self.oldtime = 0
00196
00197
              self.trigtime = []
00198
              self.waited = []
00199
              self.read = []
00200
              self.unstringed = []
00201
              self.parsedtime = []
00202
              self.savetime = []
00203
              self.cycle = []
00204
              self.boardname = boardname.lower()
00205
              self.timeout = timeout
00206
               # TODO: parse boardname, etc. in separate method
              if self.boardname in ["llnlv1", "v1", "1", 1]:
    self.boardname = "llnl_v1"
00207
00208
00209
              if self.boardname in ["llnlv4", "v4", "4", 4]:
00210
                  self.boardname = "llnl_v4"
00211
              self.commname = commname.lower()
00212
              if self.commname[0] == "g" or self.commname[0] == "e":
00213
                  self.commname = "gige"
00214
              if self.commname[0] == "r":
00215
                  self.commname = "rs422"
              self.sensorname = sensorname.lower()
00216
00217
              if self.sensorname in ["i1", "ic1", "icarus1"]:
                  self.sensorname = "icarus"
00218
00219
              if self.sensorname in ["i2", "ic2"]:
00220
                  self.sensorname = "icarus2"
              if self.sensorname == "d":
00221
                  self.sensorname = "daedalus"
00223
              self.verbose = int(verbose)
00224
              self.port = port
              self.python, self.pyth1, self.pyth2, _, _ = sys.version_info
00225
              self.PY3 = self.python >= 3
00226
00227
              self.platform = platform.system()
              self.arch, _ = platform.architecture()
00228
00229
00230
              self.FPGAVersion = ""
00231
              self.FPGANum = ""
00232
              # FPGA information here and below populated during initialization using
00233
                  getBoardInfo
              self.FPGAboardtype = ""
00234
00235
              self.FPGArad = False
00236
              self.FPGAsensor = ""
00237
              self.FPGAinterfaces = []
00238
```

```
00239
               # indicates invalid FPGA information in register# (0x80000001 accepted as valid)
00240
               self.FPGAinvalid = False
00241
00242
               self.iplist = None
00243
               self.packageroot = os.path.dirname(inspect.getfile(CameraAssembler))
00244
               self.armed = False
00245
00246
               # only one of these collections (senstiming, sensmanual) should be nonempty at
00247
                   any given time
               self.senstiming = {} # preserve HST setting against possible power failure
00248
00249
               self.sensmanual = [] # preserve manual timing
00250
               self.inittime = 0
               self.padToFull = False
00251
00252
               self.abort = False
00253
00254
               self.verbmap = {
00255
                   0: 99,
00256
                   1: logging.CRITICAL,
00257
                   2: logging.ERROR,
00258
                   3: logging.WARNING,
00259
                   4: logging.INFO,
00260
                   5: logging.DEBUG,
00261
00262
               if logtag is None:
00263
                   logtag = ""
               self.logtag = logtag
00264
00265
00266
               self.logcritbase = "CRITICAL {logtag}: ".format(logtag=logtag)
               self.logerrbase = "ERROR {logtag}: ".format(logtag=logtag)
00267
               self.logwarnbase = "WARNING {logtag}: ".format(logtag-logtag)
self.loginfobase = "INFO {logtag}: ".format(logtag=logtag)
00268
00269
               self.logdebugbase = "DEBUG {logtag}: ".format(logtag=logtag)
00270
00271
               self.logcrit = "{lb}[CA]".format(lb=self.logcritbase)
self.logerr = "{lb}[CA]".format(lb=self.logerrbase)
00272
00273
               self.logwarn = "{lb}[CA]".format(lb=self.logwarnbase)
00274
               self.loginfo = "{lb}[CA]".format(lb=self.loginfobase)
self.logdebug = "{lb}[CA]".format(lb=self.logdebugbase)
00275
00276
00277
00278
               self.verblevel = self.verbmap.get(verbose, 5) # defaults to 5 for invalid entry
00279
00280
               if logfile:
                   logging.basicConfig(format="%(message)s", filename=logfile)
00281
00282
               else:
00283
                   logging.basicConfig(format="%(message)s")
00284
               logging.getLogger().setLevel(self.verblevel)
00285
               logging.getLogger("matplotlib.font_manager").disabled = True
00286
               logging.debug(
00287
                   "{logdebug}CameraAssembler: boardname = {boardname}; commname = {commname};"
00288
                   " sensorname = {sensorname}; verbose = {verbose}; port = {port}; ip = {ip};"
                   " logfile = {logfile}; logtag = {logtag}".format(
00289
00290
                        logdebug=self.logdebug,
00291
                        boardname=boardname,
00292
                        commname=commname,
00293
                        sensorname=sensorname,
00294
                        verbose=verbose,
00295
                        port=port,
00296
                        ip=ip,
00297
                        logfile=logfile,
00298
                        logtag=logtag,
00299
00300
               )
00301
00302
               if ip:
00303
                   try:
00304
                       iphex = socket.inet_aton(ip)
                   except socket.error:
00305
00306
                        logging.critical(
00307
                            "{logcrit}CameraAssembler: invalid IP provided".format(
                                logcrit=self.logcrit
00308
00309
                            )
00310
                        )
00311
                        sys.exit(1)
                   ipnum = [0, 0, 0, 0]
00312
                   for i in range(4):
00313
00314
                        if self.PY3:
00315
                           ipnum[i] = iphex[i]
00316
                        else:
00317
                           ipnum[i] = int(iphex[i].encode("hex"), 16)
00318
                   self.iplist = ipnum
00319
```

```
00320
              self.payloaderror = False
00321
00322
              # code pulled out of __init__ to facilitate reinitialization of the board
                  without needing to instantiate a new CameraAssembler object
00323
00324
              self.initialize()
00325
00326
00328
00329
          def initBoard(self):
00330
              return self.board.initBoard()
00331
00332
          def initPots(self):
00333
              return self.board.initPots()
00334
00335
          def latchPots(self):
00336
              return self.board.latchPots()
00337
00338
          def initSensor(self):
00339
              return self.board.initSensor()
00340
00341
          def configADCs(self):
00342
              return self.board.configADCs()
00343
00344
          def disarm(self):
00345
              return self.board.disarm()
00346
00347
          def startCapture(self, mode):
00348
              return self.board.startCapture(mode)
00349
00350
          def readSRAM(self):
00351
              return self.board.readSRAM()
00352
          def waitForSRAM(self, timeout=None):
00353
00354
              return self.board.waitForSRAM(timeout)
00355
00356
          def getTimer(self):
00357
              return self.board.getTimer()
00358
00359
          def resetTimer(self):
00360
              return self.board.resetTimer()
00361
00362
          def enableLED(self, status=1):
00363
              return self.board.enableLED(status)
00364
00365
          def setLED(self, LED=1, status=1):
00366
              return self.board.setLED(LED, status)
00367
00368
          def setPowerSave(self, status=1):
00369
              return self.board.setPowerSave(status)
00370
00371
          def setPPER(self, pollperiod=None):
00372
              return self.board.setPPER(pollperiod)
00373
00374
          def getTemp(self, scale=None):
00375
              return self.board.getTemp(scale)
00376
00377
          def getPressure(self, offset=None, sensitivity=None, units=None):
00378
              return self.board.getPressure(offset, sensitivity, units)
00379
00380
          def clearStatus(self):
00381
              return self.board.clearStatus()
00382
00383
          def checkStatus(self):
00384
              return self.board.checkStatus()
00385
          def checkStatus2(self):
00387
              return self.board.checkStatus2()
00388
00389
          def reportStatus(self):
00390
              return self.board.reportStatus()
00391
00392
          def reportEdgeDetects(self):
00393
              return self.board.reportEdgeDetects()
00394
00395
          def dumpStatus(self):
00396
              return self.board.dumpStatus()
00397
00398
          def checkSensorVoltStat(self):
00399
              return self.sensor.checkSensorVoltStat()
00400
00401
          def setTiming(self, side="AB", sequence=None, delay=0):
```

```
00402
              return self.sensor.setTiming(side, sequence, delay)
00403
00404
          def setArbTiming(self, side="AB", sequence=None):
00405
              return self.sensor.setArbTiming(side, sequence)
00406
00407
          def getTiming(self, side=None, actual=None):
00408
              return self.sensor.getTiming(side, actual)
00409
00410
          def setManualShutters(self, timing=None):
00411
              return self.sensor.setManualTiming(timing)
00412
00413
          def setManualTiming(self, timing=None):
00414
              return self.sensor.setManualTiming(timing)
00416
          def getManualTiming(self):
00417
              return self.sensor.getManualTiming()
00418
00419
          def getSensTemp(self, scale=None, offset=None, slope=None, dec=1):
00420
              return self.sensor.getSensTemp(scale, offset, slope, dec)
00421
00422
          def sensorSpecific(self):
00423
              return self.sensor.sensorSpecific()
00424
00425
          def selectOscillator(self, osc=None):
00426
              return self.sensor.selectOscillator(osc)
00427
          def setInterlacing(self, ifactor=None, side=None):
00428
00429
              return self.sensor.setInterlacing(ifactor, side)
00430
          def setHighFullWell(self, flag=True):
00431
00432
              return self.sensor.setHighFullWell(flag)
00433
00434
          def setZeroDeadTime(self, flag=True, side=None):
00435
              return self.sensor.setZeroDeadTime(flag, side)
00436
00437
          def setTriggerDelay(self, delay=0):
00438
              return self.sensor.setTriggerDelay(delay)
00439
00440
          def setPhiDelay(self, side=None, delay=0):
00441
              return self.sensor.setPhiDelay(side, delay)
00442
00443
          def setExtClk(self, dilation=None, frequency=None):
00444
              return self.sensor.setExtClk(dilation, frequency)
00445
00446
          def parseReadoff(self, frames, columns=1):
00447
               return self.sensor.parseReadoff(frames, columns)
00448
00449
          def sendCMD(self, pkt):
00450
              return self.comms.sendCMD(pkt)
00451
00452
          def arm(self, mode=None):
00453
              return self.comms.arm(mode)
00454
00455
          def readFrames(self, waitOnSRAM=None, timeout=0, fast=False, columns=1):
00456
              frames, _, _ = self.comms.readoff(waitOnSRAM, timeout, fast, columns)
00457
00458
00459
          def readoff(self, waitOnSRAM=None, timeout=0, fast=None, columns=1):
00460
              return self.comms.readoff(waitOnSRAM, timeout, fast, columns)
00461
00462
          def writeSerial(self, cmd, timeout=None):
00463
              return self.comms.writeSerial(cmd, timeout)
00464
00465
          def readSerial(self, size, timeout=None):
00466
              return self.comms.readSerial(size, timeout)
00468
          def closeDevice(self):
00469
              return self.comms.closeDevice()
00470
00471
          def saveTiffs(self, frames, path=None, filename="Frame", prefix=None, index=None):
00472
              return saveTiffs(self, frames, path, filename, prefix, index)
00473
00474
          def plotFrames(self, frames, index=None):
00475
              return plotFrames(self, frames, index)
00476
00477
          def getEnter(self, text):
00478
              return getEnter(text)
00479
00480
          def checkCRC(self, rval):
00481
              return checkCRC(rval)
00482
```

```
00483
          def str2bytes(self, astring):
00484
              return str2bytes(astring)
00485
00486
          def bytes2str(self, bytesequence):
00487
              return bytes2str(bytesequence)
00488
00489
          def str2nparray(self, valstring):
00490
              return str2nparray(valstring)
00491
00492
          def flattenlist(self, mylist):
00493
              return flattenlist (mylist)
00494
00495
          def partition(self, frames, columns):
00496
               return partition(self, frames, columns)
00497
00498
00499
          def initialize(self):
00500
00501
00502
               Initialize board registers and set pots
00503
00504
               # TODO: automate sensor and board selection from firmware info
00505
00507
00508
              # get sensor
00509
               # TODO: pull sensor, board, comm id out to separate methods
              if self.sensorname == "icarus":
00510
00511
                  import nsCamera.sensors.icarus as snsr
              elif self.sensorname == "icarus2":
00512
00513
                  import nsCamera.sensors.icarus2 as snsr
00514
              elif self.sensorname == "daedalus":
00515
                  import nsCamera.sensors.daedalus as snsr
00516
              else: # catch-all for added sensors to attempt object encapsulation
                  sensormodname = ".sensors." + self.sensorname
00517
00518
                      sensormod = importlib.import_module(sensormodname, "nsCamera")
00519
00520
                   except ImportError:
                       logging.critical(self.logcrit + "invalid sensor name")
00521
00522
                      sys.exit(1)
00523
                   snsr = getattr(sensormod, self.sensorname)
00524
              self.sensor = snsr(self)
00525
              # kill existing connections (for reinitialize)
if hasattr(self, "comms"):
00526
00527
00528
                   self.closeDevice()
00529
              # get communications interface
if self.commname == "rs422":
00530
00531
00532
                   import nsCamera.comms.RS422 as comms
00533
              elif self.commname == "gige":
00534
                   import nsCamera.comms.GigE as comms
00535
              else:
00536
                   commsmodname = ".comms." + self.commname
00537
                   try:
00538
                      commsmod = importlib.import_module(commsmodname, "nsCamera")
00539
                   except ImportError:
00540
                       logging.critical(self.logcrit + "invalid comms name")
00541
                       sys.exit(1)
                   comms = getattr(commsmod, self.commname)
00542
00543
              self.comms = comms(self)
00544
00545
              # get board
00546
              if self.boardname == "llnl_v1":
00547
                  import nsCamera.boards.LLNL_v1 as brd
00548
00549
                  self.board = brd.llnl_v1(self)
00550
              elif self.boardname == "llnl_v4":
00551
                  import nsCamera.boards.LLNL_v4 as brd
00552
00553
                   self.board = brd.llnl_v4(self)
00554
              else:
00555
                  boardmodname = ".board." + self.boardname
00556
00557
                       boardmod = importlib.import_module(boardmodname, "nsCamera")
00558
                   except ImportError:
00559
                       logging.critical(self.logcrit + "invalid board name")
00560
                       svs.exit(1)
                  boardobj = getattr(boardmod, self.boardname)
self.board = boardobj(self)
00561
00562
00563
00564
               # Now that board exists, initialize board-specific aliases for sensors
```

```
00565
              self.sensor.init_board_specific()
00566
00567
00568
00569
              # TODO: make cython the standard version
00570
              # ##############
00571
              # # For cython version
00572
00573
              # # get sensor
00574
              # if self.sensorname == "icarus":
00575
                    import nsCamera.sensors.icarus as snsr
00576
                    self.sensor = snsr.icarus(self)
00577
              # elif self.sensorname == "icarus2":
00578
                import nsCamera.sensors.icarus2 as snsr
00579
                    self.sensor = snsr.icarus2(self)
00580
              # elif self.sensorname == "daedalus":
00581
                   import nsCamera.sensors.daedalus as snsr
00582
                    self.sensor = snsr.daedalus(self)
00583
00584
              # # kill existing connections (for reinitialize)
              # if hasattr(self, "comms"):
00585
00586
                    self.closeDevice()
00587
00588
              # # get communications interface
00589
              # if self.commname == "rs422":
00590
                    import nsCamera.comms.RS422 as comms
00591
                    self.comms = comms.RS422(self)
              # elif self.commname == "gige":
00592
00593
                    import nsCamera.comms.GigE as comms
00594
                    self.comms = comms.GigE(self)
00595
00596
              # # get board
              # if self.boardname == "llnl_v1":
00597
00598
                    import nsCamera.boards.LLNL v1 as brd
00599
                    self.board = brd.llnl_v1(self)
              # elif self.boardname == "llnl_v4":
00600
00601
                    import nsCamera.boards.LLNL_v4 as brd
00602
                    self.board = brd.llnl_v4(self)
              # ################
00603
00604
              err, rval = self.getRegister("FPGA_NUM")
if err or rval == "":
00605
00606
                  err, rval = self.getRegister("FPGA_NUM")
if err or rval == "":
00607
00608
00609
                      logging.critical(
                           self.logcrit + "Initialization failed: unable to communicate with"
00610
00611
                           " board. "
00612
                      )
00613
                  sys.exit(1)
00614
00615
              self.initBoard()
00616
              self.initPots()
00617
              self.initSensor()
00618
              self.initPowerCheck()
00619
              self.getBoardInfo()
00620
              self.printBoardInfo()
00621
00622
          def reinitialize(self):
00623
00624
              Reinitialize board registers and pots, reinitialize sensor timing (if
              previously set)
00625
00626
00627
              logging.info(self.loginfo + "reinitializing")
00628
              self.initialize()
00629
00630
              for side in self.senstiming:
00631
                  self.setTiming(side, self.senstiming[side][0], self.senstiming[side][1])
00632
00633
              if self.sensmanual: # should be mutually exclusive with anything in senstiming
00634
                  self.setManualShutters(self.sensmanual)
00635
00636
          def reboot(self):
00637
00638
              Perform soft reboot on board and reinitialize
00639
00640
              self.board.softReboot()
00641
              self.reinitialize()
00642
00643
          def getBoardInfo(self):
00644
00645
              Get board info from FPGA NUM register. Returns error flag if register contents
```

```
00646
                are invalid and tuple (board version number, rad tolerance flag, sensor name)
00647
              Returns:
00648
              tuple (errorFlag, (board version, rad tolerance flag, sensor name))
00649
00650
00651
              invalidFPGANum = False
00652
              interfaces = []
00653
00654
              # TODO: move to new method (combine with parsing from initialize)
              if int(self.FPGANum[0], 16) & 8:
                  if self.FPGANum[1] == "1":
00656
                      boardtype = "LLNLv1"
00657
00658
                  elif self.FPGANum[1] == "4":
                      boardtype = "LLNLv4"
00659
00660
                  else:
                      boardtype = "LLNLv?"
00661
                      invalidFPGANum = True
00662
00663
              else:
00664
                  boardtype = "SNLrevC"
00665
                  logging.warning(
00666
                      self.logwarn + "FPGA self-identifies as SNLrevC, which is not"
00667
                      " supported by this software
00668
00669
                  invalidFPGANum = True
00670
              self.FPGAboardtype = boardtype
00671
00672
              if int(self.FPGANum[6], 16) & 1:
00673
                  rad = True
00674
              else:
00675
                 rad = False
00676
              self.FPGArad = rad
00677
              if self.FPGANum[7] == "1":
00678
                  sensor = "Icarus"
00679
              elif self.FPGANum[7] == "2":
00680
                 sensor = "Daedalus"
00681
00682
              else:
                  sensor = "Undefined"
00683
                  invalidFPGANum = True
00684
00685
              self.FPGAsensor = sensor
00686
00687
              if int(self.FPGANum[5], 16) & 1:
00688
                  interfaces.append("RS422")
00689
              if int(self.FPGANum[5], 16) & 2:
00690
                  interfaces.append("GigE")
00691
              self.FPGAinterfaces = interfaces
00692
00693
              if invalidFPGANum:
00694
                  if self.FPGANum == "80000001":
00695
                      invalidFPGANum = False
00696
00697
                      logging.warning(self.logwarn + "FPGA self-identification is invalid")
00698
              self.FPGAinvalid = invalidFPGANum
00699
00700
              return invalidFPGANum, (boardtype, rad, sensor)
00701
00702
          def getRegister(self, regname):
00703
00704
              Retrieves contents of named register as hexadecimal string without '0x'
00705
00706
00707
                  regname: name of register as given in ICD
00708
00709
              Returns:
              tuple: (error string, register contents as hexadecimal string without '0x')
00710
00711
00712
              # logging.debug(self.logdebug + "getRegister: regname = " + str(regname))
00713
              logging.debug(
00714
                  "{logdebug}getRegister: regname = {regname}".format(
00715
                      logdebug=self.logdebug, regname=regname
00716
00717
              )
00718
00719
              regname = regname.upper()
00720
              if regname not in self.board.registers:
00721
                  err = "{logerr}getRegister: Invalid register name: {regname}; returning"
                  " zeros".format(logerr=self.logerr, regname=regname)
00722
00723
                  logging.error(err)
                  return err, "00000000"
00724
00725
              sendpkt = Packet(cmd="1", addr=self.board.registers[regname])
00726
              err, rval = self.comms.sendCMD(sendpkt)
```

```
00727
              if err:
00728
                  logging.error(
00729
                      "{logerr}getRegister: {regname}; {err}".format(
00730
                          logerr=self.logerr, regname=regname, err=err
00731
00732
00733
00734
              retval = rval[8:16]
00735
              logging.debug(
00736
                  "{logdebug}getRegister: retval = {retval}".format(
00737
                      logdebug=self.logdebug, retval=retval
00738
00739
              )
00740
00741
              return err, retval
00742
00743
         def setRegister(self, regname, regval):
00744
00745
              Sets named register to given value
00746
00747
              Args:
00748
                  regname: name of register as given in ICD
                  regval: value to assign to register, as integer or hexadecimal string
00749
00750
                   with or without '0x'
00751
00752
              Returns:
              tuple: (error string, response string)
00753
00754
00755
              logging.debug(
00756
                   "{logdebug}setRegister: regname = {regname}; regval = {regval}".format(
00757
                      {\tt logdebug=self.logdebug, regname=regname, regval=regval}
00758
00759
              )
00760
00761
              regname = regname.upper()
00762
              if regname not in self.board.registers:
                  err = "{logerr}setRegister: Invalid register name: {regname}".format(
00763
00764
                      logerr=self.logerr, regname=regname
00765
00766
                  logging.error(err)
00767
                  return err, "00000000"
00768
              if isinstance(regval, int):
00769
                  regval = hex(regval)
00770
00771
                  if regval[0:2] == "0x":
00772
                      regval = regval[2:]
00773
              except TypeError:
00774
                  err = "{logerr}setRegister: invalid register value parameter".format(
00775
                     logerr=self.logerr
00776
00777
                  logging.error(err)
00778
                  return err, "00000000"
00779
              pkt = Packet(addr=self.board.registers[regname], data=regval)
00780
              err, rval = self.comms.sendCMD(pkt)
00781
              if err:
00782
                  logging.error(
00783
                      "{logerr}setRegister: {regname}: {err}".format(
00784
                          logerr=self.logerr, regname=regname, err=err
00785
00786
00787
              if len(rval) < 32:</pre>
00788
                  logging.debug(
00789
                      "{logdebug}SetRegister: rval = {rval}".format(
00790
                          logdebug=self.logdebug, rval=rval
00791
00792
                  )
00793
              else:
00794
                  logging.debug(
00795
                      "{logdebug}SetRegister: rval (truncated) = {rval}".format(
00796
                          logdebug=self.logdebug, rval=rval[0:32]
00797
00798
                  )
00799
              return err, rval
00800
00801
          def resolveSubreg(self, srname):
00802
00803
              Resolves subregister name or alias, returns object associated with subregister
00804
                and flag indicating writability
00805
00806
              Aras:
00807
                  srname: name or alias of subregister
```

```
00808
00809
              tuple(subregister name string, associated object, writable flag) _{\tt n\,u\,u}
00810
00811
00812
              logging.debug(
00813
                  "{logdebug}resolveSubreg: srname = {srname}".format(
00814
                      logdebug=self.logdebug,
00815
                      srname=srname,
00816
                  )
00817
00818
              writable = False
00819
              srname = srname.upper()
00820
              if srname in self.board.subreg_aliases:
00821
                  srname = self.board.subreg_aliases[srname].upper()
00822
              if srname in self.board.subreglist:
00823
                  srobj = getattr(self.board, srname)
00824
                  writable = getattr(self.board, srname).writable
00825
              else:
00826
                  # No-object error is handled by calling function
00827
                  srobj = None
00828
              logging.debug(
00829
                   "{logdebug}resolveSubreg: srobj = {srobj}, writable={writable}".format(
00830
                      logdebug=self.logdebug, srobj=srobj, writable=writable
00831
00832
              )
00833
00834
              return srname, srobj, writable
00835
00836
          def getSubregister(self, subregname):
00837
00838
              Returns substring of register identified in board attribute 'subregname'
00839
00840
00841
                  subregname: listed in board.subreg_aliases or defined in board.subregisters
00842
00843
              Returns:
              tuple: (error string, contents of subregister as binary string without '0b')
00844
00845
00846
              logging.debug(
                   "{logdebug}getSubegister: subregname = {subregname}".format(
00847
00848
                      logdebug=self.logdebug,
00849
                      subregname=subregname,
00850
00851
              )
00852
00853
              subregname, subregobj, _ = self.resolveSubreg(subregname)
00854
              if not subregobj:
00855
                  err = "{logerr}getSubregister: invalid lookup: {subregname}; returning"
00856
                  " string of zeroes".format(logerr=self.logerr, subregname=subregname)
00857
00858
                  logging.error(err)
                  return err, "".zfill(8)
00859
00860
              err, resp = self.getRegister(subregobj.register)
00861
00862
                  logging.error(
00863
                       "{logerr}getSubregister: unable to retrieve register setting: \
00864
                      {subregname}; returning '0' string".format(
00865
                          logerr=self.logerr, subregname=subregname
00866
00867
                  )
00868
00869
                  return err, "".zfill(8)
00870
              hex_str = "0x" + resp # this should be a hexadecimalstring
              b_reg_value = "{0:0=32b}".format(int(hex_str, 16))  # convert to binary string
00871
00872
              # list indexing is reversed from bit string; the last bit of the string is at
00873
                  index 0 in the list (thus bit 0 is at index 0)
00874
              startindex = 31 - subregobj.start_bit
00875
              retval = b_req_value[startindex : startindex + subregobj.width]
00876
              logging.debug(
00877
                   "{logdebug}getSubregister: retval = {retval}".format(
00878
                      logdebug=self.logdebug, retval=retval
00879
00880
00881
              return "", retval
00882
00883
          def setSubregister(self, subregname, valstring):
00884
00885
              Sets substring of register identified in board attribute 'subregname' to
               valstring if subregister is writable
00886
00887
00888
              Args:
```

```
00889
                   subregname: listed in board.subreg_aliases or defined in board.subregisters
00890
                   valstring: integer or binary string with or without '0b'
00891
00892
              tuple: (error, packet response string) from setRegister
00893
00894
00895
              logging.debug(
00896
                   "{logdebug}setSubegister: subregname = {subregname}; valstring ="
00897
                   " {valstring}".format(
00898
                       logdebug=self.logdebug, subregname=subregname, valstring=valstring
00899
00900
00901
00902
              subregname, subregobj, writable = self.resolveSubreg(subregname)
00903
              if not subregobj:
00904
                  err = "{logerr}getSubregister: invalid lookup: {subregname}".format(
00905
                       logerr=self.logerr, subregname=subregname
00906
00907
00908
                  logging.error(err)
00909
                  return err, "0'
00910
              if not writable:
00911
                  err = "{logerr}getSubregister: not a writable subregister: {subregname}"
00912
                   "".format(logerr=self.logerr, subregname=subregname)
00913
                  logging.error(err)
00914
                   return err, "0"
00915
              if isinstance(valstring, int):
00916
                  valstring = bin(valstring)[2:]
00917
                  if valstring[0:2] == "0b":
00918
00919
                      valstring = valstring[2:]
00920
              except TypeError:
00921
                  err = "{logerr}getSubregister: invalid subregister value parameter".format(
00922
                      logerr=self.logerr
00923
00924
00925
                  logging.error(err)
return err, "0"
00926
              if len(str(valstring)) > subregobj.width:
00927
00928
                  err = "{logerr}getSubregister: ialue string is too long".format(
00929
                      logerr=self.logerr
00930
                  )
00931
00932
                  logging.error(err)
              return err, "0"
# read current value of register data
00933
00934
00935
              err, resp = self.getRegister(subregobj.register)
00936
              if err:
00937
                  logging.error(
                       "{logerr}getSubregister: unable to retrieve register setting; setting"
" of {subregname} likely failed ".format(
00938
00939
00940
                           logerr=self.logerr, subregname=subregname
00941
00942
                  )
00943
00944
                  return err, "0"
00945
              hex_str = "0x" + resp
00946
              b_reg_value = "{0:0=32b}".format(int(hex_str, 16)) # convert to binary
00947
               # list indexing is reversed from bit string; the last bit of the string is at
00948
                  index 0 in the list (thus bit 0 is at index 0)
              startindex = 31 - subregobj.start_bit
00949
00950
              valstringpadded = str(valstring).zfill(subregobj.width)
00951
              fullreg = list(b_reg_value)
00952
              fullreg[startindex : startindex + subregobj.width] = valstringpadded
               # convert binary string back to hexadecimal string for writing
00953
                               "".join(fullreg)
00954
              new_reg_value =
              h_reg_value = "{num:{fill}{width}x}".format(
00955
00956
                  num=int(new_reg_value, 2), fill="0", width=8
00957
00958
              err, retval = self.setRegister(subregobj.register, h_reg_value)
00959
              # logging.debug(self.logdebug + "retval = " + str(retval))
00960
              if len(retval) < 32:</pre>
                   logging.debug(
00961
00962
                       "{logdebug}setSubregister: retval = {retval}".format(
00963
                           logdebug=self.logdebug, retval=retval
00964
00965
                  )
00966
              else:
00967
                  logging.debug(
00968
                       "{logdebug}setSubregister: retval (truncated) = {retval}".format(
00969
                           logdebug=self.logdebug, retval=retval[0:32]
```

```
00970
                      )
00971
00972
00973
              return err, retval
00974
00975
          def submitMessages(self, messages, errorstring="Error"):
00976
00977
              Serially set multiple register / subregister values
00978
00979
              Args:
00980
                  messages: list of tuples (register name, integer or hexadecimal string with
                   or without '0x') and/or (subregister name, integer or binary string with
00981
00982
                    or without '0b')
00983
                  errorstring: error message to print in case of failure
00984
00985
              tuple (accumulated error string, response string of final message) _{\tt mum}
00986
00987
00988
              logging.debug(
00989
                  "{logdebug}submitMessages: messages = {messages}; errorstring ="
00990
                  " {errorstring}".format(
00991
                      logdebug=self.logdebug, messages=messages, errorstring=errorstring
00992
00993
              )
00994
00995
              errs = ""
              err = ""
00996
              rval = ""
00997
00998
              for m in messages:
                  if m[0].upper() in self.board.registers:
00999
01000
                      err, rval = self.setRegister(m[0].upper(), m[1])
01001
                  elif m[0].upper() in self.board.subreglist:
                      err, rval = self.setSubregister(m[0].upper(), m[1])
01002
01003
                  else:
                      err = "{logerr}submitMessages: Invalid register/subregister:"
01004
                       " {errorstring}:{m0}; ".format(
01005
01006
                          logerr=self.logerr, errorstring=errorstring, m0=m[0]
01007
01008
01009
                     logging.error(err)
01010
                  errs = errs + err
01011
              return err, rval
01012
01013
          def getPot(self, potname, errflag=False):
01014
01015
              Retrieves value of pot or ADC monitor subregister, scaled to [0,1).
01016
01017
              Args:
01018
                 potname: name of pot or monitor, e.g., VRST or MON_CH2 found in
01019
                    board.subreg_aliases or defined in board.subregisters
01020
                  errflag: if True, return tuple with error string
01021
01022
              Returns:
01023
                  if errflag:
01024
                      tuple: (error string, float value of subregister, scaled to [0,1) )
01025
01026
                     float value of subregister, scaled to [0,1)
01027
01028
              logging.debug(
01029
                  "{logdebug}getPot: potname = {potname}; errflag = {errflag}".format(
01030
                      logdebug=self.logdebug, potname=potname, errflag=errflag
01031
01032
01033
01034
              potname, potobj, _ = self.resolveSubreg(potname)
01035
              if not potobj:
                  err = "{logerr}getPot: invalid lookup: {potname}; returning 0".format(
01036
                      logerr=self.logerr, potname=potname
01037
01038
01039
01040
                  logging.error(err)
01041
                  if errflag:
01042
                     return err. 0
01043
                  return 0
01044
              err, b_pot_value = self.getSubregister(potname)
01045
              if err:
                  err = "{logerr}getPot: unable to read subregister: {potname}".format(
01046
01047
                      {\tt logerr=self.logerr,\ potname=potname}
01048
01049
01050
              # convert binary string back to decimal
```

```
01051
               f_reg_value = 1.0 * int(b_pot_value, 2)
               value = (f_reg_value - potobj.min) / (potobj.max - potobj.min)
# logging.debug(self.logdebug + "getpot: value = " + str(value))
01052
01053
01054
01055
               logging.debug(
01056
                   "{logdebug}getpot: value = {value}".format(
01057
                       logdebug=self.logdebug, value=value
01058
01059
01060
01061
               if errflag:
01062
                  return err, value
01063
               return value
01064
01065
          def setPot(self, potname, value=1.0, errflag=False):
01066
01067
               Sets value of pot to value, normalized so that '1.0' corresponds with the
01068
                fixed point maximum value of pot.
01069
01070
               Args:
01071
                  potname: common name of pot, e.g., VRST found in board.subreg_aliases or
01072
                    defined in board.subregisters
01073
                   value: float between 0 and 1
01074
                   errflag: if True, return tuple with error string
01075
01076
               Returns:
01077
                   if errflag:
                       tuple: (error string, response packet as string)
01078
01079
                   else:
01080
                       response packet as string
01081
01082
               logging.debug(
01083
                   "{logdebug}setPot: potname = {potname}; value={value} errflag = {errflag}"
                   "".format(
01084
                       logdebug=self.logdebug, potname=potname, value=value, errflag=errflag
01085
01086
01087
              )
01088
01089
               if value < 0:</pre>
                  value = 0.0
01090
              if value > 1:
01091
01092
                   value = 1.0
01093
01094
               potname, potobj, writable = self.resolveSubreg(potname)
01095
               if not potobj:
01096
                   err = "{logerr}setPot: invalid lookup: {potname}; returning '0'".format(
01097
                       logerr=self.logerr, potname=potname
01098
01099
01100
                   logging.error(err)
01101
                   if errflag:
01102
                       return err, 0
01103
                   return 0
01104
               if not writable:
01105
                   err = "{logerr}setPot: not a writable subregister: {potname}; returning '0'"
01106
                   "".format(logerr=self.logerr, potname=potname)
01107
                   logging.error(err)
01108
                   if errflag:
01109
                       return err, "0"
01110
                   return "0"
               setpoint = int(round(value * potobj.max_value))
01111
01112
               setpointpadded = "{num:{fill}{width}b}".format(
01113
                  num=setpoint, fill="0", width=potobj.width
01114
01115
               logging.debug(
                   "{logdebug}setpot: setpointpadded = {setpointpadded}".format(
01116
                       logdebug=self.logdebug, setpointpadded=setpointpadded
01117
01118
01119
              )
01120
01121
               err, rval = self.setSubregister(potname, setpointpadded)
01122
               if err:
01123
                   logging.error(
                       err="{logerr}setPot: unable to confirm setting of subregister:"
01124
                         {potname} ".format(logerr=self.logerr, potname=potname)
01125
01126
              ident = potname[3:]
01127
               if ident[0].isdigit(): # numbered pot scheme
  potnumlatch = int(ident) * 2 + 1
01128
01129
                   potnumlatchstring = "{num:{fill}{width}x}".format(
    num=potnumlatch, fill="0", width=8
01130
01131
```

```
01132
01133
                   err1, resp = self.setRegister("POT_CTL", potnumlatchstring)
01134
              else: # alphabetical DAC scheme
                 ident = ident.upper() # expects single character, e.g. 'A' from 'DACA'
identnum = ord(ident) - ord("A") # DACA -> 0
01135
01136
01137
                   potnumlatch = int(identnum) * 2 + 1
01138
                  potnumlatchstring = "{num:{fill}{width}x}".format(
01139
                      num=potnumlatch, fill="0", width=8
01140
01141
                  err1, resp = self.setRegister("DAC_CTL", potnumlatchstring)
01142
              if err1:
                  # logging.error(self.logerr + "setPot: unable to latch register")
01143
01144
01145
                  logging.error(
01146
                       err="{logerr}setPot: unable to latch register".format(
01147
                          logerr=self.logerr
01148
01149
01150
01151
              if errflag:
01152
                  return err + err1, rval
01153
              return rval
01154
01155
          def getPotV(self, potname, errflag=False):
01156
01157
              Reads voltage setting (not actual voltage) of specified pot
01158
01159
01160
                  potname: name of pot or monitor, e.g., VRST or MON_CH2 found in
01161
                    board.subreq_aliases or defined in board.subregisters
01162
                  errflag: if True, return tuple with error string
01163
01164
              Returns:
                  if errflag:
01165
01166
                      tuple: (error string, float value of pot voltage)
01167
                   else:
01168
                       float value of pot voltage
01169
01170
              logging.debug(
01171
                  self.logdebug
01172
                   + "getPotV: potname = "
01173
                  + str(potname)
                  + "; errflag = "
01174
01175
                   + str(errflag)
01176
01177
              potname, potobj, _ = self.resolveSubreg(potname)
01178
              if not potobj:
01179
                  err = (
01180
                      self.logerr + "getPotV: invalid lookup: " + potname + " , returning 0 "
01181
01182
                   logging.error(err)
01183
                   if errflag:
01184
                      return err, 0
01185
                  return 0
              err, val = self.getPot(potname, errflag=True)
logging.debug(self.logdebug + "getPotV: val = " + str(val))
01186
01187
01188
01189
                   logging.error(self.logerr + "getPotV: unable to read pot " + potname)
01190
              minV = potobj.minV
01191
              maxV = potobj.maxV
01192
              if errflag:
01193
                  return err, val * (maxV - minV)
01194
              return val * (maxV - minV)
01195
01196
          # TODO: optimize tuning speed for DACs
01197
          def setPotV(
01198
              self.
01199
              potname,
01200
              voltage,
01201
              tune=False,
01202
              accuracy=0.01,
01203
              iterations=20,
01204
              approach=0.75,
01205
              errflag=False,
01206
         ):
01207
01208
              Sets pot to specified voltage. If tune=True, uses monitor to adjust pot to
01209
                correct voltage. Tuning will attempt to tune to closest LSB on pot; if
                 'accuracy' > LSB resolution, will only complain if tuning is unable to get
01210
                the voltage within 'accuracy'
01211
01212
```

```
01213
              Args:
01214
                  potname: common name of pot, e.g., VRST found in board.subreg_aliases or
01215
                    defined in board.subregisters
01216
                  voltage: voltage bound by pot max and min (set in board object)
01217
                  tune: if True, iterate with monitor to correct voltage
01218
                  accuracy: acceptable error in volts (if None, attempts to find the closest
01219
                    possible pot setting and warns if last iteration does not reduce error
01220
                    below the resolution of the pot)
01221
                  iterations: number of iteration attempts
01222
                  approach: approximation parameter (>1 may cause overshoot)
01223
                  errflag: if True, return tuple with error string
01224
01225
              Returns:
01226
                  if errflag:
01227
                      tuple: (error string, response string)
01228
01229
                      response string
01230
01231
              logging.debug(
01232
                  self.logdebug
                  + "setPotV: potname = "
01233
01234
                  + str(potname)
01235
                  + "; voltage =
01236
                  + str(voltage)
01237
                  + "; tune =
                  + str(tune)
01238
                  + "; accuracy =
01239
01240
                  + str(accuracy)
                  + ": iterations = "
01241
01242
                  + str(iterations)
01243
                  + "; approach =
01244
                  + str(approach)
                  + "; errflag = "
01245
01246
                  + str(errflag)
01247
             )
              potname, potobj, writable = self.resolveSubreg(potname)
01248
01249
              if not potobj:
01250
                  err = (
                     self.logerr
01251
01252
                      + "setPotV: invalid lookup: "
01253
                      + potname
                      + " , returning zero "
01254
01255
01256
                  logging.error(err)
01257
                  if errflag:
                      return err, 0
01258
01259
                  return 0
01260
              if not writable:
01261
                  err = (
01262
                      self.logerr
01263
                      + "setPotV: not a writable subregister: "
01264
                      + potname
01265
                      + "; returning zero"
01266
01267
                  logging.error(err)
01268
                  if errflag:
01269
                      return err, 0
01270
                  return 0
01271
              if voltage < potobj.minV:</pre>
01272
                  voltage = potobj.minV
01273
              if voltage > potobj.maxV:
01274
                  voltage = potobj.maxV
01275
              setting = (voltage - potobj.minV) / (potobj.maxV - potobj.minV)
              logging.debug(self.logdebug + "setPotV: setting = " + str(setting))
01276
01277
              err, rval = self.setPot(potname, setting, errflag=True)
01278
              time.sleep(0.1)
01279
              # TODO: refactor tuning to separate method
01280
              if tune:
                  logging.debug(self.logdebug + "setPotV: beginning tuning")
01281
01282
                  if potname not in self.board.monitor_controls.values():
01283
                      err = (
01284
                          self.logerr
01285
                          + "setPotV: pot '"
01286
                          + potname
                          + "' does not have a corresponding monitor"
01287
01288
01289
                      logging.error(err)
01290
                      if errflag:
01291
                         return err, rval
                      return rval
01292
                  self.setPot(potname, 0.65)
01293
```

```
01294
                   time.sleep(0.2)
01295
                   err1, mon65 = self.getMonV(potname, errflag=True)
01296
                   self.setPot(potname, 0.35)
                   time.sleep(0.2)
01297
01298
                   err2, mon35 = self.getMonV(potname, errflag=True)
01299
                   # theoretical voltage range assuming linearity
01300
                   potrange = (mon65 - mon35) / 0.3
01301
                   stepsize = potrange / (potobj.max_value + 1)
01302
                   err += err1 + err2
01303
                   if err or potrange < 1:</pre>
                       err += self.logerr + "setPotV: unable to tune pot " + potname
01304
                       if potrange < 1: # potrange should be on the order of 3.3 or 5 volts
01305
01306
                           err += "; monitor shows insufficient change with pot variation;
01307
                           "retrying setPotV with tune=False"
01308
                       logging.warning(err)
01309
                       err, rval = self.setPotV(
01310
                           potname=potname, voltage=voltage, tune=False, errflag=True
01311
01312
                       if errflag:
01313
                          return err, rval
01314
                       return rval
01315
                   potzero = 0.35 - (mon35 / potrange)
                   potone = 1.65 - (mon65 / potrange)
01316
01317
                   if potzero < 0:</pre>
01318
                       potzero = 0
                   if potone > 1:
01319
01320
                       potone = 1
01321
                  if accuracy > stepsize:
    mindiff = accuracy
01322
01323
01324
                   else:
                      mindiff = stepsize
01325
01326
                   setting = potzero + (voltage / potone)
                   self.setPot(potname, setting)
01327
                   lastdiff = 0
01328
01329
                   smalladjust = 0
                   err3 = ""
01330
                   for _ in range(iterations):
01331
01332
                       err3i, measured = self.getMonV(potname, errflag=True)
01333
                       if err3i:
01334
                           err3 = err3 + err3i + " "
01335
                       diff = voltage - measured
01336
                       if abs(diff - lastdiff) < stepsize / 2:</pre>
01337
                           if (
01338
                                smalladjust > 12
01339
                               # magic number for now; if it doesn't converge after several
01340
                                    tries, it never will, usually because the setting is pinned
01341
                                    to 0 or 1 and adjust can't change it
01342
                                logging.warning(
01343
                                    self.logwarn
01344
                                    + "setPotV: Tuning converged too slowly: pot "
01345
                                    + potname
01346
                                    + " set to "
01347
                                    + str(voltage)
01348
                                    + "V, monitor returns "
01349
                                    + str(measured)
01350
                                    + "V; if this value is incorrect, consider trying "
01351
                                    + "tune=False"
01352
01353
                                logging.debug(self.logdebug + "setPotV: tuning complete")
                               if errflag:
    return "", rval
01354
01355
01356
                                return rval
01357
                           smalladjust += 1
                       if not int (2 * diff / stepsize):
01358
                           # TODO: is this check redundant with the first one?
01359
                           logging.debug(self.logdebug + "setPotV: tuning complete")
01360
01361
                           if errflag:
                               return "", rval
01362
01363
                           return rval
01364
                       adjust = approach * (diff / potrange)
01365
                       setting += adjust
01366
                       if setting > 1:
01367
                           setting = 1
01368
                       elif setting < 0:</pre>
                          setting = 0
01369
                       err1, rval = self.setPot(potname, setting, True)
01370
01371
                       lastdiff = diff
01372
                       time.sleep(0.2)
                  err4, measured = self.getMonV(potname, errflag=True)
diff = voltage - measured
01373
01374
```

```
01375
                   # code will try to get to within one stepsize, but will only complain if it
01376
                     doesn't get within mindiff
01377
                   if int(diff / mindiff):
                       logging.warning(
01378
01379
                           self.logwarn
01380
                           + "setPotV: pot "
01381
                           + potname
01382
                            + " set to "
01383
                           + str(voltage)
01384
                            + "V, monitor returns "
01385
                           + str(measured)
01386
                            + "17"
01387
                   err += err1 + err2 + err3 + err4
01388
01389
               if err:
01390
                  logging.error(self.logerr + "setPotV: errors occurred: " + err)
01391
               if errflag:
01392
                  return err, rval
01393
              logging.debug(self.logdebug + "setPotV: tuning complete")
01394
              return rval
01395
01396
          def getMonV(self, monname, errflag=False):
01397
01398
              Reads voltage from monitor named or that associated with the pot named 'monname'
01399
01400
              Args:
                  monname: name of pot or monitor, e.g., VRST or MON_CH2 found in
01401
01402
                    board.subreg_aliases or defined in board.subregisters
01403
                   errflag: if True, return tuple with error string
01404
01405
              Returns:
01406
                  if errflag:
01407
                      tuple: (error string, float value of voltage measured by monitor)
01408
                   else:
                       float value of voltage measured by monitor
01409
01410
01411
              logging.debug(
01412
                   self.logdebug
                   + "getMonV: monname = "
01413
01414
                  + str(monname)
                  + "; errflag = "
01415
01416
                  + str(errflag)
01417
              monname = monname.upper()
01418
01419
              if monname in self.board.subreg_aliases:
01420
                   monname = self.board.subreg_aliases[monname].upper()
01421
              # else:
01422
              for key, value in self.board.monitor_controls.items():
                  if value == monname:
01423
01424
                      monname = key
01425
              if monname not in self.board.monitor_controls:
01426
                  if monname in self.board.subreglist:
01427
                      pass # no change necessary
01428
                   else:
01429
01430
                           self.logerr + "getMonV: invalid lookup " + monname + ", returning 0"
01431
01432
                       logging.error(err)
01433
                       if errflag:
01434
                          return err, 0
01435
                       return 0
01436
              err, monval = self.getPot(monname, errflag=True)
              logging.debug(self.logdebug + "getMonV: monval = " + str(monval))
01437
01438
              if err:
01439
                  logging.error(
01440
                       self.logerr + "getMonV: unable to read monitor value for " + monname
01441
01442
              if self.board.ADC5_bipolar:
01443
                  if monval >= 0.5:
01444
                       monval -= 1 # handle negative measurements (two's complement)
01445
                   if errflag:
                  return err, 2 * self.board.ADC5_mult * monval * self.board.VREF
return 2 * self.board.ADC5_mult * monval * self.board.VREF
01446
01447
01448
              else:
01449
                   return err, self.board.ADC5_mult * monval * self.board.VREF
return self.board.ADC5_mult * monval * self.board.VREF
01450
01451
01452
01453
          def readImgs(self, waitOnSRAM=True, mode="Hardware"):
01454
              Combines arm() and readoff() functions
01455
```

```
01456
01457
              Returns:
01458
                  tuple (list of numpy arrays, length of downloaded payload, payload error
01459
                    flag) returned by readoff
01460
01461
              logging.info(self.loginfo + "readImgs")
01462
              self.arm(mode)
01463
              return self.readoff(waitOnSRAM)
01464
01465
          def saveFrames(self, frames, path=None, filename="frames", prefix=None):
01466
              Save list of numpy arrays to disk. If passed an unprocessed text string, saves
01467
01468
                it directly to disk for postprocessing. Use 'prefix=""' for no prefix
01469
01470
              Args:
01471
                  frames: numpy array or list of numpy arrays OR text string
01472
                  path: save path, defaults to './output'
                   filename: defaults to 'frames.bin'
01473
01474
                  prefix: prepended to filename, defaults to time/date (e.g. '160830-124704_')
01475
01476
              Returns:
              Error string
01477
01478
01479
              logging.debug(
01480
                  self.logdebug
                   + "saveFrames: path = "
01481
01482
                  + str(path)
                  + "; filename = "
01483
                  + str(filename)
01484
01485
                  + "; prefix =
01486
                  + str(prefix)
01487
              logging.info(self.loginfo + "saveFrames")
01488
01489
              err = "
              if path is None:
01490
                  path = os.path.join(os.getcwd(), "output")
01491
              if prefix is None:
01492
                  prefix = datetime.now().strftime("%y%m%d-%H%M%S%f")[:-5] + "_"
01493
              if not os.path.exists(path):
01494
01495
                  os.makedirs(path)
01496
01497
              # TODO catch save file exceptions
01498
              if isinstance(frames[0], str):
01499
                  logging.debug(self.logdebug + "saveFrames: saving text frames")
                  filename = filename + ".txt" savefile = open(os.path.join(path, prefix + filename), "w+")
01500
01501
01502
                  savefile.write(frames)
01503
              else:
01504
                  logging.debug(self.logdebug + "saveFrames: saving numerical frames")
01505
                   filename = filename + ".bin"
01506
                   stacked = np.stack(frames)
01507
01508
                       stacked = stacked.reshape(
01509
                           (
01510
                               self.sensor.nframes,
01511
                               self.sensor.height // (self.sensor.interlacing + 1),
01512
                               self.sensor.width,
01513
01514
01515
                  except Exception as e:
01516
                      err = self.logerr + "saveFrames: unable to save frames: " + str(e)
01517
                      logging.error(err)
01518
01519
                  stacked.tofile(os.path.join(path, prefix + filename))
01520
              return err
01521
01522
          def saveNumpys(
01523
              self,
01524
              frames,
01525
              path=None,
              filename="Frame",
01526
01527
              prefix=None,
01528
              index=None,
01529
         ):
01530
01531
              Save numpy array or list of numpy arrays to disk as individual numpy data files,
01532
                with frame number appended to filename.
01533
01534
              Args:
01535
                  frames: numpy array or list of numpy arrays or single numpy array
                  path: save path, defaults to './output'
01536
```

```
01537
                  filename: defaults to 'Frame' followed by frame number
                  prefix: prepended to 'filename', defaults to time/date
01538
01539
                    (e.g. '160830-124704_')
01540
                  index: number to start frame numbering
01541
01542
              Returns:
              Error string
01543
01544
01545
              logging.info(self.loginfo + "saveNumpys")
01546
              logging.debug(
01547
                  self.logdebug
01548
                  + "saveNumpys: path = "
01549
                  + str(path)
01550
                  + "; filename = "
01551
                  + str(filename)
01552
                  + "; prefix =
01553
                  + str(prefix)
01554
                  + "; index =
01555
                  + str(index)
01556
             )
              err = ""
01557
              if path is None:
01558
01559
                  path = os.path.join(os.getcwd(), "output")
01560
              if prefix is None:
01561
                  prefix = datetime.now().strftime("%y%m%d-%H%M%S%f")[:-5] + "_"
              if not os.path.exists(path):
01562
01563
                  os.makedirs(path)
01564
              if index is None:
01565
                  firstnum = self.sensor.firstframe
01566
              else:
01567
                  firstnum = index
01568
              if not isinstance(frames, list):
01569
                  frames = [frames]
01570
              # if this is a text string from fast readoff, do the numpy conversion now
01571
01572
              if isinstance(frames[0], str):
01573
                  frames = generateFrames(frames)
01574
01575
              framestemp = np.copy(frames)
01576
01577
              for idx, frame in enumerate(framestemp):
01578
                  if idx < len(framestemp) / 2:</pre>
01579
                      interlacing = self.sensor.interlacing[0]
01580
                  else:
01581
                      interlacing = self.sensor.interlacing[1]
01582
01583
                      if self.padToFull:
01584
                          frame = np.reshape(
01585
                              frame, (self.sensor.maxheight // (interlacing + 1), -1)
01586
                          )
01587
                       else:
01588
                           frame = np.reshape(
01589
                              frame,
01590
01591
                                   (self.sensor.lastrow - self.sensor.firstrow + 1)
01592
                                  // (interlacing + 1),
01593
                                   -1,
01594
                              ),
01595
                          )
01596
                      namenum = filename + "_%d" % firstnum
01597
                      nppath = os.path.join(path, prefix + namenum + ".npy")
01598
                      np.save(nppath, frame)
01599
                      firstnum += 1
01600
                  except SystemExit:
01601
01602
                  except KeyboardInterrupt:
01603
                     raise
01604
                  except Exception:
01605
                     err = self.logerr + "saveNumpys: unable to save arrays"
01606
                      logging.error(err)
01607
01608
              return err
01609
          def dumpNumpy(
01610
01611
              self,
01612
              datastream.
              path=None,
01613
              filename="Dump",
01614
01615
              prefix=None,
01616
         ):
01617
```

```
01618
              Datastream is converted directly to numpy array and saved to disk. No attempt to
01619
                parse headers or separate into individual frames is made. The packet header is
01620
                removed before saving
01621
01622
              Args:
01623
                  datastream: string to be saved
01624
                  path: save path, defaults to './output'
01625
                  filename: defaults to 'Dump'
                  prefix: prepended to 'filename', defaults to time/date
01626
                    (e.g. '160830-124704_')
01627
01628
01629
              Returns:
              Error string
01630
01631
01632
              logging.info(self.loginfo + "dumpNumpy")
01633
              logging.debug(
01634
                  self.logdebug
01635
                  + "dumpNumpy: path = "
01636
                  + str(path)
01637
                  + "; filename = "
                  + str(filename)
01638
01639
                  + "; prefix =
01640
                  + str(prefix)
01641
              )
              err = ""
01642
              if path is None:
01643
01644
                  path = os.path.join(os.getcwd(), "output")
01645
              if prefix is None:
                  prefix = time.strftime("%y%m%d-%H%M%S_", time.localtime())
01646
01647
              if not os.path.exists(path):
01648
                  os.makedirs(path)
              npdata = str2nparray(datastream[36:])
01649
01650
              try:
                  nppath = os.path.join(path, prefix + filename + ".npy")
01651
01652
                  np.save(nppath, npdata)
01653
              except SystemExit:
01654
01655
              except KeyboardInterrupt:
01656
01657
              except Exception:
01658
                  err = self.logerr + "dumpNumpy: unable to save data stream"
01659
                  logging.error(err)
01660
              return err
01661
01662
          def checkRegSet(self, regname, teststring):
01663
01664
              Quick check to confirm that data read from register matches data write
01665
01666
01667
                  regname: register to test
01668
                  teststring: value to assign to register, as integer or hexadecimal string
01669
                    with or without '0x'
01670
01671
              boolean, True if read and write values match _{\rm m\,m\,m}
01672
01673
01674
              self.setRegister(regname, teststring)
01675
              # tell board to send data; wait to clear before interrogating register contents
01676
              if regname == "SRAM_CTL":
01677
                  time.sleep(2)
01678
                  if self.commname == "rs422":
01679
                       logging.info(
01680
                          self.loginfo + "skipping 'SRAM_CTL' register check for RS422"
01681
01682
                      return True
01683
              else:
01684
                 time.sleep(0.1)
01685
              temp = self.getRegister(regname)
              resp = temp[1].upper()
01686
01687
              if resp != teststring.upper():
01688
                  logging.error(
01689
                      self.logerr
                       + "checkRegSet failure: "
01690
01691
                      + regname
                       + " ; set: "
01692
01693
                      + teststring
                       + " ; read:
01694
01695
                      + resp
01696
01697
                  return False
01698
              return True
```

```
01699
01700
          def initPowerCheck(self):
01701
              Reset software and board timers for monitoring power status
01702
01703
01704
               self.inittime = time.time()
01705
               logging.info(self.loginfo + "resetting timer for power check function")
01706
               self.resetTimer()
01707
01708
          def powerCheck(self, delta=10):
01709
01710
               Check to see if board power has persisted since powerCheck was last initialized.
01711
                 Compares time elapsed since initialization against board's timer. If the
01712
                 difference is greater than 'delta,' flag as False (power has likely failed)
01713
01714
01715
                  delta: difference in seconds permitted between software and board timers
01716
01717
              Returns:
01718
                  boolean, 'True' means timer difference is less than 'delta' parameter;
                            'False' indicates power failure
01719
01720
01721
              elapsed = time.time() - self.inittime
01722
              logging.debug(self.logdebug + "powerCheck: elapsed time = " + str(elapsed))
01723
              difference = abs(elapsed - self.getTimer())
              if difference > delta:
01724
01725
                  logging.warning(
                       self.logwarn + "powerCheck function has failed; may indicate current "
01726
01727
                       "or recent power failure "
01728
01729
              return difference < delta
01730
01731
          def printBoardInfo(self):
01732
               # TODO: add override option if logging level is above info
01733
              logging.info(
01734
                  self.loginfo
01735
                   + "Python version: "
01736
                  + str(self.python)
01737
01738
                  + str(self.pyth1)
01739
01740
                  + str(self.pyth2)
01741
              logging.info(self.loginfo + "nsCamera software version: " + self.version)
01742
              logging.info(self.loginfo + "FPGA firmware version: " + self.FPGAVersion) logging.info(self.loginfo + "FPGA implementation: " + self.FPGANum)
01743
01744
01745
              if self.FPGAinvalid:
                   logging.info(self.loginfo + "FPGA information unavailable")
01746
01747
01748
                   logging.info(self.loginfo + "Board type: " + self.FPGAboardtype)
                  logging.info(self.loginfo + "Rad-Tolerant: " + str!.FrGADatatype)
logging.info(self.loginfo + "Rad-Tolerant: " + str!.FrGADatatype)
01749
01750
                   logging.info(self.loginfo + "Sensor label: " + self.sensor.loglabel)
01751
01752
                   logging.info(
01753
                       self.loginfo + "Available interfaces: " + ", ".join(self.FPGAinterfaces)
01754
01755
               if self.commname == "gige":
01756
                   ci = self.comms.CardInfoP.contents
                   ip = ".".join(str(e) for e in [b for b in ci.IPAddr])
01757
01758
                   logging.info(
01759
                       self.loginfo + "GigE connected to " + ip + ":" + str(self.port)
01760
01761
               elif self.commname == "rs422":
                   logging.info(self.loginfo + "RS422 connected to " + self.comms.port)
01762
01763
01764
          def dumpRegisters(self):
01765
01766
               *DEPRECATED* use dumpStatus() instead
01767
01768
              List contents of all registers in board.registers. *WARNING* some status flags
01769
                will reset when read.
01770
01771
              Returns:
01772
                  Sorted list: [register name (register address) : register contents as
01773
                    hexadecimal string without '0x']
01774
01775
              dump = { } { }
01776
              for key in self.board.registers.keys():
01777
                  err, rval = self.getRegister(key)
01778
                   dump[key] = rval
01779
              reglistmax = int(max(self.board.registers.values()), 16)
```

```
01780
               dumplist = [0] * (reglistmax + 1)
01781
              for k, v in dump.items():
01782
                   regnum = self.board.registers[k]
01783
                   dumplist[int(regnum, 16)] = (
01784
                       "(" + regnum + ") {0:<24} {1}".format(k, v.upper())
01785
01786
               reglist = [a for a in dumplist if a]
01787
              return reglist
01788
01789
          def dumpSubregisters(self):
01790
01791
               *DEPRECATED* use dumpStatus() instead
01792
01793
              List contents of all subregisters in board.channel_lookups and
01794
                board.monitor_lookups.
01795
               \starWARNING\star some registers will reset when read; only the first subregister from
01796
                such a register will return the correct value, the remainder will return zeros
01797
01798
              Returns:
01799
                  dictionary {subregister name : subregister contents as binary string
01800
                    without initial '0b'}
01801
01802
              dump = \{ \}
01803
              for sub in self.board.subreglist:
01804
                  kev = sub.name
01805
                   err, resp = self.getSubregister(key)
01806
                   if err:
01807
                      logging.warning(
                           self.logwarn + "dumpSubregisters: unable to read subregister " + key
01808
01809
01810
                  val = hex(int(resp, 2))
                  dump[key] = val
01811
01812
              return dump
01813
01814
          def mmReadoff(self, waitOnSRAM, variation=None):
01815
01816
              Convenience function for parsing frames for use by MicroManager plugin
01817
              Args:
01818
                   waitOnSRAM: readoff wait flag
                   variation: format of frames generated from readoff
01819
01820
                       default - return first frame only
                       "LastFrame" - return last frame only
01821
                       "Average" - provide average of frames as single frame
"Landscape" - stitch frames together horizontally into single wide frame
01822
01823
01824
01825
              Returns:
               ndarray - single image frame
01826
01827
              frames, datalen, data_err = self.readoff(waitOnSRAM)
if variation == "LastFrame":
01828
01829
01830
                   return frames[self.sensor.nframes - 1]
01831
              elif variation == "Average":
01832
                   return np.sum(frames, axis=0) // self.sensor.nframes
01833
              elif variation == "Landscape":
01834
                   shaped = [
01835
                       np.reshape(frame, (self.sensor.maxheight, self.sensor.maxwidth))
01836
                       for frame in frames
01837
01838
                  return np.concatenate(shaped, axis=1)
01839
01840
                  return frames[0]
01841
01842
          def setFrames(self, minframe=None, maxframe=None):
01843
01844
              Sets bounds on frames returned by board, inclusive (e.g., 0,3 returns four
01845
              frames). If called without parameters, resets to full set of frames.
01846
01847
              Args:
01848
                  minframe: first frame to read from board
01849
                  maxframe: last frame to read from board
01850
01851
              Returns:
              Error string
01852
01853
01854
              logging.debug(
01855
                  self.logdebug
                   + "setFrames: minframe = "
01856
01857
                   + str(minframe)
                   + "; maxframe =
01858
01859
                   + str(maxframe)
01860
```

```
01861
              if minframe is None:
01862
                 minframe = self.sensor.minframe
01863
              if maxframe is None:
01864
                  maxframe = self.sensor.maxframe
01865
              if (
01866
                  not isinstance(minframe, int)
01867
                  or minframe < self.sensor.minframe</pre>
01868
                  or minframe > maxframe
01869
                  or not isinstance(maxframe, int)
01870
                  or maxframe > self.sensor.maxframe
01871
             ):
01872
01873
                      self.logerr + "setFrames: invalid frame limits submitted. Frame "
01874
                      "selection remains unchanged. "
01875
01876
                  logging.error(err)
01877
                  return err
01878
01879
              initframe = hex(minframe)[2:].zfill(8)
01880
              finframe = hex(maxframe)[2:].zfill(8)
              err1, _ = self.setRegister("FPA_FRAME_INITIAL", initframe)
01881
01882
              err2, _ = self.setRegister("FPA_FRAME_FINAL", finframe)
01883
              self.sensor.firstframe = minframe
01884
              self.sensor.lastframe = maxframe
01885
              self.sensor.nframes = maxframe - minframe + 1
01886
              self.comms.payloadsize = (
01887
                  self.sensor.width
01888
                  * self.sensor.height
01889
                  * self.sensor.nframes
01890
                  * self.sensor.bytesperpixel
01891
              plural = ""
01892
              if self.sensor.nframes > 1:
    plural = "s"
01893
01894
              logging.info(
01895
01896
                 self.loginfo
01897
                  + "Readoff set to "
01898
                  + str(self.sensor.nframes)
                  + " frame"
01899
01900
                  + plural
01901
                  + " ("
01902
                  + str(minframe)
01903
01904
                  + str(maxframe)
01905
                  + ")"
01906
01907
              err = err1 + err2
01908
              if err:
01909
                  logging.error(
01910
                      self.logerr + "setFrames may not have functioned properly: " + err
01911
01912
              return err
01913
01914
          def setRows(self, minrow=0, maxrow=None, padToFull=False):
01915
01916
              Sets bounds on rows returned by board, inclusive (e.g., 0,1023 returns all 1024
01917
               rows). If called without parameters, resets to full image size.
01918
01919
01920
                 minrow: first row to return from board
                  maxrow: last row to return from board
01921
01922
                  padToFull: if True, generate full size frames, padding collected rows with
01923
                   zeroes if necessary
01924
01925
              logging.debug(
01926
                  self.logdebug
01927
                  + "setRows: minrow = "
01928
                  + str(minrow)
01929
                  + "; maxrow =
01930
                  + str(maxrow)
01931
                  + "; padToFull = "
01932
                  + str(padToFull)
01933
01934
              if maxrow is None:
01935
                  maxrow = self.sensor.maxheight - 1
01936
              if (
01937
                  not isinstance(minrow, int)
01938
                  or minrow < 0
01939
                  or minrow > maxrow
01940
                  or not isinstance (maxrow, int)
01941
                  or maxrow >= self.sensor.maxheight
```

```
01942
01943
01944
                       self.logerr + "setRows: invalid row arguments submitted. Frame size"
                       " remains unchanged. "
01945
01946
01947
                   logging.error(err)
01948
                   return err
01949
01950
               initrow = hex(minrow)[2:].zfill(8)
              finrow = hex(maxrow)[2:].zfill(8)
01951
              err1, _ = self.setRegister("FPA_ROW_INITIAL", initrow)
01952
01953
              err2, _ = self.setRegister("FPA_ROW_FINAL", finrow)
01954
              self.sensor.firstrow = minrow
              self.sensor.lastrow = maxrow
01955
01956
              self.sensor.height = maxrow - minrow + 1
01957
              self.comms.payloadsize = (
01958
                  self.sensor.width
01959
                   * self.sensor.height
01960
                   * self.sensor.nframes
01961
                   * self.sensor.bytesperpixel
01962
              )
01963
01964
              if self.commname == "rs422":
01965
                  self.comms.datatimeout = (
01966
                       (1.0 * self.sensor.height / self.sensor.maxheight)
                       * 5e7
01967
                       * self.sensor.nframes
01968
                       / self.comms.baud
01969
01970
01971
              self.padToFull = padToFull
01972
              logging.info(
01973
                  self.loginfo
01974
                   + "Readoff set to "
01975
                  + str(self.sensor.height)
                  + " rows ("
01976
01977
                  + str(minrow)
01978
                  + ", "
01979
                  + str(maxrow)
01980
                   + ")"
01981
              )
01982
              err = err1 + err2
01983
              if err:
01984
                  logging.error(
01985
                       self.logerr + "setRows may not have functioned properly: " + err
01986
01987
              return err
01988
01989
          def abortReadoff(self, flag=True):
01990
01991
              Simple abort command for readoff in waiting mode--does not interrupt download in
01992
                  progress. Requires external threading to function. *WARNING* if not
01993
                  intercepted by active readoff command, will terminate next readoff command
01994
                  immediately at inception.
01995
01996
                  flag: Sets passive abort flag read by readoff command
01997
              Returns:
               boolean: updated setting of flag
01998
01999
02000
              logging.info(self.loginfo + "abortReadoff")
02001
              self.abort = flag
              return flag
02002
02003
02004
          def batchAcquire(
02005
              self,
02006
              sets=1,
               trig="Hardware",
02007
              path=None,
02008
               filename="Frame",
02009
02010
              prefix=None,
02011
              showProgress=0,
02012
         ):
02013
02014
              Acquire a series of images as fast as possible, then process and save to disk.
02015
               *WARNING* This method stores images in RAM, so the number of sets that can be
02016
                acquired in a single call is limited by available memory.
02017
02018
              Args:
                  sets: Number of acquisitions to perform trig: trigger type; 'hardware', 'software', or 'dual'
02019
02020
                  path: save path, defaults to './output' filename: defaults to 'frames.bin'
02021
02022
```

```
02023
                   prefix: prepended to filename, defaults to time/date (e.g. '160830-124704_')
02024
                     DO NOT USE unless providing a varying value (a fixed prefix will cause
02025
                     overwriting)
                   showProgress: if non-zero, show notice every 'showProgress' acquisitions and
02026
02027
                    print total acquisition time
02028
02029
              Time taken for acquisition (seconds)
02030
02031
02032
              logging.debug(
02033
                  self.logdebug
02034
                   + "batchAcquire: sets = "
02035
                   + str(sets)
                   + "; trig = "
02036
02037
                  + str(trig)
02038
                   + "; path = "
02039
                   + str(path)
02040
                  + "; filename =
02041
                   + str(filename)
02042
                   + "; prefix =
02043
                   + str(prefix)
02044
                   + "; showProgress = "
02045
                   + str(showProgress)
02046
02047
              datalist = ["0"] * sets
              timelist = [datetime.now()] * sets
02048
02049
              logging.info(
02050
                  self.loginfo
02051
                   + "batchAcquire: temporarily disabling warning and information logging "
02052
02053
              logging.getLogger().setLevel(self.verbmap.get(2))
02054
              beforeread = time.time()
02055
              for i in range (sets):
                  if showProgress and not (i + 1) % showProgress:
    print(self.loginfo + "batchAcquire: Acquiring set " + str(i + 1))
02056
02057
02058
                   self.arm(trig)
02059
                   data, datalen, data_err = self.readoff(fast=True)
02060
                  datalist[i] = data
timelist[i] = datetime.now()
02061
02062
              afterread = time.time()
02063
              if showProgress:
02064
                   print(
                       self.loginfo
02065
                       + "batchAcquire: "
02066
                       + str(afterread - beforeread)
+ " seconds for "
02067
02068
02069
                       + str(sets)
+ " sets"
02070
02071
                  )
02072
              setnum = 0
              if path is None:
02073
02074
                   path = os.path.join(os.getcwd(), "output")
02075
              for imset, imtime in zip(datalist, timelist):
02076
                  setnum = setnum + 1
02077
                   if showProgress and not setnum % showProgress:
02078
                      print(self.loginfo + "batchAcquire: Saving set " + str(setnum))
02079
                   parsed = generateFrames(self, imset)
02080
                   if prefix is None:
02081
                       setprefix = imtime.strftime("%y%m%d-%H%M%S%f")[:-2] + "_"
02082
                   else:
02083
                      setprefix = prefix
02084
                   self.saveTiffs(parsed, path, filename, prefix=setprefix)
02085
               logging.getLogger().setLevel(self.verblevel)
              logging.info(self.loginfo + "batchAcquire: re-enabling logging")
02086
02087
              return afterread - beforeread
02088
          # TODO: should this be just a flag for readoff instead of a distinct method?
02089
02090
          # TODO: make sure this handles single frames (list made already?), text frames
02091
          # TODO: add documentation
02092
          def saveHDF(
02093
              self,
02094
              frames,
              path=None,
02095
               filename="Acquisition",
02096
              prefix=None,
02097
02098
          ):
              ... ...
02099
02100
              logging.info(self.loginfo + ": saveHDF")
02101
              err =
              if path is None:
02102
02103
                   path = os.path.join(os.getcwd(), "output")
```

```
if prefix is None:
02104
                 prefix = datetime.now().strftime("%y%m%d-%H%M%S%f")[:-5] + "_"
              if not os.path.exists(path):
02107
                 os.makedirs(path)
02108
02109
             h5file = os.path.join(path, prefix + filename + ".hdf5")
02110
             with h5py.File(h5file, "w") as f:
02111
                # shotgrp = f.create_group("DATA/SHOT")
                 frame_index = 0
02112
02113
                 for frame in frames:
                     grp = f.create_group("DATA/SHOT/FRAME_0" + str(frame_index))
                      data = grp.create_dataset(
02115
02116
                          "DATA", (self.sensor.height, self.sensor.width), data=frame
02117
02118
                      frame_index += 1
02119
02120
02121 """
02122 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
02123 LLNL-CODE-838080
02124
02125 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
02126 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
02127 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
02128 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
02129 be made under this license.
02130 """
```

# 7.17 C:/Users/hill35/git/camera\_python/nsCamera/comms/GigE.py File Reference

### Classes

- · class nsCamera.comms.GigE.GigE
- class nsCamera.comms.GigE.GigE.ZESTETM1\_CARD\_INFO

### **Namespaces**

- namespace nsCamera
- namespace nsCamera.comms
- namespace nsCamera.comms.GigE

## 7.18 **GigE.py**

#### Go to the documentation of this file.

```
00001 # -*- coding: utf-8 -*-
00002 """
00003 Gigabit Ethernet interface for nsCamera.
00004
00005 Author: Jeremy Martin Hill (jerhill@llnl.gov)
00006
00007 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00008 LLNL-CODE-838080
00009
00010 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00011 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00012 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00013 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00015
00016 Version: 2.1.2 (February 2025)
```

```
00017 """
00018
00019 import ctypes as C
00020 import logging
00021 import os.path
00022 import sys
00023
00024 from nsCamera.utils.misc import generateFrames, str2bytes, bytes2str
00025
00026
00027 class GigE:
00028
00029
          Code to manage Gigabit Ethernet connection to board. Each GigE object manages a
            single OT card; to use multiple cards, instantiate multiple cameraAssembler
00031
            objects, each specifying the unique IPs of the corresponding OT card.
00032
00033
          Note: Orange Tree card must be configured before use. See the README for details
00034
00035
          Exposed methods:
00036
               arm() - puts camera into wait state for external trigger
00037
               readFrames() - waits for data ready register flag, then copies camera image data
00038
                into numpy arrays
00039
               readoff() - waits for data ready register flag, then copies camera image data
00040
                into numpy arrays; returns payload, payload size, and error message
00041
               sendCMD(pkt) - sends packet object via serial port
readSerial(size, timeout) - read 'size' bytes from connection
00042
00043
               writeSerial(outstring) - submits string 'outstring' over connection
00044
               closeDevice() - close connections and free resources
00045
               getCardIP() - returns IP address of OT card
              getCardInfo() - prints report of details of OT card and connection
00046
00047
00048
00049
          def __init__(self, camassem):
00050
               Aras:
00051
               camassem: parent cameraAssembler object
00052
00053
00054
               self.ca = camassem
               self.logcrit = self.ca.logcritbase + "[GigE] "
self.logerr = self.ca.logerrbase + "[GigE] "
00055
00056
               self.logwarn = self.ca.logwarnbase + "[GigE] "
self.loginfo = self.ca.loginfobase + "[GigE] "
00057
00058
               self.logdebug = self.ca.logdebugbase + "[GigE] "
00059
00060
               logging.info(self.loginfo + "Initializing GigE comms object")
00061
               self.mode = 1
00062
               self.writeTimeout = 10000
00063
               self.readTimeout = 10000
00064
               self.payloadsize = (
00065
                   self.ca.sensor.width
00066
                   * self.ca.sensor.height
00067
                   * self.ca.sensor.nframes
00068
                   * self.ca.sensor.bytesperpixel
00069
00070
00071
                   self.logdebug + "Payload size: " + str(self.payloadsize) + " bytes"
00072
00073
               self.skipError = False
00074
00075
               self.ZErrorDict = {
00076
                   0x8000: "Socket Error",
00077
                   0x8001: "Internal Error",
00078
                   0x8002: "Illegal Status Code",
00079
                   0x8003: "Null Parameter",
                   0x8004: "Out of Memory",
00080
                   0x8005: "Invalid Connection Type",
                   0x8006: "Illegal Connection",
00082
                   0x8007: "Socket Closed Unexpectedly",
00083
00084
                   0x8008: "Timeout",
                   0x8009: "Illegal Parameter",
00085
00086
               }
00087
00088
               if self.ca.port:
00089
                   logging.debug(
00090
                       self.logdebug + "Port supplied to GigE.py: " + str(self.ca.port)
00091
00092
                   if isinstance(self.ca.port, int) and 0 < self.ca.port < 65536:</pre>
00093
                       self.dport = self.ca.port
00094
                   else:
00095
                       logging.error(
                            self.logerr + "Invalid port number supplied, defaulting to "
00096
00097
```

7.18 GigE.py 289

```
00098
00099
                      self.dport = 20482
00100
              else:
                  self.dport = 20482 # default
00101
00102
00103
              self.ca.port = self.dport
00104
              logging.debug(self.logdebug + "Port used by GigE.py: " + str(self.dport))
00105
00106
              logging.debug(self.logdebug + "CPU architecture: " + str(self.ca.arch))
              if self.ca.arch == "64bit":
00107
00108
                  arch = "64"
              else:
00109
00110
                  arch = "32"
00111
00112
              logging.debug(self.logdebug + "Operating system: " + str(self.ca.platform))
00113
              if self.ca.platform == "Windows":
                  lib_name = "ZestETM1.dll"
00114
              elif self.ca.platform == "Linux" or self.ca.platform == "Darwin":
00115
                  lib_name = "libZestETM1.so"
00116
00117
              else:
00118
                  logging.warning(
                      self.logwarn + "System does not self-identify as Linux, Windows, "
00119
00120
                       "or Mac. Assuming posix-style libraries
00121
00122
                  lib name = "libZestETM1.so"
00123
              self.closecard = False
00124
00125
              libpath = os.path.join(self.ca.packageroot, "comms", "ZestETM1", arch, lib_name)
00126
00127
              self._zest = C.CDLL(libpath)
00128
              self.CardInfo = self.ZESTETM1_CARD_INFOZESTETM1_CARD_INFO()
00129
00130
              self.CardInfoP = C.pointer(self.CardInfo)
00131
00132
              # functions
              self.ZCountCards = self._zest.ZestETM1CountCards
00133
              self.ZCountCards.argtypes = [
00134
00135
                  C.POINTER(C.c_ulong),
                  C.POINTER(C.POINTER(self.ZESTETM1_CARD_INFOZESTETM1_CARD_INFO)),
00136
00137
                  C.c_int,
00138
00139
00140
              self.ZOpenConnection = self._zest.ZestETM1OpenConnection
00141
              self.ZOpenConnection.argtypes = [
00142
                  C.POINTER(self.ZESTETM1_CARD_INFOZESTETM1_CARD_INFO),
00143
                  C.c_int,
00144
                  C.c_ushort,
00145
                  C.c_ushort,
00146
                  C.POINTER(C.c_void_p),
00147
00148
00149
              self.ZWriteData = self._zest.ZestETM1WriteData
00150
              self.ZWriteData.argtypes = [
00151
                  C.c_void_p,
00152
                  C.c_void_p,
00153
                  C.c_ulong,
                  C.POINTER(C.c_ulong),
00154
00155
                  C.c_ulong,
00156
00157
00158
              self.ZReadData = self._zest.ZestETM1ReadData
00159
              self.ZReadData.argtypes = [
00160
                  C.c_void_p,
00161
                  C.c_void_p,
00162
                  C.c_ulong,
00163
                  C.POINTER(C.c_ulong),
00164
                  C.c_ulong,
00165
              ]
00166
00167
              self.Connection = C.c_void_p()
00168
              self.openDevice(self.ca.timeout)
00169
00170
          def sendCMD(self, pkt):
00171
00172
              Submit packet and verify the response packet.
00173
              Packet communications with FPGA omit CRC suffix, so adds fake CRC bytes to
00174
               response
00175
00176
              Args:
                  pkt: Packet object
00177
00178
```

```
00179
              tuple (error, response string)
00180
00181
00182
00183
              pktStr = pkt.pktStr()[0:16]
00184
               logging.debug(self.logdebug + "sendCMD packet: " + str(pktStr))
00185
00186
              self.ca.writeSerial(pktStr)
00187
              if (
                  hasattr(self.ca, "board")
and pktStr[4] == "0"
00188
00189
00190
                  and pktStr[5:8] == self.ca.board.registers["SRAM_CTL"]
00191
00192
                  bufsize = self.payloadsize + 16
00193
                  resptext = self.readSerial(bufsize)
00194
00195
                  if len(resptext) < 32:</pre>
00196
                      logging.debug(self.logdebug + "sendCMD resptext = " + str(resptext))
                  else:
00197
00198
                      logging.debug(
00199
                           self.logdebug
00200
                           + "sendCMD resptext (truncated) = "
00201
                           + str(resptext)[0:32]
00202
                      )
00203
00204
                  if len(resptext) < bufsize + 16:</pre>
00205
                      err += (
00206
                          self.logerr + "sendCMD- packet too small, payload may be incomplete"
00207
00208
                      logging.error(err)
00209
              else:
00210
                  # workaround for initial setup before board object has been initialized
00211
                  resp = self.readSerial(8)
                  logging.debug(self.logdebug + "sendCMD response: " + str(resp))
00212
00213
                  if len(resp) < 8:
                      err += self.logerr + "sendCMD- response too small, returning zeros" resptext = "000000000000000000000"
00214
00215
00216
                      logging.error(err)
00217
                  else:
00218
                      resptext = resp + "0000"
00219
              return err, resptext
00220
00221
          def arm(self, mode):
00222
00223
              Puts camera into wait state for trigger. Mode determines source; defaults to
00224
                'Hardware'
00225
00226
              Args:
               mode:
00227
                          'Software' | 'S' activates software, disables hardware triggering
00228
                           'Hardware' | 'H' activates hardware, disables software triggering
00229
                            Hardware is the default
00230
00231
              Returns:
              tuple (error, response string)
00232
00233
00234
              if not mode:
00235
                  mode = "Hardware"
00236
                  logging.info(self.loginfo + "arm")
00237
              logging.debug(self.logdebug + "arming mode: " + str(mode))
00238
              self.ca.clearStatus()
00239
              self.ca.latchPots()
00240
              err, resp = self.ca.startCapture(mode)
00241
              if err:
00242
                  logging.error(self.logerr + "unable to arm camera")
00243
              else:
00244
                  self.ca.armed = True
00245
                  self.skipError = True
00246
              return err, resp
00247
00248
          def readFrames(self, waitOnSRAM, timeout=0, fast=False, columns=1):
00249
00250
              Copies image data from board into numpy arrays.
00251
00252
              Args:
00253
                  waitOnSRAM: if True, wait until SRAM_READY flag is asserted to begin copying
00254
                    data
                  timeout: passed to waitForSRAM; after this many seconds begin copying data
00255
                     irrespective of SRAM_READY status; 'zero' means wait indefinitely
00256
00257
                     WARNING: If acquisition fails, the SRAM will not contain a current image,
00258
                      but the code will copy the data anyway
00259
                  fast: if False, parse and convert frames to numpy arrays; if True, return
```

7.18 GigE.py 291

```
00260
                     unprocessed text stream
00261
                   columns: 1 for single image per frame, 2 for separate hemisphere images
00262
00263
              Returns:
00264
                  list of numpy arrays OR raw text stream
00265
00266
00267
              frames, _, _ = self.readoff(waitOnSRAM, timeout, fast, columns)
00268
              return frames
00269
00270
          def readoff(self, waitOnSRAM, timeout=0, fast=False, columns=1):
00271
              Copies image data from board into numpy arrays; returns data, length of data, and error messages. Use 'readFrames()' unless you require this additional
00272
00273
00274
              information
00275
00276
              Args:
00277
                  waitOnSRAM: if True, wait until SRAM_READY flag is asserted to begin copying
00278
                    data
00279
                   timeout: passed to waitForSRAM; after this many seconds begin copying data
00280
                     irrespective of SRAM_READY status; 'zero' means wait indefinitely
00281
                     WARNING: If acquisition fails, the SRAM will not contain a current image,
00282
                      but the code will copy the data anyway
00283
                   fast: if False, parse and convert frames to numpy arrays; if True, return
00284
                    unprocessed text stream
00285
                  columns: 1 for single image per frame, 2 for separate hemisphere images
00286
00287
              Returns:
00288
                  tuple (list of numpy arrays OR raw text stream, length of downloaded payload
00289
                    in bytes, payload error flag) since CRC check is handled by TCP/IP,
00290
                    payload error flag is always False for GigE
00291
00292
              logging.info(self.loginfo + "readoff")
00293
              logging.debug(
00294
                  self.logdebug
                   + "readoff: waitonSRAM = "
00295
00296
                  + str(waitOnSRAM)
00297
                  + "; timeout =
00298
                  + str(timeout)
00299
                  + "; fast =
00300
                   + str(fast)
00301
00302
              \ensuremath{\text{\#}} Wait for data to be ready on board
               # Skip wait only if explicitly tagged 'False' ('None' defaults to True)
00303
00304
               if waitOnSRAM is not False:
00305
                   self.ca.waitForSRAM(timeout)
00306
              self.skipError = False
00307
              err, rval = self.ca.readSRAM()
00308
              if err:
00309
                   logging.error(self.logerr + "Error detected in readSRAM")
00310
              elif self.ca.boardname == "llnl_v4":
                  # self.ca.setSubregister('SWACK','1')
00311
00312
00313
               # extract the data. Remove header; the FPGA returns a packet without the CRC
00314
                  suffix
00315
              logging.debug(self.logdebug + "readoff: first 64 chars: " + str(rval[0:64]))
00316
              data = rval[32:]
00317
               if fast:
00318
                  return data, len(data) // 2, bool(err)
00319
               else:
00320
                  parsed = generateFrames(self.ca, data, columns)
00321
                   return parsed, len(data) // 2, bool(err)
00322
          def writeSerial(self, outstring, timeout):
00323
00324
00325
              Transmit string to board
00326
00327
              Args:
00328
                  outstring: string to write
00329
                   timeout: serial timeout in sec (defaults to self.writeTimeout)
00330
00331
              Returns:
              integer number of bytes written
00332
00333
00334
              logging.debug(
00335
                  self.logdebug
                   + "writeSerial: outstring = "
00336
00337
                   + str(outstring)
                   + "; timeout =
00338
                   + str(timeout)
00339
00340
              )
```

```
00341
              if not timeout:
00342
                  timeout = self.writeTimeout
00343
              outstring = str2bytes(outstring)
              outbuff = C.create_string_buffer(outstring)
00344
              outbuffp = C.pointer(outbuff)
00345
00346
              outbufflen = len(outstring)
              writelen = C.c_ulong(0)
00347
00348
              err = self.ZWriteData(
00349
                  self.Connection, outbuffp, outbufflen, C.byref(writelen), timeout
00350
00351
              if err:
                  if err == 0x4000:
00352
00353
                      logging.warning(
                          self.logerr + "OT Card emitted an undefined warning message"
00354
00355
00356
                  else:
00357
                      logging.error(
00358
                          self.logerr + "writeSerial error: " + self.ZErrorDict[err]
00359
00360
              logging.debug(self.logdebug + "writeSerial: writelen = " + str(writelen))
00361
              return writelen
00362
00363
          def readSerial(self, size, timeout=None):
00364
00365
              Read bytes from the serial port. Does not verify packets.
00366
00367
              Aras:
00368
                 size: number of bytes to read
00369
                 timeout: serial timeout in sec (defaults to self.readTimeout)
00370
00371
              Returns:
              tuple (error string, string read from serial port) """
00372
00373
00374
              logging.debug(
00375
                  self.logdebug
00376
                  + "readSerial: size = "
00377
                  + str(size)
                  + "; timeout = "
00378
00379
                  + str(timeout)
00380
00381
              if not timeout:
00382
                  timeout = self.readTimeout
              inbuff = C.create_string_buffer(size + 1)
00383
00384
              inbuffp = C.pointer(inbuff)
              readlen = C.c_ulong(0)
00385
00386
              err = self.ZReadData(self.Connection, inbuffp, size, C.byref(readlen), timeout)
00387
              if err:
00388
                  if self.skipError:
00389
                      logging.debug(
00390
                           self.logdebug + "readSerial: skipped error: " + self.ZErrorDict[err]
00391
00392
                      self.skipError = False
00393
                  elif err == 0x4000:
00394
                       logging.warning(
00395
                          self.logerr + "OT Card emitted an undefined warning message"
00396
00397
00398
                      logging.error(self.logerr + "readSerial error: " + self.ZErrorDict[err])
              return bytes2str(inbuff.raw)[:-2]
00399
00400
00401
          # TODO: check for valid timeout, probably in init in CameraAssembler
00402
          def openDevice(self, timeout=30):
00403
00404
              Find Orange Tree card and open a connection; if IP is supplied as parameter for
00405
                the CameraAssembler, bypass network search and connect directly to indicated
00406
                IP address
00407
00408
              Args:
              timeout: timeout in seconds for attempting to connect to a card \ensuremath{\tt^{muu}}
00409
00410
00411
              err = self._zest.ZestETM1Init()
00412
              if err:
00413
                  logging.critical(self.logcrit + "ZestETM1Init failure")
00414
                  svs.exit(1)
              logging.info(self.loginfo + "searching for Orange Tree cards")
00415
00416
              NumCards = C.c_ulong(0)
00417
00418
              if self.ca.iplist:
                  ubyte4 = C.c_ubyte * 4
00419
                  self.CardInfo.IPAddr = ubyte4(*self.ca.iplist)
00420
00421
                  self.CardInfo.ControlPort = C.c ushort(self.dport)
```

7.18 GigE.py 293

```
00422
                  self.CardInfo.Timeout = C.c_ulong(self.writeTimeout)
00423
                  self.closecard = False
00424
              else:
                  wait = 0
00425
00426
                  while True:
00427
                      logging.debug(
00428
                           self.logdebug + "openDevice: connection wait = " + str(wait)
00429
00430
                      err = self.ZCountCards(C.byref(NumCards), C.byref(self.CardInfoP), 1000)
                      if err:
00431
00432
                           logging.critical(self.logcrit + "CountCards failure")
00433
                           sys.exit(1)
00434
                       if NumCards.value > 0:
00435
                          break
00436
                       if wait == timeout:
00437
                           logging.critical(
00438
                               "{}No Orange Tree cards found in {} seconds".format(
00439
                                   self.logcrit, timeout
00440
00441
                           )
00442
                          sys.exit(1)
00443
                       elif not wait % 5:
00444
                           logging.info(
00445
                               "{}Still trying to connect after {} seconds...".format(
00446
                                   self.loginfo, wait
00447
00448
                          )
00449
                      wait += 1
00450
                  logging.info(
                      self.loginfo + str(NumCards.value) + " Orange Tree card(s) found"
00451
                  ) # TODO: add check for GigE bit in board description
00452
00453
              err = self.ZOpenConnection(
00454
                  self.CardInfoP, 0, self.dport, 0, C.byref(self.Connection)
00455
00456
              if err:
                  if err == 0x4000:
00457
00458
                      logging.warning(
                           self.logerr + "OT Card emitted an undefined warning message"
00459
00460
00461
                  else:
00462
                      logging.critical(
                           self.logcrit + "OpenConnection failure: " + self.ZErrorDict[err]
00463
00464
00465
                  svs.exit(1)
00466
00467
          def closeDevice(self):
00468
00469
              Close connection to Orange Tree card and free resources
00470
00471
              logging.debug(self.logdebug + "Closing connection to Orange Tree card")
00472
              self._zest.ZestETM1CloseConnection(self.Connection)
00473
              if self.closecard:
00474
                  try:
00475
                      self._zest.ZestETM1FreeCards(self.CardInfoP)
00476
                  except SystemExit:
00477
00478
                  except KeyboardInterrupt:
00479
                      raise
00480
                  except Exception:
00481
                      logging.error(self.logerr + "Error reported in OT card closure")
00482
              self._zest.ZestETM1Close()
00483
00484
          def getCardIP(self):
00485
00486
              Query IP address of OT card
00487
00488
              Returns: address of OT card as list of bytes
00489
00490
              return self.CardInfo.IPAddr
00491
00492
          # TODO: use logging.info, with override option?
          def getCardInfo(self):
00493
00494
00495
              Prints status message with information returned by OT card
00496
00497
              ci = self.CardInfoP.contents
              print("GigE Card Status:")
00498
              print("----
00499
              print("IP: " + ".".join(str(e) for e in [b for b in ci.IPAddr]))
00500
              print("ControlPort: " + str(ci.ControlPort))
print("Timeout: " + str(ci.Timeout))
00501
00502
```

```
print("HTTPPort: " + str(ci.HTTPPort))
print("MACAddr: " + ".".join(format(e, "02X") for e in [b for b in ci.MACAddr]))
print("SubNet: " + ".".join(str(e) for e in [b for b in ci.SubNet]))
print("Gateway: " + ".".join(str(e) for e in [b for b in ci.Gateway]))
print("SerialNumber: " + str(ci.SerialNumber))
print("FirmwareVersion: " + str(ci.FirmwareVersion))
00505
00506
00507
00508
                   print("HardwareVersion: " + str(ci.HardwareVersion))
00509
00510
00511
          class ZESTETM1_CARD_INFO(C.Structure):
00513
                 ubyte4 = C.c_ubyte * 4
00514
                  ubyte6 = C.c_ubyte * 6
                  _fields_ = [
                       ("IPAddr", ubyte4),
                        ("ControlPort", C.c_ushort),
00517
                       ("Timeout", C.c_ulong),
("HTTPPort", C.c_ushort),
00519
                       ("MACAddr", ubyte6),
("SubNet", ubyte4),
("Gateway", ubyte4),
00520
00521
00522
                        ("SerialNumber", C.c_ulong),
("FirmwareVersion", C.c_ulong),
00523
00524
00525
                        ("HardwareVersion", C.c_ulong),
00526
00527
00528
00529 """
00530 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00531 LLNL-CODE-838080
00533 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00534 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00535 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00536 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00537 be made under this license.
00538 """
```

# 7.19 C:/Users/hill35/git/camera\_python/nsCamera/comms/RS422.py File Reference

### **Classes**

class nsCamera.comms.RS422.RS422

### **Namespaces**

- namespace nsCamera
- namespace nsCamera.comms
- namespace nsCamera.comms.RS422

## 7.20 RS422.py

#### Go to the documentation of this file.

```
00001 # -*- coding: utf-8 -*-
00002 """
00003 RS422 driver for nsCamera
00004
00005 Author: Brad Funsten (funstenl@llnl.gov)
00006 Author: Jeremy Martin Hill (jerhill@llnl.gov)
00007
00008 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
```

7.20 RS422.py 295

```
00010
00011 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00012 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00013 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00014 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00015 be made under this license.
00016
00017 Version: 2.1.2 (February 2025)
00018 """
00019
00020 import logging
00021 import sys
00022 import time # to time the script
00024 import serial
00025 import serial.tools.list_ports # for RS422 serial link setup
00027 from nsCamera.utils.misc import generateFrames, str2bytes, bytes2str, checkCRC
00028
00029
00030 class RS422:
00031
00032
          Code to manage RS422 connection. Will automatically query available COM interfaces
            until a board is found. Use the 'port=x' parameter in cameraAssembler call to
00033
00034
            specify a particular COM interface.
00035
00036
          Exposed methods:
00037
              arm() - Puts camera into wait state for external trigger
00038
              readFrames() - waits for data ready register flag, then copies camera image data
00039
                into numpy arrays
00040
              readoff() - waits for data ready register flag, then copies camera image data
00041
                into numpy arrays; returns payload, payload size, and error message
              sendCMD(pkt) - sends packet object via serial port
readSerial(size, timeout) - read 'size' bytes from serial port
00042
00043
              writeSerial(cmd) - submits string 'cmd' (assumes string is preformed packet)
00044
00045
              closeDevice() - close serial connections
00046
00047
00048
          def __init__(self, camassem, baud=921600, par="0", stop=1):
00049
00050
              Aras:
00051
                  camassem: parent cameraAssembler object
00052
                  baud: bits per second
00053
                  par: parity type
              stop: number of stop bits
00054
00055
00056
              self.ca = camassem
              self.logcrit = self.ca.logcritbase + "[RS422] "
00057
00058
              self.logerr = self.ca.logerrbase + "[RS422] "
              self.logwarn = self.ca.logwarnbase + "[RS422] "
self.loginfo = self.ca.loginfobase + "[RS422] "
00059
00060
              self.logdebug = self.ca.logdebugbase + "[RS422] "
00061
00062
              logging.info(self.loginfo + "initializing RS422 comms object")
00063
              logging.debug(
00064
                  self.logdebug
00065
                   + "Init: baud = "
00066
                   + str(baud)
00067
                  + "; par = '
00068
                   + str(par)
00069
                   + "; stop = "
00070
                   + str(stop)
00071
00072
              self.mode = 0
00073
              self.baud = baud # Baud rate (bits/second)
00074
              self.par = par # Parity bit
              self.stop = stop # Number of stop bits
00075
00076
              self.read_timeout = 1 # default timeout for ordinary packets
00077
              self.write_timeout = 1
00078
              # TODO: make datatimeout a cameraAssembler parameter
              self.datatimeout = 60 # timeout for data read
00080
              logging.debug(
00081
                  self.logdebug + "Data timeout = " + str(self.datatimeout) + " seconds"
00082
00083
              self.PY3 = sys.version_info > (3,)
00084
              self.skipError = False
              port = ""
00085
              ports = list(serial.tools.list_ports.comports())
00086
00087
              logging.debug(self.logdebug + "Comports: " + str(ports))
00088
              for p, desc, add in ports:
00089
                   if self.ca.port is None or p == "COM" + str(self.ca.port):
                       logging.info(self.loginfo + "found comm port " + p)
00090
```

```
00091
                      try:
00092
                           with serial.Serial(
00093
                              p,
                               self.baud,
00094
00095
                              parity=self.par,
00096
                               timeout=0.01,
00097
                               write_timeout=0.01,
00098
                           ) as ser:
00099
                              ser.write(str2bytes("aaaa1000000000001a84"))
00100
                              time.sleep(1)
00101
                               s = ser.read(10)
00102
                              resp = bytes2str(s)
00103
                              logging.debug(self.logdebug + "Init response: " + str(resp))
00104
                              if (
00105
                                  resp[0:5].lower() == "aaaa9"
00106
                              ): # TODO: add check for RS422 bit in board description
00107
                                  boardid = resp[8:10]
                                   if boardid == "00":
00108
00109
                                       logging.critical(
                                          self.logcrit + "SNLrevC board detected - not "
00110
                                           "compatible with nsCamera >= 2.0"
00111
00112
00113
                                       svs.exit(1)
                                   elif boardid == "81":
00114
00115
                                       logging.info(self.loginfo + "LLNLv1 board detected")
                                   elif boardid == "84":
00116
                                      logging.info(self.loginfo + "LLNLv4 board detected")
00117
00118
                                   else:
00119
                                       logging.info(
                                           self.loginfo + "unidentified board detected"
00120
00121
00122
                                   logging.info(self.loginfo + "connected to " + p)
00123
                                   port = p
00124
                                   ser.reset_input_buffer()
00125
                                   ser.reset_output_buffer()
00126
00127
                      except Exception as e:
                          logging.error(self.logerr + "port identification: " + str(e))
00128
              if port == "":
00129
                  if self.ca.port:
00130
00131
                      logging.critical(
                           self.logcrit + "No usable board found at port " + str(self.ca.port)
00132
00133
00134
                      sys.exit(1)
00135
                  else:
00136
                      logging.critical(self.logcrit + "No usable board found")
00137
                      sys.exit(1)
00138
              self.port = port # COM port to use for RS422 link
00139
              self.ca.port = port[3:] # re-extract port number from com name
00140
00141
              self._ser = serial.Serial( # Class RS422
00142
                  port=self.port,
00143
                  baudrate=self.baud,
00144
                  parity=self.par,
00145
                  stopbits=self.stop,
00146
                  timeout=self.read_timeout, # timeout for serial read
00147
                  bytesize=serial.EIGHTBITS,
00148
00149
              self.payloadsize = (
00150
                  self.ca.sensor.width
00151
                  * self.ca.sensor.height
                  * self.ca.sensor.nframes
00152
00153
                  * self.ca.sensor.bytesperpixel
00154
00155
              logging.debug(
                  self.logdebug + "Payload size: " + str(self.payloadsize) + " bytes"
00156
00157
00158
              self._ser.flushInput()
00159
              if not self._ser.is_open:
                  logging.critical(self.logcrit + "Unable to open serial connection")
00160
00161
                  sys.exit(1)
00162
00163
          def serialClose(self):
00164
00165
              Close serial interface
00166
              logging.debug(self.logdebug + "serialclose")
00167
              self._ser.close() # close serial interface COM port
00168
00169
          def sendCMD(self, pkt):
00170
00171
```

7.20 RS422.pv 297

```
00172
              Submit packet and verify response packet. Recognizes readoff packet and adjusts
00173
              read size and timeout appropriately
00174
00175
00176
                  pkt: Packet object
00177
00178
              Returns:
              tuple (error, response string)
00179
00180
00181
              pktStr = pkt.pktStr()
00182
              logging.debug(self.logdebug + "sendCMD packet: " + str(pktStr))
00183
              self._ser.flushInput()
00184
              time.sleep(0.01) # wait 10 ms in between flushing input and output buffers
              self._ser.flushOutput()
00186
              self.ca.writeSerial(pktStr)
              err0 = ""
00187
              err = ""
00188
              resp = ""
00189
00190
              tries = 3 # TODO: make a function parameter?
00191
00192
00193
                  hasattr(self.ca, "board")
                  and pktStr[4] == "0"
00194
00195
                  and pktStr[5:8] == self.ca.board.registers["SRAM_CTL"]
00196
                   # download data payload
00197
00198
                  logging.info(
00199
                      self.loginfo + "Payload size (bytes) = " + str(self.payloadsize)
00200
                  crcresp0 = ""
00201
                  crcresp1 = ""
00202
                  smallresp = ""
00203
00204
                  emptyResponse = False
00205
                  wrongSize = False
00206
                  \ensuremath{\mathtt{\#}} TODO: refactor payload error management to another method
00207
                  for i in range(tries):
00208
                      err, resp = self.readSerial(
                           self.payloadsize + 20, timeout=self.datatimeout
00209
00210
00211
                      if err:
00212
                          logging.error(
                               self.logerr + "sendCMD: read payload failed " + pktStr + err
00213
00214
00215
                           self.ca.payloaderror = True
00216
                      else:
00217
                           if not len(resp):
00218
                               err0 = self.logerr + "sendCMD: empty response from board"
00219
                               logging.error(err0)
00220
                               emptyResponse = True
00221
                               self.ca.payloaderror = True
00222
                           elif len(resp) != 2 * (self.payloadsize + 20):
00223
                               err0 = (
00224
                                   self.logerr
00225
                                   + "sendCMD: incorrect response; expected "
00226
                                   + str(self.payloadsize + 20)
00227
                                   + " bytes, received "
                                   + str(len(resp) // 2)
00228
00229
00230
                               logging.error(err0)
00231
                               wrongSize = True
00232
                               smallresp = resp
00233
                               self.ca.payloaderror = True
00234
                           elif not checkCRC(resp[4:20]):
00235
                               err0 = (
00236
                                  self.logerr
00237
                                   + "sendCMD:
00238
                                   + pktStr
                                   + " - payload preface CRC fail"
00239
00240
00241
                               logging.error(err0)
00242
                               self.ca.payloaderror = True
00243
                               crcresp1 = resp
00244
                           elif not checkCRC(resp[24:]):
00245
                               err0 = (
                                   self.logerr + "sendCMD: " + pktStr + " - payload CRC fail"
00246
00247
00248
                               logging.error(err0)
00249
                               self.ca.payloaderror = True
00250
                               crcresp0 = resp
                           err += err0
00251
00252
                      time.sleep(5)
```

```
00253
                       if self.ca.payloaderror:
00254
                           # keep best results over multiple tries; e.g., if first try is
00255
                               bad CRC and second try is an incomplete payload, use the
00256
                               first payload
00257
                           if i == tries - 1:
00258
                                if crcresp0:
00259
                                    logging.error(
00260
                                       self.logerr + "sendCMD: Unable to acquire "
00261
                                        "CRC-confirmed payload after '
00262
                                        + str(tries)
00263
                                        + " attempts. Continuing with unconfirmed payload"
00264
00265
                                    resp = crcresp0
                                elif crcresp1:
00266
00267
                                    logging.error(
00268
                                       self.logerr + "sendCMD: Unable to acquire "
00269
                                         "CRC-confirmed readoff after '
00270
                                        + str(tries)
00271
                                        + " attempts. Continuing with unconfirmed payload"
00272
                                    )
00273
                                    resp = crcresp1
00274
                                elif wrongSize:
00275
                                    logging.error(
00276
                                        self.logerr + "sendCMD: Unable to acquire complete "
00277
                                         "pavload after '
00278
                                        + str(tries)
00279
                                        + " attempts. Dumping datastream to file."
00280
00281
                                    resp = smallresp
00282
                                    self.ca.dumpNumpy(resp)
00283
                                elif emptyResponse:
00284
                                    logging.error(
                                        self.logerr + "sendCMD: Unable to acquire any "
"payload after " + str(tries) + " attempts."
00285
00286
00287
00288
                           else:
00289
                                logging.info(
                                    self.loginfo + "Retrying download, attempt #" + str(i + 1)
00290
00291
00292
                                err = ""
                                err0 = ""
00293
00294
                                self.ca.payloaderror = False
00295
                                self.ca.writeSerial(pktStr)
00296
                       else:
                           logging.info(self.loginfo + "Download successful")
if self.ca.boardname == "llnl_v4":
00297
00298
00299
                                # self.ca.setSubregister('SWACK','1')
00300
00301
                           break
00302
00303
              else:
00304
                   # non-payload messages and workaround for initial setup before board object
00305
                   # has been initialized
00306
                   time.sleep(0.03)
00307
                   self.\_ser.timeout = 0.02
00308
                   err, resp = self.readSerial(10)
00309
                   logging.debug(self.logdebug + "sendCMD response: " + str(resp))
00310
                   if err:
00311
                       logging.error(
00312
                           self.logerr + "sendCMD: readSerial failed (regular packet) " + err
00313
                   elif not checkCRC(resp[4:20]):
00314
00315
                       err = self.logerr + "sendCMD- regular packet CRC fail: " + resp
00316
                       logging.error(err)
00317
              return err, resp
00318
00319
          def arm(self, mode):
00320
00321
              Puts camera into wait state for trigger. Mode determines source; defaults to
00322
                 'Hardware'
00323
00324
              Args:
00325
                           'Software' | 'S' activates software, disables hardware triggering
                  mode:
                           'Hardware' | 'H' activates hardware, disables software triggering
00326
00327
                             Hardware is the default
00328
00329
              Returns:
              tuple (error, response string)
00330
00331
              if not mode:
   mode = "Hardware"
00332
00333
```

7.20 RS422.pv 299

```
00334
               logging.info(self.loginfo + "arm")
00335
               logging.debug(self.logdebug + "arming mode: " + str(mode))
00336
               self.ca.clearStatus()
00337
               self.ca.latchPots()
00338
               err, resp = self.ca.startCapture(mode)
00339
00340
                   logging.error(self.logerr + "unable to arm camera")
00341
00342
                   self.ca.armed = True
00343
                   self.skipError = True
00344
               return err, resp
00345
00346
          def readFrames(self, waitOnSRAM, timeout=0, fast=False, columns=1):
00347
00348
               Copies image data from board into numpy arrays.
00349
00350
               Aras:
00351
                   waitOnSRAM: if True, wait until SRAM_READY flag is asserted to begin copying
00352
                     data
00353
                   timeout: passed to waitForSRAM; after this many seconds begin copying data
00354
                     irrespective of SRAM_READY status; 'zero' means wait indefinitely
00355
                     WARNING: If acquisition fails, the SRAM will not contain a current image,
00356
                   but the code will copy the data anyway fast: if False, parse and convert frames to numpy arrays; if True, return
00357
00358
                    unprocessed text stream
00359
                   columns: 1 for single image per frame, 2 for separate hemisphere images
00360
00361
               Returns:
00362
                   list of numpy arrays OR raw text stream
00363
00364
00365
               frames, _, _ = self.readoff(waitOnSRAM, timeout, fast, columns)
00366
               return frames
00367
00368
          def readoff(self, waitOnSRAM, timeout, fast, columns=1):
00369
               Copies image data from board into numpy arrays; returns data, length of data, and error messages. Use 'readFrames()' unless you require this additional
00370
00371
00372
               information
00373
00374
00375
                   waitOnSRAM: if True, wait until SRAM_READY flag is asserted to begin copying
00376
                     data
00377
                   timeout: passed to waitForSRAM; after this many seconds begin copying data
00378
                     irrespective of SRAM_READY status; ^\prime \, {\sf zero'} \, means wait indefinitely
00379
                     WARNING: If acquisition fails, the SRAM will not contain a current image,
00380
                       but the code will copy the data anyway
00381
                   fast: if False, parse and convert frames to numpy arrays; if True, return
00382
                     unprocessed text stream
00383
                   columns: 1 for single image per frame, 2 for separate hemisphere images
00384
00385
               Returns:
00386
                   tuple (list of numpy arrays OR raw text stream, length of downloaded payload
00387
                     in bytes, payload error flag)
00388
                   NOTE: This reduces readoff by <1 second, so will have no noticeable impact
00389
                     when using RS422
00390
00391
               logging.info(self.loginfo + "readoff")
00392
               logging.debug(
00393
                   self.logdebug
00394
                   + "readoff: waitonSRAM = "
00395
                   + str(waitOnSRAM)
00396
                   + "; timeout =
00397
                   + str(timeout)
00398
                   + "; fast =
00399
                   + str(fast)
00400
00401
               errortemp = False
00402
00403
               # Wait for data to be ready on board, turns off error messaging
               # Skip wait only if explicitly tagged 'False' ('None' defaults to True)
00404
00405
               if waitOnSRAM is not False:
00406
                   logging.getLogger().setLevel(logging.CRITICAL)
00407
                   self.ca.waitForSRAM(timeout)
00408
                   logging.getLogger().setLevel(self.ca.verblevel)
00409
00410
               # Retrieve data
00411
               err, rval = self.ca.readSRAM()
00412
               if err:
                   logging.error(self.logerr + "Error detected in readSRAM")
00413
00414
               time.sleep(0.3)
```

```
00415
              logging.debug(self.logdebug + "readoff: first 64 chars: " + str(rval[0:64]))
00416
              # extract only the read burst data. Remove header & CRC footer
00417
              read_burst_data = rval[36:-4]
00418
00419
              # Payload size as string implied by provided parameters
              expectedlength = (
00420
00421
00422
                  * (self.ca.sensor.lastframe - self.ca.sensor.firstframe + 1)
00423
                  * (self.ca.sensor.lastrow - self.ca.sensor.firstrow + 1)
00424
                  * self.ca.sensor.width
00425
00426
              padding = expectedlength - len(read_burst_data)
00427
              if padding:
00428
                  logging.warning(
00429
                       "{logwarn}readoff: Payload is shorter than expected."
00430
                      " Padding with '0's".format(logwarn=self.logwarn)
00431
00432
                  read_burst_data = read_burst_data.ljust(expectedlength, "0")
00433
00434
              if fast:
00435
                  return read_burst_data, len(read_burst_data) // 2, errortemp
00436
              else:
00437
                  parsed = generateFrames(self.ca, read_burst_data, columns)
00438
                  return parsed, len(read_burst_data) // 2, errortemp
00439
          def writeSerial(self, outstring, timeout):
00440
00441
00442
              Transmit string to board
00443
00444
              Args:
00445
                  outstring: string to write
00446
                  timeout: serial timeout in sec
00447
              Returns:
              integer length of string written to serial port
00448
00449
00450
              logging.debug(
00451
                  self.logdebug
00452
                  + "writeSerial: outstring = "
00453
                  + str(outstring)
00454
                  + "; timeout =
00455
                  + str(timeout)
00456
00457
              if timeout:
00458
                  self._ser.timeout = timeout
00459
00460
                  self._ser.timeout = self.write_timeout
00461
              lengthwritten = self._ser.write(str2bytes(outstring))
00462
              self._ser.timeout = self.read_timeout # reset if changed above
00463
              return lengthwritten
00464
00465
          def readSerial(self, size, timeout=None):
00466
00467
              Read bytes from the serial port. Does not verify packets.
00468
00469
00470
                 size: number of bytes to read
00471
                 timeout: serial timeout in sec
00472
00473
              tuple (error string, string read from serial port)
00474
00475
00476
              logging.debug(
00477
                 self.logdebug
00478
                  + "readSerial: size = "
00479
                  + str(size)
                  + "; timeout = "
00480
00481
                  + str(timeout)
00482
              err = ""
00483
00484
              if timeout:
00485
                  self._ser.timeout = timeout
00486
              else:
00487
                 self._ser.timeout = self.read_timeout
00488
              resp = self._ser.read(size)
00489
              if len(resp) < 10: # bytes
00490
                  err += (
00491
                     self.logerr + "readSerial : packet too small: '" + bytes2str(resp) + "'"
00492
00493
                  logging.error(err)
00494
              return err, bytes2str(resp)
00495
```

```
00496
         def closeDevice(self):
00497
00498
             Close primary serial interface
00499
00500
              logging.debug(self.logdebug + "Closing RS422 connection")
00501
00502
00503
00504 """
00505 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00508 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00509 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00510 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00511 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00512 be made under this license.
00513 """
```

# 7.21 C:/Users/hill35/git/camera\_python/nsCamera/comms/ZestETM1/Data.c File Reference

```
#include <stdint.h>
#include <stdlib.h>
#include <stdio.h>
#include <memory.h>
#include "ZestETM1.h"
#include "Private.h"
```

#### Classes

- struct ZESTETM1\_WRITE\_REG\_CMD
- struct ZESTETM1\_WRITE\_REG\_RESPONSE
- struct ZESTETM1 READ REG CMD
- struct ZESTETM1\_READ\_REG\_RESPONSE
- struct ZESTETM1\_MAILBOX\_INT\_CMD
- struct ZESTETM1\_MAILBOX\_INT\_RESPONSE

#### **Macros**

- #define CRT SECURE NO WARNINGS
- #define ZESTETM1\_COMMAND\_SPI 0xee
- #define ZESTETM1 COMMAND WRITE REG 0xf6
- #define ZESTETM1\_COMMAND\_READ\_REG 0xf7
- #define ZESTETM1\_COMMAND\_MAILBOX\_INT 0xf8

#### **Functions**

ZESTETM1\_STATUS ZestETM1\_OpenConnection (ZESTETM1\_CARD\_INFO \*CardInfo, ZESTETM1\_CONNECTION\_TYPE
 Type, uint16\_t Port, uint16\_t LocalPort, ZESTETM1\_CONNECTION \*Connection)

- ZESTETM1\_STATUS ZestETM1\_CloseConnection (ZESTETM1\_CONNECTION Connection)
- static ZESTETM1\_STATUS ZestETM1\_WriteData (ZESTETM1\_CONNECTION Connection, void \*Buffer, uint32 t Length, unsigned long \*Written, uint32 t Timeout)
- static ZESTETM1\_STATUS ZestETM1\_ReadData (ZESTETM1\_CONNECTION Connection, void \*Buffer, uint32 t Length, unsigned long \*Read, uint32 t Timeout)
- ZESTETM1\_STATUS ZestETM1\_SendCommand (ZESTETM1\_CARD\_INFO \*CardInfo, ZESTETM1\_CONNECTION Connection, void \*WriteBuffer, uint32 t WriteLen, void \*ReadBuffer, uint32 t ReadLen, int WaitForAck)
- ZESTETM1\_STATUS ZestETM1\_SPIReadWrite (ZESTETM1\_CARD\_INFO \*CardInfo, ZESTETM1\_CONNECTION Connection, int Device, int WordLen, uint32\_t \*WriteData, uint32\_t \*ReadData, uint32\_t Length, int ReleaseCS, int WaitForAck)
- ZESTETM1\_STATUS ZestETM1OpenConnection (ZESTETM1\_CARD\_INFO \*CardInfo, ZESTETM1\_CONNECTION\_TYPE Type, uint16\_t Port, uint16\_t LocalPort, ZESTETM1\_CONNECTION \*Connection)
- ZESTETM1 STATUS ZestETM1CloseConnection (ZESTETM1 CONNECTION Connection)
- ZESTETM1\_STATUS ZestETM1WriteData (ZESTETM1\_CONNECTION Connection, void \*Buffer, unsigned long Length, unsigned long \*Written, unsigned long Timeout)
- ZESTETM1\_STATUS ZestETM1ReadData (ZESTETM1\_CONNECTION Connection, void \*Buffer, unsigned long Length, unsigned long \*Read, unsigned long Timeout)
- ZESTETM1\_STATUS ZestETM1SPIReadWrite (ZESTETM1\_CARD\_INFO \*CardInfo, ZESTETM1\_SPI\_RATE Rate, int WordLen, void \*WriteData, void \*ReadData, unsigned long Length, int ReleaseCS)
- ZESTETM1\_STATUS ZestETM1WriteRegister (ZESTETM1\_CARD\_INFO \*CardInfo, unsigned long Addr, unsigned short Data)
- ZESTETM1\_STATUS ZestETM1ReadRegister (ZESTETM1\_CARD\_INFO \*CardInfo, unsigned long Addr, unsigned short \*Data)
- ZESTETM1 STATUS ZestETM1SetInterrupt (ZESTETM1 CARD INFO \*CardInfo)

## 7.21.1 Class Documentation

## 7.21.1.1 struct ZESTETM1\_WRITE\_REG\_CMD

Definition at line 58 of file Data.c.

#### **Class Members**

uint8_t	Addr	
uint8_t	Command	
uint16_t	Data	

## 7.21.1.2 struct ZESTETM1\_WRITE\_REG\_RESPONSE

Definition at line 64 of file Data.c.

**Class Members** 

uint8_t	Command	
---------	---------	--

## **Class Members**

uint8_t	Dummy1[2]	
uint8_t	Status	

# 7.21.1.3 struct ZESTETM1\_READ\_REG\_CMD

Definition at line 70 of file Data.c.

## **Class Members**

uint8_t	Addr	
uint8_t	Command	
uint8_t	Dummy[2]	

# 7.21.1.4 struct ZESTETM1\_READ\_REG\_RESPONSE

Definition at line 76 of file Data.c.

## **Class Members**

uint8_t	Command	
uint8_t	Status	
uint16_t	Value	

# 7.21.1.5 struct ZESTETM1\_MAILBOX\_INT\_CMD

Definition at line 82 of file Data.c.

## **Class Members**

uint8_t	Command	
uint8_t	Dummy[3]	

# 7.21.1.6 struct ZESTETM1\_MAILBOX\_INT\_RESPONSE

Definition at line 87 of file Data.c.

# **Class Members**

#### **Class Members**

uint16_t	Dummy	
uint8_t	Status	

# 7.21.2 Macro Definition Documentation

## 7.21.2.1 \_CRT\_SECURE\_NO\_WARNINGS

#define \_CRT\_SECURE\_NO\_WARNINGS

Definition at line 45 of file Data.c.

# 7.21.2.2 ZESTETM1\_COMMAND\_MAILBOX\_INT

#define ZESTETM1\_COMMAND\_MAILBOX\_INT 0xf8

Definition at line 96 of file Data.c.

# 7.21.2.3 ZESTETM1\_COMMAND\_READ\_REG

#define ZESTETM1\_COMMAND\_READ\_REG 0xf7

Definition at line 95 of file Data.c.

## 7.21.2.4 ZESTETM1\_COMMAND\_SPI

#define ZESTETM1\_COMMAND\_SPI 0xee

Definition at line 93 of file Data.c.

# 7.21.2.5 ZESTETM1\_COMMAND\_WRITE\_REG

#define ZESTETM1\_COMMAND\_WRITE\_REG 0xf6

Definition at line 94 of file Data.c.

## 7.21.3 Function Documentation

#### 7.21.3.1 ZestETM1\_CloseConnection()

```
ZESTETM1_STATUS ZestETM1_CloseConnection (
              ZESTETM1_CONNECTION Connection )
Definition at line 225 of file Data.c.
00226 {
00227
          ZESTETM1_CONNECTION_STRUCT *Conn = (ZESTETM1_CONNECTION_STRUCT *)Connection;
00228
00229
00230
00231
              return ZESTETM1_NULL_PARAMETER;
00232
00233
         if (Conn->Magic!=ZESTETM1_CONNECTION_HANDLE_MAGIC)
00234
00235
              return ZESTETM1_ILLEGAL_CONNECTION;
00236
00237
00238
         // Cleanup
00239
         closesocket(Conn->Socket);
00240
          Conn->Magic = 0;
00241
         free (Conn);
00242
00243
         return ZESTETM1_SUCCESS;
00244 }
```

## 7.21.3.2 ZestETM1\_OpenConnection()

```
ZESTETM1_STATUS ZestETM1_OpenConnection (
              ZESTETM1_CARD_INFO * CardInfo,
              ZESTETM1_CONNECTION_TYPE Type,
              uint16_t Port,
              uint16_t LocalPort,
              ZESTETM1_CONNECTION * Connection )
Definition at line 101 of file Data.c.
00106 {
00107
          ZESTETM1_CONNECTION_STRUCT *NewStruct;
00108
          SOCKET Socket = -1;
00109
         char AddrBuffer[32];
00110
         char PortBuffer[32];
00111
00112
         if (Connection==NULL || CardInfo==NULL)
00113
         {
00114
              return ZESTETM1_NULL_PARAMETER;
00115
00116
00117
          // Allocate data structure
00118
         NewStruct = malloc(sizeof(ZESTETM1_CONNECTION_STRUCT));
00119
         if (NewStruct==NULL)
00120
00121
              return ZESTETM1_OUT_OF_MEMORY;
00122
00123
00124
         // Build target addresses
         sprintf(AddrBuffer, "%d.%d.%d.%d", CardInfo->IPAddr[0], CardInfo->IPAddr[1],
00125
00126
                         CardInfo->IPAddr[2], CardInfo->IPAddr[3]);
00127
         sprintf(PortBuffer, "%d", Port);
00128
00129
          if (Type==ZESTETM1_TYPE_UDP)
00130
00131
              // Open UDP connection
             struct sockaddr_in SourceIP;
00132
00133
              int SourceLen = (int)sizeof(struct sockaddr in);
             Socket = socket(AF_INET, SOCK_DGRAM, IPPROTO_UDP);
00134
00135
             if (Socket<0)</pre>
```

```
00136
                   return ZESTETM1_SOCKET_ERROR;
00137
              SourceIP.sin_family = AF_INET;
00138
               SourceIP.sin_addr.s_addr = htonl(INADDR_ANY);
               SourceIP.sin_port = htons(LocalPort);
00139
00140
              bind(Socket, (const struct sockaddr *)&SourceIP, SourceLen);
00141
00142
              NewStruct->Target.sin_family = AF_INET;
00143
              NewStruct->Target.sin_addr.s_addr = inet_addr(AddrBuffer);
00144
              NewStruct->Target.sin_port = htons(atoi(PortBuffer));
00145
00146
          else if (Type==ZESTETM1_TYPE_TCP)
00147
00148
               // Open TCP connection
              struct addrinfo *AddrResult = NULL,
00149
00150
                               *Ptr = NULL,
00151
                               Hints;
00152
              int Result;
00153
              struct sockaddr_in SourceIP;
00154
              int SourceLen = (int)sizeof(struct sockaddr_in);
00155
              memset(&Hints, 0, sizeof(Hints));
Hints.ai_family = AF_UNSPEC;
00156
00157
00158
              Hints.ai_socktype = SOCK_STREAM;
00159
              Hints.ai_protocol = IPPROTO_TCP;
00160
00161
               // Resolve the server address and port
              Result = getaddrinfo(AddrBuffer, PortBuffer, &Hints, &AddrResult);
00162
00163
              if (Result!=0)
00164
00165
                   free (NewStruct);
                   return ZESTETM1_SOCKET_ERROR;
00166
00167
00168
00169
              \ensuremath{//} Attempt to connect to an address until one succeeds
              for (Ptr=AddrResult; Ptr!=NULL; Ptr=Ptr->ai_next)
00170
00171
00172
                   // Create a SOCKET for connecting to server
00173
                   Socket = socket(Ptr->ai_family, Ptr->ai_socktype,
00174
                                   Ptr->ai_protocol);
00175
                   if (Socket<0)
00176
00177
                       freeaddrinfo(AddrResult);
00178
                       free (NewStruct);
                       return ZESTETM1_SOCKET_ERROR;
00179
00180
00181
00182
                   // Connect to ZestETM1
00183
                   Result = connect(Socket, Ptr->ai_addr, (int)Ptr->ai_addrlen);
00184
                   if (Result<0)</pre>
00185
00186
                       closesocket(Socket);
00187
                       Socket = -1;
00188
                       continue;
00189
00190
                   break;
00191
00192
00193
               SourceIP.sin_family = AF_INET;
00194
               SourceIP.sin_addr.s_addr = htonl(INADDR_ANY);
00195
               SourceIP.sin_port = 0;
00196
              bind(Socket, (const struct sockaddr *)&SourceIP, SourceLen);
00197
               freeaddrinfo(AddrResult);
00198
00199
          else
00200
00201
               free (NewStruct);
00202
              return ZESTETM1_INVALID_CONNECTION_TYPE;
00203
00204
00205
          if (Socket==-1)
00206
          {
00207
              free (NewStruct);
              return ZESTETM1_SOCKET_ERROR;
00208
00209
00210
00211
          NewStruct->Magic = ZESTETM1_CONNECTION_HANDLE_MAGIC;
00212
          NewStruct->Type = Type;
          NewStruct->Port = Port;
00213
00214
          NewStruct->LocalPort = LocalPort;
          NewStruct->Socket = Socket;
00215
00216
          NewStruct->CardInfo = CardInfo;
```

## 7.21.3.3 ZestETM1\_ReadData()

```
00366
00367
00368
          int Result;
00369
          struct sockaddr_in Target;
00370
          int TargetLen = (int)sizeof(struct sockaddr_in);
00371
00372
          if (Conn==NULL)
00373
00374
              if (Read!=NULL) *Read = 0;
00375
              return ZESTETM1_NULL_PARAMETER;
00376
00377
          if (Conn->Magic!=ZESTETM1 CONNECTION HANDLE MAGIC)
00378
00379
              if (Read!=NULL) *Read = 0;
00380
              return ZESTETM1_ILLEGAL_CONNECTION;
00381
00382
00383
00384
00385
              int OK = 1;
00386
00387
              // Wait for socket to become ready
00388
00389
                  int CurTime = 0;
00390
                  while (CurTime!=Timeout)
00391
00392
                      int T = (Timeout-CurTime) < 1000 ? Timeout-CurTime : 1000;</pre>
00393
                      FD_ZERO(&ReadFDS);
00394
                      FD_SET(Conn->Socket, &ReadFDS);
00395
                      Time.tv_sec = T/1000;
00396
                      Time.tv_usec = (T%1000) *1000;
00397
                      Result = select((int)Conn->Socket+1, &ReadFDS, NULL, NULL, &Time);
00398
                       if (Result>0) break;
00399
                      CurTime+=T;
00400
00401
00402
              if (Result<0 || !FD_ISSET(Conn->Socket, &ReadFDS))
             {
00404
                  if (Read!=NULL) *Read = BufferPtr;
00405
                  return ZESTETM1_TIMEOUT;//FIXME: Any other errors?
00406
             }
00407
00408
              if (Conn->Type==ZESTETM1_TYPE_UDP)
00409
00410
                  Result = recvfrom(Conn->Socket, (char *)Buffer+BufferPtr,
                                     Length-BufferPtr, 0,
00411
00412
                                     (struct sockaddr *) & Target, & TargetLen);
00413
                  if (Conn->LocalPort!=0 && Target.sin_port!=htons(Conn->LocalPort))
00414
                      OK = 0;
00415
              else if (Conn->Type==ZESTETM1_TYPE_TCP)
00416
00417
00418
                  Result = recv(Conn->Socket, (char *)Buffer+BufferPtr,
00419
                                Length-BufferPtr, 0);
00420
```

```
00421
              else
00422
              {
00423
                  if (Read!=NULL) *Read = BufferPtr;
00424
                  return ZESTETM1_ILLEGAL_CONNECTION;
00425
00426
00427
              // Update counters
00428
              if (Result>0 && OK==1)
00429
              {
00430
                  BufferPtr += Result;
00431
00432
              else if (Result==0)
00433
                  // Connection closed
                  if (Read!=NULL) *Read = BufferPtr;
00435
00436
                  return ZESTETM1_SOCKET_CLOSED;
00437
00438
              else if (Result<0)
00439
00440
                  // Socket error
00441
                  if (Read!=NULL) *Read = BufferPtr;
00442
                  return ZESTETM1_SOCKET_ERROR;
00443
00444
         } while (Result>0 && BufferPtr<Length);
00445
          if (Read!=NULL) *Read = BufferPtr;
00446
00447
          return ZESTETM1_SUCCESS;
00448 }
```

## 7.21.3.4 ZestETM1\_SendCommand()

ZESTETM1\_STATUS ZestETM1\_SendCommand (

```
ZESTETM1_CARD_INFO * CardInfo,
               ZESTETM1_CONNECTION Connection,
               void * WriteBuffer,
               uint32_t WriteLen,
               void * ReadBuffer,
               uint32_t ReadLen,
               int WaitForAck )
Definition at line 453 of file Data.c.
00458 {
00459
          ZESTETM1_STATUS Result;
00460
          unsigned long Written;
00461
          unsigned long Received;
00462
00463
          // Send/receive data
00464
          Result = ZestETM1_WriteData(Connection, WriteBuffer, WriteLen, &Written,
00465
                                     CardInfo->Timeout);
00466
          if (Result!=ZESTETM1_SUCCESS)
00467
00468
              return Result;
00469
00470
          if (Written!=WriteLen)
00471
00472
              return ZESTETM1_INTERNAL_ERROR;
00473
00474
          *((uint8_t *)ReadBuffer) = 0;
00475
          if (WaitForAck==1)
00476
00477
              Result = ZestETM1_ReadData(Connection, ReadBuffer, ReadLen, &Received,
00478
                                        CardInfo->Timeout);
00479
              if (Result!=ZESTETM1_SUCCESS)
00480
00481
                  return Result:
00482
00483
              if (Received!=ReadLen)
00484
00485
                  return ZESTETM1 INTERNAL ERROR;
00486
00487
00488
00489
          return ZESTETM1_SUCCESS;
00490 }
```

## 7.21.3.5 ZestETM1\_SPIReadWrite()

```
ZESTETM1_STATUS ZestETM1_SPIReadWrite (
               ZESTETM1_CARD_INFO * CardInfo,
               ZESTETM1_CONNECTION Connection,
               int Device,
              int WordLen.
               uint32_t * WriteData,
              uint32_t * ReadData,
              uint32_t Length,
               int ReleaseCS,
               int WaitForAck )
Definition at line 495 of file Data.c.
00501 {
00502
          uint8_t Buffer[65536];
00503
          uint32_t *BufPtr;
00504
          uint32_t i;
00505
          ZESTETM1_STATUS Result;
00507
          // Build command
00508
         Buffer[0] = ZESTETM1_COMMAND_SPI;// Command byte
                                 // SPI device
00509
          Buffer[1] = Device;
00510
          Buffer[2] = WordLen;
                                          // Word length
          Buffer[3] = ReleaseCS;
00511
                                          // Release CS on completion
00512
          BufPtr = (uint32_t *) (Buffer+4);
00513
          BufPtr[0] = WriteData==NULL ? 0 : ZESTETM1_REVERSE(Length);
00514
                                                                          // Num words
00515
          BufPtr[1] = ReadData==NULL ? 0 : ZESTETM1_REVERSE(Length);
00516
          if (WriteData!=NULL)
00517
00518
              for (i=0; i<Length; i++)</pre>
00519
                  BufPtr[2+i] = ZESTETM1_REVERSE(WriteData[i]);
00520
00521
00522
00523
         else
00524
              memset (BufPtr+2, 0, Length *4);
00525
00526
          // Send command
00527
         Result = ZestETM1_SendCommand(CardInfo, Connection, Buffer,
00528
                                       WriteData == NULL ? 12 : 12 + Length * 4,
                                       Buffer, ReadData==NULL ? 4 : 4+Length*4,
00529
00530
                                       WaitForAck);
         if (WaitForAck==0)
00531
00532
              return Result;
00533
          // Extract results
00534
         if (Buffer[0]!=ZESTETM1_COMMAND_SPI || Buffer[1]!=0)
00535
00536
00537
              return ZESTETM1_INTERNAL_ERROR;
00538
00539
00540
         BufPtr = (uint32_t *)(Buffer+4);
00541
          if (ReadData!=NULL)
00542
00543
              for (i=0; i<Length; i++)</pre>
00544
00545
                  ReadData[i] = ZESTETM1_REVERSE(BufPtr[i]);
00546
00547
00548
00549
          return ZESTETM1_SUCCESS;
00550 }
```

# 7.21.3.6 ZestETM1\_WriteData()

```
void * Buffer,
uint32_t Length,
unsigned long * Written,
uint32_t Timeout ) [static]
```

#### Definition at line 249 of file Data.c.

```
00255
          ZESTETM1_CONNECTION_STRUCT *Conn = (ZESTETM1_CONNECTION_STRUCT *)Connection;
00256
          fd_set WriteFDS;
00257
          struct timeval Time;
00258
          uint32_t BufferPtr = 0;
00259
          int Result;
          int TargetLen = (int)sizeof(struct sockaddr_in);
00260
00261
          unsigned int MaxSize;
00262
00263
          if (Conn==NULL)
00264
          {
00265
              if (Written!=NULL) *Written = 0;
00266
              return ZESTETM1_NULL_PARAMETER;
00267
00268
          if (Conn->Magic!=ZESTETM1_CONNECTION_HANDLE_MAGIC)
00269
         {
00270
              if (Written!=NULL) *Written = 0;
00271
              return ZESTETM1_ILLEGAL_CONNECTION;
00272
00273
          if (Conn->Type==ZESTETM1_TYPE_UDP)
00274
00275
          {
00276
              // Get maximum transfer size
00277 #if defined(MSVC) || defined(WINGCC)
00278
              int MaxSizeLen = sizeof(MaxSize);
00279
              getsockopt(Conn->Socket, SOL_SOCKET, SO_MAX_MSG_SIZE,
00280
                         (char *)&MaxSize, &MaxSizeLen);
00281 #else
              //FIXME: Linux doesn't support SO_MAX_MSG_SIZE
00282
00283
              MaxSize = 65507;
00284 #endif
00285
         }
00286
00287
          do
00288
          {
00289
              uint32_t Bytes = Length-BufferPtr;
00290
00291
              if (Conn->Type==ZESTETM1_TYPE_UDP && Bytes>MaxSize)
00292
                  Bytes = MaxSize;
00293
00294
              // Wait for socket to become ready
00295
00296
                  int CurTime = 0;
00297
                  while (CurTime!=Timeout)
00298
00299
                       int T = (Timeout-CurTime) < 1000 ? Timeout-CurTime : 1000;
00300
                      FD_ZERO(&WriteFDS);
00301
                      FD_SET(Conn->Socket, &WriteFDS);
00302
                      Time.tv\_sec = T/1000;
00303
                      Time.tv\_usec = (T%1000)*1000;
00304
                      Result = select((int)Conn->Socket+1, NULL, &WriteFDS, NULL, &Time);
00305
                       if (Result>0) break;
00306
                      CurTime+=T;
00307
00308
00309
              if (Result<0 || !FD_ISSET(Conn->Socket, &WriteFDS))
00310
              {
00311
                  if (Written!=NULL) *Written = BufferPtr;
00312
                  return ZESTETM1_TIMEOUT;//FIXME: Any other errors?
00313
00314
00315
              if (Conn->Type==ZESTETM1_TYPE_UDP)
00316
00317
                  Result = sendto(Conn->Socket, (char *)Buffer+BufferPtr,
00318
                                   Bytes, 0,
00319
                                   (struct sockaddr *) &Conn->Target, TargetLen);
00320
              else if (Conn->Type==ZESTETM1_TYPE_TCP)
00321
00322
00323
                  Result = send(Conn->Socket, (char *)Buffer+BufferPtr,
00324
                                Bytes, 0);
00325
00326
              else
```

```
00327
              {
00328
                  if (Written!=NULL) *Written = BufferPtr;
00329
                  return ZESTETM1_ILLEGAL_CONNECTION;
00330
00331
00332
             // Update counters
00333
              if (Result>0)
00334
00335
                  BufferPtr += Result;
00336
00337
             else if (Result==0)
00338
00339
                  // Connection closed
00340
                  if (Written!=NULL) *Written = BufferPtr;
00341
                 return ZESTETM1_SOCKET_CLOSED;
00342
00343
             else
00344
             {
00345
                  // Socket error
00346
                  if (Written!=NULL) *Written = BufferPtr;
00347
                 return ZESTETM1_SOCKET_ERROR;
00348
00349
        } while (Result>0 && BufferPtr<Length);
00350
00351
         if (Written!=NULL) *Written = BufferPtr;
00352
          return ZESTETM1 SUCCESS;
00353 }
```

## 7.21.3.7 ZestETM1CloseConnection()

#### Definition at line 575 of file Data.c.

```
00576 {
00577
          ZESTETM1_STATUS Result;
          ZESTETM1_CONNECTION_STRUCT *Conn = (ZESTETM1_CONNECTION_STRUCT*)Connection;
00578
00579
00580
          Result = ZestETM1_CloseConnection(Connection);
00581
          if (Result!=ZESTETM1_SUCCESS)
00582
00583
              ZESTETM1_ERROR_CONN("ZestETM1CloseConnection", Result);
00584
00585
00586
          return ZESTETM1_SUCCESS;
00587 }
```

## 7.21.3.8 ZestETM1OpenConnection()

```
ZESTETM1_STATUS ZestETM1OpenConnection (
              ZESTETM1_CARD_INFO * CardInfo,
              ZESTETM1_CONNECTION_TYPE Type,
              uint16_t Port,
              uint16_t LocalPort,
              ZESTETM1_CONNECTION * Connection )
Definition at line 555 of file Data.c.
00560 {
00561
         ZESTETM1_STATUS Result;
00562
00563
         Result = ZestETM1_OpenConnection(CardInfo, Type, Port, LocalPort, Connection);
00564
         if (Result!=ZESTETM1_SUCCESS)
00565
             ZESTETM1_ERROR("ZestETM1OpenConnection", Result);
00566
00567
00568
00569
         return ZESTETM1_SUCCESS;
```

00570 }

## 7.21.3.9 ZestETM1ReadData()

```
ZESTETM1 STATUS ZestETM1ReadData (
               ZESTETM1_CONNECTION Connection,
               void * Buffer,
               unsigned long Length,
               unsigned long * Read,
               unsigned long Timeout )
Definition at line 613 of file Data.c.
00618 {
00619
          ZESTETM1_STATUS Result;
00620
          ZESTETM1_CONNECTION_STRUCT *Conn = (ZESTETM1_CONNECTION_STRUCT*)Connection;
00621
00622
          Result = ZestETM1_ReadData(Connection, Buffer, Length, Read, Timeout);
          if (Result!=ZESTETM1_SUCCESS)
00623
00624
00625
              ZESTETM1_ERROR_CONN("ZestETM1ReadData", Result);
00626
00627
00628
          return ZESTETM1_SUCCESS;
00629 }
7.21.3.10 ZestETM1ReadRegister()
ZESTETM1_STATUS ZestETM1ReadRegister (
               ZESTETM1_CARD_INFO * CardInfo,
               unsigned long Addr,
               unsigned short * Data )
Definition at line 728 of file Data.c.
00729 {
00730
          ZESTETM1 STATUS Result:
00731
          ZESTETM1_CONNECTION Connection;
00732
          ZESTETM1_READ_REG_CMD Cmd;
00733
          ZESTETM1_READ_REG_RESPONSE Response;
00734
00735
          if (CardInfo==NULL || Data==NULL)
00736
00737
              ZESTETM1_ERROR("ZestETM1ReadRegister", ZESTETM1_NULL_PARAMETER);
00738
00739
00740
00741
              ZESTETM1_ERROR("ZestETM1ReadRegister", ZESTETM1_ILLEGAL_PARAMETER);
00742
00743
00744
          Result = ZestETM1_OpenConnection(CardInfo, ZESTETM1_TYPE_TCP, CardInfo->ControlPort, 0, &Connection);
00745
          if (Result!=ZESTETM1_SUCCESS)
00746
00747
              ZESTETM1_ERROR("ZestETM1ReadRegister", Result);
00748
00749
00750
          // Read register from the device
00751
          Cmd.Command = ZESTETM1_COMMAND_READ_REG;
00752
          Cmd.Addr = (uint8_t)Addr;
00753
          Result = ZestETM1_SendCommand(CardInfo, Connection,
00754
                                       &Cmd, sizeof(Cmd),
00755
                                       &Response, sizeof(Response), 1);
00756
          if (Result!=ZESTETM1_SUCCESS)
00757
          {
00758
              ZestETM1_CloseConnection(Connection);
00759
              ZESTETM1_ERROR("ZestETM1ReadRegister", Result);
00760
00761
          if (Response.Command!=Cmd.Command || Response.Status!=0)
00762
00763
              ZestETM1 CloseConnection(Connection);
00764
              ZESTETM1_ERROR("ZestETM1ReadRegister", ZESTETM1_INTERNAL_ERROR);
00765
00766
          *Data = ((Response.Value»8)&0xff) | ((Response.Value&0xff) <8);
```

```
00767
00768    Result = ZestETM1_CloseConnection(Connection);
00769    if (Result!=ZESTETM1_SUCCESS)
00770    {
        ZESTETM1_ERROR("ZestETM1ReadRegister", Result);
00772    }
00773    return ZESTETM1_SUCCESS;
00774    return ZESTETM1_SUCCESS;
```

#### 7.21.3.11 ZestETM1SetInterrupt()

```
ZESTETM1_STATUS ZestETM1SetInterrupt (
               ZESTETM1_CARD_INFO * CardInfo )
Definition at line 780 of file Data.c.
00781 {
00782
          ZESTETM1_STATUS Result;
00783
          ZESTETM1_CONNECTION Connection;
00784
          ZESTETM1_MAILBOX_INT_CMD Cmd;
00785
          ZESTETM1_MAILBOX_INT_RESPONSE Response;
00786
00787
          if (CardInfo==NULL)
00788
00789
              ZESTETM1 ERROR ("ZestETM1SetInterrupt", ZESTETM1 NULL PARAMETER);
00790
00791
00792
          Result = ZestETM1_OpenConnection(CardInfo, ZESTETM1_TYPE_TCP, CardInfo->ControlPort, 0, &Connection);
00793
          if (Result!=ZESTETM1_SUCCESS)
00794
          {
00795
              ZESTETM1_ERROR("ZestETM1SetInterrupt", Result);
00796
00797
00798
          \ensuremath{//} Send command to set interrupt
00799
          Cmd.Command = ZESTETM1 COMMAND MAILBOX INT;
00800
          Result = ZestETM1_SendCommand(CardInfo, Connection,
00801
                                        &Cmd, sizeof(Cmd),
00802
                                        &Response, sizeof(Response), 1);
00803
          if (Result!=ZESTETM1_SUCCESS)
00804
          {
              ZestETM1_CloseConnection(Connection);
00805
00806
              ZESTETM1_ERROR("ZestETM1SetInterrupt", Result);
00807
00808
          if (Response.Command!=Cmd.Command || Response.Status!=0)
00809
00810
              ZestETM1_CloseConnection(Connection);
00811
              ZESTETM1_ERROR("ZestETM1SetInterrupt", ZESTETM1_INTERNAL_ERROR);
00812
00813
00814
          Result = ZestETM1_CloseConnection(Connection);
00815
          if (Result!=ZESTETM1_SUCCESS)
00816
          {
00817
              ZESTETM1_ERROR("ZestETM1SetInterrupt", Result);
00818
00819
00820
          return ZESTETM1_SUCCESS;
00821 }
```

#### 7.21.3.12 ZestETM1SPIReadWrite()

```
unsigned long Length,
int ReleaseCS )
```

```
Definition at line 635 of file Data.c.
```

```
00639 {
          ZESTETM1_STATUS Result;
00640
00641
          ZESTETM1_CONNECTION Connection;
          unsigned long RateVal = Rate==ZESTETM1_SPI_RATE_35 ? ZESTETM1_RATE_40MHz :
00642
00643
                                  Rate==ZESTETM1_SPI_RATE_17_5 ? ZESTETM1_RATE_20MHz : ZESTETM1_RATE_10MHz;
00644
00645
          if (CardInfo==NULL || (WriteData==NULL && ReadData==NULL))
00646
00647
              ZESTETM1_ERROR("ZestETM1SPIReadWrite", ZESTETM1_NULL_PARAMETER);
00648
00649
          if (WordLen<1 || WordLen>32 || Length>16384)
00650
00651
              ZESTETM1_ERROR("ZestETM1SPIReadWrite", ZESTETM1_ILLEGAL_PARAMETER);
00652
00653
00654
          Result = ZestETM1_OpenConnection(CardInfo, ZESTETM1_TYPE_TCP, CardInfo->ControlPort, 0, &Connection);
00655
          if (Result!=ZESTETM1_SUCCESS)
00656
00657
              ZESTETM1_ERROR("ZestETM1SPIReadWrite", Result);
00658
00659
00660
          Result = ZestETM1_SPIReadWrite(CardInfo, Connection, ZESTETM1_USER_DEVICE_ID|RateVal, WordLen,
      (uint32_t *)WriteData, (uint32_t *)ReadData, Length, ReleaseCS, 1);
00661
          if (Result!=ZESTETM1_SUCCESS)
00662
          {
00663
              ZESTETM1_ERROR("ZestETM1SPIReadWrite", Result);
00664
00665
00666
          Result = ZestETM1_CloseConnection(Connection);
00667
          if (Result!=ZESTETM1_SUCCESS)
00668
00669
              ZESTETM1_ERROR("ZestETM1SPIReadWrite", Result);
00670
00671
00672
          return ZESTETM1_SUCCESS;
00673 }
```

## 7.21.3.13 ZestETM1WriteData()

```
ZESTETM1_STATUS ZestETM1WriteData (
    ZESTETM1_CONNECTION Connection,
    void * Buffer,
    unsigned long Length,
    unsigned long * Written,
    unsigned long Timeout )
```

#### Definition at line 592 of file Data.c.

```
00597 {
00598
          ZESTETM1_STATUS Result;
00599
          ZESTETM1_CONNECTION_STRUCT *Conn = (ZESTETM1_CONNECTION_STRUCT*)Connection;
00600
00601
          Result = ZestETM1_WriteData(Connection, Buffer, Length, Written, Timeout);
00602
          if (Result!=ZESTETM1_SUCCESS)
00603
          {
00604
              ZESTETM1_ERROR_CONN("ZestETM1WriteData", Result);
00605
00606
00607
          return ZESTETM1_SUCCESS;
00608 }
```

#### 7.21.3.14 ZestETM1WriteRegister()

7.22 Data.c 315

```
unsigned long Addr,
unsigned short Data )
```

## Definition at line 678 of file Data.c.

```
00679 {
00680
          ZESTETM1_STATUS Result;
00681
          ZESTETM1_CONNECTION Connection;
00682
          ZESTETM1_WRITE_REG_CMD Cmd;
00683
          ZESTETM1_WRITE_REG_RESPONSE Response;
00684
00685
          if (CardInfo==NULL)
00686
00687
              ZESTETM1_ERROR("ZestETM1WriteRegister", ZESTETM1_NULL_PARAMETER);
00688
00689
00690
          {
00691
              ZESTETM1_ERROR("ZestETM1WriteRegister", ZESTETM1_ILLEGAL_PARAMETER);
00692
00693
00694
          Result = ZestETM1_OpenConnection(CardInfo, ZESTETM1_TYPE_TCP, CardInfo->ControlPort, 0, &Connection);
00695
          if (Result!=ZESTETM1_SUCCESS)
00696
00697
              ZESTETM1_ERROR("ZestETM1WriteRegister", Result);
00698
00699
00700
          // Write register to the device
00701
          //FIXME: Do we want to be able to copy this value to flash?
00702
          Cmd.Command = ZESTETM1_COMMAND_WRITE_REG;
00703
          Cmd.Addr = (uint8 t)Addr;
00704
          Cmd.Data = ((Data>8)&0xff) | (Data&0xff);
00705
          Result = ZestETM1_SendCommand(CardInfo, Connection,
00706
                                        &Cmd, sizeof(Cmd),
00707
                                        &Response, sizeof(Response), 1);
00708
          if (Result!=ZESTETM1 SUCCESS)
00709
00710
              ZestETM1 CloseConnection(Connection);
00711
              ZESTETM1_ERROR("ZestETM1WriteRegister", Result);
00712
00713
          if (Response.Command!=Cmd.Command || Response.Status!=0)
00714
          {
00715
              ZestETM1 CloseConnection(Connection);
00716
              ZESTETM1_ERROR("ZestETM1WriteRegister", ZESTETM1_INTERNAL_ERROR);
00717
00718
          Result = ZestETM1_CloseConnection(Connection);
00719
          if (Result!=ZESTETM1_SUCCESS)
00720
00721
          {
00722
              ZESTETM1_ERROR("ZestETM1WriteRegister", Result);
00723
00724
00725
          return ZESTETM1_SUCCESS;
00726 }
```

## 7.22 Data.c

## Go to the documentation of this file.

```
00001 // File:
                    Data.c
00002 //
00003 // Purpose:
00004 //
           ZestETM1 Host Library
00005 //
            Data transfer functions
00006 //
00007 // Version: 1.00
00008 // Date: 11/10/12
00009
00010 // Copyright (C) 2012 Orange Tree Technologies Ltd. All rights reserved.
00011 // Orange Tree Technologies grants the purchaser of a ZestETM1 the right to use and
00012 // modify this source code in any form in designs that target the ZestETM1.
00013 // Orange Tree Technologies prohibits the use of this source code or any modification of
00014 // it in any form in designs that target any other hardware unless the purchaser of the
00015 // ZestETM1 has purchased the appropriate licence from Orange Tree Technologies.
00016 // Contact Orange Tree Technologies if you want to purchase such a licence.
00017
```

```
00020 //**
           Disclaimer: LIMITED WARRANTY AND DISCLAIMER. These designs are
                     provided to you "as is". Orange Tree Technologies and its licensors
00021 //**
00022 //**
                      make and you receive no warranties or conditions, express, implied,
00023 //**
                      statutory or otherwise, and Orange Tree Technologies specifically
00024 //**
                      disclaims any implied warranties of merchantability, non-infringement,
00025 //**
                      or fitness for a particular purpose. Orange Tree Technologies does not
00026 //**
                      warrant that the functions contained in these designs will meet your
                     requirements, or that the operation of these designs will be
                      uninterrupted or error free, or that defects in the Designs will be
00029 //**
                      corrected. Furthermore, Orange Tree Technologies does not warrant or
00030 //**
                      make any representations regarding use or the results of the use of the
00031 //**
                     designs in terms of correctness, accuracy, reliability, or otherwise.
00032 //**
00033 //**
                     LIMITATION OF LIABILITY. In no event will Orange Tree Technologies
00034 //**
                      or its licensors be liable for any loss of data, lost profits, cost or
00035 //**
                      procurement of substitute goods or services, or for any special,
00036 //**
                      incidental, consequential, or indirect damages arising from the use or
00037 //**
                      operation of the designs or accompanying documentation, however caused
00038 //**
                      and on any theory of liability. This limitation will apply even if
00039 //**
                      Orange Tree Technologies has been advised of the possibility of such
00040 //**
                      damage. This limitation shall apply notwithstanding the failure of the
00041 //**
                      essential purpose of any limited remedies herein.
00042 //**
00043 //****************************
00044
00045 #define CRT SECURE NO WARNINGS
00046 #ifdef WINGCC
00047 #define __USE_W32_SOCKETS
00048 #endif
00049
00050 #include <stdint.h>
00051 #include <stdlib.h>
00052 #include <stdio.h>
00053 #include <memory.h>
00054 #include "ZestETM1.h"
00055 #include "Private.h"
00056
00057 // Read/write register command structures
00058 typedef struct
00059 {
00060
         uint8_t Command;
00061
         uint8_t Addr;
00062
        uint16_t Data;
00063 } ZESTETM1_WRITE_REG_CMD;
00064 typedef struct
00065 {
00066
         uint8_t Command;
00067
         uint8_t Status;
00068
         uint8_t Dummy1[2];
00069 } ZESTETM1_WRITE_REG_RESPONSE;
00070 typedef struct
00071 {
00072
         uint8_t Command;
00073
         uint8_t Addr;
00074
         uint8_t Dummy[2];
00075 } ZESTETM1_READ_REG_CMD;
00076 typedef struct
00077 {
00078
         uint8_t Command;
00079
         uint8_t Status;
08000
         uint16_t Value;
00081 } ZESTETM1_READ_REG_RESPONSE;
00082 typedef struct
00083 {
         uint8_t Command;
00085
         uint8_t Dummy[3];
00086 } ZESTETM1_MAILBOX_INT_CMD;
00087 typedef struct
00088 {
00089
         uint8 t Command;
00090
         uint8_t Status;
00091    uint16_t Dummy;
00092 } ZESTETM1_MAILBOX_INT_RESPONSE;
00093 #define ZESTETM1_COMMAND_SPI 0xee
00094 #define ZESTETM1_COMMAND_WRITE_REG 0xf6
00095 #define ZESTETM1_COMMAND_READ_REG 0xf7
00096 #define ZESTETM1_COMMAND_MAILBOX_INT 0xf8
00097
00098 /***********************************
```

7.22 Data.c 317

```
00099 \star Open a connection to a ZestETM1 for data transfer (internal version) \star
00101 ZESTETM1_STATUS ZestETM1_OpenConnection(ZESTETM1_CARD_INFO *CardInfo,
                                               ZESTETM1_CONNECTION_TYPE Type,
00102
00103
                                                uint16_t Port,
00104
                                                uint16_t LocalPort,
00105
                                                ZESTETM1_CONNECTION *Connection)
00106 {
00107
          ZESTETM1_CONNECTION_STRUCT *NewStruct;
          SOCKET Socket = -1;
00108
          char AddrBuffer[32];
00109
00110
          char PortBuffer[32];
00111
00112
          if (Connection==NULL || CardInfo==NULL)
00113
00114
              return ZESTETM1_NULL_PARAMETER;
00115
00116
00117
          // Allocate data structure
00118
          NewStruct = malloc(sizeof(ZESTETM1_CONNECTION_STRUCT));
00119
          if (NewStruct==NULL)
00120
00121
              return ZESTETM1 OUT OF MEMORY;
00122
00123
00124
          // Build target addresses
          sprintf(AddrBuffer, "%d.%d.%d", CardInfo->IPAddr[0], CardInfo->IPAddr[1],
00125
00126
                           CardInfo->IPAddr[2], CardInfo->IPAddr[3]);
00127
          sprintf(PortBuffer, "%d", Port);
00128
00129
          if (Type==ZESTETM1_TYPE_UDP)
00130
00131
              // Open UDP connection
00132
              struct sockaddr_in SourceIP;
              int SourceLen = (int)sizeof(struct sockaddr_in);
00133
00134
              Socket = socket(AF_INET, SOCK_DGRAM, IPPROTO_UDP);
00135
              if (Socket<0)
00136
                  return ZESTETM1_SOCKET_ERROR;
00137
              SourceIP.sin_family = AF_INET;
              SourceIP.sin_addr.s_addr = htonl(INADDR_ANY);
00138
00139
              SourceIP.sin_port = htons(LocalPort);
00140
              bind(Socket, (const struct sockaddr *)&SourceIP, SourceLen);
00141
00142
              NewStruct->Target.sin_family = AF_INET;
00143
              NewStruct->Target.sin_addr.s_addr = inet_addr(AddrBuffer);
00144
              NewStruct->Target.sin_port = htons(atoi(PortBuffer));
00145
00146
          else if (Type==ZESTETM1_TYPE_TCP)
00147
00148
              // Open TCP connection
00149
              struct addrinfo *AddrResult = NULL,
00150
                               \starPtr = NULL,
00151
                               Hints;
00152
              int Result;
00153
              struct sockaddr_in SourceIP;
00154
              int SourceLen = (int)sizeof(struct sockaddr_in);
00155
00156
              memset(&Hints, 0, sizeof(Hints));
00157
              Hints.ai_family = AF_UNSPEC;
00158
              Hints.ai_socktype = SOCK_STREAM;
              Hints.ai_protocol = IPPROTO_TCP;
00159
00160
00161
              // Resolve the server address and port
00162
              Result = getaddrinfo(AddrBuffer, PortBuffer, &Hints, &AddrResult);
00163
              if (Result!=0)
00164
              {
00165
                  free (NewStruct);
00166
                  return ZESTETM1_SOCKET_ERROR;
00167
00168
00169
              // Attempt to connect to an address until one succeeds
00170
              for (Ptr=AddrResult; Ptr!=NULL; Ptr=Ptr->ai_next)
00171
00172
                   // Create a SOCKET for connecting to server
00173
                  Socket = socket(Ptr->ai_family, Ptr->ai_socktype,
00174
                                   Ptr->ai_protocol);
00175
                   if (Socket<0)
00176
00177
                       freeaddrinfo(AddrResult);
00178
                      free (NewStruct);
                       return ZESTETM1_SOCKET_ERROR;
00179
```

```
00180
00181
00182
                 // Connect to ZestETM1
00183
                 Result = connect(Socket, Ptr->ai_addr, (int)Ptr->ai_addrlen);
00184
                 if (Result<0)</pre>
00185
00186
                     closesocket(Socket);
00187
                     Socket = -1;
00188
                     continue;
00189
00190
                 break:
00191
00192
00193
             SourceIP.sin_family = AF_INET;
00194
             SourceIP.sin_addr.s_addr = htonl(INADDR_ANY);
00195
             SourceIP.sin_port = 0;
00196
             bind(Socket, (const struct sockaddr *)&SourceIP, SourceLen);
00197
             freeaddrinfo(AddrResult);
00198
00199
         else
00200
         {
00201
             free (NewStruct);
00202
             return ZESTETM1_INVALID_CONNECTION_TYPE;
00203
00204
00205
         if (Socket==-1)
00206
00207
             free (NewStruct);
             return ZESTETM1_SOCKET_ERROR;
00208
00209
00210
00211
         NewStruct->Magic = ZESTETM1_CONNECTION_HANDLE_MAGIC;
00212
         NewStruct->Type = Type;
         NewStruct->Port = Port;
00213
         NewStruct->LocalPort = LocalPort;
00214
         NewStruct->Socket = Socket;
00215
00216
         NewStruct->CardInfo = CardInfo;
00217
         *Connection = NewStruct;
00218
         return ZESTETM1 SUCCESS:
00219
00220 }
00221
00222 /****************************
00223 \star Close a connection to a ZestETM1 (internal version \star
00225 ZESTETM1_STATUS ZestETM1_CloseConnection(ZESTETM1_CONNECTION Connection)
00226 {
00227
         ZESTETM1_CONNECTION_STRUCT *Conn = (ZESTETM1_CONNECTION_STRUCT *)Connection;
00228
00229
         if (Conn==NULL)
00230
00231
             return ZESTETM1_NULL_PARAMETER;
00232
00233
         if (Conn->Magic!=ZESTETM1_CONNECTION_HANDLE_MAGIC)
00234
        {
00235
             return ZESTETM1_ILLEGAL_CONNECTION;
00236
00237
         // Cleanup
00238
00239
         closesocket(Conn->Socket);
00240
         Conn->Magic = 0;
00241
         free (Conn);
00242
00243
         return ZESTETM1_SUCCESS;
00244 }
00247 * Write data to ZestETM1 connection (internal version) *
00249 static ZESTETM1_STATUS ZestETM1_WriteData(ZESTETM1_CONNECTION Connection,
00250
                                            void *Buffer,
00251
                                            uint32 t Length,
00252
                                            unsigned long *Written,
00253
                                            uint32_t Timeout)
00254 {
         ZESTETM1_CONNECTION_STRUCT *Conn = (ZESTETM1_CONNECTION_STRUCT *)Connection;
00255
00256
         fd set WriteFDS;
         struct timeval Time;
00257
00258
         uint32_t BufferPtr = 0;
00259
         int Result;
00260
         int TargetLen = (int)sizeof(struct sockaddr_in);
```

7.22 Data.c 319

```
00261
          unsigned int MaxSize;
00262
00263
          if (Conn==NULL)
00264
00265
               if (Written!=NULL) *Written = 0;
00266
              return ZESTETM1_NULL_PARAMETER;
00267
00268
          if (Conn->Magic!=ZESTETM1_CONNECTION_HANDLE_MAGIC)
00269
00270
              if (Written!=NULL) *Written = 0;
00271
              return ZESTETM1_ILLEGAL_CONNECTION;
00272
00273
00274
          if (Conn->Type==ZESTETM1_TYPE_UDP)
00275
          {
00276
              // Get maximum transfer size
00277 #if defined(MSVC) || defined(WINGCC)
00278
              int MaxSizeLen = sizeof(MaxSize);
              getsockopt(Conn->Socket, SOL_SOCKET, SO_MAX_MSG_SIZE,
00279
                         (char *) &MaxSize, &MaxSizeLen);
00280
00281 #else
00282
              //FIXME: Linux doesn't support SO_MAX_MSG_SIZE
00283
              MaxSize = 65507;
00284 #endif
00285
         }
00286
00287
          do
00288
          {
00289
              uint32_t Bytes = Length-BufferPtr;
00290
              if (Conn->Type==ZESTETM1_TYPE_UDP && Bytes>MaxSize)
00291
00292
                  Bytes = MaxSize;
00293
00294
              // Wait for socket to become ready
00295
                   int CurTime = 0;
00296
00297
                  while (CurTime!=Timeout)
00298
00299
                       int T = (Timeout-CurTime) < 1000 ? Timeout-CurTime : 1000;</pre>
                       FD_ZERO(&WriteFDS);
00300
00301
                       FD_SET(Conn->Socket, &WriteFDS);
00302
                       Time.tv_sec = T/1000;
                       Time.tv_usec = (T%1000) *1000;
00303
00304
                       Result = select((int)Conn->Socket+1, NULL, &WriteFDS, NULL, &Time);
00305
                       if (Result>0) break;
00306
                       CurTime+=T;
00307
00308
00309
              if (Result<0 || !FD_ISSET(Conn->Socket, &WriteFDS))
00310
00311
                   if (Written!=NULL) *Written = BufferPtr;
00312
                   return ZESTETM1_TIMEOUT;//FIXME: Any other errors?
00313
00314
00315
              if (Conn->Type==ZESTETM1_TYPE_UDP)
00316
00317
                   Result = sendto(Conn->Socket, (char *)Buffer+BufferPtr,
00318
                                   Bytes, 0,
00319
                                   (struct sockaddr *) &Conn->Target, TargetLen);
00320
              else if (Conn->Type==ZESTETM1_TYPE_TCP)
00321
00322
00323
                  Result = send(Conn->Socket, (char *)Buffer+BufferPtr,
00324
                                 Bytes, 0);
00325
              }
00326
              else
00327
              {
00328
                   if (Written!=NULL) *Written = BufferPtr;
00329
                  return ZESTETM1_ILLEGAL_CONNECTION;
00330
              }
00331
00332
              // Update counters
00333
              if (Result>0)
00334
              {
00335
                  BufferPtr += Result;
00336
              else if (Result==0)
00337
00338
00339
                   // Connection closed
                  if (Written!=NULL) *Written = BufferPtr;
return ZESTETM1_SOCKET_CLOSED;
00340
00341
```

```
00342
00343
              else
00344
                  // Socket error
00345
00346
                  if (Written!=NULL) *Written = BufferPtr;
00347
                 return ZESTETM1_SOCKET_ERROR;
00348
00349
         } while (Result>0 && BufferPtr<Length);
00350
          if (Written!=NULL) *Written = BufferPtr;
00351
00352
          return ZESTETM1_SUCCESS;
00353 }
00354
00356 * Read data from a ZestETM1 connection (internal version) *
00358 static ZESTETM1_STATUS ZestETM1_ReadData(ZESTETM1_CONNECTION Connection,
00359
                                             void *Buffer,
00360
                                             uint32_t Length,
00361
                                             unsigned long *Read,
00362
                                             uint32_t Timeout)
00363 {
00364
          ZESTETM1_CONNECTION_STRUCT *Conn = (ZESTETM1_CONNECTION_STRUCT *)Connection;
00365
          fd_set ReadFDS;
00366
          struct timeval Time:
00367
          uint32 t BufferPtr = 0;
00368
          int Result:
00369
          struct sockaddr_in Target;
00370
          int TargetLen = (int)sizeof(struct sockaddr_in);
00371
00372
          if (Conn==NULL)
00373
              if (Read!=NULL) *Read = 0;
00374
00375
              return ZESTETM1 NULL PARAMETER;
00376
00377
         if (Conn->Magic!=ZESTETM1_CONNECTION_HANDLE_MAGIC)
00378
              if (Read!=NULL) *Read = 0;
00379
00380
              return ZESTETM1_ILLEGAL_CONNECTION;
00381
00382
00383
00384
         {
              int OK = 1;
00385
00386
00387
              // Wait for socket to become ready
00388
00389
                  int CurTime = 0;
00390
                  while (CurTime!=Timeout)
00391
00392
                      int T = (Timeout-CurTime) < 1000 ? Timeout-CurTime : 1000;</pre>
00393
                      FD_ZERO(&ReadFDS);
00394
                      FD_SET(Conn->Socket, &ReadFDS);
00395
                      Time.tv_sec = T/1000;
00396
                      Time.tv\_usec = (T%1000)*1000;
00397
                      Result = select((int)Conn->Socket+1, &ReadFDS, NULL, NULL, &Time);
00398
                      if (Result>0) break;
00399
                      CurTime+=T;
00400
00401
00402
              if (Result<0 || !FD_ISSET(Conn->Socket, &ReadFDS))
00403
00404
                  if (Read!=NULL) *Read = BufferPtr;
00405
                  return ZESTETM1_TIMEOUT;//FIXME: Any other errors?
00406
00407
00408
              if (Conn->Type==ZESTETM1_TYPE_UDP)
00409
              {
00410
                  Result = recvfrom(Conn->Socket, (char *)Buffer+BufferPtr,
                                    Length-BufferPtr, 0,
00411
00412
                                    (struct sockaddr *) & Target, & TargetLen);
00413
                  if (Conn->LocalPort!=0 && Target.sin_port!=htons(Conn->LocalPort))
00414
                      OK = 0;
00415
00416
              else if (Conn->Type==ZESTETM1_TYPE_TCP)
00417
                  Result = recv(Conn->Socket, (char *)Buffer+BufferPtr,
00418
00419
                               Length-BufferPtr, 0);
00420
              }
00421
              else
00422
```

7.22 Data.c 321

```
00423
                  if (Read!=NULL) *Read = BufferPtr;
00424
                 return ZESTETM1_ILLEGAL_CONNECTION;
00425
00426
00427
             // Update counters
00428
             if (Result>0 && OK==1)
00429
             {
00430
                  BufferPtr += Result;
00431
00432
             else if (Result==0)
00433
00434
                  // Connection closed
00435
                  if (Read!=NULL) *Read = BufferPtr;
                 return ZESTETM1_SOCKET_CLOSED;
00436
00437
00438
             else if (Result<0)</pre>
00439
00440
                  // Socket error
                  if (Read!=NULL) *Read = BufferPtr;
00441
                 return ZESTETM1_SOCKET_ERROR;
00442
00443
00444
         } while (Result>0 && BufferPtr<Length);</pre>
00445
00446
         if (Read!=NULL) *Read = BufferPtr;
         return ZESTETM1_SUCCESS;
00447
00448 }
00449
00450 /**************************
00451 \star Send a control command to GigEx and get response \star
00453 ZESTETM1_STATUS ZestETM1_SendCommand(ZESTETM1_CARD_INFO *CardInfo,
00454
                                        ZESTETM1 CONNECTION Connection.
00455
                                         void *WriteBuffer, uint32_t WriteLen,
00456
                                         void *ReadBuffer, uint32_t ReadLen,
00457
                                         int WaitForAck)
00458 {
00459
          ZESTETM1_STATUS Result;
00460
          unsigned long Written;
00461
         unsigned long Received;
00462
00463
          // Send/receive data
00464
          Result = ZestETM1_WriteData(Connection, WriteBuffer, WriteLen, &Written,
00465
                                     CardInfo->Timeout);
00466
          if (Result!=ZESTETM1_SUCCESS)
00467
00468
              return Result:
00469
00470
          if (Written!=WriteLen)
00471
00472
              return ZESTETM1_INTERNAL_ERROR;
00473
00474
          *((uint8_t *)ReadBuffer) = 0;
00475
          if (WaitForAck==1)
00476
00477
              Result = ZestETM1_ReadData(Connection, ReadBuffer, ReadLen, &Received,
00478
                                       CardInfo->Timeout);
00479
              if (Result!=ZESTETM1_SUCCESS)
00480
              {
00481
                  return Result;
00482
00483
              if (Received!=ReadLen)
00484
00485
                  return ZESTETM1_INTERNAL_ERROR;
00486
00487
00488
          return ZESTETM1_SUCCESS;
00489
00490 }
00491
00493 * Send data to/from SPI port on GigEx device *
00494 ******************************
00495 ZESTETM1_STATUS ZestETM1_SPIReadWrite(ZESTETM1_CARD_INFO *CardInfo,
00496
                                          ZESTETM1_CONNECTION Connection,
00497
                                          int Device.
00498
                                          int WordLen, uint32_t *WriteData,
00499
                                          uint32_t *ReadData, uint32_t Length,
int ReleaseCS, int WaitForAck)
00500
00501 {
00502
          uint8 t Buffer[65536];
00503
         uint32 t *BufPtr;
```

```
uint32_t i;
00504
00505
          ZESTETM1_STATUS Result;
00506
          // Build command
00507
00508
          Buffer[0] = ZESTETM1_COMMAND_SPI;// Command byte
                                      // SPI device
00509
          Buffer[1] = Device;
00510
          Buffer[2] = WordLen;
                                         // Word length
00511
          Buffer[3] = ReleaseCS;
                                         // Release CS on completion
00512
          BufPtr = (uint32_t *) (Buffer+4);
00513
00514
          BufPtr[0] = WriteData==NULL ? 0 : ZESTETM1_REVERSE(Length);
                                                                        // Num words
00515
          BufPtr[1] = ReadData==NULL ? 0 : ZESTETM1_REVERSE(Length);
00516
          if (WriteData!=NULL)
00517
00518
              for (i=0; i<Length; i++)</pre>
00519
00520
                 BufPtr[2+i] = ZESTETM1_REVERSE(WriteData[i]);
00521
00522
00523
         else
00524
             memset (BufPtr+2, 0, Length * 4);
00525
00526
          // Send command
00527
          Result = ZestETM1_SendCommand(CardInfo, Connection, Buffer,
00528
                                      WriteData == NULL ? 12 : 12 + Length * 4,
00529
                                      Buffer, ReadData == NULL ? 4 : 4 + Length * 4,
00530
                                      WaitForAck):
00531
         if (WaitForAck==0)
00532
             return Result;
00533
00534
         // Extract results
         if (Buffer[0]!=ZESTETM1_COMMAND_SPI || Buffer[1]!=0)
00535
00536
00537
             return ZESTETM1 INTERNAL ERROR;
00538
         }
00539
         BufPtr = (uint32_t *)(Buffer+4);
00540
00541
         if (ReadData!=NULL)
00542
00543
             for (i=0; i<Length; i++)</pre>
00544
                 ReadData[i] = ZESTETM1_REVERSE(BufPtr[i]);
00545
00546
00547
         }
00548
00549
          return ZESTETM1_SUCCESS;
00550 }
00551
00552 /****************************
00553 * Open a connection to a ZestETM1 for data transfer *
00555 ZESTETM1_STATUS ZestETM1OpenConnection(ZESTETM1_CARD_INFO *CardInfo,
00556
                                          ZESTETM1_CONNECTION_TYPE Type,
                                          uint16_t Port,
00557
00558
                                          uint16_t LocalPort,
00559
                                          ZESTETM1_CONNECTION *Connection)
00560 {
00561
          ZESTETM1_STATUS Result;
00562
00563
          Result = ZestETM1_OpenConnection(CardInfo, Type, Port, LocalPort, Connection);
00564
          if (Result!=ZESTETM1_SUCCESS)
00565
00566
             ZESTETM1_ERROR("ZestETM1OpenConnection", Result);
00567
00568
00569
          return ZESTETM1_SUCCESS;
00570 }
00571
00572 /*******************
00573 * Close a connection to a ZestETM1 *
0.0574 *******************************
00575 ZESTETM1_STATUS ZestETM1CloseConnection(ZESTETM1_CONNECTION Connection)
00576 {
00577
          ZESTETM1_STATUS Result;
          ZESTETM1_CONNECTION_STRUCT *Conn = (ZESTETM1_CONNECTION_STRUCT*)Connection;
00578
00579
          Result = ZestETM1 CloseConnection(Connection);
00580
00581
          if (Result!=ZESTETM1 SUCCESS)
00582
00583
             ZESTETM1 ERROR CONN("ZestETM1CloseConnection", Result);
00584
```

7.22 Data.c 323

```
00585
00586
          return ZESTETM1_SUCCESS;
00587 }
00588
00589 /*****************
00590 * Write data to ZestETM1 connection *
00591 ****************************
00592 ZESTETM1_STATUS ZestETM1WriteData(ZESTETM1_CONNECTION Connection,
00593
                                      void *Buffer,
00594
                                      unsigned long Length,
00595
                                      unsigned long *Written,
00596
                                      unsigned long Timeout)
00597 {
00598
          ZESTETM1_STATUS Result;
00599
          ZESTETM1_CONNECTION_STRUCT *Conn = (ZESTETM1_CONNECTION_STRUCT*)Connection;
00600
00601
          Result = ZestETM1_WriteData(Connection, Buffer, Length, Written, Timeout);
00602
         if (Result!=ZESTETM1_SUCCESS)
00603
00604
              ZESTETM1_ERROR_CONN("ZestETM1WriteData", Result);
00605
00606
00607
          return ZESTETM1 SUCCESS;
00608 }
00609
00610 /******************
00611 * Read data from a ZestETM1 connection *
00612 ******************************
00613 ZESTETM1_STATUS ZestETM1ReadData(ZESTETM1_CONNECTION Connection,
00614
                                     void *Buffer.
00615
                                     unsigned long Length,
00616
                                     unsigned long *Read,
00617
                                     unsigned long Timeout)
00618 {
          ZESTETM1 STATUS Result:
00619
          ZESTETM1_CONNECTION_STRUCT *Conn = (ZESTETM1_CONNECTION_STRUCT*)Connection;
00620
00621
00622
          Result = ZestETM1_ReadData(Connection, Buffer, Length, Read, Timeout);
00623
          if (Result!=ZESTETM1_SUCCESS)
00624
          {
00625
              ZESTETM1_ERROR_CONN("ZestETM1ReadData", Result);
00626
00627
00628
          return ZESTETM1 SUCCESS;
00629 }
00630
00631
00632 /*****************
00633 * Read data from master SPI port *
00634 ****************************
00635 ZESTETM1_STATUS ZestETM1SPIReadWrite(ZESTETM1_CARD_INFO *CardInfo, ZESTETM1_SPI_RATE Rate,
00636
                                            int WordLen, void *WriteData,
00637
                                            void *ReadData, unsigned long Length,
00638
                                           int ReleaseCS)
00639 {
00640
          ZESTETM1_STATUS Result;
00641
          ZESTETM1_CONNECTION Connection;
00642
          unsigned long RateVal = Rate==ZESTETM1_SPI_RATE_35 ? ZESTETM1_RATE_40MHz :
00643
                                  Rate==ZESTETM1_SPI_RATE_17_5 ? ZESTETM1_RATE_20MHz : ZESTETM1_RATE_10MHz;
00644
00645
          if (CardInfo==NULL || (WriteData==NULL && ReadData==NULL))
00646
00647
              ZESTETM1_ERROR("ZestETM1SPIReadWrite", ZESTETM1_NULL_PARAMETER);
00648
00649
          if (WordLen<1 || WordLen>32 || Length>16384)
00650
         {
              ZESTETM1_ERROR("ZestETM1SPIReadWrite", ZESTETM1_ILLEGAL_PARAMETER);
00651
00652
         }
00653
00654
          Result = ZestETM1_OpenConnection(CardInfo, ZESTETM1_TYPE_TCP, CardInfo->ControlPort, 0, &Connection);
00655
          if (Result!=ZESTETM1_SUCCESS)
00656
          {
00657
              ZESTETM1_ERROR("ZestETM1SPIReadWrite", Result);
00658
00659
      Result = ZestETM1_SPIReadWrite(CardInfo, Connection, ZESTETM1_USER_DEVICE_ID|RateVal, WordLen,
(uint32_t *)WriteData, (uint32_t *)ReadData, Length, ReleaseCS, 1);
00660
00661
          if (Result!=ZESTETM1_SUCCESS)
00662
          {
00663
              ZESTETM1 ERROR("ZestETM1SPIReadWrite", Result);
00664
```

```
00665
00666
          Result = ZestETM1_CloseConnection(Connection);
00667
          if (Result!=ZESTETM1_SUCCESS)
00668
00669
              ZESTETM1_ERROR("ZestETM1SPIReadWrite", Result);
00670
00671
00672
          return ZESTETM1_SUCCESS;
00673 }
00674
00675 /*******************
00676 * Read/Write user interface registers *
00678 ZESTETM1_STATUS ZestETM1WriteRegister(ZESTETM1_CARD_INFO *CardInfo, unsigned long Addr, unsigned short
     Data)
00679 {
00680
          ZESTETM1_STATUS Result;
00681
          ZESTETM1_CONNECTION Connection;
          ZESTETM1_WRITE_REG_CMD Cmd;
00682
          ZESTETM1_WRITE_REG_RESPONSE Response;
00683
00684
00685
          if (CardInfo==NULL)
00686
         {
00687
              ZESTETM1_ERROR("ZestETM1WriteRegister", ZESTETM1_NULL_PARAMETER);
00688
00689
         if (Addr>127)
00690
         {
00691
              ZESTETM1_ERROR("ZestETM1WriteRegister", ZESTETM1_ILLEGAL_PARAMETER);
00692
00693
00694
          Result = ZestETM1_OpenConnection(CardInfo, ZESTETM1_TYPE_TCP, CardInfo->ControlPort, 0, &Connection);
00695
          if (Result!=ZESTETM1_SUCCESS)
00696
              ZESTETM1_ERROR("ZestETM1WriteRegister", Result);
00697
00698
00699
00700
          // Write register to the device
00701
          //FIXME: Do we want to be able to copy this value to flash?
          Cmd.Command = ZESTETM1_COMMAND_WRITE_REG;
00702
00703
          Cmd.Addr = (uint8 t)Addr;
00704
          Cmd.Data = ((Data *8) & 0xff) | (Data & 0xff);
00705
          Result = ZestETM1_SendCommand(CardInfo, Connection,
00706
                                       &Cmd, sizeof(Cmd),
00707
                                       &Response, sizeof(Response), 1);
00708
          if (Result!=ZESTETM1_SUCCESS)
00709
              ZestETM1_CloseConnection(Connection);
ZESTETM1_ERROR("ZestETM1WriteRegister", Result);
00710
00711
00712
00713
          if (Response.Command!=Cmd.Command || Response.Status!=0)
00714
00715
              ZestETM1_CloseConnection(Connection);
00716
              ZESTETM1_ERROR("ZestETM1WriteRegister", ZESTETM1_INTERNAL_ERROR);
00717
00718
00719
          Result = ZestETM1_CloseConnection(Connection);
00720
          if (Result!=ZESTETM1_SUCCESS)
00721
          {
00722
              ZESTETM1_ERROR("ZestETM1WriteRegister", Result);
00723
00724
00725
          return ZESTETM1_SUCCESS;
00726 }
00727
00728 ZESTETM1_STATUS ZestETM1ReadRegister(ZESTETM1_CARD_INFO *CardInfo, unsigned long Addr, unsigned short
00729 {
00730
          ZESTETM1_STATUS Result;
00731
          ZESTETM1_CONNECTION Connection;
          ZESTETM1_READ_REG_CMD Cmd;
00732
00733
          ZESTETM1_READ_REG_RESPONSE Response;
00734
00735
          if (CardInfo==NULL || Data==NULL)
00736
          {
00737
              ZESTETM1_ERROR("ZestETM1ReadRegister", ZESTETM1_NULL_PARAMETER);
00738
00739
          if (Addr>127)
00740
00741
              ZESTETM1 ERROR ("ZestETM1ReadRegister", ZESTETM1 ILLEGAL PARAMETER);
00742
00743
```

7.22 Data.c 325

```
00744
          Result = ZestETM1_OpenConnection(CardInfo, ZESTETM1_TYPE_TCP, CardInfo->ControlPort, 0, &Connection);
00745
          if (Result!=ZESTETM1_SUCCESS)
00746
00747
              ZESTETM1_ERROR("ZestETM1ReadRegister", Result);
00748
00749
00750
          // Read register from the device
00751
          Cmd.Command = ZESTETM1_COMMAND_READ_REG;
00752
          Cmd.Addr = (uint8_t)Addr;
00753
          Result = ZestETM1_SendCommand(CardInfo, Connection,
00754
                                        &Cmd, sizeof(Cmd),
00755
                                        &Response, sizeof(Response), 1);
00756
          if (Result!=ZESTETM1_SUCCESS)
00757
00758
              ZestETM1_CloseConnection(Connection);
00759
              ZESTETM1_ERROR("ZestETM1ReadRegister", Result);
00760
00761
          if (Response.Command!=Cmd.Command || Response.Status!=0)
00762
00763
              ZestETM1_CloseConnection(Connection);
00764
              ZESTETM1_ERROR("ZestETM1ReadRegister", ZESTETM1_INTERNAL_ERROR);
00765
00766
          *Data = ((Response.Value»8)&0xff) | ((Response.Value&0xff) < 8);
00767
00768
          Result = ZestETM1 CloseConnection(Connection);
00769
          if (Result!=ZESTETM1 SUCCESS)
00770
00771
              ZESTETM1 ERROR("ZestETM1ReadRegister", Result);
00772
00773
00774
          return ZESTETM1 SUCCESS;
00775 }
00776
00777 /**************
00778 * Set mailbox interrupt *
00779 ****************/
00780 ZESTETM1_STATUS ZestETM1SetInterrupt(ZESTETM1_CARD_INFO *CardInfo)
00781 {
00782
          ZESTETM1_STATUS Result;
00783
          ZESTETM1 CONNECTION Connection:
          ZESTETM1_MAILBOX_INT_CMD Cmd;
00784
00785
          ZESTETM1_MAILBOX_INT_RESPONSE Response;
00786
00787
          if (CardInfo==NULL)
00788
00789
              ZESTETM1_ERROR("ZestETM1SetInterrupt", ZESTETM1_NULL_PARAMETER);
00790
00791
00792
          Result = ZestETM1_OpenConnection(CardInfo, ZESTETM1_TYPE_TCP, CardInfo->ControlPort, 0, &Connection);
00793
          if (Result!=ZESTETM1_SUCCESS)
00794
00795
              ZESTETM1_ERROR("ZestETM1SetInterrupt", Result);
00796
00797
00798
          // Send command to set interrupt
00799
          Cmd.Command = ZESTETM1_COMMAND_MAILBOX_INT;
00800
          Result = ZestETM1_SendCommand(CardInfo, Connection,
00801
                                        &Cmd, sizeof(Cmd),
                                        &Response, sizeof(Response), 1);
00802
00803
          if (Result!=ZESTETM1_SUCCESS)
00804
00805
              ZestETM1_CloseConnection(Connection);
00806
              ZESTETM1_ERROR("ZestETM1SetInterrupt", Result);
00807
00808
          if (Response.Command!=Cmd.Command || Response.Status!=0)
00809
00810
              ZestETM1_CloseConnection(Connection);
00811
              ZESTETM1_ERROR("ZestETM1SetInterrupt", ZESTETM1_INTERNAL_ERROR);
00812
00813
00814
          Result = ZestETM1_CloseConnection(Connection);
00815
          if (Result!=ZESTETM1 SUCCESS)
00816
00817
              ZESTETM1 ERROR("ZestETM1SetInterrupt", Result);
00818
00819
00820
          return ZESTETM1_SUCCESS;
00821 }
```

# 7.23 C:/Users/hill35/git/camera\_python/nsCamera/comms/ZestETM1/Error.c File Reference

```
#include <stdint.h>
#include "ZestETM1.h"
#include "Private.h"
```

#### **Functions**

- ZESTETM1 STATUS ZestETM1RegisterErrorHandler (ZESTETM1 ERROR FUNC Function)
- ZESTETM1\_STATUS ZestETM1GetErrorMessage (ZESTETM1\_STATUS Status, char \*\*Buffer)

#### **Variables**

00078 }

- char \* ZestETM1\_ErrorStrings []
- ZESTETM1\_ERROR\_FUNC ZestETM1\_ErrorHandler

## 7.23.1 Function Documentation

## 7.23.1.1 ZestETM1GetErrorMessage()

```
ZESTETM1_STATUS ZestETM1GetErrorMessage (
              ZESTETM1_STATUS Status,
              char ** Buffer )
Definition at line 84 of file Error.c.
00086 {
00087
          if (Status>ZESTETM1_MAX_ERROR | |
00088
              (Status<ZESTETM1_ERROR_BASE && Status>=ZESTETM1_MAX_WARNING) ||
00089
              (Status<ZESTETM1_WARNING_BASE && Status>=ZESTETM1_MAX_INFO))
00090
00091
             return ZESTETM1_ILLEGAL_STATUS_CODE;
00092
00093
00094
         *Buffer = ZESTETM1_ERROR_STRING(Status);
00095
         return ZESTETM1_SUCCESS;
00096 }
```

## 7.23.1.2 ZestETM1RegisterErrorHandler()

7.24 Error.c 327

## 7.23.2 Variable Documentation

#### 7.23.2.1 ZestETM1\_ErrorHandler

```
{\tt ZESTETM1\_ERROR\_FUNC} \  \, {\tt ZestETM1\_ErrorHandler}
```

Definition at line 67 of file Error.c.

## 7.23.2.2 ZestETM1\_ErrorStrings

```
char* ZestETM1_ErrorStrings[]
Initial value:
    "Success (no error)",
    "Error communicating with socket",
    "An unspecified internal error occurred",
    "Status code is out of range",
    "NULL was used illegally as one of the parameter values",
    "Not enough memory to complete the requested operation",
    "The requested connection type is invalid",
    "The requested connection is invalid",
    "The connection was closed unexpectedly",
    "Operation timed out",
    "One of the parameters has an illegal value",
}
Definition at line 52 of file Error.c.
00053 {
00054
          "Success (no error)",
          "Error communicating with socket",
00055
00056
          "An unspecified internal error occurred",
          "Status code is out of range",
00057
          "NULL was used illegally as one of the parameter values",
00058
          "Not enough memory to complete the requested operation",
00059
00060
          "The requested connection type is invalid",
00061
          "The requested connection is invalid",
          "The connection was closed unexpectedly",
00062
          "Operation timed out",
00063
          "One of the parameters has an illegal value",
00064
00065
00066 };
```

## 7.24 Error.c

#### Go to the documentation of this file.

```
00001 // File:
                      Error.c
00002 //
00003 // Purpose:
00004 //
             ZestETM1 Host Library
00005 //
             Error functions
00006 //
00007 // Version: 1.00
00008 // Date: 11/10/12
00010 // Copyright (C) 2012 Orange Tree Technologies Ltd. All rights reserved.
00011 // Orange Tree Technologies grants the purchaser of a ZestETM1 the right to use and
00012 // modify this source code in any form in designs that target the ZestETM1.
00013 // Orange Tree Technologies prohibits the use of this source code or any modification of
00014 // it in any form in designs that target any other hardware unless the purchaser of the 00015 // ZestETM1 has purchased the appropriate licence from Orange Tree Technologies.
00016 // Contact Orange Tree Technologies if you want to purchase such a licence.
```

```
00017
00019 //**
00020 //**
          Disclaimer: LIMITED WARRANTY AND DISCLAIMER. These designs are
              provided to you "as is". Orange Tree Technologies and its licensors
00021 //**
00022 //**
                     make and you receive no warranties or conditions, express, implied,
00023 //**
                     statutory or otherwise, and Orange Tree Technologies specifically
00024 //**
                    disclaims any implied warranties of merchantability, non-infringement,
00025 //**
                    or fitness for a particular purpose. Orange Tree Technologies does not
00026 //**
                    warrant that the functions contained in these designs will meet your
00027 //**
                    requirements, or that the operation of these designs will be
00028 //**
                    uninterrupted or error free, or that defects in the Designs will be
00029 //**
                    corrected. Furthermore, Orange Tree Technologies does not warrant or
00030 //**
                    make any representations regarding use or the results of the use of the
00031 //**
                    designs in terms of correctness, accuracy, reliability, or otherwise.
00032 //**
00033 //**
                    LIMITATION OF LIABILITY. In no event will Orange Tree Technologies
                     or its licensors be liable for any loss of data, lost profits, cost or
00034 //**
00035 //**
                     procurement of substitute goods or services, or for any special,
00036 //**
                     incidental, consequential, or indirect damages arising from the use or
                     operation of the designs or accompanying documentation, however caused
00037 //**
00038 //**
                     and on any theory of liability. This limitation will apply even if
00039 //**
                     Orange Tree Technologies has been advised of the possibility of such
00040 //**
                    damage. This limitation shall apply notwithstanding the failure of the
00041 //**
                    essential purpose of any limited remedies herein.
00042 //**
00043 //***************************
00044
00045 #include <stdint.h>
00046 #include "ZestETM1.h"
00047 #include "Private.h"
00048
00049 /*****************************
00050 * Globals
00052 char *ZestETM1_ErrorStrings[] =
00053 {
        "Success (no error)",
00054
        "Error communicating with socket", \ 
00055
00056
         "An unspecified internal error occurred",
00057
        "Status code is out of range",
00058
         "NULL was used illegally as one of the parameter values",
00059
         "Not enough memory to complete the requested operation",
00060
         "The requested connection type is invalid",
00061
        "The requested connection is invalid",
00062
         "The connection was closed unexpectedly",
00063
        "Operation timed out",
00064
        "One of the parameters has an illegal value",
00065
00066 };
00067 ZESTETM1_ERROR_FUNC ZestETM1_ErrorHandler;
00068
00069
00070 /****************************
00071 * Register a user error handling function to be called
00072 \star Set to NULL to disable error callbacks
00074 ZESTETM1_STATUS ZestETM1RegisterErrorHandler(ZESTETM1_ERROR_FUNC Function)
00075 {
00076
        ZestETM1_ErrorHandler = Function;
00077
        return ZESTETM1_SUCCESS;
00078 }
00079
00082 * Get a human-readable error string for a status code
00083 *****************************
00084 ZESTETM1_STATUS ZestETM1GetErrorMessage(ZESTETM1_STATUS Status,
00085
                                      char **Buffer)
00086 {
00087
        if (Status>ZESTETM1_MAX_ERROR | |
00088
            (Status<ZESTETM1_ERROR_BASE && Status>=ZESTETM1_MAX_WARNING) ||
00089
            (Status<ZESTETM1_WARNING_BASE && Status>=ZESTETM1_MAX_INFO))
00090
00091
            return ZESTETM1_ILLEGAL_STATUS_CODE;
00092
        }
00093
00094
        *Buffer = ZESTETM1_ERROR_STRING(Status);
00095
        return ZESTETM1 SUCCESS;
00096 }
00097
```

00098

# 7.25 C:/Users/hill35/git/camera\_python/nsCamera/comms/ZestETM1/Main.c File Reference

```
#include <stdint.h>
#include <sys/socket.h>
#include "ZestETM1.h"
#include "Private.h"
```

#### **Functions**

- ZESTETM1\_STATUS ZestETM1Init (void)
- ZESTETM1\_STATUS ZestETM1Close (void)

## 7.25.1 Function Documentation

## 7.25.1.1 ZestETM1Close()

# 7.25.1.2 ZestETM1Init()

# Definition at line 58 of file Main.c.

```
00060 #if defined(MSVC) || defined(WINGCC)
         WORD VersionRequested;
00061
       WSADATA WSAData;
00062
00063
         int Error;
00064
00065
         VersionRequested = MAKEWORD(2, 2);
00066
        Error = WSAStartup(VersionRequested, &WSAData);
        if (Error!=0)
00067
00068
00069
             ZESTETM1_ERROR_GENERAL("ZestETM1Init", ZESTETM1_SOCKET_ERROR);
00070
00071
00072
         // Confirm that the WinSock DLL supports 2.2.
00073
         // Note that if the DLL supports versions greater
```

```
00074
          // than 2.2 in addition to 2.2, it will still return
00075
          // 2.2 in Version since that is the version we
00076
          // requested.
00077
          if (LOBYTE(WSAData.wVersion)!=2 ||
00078
              HIBYTE(WSAData.wVersion)!=2)
00079
08000
              WSACleanup();
00081
              ZESTETM1_ERROR_GENERAL("ZestETM1Init", ZESTETM1_SOCKET_ERROR);
00082
00083 #endif
00084
          return ZESTETM1_SUCCESS;
00086 }
```

## 7.26 Main.c

## Go to the documentation of this file.

```
00001 // File:
                  Main.c
00002 //
00003 // Purpose:
00004 //
         ZestETM1 Host Library
00005 //
          Main functions
00006 //
00007 // Version: 1.00
00008 // Date: 11/10/12
00009
00010 // Copyright (C) 2012 Orange Tree Technologies Ltd. All rights reserved.
00011 // Orange Tree Technologies grants the purchaser of a ZestETM1 the right to use and
00012 // modify this source code in any form in designs that target the {\tt ZestETM1}.
00013 // Orange Tree Technologies prohibits the use of this source code or any modification of
00014 // it in any form in designs that target any other hardware unless the purchaser of the
00015 // ZestETM1 has purchased the appropriate licence from Orange Tree Technologies.
00016 // Contact Orange Tree Technologies if you want to purchase such a licence.
00017
00019 //**
00020 //** Disclaimer: LIMITED WARRANTY AND DISCLAIMER. These designs are
00021 //**
                      provided to you "as is". Orange Tree Technologies and its licensors
00022 //**
                      make and you receive no warranties or conditions, express, implied,
00023 //**
                      statutory or otherwise, and Orange Tree Technologies specifically
00024 //**
                     disclaims any implied warranties of merchantability, non-infringement,
00025 //**
                     or fitness for a particular purpose. Orange Tree Technologies does not
00026 //**
                      warrant that the functions contained in these designs will meet your
00027 //**
                      requirements, or that the operation of these designs will be
00028 //**
                      uninterrupted or error free, or that defects in the Designs will be
00029 //**
                      corrected. Furthermore, Orange Tree Technologies does not warrant or
00030 //**
                      make any representations regarding use or the results of the use of the
00031 //**
                      designs in terms of correctness, accuracy, reliability, or otherwise.
00032 //**
00033 //**
                      LIMITATION OF LIABILITY. In no event will Orange Tree Technologies
00034 //**
                      or its licensors be liable for any loss of data, lost profits, cost or
00035 //**
                      procurement of substitute goods or services, or for any special,
00036 //**
                      incidental, consequential, or indirect damages arising from the use or
00037 //**
                      operation of the designs or accompanying documentation, however caused
00038 //**
                      and on any theory of liability. This limitation will apply even if
                      Orange Tree Technologies has been advised of the possibility of such
00039 //**
00040 //**
                      damage. This limitation shall apply notwithstanding the failure of the
00041 //**
                      essential purpose of any limited remedies herein.
00042 //**
00043 //***
00045 #include <stdint.h>
00046 #if defined(MSVC) || defined(WINGCC)
00047 #include <winsock2.h>
00048 #else
00049 #include <sys/socket.h>
00050 #endif
00051 #include "ZestETM1.h"
00052 #include "Private.h"
00053
00054 /***********************
00055 * Main initialisation function.
00056 * Must be called before other ZestETM1 functions. *
00058 ZESTETM1 STATUS ZestETM1Init(void)
```

```
00059 {
00060 #if defined(MSVC) || defined(WINGCC)
         WORD VersionRequested;
        WSADATA WSAData;
00063
        int Error;
00064
00065
        VersionRequested = MAKEWORD(2, 2);
00066
      Error = WSAStartup(VersionRequested, &WSAData);
00067
       {
00069
             ZESTETM1_ERROR_GENERAL("ZestETM1Init", ZESTETM1_SOCKET_ERROR);
00070
00071
        // Confirm that the WinSock DLL supports 2.2.
         // Note that if the DLL supports versions greater
        // than 2.2 in addition to 2.2, it will still return
00075
         // 2.2 in Version since that is the version we
        // requested.
        if (LOBYTE (WSAData.wVersion)!=2 ||
00077
00078
            HIBYTE(WSAData.wVersion)!=2)
00079
08000
             WSACleanup();
00081
             ZESTETM1_ERROR_GENERAL("ZestETM1Init", ZESTETM1_SOCKET_ERROR);
00082
00083 #endif
00084
         return ZESTETM1_SUCCESS;
00085
00086 }
00087
00088 /***************************
00089 * Main clean up function.
00090 * Must be called after other ZestETM1 functions. *
00092 ZESTETM1_STATUS ZestETM1Close(void)
00093 {
00094 #if defined(MSVC) || defined(WINGCC)
00095
       WSACleanup();
00096 #endif
00097
00098
         return ZESTETM1_SUCCESS;
00099 }
```

# 7.27 C:/Users/hill35/git/camera\_python/nsCamera/comms/ZestETM1/ Private.h File Reference

```
#include <stdint.h>
#include <unistd.h>
#include <sys/socket.h>
#include <netinet/in.h>
#include <netdb.h>
#include <arpa/inet.h>
```

## Classes

struct ZESTETM1 CONNECTION STRUCT

#### **Macros**

- #define SD\_BOTH SHUT\_RDWR
- · #define closesocket close
- #define ZESTETM1 ERROR(f, x)

- #define ZESTETM1\_ERROR\_GENERAL(f, x)
- #define ZESTETM1\_ERROR\_CONN(f, x)
- #define ZESTETM1\_ERROR\_STRING(x)
- #define ZESTETM1 CONNECTION\_HANDLE\_MAGIC 0xdeadbed1
- #define ZESTETM1 DEFAULT TIMEOUT 10000
- #define ZESTETM1 RATE 40MHz (0<<4)
- #define ZESTETM1 RATE 20MHz (1<<4)
- #define ZESTETM1 RATE 10MHz (2<<4)
- #define ZESTETM1\_USER\_DEVICE\_ID (1)
- #define ZESTETM1\_REVERSE(x) ((((x)&0xff)<<24) | (((x)&0xff00)<<8) | (((x)&0xff0000)>>8) | (((x)&0xff000000)>>24))

#### **Typedefs**

typedef int SOCKET

#### **Functions**

- ZESTETM1\_STATUS ZestETM1\_OpenConnection (ZESTETM1\_CARD\_INFO \*CardInfo, ZESTETM1\_CONNECTION\_TYPE Type, uint16\_t Port, uint16\_t LocalPort, ZESTETM1\_CONNECTION \*Connection)
- ZESTETM1 STATUS ZestETM1 CloseConnection (ZESTETM1 CONNECTION Connection)
- ZESTETM1\_STATUS ZestETM1\_SendCommand (ZESTETM1\_CARD\_INFO \*CardInfo, ZESTETM1\_CONNECTION Connection, void \*WriteBuffer, uint32\_t WriteLen, void \*ReadBuffer, uint32\_t ReadLen, int WaitForAck)
- ZESTETM1\_STATUS ZestETM1\_SPIReadWrite (ZESTETM1\_CARD\_INFO \*CardInfo, ZESTETM1\_CONNECTION Connection, int Device, int WordLen, uint32\_t \*WriteData, uint32\_t \*ReadData, uint32\_t Length, int ReleaseCS, int WaitForAck)
- ZESTETM1\_STATUS ZestETM1\_WriteFlash (ZESTETM1\_CARD\_INFO \*CardInfo, uint32\_t Address, void \*Buffer, uint32\_t Length)
- ZESTETM1\_STATUS ZestETM1\_EraseFlashSector (ZESTETM1\_CARD\_INFO \*CardInfo, ZESTETM1\_CONNECTION Connection, uint32\_t Address)

## **Variables**

- ZESTETM1 ERROR FUNC ZestETM1 ErrorHandler
- char \* ZestETM1\_ErrorStrings []

#### 7.27.1 Class Documentation

#### 7.27.1.1 struct ZESTETM1 CONNECTION STRUCT

Definition at line 93 of file Private.h.

#### **Class Members**

ZESTETM1_CARD_INFO *	CardInfo	
uint16_t	LocalPort	
uint32_t	Magic	
uint16_t	Port	
SOCKET	Socket	
struct sockaddr_in	Target	
ZESTETM1_CONNECTION_TYPE	Туре	

Generated by Doxygen

# 7.27.2 Macro Definition Documentation

## 7.27.2.1 closesocket

```
#define closesocket close
```

Definition at line 58 of file Private.h.

## 7.27.2.2 SD\_BOTH

```
#define SD_BOTH SHUT_RDWR
```

Definition at line 57 of file Private.h.

## 7.27.2.3 ZESTETM1\_CONNECTION\_HANDLE\_MAGIC

```
#define ZESTETM1_CONNECTION_HANDLE_MAGIC 0xdeadbed1
```

Definition at line 92 of file Private.h.

## 7.27.2.4 ZESTETM1\_DEFAULT\_TIMEOUT

```
#define ZESTETM1_DEFAULT_TIMEOUT 10000
```

Definition at line 107 of file Private.h.

#### 7.27.2.5 ZESTETM1 ERROR

## 7.27.2.6 ZESTETM1\_ERROR\_CONN

#define ZESTETM1\_ERROR\_CONN(

```
f,
                 x )
Value:
        if (ZestETM1_ErrorHandler!=NULL) \
        ZestETM1_ErrorHandler(f, (Conn!=NULL ? Conn->CardInfo : NULL), x, ZESTETM1_ERROR_STRING(x)); \
        return (x); \
    }
Definition at line 77 of file Private.h.
00077 #define ZESTETM1_ERROR_CONN(f, x) \
00078
00079
               if (ZestETM1 ErrorHandler!=NULL) \
              ZestETM1_ErrorHandler(f, (Conn!=NULL ? Conn->CardInfo : NULL), x, ZESTETM1_ERROR_STRING(x)); \
00080
00081
               return (x); \
00082
7.27.2.7 ZESTETM1 ERROR GENERAL
#define ZESTETM1_ERROR_GENERAL(
                 f,
                 x )
Value:
        if (ZestETM1_ErrorHandler!=NULL) \setminus
            {\tt ZestETM1\_ErrorHandler(f, NULL, x, ZESTETM1\_ERROR\_STRING(x));} \setminus \\
        return (x); \
    }
Definition at line 71 of file Private.h.
00071 #define ZESTETM1_ERROR_GENERAL(f, x) \
00072
               if (ZestETM1_ErrorHandler!=NULL) \ ZestETM1_ErrorHandler(f, NULL, x, ZESTETM1_ERROR_STRING(x)); \
00073
00074
              return (x); \
00075
00076
7.27.2.8 ZESTETM1_ERROR_STRING
#define ZESTETM1_ERROR_STRING(
                 X)
Value:
    ZestETM1_ErrorStrings[(x)>=ZESTETM1_ERROR_BASE ? \
                              (x)-ZESTETM1_ERROR_BASE+(ZESTETM1_MAX_INFO-ZESTETM1_INFO_BASE)+(ZESTETM1_MAX_WARNING-ZESTETM1_WARNING
                         ((x)>=ZESTETM1_WARNING_BASE ?
      (x) -ZESTETM1_WARNING_BASE+(ZESTETM1_MAX_INFO-ZESTETM1_INFO_BASE) : (x) -ZESTETM1_INFO_BASE)]
Definition at line 83 of file Private.h.
00083 #define ZESTETM1_ERROR_STRING(x) \
          {\tt ZestETM1\_ErrorStrings[(x)>=ZESTETM1\_ERROR\_BASE~?~\backslash}
00084
00085
      (x) -ZESTETM1_ERROR_BASE+ (ZESTETM1_MAX_INFO-ZESTETM1_INFO_BASE) + (ZESTETM1_MAX_WARNING-ZESTETM1_WARNING_BASE)
                                ((x)>=ZESTETM1_WARNING_BASE ?
00086
      (x) -ZESTETM1_WARNING_BASE+(ZESTETM1_MAX_INFO-ZESTETM1_INFO_BASE) : (x) -ZESTETM1_INFO_BASE)]
```

## 7.27.2.9 ZESTETM1\_RATE\_10MHz

```
#define ZESTETM1_RATE_10MHz (2<<4)</pre>
```

Definition at line 112 of file Private.h.

## 7.27.2.10 ZESTETM1\_RATE\_20MHz

```
#define ZESTETM1_RATE_20MHz (1<<4)</pre>
```

Definition at line 111 of file Private.h.

# 7.27.2.11 ZESTETM1\_RATE\_40MHz

```
#define ZESTETM1_RATE_40MHz (0<<4)</pre>
```

Definition at line 110 of file Private.h.

# 7.27.2.12 ZESTETM1\_REVERSE

Definition at line 116 of file Private.h.

## 7.27.2.13 ZESTETM1\_USER\_DEVICE\_ID

```
#define ZESTETM1_USER_DEVICE_ID (1)
```

Definition at line 113 of file Private.h.

# 7.27.3 Typedef Documentation

#### 7.27.3.1 SOCKET

```
typedef int SOCKET
```

Definition at line 56 of file Private.h.

## 7.27.4 Function Documentation

#### 7.27.4.1 ZestETM1\_CloseConnection()

```
ZESTETM1_STATUS ZestETM1_CloseConnection (
              ZESTETM1_CONNECTION Connection )
Definition at line 225 of file Data.c.
00226 {
00227
          ZESTETM1_CONNECTION_STRUCT *Conn = (ZESTETM1_CONNECTION_STRUCT *)Connection;
00228
00229
          if (Conn==NULL)
00230
          {
              return ZESTETM1_NULL_PARAMETER;
00231
00232
00233
          if (Conn->Magic!=ZESTETM1_CONNECTION_HANDLE_MAGIC)
00234
         {
              return ZESTETM1_ILLEGAL_CONNECTION;
00235
00236
00237
          // Cleanup
00238
          closesocket(Conn->Socket);
00239
00240
          Conn->Magic = 0;
00241
          free (Conn);
00242
00243
          return ZESTETM1_SUCCESS;
00244 }
```

## 7.27.4.2 ZestETM1 EraseFlashSector()

## 7.27.4.3 ZestETM1\_OpenConnection()

## Definition at line 101 of file Data.c.

```
00106 {
          ZESTETM1_CONNECTION_STRUCT *NewStruct;
00108
          SOCKET Socket = -1;
00109
          char AddrBuffer[32];
00110
          char PortBuffer[32];
00111
00112
          if (Connection==NULL || CardInfo==NULL)
00113
         {
00114
              return ZESTETM1_NULL_PARAMETER;
00115
00116
00117
          // Allocate data structure
          NewStruct = malloc(sizeof(ZESTETM1_CONNECTION_STRUCT));
00118
          if (NewStruct==NULL)
00119
00120
00121
              return ZESTETM1_OUT_OF_MEMORY;
00122
```

```
00123
00124
          // Build target addresses
00125
          sprintf(AddrBuffer, "%d.%d.%d.%d", CardInfo->IPAddr[0], CardInfo->IPAddr[1],
00126
                            CardInfo->IPAddr[2], CardInfo->IPAddr[3]);
00127
          sprintf(PortBuffer, "%d", Port);
00128
00129
          if (Type==ZESTETM1_TYPE_UDP)
00130
00131
               // Open UDP connection
00132
              struct sockaddr_in SourceIP;
00133
               int SourceLen = (int)sizeof(struct sockaddr_in);
              Socket = socket(AF_INET, SOCK_DGRAM, IPPROTO_UDP);
00134
00135
              if (Socket<0)</pre>
                   return ZESTETM1_SOCKET_ERROR;
00137
              SourceIP.sin_family = AF_INET;
00138
              SourceIP.sin_addr.s_addr = htonl(INADDR_ANY);
00139
              SourceIP.sin_port = htons(LocalPort);
00140
              bind(Socket, (const struct sockaddr *)&SourceIP, SourceLen);
00141
00142
              NewStruct->Target.sin_family = AF_INET;
              NewStruct->Target.sin_addr.s_addr = inet_addr(AddrBuffer);
NewStruct->Target.sin_port = htons(atoi(PortBuffer));
00143
00144
00145
00146
          else if (Type==ZESTETM1_TYPE_TCP)
00147
00148
               // Open TCP connection
00149
              struct addrinfo *AddrResult = NULL,
00150
                                *Ptr = NULL,
00151
                               Hints:
00152
              int Result:
              struct sockaddr_in SourceIP;
00153
00154
              int SourceLen = (int)sizeof(struct sockaddr_in);
00155
00156
              memset(&Hints, 0, sizeof(Hints));
              Hints.ai_family = AF_UNSPEC;
00157
              Hints.ai_socktype = SOCK_STREAM;
00158
              Hints.ai_protocol = IPPROTO_TCP;
00159
00160
00161
               // Resolve the server address and port
              Result = getaddrinfo(AddrBuffer, PortBuffer, &Hints, &AddrResult);
00162
00163
              if (Result!=0)
00164
               {
00165
                   free (NewStruct);
00166
                   return ZESTETM1_SOCKET_ERROR;
00167
00168
00169
              \ensuremath{//} Attempt to connect to an address until one succeeds
00170
              for (Ptr=AddrResult; Ptr!=NULL; Ptr=Ptr->ai_next)
00171
00172
                   // Create a SOCKET for connecting to server
00173
                   Socket = socket(Ptr->ai_family, Ptr->ai_socktype,
00174
                                    Ptr->ai_protocol);
00175
                   if (Socket<0)
00176
                   {
00177
                       freeaddrinfo(AddrResult);
00178
                       free (NewStruct);
00179
                       return ZESTETM1_SOCKET_ERROR;
00180
00181
00182
                   // Connect to ZestETM1
00183
                   Result = connect(Socket, Ptr->ai_addr, (int)Ptr->ai_addrlen);
00184
                   if (Result<0)</pre>
00185
                   {
00186
                       closesocket(Socket);
                       Socket = -1;
00188
                       continue;
00189
00190
                   break;
00191
              }
00192
00193
              SourceIP.sin_family = AF_INET;
00194
              SourceIP.sin_addr.s_addr = htonl(INADDR_ANY);
00195
              SourceIP.sin_port = 0;
00196
              bind(Socket, (const struct sockaddr *)&SourceIP, SourceLen);
              freeaddrinfo(AddrResult);
00197
00198
00199
          else
00200
          {
00201
              free (NewStruct);
              return ZESTETM1_INVALID_CONNECTION_TYPE;
00202
00203
```

```
00204
00205
          if (Socket==-1)
00206
00207
              free (NewStruct);
00208
              return ZESTETM1_SOCKET_ERROR;
00209
00210
00211
          NewStruct->Magic = ZESTETM1_CONNECTION_HANDLE_MAGIC;
00212
          NewStruct->Type = Type;
00213
          NewStruct->Port = Port;
00214
          NewStruct->LocalPort = LocalPort;
00215
          NewStruct->Socket = Socket;
00216
          NewStruct->CardInfo = CardInfo;
          *Connection = NewStruct;
00218
00219
          return ZESTETM1_SUCCESS;
00220 }
```

# 7.27.4.4 ZestETM1\_SendCommand()

```
ZESTETM1_STATUS ZestETM1_SendCommand (
              ZESTETM1_CARD_INFO * CardInfo,
              ZESTETM1_CONNECTION Connection,
              void * WriteBuffer,
              uint32_t WriteLen,
              void * ReadBuffer,
              uint32 t ReadLen,
              int WaitForAck )
Definition at line 453 of file Data.c.
00458 {
          ZESTETM1_STATUS Result;
00459
00460
          unsigned long Written;
00461
          unsigned long Received;
00462
00463
          // Send/receive data
00464
          Result = ZestETM1_WriteData(Connection, WriteBuffer, WriteLen, &Written,
00465
                                     CardInfo->Timeout);
00466
          if (Result!=ZESTETM1_SUCCESS)
00467
00468
              return Result:
00469
          if (Written!=WriteLen)
00470
00471
              return ZESTETM1_INTERNAL_ERROR;
00472
00473
          *((uint8_t *)ReadBuffer) = 0;
00474
00475
          if (WaitForAck==1)
00476
00477
              Result = ZestETM1_ReadData(Connection, ReadBuffer, ReadLen, &Received,
00478
                                        CardInfo->Timeout);
00479
              if (Result!=ZESTETM1_SUCCESS)
00480
00481
                  return Result;
00482
00483
              if (Received!=ReadLen)
00484
00485
                  return ZESTETM1_INTERNAL_ERROR;
00486
00487
00488
00489
          return ZESTETM1_SUCCESS;
```

# 7.27.4.5 ZestETM1\_SPIReadWrite()

00490 }

```
ZESTETM1_CONNECTION Connection,
               int Device,
               int WordLen,
               uint32_t * WriteData,
               uint32_t * ReadData,
               uint32_t Length,
               int ReleaseCS,
               int WaitForAck )
Definition at line 495 of file Data.c.
00502
          uint8_t Buffer[65536];
00503
          uint32_t *BufPtr;
00504
          uint32_t i;
00505
          ZESTETM1_STATUS Result;
00506
00507
          // Build command
          Buffer[0] = ZESTETM1_COMMAND_SPI;// Command byte
00508
          Buffer[1] = Device;
                                  // SPI device
          Buffer[2] = WordLen;
                                           // Word length
00510
          Buffer[3] = ReleaseCS;
00511
                                           // Release CS on completion
00512
          BufPtr = (uint32_t *) (Buffer+4);
00513
00514
          BufPtr[0] = WriteData==NULL ? 0 : ZESTETM1_REVERSE(Length);
                                                                            // Num words
00515
          BufPtr[1] = ReadData==NULL ? 0 : ZESTETM1_REVERSE(Length);
00516
          if (WriteData!=NULL)
00517
00518
              for (i=0; i<Length; i++)</pre>
00519
00520
                  BufPtr[2+i] = ZESTETM1_REVERSE(WriteData[i]);
00521
00522
00523
          else
              memset (BufPtr+2, 0, Length *4);
00524
00525
00526
          // Send command
         Result = ZestETM1_SendCommand(CardInfo, Connection, Buffer, WriteData==NULL ? 12 : 12+Length*4,
00527
00528
                                        Buffer, ReadData==NULL ? 4 : 4+Length*4,
00529
00530
                                        WaitForAck);
          if (WaitForAck==0)
00531
00532
              return Result;
00533
00534
          // Extract results
00535
          if (Buffer[0]!=ZESTETM1_COMMAND_SPI || Buffer[1]!=0)
00536
00537
              return ZESTETM1_INTERNAL_ERROR;
00538
00539
00540
          BufPtr = (uint32_t *)(Buffer+4);
00541
          if (ReadData!=NULL)
00542
00543
              for (i=0; i<Length; i++)</pre>
00544
00545
                  ReadData[i] = ZESTETM1_REVERSE(BufPtr[i]);
00546
00547
00548
00549
          return ZESTETM1_SUCCESS;
00550 }
```

#### 7.27.4.6 ZestETM1\_WriteFlash()

# 7.27.5 Variable Documentation

#### 7.27.5.1 ZestETM1\_ErrorHandler

```
ZESTETM1_ERROR_FUNC ZestETM1_ErrorHandler [extern]

Definition at line 67 of file Error.c.
```

# 7.27.5.2 ZestETM1\_ErrorStrings

```
char* ZestETM1_ErrorStrings[] [extern]
Definition at line 52 of file Error.c.
          "Success (no error)",
00055
          "Error communicating with socket",
00056
          "An unspecified internal error occurred",
00057
          "Status code is out of range",
00058
          "NULL was used illegally as one of the parameter values",
00059
          "Not enough memory to complete the requested operation",
          "The requested connection type is invalid",
00060
00061
          "The requested connection is invalid",
          "The connection was closed unexpectedly"
00062
00063
          "Operation timed out",
00064
          "One of the parameters has an illegal value",
00065
00066 };
```

# 7.28 Private.h

#### Go to the documentation of this file.

```
00001 // File:
                   Private.h
00002 //
00003 // Purpose:
00004 //
           ZestETM1 Host Library
00005 //
           Internal header file
00006 //
00007 // Version: 1.00
00008 // Date: 11/10/12
00010 // Copyright (C) 2012 Orange Tree Technologies Ltd. All rights reserved.
00011 // Orange Tree Technologies grants the purchaser of a ZestETM1 the right to use and
00012 // modify this source code in any form in designs that target the ZestETM1.
00013 // Orange Tree Technologies prohibits the use of this source code or any modification of
00014 // it in any form in designs that target any other hardware unless the purchaser of the
00015 // ZestETM1 has purchased the appropriate licence from Orange Tree Technologies.
00016 // Contact Orange Tree Technologies if you want to purchase such a licence.
00019 //**
00020 //** Disclaimer: LIMITED WARRANTY AND DISCLAIMER. These designs are
                     provided to you "as is". Orange Tree Technologies and its licensors
00021 //**
00022 //**
                       make and you receive no warranties or conditions, express, implied,
00023 //**
                       statutory or otherwise, and Orange Tree Technologies specifically
00024 //**
                       disclaims any implied warranties of merchantability, non-infringement,
                       or fitness for a particular purpose. Orange Tree Technologies does not
00025 //**
00026 //**
                       warrant that the functions contained in these designs will meet your
00027 //**
                       requirements, or that the operation of these designs will be
00028 //**
                       uninterrupted or error free, or that defects in the Designs will be
                       corrected. Furthermore, Orange Tree Technologies does not warrant or
00029 //**
00030 //**
                       make any representations regarding use or the results of the use of the
00031 //**
                       designs in terms of correctness, accuracy, reliability, or otherwise.
00032 //**
00033 //**
                       LIMITATION OF LIABILITY. In no event will Orange Tree Technologies
00034 //**
                       or its licensors be liable for any loss of data, lost profits, cost or
```

7.28 Private.h 341

```
00035 //**
                       procurement of substitute goods or services, or for any special,
00036 //**
                       incidental, consequential, or indirect damages arising from the use or
00037 //**
                       operation of the designs or accompanying documentation, however caused
00038 //**
                       and on any theory of liability. This limitation will apply even if
00039 //**
                       Orange Tree Technologies has been advised of the possibility of such
00040 //**
                       damage. This limitation shall apply notwithstanding the failure of the
00041 //**
                       essential purpose of any limited remedies herein.
00042 //**
00044
00045 #include <stdint.h>
00047 #if defined(MSVC) || defined(WINGCC)
00048 #include "Winsock2.h"
00049 #include "Ws2tcpip.h"
00050 #else
00051 #include <unistd.h>
00052 #include <sys/socket.h>
00053 #include <netinet/in.h>
00054 #include <netdb.h>
00055 #include <arpa/inet.h>
00056 typedef int SOCKET;
00057 #define SD_BOTH SHUT_RDWR
00058 #define closesocket close
00059 #endif
00060
00061 /********
00062 * Error macro *
00063 ***********
00064 extern ZESTETM1_ERROR_FUNC ZestETM1_ErrorHandler;
00065 #define ZESTETM1_ERROR(f, x) \
00066
00067
             if (ZestETM1_ErrorHandler!=NULL) \
                ZestETM1_ErrorHandler(f, CardInfo, x, ZESTETM1_ERROR_STRING(x)); \
00068
             return (x); \
00069
00070
00071 #define ZESTETM1_ERROR_GENERAL(f, x) \
00072
00073
             if (ZestETM1_ErrorHandler!=NULL) \
                 ZestETM1_ErrorHandler(f, NULL, x, ZESTETM1_ERROR_STRING(x)); \
00074
00075
             return (x); \
00076
00077 #define ZESTETM1_ERROR_CONN(f, x) \
00078
        { \
00079
             if (ZestETM1_ErrorHandler!=NULL) \
00080
             ZestETM1_ErrorHandler(f, (Conn!=NULL ? Conn->CardInfo : NULL), x, ZESTETM1_ERROR_STRING(x)); \
00081
             return (x); \
00082
00083 #define ZESTETM1_ERROR_STRING(x) \
00084
         ZestETM1_ErrorStrings[(x)>=ZESTETM1_ERROR_BASE ? \
00085
      (x) -ZESTETM1_ERROR_BASE+(ZESTETM1_MAX_INFO-ZESTETM1_INFO_BASE)+(ZESTETM1_MAX_WARNING-ZESTETM1_WARNING_BASE)
00086
                             ((x)>=ZESTETM1_WARNING_BASE ?
     (x)-ZESTETM1_WARNING_BASE+(ZESTETM1_MAX_INFO-ZESTETM1_INFO_BASE) : (x)-ZESTETM1_INFO_BASE)]
00087 extern char *ZestETM1_ErrorStrings[];
00088
00089 /*******************
00090 * Network connection descriptor structure
00091 ******************************
00092 #define ZESTETM1_CONNECTION_HANDLE_MAGIC 0xdeadbed1
00093 typedef struct
00094 {
00095
         uint32_t Magic;
         ZESTETM1_CARD_INFO *CardInfo;
00096
         ZESTETM1_CONNECTION_TYPE Type;
00097
00098
        struct sockaddr_in Target;
00099
         uint16_t Port;
        uint16_t LocalPort;
SOCKET Socket;
00100
00101
00102 } ZESTETM1_CONNECTION_STRUCT;
00103
00104 /********
00105 * Constants *
00106 *********/
00107 #define ZESTETM1 DEFAULT TIMEOUT 10000
00108
00109 // SPI Device ID and clock
00110 #define ZESTETM1 RATE 40MHz (0«4)
00111 #define ZESTETM1_RATE_20MHz (1«4)
00112 #define ZESTETM1_RATE_10MHz (2«4)
```

```
00113 #define ZESTETM1_USER_DEVICE_ID
00115 // Reverse bytes in 32 bit word
00116 #define ZESTETM1_REVERSE(x) ((((x)&0xff) < 24) | (((x)&0xff00) < 8) | (((x)&0xff0000) > 8) |
      (((x)&0xff000000)»24))
00117
00119 * Local functions *
00120 *************
00121 ZESTETM1_STATUS ZestETM1_OpenConnection(ZESTETM1_CARD_INFO *CardInfo,
                                               ZESTETM1_CONNECTION_TYPE Type,
00123
                                               uint16_t Port,
00124
                                               uint16_t LocalPort,
                                               ZESTETM1_CONNECTION *Connection);
00125
00126 ZESTETM1_STATUS ZestETM1_CloseConnection(ZESTETM1_CONNECTION Connection);
00127 ZESTETM1_STATUS ZestETM1_SendCommand(ZESTETM1_CARD_INFO *CardInfo,
00128
                                           ZESTETM1_CONNECTION Connection,
00129
                                            void *WriteBuffer, uint32_t WriteLen,
00130
                                            void *ReadBuffer, uint32_t ReadLen,
00131
                                            int WaitForAck);
00132 ZESTETM1_STATUS ZestETM1_SPIReadWrite(ZESTETM1_CARD_INFO *CardInfo,
                                             ZESTETM1_CONNECTION Connection,
00133
00134
                                             int Device,
00135
                                             int WordLen, uint32 t *WriteData,
                                             uint32_t *ReadData, uint32_t Length,
00136
                                             int ReleaseCS, int WaitForAck);
00137
00138 ZESTETM1_STATUS ZestETM1_WriteFlash(ZESTETM1_CARD_INFO *CardInfo,
00139
                                          uint32 t Address,
00140
                                           void *Buffer,
00141
                                          uint32_t Length);
00142 ZESTETM1_STATUS ZestETM1_EraseFlashSector(ZESTETM1_CARD_INFO *CardInfo,
                                                 ZESTETM1_CONNECTION Connection,
00143
00144
                                                 uint32_t Address);
00145
00146
```

# 7.29 C:/Users/hill35/git/camera\_python/nsCamera/comms/ZestETM1/ UPnP.c File Reference

```
#include <stdint.h>
#include <stdib.h>
#include <stdio.h>
#include <memory.h>
#include "ZestETM1.h"
#include "Private.h"
#include <sys/types.h>
#include <sys/ioctl.h>
#include <arpa/inet.h>
#include <ifaddrs.h>
```

# Classes

- struct ZESTETM1 GET SETTINGS CMD
- struct ZESTETM1\_GET\_SETTINGS\_RESPONSE

#### **Macros**

- #define CRT SECURE NO WARNINGS
- #define \_strnicmp strncasecmp
- #define ZESTETM1\_UPNP\_ADDR "239.255.255.250"
- #define ZESTETM1 UPNP PORT 1900
- #define ZESTETM1 COMMAND GET SETTINGS 0xf0

#### **Functions**

- static int ZestETM1\_HTTPGet (char \*IPAddr, char \*Port, char \*FileName, void \*Buffer, int BufferLength, int Wait)
- static ZESTETM1\_STATUS ZestETM1\_ReadSettings (ZESTETM1\_CARD\_INFO \*CardInfo)
- static void ZestETM1\_GetCardInfo (char \*Location, uint32\_t \*NumCards, ZESTETM1\_CARD\_INFO \*\*CardInfo, int Wait)
- static int ZestETM1\_JoinGroup (SOCKET sd, uint32\_t grpaddr, uint32\_t iaddr)
- static int ZestETM1 LeaveGroup (SOCKET sd, uint32 t grpaddr, uint32 t iaddr)
- static ZESTETM1\_STATUS ZestETM1\_GetAllAdapters (uint32\_t \*NumAdapters, struct sockaddr\_in \*\*Adapters)
- ZESTETM1\_STATUS ZestETM1CountCards (unsigned long \*NumCards, ZESTETM1\_CARD\_INFO \*\*CardInfo, unsigned long Wait)
- ZESTETM1\_STATUS ZestETM1FreeCards (ZESTETM1\_CARD\_INFO \*CardInfo)
- ZESTETM1\_STATUS ZestETM1GetCardInfo (ZESTETM1\_CARD\_INFO \*CardInfo)

# **Variables**

static char \* ZestETM1 SearchReq

# 7.29.1 Class Documentation

# 7.29.1.1 struct ZESTETM1\_GET\_SETTINGS\_CMD

Definition at line 82 of file UPnP.c.

#### **Class Members**

uint8_t	Command	
uint8_t	Dummy[3]	

# 7.29.1.2 struct ZESTETM1\_GET\_SETTINGS\_RESPONSE

Definition at line 87 of file UPnP.c.

#### **Class Members**

uint8_t	Command	
uint16_t	ControlPort	
uint8_t	Dummy1[2]	
uint8_t	Dummy2[2]	
uint32_t	Gateway	
uint16_t	HardwareVersion	
uint16_t	HTTPPort	
uint32_t	IPAddr	
uint8_t	MACAddr[6]	
uint32_t	SerialNumber	

# **Class Members**

uint16_t	SoftwareVersion	
uint8_t	Status	
uint32_t	SubNet	

# 7.29.2 Macro Definition Documentation

# 7.29.2.1 \_CRT\_SECURE\_NO\_WARNINGS

#define \_CRT\_SECURE\_NO\_WARNINGS

Definition at line 45 of file UPnP.c.

# 7.29.2.2 \_strnicmp

#define \_strnicmp strncasecmp

Definition at line 61 of file UPnP.c.

# 7.29.2.3 ZESTETM1\_COMMAND\_GET\_SETTINGS

#define ZESTETM1\_COMMAND\_GET\_SETTINGS 0xf0

Definition at line 103 of file UPnP.c.

# 7.29.2.4 ZESTETM1\_UPNP\_ADDR

#define ZESTETM1\_UPNP\_ADDR "239.255.255.250"

Definition at line 78 of file UPnP.c.

# 7.29.2.5 ZESTETM1\_UPNP\_PORT

#define ZESTETM1\_UPNP\_PORT 1900

Definition at line 79 of file UPnP.c.

# 7.29.3 Function Documentation

# 7.29.3.1 ZestETM1\_GetAllAdapters()

```
static ZESTETM1_STATUS ZestETM1_GetAllAdapters (
               uint32_t * NumAdapters,
               struct sockaddr_in ** Adapters ) [static]
Definition at line 431 of file UPnP.c.
00432
00433 #if defined(MSVC) || defined(WINGCC)
          SOCKET Socket;
00434
          SOCKET_ADDRESS_LIST *AddressListPtr;
00435
00436
          DWORD BytesRequired;
00437
          int i;
00438
          int Count = 0;
00439
00440
          *NumAdapters = 0;
00441
          *Adapters = NULL;
00442
00443
          Socket = socket(AF_INET, SOCK_DGRAM, IPPROTO_UDP);
00444
          if (Socket<0)
00445
             return ZESTETM1_INTERNAL_ERROR;
00446
00447
         WSAIoctl(Socket, SIO_ADDRESS_LIST_QUERY, NULL, 0,
          NULL, 0, (LPDWORD)&BytesRequired, NULL, NULL);
AddressListPtr = (SOCKET_ADDRESS_LIST *)malloc(BytesRequired);
00448
00449
00450
          if (AddressListPtr==NULL)
00451
              closesocket(Socket);
00452
00453
              return ZESTETM1_OUT_OF_MEMORY;
00454
00455
         if (WSAIoctl(Socket, SIO_ADDRESS_LIST_QUERY, NULL, 0,
00456
                       AddressListPtr, BytesRequired, &BytesRequired, NULL, NULL) ==SOCKET_ERROR)
00457
00458
              free (AddressListPtr);
00459
              closesocket(Socket);
00460
              return ZESTETM1_INTERNAL_ERROR;
00461
00462
00463
         for (i=0; i<AddressListPtr->iAddressCount; i++)
00464
              if (AddressListPtr->Address[i].iSockaddrLength==sizeof(struct sockaddr_in))
              {
00467
                  Count++;
00468
                  (*Adapters) = (struct sockaddr_in *)realloc(*Adapters, Count*sizeof(struct sockaddr_in));
00469
                  if ((*Adapters) ==NULL)
00470
00471
                      free (AddressListPtr);
00472
                      closesocket(Socket);
00473
                      return ZESTETM1_OUT_OF_MEMORY;
00474
00475
                  memcpy(&(*Adapters)[Count-1], AddressListPtr->Address[i].lpSockaddr, sizeof(struct
      sockaddr in));
00476
00477
00478
          *NumAdapters = Count;
00479
          closesocket (Socket):
00480
          free (AddressListPtr):
00481
00482
          return ZESTETM1_SUCCESS;
00483 #else
         struct ifaddrs *Interfaces;
00484
00485
          struct ifaddrs *Ptr;
00486
          uint32_t Count = 0;
00487
00488
          if (getifaddrs(&Interfaces)!=0)
00489
              return ZESTETM1_INTERNAL_ERROR;
00490
00491
          *NumAdapters = 0;
00492
          *Adapters = NULL;
00493
00494
          Ptr = Interfaces;
00495
          while (Ptr!=NULL)
00496
```

```
00497
              if (Ptr->ifa_addr!=NULL)
00498
00499
00500
                   (*Adapters) = (struct sockaddr_in *)realloc(*Adapters, Count*sizeof(struct sockaddr_in));
00501
                   if ((*Adapters)==NULL)
00502
00503
                      freeifaddrs(Interfaces);
00504
                      return ZESTETM1_OUT_OF_MEMORY;
00505
00506
                  memcpy(&(*Adapters)[Count-1], Ptr->ifa_addr, sizeof(struct sockaddr_in));
00507
00508
              Ptr = Ptr->ifa_next;
00509
00510
          *NumAdapters = Count;
00511
          freeifaddrs (Interfaces);
00512
00513
          return ZESTETM1_SUCCESS;
00514 #endif
00515 }
```

#### 7.29.3.2 ZestETM1 GetCardInfo()

#### Definition at line 279 of file UPnP.c.

```
00283 {
00284
          static char *Prefix = "http://";
00285
          char Buffer[65536];
00286
          char *IPAddr = NULL;
00287
          char *FileName = NULL;
00288
          char *Port = "80";
00289
          int i;
00290
00291
          // Extract address and port
00292
          for (i=0; Location[i]!=0 && Location[i]!='\r' && Location[i]!='\n'; i++)
00293
00294
              if (IPAddr==NULL && Prefix[i]==0) IPAddr = Location+i;
00295
              if (IPAddr==NULL && Location[i]!=Prefix[i])
00296
                   return;
00297
              if (IPAddr!=NULL)
00298
              {
                   if (Location[i]=='/' || Location[i]==0 ||
    Location[i]=='\r' || Location[i]=='\r')
00299
00300
00301
00302
                      Location[i]=0;
00303
                      i++;
break;
00304
00305
                  if (Location[i]==':')
00306
00307
00308
                      00309
                              Location[j]!='\n' && Location[j]!='/') j++;
00310
                      Location[i]=0;
00311
00312
                      Location[j]=0;
00313
                      Port = Location+i+1;
                      i = j+1;
00314
                      break;
00315
00316
                  }
00317
              }
00318
          if (IPAddr==NULL || Location[i]==0 || Location[i]=='\r' || Location[i]=='\n')
00319
00320
              return;
00321
00322
          // Extract XML filename
00323
          FileName = Location+i;
          while (Location[i]!=0 && Location[i]!='\r' && Location[i]!='\n') i++;
00324
00325
          Location[i] = 0;
00326
```

```
00327
          // Get XML file
00328
          memset(Buffer, 0, sizeof(Buffer));
00329
          if (ZestETM1_HTTPGet(IPAddr, Port, FileName,
00330
                               Buffer, sizeof(Buffer), Wait)>0)
00331
          {
00332
              uint8_t IPAddr[4];
00333
              uint32_t i;
00334
              uint16_t Port;
00335
00336
              // Parse XML for service description URL
00337
              char *ControlURL = strstr(Buffer, "<controlURL>");
00338
              if (ControlURL==NULL)
00339
                  return;
00340
00341
              ControlURL += 12;
00342
              for (i=0; i<4; i++)</pre>
00343
00344
                   char *EndPtr;
00345
                   IPAddr[i] = (uint8_t) strtoul(ControlURL, &EndPtr, 10);
                  if ((i!=3 && *EndPtr!='.') || (i==3 && *EndPtr!=':'))
00346
00347
                      break;
00348
                  ControlURL = EndPtr+1;
00349
00350
              if (i!=4)
00351
                  return:
00352
00353
              Port = atoi(ControlURL);
00354
00355
              // Make sure only unique devices are added to the list!
00356
              if (*CardInfo!=NULL)
00357
00358
                   for (i=0; i<*NumCards; i++)</pre>
00359
                       if ((*CardInfo)[i].ControlPort==Port &&
00360
                           \texttt{memcmp((*CardInfo)[i].IPAddr, IPAddr, sizeof(IPAddr))==0)}
00361
00362
00363
                           break;
00364
00365
                   }
00366
00367
              if (*CardInfo==NULL || i==*NumCards)
00368
00369
                  uint32_t Index = (*NumCards);
                  ZESTETM1_CARD_INFO *NewBuffer;
00370
00371
00372
                   // Allocate space for new card info structure
00373
                   if ((*CardInfo) ==NULL)
00374
                      NewBuffer = malloc(sizeof(ZESTETM1_CARD_INFO)*(Index+1));
00375
00376
                      NewBuffer = realloc(*CardInfo, sizeof(ZESTETM1_CARD_INFO)*(Index+1));
00377
00378
                   // Get new card settings
00379
                  NewBuffer[Index].ControlPort = Port;
00380
                  memcpy(NewBuffer[Index].IPAddr, IPAddr, sizeof(IPAddr));
00381
                  NewBuffer[Index].Timeout = Wait;
00382
                   if (ZestETM1_ReadSettings(&(NewBuffer[Index]))!=ZESTETM1_SUCCESS)
00383
00384
                       if (*CardInfo==NULL)
00385
                           free (NewBuffer);
00386
                       else
00387
00388
                           (*CardInfo) = NewBuffer;
00389
                           memset(NewBuffer[Index].IPAddr, 0, sizeof(NewBuffer[Index].IPAddr));
00390
                           NewBuffer[Index].ControlPort = 0;
00391
00392
00393
                  else
00394
                   {
00395
                      NewBuffer[Index].Timeout = ZESTETM1_DEFAULT_TIMEOUT;
00396
                       (*CardInfo) = NewBuffer;
00397
                       (*NumCards)++;
00398
                  }
00399
              }
00400
          }
00401 }
```

# 7.29.3.3 ZestETM1\_HTTPGet()

```
static int ZestETM1_HTTPGet (
             char * IPAddr,
             char * Port,
             char * FileName.
             void * Buffer,
             int BufferLength,
             int Wait ) [static]
Definition at line 108 of file UPnP.c.
```

```
00111 {
          struct addrinfo *Addr = NULL,
00112
                           *Ptr = NULL,
00113
00114
                          Hints;
          int Result;
00115
          SOCKET Socket;
00116
00117
          int Offset = 0;
          char Req[1024];
00118
00119
          char *HdrEnd;
00120
          fd_set ReadFDS;
00121
          struct timeval Timeout;
00122
00123
          Timeout.tv_sec = Wait/1000;
00124
          Timeout.tv_usec = (Wait%1000) *1000;
00125
00126
          \ensuremath{//} Attempt to connect to the address
00127
          memset(&Hints, 0, sizeof(Hints));
00128
          Hints.ai_family = AF_UNSPEC;
00129
          Hints.ai_socktype = SOCK_STREAM;
00130
          Hints.ai_protocol = IPPROTO_TCP;
00131
00132
          // Resolve the server address and port
00133
          Result = getaddrinfo(IPAddr, Port, &Hints, &Addr);
00134
          if (Result!=0)
00135
              return -1;
00136
00137
          for (Ptr=Addr; Ptr!=NULL; Ptr=Ptr->ai_next)
00138
00139
               // Create a SOCKET for connecting to server
00140
              Socket = socket(Ptr->ai_family, Ptr->ai_socktype,
00141
                               Ptr->ai_protocol);
00142
               if (Socket<0)
00143
              {
00144
                  freeaddrinfo(Addr);
00145
                  return -1;
00146
00147
00148
              // Connect to ZestETM1
00149
              Result = connect(Socket, Ptr->ai_addr, (int)Ptr->ai_addrlen);
00150
              if (Result<0)
00151
                  closesocket(Socket);
00152
00153
                  Socket = -1;
00154
                  continue;
00155
00156
              break:
00157
          freeaddrinfo(Addr);
00158
00159
          // Send GET request
00160
          sprintf(Req, "GET /%s HTTP/1.1\r\nHOST: %s:%s\r\nContent-length: 0\r\n\r\n", FileName, IPAddr, Port);
00161
00162
          Result = send(Socket, Req, (int)strlen(Req), 0);
          if (Result!=strlen(Req))
00163
00164
          {
00165
              closesocket(Socket);
              return -1;
00166
00167
          }
00168
00169
          // Get response
00170
          ((char*)Buffer)[0] = 0;
00171
00172
          {
00173
              FD_ZERO(&ReadFDS);
00174
              FD_SET(Socket, &ReadFDS);
```

```
00175
              Result = select((int)Socket+1, &ReadFDS, NULL, NULL, &Timeout);
00176
              if (Result<0)</pre>
00177
              {
00178
                  closesocket(Socket);
00179
                 return -1;
00180
00181
              if (!FD_ISSET(Socket, &ReadFDS)) break;
00182
              Result = recv(Socket, (char *)Buffer+Offset, BufferLength-Offset, 0);
00183
              if (Result<0)</pre>
00184
             {
00185
                 closesocket(Socket);
00186
                 return -1;
00187
             Offset+=Result;
00188
00189
              if (Offset==BufferLength)
00190
00191
                 break:
00192
00193
         } while (Result!=0);
00194
00195
         // Check status response
         if (_strnicmp("HTTP/1.1 200 OK", Buffer, 15)!=0)
00196
00197
         {
00198
              closesocket(Socket);
00199
              return -1;
00200
00201
00202
         // Remove HTTP header
         00203
00204
         if (HdrEnd==NULL)
00205
00206
              Offset = 0;
00207
00208
         else
00209
        {
              Offset -= (int) (HdrEnd+4-(char*)Buffer);
00210
00211
              memcpy(Buffer, HdrEnd+4, Offset);
00212
00213
00214
         closesocket (Socket);
00215
         return Offset;
00216 }
```

# 7.29.3.4 ZestETM1\_JoinGroup()

# Definition at line 407 of file UPnP.c.

#### 7.29.3.5 ZestETM1 LeaveGroup()

# Definition at line 417 of file UPnP.c.

# 7.29.3.6 ZestETM1\_ReadSettings()

```
static ZESTETM1_STATUS ZestETM1_ReadSettings (
               ZESTETM1_CARD_INFO * CardInfo ) [static]
Definition at line 221 of file UPnP.c.
00222 {
00223
          ZESTETM1 STATUS Result:
00224
          ZESTETM1_CONNECTION Connection;
          ZESTETM1_GET_SETTINGS_CMD Cmd;
00225
00226
          ZESTETM1_GET_SETTINGS_RESPONSE Response;
00227
00228
          if (CardInfo==NULL)
00229
          {
00230
              return ZESTETM1 NULL PARAMETER:
00231
00232
00233
          // Open control connection
          Result = ZestETM1_OpenConnection(CardInfo, ZESTETM1_TYPE_TCP,
00234
00235
                                            CardInfo->ControlPort, 0, &Connection);
00236
          if (Result!=ZESTETM1_SUCCESS)
00237
00238
              return Result:
00239
00240
00241
          // Get the settings from the device
00242
          Cmd.Command = ZESTETM1_COMMAND_GET_SETTINGS;
00243
          Result = ZestETM1_SendCommand(CardInfo, Connection,
00244
                                        &Cmd, sizeof(Cmd),
00245
                                         &Response, sizeof(Response), 1);
          if (Result!=ZESTETM1_SUCCESS)
00246
00247
00248
              ZestETM1_CloseConnection(Connection);
00249
              return Result;
00250
00251
00252
          // Result values (including endian conversion)
00253
          CardInfo->FirmwareVersion = ((Response.SoftwareVersion»8)&0xff)|((Response.SoftwareVersion&0xff) < 8);
00254
           {\tt CardInfo-> Hardware Version = ((Response. Hardware Version * 8) \& 0xff) | ((Response. Hardware Version \& 0xff) ( 8); } \\
00255
          CardInfo->Gateway[0] = (uint8_t)((Response.Gateway>0)&0xff);
00256
          CardInfo->Gateway[1] = (uint8_t) ((Response.Gateway»8)&0xff);
00257
          CardInfo->Gateway[2] = (uint8_t)((Response.Gateway»16)&0xff);
00258
          CardInfo->Gateway[3] = (uint8_t)((Response.Gateway»24)&0xff);
          CardInfo->SubNet[0] = (uint8_t) ((Response.SubNet»0)&0xff);
00259
          CardInfo->SubNet[1] = (uint8_t) ((Response.SubNet>8)&0xff);
00260
          CardInfo->SubNet[2] = (uint8_t) ((Response.SubNet**)16) &0xff);
00261
          CardInfo->SubNet[3] = (uint8_t) ((Response.SubNet»24)&0xff);
00262
          memcpy(CardInfo->MACAddr, Response.MACAddr, 6);
00263
00264
          CardInfo->ControlPort = ((Response.ControlPort*8)&Oxff) | ((Response.ControlPort&Oxff) < 8);</pre>
          CardInfo->HTTPPort = ((Response.HTTPPort%0xff) | ((Response.HTTPPort&0xff) «8);
00265
00266
          CardInfo->SerialNumber = ZESTETM1_REVERSE(Response.SerialNumber);
00267
00268
          Result = ZestETM1_CloseConnection(Connection);
00269
          if (Result!=ZESTETM1_SUCCESS)
00270
              return Result;
00271
00272
          return ZESTETM1 SUCCESS;
```

# 7.29.3.7 ZestETM1CountCards()

00273 }

```
ZESTETM1_CARD_INFO ** CardInfo,
               unsigned long Wait )
Definition at line 521 of file UPnP.c.
00524
          SOCKET Socket;
00525
          struct sockaddr_in DestIP;
          struct sockaddr_in SourceIP;
00527
          int SourceIPLength;
          int Flag = 1;
00529
          int Result;
         char Req[1024];
00530
00531
          char Response[1024];
          ZESTETM1_CARD_INFO *Cards = NULL;
00533
          uint32_t CardCount = 0;
00534
          int i;
00535
          struct timeval Timeout;
00536
         fd_set ReadFDS;
00537
          uint32 t Interface;
00538
         uint32_t NumAdapters;
         struct sockaddr_in *Adapters;
ZESTETM1_STATUS Status;
00539
00540
00541
00542
          *NumCards = 0:
00543
          *CardInfo = NULL;
00544
00545
          \ensuremath{//} Get a list of all adapters
          Status = ZestETM1_GetAllAdapters(&NumAdapters, &Adapters);
00546
          if (Status!=ZESTETM1_SUCCESS)
00547
              ZESTETM1_ERROR_GENERAL("ZestETM1CountCards", Status);
00548
00549
          if (NumAdapters==0)
00550
         {
00551
              *NumCards = 0:
              return ZESTETM1_SUCCESS;
00552
00553
00554
          // Send queries on all interfaces
00555
00556
          for (Interface=0; Interface<NumAdapters; Interface++)</pre>
00557
00558
              if (Adapters[Interface].sin_family!=AF_INET)
00559
                  continue;
00560
00561
              \ensuremath{//} Open socket for search requests
00562
              Socket = socket(AF_INET, SOCK_DGRAM, IPPROTO_UDP);
00563
              if (Socket<0)</pre>
00564
00565
                   continue;
00566
00567
00568
              \ensuremath{//} Set reuse port to on to allow multiple binds per host
00569
             if (setsockopt(Socket, SOL_SOCKET, SO_REUSEADDR, (char *)&Flag,
00570
                              sizeof(Flag))<0)
00571
00572
                  closesocket(Socket);
00573
                  continue:
00574
00575
00576
              // Bind to port for receiving responses
00577
              SourceIP.sin_family = AF_INET;
00578
              SourceIP.sin_addr.s_addr = htonl(INADDR_ANY);
00579
              SourceIP.sin_port = 0;
00580
              Result = bind(Socket, (const struct sockaddr *)(&Adapters[Interface]), sizeof(struct
     sockaddr_in));
00581
         if (Result<0)
00582
              {
00583
                  closesocket(Socket);
00584
                  continue;
00585
00586
00587
              // Join multicast group
00588
              if (ZestETM1_JoinGroup(Socket, inet_addr(ZESTETM1_UPNP_ADDR),
00589
                                     htonl(INADDR ANY))<0)
00590
00591
                  closesocket(Socket);
00592
                  continue;
00593
              }
00594
              // Send M-SEARCH request
00595
```

// Send more than once as UDP is unreliable

DestIP.sin\_family = AF\_INET;

00596

```
00598
              DestIP.sin_addr.s_addr = inet_addr(ZESTETM1_UPNP_ADDR);
00599
              DestIP.sin_port = htons(1900);
00600
              sprintf(Req, ZestETM1_SearchReq, (Wait+999)/1000);
00601
               for (i=0; i<3; i++)</pre>
00602
               {
00603
                   Result = sendto(Socket, ZestETM1_SearchReq,
00604
                                   (int) strlen (ZestETM1_SearchReq),
00605
                                   0, (struct sockaddr *)&DestIP, sizeof(DestIP));
00606
                   if (Result!=strlen(ZestETM1_SearchReq))
00607
                   {
00608
                       ZestETM1_LeaveGroup(Socket, inet_addr(ZESTETM1_UPNP_ADDR),
00609
                                          htonl(INADDR_ANY));
00610
                       closesocket(Socket);
00611
                       continue;
00612
                   }
00613
              }
00614
00615
               // Read responses
00616
              Timeout.tv_sec = Wait/1000;
              Timeout.tv_usec = (Wait%1000) *1000;
00617
00618
00619
              {
00620
                  FD_ZERO(&ReadFDS);
00621
                  FD_SET(Socket, &ReadFDS);
00622
                  Result = select((int)Socket+1, &ReadFDS, NULL, NULL, &Timeout);
00623
                   if (Result<0)
00624
00625
                       break:
00626
00627
                   if (!FD_ISSET(Socket, &ReadFDS)) break;
00628
                  SourceIPLength = sizeof(SourceIP);
00629
00630
                  Result = recvfrom(Socket, Response, sizeof(Response),
00631
                                     0, (struct sockaddr *) & SourceIP, & SourceIPLength);
00632
                   if (Result<0)
00633
                   {
00634
                       // Error!
00635
                       break;
00636
                  else if (Result==0)
00637
00638
00639
                       // Clean shut down
00640
                       break;
00641
                   }
00642
                  else
00643
00644
                       // Parse results
                       if (_strnicmp("NOTIFY", Response, 6) == 0 ||
00645
00646
                           _strnicmp("HTTP/1.1 200 OK", Response, 15) == 0)
00647
00648
                           // Check its a {\tt GigExpedite} and request XML description
00649
                           char *Server = strstr(Response, "GigExpedite2");
00650
                           char *Location = strstr(Response, "LOCATION");
00651
                           if (Server!=NULL && Location!=NULL)
00652
00653
00654
                               for (c=Location+10; *c!=0 && *c!='\n' && *c!='\r'; c++);
00655
00656
                               ZestETM1_GetCardInfo(Location+10, &CardCount, &Cards, Wait);
00657
00658
00659
00660
              } while(1);
00661
00662
               // Leave multicast group
              ZestETM1_LeaveGroup(Socket, inet_addr(ZESTETM1_UPNP_ADDR),
00663
                                  htonl(INADDR_ANY));
00664
00665
00666
              // Close socket
00667
              closesocket(Socket);
00668
00669
00670
          *NumCards = CardCount;
          *CardInfo = Cards;
00671
00672
          free (Adapters):
00673
00674
          return ZESTETM1 SUCCESS;
00675 }
```

# 7.29.3.8 ZestETM1FreeCards()

#### 7.29.3.9 ZestETM1GetCardInfo()

```
{\tt ZESTETM1\_STATUS} \  \  {\tt ZESTETM1\_CARD\_INFO} \  \  * \  \  {\tt CardInfo} \  \  )
```

# Definition at line 691 of file UPnP.c.

```
00692 {
00693     ZESTETM1_STATUS Status;
00694
00695     Status = ZestETM1_ReadSettings(CardInfo);
00696     if (Status!=ZESTETM1_SUCCESS)
00697     {
          ZESTETM1_ERROR("ZestETM1GetCardInfo", Status);
00699     }
00700     return ZESTETM1_SUCCESS;
00702 }
```

# 7.29.4 Variable Documentation

# 7.29.4.1 ZestETM1\_SearchReq

```
char* ZestETM1_SearchReq [static]
```

#### Initial value:

```
=
"M-SEARCH * HTTP/1.1\r\n"
"ST: upnp:rootdevice\r\n"
"MX: %d\r\n"
"MAN: \"ssdp:discover\"\r\n"
"HOST: 239.255.255.250:1900\r\n"
```

Definition at line 69 of file UPnP.c.

# 7.30 **UPnP.c**

00076

#### Go to the documentation of this file. 00001 // File: 00002 // 00003 // Purpose: ZestETM1 Host Library 00004 // 00005 // UPnP board discovery functions 00006 // 00007 // Version: 1.00 00008 // Date: 11/10/12 00010 // Copyright (C) 2012 Orange Tree Technologies Ltd. All rights reserved. 00011 // Orange Tree Technologies grants the purchaser of a ZestETM1 the right to use and 00012 // modify this source code in any form in designs that target the ZestETM1. 00013 // Orange Tree Technologies prohibits the use of this source code or any modification of 00014 // it in any form in designs that target any other hardware unless the purchaser of the 00015 // ZestETM1 has purchased the appropriate licence from Orange Tree Technologies. 00016 // Contact Orange Tree Technologies if you want to purchase such a licence. 00017 00018 //\* 00019 //\*\* 00020 //\*\* Disclaimer: LIMITED WARRANTY AND DISCLAIMER. These designs are provided to you "as is". Orange Tree Technologies and its licensors 00021 //\*\* 00022 //\*\* make and you receive no warranties or conditions, express, implied, 00023 //\*\* statutory or otherwise, and Orange Tree Technologies specifically 00024 //\*\* disclaims any implied warranties of merchantability, non-infringement, or fitness for a particular purpose. Orange Tree Technologies does not 00025 //\*\* 00026 //\*\* warrant that the functions contained in these designs will meet your 00027 //\*\* requirements, or that the operation of these designs will be 00028 //\*\* uninterrupted or error free, or that defects in the Designs will be 00029 //\*\* corrected. Furthermore, Orange Tree Technologies does not warrant or 00030 //\*\* make any representations regarding use or the results of the use of the 00031 //\*\* designs in terms of correctness, accuracy, reliability, or otherwise. 00032 //\*\* 00033 //\*\* LIMITATION OF LIABILITY. In no event will Orange Tree Technologies 00034 //\*\* or its licensors be liable for any loss of data, lost profits, cost or 00035 //\*\* procurement of substitute goods or services, or for any special, 00036 //\*\* incidental, consequential, or indirect damages arising from the use or 00037 //\*\* operation of the designs or accompanying documentation, however caused 00038 //\*\* and on any theory of liability. This limitation will apply even if 00039 //\*\* Orange Tree Technologies has been advised of the possibility of such damage. This limitation shall apply notwithstanding the failure of the 00040 //\*\* 00041 //\*\* essential purpose of any limited remedies herein. 00042 //\*\* 00045 #define \_CRT\_SECURE\_NO\_WARNINGS 00046 #ifdef WINGCC 00047 #define \_\_USE\_W32\_SOCKETS 00048 #endif 00049 00050 #include <stdint.h> 00051 #include <stdlib.h> 00052 #include <stdio.h> 00053 #include <memory.h> 00054 #include "ZestETM1.h" 00055 #include "Private.h" 00056 00057 #ifdef WINGCC 00058 #define \_strnicmp strncasecmp 00059 #endif 00060 #if !defined(MSVC) && !defined(WINGCC) 00061 #define \_strnicmp strncasecmp 00062 #include <sys/types.h> 00063 #include <sys/ioctl.h> 00064 #include <arpa/inet.h> 00065 #include <ifaddrs.h> 00066 #endif 00068 // UPnP search string 00069 static char \*ZestETM1\_SearchReg = 00070 "M-SEARCH \* HTTP/1.1\r\n" 00071 "ST: upnp:rootdevice\r\n" 00072 "MX: %d\r\n" 00073 "MAN: \"ssdp:discover\"\r\n" 00074 "HOST: 239.255.255.250:1900\r\n" 00075;

7.30 UPnP.c 355

```
00077 // UPnP broadcast address and port
00078 #define ZESTETM1_UPNP_ADDR "239.255.255.250"
00079 #define ZESTETM1_UPNP_PORT 1900
08000
00081 // Get settings command structure
00082 typedef struct
00083 {
00084
          uint8_t Command;
00085
          uint8_t Dummy[3];
00086 } ZESTETM1_GET_SETTINGS_CMD;
00087 typedef struct
00089
          uint8_t Command;
00090
          uint8_t Status;
00091
          uint8_t Dummy1[2];
00092
          uint16_t SoftwareVersion;
00093
          uint16_t HardwareVersion;
00094
          uint32_t SerialNumber;
          uint32_t IPAddr;
00095
          uint32_t Gateway;
00096
00097
          uint32_t SubNet;
00098
          uint16_t HTTPPort;
          uint16_t ControlPort;
00099
00100
          uint8_t MACAddr[6];
00101    uint8_t Dummy2[2];
00102 } ZESTETM1_GET_SETTINGS_RESPONSE;
00103 #define ZESTETM1_COMMAND_GET_SETTINGS 0xf0
00104
00105 /*************
00106 * Issue HTTP GET command *
0.0107 ******************
00108 static int ZestETM1_HTTPGet(char *IPAddr, char *Port,
00109
                                  char *FileName, void *Buffer,
00110
                                  int BufferLength, int Wait)
00111 {
          struct addrinfo *Addr = NULL,
00112
00113
                           *Ptr = NULL,
00114
                           Hints:
00115
          int Result;
00116
          SOCKET Socket:
00117
          int Offset = 0;
00118
          char Req[1024];
00119
          char *HdrEnd;
00120
          fd set ReadFDS;
00121
          struct timeval Timeout;
00122
00123
          Timeout.tv_sec = Wait/1000;
00124
          Timeout.tv\_usec = (Wait%1000) *1000;
00125
00126
          // Attempt to connect to the address
00127
          memset(&Hints, 0, sizeof(Hints));
00128
          Hints.ai_family = AF_UNSPEC;
          Hints.ai_socktype = SOCK_STREAM;
Hints.ai_protocol = IPPROTO_TCP;
00129
00130
00131
00132
          // Resolve the server address and port
00133
          Result = getaddrinfo(IPAddr, Port, &Hints, &Addr);
00134
          if (Result!=0)
00135
00136
00137
          for (Ptr=Addr; Ptr!=NULL; Ptr=Ptr->ai_next)
00138
00139
               // Create a SOCKET for connecting to server
00140
              Socket = socket(Ptr->ai_family, Ptr->ai_socktype,
00141
                               Ptr->ai_protocol);
00142
               if (Socket<0)
00143
              {
00144
                   freeaddrinfo(Addr);
00145
                   return -1;
00146
00147
               // Connect to ZestETM1
00148
00149
              Result = connect(Socket, Ptr->ai_addr, (int)Ptr->ai_addrlen);
00150
               if (Result<0)
00151
00152
                   closesocket (Socket);
                   Socket = -1;
00153
00154
                   continue:
00155
00156
              break:
00157
```

```
00158
          freeaddrinfo(Addr);
00159
          // Send GET request
00160
00161
          sprintf(Req, "GET / \$s \ HTTP/1.1\ r\ hHOST: \$s:\$s\ r\ nContent-length: 0\ r\ n', r', n'', FileName, IPAddr, Port);
00162
          Result = send(Socket, Req, (int)strlen(Req), 0);
00163
          if (Result!=strlen(Req))
00164
          {
00165
              closesocket(Socket);
00166
              return -1;
00167
00168
00169
          // Get response
00170
          ((char*)Buffer)[0] = 0;
00171
00172
          {
00173
              FD_ZERO(&ReadFDS);
00174
              FD_SET(Socket, &ReadFDS);
00175
              Result = select((int)Socket+1, &ReadFDS, NULL, NULL, &Timeout);
00176
              if (Result<0)
00177
              {
00178
                  closesocket(Socket);
00179
                  return -1;
00180
00181
              if (!FD_ISSET(Socket, &ReadFDS)) break;
00182
              Result = recv(Socket, (char *)Buffer+Offset, BufferLength-Offset, 0);
00183
              if (Result<0)
00184
              {
00185
                  closesocket(Socket);
00186
                  return -1;
00187
00188
              Offset += Result:
00189
              if (Offset==BufferLength)
00190
00191
                  break;
00192
00193
          } while (Result!=0);
00194
00195
          // Check status response
          if (_strnicmp("HTTP/1.1 200 OK", Buffer, 15)!=0)
00196
00197
00198
              closesocket(Socket);
00199
              return -1;
00200
00201
00202
          00203
          HdrEnd = strstr(Buffer, "\r\n\r\n");
00204
          if (HdrEnd==NULL)
00205
00206
              Offset = 0;
00207
00208
         else
00209
          {
00210
              Offset -= (int) (HdrEnd+4-(char*)Buffer);
00211
              memcpy(Buffer, HdrEnd+4, Offset);
00212
00213
00214
          closesocket(Socket);
00215
          return Offset;
00216 }
00217
00218 /****************
00219 * Read settings from ETM1 flash *
00221 static ZESTETM1_STATUS ZestETM1_ReadSettings(ZESTETM1_CARD_INFO *CardInfo)
00222 {
00223
          ZESTETM1_STATUS Result;
          ZESTETM1_CONNECTION Connection;
          ZESTETM1_GET_SETTINGS_CMD Cmd;
00225
          ZESTETM1_GET_SETTINGS_RESPONSE Response;
00226
00227
00228
          if (CardInfo==NULL)
00229
         {
00230
              return ZESTETM1_NULL_PARAMETER;
00231
00232
00233
          // Open control connection
00234
          Result = ZestETM1_OpenConnection(CardInfo, ZESTETM1_TYPE_TCP,
00235
                                           CardInfo->ControlPort, 0, &Connection);
00236
          if (Result!=ZESTETM1_SUCCESS)
00237
00238
              return Result:
```

7.30 UPnP.c 357

```
00239
00240
00241
            // Get the settings from the device
00242
           Cmd.Command = ZESTETM1_COMMAND_GET_SETTINGS;
00243
           Result = ZestETM1_SendCommand(CardInfo, Connection,
00244
                                              &Cmd, sizeof(Cmd),
00245
                                              &Response, sizeof(Response), 1);
00246
           if (Result!=ZESTETM1_SUCCESS)
00247
           {
00248
                ZestETM1_CloseConnection(Connection);
00249
                return Result;
00250
00251
00252
            // Result values (including endian conversion)
           CardInfo->FirmwareVersion = ((Response.SoftwareVersion»8)&0xff)|((Response.SoftwareVersion&0xff) <8);
CardInfo->HardwareVersion = ((Response.HardwareVersion»8)&0xff)|((Response.HardwareVersion&0xff) <8);
00253
00254
00255
           CardInfo->Gateway[0] = (uint8_t)((Response.Gateway>0)&0xff);
           CardInfo->Gateway[1] = (uint8_t) ((Response.Gateway>8)&0xff);
00256
00257
           CardInfo->Gateway[2] = (uint8_t)((Response.Gateway»16)&0xff);
           CardInfo->Gateway[3] = (uint8_t)((Response.Gateway>24)&0xff);
00258
           CardInfo->SubNet[0] = (uint8_t)((Response.SubNet**)0 & 0 x ff);
CardInfo->SubNet[1] = (uint8_t)((Response.SubNet**)8) & 0 x ff);
00259
00260
           CardInfo->SubNet[2] = (uint8_t) ((Response.SubNet*)6) &Oxff);
CardInfo->SubNet[3] = (uint8_t) ((Response.SubNet*)24) &Oxff);
00261
00262
00263
           memcpy(CardInfo->MACAddr, Response.MACAddr, 6);
CardInfo->ControlPort = ((Response.ControlPort&0xff) ((Response.ControlPort&0xff) (%);
00264
           CardInfo->HTTPPort = ((Response.HTTPPortx8)&Oxff)|((Response.HTTPPort&Oxff) x8);
CardInfo->SerialNumber = ZESTETM1_REVERSE(Response.SerialNumber);
00265
00266
00267
           Result = ZestETM1_CloseConnection(Connection);
00268
00269
           if (Result!=ZESTETM1_SUCCESS)
00270
                return Result:
00271
           return ZESTETM1_SUCCESS;
00272
00273 }
00274
00275
00276 /*****************
00277 * Read information about a card *
00278 *******************/
00279 static void ZestETM1_GetCardInfo(char *Location,
00280
                                             uint32_t *NumCards,
00281
                                             ZESTETM1_CARD_INFO **CardInfo,
00282
                                             int Wait)
00283 {
00284
           static char *Prefix = "http://";
00285
           char Buffer[65536];
00286
           char *IPAddr = NULL;
00287
           char *FileName = NULL;
00288
           char *Port = "80";
00289
           int i;
00290
00291
           // Extract address and port
00292
           for (i=0; Location[i]!=0 && Location[i]!='\r' && Location[i]!='\r' i++)
00293
00294
                 if (IPAddr==NULL && Prefix[i]==0) IPAddr = Location+i;
00295
                if (IPAddr==NULL && Location[i]!=Prefix[i])
00296
                     return;
00297
                if (IPAddr!=NULL)
00298
                {
                     if (Location[i]=='/' || Location[i]==0 ||
    Location[i]=='\r' || Location[i]=='\r')
00299
00300
00301
00302
                          Location[i]=0;
00303
                          i++;
00304
                          break;
00305
00306
                     if (Location[i]==':')
00307
00308
                          int j=i;
00309
                          while (Location[j]!=0 && Location[j]!='\r' &&
00310
                                  Location[j]!='\n' && Location[j]!='/') j++;
00311
                          Location[i]=0;
00312
                          Location[i]=0;
00313
                          Port = Location+i+1;
00314
                          i = j+1;
00315
                          break:
00316
                     }
00317
00318
00319
           if (IPAddr==NULL || Location[i]==0 || Location[i]=='\r' || Location[i]=='\n')
```

```
00320
              return;
00321
00322
          // Extract XML filename
00323
          FileName = Location+i;
00324
          while (Location[i]!=0 && Location[i]!='\r' && Location[i]!='\n') i++;
00325
          Location[i] = 0;
00326
00327
          // Get XML file
00328
          memset(Buffer, 0, sizeof(Buffer));
          if (ZestETM1_HTTPGet(IPAddr, Port, FileName,
00329
00330
                               Buffer, sizeof(Buffer), Wait)>0)
00331
00332
              uint8_t IPAddr[4];
              uint32_t i;
00333
00334
              uint16_t Port;
00335
00336
              // Parse XML for service description URL
              char *ControlURL = strstr(Buffer, "<controlURL>");
00337
00338
              if (ControlURL==NULL)
00339
                  return;
00340
00341
              ControlURL += 12;
              for (i=0; i<4; i++)</pre>
00342
00343
00344
                  char *EndPtr:
                  IPAddr[i] = (uint8 t)strtoul(ControlURL, &EndPtr, 10);
00345
                  if ((i!=3 && *EndPtr!='.') || (i==3 && *EndPtr!=':'))
00346
00347
                      break;
00348
                  ControlURL = EndPtr+1;
00349
00350
              if (i!=4)
00351
                  return:
00352
00353
              Port = atoi(ControlURL);
00354
              // Make sure only unique devices are added to the list!
00355
00356
              if (*CardInfo!=NULL)
00357
00358
                   for (i=0; i<*NumCards; i++)</pre>
00359
00360
                       if ((*CardInfo)[i].ControlPort==Port &&
00361
                           memcmp((*CardInfo)[i].IPAddr, IPAddr, sizeof(IPAddr))==0)
00362
00363
                           break;
00364
00365
00366
00367
              if (*CardInfo==NULL || i==*NumCards)
00368
00369
                  uint32_t Index = (*NumCards);
00370
                  ZESTETM1_CARD_INFO *NewBuffer;
00371
00372
                   // Allocate space for new card info structure
00373
                   if ((*CardInfo) ==NULL)
00374
                      NewBuffer = malloc(sizeof(ZESTETM1_CARD_INFO)*(Index+1));
00375
00376
                      NewBuffer = realloc(*CardInfo, sizeof(ZESTETM1_CARD_INFO)*(Index+1));
00377
00378
                   // Get new card settings
00379
                  NewBuffer[Index].ControlPort = Port;
00380
                  memcpy(NewBuffer[Index].IPAddr, IPAddr, sizeof(IPAddr));
00381
                  NewBuffer[Index].Timeout = Wait;
00382
                  if (ZestETM1_ReadSettings(&(NewBuffer[Index]))!=ZESTETM1_SUCCESS)
00383
                  {
00384
                       if (*CardInfo==NULL)
00385
                          free (NewBuffer);
00386
                       else
00387
                       {
00388
                           (*CardInfo) = NewBuffer;
                           memset(NewBuffer[Index].IPAddr, 0, sizeof(NewBuffer[Index].IPAddr));
00389
00390
                           NewBuffer[Index].ControlPort = 0;
00391
00392
                  }
00393
                  else
00394
00395
                      NewBuffer[Index].Timeout = ZESTETM1_DEFAULT_TIMEOUT;
00396
                       (*CardInfo) = NewBuffer;
00397
                       (*NumCards)++;
00398
00399
00400
          }
```

7.30 UPnP.c 359

```
00401 }
00402
00403
00404 /***********************
00405 \star Multicasting functions to join and leave a group \star
00407 static int ZestETM1_JoinGroup(SOCKET sd, uint32_t grpaddr,
00408
                                  uint32_t iaddr)
00409 {
00410
         struct ip_mreq imr;
00411
00412
          imr.imr_multiaddr.s_addr = grpaddr;
00413
         imr.imr_interface.s_addr = iaddr;
         return setsockopt(sd, IPPROTO_IP, IP_ADD_MEMBERSHIP,
00415
                           (const char *)&imr, sizeof(imr));
00416 }
00417 static int ZestETM1_LeaveGroup(SOCKET sd, uint32_t grpaddr,
00418
                                    uint32 t iaddr)
00419 {
00420
         struct ip_mreq imr;
00421
         imr.imr_multiaddr.s_addr = grpaddr;
00422
         imr.imr_interface.s_addr = iaddr;
00423
00424
         return setsockopt(sd, IPPROTO_IP, IP_DROP_MEMBERSHIP,
00425
                           (const char *)&imr, sizeof(imr));
00426 }
00427
00428 /**********************
00429 \star Get an array with all local IP addresses \star
00430 *****************************
00431 static ZESTETM1_STATUS ZestETM1_GetAllAdapters(uint32_t *NumAdapters, struct sockaddr_in **Adapters)
00432 {
00433 #if defined(MSVC) || defined(WINGCC)
          SOCKET Socket;
00434
          SOCKET_ADDRESS_LIST *AddressListPtr;
00435
00436
         DWORD BytesRequired;
00437
          int i;
00438
          int Count = 0;
00439
00440
         *NumAdapters = 0;
00441
         *Adapters = NULL;
00442
00443
          Socket = socket(AF_INET, SOCK_DGRAM, IPPROTO_UDP);
         if (Socket<0)</pre>
00444
00445
             return ZESTETM1_INTERNAL_ERROR;
00446
00447
         WSAIoctl(Socket, SIO_ADDRESS_LIST_QUERY, NULL, 0,
00448
                  NULL, 0, (LPDWORD) & BytesRequired, NULL, NULL);
          AddressListPtr = (SOCKET_ADDRESS_LIST *) malloc(BytesRequired);
00449
00450
          if (AddressListPtr==NULL)
00451
00452
              closesocket(Socket);
00453
              return ZESTETM1_OUT_OF_MEMORY;
00454
00455
         if (WSAIoctl(Socket, SIO_ADDRESS_LIST_QUERY, NULL, 0,
00456
                       AddressListPtr, BytesRequired, &BytesRequired, NULL, NULL) ==SOCKET_ERROR)
00457
00458
              free (AddressListPtr);
00459
              closesocket(Socket);
00460
              return ZESTETM1_INTERNAL_ERROR;
00461
00462
00463
         for (i=0; i<AddressListPtr->iAddressCount; i++)
00464
00465
              if (AddressListPtr->Address[i].iSockaddrLength==sizeof(struct sockaddr_in))
00466
              {
00467
00468
                  (*Adapters) = (struct sockaddr_in *)realloc(*Adapters, Count*sizeof(struct sockaddr_in));
00469
                  if ((*Adapters) ==NULL)
00470
00471
                      free (AddressListPtr);
00472
                      closesocket(Socket);
00473
                      return ZESTETM1_OUT_OF_MEMORY;
00474
00475
                 memcpy(&(*Adapters)[Count-1], AddressListPtr->Address[i].lpSockaddr, sizeof(struct
     sockaddr in));
00476
00477
00478
          *NumAdapters = Count;
00479
          closesocket (Socket):
00480
         free (AddressListPtr):
```

```
00481
00482
         return ZESTETM1_SUCCESS;
00483 #else
00484
        struct ifaddrs *Interfaces;
00485
         struct ifaddrs *Ptr;
00486
         uint32_t Count = 0;
00487
00488
         if (getifaddrs(&Interfaces)!=0)
00489
            return ZESTETM1_INTERNAL_ERROR;
00490
00491
         *NumAdapters = 0;
00492
         *Adapters = NULL;
00493
00494
         Ptr = Interfaces;
00495
         while (Ptr!=NULL)
00496
00497
             if (Ptr->ifa_addr!=NULL)
00498
             {
00499
                 Count++;
00500
                 (*Adapters) = (struct sockaddr_in *)realloc(*Adapters, Count*sizeof(struct sockaddr_in));
00501
                 if ((*Adapters)==NULL)
00502
                 {
00503
                     freeifaddrs (Interfaces);
00504
                     return ZESTETM1_OUT_OF_MEMORY;
00505
00506
                 memcpy(&(*Adapters)[Count-1], Ptr->ifa_addr, sizeof(struct sockaddr_in));
00507
00508
             Ptr = Ptr->ifa next;
00509
         *NumAdapters = Count;
00510
00511
         freeifaddrs (Interfaces);
00512
00513
         return ZESTETM1 SUCCESS:
00514 #endif
00515 }
00516
00518 * Scan networks for ZestETM1 cards and return the number of attached devices
00519 \star and details about each one
00521 ZESTETM1_STATUS ZestETM1CountCards(unsigned long *NumCards,
00522
                                     ZESTETM1_CARD_INFO **CardInfo, unsigned long Wait)
00523 {
00524
         SOCKET Socket;
00525
         struct sockaddr_in DestIP;
00526
         struct sockaddr_in SourceIP;
00527
         int SourceIPLength;
00528
         int Flag = 1;
00529
         int Result;
00530
         char Req[1024];
00531
         char Response[1024];
00532
         ZESTETM1_CARD_INFO *Cards = NULL;
00533
         uint32_t CardCount = 0;
00534
         int i;
00535
         struct timeval Timeout;
00536
         fd_set ReadFDS;
00537
         uint32_t Interface;
00538
         uint32_t NumAdapters;
00539
         struct sockaddr_in *Adapters;
00540
         ZESTETM1_STATUS Status;
00541
00542
         *NumCards = 0;
00543
         *CardInfo = NULL;
00544
00545
         // Get a list of all adapters
         Status = ZestETM1_GetAllAdapters(&NumAdapters, &Adapters);
00546
00547
         if (Status!=ZESTETM1_SUCCESS)
             ZESTETM1_ERROR_GENERAL("ZestETM1CountCards", Status);
00548
00549
         if (NumAdapters==0)
00550
         {
00551
             *NumCards = 0;
00552
             return ZESTETM1_SUCCESS;
00553
00554
00555
         // Send queries on all interfaces
00556
         for (Interface=0; Interface<NumAdapters; Interface++)</pre>
00557
00558
             if (Adapters[Interface].sin_family!=AF_INET)
00559
                 continue;
00560
00561
             // Open socket for search requests
```

7.30 UPnP.c 361

```
00562
              Socket = socket(AF_INET, SOCK_DGRAM, IPPROTO_UDP);
00563
               if (Socket<0)
00564
00565
                  continue;
00566
00567
00568
              // Set reuse port to on to allow multiple binds per host
00569
              if (setsockopt(Socket, SOL_SOCKET, SO_REUSEADDR, (char *)&Flag,
00570
                              sizeof(Flag))<0)
00571
00572
                  closesocket(Socket);
00573
                  continue;
00574
              }
00575
00576
               // Bind to port for receiving responses
00577
              SourceIP.sin_family = AF_INET;
00578
              SourceIP.sin_addr.s_addr = htonl(INADDR_ANY);
00579
              SourceIP.sin_port = 0;
00580
              Result = bind(Socket, (const struct sockaddr *)(&Adapters[Interface]), sizeof(struct
      sockaddr in));
00581
              if (Result<0)
00582
              {
00583
                  closesocket(Socket);
00584
                  continue:
00585
00586
00587
              // Join multicast group
              if (ZestETM1_JoinGroup(Socket, inet_addr(ZESTETM1_UPNP_ADDR),
00588
00589
                                    htonl(INADDR_ANY))<0)
00590
00591
                  closesocket(Socket);
00592
                  continue;
00593
              }
00594
              // Send M-SEARCH request
00595
              // Send more than once as UDP is unreliable
00596
00597
              DestIP.sin_family = AF_INET;
00598
              DestIP.sin_addr.s_addr = inet_addr(ZESTETM1_UPNP_ADDR);
              DestIP.sin_port = htons(1900);
00599
00600
              sprintf(Req, ZestETM1_SearchReq, (Wait+999)/1000);
00601
              for (i=0; i<3; i++)</pre>
00602
               {
00603
                  Result = sendto(Socket, ZestETM1_SearchReq,
00604
                                   (int)strlen(ZestETM1_SearchReq),
00605
                                   0, (struct sockaddr *)&DestIP, sizeof(DestIP));
00606
                  if (Result!=strlen(ZestETM1_SearchReq))
00607
00608
                       ZestETM1_LeaveGroup(Socket, inet_addr(ZESTETM1_UPNP_ADDR),
00609
                                          htonl(INADDR_ANY));
00610
                       closesocket(Socket);
00611
                       continue;
00612
00613
00614
00615
               // Read responses
00616
              Timeout.tv_sec = Wait/1000;
00617
              Timeout.tv_usec = (Wait%1000)*1000;
00618
00619
00620
                  FD_ZERO(&ReadFDS);
00621
                  FD_SET(Socket, &ReadFDS);
00622
                  Result = select((int)Socket+1, &ReadFDS, NULL, NULL, &Timeout);
00623
                  if (Result<0)</pre>
00624
                  {
00625
00626
00627
                  if (!FD_ISSET(Socket, &ReadFDS)) break;
00628
00629
                  SourceIPLength = sizeof(SourceIP);
00630
                  Result = recvfrom(Socket, Response, sizeof(Response),
00631
                                     0, (struct sockaddr *) &SourceIP, &SourceIPLength);
00632
                   if (Result<0)
00633
00634
                       // Error!
00635
                      break:
00636
00637
                  else if (Result==0)
00638
00639
                       // Clean shut down
00640
                      break:
00641
                  }
```

```
00642
00644
                     // Parse results
                    if (_strnicmp("NOTIFY", Response, 6) == 0 ||
00645
                        _strnicmp("HTTP/1.1 200 OK", Response, 15) == 0)
00647
                        // Check its a GigExpedite and request XML description
00649
                        char *Server = strstr(Response, "GigExpedite2");
00650
                        char *Location = strstr(Response, "LOCATION");
                        if (Server!=NULL && Location!=NULL)
00652
00654
                            for (c=Location+10; *c!=0 && *c!='\n' && *c!='\r'; c++);
00656
                            ZestETM1_GetCardInfo(Location+10, &CardCount, &Cards, Wait);
00658
00659
             } while(1);
00660
00661
00662
             // Leave multicast group
00663
             ZestETM1_LeaveGroup(Socket, inet_addr(ZESTETM1_UPNP_ADDR),
00664
                               htonl(INADDR ANY));
00665
00666
             // Close socket
00667
             closesocket(Socket);
        }
00668
00669
00670
        *NumCards = CardCount;
*CardInfo = Cards;
00671
00672
        free (Adapters);
00673
00674
         return ZESTETM1 SUCCESS:
00675 }
00676
00678 * Free data structures returned by ZestETM1CountCards *
00680 ZESTETM1_STATUS ZestETM1FreeCards(ZESTETM1_CARD_INFO *CardInfo)
00681 {
         if (CardInfo!=NULL)
00682
00683
             free (CardInfo);
00684
00685
         return ZESTETM1 SUCCESS;
00686 }
00687
00688 /*****************
00689 \star Fill in card information fields \star
00690 ****************************
00691 ZESTETM1_STATUS ZestETM1GetCardInfo(ZESTETM1_CARD_INFO *CardInfo)
00692 {
00693
         ZESTETM1_STATUS Status;
00694
00695
         Status = ZestETM1_ReadSettings(CardInfo);
00696
         if (Status!=ZESTETM1_SUCCESS)
00697
00698
             ZESTETM1_ERROR("ZestETM1GetCardInfo", Status);
00699
00700
00701
         return ZESTETM1_SUCCESS;
00702 }
00703
```

# 7.31 C:/Users/hill35/git/camera\_python/nsCamera/comms/ZestETM1/Zest⊸ ETM1.h File Reference

#### Classes

• struct ZESTETM1 CARD INFO

#### **Macros**

- #define ZESTETM1 VERSION FALLBACK 0x8000
- #define ZESTETM1 INFO BASE 0
- #define ZESTETM1 WARNING BASE 0x4000
- #define ZESTETM1\_ERROR\_BASE 0x8000

# **Typedefs**

- typedef void \* ZESTETM1\_HANDLE
- typedef void \* ZESTETM1\_CONNECTION
- typedef void(\* ZESTETM1\_ERROR\_FUNC) (const char \*Function, ZESTETM1\_CARD\_INFO \*CardInfo, ZESTETM1\_STATUS Status, const char \*Msg)

#### **Enumerations**

- enum ZESTETM1 CONNECTION TYPE { ZESTETM1 TYPE TCP , ZESTETM1 TYPE UDP }
- enum ZESTETM1\_SPI\_RATE { ZESTETM1\_SPI\_RATE\_35 , ZESTETM1\_SPI\_RATE\_17\_5 , ZESTETM1\_SPI\_RATE\_8\_75 }
- enum ZESTETM1\_STATUS {
   ZESTETM1\_SUCCESS = ZESTETM1\_INFO\_BASE , ZESTETM1\_MAX\_INFO , ZESTETM1\_MAX\_WARNING
   = ZESTETM1\_WARNING\_BASE , ZESTETM1\_SOCKET\_ERROR = ZESTETM1\_ERROR\_BASE ,
   ZESTETM1\_INTERNAL\_ERROR , ZESTETM1\_ILLEGAL\_STATUS\_CODE , ZESTETM1\_NULL\_PARAMETER ,
   ZESTETM1\_OUT\_OF\_MEMORY ,
   ZESTETM1\_INVALID\_CONNECTION\_TYPE , ZESTETM1\_ILLEGAL\_CONNECTION , ZESTETM1\_SOCKET\_CLOSED , ZESTETM1\_TIMEOUT ,
   ZESTETM1\_ILLEGAL\_PARAMETER , ZESTETM1\_MAX\_ERROR }

# **Functions**

- ZESTETM1 STATUS ZestETM1Init (void)
- ZESTETM1 STATUS ZestETM1Close (void)
- ZESTETM1\_STATUS ZestETM1CountCards (unsigned long \*NumCards, ZESTETM1\_CARD\_INFO \*\*CardInfo, unsigned long Wait)
- ZESTETM1\_STATUS ZestETM1GetCardInfo (ZESTETM1\_CARD\_INFO \*CardInfo)
- ZESTETM1\_STATUS ZestETM1FreeCards (ZESTETM1\_CARD\_INFO \*CardInfo)
- ZESTETM1\_STATUS ZestETM1RegisterErrorHandler (ZESTETM1\_ERROR\_FUNC Function)
- ZESTETM1\_STATUS ZestETM1GetErrorMessage (ZESTETM1\_STATUS Status, char \*\*Buffer)
- ZESTETM1\_STATUS ZestETM1OpenConnection (ZESTETM1\_CARD\_INFO \*CardInfo, ZESTETM1\_CONNECTION\_TYPE
  Type, unsigned short Port, unsigned short LocalPort, ZESTETM1\_CONNECTION \*Connection)
- ZESTETM1\_STATUS ZestETM1CloseConnection (ZESTETM1\_CONNECTION Connection)
- ZESTETM1\_STATUS ZestETM1WriteData (ZESTETM1\_CONNECTION Connection, void \*Buffer, unsigned long Length, unsigned long \*Written, unsigned long Timeout)
- ZESTETM1\_STATUS ZestETM1ReadData (ZESTETM1\_CONNECTION Connection, void \*Buffer, unsigned long Length, unsigned long \*Read, unsigned long Timeout)
- ZESTETM1\_STATUS ZestETM1SPIReadWrite (ZESTETM1\_CARD\_INFO \*CardInfo, ZESTETM1\_SPI\_RATE Rate, int WordLen, void \*WriteData, void \*ReadData, unsigned long Length, int ReleaseCS)
- ZESTETM1\_STATUS ZestETM1WriteRegister (ZESTETM1\_CARD\_INFO \*CardInfo, unsigned long Addr, unsigned short Data)
- ZESTETM1\_STATUS ZestETM1ReadRegister (ZESTETM1\_CARD\_INFO \*CardInfo, unsigned long Addr, unsigned short \*Data)
- ZESTETM1 STATUS ZestETM1SetInterrupt (ZESTETM1 CARD INFO \*CardInfo)

# 7.31.1 Class Documentation

# 7.31.1.1 struct ZESTETM1\_CARD\_INFO

Definition at line 30 of file ZestETM1.h.

# **Class Members**

unsigned short	ControlPort	
unsigned long	FirmwareVersion	
unsigned char	Gateway[4]	
unsigned long	HardwareVersion	
unsigned short	HTTPPort	
unsigned char	IPAddr[4]	
unsigned char	MACAddr[6]	
unsigned long	SerialNumber	
unsigned char	SubNet[4]	
unsigned long	Timeout	

# 7.31.2 Macro Definition Documentation

# 7.31.2.1 ZESTETM1\_ERROR\_BASE

#define ZESTETM1\_ERROR\_BASE 0x8000

Definition at line 77 of file ZestETM1.h.

# 7.31.2.2 ZESTETM1\_INFO\_BASE

#define ZESTETM1\_INFO\_BASE 0

Definition at line 75 of file ZestETM1.h.

# 7.31.2.3 ZESTETM1\_VERSION\_FALLBACK

#define ZESTETM1\_VERSION\_FALLBACK 0x8000

Definition at line 50 of file ZestETM1.h.

# 7.31.2.4 ZESTETM1\_WARNING\_BASE

#define ZESTETM1\_WARNING\_BASE 0x4000

Definition at line 76 of file ZestETM1.h.

# 7.31.3 Typedef Documentation

# 7.31.3.1 ZESTETM1\_CONNECTION

```
typedef void* ZESTETM1_CONNECTION
```

Definition at line 55 of file ZestETM1.h.

# 7.31.3.2 ZESTETM1\_ERROR\_FUNC

```
typedef void(* ZESTETM1_ERROR_FUNC) (const char *Function, ZESTETM1_CARD_INFO *CardInfo, ZESTETM1_STATUS Status, const char *Msq)
```

Definition at line 98 of file ZestETM1.h.

# 7.31.3.3 ZESTETM1\_HANDLE

```
typedef void* ZESTETM1_HANDLE
```

Definition at line 24 of file ZestETM1.h.

# 7.31.4 Enumeration Type Documentation

# 7.31.4.1 ZESTETM1\_CONNECTION\_TYPE

```
enum ZESTETM1_CONNECTION_TYPE
```

# Enumerator

```
ZESTETM1_TYPE_TCP
ZESTETM1_TYPE_UDP
```

# Definition at line 56 of file ZestETM1.h.

# 7.31.4.2 ZESTETM1 SPI RATE

```
enum ZESTETM1_SPI_RATE
```

# **Enumerator**

ZESTETM1 SPI RATE 35

#### Enumerator

```
ZESTETM1_SPI_RATE_17↔
_5
ZESTETM1_SPI_RATE_8_75
```

# Definition at line 65 of file ZestETM1.h.

# 7.31.4.3 ZESTETM1\_STATUS

```
enum ZESTETM1_STATUS
```

#### Enumerator

ZESTETM1_SUCCESS	
ZESTETM1_MAX_INFO	
ZESTETM1_MAX_WARNING	
ZESTETM1_SOCKET_ERROR	
ZESTETM1_INTERNAL_ERROR	
ZESTETM1_ILLEGAL_STATUS_CODE	
ZESTETM1_NULL_PARAMETER	
ZESTETM1_OUT_OF_MEMORY	
ZESTETM1_INVALID_CONNECTION_TYPE	
ZESTETM1_ILLEGAL_CONNECTION	
ZESTETM1_SOCKET_CLOSED	
ZESTETM1_TIMEOUT	
ZESTETM1_ILLEGAL_PARAMETER	
ZESTETM1_MAX_ERROR	

#### Definition at line 78 of file ZestETM1.h.

```
00079 {
              ZESTETM1_SUCCESS = ZESTETM1_INFO_BASE,
00081
              ZESTETM1_MAX_INFO,
00082
00083
              ZESTETM1_MAX_WARNING = ZESTETM1_WARNING_BASE,
00084
00085
              ZESTETM1_SOCKET_ERROR = ZESTETM1_ERROR_BASE,
             ZESTETM1_INTERNAL_ERROR,
ZESTETM1_ILLEGAL_STATUS_CODE,
00086
00087
00088
              ZESTETM1_NULL_PARAMETER,
             ZESTETM1_OUT_OF_MEMORY,
ZESTETM1_INVALID_CONNECTION_TYPE,
00089
00090
             ZESTETMI_INVALID_CONNECTION_
ZESTETMI_SOCKET_CLOSED,
ZESTETMI_SOCKET_CLOSED,
ZESTETMI_TIMEOUT,
ZESTETMI_ILLEGAL_PARAMETER,
00091
00092
00093
00094
00095
00096
             ZESTETM1_MAX_ERROR
00097 } ZESTETM1_STATUS;
```

# 7.31.5 Function Documentation

#### 7.31.5.1 ZestETM1Close()

# 7.31.5.2 ZestETM1CloseConnection()

```
ZESTETM1_STATUS ZestETM1CloseConnection ( {\tt ZESTETM1\_CONNECTION}\ \ Connection\ )
```

#### Definition at line 575 of file Data.c.

```
ZESTETM1_STATUS Result;
00578
          ZESTETM1_CONNECTION_STRUCT *Conn = (ZESTETM1_CONNECTION_STRUCT*)Connection;
00579
00580
         Result = ZestETM1_CloseConnection(Connection);
        if (Result!=ZESTETM1_SUCCESS)
{
00581
00582
00583
              ZESTETM1_ERROR_CONN("ZestETM1CloseConnection", Result);
00584
00585
00586
          return ZESTETM1_SUCCESS;
00587 }
```

#### 7.31.5.3 ZestETM1CountCards()

```
ZESTETM1_STATUS ZestETM1CountCards (
          unsigned long * NumCards,
          ZESTETM1_CARD_INFO ** CardInfo,
          unsigned long Wait )
```

# Definition at line 521 of file UPnP.c.

```
00523 {
00524
          SOCKET Socket;
          struct sockaddr_in DestIP;
00526
          struct sockaddr_in SourceIP;
          int SourceIPLength;
00528
          int Flag = 1;
00529
         int Result;
          char Req[1024];
00530
         char Response[1024];
ZESTETM1_CARD_INFO *Cards = NULL;
00531
00532
00533
         uint32_t CardCount = 0;
00534
          int i;
00535
         struct timeval Timeout;
00536
          fd_set ReadFDS;
00537
          uint32 t Interface;
00538
          uint32_t NumAdapters;
          struct sockaddr_in *Adapters;
00539
00540
          ZESTETM1_STATUS Status;
00541
00542
          *NumCards = 0;
```

```
00543
          *CardInfo = NULL;
00544
00545
          // Get a list of all adapters
00546
          Status = ZestETM1_GetAllAdapters(&NumAdapters, &Adapters);
00547
          if (Status!=ZESTETM1_SUCCESS)
00548
              ZESTETM1_ERROR_GENERAL("ZestETM1CountCards", Status);
00549
          if (NumAdapters==0)
00550
00551
              *NumCards = 0;
00552
              return ZESTETM1_SUCCESS;
00553
00554
00555
          // Send queries on all interfaces
00556
          for (Interface=0; Interface<NumAdapters; Interface++)</pre>
00557
00558
              if (Adapters[Interface].sin_family!=AF_INET)
00559
                 continue:
00560
00561
              // Open socket for search requests
              Socket = socket(AF_INET, SOCK_DGRAM, IPPROTO_UDP);
00562
00563
              if (Socket<0)
00564
00565
                  continue;
00566
00567
              // Set reuse port to on to allow multiple binds per host
00568
             00569
00570
00571
              {
00572
                  closesocket(Socket);
00573
                 continue;
00574
00575
              // Bind to port for receiving responses
00576
00577
              SourceIP.sin_family = AF_INET;
00578
              SourceIP.sin_addr.s_addr = htonl(INADDR_ANY);
00579
              SourceIP.sin_port = 0;
             Result = bind(Socket, (const struct sockaddr *)(&Adapters[Interface]), sizeof(struct
00580
     sockaddr_in));
00581
             if (Result<0)
00582
00583
                  closesocket(Socket);
00584
                  continue;
00585
             }
00586
00587
              // Join multicast group
00588
              if (ZestETM1_JoinGroup(Socket, inet_addr(ZESTETM1_UPNP_ADDR),
00589
                                    htonl(INADDR_ANY))<0)
00590
00591
                  closesocket(Socket);
00592
                  continue;
00593
00594
              // Send M-SEARCH request
00595
00596
              // Send more than once as UDP is unreliable
00597
              DestIP.sin_family = AF_INET;
00598
              DestIP.sin_addr.s_addr = inet_addr(ZESTETM1_UPNP_ADDR);
00599
              DestIP.sin_port = htons(1900);
00600
              sprintf(Req, ZestETM1_SearchReq, (Wait+999)/1000);
00601
              for (i=0; i<3; i++)
00602
00603
                  Result = sendto(Socket, ZestETM1_SearchReq,
00604
                                  (int) strlen (ZestETM1_SearchReq),
00605
                                  0, (struct sockaddr *)&DestIP, sizeof(DestIP));
                  if (Result!=strlen(ZestETM1_SearchReq))
00606
00607
                  {
                      ZestETM1_LeaveGroup(Socket, inet_addr(ZESTETM1_UPNP_ADDR),
00609
                                         htonl(INADDR_ANY));
00610
                      closesocket(Socket);
00611
                      continue;
00612
                  }
00613
             }
00614
00615
              // Read responses
00616
              Timeout.tv_sec = Wait/1000;
00617
              Timeout.tv_usec = (Wait%1000)*1000;
00618
00619
00620
                  FD_ZERO(&ReadFDS);
                  FD_SET(Socket, &ReadFDS);
00621
                  Result = select((int)Socket+1, &ReadFDS, NULL, NULL, &Timeout);
00622
```

```
00623
                    if (Result<0)
00624
00625
                        break;
00626
00627
                    if (!FD_ISSET(Socket, &ReadFDS)) break;
                    SourceIPLength = sizeof(SourceIP);
00630
                    Result = recvfrom(Socket, Response, sizeof(Response),
00631
                                        0, (struct sockaddr *)&SourceIP, &SourceIPLength);
                    if (Result<0)
00633
00634
                         // Error!
00635
                        break;
00636
00637
                    else if (Result==0)
00638
00639
                         // Clean shut down
00640
                        break:
00641
00642
                    else
00643
                    {
00644
                         // Parse results
                        if (_strnicmp("NOTIFY", Response, 6) == 0 ||
    _strnicmp("HTTP/1.1 200 OK", Response, 15) == 0)
00645
00646
00647
                             // Check its a GigExpedite and request XML description
char *Server = strstr(Response, "GigExpedite2");
00648
00649
                             char *Location = strstr(Response, "LOCATION");
00650
00651
                             if (Server!=NULL && Location!=NULL)
00652
00653
                                 char *c:
00654
                                 for (c=Location+10; *c!=0 && *c!='\n' && *c!='\r'; c++);
00655
                                  *C = 0:
00656
                                 ZestETM1_GetCardInfo(Location+10, &CardCount, &Cards, Wait);
00657
00658
00659
00660
               } while(1);
00661
               // Leave multicast group
00662
               ZestETM1_LeaveGroup(Socket, inet_addr(ZESTETM1_UPNP_ADDR),
00663
00664
                                    htonl(INADDR_ANY));
00665
               // Close socket
00666
00667
               closesocket(Socket);
00668
00669
00670
           *NumCards = CardCount;
00671
           *CardInfo = Cards;
00672
           free (Adapters);
00673
00674
           return ZESTETM1_SUCCESS;
00675 }
```

# 7.31.5.4 ZestETM1FreeCards()

free (CardInfo);

return ZESTETM1\_SUCCESS;

# 7.31.5.5 ZestETM1GetCardInfo()

00683

00684 00685

00686 }

# Definition at line 691 of file UPnP.c.

```
00692 {
          ZESTETM1_STATUS Status;
00693
00694
00695
          Status = ZestETM1_ReadSettings(CardInfo);
00696
          if (Status!=ZESTETM1_SUCCESS)
00697
00698
              ZESTETM1_ERROR("ZestETM1GetCardInfo", Status);
00699
00700
00701
          return ZESTETM1_SUCCESS;
00702 }
```

# 7.31.5.6 ZestETM1GetErrorMessage()

```
ZESTETM1_STATUS ZestETM1GetErrorMessage (
               ZESTETM1_STATUS Status,
               char ** Buffer )
Definition at line 84 of file Error.c.
00086 {
00087
          if (Status>ZESTETM1_MAX_ERROR | |
00088
              (Status<ZESTETM1_ERROR_BASE && Status>=ZESTETM1_MAX_WARNING) ||
00089
              (Status<ZESTETM1_WARNING_BASE && Status>=ZESTETM1_MAX_INFO))
00090
00091
              return ZESTETM1_ILLEGAL_STATUS_CODE;
00092
00093
00094
          *Buffer = ZESTETM1_ERROR_STRING(Status);
```

# 7.31.5.7 ZestETM1Init()

00095

00096 }

return ZESTETM1\_SUCCESS;

# Definition at line 58 of file Main.c.

```
00060 #if defined(MSVC) || defined(WINGCC)
00061
          WORD VersionRequested;
          WSADATA WSAData;
00062
00063
          int Error;
00064
00065
          VersionRequested = MAKEWORD(2, 2);
00066
          Error = WSAStartup(VersionRequested, &WSAData);
00067
          if (Error!=0)
00068
00069
              ZESTETM1_ERROR_GENERAL("ZestETM1Init", ZESTETM1_SOCKET_ERROR);
00070
00071
00072
          // Confirm that the WinSock DLL supports 2.2.
00073
          // Note that if the DLL supports versions greater
00074
          \ensuremath{//} than 2.2 in addition to 2.2, it will still return
00075
          // 2.2 in Version since that is the version we
00076
          // requested.
00077
          if (LOBYTE(WSAData.wVersion)!=2 ||
00078
              HIBYTE(WSAData.wVersion)!=2)
00079
              WSACleanup();
00080
              ZESTETM1_ERROR_GENERAL("ZestETM1Init", ZESTETM1_SOCKET_ERROR);
00081
00082
00083 #endif
00084
00085
          return ZESTETM1_SUCCESS;
00086 }
```

# 7.31.5.8 ZestETM1OpenConnection()

# 7.31.5.9 ZestETM1ReadData()

```
ZESTETM1_STATUS ZestETM1ReadData (
    ZESTETM1_CONNECTION Connection,
    void * Buffer,
    unsigned long Length,
    unsigned long * Read,
    unsigned long Timeout )
```

# Definition at line 613 of file Data.c.

```
00618 {
00619
          ZESTETM1_STATUS Result;
00620
          ZESTETM1_CONNECTION_STRUCT *Conn = (ZESTETM1_CONNECTION_STRUCT*)Connection;
00621
00622
          Result = ZestETM1_ReadData(Connection, Buffer, Length, Read, Timeout);
00623
         if (Result!=ZESTETM1_SUCCESS)
00624
00625
              ZESTETM1_ERROR_CONN("ZestETM1ReadData", Result);
00626
00627
          return ZESTETM1_SUCCESS;
00628
00629 }
```

#### 7.31.5.10 ZestETM1ReadRegister()

#### Definition at line 728 of file Data.c.

```
00730
          ZESTETM1_STATUS Result;
00731
          ZESTETM1_CONNECTION Connection;
00732
          ZESTETM1_READ_REG_CMD Cmd;
00733
          ZESTETM1_READ_REG_RESPONSE Response;
00734
00735
          if (CardInfo==NULL || Data==NULL)
00736
00737
              ZESTETM1_ERROR("ZestETM1ReadRegister", ZESTETM1_NULL_PARAMETER);
00738
00739
         if (Addr>127)
00740
         {
00741
              ZESTETM1_ERROR("ZestETM1ReadRegister", ZESTETM1_ILLEGAL_PARAMETER);
00742
00743
00744
          Result = ZestETM1_OpenConnection(CardInfo, ZESTETM1_TYPE_TCP, CardInfo->ControlPort, 0, &Connection);
00745
          if (Result!=ZESTETM1_SUCCESS)
00746
00747
              ZESTETM1_ERROR("ZestETM1ReadRegister", Result);
00748
00749
```

```
00750
          // Read register from the device
00751
          Cmd.Command = ZESTETM1_COMMAND_READ_REG;
00752
          Cmd.Addr = (uint8_t)Addr;
00753
          Result = ZestETM1_SendCommand(CardInfo, Connection,
00754
                                        &Cmd, sizeof(Cmd),
                                        &Response, sizeof(Response), 1);
00755
00756
          if (Result!=ZESTETM1_SUCCESS)
00757
00758
              ZestETM1_CloseConnection(Connection);
00759
              ZESTETM1_ERROR("ZestETM1ReadRegister", Result);
00760
00761
          if (Response.Command!=Cmd.Command || Response.Status!=0)
00762
00763
              ZestETM1_CloseConnection(Connection);
00764
              ZESTETM1_ERROR("ZestETM1ReadRegister", ZESTETM1_INTERNAL_ERROR);
00765
00766
          *Data = ((Response.Value»8)&0xff) | ((Response.Value&0xff) < 8);
00767
00768
          Result = ZestETM1_CloseConnection(Connection);
00769
          if (Result!=ZESTETM1_SUCCESS)
00770
00771
              ZESTETM1_ERROR("ZestETM1ReadRegister", Result);
00772
00773
00774
          return ZESTETM1 SUCCESS:
00775 }
```

# 7.31.5.11 ZestETM1RegisterErrorHandler()

#### Definition at line 74 of file Error.c.

```
00076 ZestETM1_ErrorHandler = Function;
00077 return ZESTETM1_SUCCESS;
00078 }
```

# 7.31.5.12 ZestETM1SetInterrupt()

# Definition at line 780 of file Data.c.

```
00781 {
00782
          ZESTETM1_STATUS Result;
00783
          ZESTETM1_CONNECTION Connection;
00784
          ZESTETM1_MAILBOX_INT_CMD Cmd;
00785
          ZESTETM1_MAILBOX_INT_RESPONSE Response;
00786
00787
          if (CardInfo==NULL)
00788
          {
00789
              ZESTETM1_ERROR("ZestETM1SetInterrupt", ZESTETM1_NULL_PARAMETER);
00790
00791
00792
          Result = ZestETM1_OpenConnection(CardInfo, ZESTETM1_TYPE_TCP, CardInfo->ControlPort, 0, &Connection);
00793
          if (Result!=ZESTETM1_SUCCESS)
00794
00795
              ZESTETM1_ERROR("ZestETM1SetInterrupt", Result);
00796
00797
00798
          // Send command to set interrupt
00799
          Cmd.Command = ZESTETM1_COMMAND_MAILBOX_INT;
          Result = ZestETM1_SendCommand(CardInfo, Connection,
00800
00801
                                        &Cmd, sizeof(Cmd),
                                        &Response, sizeof(Response), 1);
00802
00803
          if (Result!=ZESTETM1_SUCCESS)
00804
00805
              ZestETM1 CloseConnection(Connection);
```

```
00806
              ZESTETM1_ERROR("ZestETM1SetInterrupt", Result);
00807
00808
          if (Response.Command!=Cmd.Command || Response.Status!=0)
00809
00810
              ZestETM1_CloseConnection(Connection);
00811
              ZESTETM1_ERROR("ZestETM1SetInterrupt", ZESTETM1_INTERNAL_ERROR);
00812
00813
00814
         Result = ZestETM1_CloseConnection(Connection);
00815
         if (Result!=ZESTETM1_SUCCESS)
00816
         {
00817
              ZESTETM1_ERROR("ZestETM1SetInterrupt", Result);
00818
00820
          return ZESTETM1_SUCCESS;
00821 }
```

# 7.31.5.13 ZestETM1SPIReadWrite()

#### Definition at line 635 of file Data.c.

```
00639 {
00640
          ZESTETM1_STATUS Result;
00641
          ZESTETM1_CONNECTION Connection;
00642
          unsigned long RateVal = Rate==ZESTETM1_SPI_RATE_35 ? ZESTETM1_RATE_40MHz :
00643
                                  Rate==ZESTETM1_SPI_RATE_17_5 ? ZESTETM1_RATE_20MHz : ZESTETM1_RATE_10MHz;
00644
00645
          if (CardInfo==NULL || (WriteData==NULL && ReadData==NULL))
00646
00647
              ZESTETM1_ERROR("ZestETM1SPIReadWrite", ZESTETM1_NULL_PARAMETER);
00648
00649
          if (WordLen<1 || WordLen>32 || Length>16384)
00650
00651
              ZESTETM1_ERROR("ZestETM1SPIReadWrite", ZESTETM1_ILLEGAL_PARAMETER);
00652
00653
00654
          Result = ZestETM1_OpenConnection(CardInfo, ZESTETM1_TYPE_TCP, CardInfo->ControlPort, 0, &Connection);
00655
          if (Result!=ZESTETM1_SUCCESS)
00656
00657
              ZESTETM1_ERROR("ZestETM1SPIReadWrite", Result);
00658
00659
          Result = ZestETM1_SPIReadWrite(CardInfo, Connection, ZESTETM1_USER_DEVICE_ID|RateVal, WordLen,
00660
      (uint32_t *)WriteData, (uint32_t *)ReadData, Length, ReleaseCS, 1);
00661
         if (Result!=ZESTETM1_SUCCESS)
00662
00663
              ZESTETM1_ERROR("ZestETM1SPIReadWrite", Result);
00664
00665
00666
         Result = ZestETM1_CloseConnection(Connection);
00667
          if (Result!=ZESTETM1_SUCCESS)
00668
00669
              ZESTETM1_ERROR("ZestETM1SPIReadWrite", Result);
00670
00671
00672
          return ZESTETM1 SUCCESS;
00673 }
```

# 7.31.5.14 ZestETM1WriteData()

```
void * Buffer,
               unsigned long Length,
               unsigned long * Written,
               unsigned long Timeout )
Definition at line 592 of file Data.c.
00597 {
00598
          ZESTETM1_STATUS Result;
00599
          ZESTETM1_CONNECTION_STRUCT *Conn = (ZESTETM1_CONNECTION_STRUCT*) Connection;
00600
00601
          Result = ZestETM1_WriteData(Connection, Buffer, Length, Written, Timeout);
00602
          if (Result!=ZESTETM1_SUCCESS)
00603
00604
              ZESTETM1_ERROR_CONN("ZestETM1WriteData", Result);
00605
00606
00607
          return ZESTETM1_SUCCESS;
00608 }
7.31.5.15 ZestETM1WriteRegister()
ZESTETM1 STATUS ZestETM1WriteRegister (
               ZESTETM1_CARD_INFO * CardInfo,
               unsigned long Addr,
               unsigned short Data )
Definition at line 678 of file Data.c.
00679 {
00680
          ZESTETM1_STATUS Result;
00681
          ZESTETM1_CONNECTION Connection;
          ZESTETM1 WRITE REG CMD Cmd:
00682
          ZESTETM1_WRITE_REG_RESPONSE Response;
00683
00684
00685
          if (CardInfo==NULL)
00686
00687
              ZESTETM1_ERROR("ZestETM1WriteRegister", ZESTETM1_NULL_PARAMETER);
00688
00689
          if (Addr>127)
00690
00691
              ZESTETM1_ERROR("ZestETM1WriteRegister", ZESTETM1_ILLEGAL_PARAMETER);
00692
00693
00694
          Result = ZestETM1_OpenConnection(CardInfo, ZESTETM1_TYPE_TCP, CardInfo->ControlPort, 0, &Connection);
00695
          if (Result!=ZESTETM1_SUCCESS)
00696
00697
              ZESTETM1_ERROR("ZestETM1WriteRegister", Result);
00698
00699
00700
          \ensuremath{//} Write register to the device
00701
          //FIXME: Do we want to be able to copy this value to flash?
00702
          Cmd.Command = ZESTETM1_COMMAND_WRITE_REG;
00703
          Cmd.Addr = (uint8_t)Addr;
00704
          Cmd.Data = ((Data * 8) & 0xff) | (Data & 0xff);
00705
          Result = ZestETM1_SendCommand(CardInfo, Connection,
00706
                                        &Cmd, sizeof(Cmd),
00707
                                        &Response, sizeof(Response), 1);
00708
          if (Result!=ZESTETM1_SUCCESS)
00709
          {
00710
              ZestETM1_CloseConnection(Connection);
00711
              ZESTETM1_ERROR("ZestETM1WriteRegister", Result);
00712
00713
          if (Response.Command!=Cmd.Command || Response.Status!=0)
00714
          {
00715
              ZestETM1_CloseConnection(Connection);
00716
              ZESTETM1_ERROR("ZestETM1WriteRegister", ZESTETM1_INTERNAL_ERROR);
00717
00718
00719
          Result = ZestETM1_CloseConnection(Connection);
00720
          if (Result!=ZESTETM1_SUCCESS)
00721
00722
              ZESTETM1_ERROR("ZestETM1WriteRegister", Result);
00723
00724
00725
          return ZESTETM1_SUCCESS;
```

00726 }

7.32 ZestETM1.h 375

# 7.32 ZestETM1.h

#### Go to the documentation of this file. 00003 \* (c) 2012 Orange Tree Technologies Ltd 00005 \* ZestETM1.h 00006 \* Version 1.0 00008 \* Header file for ZestETM1 Ethernet module 00011 00012 #ifndef \_\_\_ZESTETM1\_H\_ 00013 #define \_\_ZESTETM1\_H\_ 00015 #ifdef \_\_\_c 00016 extern "C" \_cplusplus 00017 { 00018 #endif 00019 00020 00021 /\*\*\*\*\*\*\*\*\*\*\*\*\*\*\*\*\* 00022 \* Handle for referencing modules \* 0.0023 \* 00024 typedef void \*ZESTETM1\_HANDLE; 00025 00026 00027 /\*\*\*\*\*\*\*\*\*\*\*\*\*\*\*\* 00028 $\star$ Card information structure $\star$ 00029 \* 00030 typedef struct 00031 { // These must be filled in before calling functions 00032 unsigned char IPAddr[4]; 00034 unsigned short ControlPort; 00035 unsigned long Timeout; 00036 00037 // These are for information purposes only 00038 unsigned short HTTPPort; 00039 unsigned char MACAddr[6]; 00040 unsigned char SubNet[4]; 00041 unsigned char Gateway[4]; 00042 unsigned long SerialNumber; 00043 unsigned long FirmwareVersion; 00044 unsigned long HardwareVersion; 00045 } ZESTETM1\_CARD\_INFO; 00046 00047 // Fallback mask 00048 // This will be set if the GigExpedite is operating in firmware version fallback mode 00049 // due to a failed upload of firmware 00050 #define ZESTETM1\_VERSION\_FALLBACK 0x8000 00053 \* Data transfer definitions \* 00054 \* 00055 typedef void \*ZESTETM1\_CONNECTION; 00056 typedef enum 00057 { 00058 ZESTETM1\_TYPE\_TCP, ZESTETM1\_TYPE\_UDP 00060 } ZESTETM1\_CONNECTION\_TYPE; 00062 /\*\*\*\*\*\*\*\*\*\*\*\*\*\* 00063 \* Master SPI clock rates \* 00064 \* 00065 typedef enum 00066 { 00067 ZESTETM1\_SPI\_RATE\_35, ZESTETM1\_SPI\_RATE\_17\_5, ZESTETM1\_SPI\_RATE\_8\_75, 00068 00069 00070 } ZESTETM1\_SPI\_RATE; 00071 00072 /\*\*\*\*\*\*\*\*\*\*\*\*\*\* 00073 \* Function return codes \* 00074 \*\*\*\*\*\*\*\*\*\*\*\*\*\*\*\* 00075 #define ZESTETM1\_INFO\_BASE 0 00076 #define ZESTETM1\_WARNING\_BASE 0x4000

```
00077 #define ZESTETM1_ERROR_BASE 0x8000
00078 typedef enum
00079 {
08000
          ZESTETM1_SUCCESS = ZESTETM1_INFO_BASE,
00081
          ZESTETM1_MAX_INFO,
00082
00083
          ZESTETM1_MAX_WARNING = ZESTETM1_WARNING_BASE,
00084
00085
          ZESTETM1_SOCKET_ERROR = ZESTETM1_ERROR_BASE,
          ZESTETM1_INTERNAL_ERROR,
00086
00087
          ZESTETM1_ILLEGAL_STATUS_CODE,
00088
          ZESTETM1_NULL_PARAMETER,
00089
          ZESTETM1_OUT_OF_MEMORY,
          ZESTETM1_INVALID_CONNECTION_TYPE,
00090
00091
          ZESTETM1_ILLEGAL_CONNECTION,
00092
          ZESTETM1_SOCKET_CLOSED,
          ZESTETM1_TIMEOUT,
00093
00094
          ZESTETM1_ILLEGAL_PARAMETER,
00095
00096
          ZESTETM1_MAX_ERROR
00097 } ZESTETM1_STATUS;
00098 typedef void (*ZESTETM1_ERROR_FUNC) (const char *Function,
00099
                                          ZESTETM1_CARD_INFO *CardInfo,
00100
                                          ZESTETM1_STATUS Status,
00101
                                          const char *Msg);
00102
00103
00104 /************
00105 * Function prototypes *
00106 ****************
00107 ZESTETM1_STATUS ZestETM1Init(void);
00108 ZESTETM1_STATUS ZestETM1Close(void);
00109 ZESTETM1_STATUS ZestETM1CountCards(unsigned long *NumCards,
                                        ZESTETM1 CARD INFO **CardInfo,
00110
                                        unsigned long Wait);
00111
00112 ZESTETM1_STATUS ZestETM1GetCardInfo(ZESTETM1_CARD_INFO *CardInfo);
00113 ZESTETM1_STATUS ZestETM1FreeCards(ZESTETM1_CARD_INFO *CardInfo);
00114
00115 ZESTETM1 STATUS ZestETM1RegisterErrorHandler(ZESTETM1 ERROR FUNC Function);
00116 ZESTETM1_STATUS ZestETM1GetErrorMessage(ZESTETM1_STATUS Status,
00117
                                             char **Buffer);
00118
00119 ZESTETM1_STATUS ZestETM1OpenConnection(ZESTETM1_CARD_INFO *CardInfo,
00120
                                            ZESTETM1_CONNECTION_TYPE Type,
00121
                                            unsigned short Port,
00122
                                            unsigned short LocalPort,
00123
                                            ZESTETM1_CONNECTION *Connection);
00124 ZESTETM1_STATUS ZestETM1CloseConnection(ZESTETM1_CONNECTION Connection);
00125 ZESTETM1_STATUS ZestETM1WriteData(ZESTETM1_CONNECTION Connection,
00126
                                       void *Buffer,
00127
                                       unsigned long Length,
00128
                                       unsigned long *Written,
00129
                                       unsigned long Timeout);
00130 ZESTETM1_STATUS ZestETM1ReadData(ZESTETM1_CONNECTION Connection,
00131
                                      void *Buffer,
00132
                                      unsigned long Length,
00133
                                      unsigned long *Read,
00134
                                      unsigned long Timeout);
00135
00136 ZESTETM1_STATUS ZestETM1SPIReadWrite(ZESTETM1_CARD_INFO *CardInfo, ZESTETM1_SPI_RATE Rate,
                                            int WordLen, void *WriteData,
                                            void *ReadData, unsigned long Length,
00138
00139
                                            int ReleaseCS);
00140
00141 ZESTETM1_STATUS ZestETM1WriteRegister(ZESTETM1_CARD_INFO *CardInfo, unsigned long Addr, unsigned short
00142 ZESTETM1_STATUS ZestETM1ReadRegister(ZESTETM1_CARD_INFO *CardInfo, unsigned long Addr, unsigned short
      *Data);
00143 ZESTETM1_STATUS ZestETM1SetInterrupt(ZESTETM1_CARD_INFO *CardInfo);
00144
00145 #ifdef __cplusplus
00146 }
00147 #endif
00148
00149 #endif // __ZESTETM1_H__
00150
```

# 7.33 C:/Users/hill35/git/camera\_python/nsCamera/sensors/daedalus.py File Reference

# Classes

class nsCamera.sensors.daedalus.daedalus

### **Namespaces**

- · namespace nsCamera
- namespace nsCamera.sensors
- · namespace nsCamera.sensors.daedalus

# 7.34 daedalus.py

#### Go to the documentation of this file.

```
00001 # -*- coding: utf-8 -*
00002 """
00003 Parameters and functions specific to the daedalus three-frame sensor
00004
00005
00006 Author: Jeremy Martin Hill (jerhill@llnl.gov)
00007
00008 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00009 LLNL-CODE-838080
00010
00011 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00012 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DDE)
00013 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00014 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00015 be made under this license.
00016
00017 Version: 2.1.2 (February 2025)
00018 """
00019
00020 import logging
00021 import numbers
00022 from collections import OrderedDict
00023
00024 import numpy as np
00026 from nsCamera.sensors.sensorBase import sensorBase
00027 from nsCamera.utils.misc import flattenlist
00029
00030 class daedalus(sensorBase):
00031 specwarn = ""
00032 minframe = 0 # fixed value for sensor
         maxframe = 2 # fixed value for sensor
       maxwidth = 512 # fixed value for sensor
         maxheight = 1024 # fixed value for sensor
        bytesperpixel = 2 fpganumID = 2 # last nybble of FPGA_NUM
00036
00037
00038
        detect = "DAEDALUS_DET'
          sensfam = "Daedalus"
00039
         loglabel = "[Daedalus] "
00040
00041
          ZDT = False
00042
         HFW = False
00043
          firstframe = 0
00044
          lastframe = 2
00045
          nframes = 3
          width = 512
00046
          height = 1024
00047
00048
         firstrow = 0
          lastrow = 1023
00049
```

```
interlacing = [0, 0]
00050
00051
           columns = 1
00052
           padToFull = True
           toffset = -165.76 # default temperature sensor offset
00053
00054
           tslope = 81.36  # default temperature sensor slope
00055
00056
                 _init__(self, ca):
00057
               self.caca = ca
00058
               super(daedalus, self).__init__(ca)
00059
00060
               self.sens_registers = OrderedDict(
00061
00062
                        "HST_READBACK_A_LO": "018",
                        "HST_READBACK_A_HI": "019",
00063
                        "HST_READBACK_B_LO": "01A",
"HST_READBACK_B_HI": "01B",
00064
00065
                        "HSTALLWEN_WAIT_TIME": "03F",
00066
                        "VRESET_HIGH_VALUE": "04A",
00067
                        "FRAME_ORDER_SEL": "04B",
00068
                        "EXT_PHI_CLK_SHO_ON": "050"
00069
00070
                        "EXT_PHI_CLK_SH0_OFF": "051",
00071
                        "EXT_PHI_CLK_SH1_ON": "052",
00072
                        "EXT_PHI_CLK_SH1_OFF": "053",
00073
                        "EXT_PHI_CLK_SH2_ON": "054",
                        "HST_TRIGGER_DELAY_DATA_LO": "120",
"HST_TRIGGER_DELAY_DATA_HI": "121",
00074
00075
00076
                        "HST_PHI_DELAY_DATA": "122",
                        "HST_EXT_CLK_HALF_PER": "129",
00077
                        "HST_COUNT_TRIG": "130",
"HST_DELAY_EN": "131",
00078
00079
                        "RSL_HFW_MODE_EN": "133"
00080
                        "RSL_ZDT_MODE_B_EN": "135",
00081
                        "RSL_ZDT_MODE_A_EN": "136",
"BGTRIMA": "137",
"BGTRIMB": "138",
00082
00083
00084
                        "COLUMN_TEST_EN": "139",
00085
                        "RSL_CONFIG_DATA_BO": "140",
"RSL_CONFIG_DATA_B1": "141",
00086
00087
                        "RSL_CONFIG_DATA_B2": "142"
00088
                        "RSL_CONFIG_DATA_B3": "143"
00089
                        "RSL_CONFIG_DATA_B4": "144",
00090
                        "RSL_CONFIG_DATA_B5": "145"
00091
                        "RSL_CONFIG_DATA_B6": "146"
00092
                        "RSL_CONFIG_DATA_B7": "147"
00093
                        "RSL_CONFIG_DATA_B8": "148"
00094
                        "RSL_CONFIG_DATA_B9": "149"
00095
00096
                        "RSL_CONFIG_DATA_B10": "14A"
                        "RSL_CONFIG_DATA_B11": "14B"
00097
                        "RSL_CONFIG_DATA_B12": "14C",
00098
00099
                        "RSL_CONFIG_DATA_B13": "14D"
                        "RSL_CONFIG_DATA_B14": "14E",
00100
00101
                        "RSL_CONFIG_DATA_B15": "14F",
00102
                        "RSL_CONFIG_DATA_B16": "150",
00103
                        "RSL_CONFIG_DATA_B17": "151",
00104
                        "RSL_CONFIG_DATA_B18": "152",
00105
                        "RSL_CONFIG_DATA_B19": "153",
00106
                        "RSL_CONFIG_DATA_B20": "154",
00107
                        "RSL_CONFIG_DATA_B21": "155",
00108
                        "RSL_CONFIG_DATA_B22": "156",
00109
                        "RSL_CONFIG_DATA_B23": "157",
                        "RSL_CONFIG_DATA_B24": "158",
00110
00111
                        "RSL_CONFIG_DATA_B25": "159",
00112
                        "RSL_CONFIG_DATA_B26": "15A",
                        "RSL_CONFIG_DATA_B27": "15B",
00113
                        "RSL_CONFIG_DATA_B28": "15C",
00114
                        "RSL_CONFIG_DATA_B29": "15D",
00115
                        "RSL_CONFIG_DATA_B30": "15E",
00116
00117
                        "RSL_CONFIG_DATA_B31": "15F",
                        "RSL_CONFIG_DATA_A0": "160",
00118
                        "RSL_CONFIG_DATA_A1": "161",
00119
00120
                        "RSL_CONFIG_DATA_A2": "162",
                        "RSL_CONFIG_DATA_A3": "163",
00121
                        "RSL_CONFIG_DATA_A4": "164",
00122
                        "RSL_CONFIG_DATA_A5": "165",
00123
                        "RSL_CONFIG_DATA_A6": "166",
00124
                        "RSL_CONFIG_DATA_A7": "167",
00125
                        "RSL CONFIG DATA A8": "168",
00126
                        "RSL_CONFIG_DATA_A9": "169",
00127
                        "RSL_CONFIG_DATA_A10": "16A"
00128
                        "RSL_CONFIG_DATA_A11": "16B",
00129
                        "RSL_CONFIG_DATA_A12": "16C",
00130
```

7.34 daedalus.py 379

```
00131
                                "RSL_CONFIG_DATA_A13": "16D",
                                "RSL_CONFIG_DATA_A14": "16E",
00132
00133
                                "RSL_CONFIG_DATA_A15": "16F",
00134
                                "RSL_CONFIG_DATA_A16": "170",
                                "RSL_CONFIG_DATA_A17": "171",
00135
00136
                                "RSL_CONFIG_DATA_A18": "172",
00137
                                "RSL_CONFIG_DATA_A19": "173",
00138
                                "RSL_CONFIG_DATA_A20": "174",
                                "RSL_CONFIG_DATA_A21": "175",
00139
                                "RSL_CONFIG_DATA_A22": "176",
00140
                                "RSL_CONFIG_DATA_A23": "177"
00141
                                "RSL_CONFIG_DATA_A24": "178",
00142
00143
                                "RSL_CONFIG_DATA_A25": "179",
                                "RSL_CONFIG_DATA_A26": "17A",
00144
                                "RSL_CONFIG_DATA_A27": "17B",
00145
00146
                                "RSL_CONFIG_DATA_A28": "17C",
                                "RSL_CONFIG_DATA_A29": "17D",
00147
                                "RSL_CONFIG_DATA_A30": "17E",
00148
00149
                                "RSL_CONFIG_DATA_A31": "17F",
00150
                          }
00151
                    )
00152
00153
                    self.sens subregisters = [
00154
                          ("HST_MODE", "HS_TIMING_CTL", 0, 1, True), ("SLOWREADOFF_O", "CTRL_REG", 4, 1, True), ("SLOWREADOFF_1", "CTRL_REG", 5, 1, True), ("MANSHUT_MODE", "CTRL_REG", 8, 1, True),
00156
00157
00158
00159
                           ("INTERLACING_EN", "CTRL_REG", 9, 1, True),
00160
                          ("HFW", "RSL_HFW_MODE_EN", 0, 1, True), ("2DT_A", "RSL_ZDT_MODE_A_EN", 0, 1, True), ("ZDT_B", "RSL_ZDT_MODE_B_EN", 0, 1, True), ("ZDT_B", "RSL_ZDT_MODE_B_EN", 0, 1, True),
00161
00162
00163
                          ("HST_DEL_EN", "HST_DELAY_EN", 0, 1, True),
("PHI_DELAY_A", "HST_PHI_DELAY_DATA", 9, 10, True),
("PHI_DELAY_B", "HST_PHI_DELAY_DATA", 29, 10, True),
00164
00165
00166
00167
                          # Assume that daedalus is not to be used with v1 board
00168
                          ("VRESET_HIGH", "VRESET_HIGH_VALUE", 15, 16, True),
00169
                          ("STAT_SHORISEUR", "STAT_REG", 3, 1, False), ("STAT_SHOFALLUR", "STAT_REG", 4, 1, False),
00172
00173
                           ("STAT_RSLNALLWENA", "STAT_REG", 12, 1, False), ("STAT_RSLNALLWENB", "STAT_REG", 15, 1, False),
00174
00175
                           # ("STAT_CONFIGHSTDONE", "STAT_REG", 16, 1, False),
00176
00177
00178
00179
              # TODO: add warning if daedalus and v1 board are together
00180
              def sensorSpecific(self):
00181
00182
                    Returns:
                    list of tuples, (Sensor-specific register, default setting)
00183
00184
00185
                    return [
                        ("FPA_FRAME_INITIAL", "00000000"),
("FPA_FRAME_FINAL", "00000002"),
("FPA_ROW_INITIAL", "00000000"),
("FPA_ROW_FINAL", "000003FF"),
00186
00187
00188
00189
                          ("HS_TIMING_DATA_ALO", "00006666"), # 0db6 = 2-1; 6666 = 2-2 ("HS_TIMING_DATA_AHI", "00000000"), ("HS_TIMING_DATA_BLO", "00006666"),
00190
00191
00192
                           ("HS_TIMING_DATA_BHI", "00000000"),
00193
                          ("FRAME_ORDER_SEL", "00000000"),
("RSL_HFW_MODE_EN", "00000000"),
00194
00195
                          ("RSL_ZDT_MODE_B_EN", "00000000"),
00196
                          ("RSL_ZDT_MODE_A_EN", "00000000"),
("RSL_CONFIG_DATA_BO", "00000000"),
00197
00198
                          ("RSL_CONFIG_DATA_B1", "00000000"), ("RSL_CONFIG_DATA_B2", "00000000"),
00199
00200
00201
                          ("RSL_CONFIG_DATA_B3", "00000000"),
                          ("RSL_CONFIG_DATA_B4", "00000000"), ("RSL_CONFIG_DATA_B5", "00000000"),
00202
00203
00204
                          ("RSL_CONFIG_DATA_B6", "00000000"),
                          ("RSL_CONFIG_DATA_B7", "00000000"), ("RSL_CONFIG_DATA_B8", "00000000"),
00205
00206
                          ("RSL_CONFIG_DATA_B9", "00000000"), ("RSL_CONFIG_DATA_B10", "000000000"),
00207
00208
                          ("RSL_CONFIG_DATA_B11", "00000000"),
("RSL_CONFIG_DATA_B12", "00000000"),
00209
00210
                          ("RSL_CONFIG_DATA_B13", "00000000"), ("RSL_CONFIG_DATA_B14", "00000000"),
00211
00212
                          ("RSL_CONFIG_DATA_B15", "00000000"), ("RSL_CONFIG_DATA_B16", "00000000"),
00213
00214
```

```
00215
                         ("RSL_CONFIG_DATA_B17", "00000000"),
                         ("RSL_CONFIG_DATA_B18", "00000000"), ("RSL_CONFIG_DATA_B19", "00000000"), ("RSL_CONFIG_DATA_B20", "00000000"),
00216
00217
00218
                         ("RSL_CONFIG_DATA_B21", "00000000"),
00219
                         ("RSL_CONFIG_DATA_B22", "00000000"),
00220
00221
                         ("RSL_CONFIG_DATA_B23", "00000000"),
00222
                         ("RSL_CONFIG_DATA_B24", "00000000"),
                         ("RSL_CONFIG_DATA_B25", "00000000"),
00223
                         ("RSL_CONFIG_DATA_B26", "00000000"),
("RSL_CONFIG_DATA_B26", "00000000"),
00225
                         ("RSL_CONFIG_DATA_B28", "00000000"),
00226
                         ("RSL_CONFIG_DATA_B29", "00000000"), ("RSL_CONFIG_DATA_B30", "00000000"),
00227
                         ("RSL_CONFIG_DATA_B31", "00000000"), ("RSL_CONFIG_DATA_A0", "00000000"),
00229
00230
00231
                         ("RSL_CONFIG_DATA_A1", "00000000"),
00232
                         ("RSL_CONFIG_DATA_A2", "00000000"),
00233
                         ("RSL_CONFIG_DATA_A3", "00000000"),
00234
                         ("RSL_CONFIG_DATA_A4", "00000000"),
                         ("RSL_CONFIG_DATA_A5", "00000000"), ("RSL_CONFIG_DATA_A6", "00000000"),
00235
00236
                         ("RSL_CONFIG_DATA_A7",
00237
                                                          "00000000"),
00238
                         ("RSL_CONFIG_DATA_A8",
                                                         "00000000"),
                         ("RSL_CONFIG_DATA_A9", "000000000"), ("RSL_CONFIG_DATA_A10", "00000000"), ("RSL_CONFIG_DATA_A11", "00000000"),
00239
00240
00241
                         ("RSL_CONFIG_DATA_A12", "00000000"),
00242
                         ("RSL_CONFIG_DATA_A13", "00000000"), ("RSL_CONFIG_DATA_A14", "00000000"),
00243
00244
                         ("RSL_CONFIG_DATA_A15", "00000000"), ("RSL_CONFIG_DATA_A16", "00000000"),
00245
00246
                         ("RSL_CONFIG_DATA_A17", "00000000"), ("RSL_CONFIG_DATA_A18", "000000000"),
00247
00248
                         ("RSL_CONFIG_DATA_A19", "00000000"), ("RSL_CONFIG_DATA_A20", "00000000"),
00249
00250
                         ("RSL_CONFIG_DATA_A21", "00000000"),
("RSL_CONFIG_DATA_A22", "00000000"),
00251
00252
                         ("RSL_CONFIG_DATA_A23", "00000000"), ("RSL_CONFIG_DATA_A24", "00000000"),
00253
00254
                         ("RSL_CONFIG_DATA_A24", "00000000"),
("RSL_CONFIG_DATA_A25", "00000000"),
("RSL_CONFIG_DATA_A26", "00000000"),
("RSL_CONFIG_DATA_A26", "00000000"),
("RSL_CONFIG_DATA_A28", "00000000"),
("RSL_CONFIG_DATA_A28", "00000000"),
("RSL_CONFIG_DATA_A30", "00000000"),
("RSL_CONFIG_DATA_A31", "00000000"),
("HST_TDIGGED_DETAY_DATA_10", "0000000"),
00255
00256
00257
00258
00259
00260
00261
                         ("HST_TRIGGER_DELAY_DATA_LO", "00000000"), ("HST_TRIGGER_DELAY_DATA_HI", "00000000"),
00262
00263
00264
                         ("HST_PHI_DELAY_DATA", "00000000"),
                         ("SLOWREADOFF_0", "0"), ("SLOWREADOFF_1", "0"),
00265
00266
00267
00268
00269
             def setInterlacing(self, ifactor=None, side=None):
00270
00271
                   Sets interlacing factor. NOTE: if called directly when HFW or ZDT mode is active,
00272
                   this will disengage those modes automatically. If hemispheres have different
00273
                   factors when the image is acquired, the resulting frames are separated into
00274
                   half-width images
00275
00276
                   Args:
00277
                        ifactor: number of interlaced lines (generates ifactor + 1 images per frame)
00278
                           defaults to 0 (no interlacing)
00279
                         side: identify particular hemisphere (A or B) to control. If left blank,
                           control both hemispheres
00282
                   Returns:
                   integer: active interlacing factor (unchanged if error) _{\tt mum}
00283
00284
00285
                   logging.debug(self.logdebug + "setInterlacing; ifactor = " + str(ifactor))
00286
                   if ifactor is None:
00287
                         ifactor = 0
00288
                   if (
00289
                        not isinstance(ifactor, int)
00290
                        or ifactor < 0
                        or ifactor > (self.maxheight - 1)
00291
00292
00293
                               self.logerr + "invalid interlacing factor submitted. "
00294
                               "Interlacing remains unchanged. "
00295
```

7.34 daedalus.py 381

```
00296
00297
                   logging.error(err)
00298
                   return self.interlacinginterlacing
00299
               if self.HFWHFW:
00300
                   logging.warning(
00301
                       self.logwarn + "HFW mode will be disengaged because of new "
00302
                        "interlacing setting '
00303
00304
                   self.setHighFullWellsetHighFullWell(False)
00305
               if self.ZDTZDT:
00306
                   logging.warning(
                       self.logwarn + "ZDT mode will be disengaged because of new "
00307
00308
                        "interlacing setting
00309
00310
                   self.setZeroDeadTimesetZeroDeadTime(False)
00311
               if ifactor == 0:
00312
                   bitscheme = self.maxheight * [0]
00313
                   # deactivating one side shouldn't turn off enable for both sides
00314
                   # TODO: is it a problem if sides are set separately, so interlacing is zero
00315
                      but still enabled?
00316
                   if side is None:
00317
                       self.caca.setSubregister("INTERLACING_EN", "0")
00318
               else:
00319
                  pattern = [0] + ifactor * [1]
00320
                   reps = 1 + self.maxheight // (ifactor + 1)
                   bitscheme = (reps * pattern)[0 : self.maxheight]
00321
                   self.caca.setSubregister("INTERLACING_EN", "1")
00322
               err = ""
00323
00324
              for regnum in range(32):
                   regbits = bitscheme[32 * regnum : 32 * (regnum + 1)]
00325
00326
                   logging.debug(self.logdebug + "regbits = " + str(regbits))
00327
                   \# generated pattern is reverse order from placement in register (element 0
                      of the list is the LSB of the register)
00328
00329
                   bitsrev = regbits[::-1]
                   b = "".join(s) # assemble as binary number for processing
00330
00331
                   hexval = "%x" % int(b, 2)
00332
00333
                   val = hexval.zfill(8)
00334
                   err0 = ""
                   err1 = ""
00335
                   if side is None or side.lower() == "a":
   lname = "RSL_CONFIG_DATA_A" + str(regnum)
00336
00337
00338
                       err1, _ = self.caca.setRegister(lname, val)
00339
                       self.interlacinginterlacing[1] = ifactor
                   if side is None or side.lower() == "b":
    rname = "RSL_CONFIG_DATA_B" + str(regnum)
00340
00341
00342
                       err0, _ = self.caca.setRegister(rname, val)
00343
                       self.interlacinginterlacing[0] = ifactor
00344
                   err = err + err0 + err1
00345
               logging.error(self.logerr + "interlacing may not be set correctly: " + err) logging.info(self.loginfo + "Interlacing set to " + str(self.interlacinginterlacing))
00346
00347
00348
               if self.interlacinginterlacing[0] == self.interlacinginterlacing[1]:
00349
                   self.columnscolumns = 1
00350
               else:
00351
                   self.columnscolumns = 2
00352
               return self.interlacinginterlacing
00353
00354
          def setHighFullWell(self, flag):
00355
               Activates High Full Well mode. All frames are acquired simultaneously. Zero Dead
00356
00357
                 Time mode and interlacing will be automatically deactivated and column number
00358
                 will be reset to 0. NOTE: after deactivating HFW, the board remains in
00359
                 uninterlaced mode (interlacing = 0)
00360
00361
               Args:
00362
                   flag: True to activate HFW mode, False to deactivate
00363
00364
               Returns:
               Error message
00365
00366
00367
               logging.debug(self.logdebug + "setHighFullWell; flag = " + str(flag))
               err0 = ""
00368
00369
               if flag:
                   if self.ZDTZDT:
00370
00371
                       logging.warning(
00372
                            self.logwarn + "ZDT mode will be disengaged because of HFW "
00373
                            "setting "
00374
00375
                       err0 = self.setZeroDeadTimesetZeroDeadTime(False)
00376
                   err1, _ = self.caca.setSubregister("HFW", "1")
```

```
00377
                  self.HFWHFW = False # preclude HFW deactivation message in setInterlacing
00378
                  self.setInterlacingsetInterlacing(0)
00379
                  self.HFWHFW = True
00380
                  logging.info(self.loginfo + "High Full Well mode active")
00381
              else:
00382
                  self.HFWHFW = False
00383
                  err1, _ = self.caca.setSubregister("HFW", "0")
00384
                  logging.info(self.loginfo + "High Full Well mode inactivate")
00385
              err = err0 + err1
00386
00387
                  logging.error(self.logerr + "HFW option may not be set correctly ")
00388
              return err
00389
00390
          def setZeroDeadTime(self, flag=True, side=None):
00391
00392
              Activates Zero Dead Time mode. Even rows follow the assigned HST schedule; odd
                rows are acquired while the 'shutter' for the even rows are closed. High Full
00393
                Well mode and interlacing will be automatically deactivated.
00394
00395
              *NOTE* after deactivating ZDT, the board reverts to uninterlaced mode
00396
                (interlacing = 0)
00397
00398
              Args:
00399
                  flag: True to activate ZDT mode, False to deactivate
00400
                  side: identify particular hemisphere (A or B) to control. If left blank,
00401
                    control both hemispheres
00402
00403
              Returns:
              Error message
00404
00405
00406
              logging.debug(self.logdebug + "setZeroDeadTime; flag = " + str(flag))
00407
              err0 = '
              err1 = ""
00408
              err2 = ""
00409
00410
              if flag:
                  if self.HFWHFW:
00411
00412
                      logging.warning(
                           self.logwarn + "HFW mode will be disengaged because of ZDT " "setting "
00413
00414
00415
00416
                      err0 = self.setHighFullWellsetHighFullWell(False)
00417
                  if side is None or side.lower() == "a":
00418
                      err2, _ = self.caca.setSubregister("ZDT_A", "1")
                  self.interlacinginterlacing[0] = 1
if side is None or side.lower() == "b":
00419
00420
00421
                      err1, _ = self.caca.setSubregister("ZDT_B", "1")
00422
                      self.interlacinginterlacing[1] = 1
00423
                  # self.ZDT = False # preclude ZDT deactivation message in setInterlacing
00424
                  # if self.interlacing != [0, 0]:
00425
                        self.setInterlacing(0)
00426
                  # TODO: need to handle flags when ZDT active for just one side
00427
                  self.ZDTZDT = True
00428
                  logging.info(
00429
                      self.loginfo + "Zero Dead Time mode active; actual interlacing = 1"
00430
00431
              else:
00432
                  self.ZDTZDT = False
00433
                  if side is None or side.lower() == "a":
00434
                      err2, _ = self.caca.setSubregister("ZDT_A", "0")
00435
                  if side is None or side.lower() == "b":
00436
                      err1, _ = self.caca.setSubregister("ZDT_B", "0")
                  self.setInterlacingsetInterlacing(0)
00437
00438
                  logging.info(self.loginfo + "Zero Dead Time mode inactivate")
00439
              err = err0 + err1 + err2
              if err:
00440
00441
                  logging.error(self.logerr + "ZDT option may not be set correctly ")
00442
              return err
00443
00444
          def selectOscillator(self, osc=None):
00445
00446
              Selects oscillator to control sensor timing
00447
              Args:
00448
                  osc: 500|100|'ring'|external', defaults to 500 MHz
00449
00450
              Returns:
              error message as string
00451
00452
              logging.info(self.loginfo + "selectOscillator; osc = " + str(osc))
00453
00454
              if osc is None:
00455
                 osc = 500
00456
              osc = str(osc)
              if osc[:3] == "500":
00457
```

7.34 daedalus.py 383

```
payload = "00"
00458
00459
               elif osc[:3] == "100":
                   payload = "01"
00460
00461
               elif osc.upper()[:3] == "RIN":
                   payload = "10"
00462
00463
               elif osc.upper()[:3] in ["EXT"]:
00464
                   payload = "11"
00465
00466
                   err = (
                        self.logerr + "selectOscillator: invalid parameter supplied. "
00467
00468
                         "Oscillator selection is unchanged."
00469
00470
                    logging.error(err)
00471
                    return err
00472
               self.caca.setSubregister("OSC_SELECT", payload)
00473
00474
           def setTriggerDelay(self, delay=0):
00475
00476
               Use trigger delay timer. Actual delay is rounded down to multiple of .15 ns, up
00477
                 to a maximum delay of 6 ns
00478
00479
00480
                   delay: trigger delay in ns
00481
00482
               Returns:
               String of errors, if any
00483
00484
00485
               logging.debug(self.logdebug + "setTriggerDelay; delay = " + str(delay))
00486
                    not (isinstance(delay, int) or isinstance(delay, float))
00487
00488
                   or delay < 0
00489
                   or delay > 6
00490
               ) •
00491
                    err = (
                        self.logerr + "invalid trigger delay submitted. Delay remains "
00492
                         "unchanged. "
00493
00494
00495
                    logging.error(err)
00496
                    return err
00497
               delayblocks = int(delay / 0.15)
00498
               if delayblocks < 0:</pre>
00499
                   delayblocks = 0
00500
               if delayblocks > 40:
00501
                   delayblocks = 40
               \begin{array}{l} \mbox{delayseq} = (40 - \mbox{delayblocks}) \; \star \; [0] \; + \; \mbox{delayblocks} \; \star \; [1] \\ \mbox{seqstr} = \; "". \mbox{join}(\mbox{str}(\mbox{x}) \; \; \mbox{for} \; \mbox{x} \; \; \mbox{in} \; \mbox{delayseq}) \end{array}
00502
00503
               seqhex = "%x" % int(seqstr, 2)
logging.debug(self.logdebug + "seqhex = " + str(seqhex))
00504
00505
00506
               highpart = seqhex[-10:-8].zfill(8)
00507
               lowpart = seqhex[-8:].zfill(8)
00508
               self.caca.setSubregister("HST_DEL_EN", "1")
               err0, _ = self.caca.setRegister("HST_TRIGGER_DELAY_DATA_LO", lowpart)
00509
               err1, _ = self.caca.setRegister("HST_TRIGGER_DELAY_DATA_HI", highpart)
00510
00511
               err2, _ = self.caca.setSubregister("HST_MODE", "1")
00512
               delayed = delayblocks * 0.15
00513
               logging.info(self.loginfo + "Actual trigger delay = " + str(delayed) + " ns")
00514
               return err0 + err1 + err2
00515
00516
           def setPhiDelay(self, side=None, delay=0):
00517
00518
               Use phi delay timer. Actual delay is rounded down to multiple of .15 ns, up to a
00519
                 maximum delay of 1.5 ns
00520
00521
                   side: hemisphere to delay; if None, delay both hemispheres
00522
                   delay: phi delay in ns
00523
00524
               String of errors, if any
00525
00526
00527
               logging.debug(self.logdebug + "setPhiDelay; delay = " + str(delay))
00528
00529
                   not (isinstance(delay, int) or isinstance(delay, float))
00530
                   or delay < 0
00531
                   or delay > 1.5
00532
00533
                    err = (
                        self.logerr + "invalid phi delay submitted. Delay remains "
00534
00535
                         "unchanged.
00536
00537
                    logging.error(err)
00538
                    return err
```

```
00539
               delayblocks = int(delay / 0.15)
00540
               if delayblocks < 0:</pre>
00541
                   delayblocks = 0
00542
               if delayblocks > 10:
00543
                   delayblocks = 10
               delayseq = (10 - delayblocks) * [0] + delayblocks * [1]
seqstr = "".join(str(x) for x in delayseq)
00544
00545
               err1 = ""
00546
00547
               err2 = ""
               if side is None or side.upper() == "A":
00548
00549
                   err1, _ = self.caca.setSubregister("PHI_DELAY_A", seqstr)
00550
               if side is None or side.upper() == "B":
               err2, _ = self.caca.setSubregister("PHI_DELAY_B", seqstr) delayed = delayblocks \star 0.15
00551
00552
00553
               logging.info(self.loginfo + "Actual phi delay = " + str(delayed) + " ns")
00554
               return err1 + err2
00555
00556
          def setExtClk(self, dilation=None, frequency=None):
00557
00558
               Override the standard board clock with the external clock.
00559
               Args:
00560
                   dilation: ratio of base frequency (500 MHz) to desired external clock
                     frequency. Default is 25. Overridden if frequency parameter is provided
00561
00562
                   frequency: Desired frequency for phi clock.
00563
               Returns:
               error message as string
00564
00565
00566
               logging.debug(
00567
                   self.logdebug
                   + "setExtClk; dilation = "
00568
00569
                   + str(dilation)
00570
                   + "; frequency = "
00571
                   + str(frequency)
00572
               if not (isinstance(frequency, int) or isinstance(frequency, float)):
00573
00574
                   err = (
00575
                       self.logerr
00576
                        + "invalid external clock frequency submitted. Clock is not "
                        + "operating"
00577
00578
00579
                   logging.error(err)
00580
                   return err
00581
               self.caca.selectOscillator("external")
00582
               if not dilation:
00583
                   dilation = 25
00584
               if not frequency:
               frequency = 5e7 / float(dilation)
count = 2e7 / float(frequency) - 1 # base phi clock is 20 MHz?
00585
00586
00587
               if count < 0:</pre>
00588
                   count = 0
00589
                   warn = (
00590
                       self.logwarn
00591
                        + "external clock frequency exceeding maximum. Frequency set to "
00592
                        + "maximum (20 MHz)"
00593
00594
                   logging.warning(warn)
00595
               if count > 0xFFFFFFFF:
00596
                   count = 0xFFFFFFF
00597
               counthex = hex(int(count))[2:].zfill(8)
00598
               self.caca.setRegister("HST_EXT_CLK_HALF_PER", counthex)
00599
00600
           # TODO: enable exponential form for all large number inputs (accept floats)
00601
00602
           def setManualShutters(self, timing=None):
00603
00604
               Legacy alias for setManualTiming()
00605
00606
               self.setManualTimingsetManualTiming(timing)
00607
00608
           def setManualTiming(self, timing=None):
00609
00610
               Manual shutter timing, five intervals given in nanoseconds, e.g., [100,50,100,50,100] for frame 0 open for 100 ns, an interframe pause of 50 ns,
00611
                  frame 1 open for 100 ns, etc. Timing is set for both hemispheres.
00612
00613
00614
               The actual timing is rounded down to the nearest multiple of 25\ \mathrm{ns.} (Each
                 count = 25 ns. e.g., a request for 140 ns rounds down to a count of '5',
00615
00616
                 which corresponds to 125 ns))
                   - Minimum timing is 75 ns
00617
                   - Maximum is 25 \times 2^30 ns (approximately 27 seconds)
00618
00619
```

7.34 daedalus.py 385

```
00620
00621
                     timing: 5-element list in nanoseconds
00622
00623
                Returns:
                tuple (error string, response string from final message)
00624
00625
00626
                if timing is None:
00627
                     logging.info(
00628
                         self.loginfo
                          + "No manual timing setting provided, defaulting to (100, 150, 100, "
00629
00630
                          " 150, 100, 150, 100) for both hemispheres"
00631
00632
                     timing = [(100, 150, 100, 150, 100)]
00633
00634
                logging.info(self.loginfo + "Manual shutter sequence: " + str(timing))
00635
                flattened = flattenlist(timing)
00636
00637
                     len(flattened) != 5
00638
                     or not all(isinstance(x, (int, float)) for x in flattened)
                     or not all (x \ge 25 \text{ for } x \text{ in flattened})
00639
00640
00641
                     err = self.logerr + "Invalid manual shutter timing list: " + str(timing)
                     logging.error(err + "; timing settings unchanged") return err, "00000000"
00642
00643
00644
00645
                timecounts = [int(a // 25) for a in flattened]
                self.caca.sensmanual = timing
self.caca.senstiming = {} # clear HST settings from ca object
00646
00647
00648
00649
                control_messages = [
                     ("MANSHUT_MODE", "1"),
00650
                     ("EXT_PHI_CLK_SH0_ON", "{0:#0{1}x}".format(timecounts[0], 10)[2:10]),
("EXT_PHI_CLK_SH0_OFF", "{0:#0{1}x}".format(timecounts[1], 10)[2:10]),
("EXT_PHI_CLK_SH1_ON", "{0:#0{1}x}".format(timecounts[2], 10)[2:10]),
("EXT_PHI_CLK_SH1_OFF", "{0:#0{1}x}".format(timecounts[3], 10)[2:10]),
("EXT_PHI_CLK_SH2_ON", "{0:#0{1}x}".format(timecounts[4], 10)[2:10]),
00651
00652
00653
00654
00655
00656
00657
                return self.caca.submitMessages(control_messages, " setManualShutters: ")
00658
           def getManualTiming(self):
00659
00660
00661
                Read off manual shutter timing settings
00662
                Returns:
                list of manual timing intervals
00663
00664
00665
                timing = []
00666
                for reg in [
00667
                     "EXT_PHI_CLK_SH0_ON",
00668
                     "EXT_PHI_CLK_SH0_OFF",
00669
                     "EXT_PHI_CLK_SH1_ON",
00670
                     "EXT_PHI_CLK_SH1_OFF",
00671
                     "EXT_PHI_CLK_SH2_ON",
00672
00673
                     _, reghex = self.caca.getRegister(reg)
00674
                     timing.append(25 * int(reghex, 16))
00675
00676
00677
           def getSensTemp(self, scale=None, offset=None, slope=None, dec=1):
00678
00679
                Read temperature sensor located on the Daedalus sensor
00680
                Args:
00681
                     scale: temperature scale to report (defaults to C, options are F and K)
00682
                     offset: offset of linear fit of sensor response (defaults to self.toffset)
00683
                     slope: slope of linear fit of sensor response (defaults to self.tslope)
                     dec: round to 'dec' digits after the decimal point
00684
00685
00686
                temperature as float on given scale, rounded to .1 degree _{\tt u,u}
00687
00688
00689
                err, rval = self.caca.getMonV("MON_TSENSE_OUT", errflag=True)
00690
                if err:
00691
                     logging.error(
00692
                         self.logerr + "unable to retrieve temperature information ("
                          'getTemp), returning "0" '
00693
00694
00695
                     return 0.0
00696
                if offset is None:
00697
                     offset = self.toffset
00698
                if slope is None:
00699
                     slope = self.tslope
00700
```

```
00701
              ctemp = offset + slope * rval
00702
              if scale == "K":
00703
                   temp = round(ctemp + 273.15, dec)
00704
               elif scale == "F":
00705
                   temp = round(1.8 \star ctemp + 32, dec)
00706
00707
                  temp = round(ctemp, dec)
00708
               return temp
00709
00710
          def parseReadoff(self, frames, columns):
00711
00712
              Parses frames from board into images
00713
              Args:
00714
                   frames: list of data arrays (frames) returned from board
00715
                   columns: 1 (full width image) or 2 (hemispheres generate distinct images)
00716
              list of data arrays (frames) reordered and deinterlaced
00717
00718
00719
              logging.debug(self.logdebug + "parseReadoff")
00720
              w = self.width
               if hasattr(self, "ca"): # TODO: this may no longer be necessary
00721
00722
                   padIt = self.caca.padToFull
00723
              else:
00724
                  padIt = self.padToFull
00725
               if padIt:
00726
                  rows = self.maxheight
00727
              else:
00728
                 rows = self.lastrow - self.firstrow + 1
00729
              parsed = []
00730
              for frame in frames:
                   current = np.zeros((rows, w), dtype=np.uint16)
mapped = np.zeros((rows, w), dtype=np.uint16)
00731
00732
00733
                   frame = frame.reshape(rows, w)
00734
                   for entry in range(int(w / 2)): col = 32 \star (entry % 8) + entry // 8 # lookup from daedlookup.xls
00735
00736
00737
                       for row in range(rows):
00738
                           current[row][col] = frame[row][2 * entry]
00739
                           current[row][col + 256] = frame[row][2 * entry + 1]
00740
00741
                   for row in range(rows):
00742
                       mapped[row][0:32] = current[row][320:352]
                       mapped[row][32:64] = current[row][352:384]
mapped[row][64:96] = current[row][192:224]
00743
00744
                       mapped[row][96:128] = current[row][160:192]
00745
00746
                       mapped[row][128:160] = current[row][256:288]
                       mapped[row][160:192] = current[row][288:320]
00747
                       mapped[row][192:224] = current[row][416:448]
00748
00749
                       mapped[row][224:256] = current[row][32:64]
00750
                       mapped[row][256:288] = current[row][128:160]
                       mapped[row][288:320] = current[row][224:256]
00751
00752
                       mapped[row][320:352] = current[row][384:416]
                       mapped[row][352:384] = current[row][448:480]
00753
00754
                       mapped[row][384:416] = current[row][480:512]
00755
                       mapped[row][416:448] = current[row][0:32]
00756
                       mapped[row][448:480] = current[row][64:96]
00757
                       mapped[row][480:512] = current[row][96:128]
00758
                   parsed.append(mapped)
00759
00760
               images = self.caca.partition(parsed, columns)
00761
               flatimages = [flattenlist(x) for x in images]
               return flatimages
00762
00763
00764
          def reportStatusSensor(self, statusbits, statusbits2):
00765
00766
              Print status messages from sensor-specific bits of status register or object
00767
                status flags
00768
00769
              Args:
00770
                   statusbits: result of checkStatus()
00771
                   statusbits2: result of checkStatus2()
00772
00773
              if int(statusbits[3]):
00774
                   print(self.loginfo + "SHO_rise_B_edge detected")
00775
              if int(statusbits[4]):
00776
                   print(self.loginfo + "SHO_fall_B_edge detected")
00777
               if int(statusbits[12]):
00778
                   print(self.loginfo + "RSLNALLWENB detected")
00779
              if int(statusbits[15]):
                   print(self.loginfo + "RSLNALLWENA detected")
00780
00781
              if self.HFWHFW:
```

```
print(self.loginfo + "High Full Well mode active")
              # TODO: handle two hemispheres for ZDT
00784
             elif self.ZDTZDT:
                 print(self.loginfo + "Zero Dead Time mode active")
00785
00786
             elif self.interlacinginterlacing != [0, 0]:
00787
00788
                      "{loginfo}Interlacing active: {interlacing}".format(
00789
                          loginfo=self.loginfo, interlacing=str(self.interlacinginterlacing)
00790
00791
00792
              if self.caca.sensmanual == []:
00793
                 print(
                      "{loginfo}High-speed timing: A:{Atiming}, B:{Btiming}".format(
00795
                          loginfo=self.loginfo,
                          Atiming=self.getTiming(side="A", actual=True),
00797
                          Btiming=self.getTiming(side="B", actual=True),
00798
                     )
00799
                 )
00800
             else:
00801
                print(
00802
                      "{loginfo}Manual timing set to {timing}".format(
                          loginfo=self.loginfo, timing=self.getManualTiminggetManualTiming()
00803
00804
                     )
00805
00806
00807
00808 """
00809 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00810 LLNL-CODE-838080
00812 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00813 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00814 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00815 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00816 be made under this license.
00817 """
```

# 7.35 C:/Users/hill35/git/camera\_python/nsCamera/sensors/icarus.py File Reference

# **Classes**

class nsCamera.sensors.icarus.icarus

# **Namespaces**

- namespace nsCamera
- namespace nsCamera.sensors
- namespace nsCamera.sensors.icarus

# 7.36 icarus.py

### Go to the documentation of this file.

```
00001 # -*- coding: utf-8 -*-
00002 """
00003 Parameters and functions specific to the icarus two-frame sensor
00004
00005 ***Do not use this file as a template for new code development***
00006
00007 Author: Jeremy Martin Hill (jerhill@llnl.gov)
00008
00009 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
```

```
00010 LLNL-CODE-838080
00011
00012 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00013 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00014 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00015 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00016 be made under this license.
00017
00018 Version: 2.1.2 (February 2025)
00019 """
00020
00021 import logging
00022 from collections import OrderedDict
00024 from nsCamera.sensors.sensorBase import sensorBase
00026
00027 class icarus(sensorBase):
00028
           specwarn = " and the use of the Icarus model 1 sensor"
           minframe = 1 # fixed value for sensor
00029
00030
           maxframe = 2 # fixed value for sensor
00031
           # WARNING: the camera will always 'acquire' four frames, but will only generate
          # images for the middle two; HST and manual shutters will manage all four # frames
00032
00033
           maxwidth = 512 # fixed value for sensor
00034
           maxheight = 1024 # fixed value for sensor
00035
           bytesperpixel = 2
00036
           icarustype = 1 # 2-frame version
00037
           fpganumID = 1  # last nybble of FPGA_NUM detect = "ICARUS_DET"
00038
00039
           sensfam = "Icarus"
00040
           loglabel = "[Icarus1] "
00041
00042
           firstframe = 1
00043
           lastframe = 2
           nframes = 2
00044
           width = 512
00045
           height = 1024
00046
00047
           firstrow = 0
00048
           lastrow = 1023
           interlacing = [0, 0] # N/A for icarus
00049
00050
           columns = 1
00051
           padToFull = True
00052
00053
           def __init__(self, ca):
00054
                self.caca = ca
00055
                super(icarus, self).__init__(ca)
00056
00057
                self.sens_registers = OrderedDict(
00058
00059
                         "VRESET_WAIT_TIME": "03E",
00060
                         "ICARUS_VER_SEL": "041",
00061
                         "VRESET_HIGH_VALUE": "04A"
                         "MISC_SENSOR_CTL": "04C",
"MANUAL_SHUTTERS_MODE": "050",
00062
00063
00064
                         "W0_INTEGRATION": "051",
00065
                         "WO_INTERFRAME": "052",
                         "W1_INTEGRATION": "053",
00066
                         "W1_INTERFRAME": "054",
"W2_INTEGRATION": "055",
00067
00068
00069
                         "W2_INTERFRAME": "056",
                         "W3_INTEGRATION": "057"
00070
00071
                         "W0_INTEGRATION_B": "058",
00072
                         "WO_INTERFRAME_B": "059",
                         "W1_INTEGRATION_B": "05A",
00073
00074
                         "W1_INTERFRAME_B": "05B",
00075
                         "W2_INTEGRATION_B": "05C",
00076
                         "W2_INTERFRAME_B": "05D"
                         "W3_INTEGRATION_B": "05E",
00077
00078
                         "TIME_ROW_DCD": "05F",
00079
                    }
00080
00081
00082
                self.sens_subregisters = [
00083
00085
                     ("MANSHUT_MODE", "MANUAL_SHUTTERS_MODE", 0, 1, True),
                     ("REVREAD", "CTRL_REG", 4, 1, True),
00086
                    ("ROWDEAD", "CIRL_REG", 4, 1, 1rue),

("PDBIAS_LOW", "CTRL_REG", 6, 1, True),

("ROWDCD_CTL", "CTRL_REG", 7, 1, True),

("ACCUMULATION_CTL", "MISC_SENSOR_CTL", 0, 1, True),

("HST_TST_ANRST_EN", "MISC_SENSOR_CTL", 1, 1, True),

("HST_TST_BNRST_EN", "MISC_SENSOR_CTL", 2, 1, True),
00087
00088
00089
00090
00091
```

7.36 icarus.py 389

```
00092
                       ("HST_TST_ANRST_IN", "MISC_SENSOR_CTL", 3, 1, True),
                      ("HST_TST_BNRST_IN", "MISC_SENSOR_CTL", 4, 1, True),
("HST_PXL_RST_EN", "MISC_SENSOR_CTL", 5, 1, True),
("HST_CONT_MODE", "MISC_SENSOR_CTL", 6, 1, True),
00093
00094
00095
                      ("COL_DCD_EN", "MISC_SENSOR_CTL", 7, 1, True), ("COL_READOUT_EN", "MISC_SENSOR_CTL", 8, 1, True),
00096
00097
00098
                       ("READOFF_DELAY_EN", "TRIGGER_CTL", 4, 1, True),
00099
                      ("STAT_W3TOPAEDGE1", "STAT_REG", 3, 1, False),
("STAT_W3TOPBEDGE1", "STAT_REG", 4, 1, False),
("STAT_HST_ALL_W_EN_DETECTED", "STAT_REG", 12, 1, False),
("PDBIAS_UNREADY", "STAT_REG2", 5, 1, False),
00102
00103
00104
00105
00106
                 1
00108
                 if self.caca.boardname == "llnl_v1":
00109
                      self.sens_subregisters.append(
00110
                           ("VRESET_HIGH", "VRESET_HIGH_VALUE", 7, 8, True)
00111
                     )
00112
                 else:
00113
                      self.sens_subregisters.extend(
00114
00115
                                ("VRESET_HIGH", "VRESET_HIGH_VALUE", 15, 16, True),
00116
                                ("READOFF_DELAY_EN", "TRIGGER_CTL", 4, 1, True),
00117
00118
00119
                      self.sens registers.update({"DELAY ASSERTION ROWDCD EN": "04F"})
00120
            def checkSensorVoltStat(self):
00121
00122
00123
                 Checks register tied to sensor select jumpers to confirm match with sensor
00124
                   object
00125
00126
                 Returns:
                 boolean, True if jumpers select for Icarus sensor
00127
00128
                 logging.debug(self.logdebug + "checkSensorVoltStat")
00129
00130
                 err, status = self.caca.getSubregister("ICARUS_DET")
00131
00132
                      logging.error(self.logerr + "unable to confirm sensor status")
00133
                      return False
00134
                 if not int(status):
00135
                      logging.error(self.logerr + "Icarus sensor not detected")
00136
                      return False
00137
                 return True
00138
00139
            def sensorSpecific(self):
00140
00141
                 Returns:
                 list of tuples, (Sensor-specific register, default setting)
00142
00143
00144
                 icarussettings = [
                      ("ICARUS_VER_SEL", "00000001"),
00145
                      ("FPA_FRAME_INITIAL", "00000001"),
("FPA_FRAME_FINAL", "00000002"),
("FPA_ROW_INITIAL", "00000000"),
("FPA_ROW_FINAL", "000003FF"),
00146
00147
00148
00149
                      ("VRESET_WAIT_TIME", "000927C0"),
("HS_TIMING_DATA_BHI", "00000000"),
("HS_TIMING_DATA_BLO", "00006666"),
("HS_TIMING_DATA_AHI", "00000000"),
("HS_TIMING_DATA_ALO", "00006666"),
00150
00151
00152
                                                                   \# 0db6 = 2-1; 6666 = 2-2
00153
00154
00155
00156
                 if self.caca.boardname == "llnl_v1":
00157
                      icarussettings.append(
                           ("VRESET_HIGH_VALUE", "000000D5") # 3.3 V (FF = 3.96)
00158
00159
00160
                 else:
00161
                     icarussettings.append(("VRESET_HIGH_VALUE", "0000FFFF"))
00162
                 return icarussettings
00163
00164
00165 """
00166 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00167 LLNL-CODE-838080
00168
00169 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00170 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00171 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00172 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00173 be made under this license.
00174 """
```

# 7.37 C:/Users/hill35/git/camera\_python/nsCamera/sensors/icarus2.py File Reference

# Classes

class nsCamera.sensors.icarus2.icarus2

### **Namespaces**

- · namespace nsCamera
- namespace nsCamera.sensors
- · namespace nsCamera.sensors.icarus2

# 7.38 icarus2.py

#### Go to the documentation of this file.

```
00001 # -*- coding: utf-8 -*
00002 """
00003 Parameters and functions specific to the four-frame icarus2 sensor
00004
00005 Author: Jeremy Martin Hill (jerhill@llnl.gov)
00006
00007 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00008 LLNL-CODE-838080
00009
00010 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00011 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00012 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00013 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00014 be made under this license.
00015
00016 Version: 2.1.2 (February 2025) 00017 """
00018
00019 from collections import OrderedDict
00020
00021 from nsCamera.sensors.sensorBase import sensorBase
00022
00023
00024 class icarus2(sensorBase):
00025 specwarn = ""
00026
         minframe = 0 # fixed value for sensor
00027
         maxframe = 3 # fixed value for sensor
       maxwidth = 512 # fixed value for sensor
00029
         maxheight = 1024 # fixed value for sensor
00030
         bytesperpixel = 2
00031
          icarustype = 0 # 4-frame version
         fpganumID = 1 # last nybble of FPGA_NUM
00033
          detect = "ICARUS_DET"
00034
          sensfam = "Icarus"
00035
          loglabel = "[Icarus2] "
          firstframe = 0
00036
          lastframe = 3
00037
00038
          nframes = 4
00039
          width = 512
          height = 1024
00040
00041
          firstrow = 0
00042
          lastrow = 1023
00043
          interlacing = [0, 0] # N/A for icarus
00044
          columns = 1
         padToFull = True
00045
00046
00047
          def __init__(self, ca):
00048
              self.caca = ca
00049
              super(icarus2, self).__init__(ca)
```

7.38 icarus2.py 391

```
00050
00051
                    self.sens_registers = OrderedDict(
00052
00053
                                "VRESET_WAIT_TIME": "03E",
00054
                                "ICARUS_VER_SEL": "041",
00055
                                "MISC_SENSOR_CTL": "04C"
00056
                                "MANUAL_SHUTTERS_MODE": "050",
00057
                                "W0_INTEGRATION": "051",
00058
                                "WO_INTERFRAME": "052",
                                "W1_INTEGRATION": "053",
                                "W1_INTERFRAME": "054",
00060
                                "W2_INTEGRATION": "055",
00061
                               "W2_INTERFRAME": "056",
"W3_INTEGRATION": "057"
00062
00063
00064
                                "WO_INTEGRATION_B": "058",
00065
                                "WO_INTERFRAME_B": "059",
00066
                                "W1_INTEGRATION_B": "05A",
                                "W1_INTERFRAME_B": "05B",
00067
                                "W2_INTEGRATION_B": "05C",
00068
                                "W2_INTERFRAME_B": "05D",
00069
                                "W3_INTEGRATION_B": "05E",
00070
00071
                                "TIME_ROW_DCD": "05F",
00072
                         }
00073
                   )
00074
00075
                    self.sens subregisters = [
00076
00078
                          ("MANSHUT_MODE", "MANUAL_SHUTTERS_MODE", 0, 1, True),
                         ("MANSHUT_MODE", "MANUAL_SHUTTERS_MODE", U, 1, True)
("REVREAD", "CTRL_REG", 4, 1, True),
("PDBIAS_LOW", "CTRL_REG", 6, 1, True),
("ROWDCD_CTL", "CTRL_REG", 7, 1, True),
("ACCUMULATION_CTL", "MISC_SENSOR_CTL", 0, 1, True),
("HST_TST_ANRST_EN", "MISC_SENSOR_CTL", 1, 1, True),
("HST_TST_BNRST_EN", "MISC_SENSOR_CTL", 2, 1, True),
("HST_TST_ANRST_IN", "MISC_SENSOR_CTL", 3, 1, True),
("HST_TST_BNRST_IN", "MISC_SENSOR_CTL", 4, 1, True),
("HST_TST_BNRST_IN", "MISC_SENSOR_CTL", 4, 1, True),
("HST_TST_BNRST_IN", "MISC_SENSOR_CTL", 5, 1, True),
00079
00080
00081
00082
00083
00084
00085
00086
                          ("HST_PXL_RST_EN", "MISC_SENSOR_CTL", 5, 1, True), ("HST_CONT_MODE", "MISC_SENSOR_CTL", 6, 1, True),
00087
00088
                          ("COL_DCD_EN", "MISC_SENSOR_CTL", 7, True), ("COL_READOUT_EN", "MISC_SENSOR_CTL", 8, 1, True),
00089
00090
00091
                          ("STAT_W3TOPAEDGE1", "STAT_REG", 3, 1, False),
("STAT_W3TOPBEDGE1", "STAT_REG", 4, 1, False),
("STAT_HST_ALL_W_EN_DETECTED", "STAT_REG", 12, 1, False),
00094
00095
00096
                          ("PDBIAS_UNREADY", "STAT_REG2", 5, 1, False),
00097
00098
00099
                    if self.caca.boardname == "llnl_v4":
00100
                         self.sens_subregisters.append(
    ("READOFF_DELAY_EN", "TRIGGER_CTL", 4, 1, True)
00101
00102
00103
00104
                          self.sens_registers.update({"DELAY_ASSERTION_ROWDCD_EN": "04F"})
00105
00106
              # TODO: clean up static methods
00107
              def sensorSpecific(self):
00108
00109
                    list of tuples, (Sensor-specific register, default setting)
00110
00111
00112
                    return [
                       ("ICARUS_VER_SEL", "00000000"),
00113
                         ("FPA_FRAME_INITIAL", "00000000"),
("FPA_FRAME_FINAL", "00000003"),
("FPA_ROW_INITIAL", "00000000"),
("FPA_ROW_FINAL", "000003FF"),
00114
00115
00116
                          ("HS_TIMING_DATA_BHI", "00000000"),
("HS_TIMING_DATA_BLO", "00006666"), # 0db6 = 2-1; 6666 = 2-2
00118
00119
                          ("HS_TIMING_DATA_AHI", "00000000"),
("HS_TIMING_DATA_ALO", "00006666"),
00120
00121
00122
                   ]
00123
00124
00125 """
00126 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00127 LLNL-CODE-838080
00128
00129 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00130 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00131 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00132 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00133 be made under this license.
```

00134 """

# 7.39 C:/Users/hill35/git/camera\_python/nsCamera/sensors/sensorBase.py File Reference

# Classes

· class nsCamera.sensors.sensorBase.sensorBase

# **Namespaces**

- namespace nsCamera
- namespace nsCamera.sensors
- namespace nsCamera.sensors.sensorBase

# 7.40 sensorBase.py

# Go to the documentation of this file.

```
00001 # -*- coding: utf-8 -*-
00003 Superclass for nsCamera sensors
00004
00005 Author: Jeremy Martin Hill (jerhill@llnl.gov)
00006
00007 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00008 LLNL-CODE-838080
00009
00010 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00011 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00012 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00013 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00014 be made under this license.
00015
00016 Version: 2.1.2 (February 2025)
00017 """
00018 import itertools
00019 import logging
00020 import numbers
00021
00022 from nsCamera.utils.misc import flattenlist, makeLogLabels
00024
00025 class sensorBase(object):
00026
00027
          Base class for sensors. 'Virtual' methods below default to Icarus behavior.
00028
          daedalus.py overrides these methods as necessary
00030
00031
          def __init__(self, camassem):
00032
              self.ca = camassem
00033
              # skip board settings if no board object exists
              if hasattr(self.ca, "board"):
00034
                  self.init_board_specific()
00035
00036
00037
                  self.logcrit,
00038
00039
                  self.logerr,
00040
                  self.logwarn.
00041
                  self.loginfo,
                  self.logdebug,
00042
00043
              ) = makeLogLabels(self.ca.logtag, self.loglabel)
00044
```

7.40 sensorBase.py 393

```
00045
              # skip assignment if no comms object exists
00046
              if hasattr(self.ca, "comms"):
00047
                  self.ca.comms.payloadsize = (
00048
                      self.width * self.height * self.nframes * self.bytesperpixel
00049
00050
00051
              logging.info(self.loginfo + "Initializing sensor object")
00052
00053
          def init_board_specific(self):
00054
              """Initialize aliases and subregisters specific to the current board and sensor."""
00055
00056
              if self.ca.sensorname == "icarus" or self.ca.sensorname == "icarus2":
00057
                  self.ca.board.subreg_aliases = self.ca.board.icarus_subreg_aliases
                  self.ca.board.monitor_controls = self.ca.board.icarus_monitor_controls
00058
00059
              else:
00060
                  self.ca.board.subreg_aliases = self.ca.board.daedalus_subreg_aliases
00061
                  self.ca.board.monitor_controls = self.ca.board.daedalus_monitor_controls
00062
          # TODO: Check if 'jumpers' still apply for newer boards
def checkSensorVoltStat(self):
00063
00064
00065
00066
              Checks register tied to sensor select jumpers to confirm match with sensor
00067
              object
00068
00069
              Returns:
              boolean, True if jumpers select for Icarus sensor _{\tt m,m}
00070
00071
00072
              logging.debug(self.logdebug + "checkSensorVoltStat")
00073
              err, status = self.ca.getSubregister(self.detect)
00074
              if err:
00075
                  logging.error(self.logerr + "Unable to confirm sensor status")
00076
                  return False
00077
              if not int(status):
00078
                  logging.error(self.logerr + self.sensfam + " sensor not detected")
00079
                  return False
00080
              return True
00081
00082
          def setInterlacing(self, ifactor):
00083
00084
              Virtual function; feature is not implemented on Icarus
00085
              Overridden in daedalus.py
00086
00087
              Returns:
              integer 0
00088
00089
00090
              if ifactor:
00091
                  logging.warning(
                      self.logwarn + "Interlacing is not supported by Icarus sensors. "
00092
00093
00094
              return 0
00095
00096
          def setHighFullWell(self, flag):
00097
00098
              Virtual function; feature is not implemented on Icarus
00099
              Overridden in daedalus.py
00100
00101
              if flag:
00102
                  logging.warning(
00103
                      self.logwarn + "HighFullWell mode is not supported by Icarus sensors. "
00104
00105
00106
          def setZeroDeadTime(self, flag):
00107
00108
              Virtual function; feature is not implemented on Icarus
00109
              Overridden in daedalus.py
00110
00111
              if flag:
00112
                  logging.warning(
00113
                      self.logwarn + "ZeroDeadTime mode is not supported by Icarus sensors. "
00114
00115
00116
          def setTriggerDelay(self, delay):
00117
00118
              Virtual function; feature is not implemented on Icarus
00119
              Overridden in daedalus.pv
00120
00121
              if delay:
00122
                  logging.warning(
00123
                      self.logwarn + "Trigger Delay is not supported by Icarus sensors. "
00124
00125
```

```
00126
          def setPhiDelay(self, delay):
00127
00128
               Virtual function; feature is not implemented on Icarus
               Overridden in daedalus.py
00129
00130
00131
               if delay:
00132
                   logging.warning(
00133
                       self.logwarn + "Phi Delay is not supported by Icarus sensors. "
00134
00135
00136
          def setExtClk(self, delay):
00137
00138
               Virtual function; feature is not implemented on Icarus
00139
               Overridden in daedalus.py
00140
00141
               if delay:
00142
                   logging.warning(
00143
                       self.logwarn + "External Phi Clock is not supported by Icarus sensors. "
00144
00145
00146
          # TODO: update docs to reflect all sensors
00147
          # TODO: remove default timing?
00148
          # TODO: double-check o+c>20 ns in doc block
          def setTiming(self, side="AB", sequence=None, delay=0):
00149
00150
               Sets timing registers based on 'sequence.' Requesting (0,0) timing will clear the
00151
00152
                 timing register.
               *WARNING* if the entire sequence does not fit into the 40-bit register space,
00153
00154
                 then the actual timings generated may differ from those requested. If the
00155
                 timing sequence fits only once into the register space (i.e., for a single
                 frame, open + closed > 20 ns), then the actual timing will be (n, 40-n), irrespective of the setting of second parameter, e.g. (35,1) will actually
00156
00157
00158
                 result in (35,5) timing.
00159
               \star \text{NOTE}\star Icarus sensors generally cannot use 1 ns timing, so all values for these
                 devices (besides the delay) should be at least 2 \ensuremath{\text{ns}}
00160
00161
00162
               Args:
                   side: Hemisphere 'A', 'B', 'AB'
00163
                   sequence: two-element tuple of timing durations in ns, e.g., '(5,2)'
00164
00165
                   delay: initial delay in ns (1 ns delay is acceptable)
00166
00167
               Returns:
00168
                   tuple (error string, 10-character hexadecimal representation of timing
00169
                     sequence)
00170
00171
               logging.info(
00172
                   "{}setTiming; side={}, sequence={}, delay={}".format(
00173
                       self.loginfo, side, sequence, delay
00174
00175
00176
               if sequence is None:
00177
                   sequence = (3, 2)
00178
               if delay is None:
00179
                   delay = 0
00180
               logging.info(
00181
                   self.loginfo
00182
                   + "HST side "
00183
                   + side.upper()
00184
                   + ": "
00185
                   + str(sequence)
00186
                   + "; delay =
00187
                   + str(delay)
00188
               err = ""
00189
00190
               if len(sequence) != 2:
00191
                   err = (
                       self.logerr
00192
00193
                        + "Invalid sequence setting for side: "
00194
                        + side
00195
                        + "; timing settings are unchanged"
00196
00197
                   logging.error(err)
00198
                   return err, "0000000000"
               if side.upper() == "AB":
00199
                   err1, _ = self.setTiming(side="A", sequence=sequence, delay=delay)
00200
00201
                   err2, full40hex = self.setTiming(side="B", sequence=sequence, delay=delay)
00202
                   return err1 + err2, full40hex
00203
               if side.upper() == "A":
00204
                   lowreg = "HS_TIMING_DATA_ALO"
              highreg = "HS_TIMING_DATA_AHI"
elif side.upper() == "B":
00205
00206
```

7.40 sensorBase.py 395

```
00207
                   lowreg = "HS_TIMING_DATA_BLO"
                   highreg = "HS_TIMING_DATA_BHI"
00208
00209
               else:
00210
                   err = (
00211
                       self.logerr
00212
                        + "setTiming: Invalid sensor side: "
00213
                        + side
00214
                        + "; timing settings unchanged"
00215
00216
                   logging.error(err)
00217
                   return err, "0000000000"
00218
               if (sequence[0] + sequence[1]) + delay > 40:
00219
                   err = (
                       self.logerr
00220
00221
                        + "setTiming: Timing sequence is too long to be implemented; "
00222
                        "timing settings unchanged "
00223
00224
                   logging.error(err)
                   return err, "0000000000"
00225
00226
00227
               self.ca.senstiming[side.upper()] = (sequence, delay)
00228
               self.ca.sensmanual = [] # clear manual settings from ca
00229
00230
               full40 = [0] * 40
00231
               bitlist = []
00232
               flag = 1
00233
               sequence = sequence[:2] # TODO: is this redundant?
00234
               for a in sequence:
                   add = [flag] * a
00235
                   bitlist += add
00236
00237
                   if flag:
00238
                       flag = 0
00239
                   else:
00240
                      flag = 1
               # automatically truncates sequence to 39 characters
00241
               logging.debug(self.logdebug + "bitlist = " + str(bitlist))
00242
              if bitlist: # skip this if timing is [0,0]
    reversedlist = bitlist[39::-1]
00243
00244
00245
                   trunclist = reversedlist[:]
00246
                   while trunclist[0] == 0:
00247
                       trunclist.pop(0)
00248
                   # fullrepeat counts open/closed cycles, doesn't include final frame
                   fullrepeats = (40 - len(trunclist) - delay) // len(reversedlist)
logging.debug(self.logdebug + "fullrepeats = " + str(fullrepeats))
00249
00250
00251
                   \# Pattern from sequence repeated to fit inside 40 bits
00252
                   repeated = trunclist + reversedlist * fullrepeats
00253
                   full40[-(len(repeated) + delay + 1) : -(delay + 1)] = repeated
00254
               else:
00255
                   logging.warning(self.logwarn + "setTiming: all-zero timing supplied")
                   fullrepeats = self.nframes
00256
               full40bin = "".join(str(x) for x in full40)
logging.debug(self.logdebug + "full40bin = " + str(full40bin))
00257
00258
00259
               full40hex = "%x" % int(full40bin, 2)
00260
               logging.debug(self.logdebug + "full40hex = " + str(full40hex))
00261
               highpart = full40hex[-10:-8].zfill(8)
00262
               lowpart = full40hex[-8:].zfill(8)
00263
               err0, _ = self.ca.setRegister(lowreg, lowpart)
00264
               errl, _ = self.ca.setRegister(highreg, highpart)
               err2, _ = self.ca.setSubregister("MANSHUT_MODE", "0")
00265
00266
               err3, _ = self.ca.setSubregister("HST_MODE", "1")
               err = err0 + err1 + err2 + err3
00267
00268
               if err:
00269
                   logging.error(
00270
                       self.logerr + "setTiming: Timing may not have been set correctly"
00271
00272
               if fullrepeats < self.nframes - 1:</pre>
00273
                   actual = self.getTiming(side, actual=True)
00274
                   if self.fpganumID == 2:
00275
                       expected = [delay] + 2 * list(sequence) + [sequence[0]]
00276
                   else:
00277
                      expected = [delay] + 3 * list(sequence) + [sequence[0]]
00278
                   if actual != expected:
00279
                       logging.warning(
00280
                           self.logwarn
                            + "setTiming: Due to sequence length"
00281
00282
                            + self.specwarn
00283
                            + ", the actual timing '
00284
                            "sequence for side
00285
                            + side
00286
                            + " will be "
00287
```

```
00288
                            + str(actual[0])
00289
                            + "}"
+ ""
00290
                            + str(actual[1 : 2 * self.nframes])
00291
00292
00293
               elif self.ca.sensorname == "icarus":
00294
                   f0delay = sequence[0] + sequence[1]
00295
                   logging.warning(
00296
                       self.logwarn + "setTiming: Due to use of the Icarus model 1 sensor, the"
00297
                        " initial delay for side "
00298
                        + " will actually be "
00299
00300
                        + str(delay + f0delay)
00301
                        + " nanoseconds"
00302
                   )
00303
               return err, full40hex
00304
00305
           # TODO: restore after power cycle?
00306
           # TODO: smart interpretation of Icarus1 timing?
00307
           # TODO: error checking like in getTiming
00308
           def setArbTiming(self, side="AB", sequence=None):
00309
00310
               Set arbitrary high-speed timing sequence.
00311
               Aras:
00312
                   side: Hemisphere 'A', 'B', 'AB'
                   sequence: list of arbitrary timing intervals, beginning with initial delay. The conventional timing (3,2) with delay = 0 would be represented by
00313
00314
                      [0,3,2,3,2,3] on icarus devices, [0,3,2,3,2,3] on daedalus. If used
00315
                     for interlacing or ZDT, you should populate the entire 40-bit register, e.g., [0,3,2,3,2,3,2,3,2,3,2,3,2,3,2,3,2]
00316
00317
00318
00319
                   *NOTE* Icarus sensors generally cannot use 1 ns timing, so should use at
00320
                     least 2 ns for frames 2 and 3 integration and interframe times (an initial
00321
                     delay of only 1 ns is acceptable)
00322
00323
                   \star NOTE \star although the Icarus model 1 only images the middle two frames, timing
00324
                      entries must be provided for all four frames; to implement frame 1 open
00325
                      for X ns, shutter closed for Y ns, and frame 2 open for Z ns, use the
00326
                     sequence [0,1,1,X,Y,Z,1,1]
00327
00328
                   *WARNING* arbitrary timings will not be restored after a board power cycle
00329
00330
               Returns:
               list: Actual timing results
00331
00332
00333
               logging.info(
00334
                   "{}setArbTiming; side={}, sequence={}".format(self.loginfo, side, sequence)
00335
00336
               if sequence is None:
                   if self.sensfam == "Daedalus":
00337
00338
                       sequence = [0, 2, 3, 4, 5, 6]
00339
                   else:
00340
                       sequence = [0, 2, 3, 4, 5, 6, 7, 8]
00341
00342
                   self.loginfo + "HST side " + side.upper() + " (arbitrary): " + str(sequence)
00343
00344
               if side.upper() == "AB":
                   err1, _ = self.setArbTiming(side="A", sequence=sequence)
err2, actual = self.setArbTiming(side="B", sequence=sequence)
00345
00346
00347
                    return err1 + err2, actual
               if side.upper() == "A":
00348
00349
                   lowreg = "HS_TIMING_DATA_ALO"
                   highreg = "HS_TIMING_DATA_AHI"
00350
               elif side.upper() == "B":
00351
                   lowreg = "HS_TIMING_DATA_BLO"
00352
                   highreg = "HS_TIMING_DATA_BHI"
00353
               else:
00354
00355
                   err = (
00356
                       self.logerr
                        + "Invalid sensor side: "
00357
00358
                        + side
00359
                        + "; timing settings unchanged"
00360
00361
                   logging.error("{}setArbTiming: {}".format(self.logerr, err))
                   return err, "0000000000"
00362
00363
00364
               full40 = [0] * 40
00365
               bitlist = []
               flag = 0 # similar to setTiming, but starts with delay
00366
00367
00368
               for a in sequence:
```

7.40 sensorBase.py 397

```
00369
                  add = [flag] * a
00370
                  bitlist += add
00371
                   if flag:
                      flag = 0
00372
00373
                  else:
00374
                      flag = 1
00375
00376
              logging.debug(self.logdebug + "bitlist = " + str(bitlist))
00377
              reversedlist = bitlist[39::-1]
00378
              full40[-(len(reversedlist) + 1) : -1] = reversedlist
00379
              full40bin = "".join(str(x) for x in full40)
              logging.debug(self.logdebug + "full40bin = " + str(full40bin))
00380
00381
              full40hex = "%x" % int(full40bin, 2)
              logging.debug(self.logdebug + "full40hex = " + str(full40hex))
00383
              highpart = full40hex[-10:-8].zfill(8)
00384
              lowpart = full40hex[-8:].zfill(8)
00385
              self.ca.setRegister(lowreg, lowpart)
00386
              self.ca.setRegister(highreg, highpart)
              # deactivates manual shutter mode if previously engaged
00387
00388
              self.ca.setSubregister("MANSHUT_MODE", "0")
00389
              self.ca.setSubregister("HST_MODE", "1")
00390
              actual = self.getTiming(side, actual=True)
00391
              f0delay = sequence[1] + sequence[2]
00392
00393
              if self.ca.sensorname == "icarus":
00394
                  if actual != sequence[:1] + sequence[3:6]:
00395
                      logging.warning(
00396
                           self.logwarn + "Due to sequence length and use of the Icarus model "
00397
                           "1 sensor, the actual timing sequence for side "
00398
                           + side
00399
                           + " will be "
                           + "{"
00400
00401
                           + str(actual[0] + f0delay)
                          + "}"
00402
00403
                           + str(actual[1 : 2 * self.nframes])
00404
00405
                      )
00406
                  else:
00407
                      logging.warning(
                           self.logwarn + "Due to use of the Icarus model 1 sensor, the actual"
00408
                           " timing sequence for side "
00409
00410
                           + side
                           + " will be " + "{"
00411
00412
00413
                           + str(actual[0] + f0delay)
00414
                          + "}"
00415
00416
                           + str(actual[1 : 2 * self.nframes])
00417
00418
              else:
00419
                  if actual != sequence:
00420
                       logging.warning(
00421
                           self.logwarn + "Due to sequence length, actual timing sequence "
00422
                           "for side "
00423
                           + side
00424
                           + " will be "
                           + "{"
00425
00426
                           + str(actual[0])
00427
                           + "}"
                           + " "
00428
00429
                           + str(actual[1 : 2 * self.nframes])
00430
00431
              return "", actual
00432
00433
          # TODO: figure out how to handle interlacing?
00434
          def getTiming(self, side, actual):
00435
00436
              actual = True: returns actual high speed intervals that will be generated by the
00437
                          FPGA as list
00438
                        False: Returns high speed timing settings as set by setTiming. Assumes
00439
                           that timing was set via the setTiming method--it will not accurately
00440
                           report arbitrary timings set by direct register sets or manual
00441
                           shutter control
00442
00443
              Aras:
00444
                  side: Hemisphere 'A' or 'B'
00445
                  actual: False: return HST settings
00446
                           True: calculate and return actual HST behavior
00447
00448
              Returns:
00449
                  actual= True: list of shutter intervals;
```

```
00450
                               icarus: [delay, open0, closed0, open1, closed1, open2, closed2,
00451
                                         open31
00452
                               daedalus: [delay, open0, closed0, open1, closed1, open2]
00453
                           False: tuple (hemisphere label,
00454
                                            'open shutter' in ns,
00455
                                            'closed shutter' in ns,
00456
                                            initial delay in ns)
00457
00458
              logging.info("{}getTiming".format(self.loginfo))
00459
              if side is None:
00460
00461
                  side = "A"
00462
              logging.info(self.loginfo + "get timing, side " + side.upper())
              if side.upper() == "A":
00463
00464
                   lowreg = "HS_TIMING_DATA_ALO"
00465
                  highreg = "HS_TIMING_DATA_AHI"
00466
              elif side.upper() == "B":
                   lowreg = "HS_TIMING_DATA_BLO"
00467
                  highreg = "HS_TIMING_DATA_BHI"
00468
00469
              else:
00470
                  logging.error(
00471
                       self.logerr
00472
                       + "Invalid sensor side: "
00473
                       + side
00474
                       + "; timing settings unchanged"
00475
00476
                   return "", 0, 0, 0
              err, lowpart = self.ca.getRegister(lowreg)
00477
00478
              err1, highpart = self.ca.getRegister(highreg)
00479
              if err or err1:
00480
                  logging.error(
                      self.logerr + "Unable to retrieve timing setting (getTiming), "
00481
00482
                       "returning zeroes "
00483
              return side.upper(), 0, 0, 0
full40hex = highpart[-2:] + lowpart.zfill(8)
00484
00485
              full40bin = "{0:0=40b}".format(int(full40bex, 16))
logging.debug(self.logdebug + "full40bin = " + str(full40bin))
00486
00487
              if actual:
00488
                   if full40bin == "0" * 40: # all-zero timing
00489
                       if self.fpganumID == 2:
00490
00491
                           times = [0] * 6
00492
                       else:
                          times = [0] * 8
00493
00494
                   else:
00495
                       full160 = 4 * full40bin
00496
                       gblist = [[k, len(list(g))] for k, g in itertools.groupby(full160)]
00497
                       if self.fpganumID == 2:
00498
                           times = [int(x[1]) for x in gblist[:-7:-1]]
00499
                       else:
00500
                           times = [int(x[1]) for x in gblist[:-9:-1]]
00501
                       times[0] = times[0] - 1
                   if self.ca.sensorname == "icarus":
00502
00503
                       \# get timing for frames 1 and 2, keep delay as offset
00504
                       # TODO: should this give a 'correct' offset from frame 0?
00505
                       times12 = [times[0]] + times[3:6]
00506
                       return times12
00507
                   return times
00508
              else:
00509
                  if full40bin == "0" * 40: # all-zero timing
00510
                      timeon, timeoff, delay = (0, 0, 0)
00511
00512
                       gblist = [[k, len(list(g))] for k, g in itertools.groupby(full40bin)]
                       delay = gblist[-1][1] - 1
00513
00514
                       timeon = gblist[-2][1]
00515
00516
                       if self.ca.sensorname == "icarus":
00517
                           if len(gblist) == 2: # 39,1 corner case
00518
                               timeoff = 1
00519
                           elif len(gblist) == 3: # sequence fits only once
00520
                              timeoff = 40 - timeon
00521
                           else:
00522
                               timeoff = qblist[-3][1]
00523
                       else:
00524
                           if len(gblist) < self.nframes: # sequence fits only once</pre>
00525
                               timeoff = 40 - timeon
00526
                           else:
                               # TODO: confirm '-3' works for daedalus
00527
                               timeoff = gblist[-3][1]
00528
00529
                   return side.upper(), timeon, timeoff, delay
00530
```

7.40 sensorBase.py 399

```
00531
             def setManualShutters(self, timing=None):
00532
                  Legacy alias for setManualTiming()
00533
00534
00535
                  self.setManualTiming(timing)
00536
00537
             def setManualTiming(self, timing=None):
00538
00539
                  Manual shutter timing, seven intervals to assign to both hemispheres, e.g.,
                     [(100,150,100,150,100,150,100)] for frame 0 open for 100 ns, an interframe
00540
00541
                     pause of 50 ns, frame 1 open for 100 ns, etc.
                  Provide two sets of seven intervals, e.g., [(100,150,100,150,100,150,100),
00542
00543
                     (200, 250, 200, 250, 200, 250, 200)] to program the A and B hemispheres
00544
                     independently
00545
00546
                  Overridden in daedalus.py
00547
00548
                  The timing list is flattened before processing; the suggested tuple structure is
00549
                     just for clarity (first tuple is A, second is B) and is optional.
00550
00551
                  The actual timing is rounded down to the nearest multiple of 25 ns. (Each
00552
                    count = 25 ns. e.g., a request for 140 ns rounds down to a count of '5',
00553
                     which corresponds to 125 ns))
00554
                       - Minimum timing is 75 ns
00555
                        - Maximum is 25 * 2^30 ns (approximately 27 seconds)
00556
00557
                  Aras:
00558
                       timing: 7- or 14-element list (substructure optional) in nanoseconds
00559
00560
                  tuple (error string, response string from final message)
00561
00562
                  if timing is None:
00563
00564
                       logging.info(
                             self.loginfo
00565
00566
                             + "No manual timing setting provided, defaulting to (100, 150, 100, "
00567
                             " 150, 100, 150, 100) for both hemispheres"
00568
                       timing = [
(100, 150, 100, 150, 100, 150, 100),
00569
00570
00571
                             (100, 150, 100, 150, 100, 150, 100),
00572
                  logging.info(self.loginfo + "Manual shutter sequence: " + str(timing))
00573
00574
                  flattened = flattenlist(timing)
00575
                   if len(flattened) == 7:
00576
                        flattened = 2 * flattened
00577
00578
                       len(flattened) != 14
00579
                       or not all(isinstance(x, numbers.Real) for x in flattened)
00580
                        or not all (x \ge 75 \text{ for } x \text{ in flattened})
00581
                        or not all(x \leq 26843545600 for x in flattened)
00582
00583
                        err = self.logerr + "Invalid manual shutter timing list: " + str(timing)
                       logging.error(err + "; timing settings unchanged")
return err, "00000000"
00584
00585
00586
00587
                  timecounts = [int(a // 25) for a in flattened]
00588
                  self.ca.sensmanual = timing
00589
                  self.ca.senstiming = {} # clear HST settings from ca object
00590
00591
                  control_messages = [
                        ("WO_INTEGRATION", "{0:#0{1}x}".format(timecounts[0], 10)[2:10]), ("WO_INTEGRATION", "{0:#0{1}x}".format(timecounts[1], 10)[2:10]), ("W1_INTEGRATION", "{0:#0{1}x}".format(timecounts[2], 10)[2:10]), ("W1_INTEGRAME", "{0:#0{1}x}".format(timecounts[3], 10)[2:10]),
00592
00593
00594
00595
                        ("W2_INTEGRATION", "{0:#0{1}x}".format(timecounts[4], 10)[2:10]), ("W2_INTEGRATION", "{0:#0{1}x}".format(timecounts[5], 10)[2:10]),
00597
                        ("W3_INTEGRATION", "{0:#0{1}x}".format(timecounts[6], 10)[2:10]),
00598
                        ("W3_INTEGRATION", "{0:#0{1}x}".format (timecounts[6], 10)[2:10]),
("W0_INTEGRATION_B", "{0:#0{1}x}".format (timecounts[7], 10)[2:10]),
("W0_INTERFRAME_B", "{0:#0{1}x}".format (timecounts[8], 10)[2:10]),
("W1_INTEGRATION_B", "{0:#0{1}x}".format (timecounts[9], 10)[2:10]),
("W2_INTEGRATION_B", "{0:#0{1}x}".format (timecounts[10], 10)[2:10]),
("W2_INTEGRATION_B", "{0:#0{1}x}".format (timecounts[11], 10)[2:10]),
("W2_INTEGRATION_B", "{0:#0{1}x}".format (timecounts[12], 10)[2:10]),
("W3_INTEGRATION_B", "{0:#0{1}x}".format (timecounts[13], 10)[2:10]),
00599
00600
00601
00602
00603
00604
00605
                        ("HST_MODE", "0"),
00606
                        ("MANSHUT_MODE", "1"),
00607
00608
00609
                  return self.ca.submitMessages(control_messages, " setManualShutters: ")
00610
             def getManualTiming(self):
00611
```

```
00612
00613
              Read off manual shutter timing settings
00614
              Overridden in daedalus.py
00615
              list of 2 lists of timing from A and B sides, respectively """
00616
00617
00618
              aside = []
00619
              bside = []
00620
              for reg in [
00621
                   "WO_INTEGRATION",
00622
                   "WO_INTERFRAME"
00623
                   "W1_INTEGRATION"
00624
                   "W1_INTERFRAME",
00625
                  "W2_INTEGRATION",
00626
                   "W2_INTERFRAME",
00627
                   "W3_INTEGRATION",
00628
              1:
00629
                   _, reghex = self.ca.getRegister(reg)
00630
                   aside.append(25 * int(reghex, 16))
00631
              for reg in [
00632
                   "WO_INTEGRATION_B",
00633
                   "WO_INTERFRAME_B"
                   "W1_INTEGRATION_B",
00634
00635
                   "W1_INTERFRAME_B",
00636
                   "W2_INTEGRATION_B",
00637
                   "W2 INTERFRAME B",
                   "W3_INTEGRATION_B",
00638
00639
              1:
00640
                    , reghex = self.ca.getRegister(reg)
                  bside.append(25 * int(reghex, 16))
00641
00642
              return [aside, bside]
00643
          def getSensTemp(self, scale=None, offset=None, slope=None, dec=None):
00644
00645
              \label{thm:continuous} \mbox{ Virtual method (Temperature sensor is not present on Icarus sensors). Returns \ \mbox{0.}
00646
00647
              Overridden by Daedalus method
00648
00649
              return 0
00650
          def selectOscillator(self, osc=None):
00651
00652
00653
              Selects oscillator to control sensor timing
00654
              Overridden in daedalus.py
00655
              Args:
00656
                  osc: 'relaxation'|'ring'|'ringnoosc'|'external', defaults to relaxation
00657
00658
              Returns:
              error message as string
00659
00660
00661
              logging.info(self.loginfo + "selectOscillator; osc = " + str(osc))
00662
              if osc is None:
00663
                  osc = "rel"
00664
              osc = str(osc)
00665
              if osc.upper()[:3] == "REL":
                  payload = "00"
00666
00667
              elif osc.upper()[:3] == "RIN":
00668
                if "NO" in osc.upper() or "0" in osc:
00669
                      payload = "10"
00670
00671
                      payload = "01"
00672
              elif osc.lower()[:3] in ["ext", "phi"]:
                  payload = "11"
00673
              else:
00674
00675
                  err = (
00676
                       self.logerr + "selectOscillator: invalid parameter supplied. "
00677
                       "Oscillator selection is unchanged."
00678
00679
                  logging.error(err)
00680
                   return err
              self.ca.setSubregister("OSC_SELECT", payload)
00681
00682
00683
          def parseReadoff(self, frames, columns):
00684
00685
              Virtual method (Order parsing is unnecessary for Icarus, continue to hemisphere
00686
                parsing.)
00687
              Overridden by Daedalus method
00688
00689
              return self.ca.partition(frames, columns)
00690
00691
          def getSensorStatus(self):
00692
```

```
00693
             Wrapper for reportSensorStatus so that the user doesn't have to query statusbits
00694
00695
             sb1 = self.ca.board.checkstatus()
00696
             sb2 = self.ca.board.checkstatus2()
00697
             self.reportStatusSensor(sb1, sb2)
00698
00699
         def reportStatusSensor(self, statusbits, statusbits2):
00700
00701
             Print status messages from sensor-specific bits of status register, default for
00702
               Icarus family sensors
00703
             Args:
00704
                statusbits: result of checkStatus()
             statusbits2: result of checkStatus2()
00705
            if int(statusbits[3]):
00708
                 print(self.loginfo + "W3_Top_A_Edge1 detected")
00709
            if int(statusbits[4]):
00710
                 print(self.loginfo + "W3_Top_B_Edge1 detected")
             if int(statusbits[12]):
                 print(self.loginfo + "HST_All_W_En detected")
00712
             if self.ca.boardname == "llnl_v4" and int(statusbits2[5]):
00713
                 print(self.loginfo + "PDBIAS Unready")
00714
00715
00716
00717 # TODO: add function to control TIME_ROW_DCD delay
00718
00719 """
00720 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00721 LLNL-CODE-838080
00723 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00724 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00725 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00726 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00727 be made under this license.
00728 """
```

# 7.41 C:/Users/hill35/git/camera\_python/nsCamera/utils/crc16pure.py File Reference

# **Namespaces**

- namespace nsCamera
- namespace nsCamera.utils
- namespace nsCamera.utils.crc16pure

### **Functions**

- nsCamera.utils.crc16pure.\_crc16 (data, crc, table)
- nsCamera.utils.crc16pure.crc16xmodem (data, crc=0)

# **Variables**

list nsCamera.utils.crc16pure.CRC16 XMODEM TABLE

# 7.42 crc16pure.py

#### Go to the documentation of this file.

```
00001 #!/usr/bin/env python
00002 """Pure python library for calculating CRC16
00003 NOTE: modified slightly to combine Python 2 and Python 3 versions in single file
00005
00006
00024
00025 import sys
00026
00027 # table for calculating CRC
00028 # this particular table was generated using pycrc v0.7.6, http://www.ttyl.net/pycrc/
00029 # using the configuration:
00030 # *
               Width
00031 # *
               Poly
                             = 0 \times 1021
00032 # *
                             = 0x0000
               XorIn
00033 # *
                            = False
               ReflectIn
00034 # *
                             = 0x0000
               XorOut
               ReflectOut
                            = False
00035 # *
00036 # *
               Algorithm
                             = table-driven
00037 # by following command:
00038 # python pycrc.py --model xmodem --algorithm table-driven --generate c 00039 CRC16_XMODEM_TABLE = [
           0x0000,
00040
00041
           0x1021,
00042
           0x2042,
00043
           0x3063,
00044
           0x4084,
00045
           0x50A5,
00046
           0x60C6,
00047
           0x70E7,
00048
           0x8108,
00049
           0x9129,
00050
           0xA14A,
00051
           0xB16B,
00052
           0xC18C,
00053
           0xD1AD,
00054
           0xE1CE,
00055
           0xF1EF,
00056
           0x1231,
           0x0210,
00057
00058
           0x3273,
00059
           0x2252,
00060
           0x52B5,
00061
           0x4294,
00062
           0x72F7,
00063
           0x62D6,
00064
           0x9339,
00065
           0x8318,
00066
           0xB37B,
00067
           0xA35A,
00068
           0xD3BD,
00069
           0xC39C,
00070
           0xF3FF,
00071
           0xE3DE,
00072
           0x2462,
00073
00074
           0x0420,
00075
           0x1401,
00076
           0x64E6,
00077
           0x74C7,
00078
           0x44A4,
00079
           0x5485,
08000
           0xA56A,
00081
           0xB54B,
00082
           0x8528,
00083
           0x9509,
00084
           0xE5EE,
00085
           0xF5CF,
00086
           0xC5AC,
00087
           0xD58D,
00088
           0x3653,
00089
           0x2672,
00090
           0x1611,
00091
           0x0630,
00092
           0x76D7,
00093
           0x66F6,
```

7.42 crc16pure.py 403

00094	0x5695,
00095	0x46B4,
00096	0xB75B,
00097	0xA77A,
00098	0x9719,
	0.0720
00099	0x8738,
00100	0xF7DF,
00101	0xE7FE,
	OAE/FE,
00102	0xD79D,
00103	0xC7BC,
00104	0x48C4,
00105	0x58E5,
00106	0x6886,
00107	0x78A7,
00108	0x0840,
00109	0x1861,
00110	0x2802,
00111	0x3823,
00112	0xC9CC,
00113	0xD9ED,
00114	0xE98E,
	0xF9AF,
00115	
00116	0x8948,
00117	0x9969,
00118	0xA90A,
00119	0xB92B,
00120	0x5AF5,
00121	0x4AD4,
00122	0x7AB7,
00123	0x6A96,
00124	0x1A71,
00125	0x0A50,
00126	0x3A33,
00127	0x2A12,
00128	0xDBFD,
00129	0xCBDC,
00130	
00131	0xEB9E,
00132	0x9B79,
00133	0x8B58,
00133	
00134	0xBB3B,
00135	0xAB1A,
00136	0x6CA6,
00137	0x7C87,
00138	0x4CE4,
00139	0x5CC5,
00140	0x2C22,
00141	0x3C03,
00142	0x0C60,
00143	0x1C41,
00144	0xEDAE,
00145	0xFD8F,
00146	0xCDEC,
00147	
00148	0xDDCD,
	0xAD2A,
00149	
	0xAD2A, 0xBD0B,
00150	0xAD2A, 0xBD0B, 0x8D68,
00150 00151	0xAD2A, 0xBD0B, 0x8D68, 0x9D49,
00150 00151 00152	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97,
00150 00151 00152 00153	0xAD2A, 0xBD0B, 0x8D68, 0x9D49,
00150 00151 00152 00153	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6,
00150 00151 00152 00153 00154	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5,
00150 00151 00152 00153 00154 00155	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4,
00150 00151 00152 00153 00154 00155 00156	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13,
00150 00151 00152 00153 00154 00155	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4,
00150 00151 00152 00153 00154 00155 00156	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x2E32,
00150 00151 00152 00153 00154 00155 00156 00157	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x2E32, 0x1E51,
00150 00151 00152 00153 00154 00155 00156 00157 00158	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x2E32, 0x1E51, 0x0E70,
00150 00151 00152 00153 00154 00155 00156 00157 00158 00159	0xAD2A, 0xBD0B, 0x8D0B, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x2E32, 0x1E51, 0x0E70, 0xFF9F,
00150 00151 00152 00153 00154 00155 00156 00157 00158 00159 00160	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x2E32, 0x1E51, 0x0E70,
00150 00151 00152 00153 00154 00155 00156 00157 00158 00159	0xAD2A, 0xBD0B, 0x8D0B, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x2E32, 0x1E51, 0x0E70, 0xFF9F,
00150 00151 00152 00153 00154 00155 00156 00157 00158 00159 00160 00161	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x2E32, 0x1E51, 0x0F70, 0xFF9F, 0xEFBE, 0xDFDD,
00150 00151 00151 00153 00154 00155 00156 00157 00158 00159 00160 00161 00162	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF1, 0x0E70, 0xFF9F, 0xEFBE, 0xDFDD, 0xCFFC,
00150 00151 00152 00153 00154 00155 00156 00157 00158 00159 00160 00161 00162 00163	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x2E32, 0x1E51, 0x0E70, 0xFF9F, 0xDFDD, 0xCFFC, 0xBF1B,
00150 00151 00152 00153 00154 00155 00156 00157 00158 00159 00160 00161 00162 00163 00164	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x2E32, 0x1E51, 0x0F70, 0xFF9F, 0xFFBE, 0xDFDD, 0xCFFC, 0xBF1B, 0xAF3A,
00150 00151 00152 00153 00154 00155 00156 00157 00158 00159 00160 00161 00162 00163	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x2E32, 0x1E51, 0x0F70, 0xFF9F, 0xFFBE, 0xDFDD, 0xCFFC, 0xBF1B, 0xAF3A,
00150 00151 00152 00153 00154 00155 00156 00157 00158 00159 00160 00161 00162 00163 00164 00165	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x2E32, 0x1E51, 0x0F70, 0xFFBE, 0xDFDD, 0xCFFC, 0xAF3A, 0x9F59,
00150 00151 00152 00153 00154 00155 00156 00157 00158 00159 00160 00161 00162 00163 00164 00165 00166	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x2E32, 0x1E51, 0x0F70, 0xFF9F, 0xFFBE, 0xDFDD, 0xCFFC, 0xBF1B, 0x9F59, 0x8F78,
00150 00151 00152 00153 00154 00155 00156 00157 00158 00159 00160 00161 00162 00163 00164 00165 00166 00167	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF1, 0x0E70, 0xFF9F, 0xEFBE, 0xDFDD, 0xCFFC, 0xBF1B, 0xAF3A, 0x9F39, 0x8F78, 0x9F88,
00150 00151 00152 00153 00154 00155 00156 00157 00158 00160 00161 00162 00163 00164 00165 00166 00167 00168	0xAD2A, 0xBD0B, 0xBD68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x0E70, 0xFF9F, 0xDFDD, 0xCFFC, 0xBF1B, 0xAF3A, 0x9F59, 0x8F78, 0x9F78, 0x9F78,
00150 00151 00152 00153 00154 00155 00156 00157 00158 00159 00160 00161 00162 00163 00164 00165 00166 00167 00168 00169	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF1, 0x0E70, 0xFF9F, 0xEFBE, 0xDFDD, 0xCFFC, 0xBF1B, 0xAF3A, 0x9F39, 0x8F78, 0x9F88,
00150 00151 00152 00153 00154 00155 00156 00157 00158 00159 00160 00161 00162 00163 00164 00165 00166 00167 00168 00169	0xAD2A, 0xBD0B, 0xBD68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x0E70, 0xFF9F, 0xDFDD, 0xCFFC, 0xBF1B, 0xAF3A, 0x9F59, 0x8F78, 0x9F78, 0x9F78,
00150 00151 00152 00153 00154 00155 00156 00157 00158 00159 00160 00161 00162 00163 00164 00165 00166 00167 00168 00169 00170	0xAD2A, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x2E32, 0x1E51, 0x0F70, 0xFFG, 0xFFG, 0xFFG, 0xBFBB, 0xOF70, 0xBFBB, 0xAF3A, 0x9F59, 0x8F78, 0x8F78, 0x8E78
00150 00151 00151 00152 00153 00154 00155 00156 00157 00158 00159 00160 00161 00162 00163 00164 00165 00166 00166 00166 00167 00168 00169 00170 00171	0xAD2A, 0xBD0B, 0x8D0B, 0x8D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x2E32, 0x1E51, 0x0F70, 0xFF9F, 0xCFFC, 0xBF1B, 0xAF3A, 0x9F59, 0x8F78, 0x9F8, 0x8F18, 0x9T8, 0x8F18, 0x9T80, 0x8F18, 0x9T80, 0x8F18, 0x9T80, 0x8F18, 0x9T80, 0x8F18, 0x8T80, 0x8T8
00150 00151 00152 00153 00154 00155 00156 00157 00158 00160 00161 00162 00163 00164 00165 00166 00167 00168 00169 00170 00171 00172	0xAD2A, 0xBD0B, 0xBD0B, 0x8D68, 0x9D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x0E70, 0xFF9F, 0xDF7D, 0xFF6C, 0xBF1B, 0xAF3A, 0x9F59, 0x8F78, 0x9B1CA, 0xA1EB, 0xB1CA, 0xA1EB, 0xB1CA, 0xA1EB, 0xB1CA, 0xA1EB,
00150 00151 00151 00152 00153 00154 00155 00156 00157 00158 00159 00160 00161 00162 00163 00164 00165 00166 00166 00166 00167 00168 00169 00170 00171	0xAD2A, 0xBD0B, 0x8D0B, 0x8D49, 0x7E97, 0x6EB6, 0x5ED5, 0x4EF4, 0x3E13, 0x2E32, 0x1E51, 0x0F70, 0xFF9F, 0xCFFC, 0xBF1B, 0xAF3A, 0x9F59, 0x8F78, 0x9F8, 0x8F18, 0x9T8, 0x8F18, 0x9T80, 0x8F18, 0x9T80, 0x8F18, 0x9T80, 0x8F18, 0x9T80, 0x8F18, 0x8T80, 0x8T8

00175	0xE16F,
00176	0x1080,
00177	0x00A1,
00178 00179	0x30C2, 0x20E3,
00179	0x5004,
00181	0x4025,
00182	0x7046,
00183	0x6067,
00184 00185	0x83B9, 0x9398,
00186	0xA3FB,
00187	0xB3DA,
00188 00189	0xC33D, 0xD31C,
00109	0xE37F,
00191	0xF35E,
00192	0x02B1,
00193	0x1290,
00194	0x22F3, 0x32D2,
00196	0x4235,
00197	0x5214,
00198	0x6277,
00199	0x7256, 0xB5EA,
00201	0xA5CB,
00202	0x95A8,
00203	0x8589,
00204 00205	0xF56E, 0xE54F,
00206	0xD52C,
00207	0xC50D,
00208	0x34E2,
00209 00210	0x24C3, 0x14A0,
00210	0x0481,
00212	0x7466,
00213	0x6447,
00214 00215	0x5424, 0x4405,
00215	0x4403,
00217	0xB7FA,
00218	0x8799,
00219 00220	0x97B8, 0xE75F,
00220	0xE73E,
00222	0xC71D,
00223	0xD73C,
00224 00225	0x26D3, 0x36F2,
00226	0x0691,
00227	0x16B0,
00228	0x6657,
00229 00230	0x7676, 0x4615,
00231	0x5634,
00232	0xD94C,
00233	0xC96D,
00234 00235	0xF90E, 0xE92F,
00236	0x99C8,
00237	0x89E9,
00238	0xB98A,
00239 00240	0xA9AB, 0x5844,
00241	0x4865,
00242	0x7806,
00243	0x6827,
00244 00245	0x18C0, 0x08E1,
00246	0x3882,
00247	0x28A3,
00248	0xCB7D,
00249 00250	0xDB5C, 0xEB3F,
00250	0xFB1E,
00252	0x8BF9,
00253	0x9BD8,
00254 00255	0xABBB, 0xBB9A,

```
00256
           0x4A75,
00257
           0x5A54,
00258
           0x6A37,
00259
           0x7A16,
00260
00261
00262
           0x2AB3,
00263
           0x3A92,
00264
           0xFD2E,
00265
           0xEDOF,
00267
           0xCD4D,
00268
           0xBDAA,
00269
           0xAD8B,
00270
           0x9DE8,
00271
           0x8DC9,
00272
           0x7C26,
00273
           0x6C07,
00274
           0x5C64,
00275
           0x4C45,
00276
           0x3CA2,
00277
           0x2C83,
00278
           0x1CE0,
00279
           0x0CC1,
00280
           0xEF1F.
00281
           0xFF3E,
00282
           0xCF5D,
00283
           0xDF7C,
00284
           0xAF9B.
00285
           0xBFBA,
00286
           0x8FD9.
00287
           0x9FF8,
00288
           0x6E17,
00289
           0x7E36,
00290
           0x4E55,
00291
           0x5E74,
00292
           0x2E93,
00293
           0x3EB2,
00294
           0x0ED1,
00295
           0x1EF0,
00296 ]
00297
00298
00299 def _crc16(data, crc, table):
00300 """Calculate CRC16 using the given table.
           `data`

    data for calculating CRC, must be a string
    initial value
    table for caclulating CRC (list of 256 integers)

00301
00302
           `crc`
          `table`
00303
00304
          Return calculated value of CRC
00305
00306
         for byte in data:
           if sys.version_info > (3,):
00307
00308
                    crc = ((crc « 8) & 0xFF00) ^ table[((crc » 8) & 0xFF) ^ byte]
00309
00310
                   crc = ((crc « 8) & 0xFF00) ^ table[((crc » 8) & 0xFF) ^ ord(byte)]
00311
00312
         return crc & 0xFFFF
00313
00314
00315 def crc16xmodem(data, crc=0):
00316
           """Calculate CRC-CCITT (XModem) variant of CRC16.
                    - data for calculating CRC, must be a string
- initial value
00318
          `crc`
00319
          Return calculated value of CRC
           return _crc16(data, crc, CRC16_XMODEM_TABLE)
```

# 7.43 C:/Users/hill35/git/camera\_python/nsCamera/utils/FlatField.py File Reference

# Namespaces

namespace nsCamera

- namespace nsCamera.utils
- · namespace nsCamera.utils.FlatField

#### **Functions**

- nsCamera.utils.FlatField.getFilenames (frame="Frame 1")
- nsCamera.utils.FlatField.getROlvector (imgfilename, roi)
- nsCamera.utils.FlatField.tslopes (x, y)
- nsCamera.utils.FlatField.generateFF (FRAMES=["Frame\_0", "Frame\_1", "Frame\_2", "Frame\_3"], roi=[0, 0, 512, 1024], directory="", ncores=-1)
- nsCamera.utils.FlatField.removeFF (filename, directory="", roi=[0, 0, 512, 1024])
- nsCamera.utils.FlatField.removeFFall (directory="", FRAMES=["Frame\_0", "Frame\_1", "Frame\_2", "Frame\_3"], roi=[0, 0, 512, 1024])

# **Variables**

- nsCamera.utils.FlatField.parser = argparse.ArgumentParser()
- nsCamera.utils.FlatField.action
- · nsCamera.utils.FlatField.dest
- · nsCamera.utils.FlatField.default
- · nsCamera.utils.FlatField.help
- · nsCamera.utils.FlatField.nargs
- nsCamera.utils.FlatField.args = parser.parse\_args()
- list nsCamera.utils.FlatField.framelist = ["Frame " + str(frame) for frame in args.frames]
- · nsCamera.utils.FlatField.directory

# 7.44 FlatField.py

# Go to the documentation of this file.

```
00001 # -*- coding: utf-8 -*
00002 """
00003 Functions for batch flat-field image corrections
00005 ***Do not use this file as a template for new code development***
00006
00007 Author: Jeremy Martin Hill (jerhill@llnl.gov)
00008 Author: Matthew Dayton (dayton5@llnl.gov)
00010 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00011 LLNL-CODE-838080
00013 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00014 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00015 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00016 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00017 be made under this license.
00018
00019 Version: 2.1.2 (February 2025)
00020 """
00021
00022 import os
00023 import re
00024
00025 import numpy as np
00026 from PIL import Image
00027 from joblib import parallel, delayed
00028 from scipy.stats import theilslopes
00029 from skimage.external.tifffile import imread
```

7.44 FlatField.pv 407

```
00030
00031
00032 def getFilenames(frame="Frame 1"):
00033
00034
          get a list of tiff filenames in current working director for frame
00035
00036
          onlyfiles = next(os.walk("./"))[2]
00037
          return [k for k in onlyfiles if frame in k and "tif" in k]
00038
00039
00040 def getROIvector(imgfilename, roi):
00041
00042
          return a numpy row vector of version of the image
00043
00044
          img = imread(imgfilename)
00045
          vroi = img[(roi[1]) : (roi[3]), (roi[0]) : (roi[2])].flattenlist()
00046
          return vroi
00047
00048
00049 def tslopes(x, y):
00050
          theilslopes implements a method for robust linear regression.
00051
00052
          It computes the slope as the median of all slopes between paired values.
00053
00054
          val = theilslopes(x, y)
00055
          return [val[0], val[1]]
00056
00057
00058 def generateFF(
          FRAMES=["Frame_0", "Frame_1", "Frame_2", "Frame_3"],
00059
          roi=[0, 0, 512, 1024], directory="",
00060
00061
00062
          ncores=-1.
00063):
00064
          # TODO: documentation
          \ensuremath{\text{\#}} use of ROI here not compatible with use of ROI in \ensuremath{\text{removeFF}}
00065
00066
00067
          if directory:
00068
              cwd = os.getcwd()
00069
              newpath = os.path.join(cwd, directory)
00070
              os.chdir(newpath)
00071
          if not FRAMES:
              print("No framelist provided, defaulting to four frames")
FRAMES = ["Frame_0", "Frame_1", "Frame_2", "Frame_3"]
00072
00073
00074
          for f in FRAMES:
00075
               files = getFilenames(frame=f)
00076
              imgslist = [getROIvector(fn, roi) for fn in files] # a list of flattened images
00077
              imgsarray = np.vstack(imgslist) # turn the list into an array
              npix = np.shape(imgsarray)[1] # total number of pixels
00078
              x = np.median(imgsarray, axis=1) # median of each image used for flat fielding
00079
00080
              y = []
00081
               for i in range(npix):
00082
                  # each member of y represents a pixel, as a list of magnitudes over all the
00083
00084
                   y.append(imgsarray[:, i])
00085
               # get pixel gain and offset for flatfield ff using Thiel-Sen slopes
00086
               ff = []
00087
               ff = parallel.Parallel(n_jobs=ncores, verbose=5, pre_dispatch="2 * n_jobs")(
00088
                   delayed(tslopes)(x, pixel) for pixel in y
00089
00090
               \# x is the dependent variable; here uses median of image as characteristic of
00091
                  noise level
00092
              m, c = zip(*ff)
                                # separate into gain and offset
00093
              m = np.array(m)
              m[m < 0.1] = 0.1 # handle outliers
00094
              m[m > 1000] = 1000 # handle outliers
00095
00096
00097
              m = m.reshape(roi[3] - roi[1], roi[2] - roi[0]) # turn into matrix
              c = np.array(c).reshape(roi[3] - roi[1], roi[2] - roi[0]) # turn into matrix
00098
00099
00100
              with open("px_gain_%s.txt" % f.replace("Frame_", "f"), "w+") as file:
00101
                  np.savetxt(file, m)
               with open("px_off_%s.txt" % f.replace("Frame_", "f"), "w+") as file:
00102
00103
                  np.savetxt(file, c)
00104
00105
00106 def removeFF(filename, directory="", roi=[0, 0, 512, 1024]):
00107
          if directory:
00108
              cwd = os.getcwd()
00109
              newpath = os.path.join(cwd, directory)
00110
              os.chdir(newpath)
```

```
00111
          framenum = re.search("Frame_(\d)", filename).group(1)
00112
          gainFilename = "px_gain_f" + framenum + ".txt"
00113
          gainall = np.loadtxt(gainFilename)
          gain = gainall[(roi[1]) : (roi[3]), (roi[0]) : (roi[2])]
00114
00115
          offFilename = "px_off_f" + framenum + ".txt
00116
          offsetall = np.loadtxt(offFilename, dtype="uint32")
00117
          offset = offsetall[(roi[1]) : (roi[3]), (roi[0]) : (roi[2])]
00118
00119
          beforeImageall = imread(filename)
          beforeImage = beforeImageall[(roi[1]) : (roi[3]), (roi[0]) : (roi[2])]
00120
00121
          imageMed = np.median(beforeImage)
00122
00123
          flat = imageMed * gain + offset
          flat = flat.clip(0)
00124
00125
          fix = beforeImage - flat
00126
          clipped = fix.clip(0)
00127
          fixinit = clipped.astype("uint16")
00128
          fiximg = Image.fromarray(fixinit)
00129
00130
          fixFilename = filename[:-4] + "ff" + filename[-4:]
          fiximg.save(fixFilename)
00131
00132
00133 def removeFFall(
00134
          directory="",
00135
          FRAMES=["Frame_0", "Frame_1", "Frame_2", "Frame_3"],
          roi=[0, 0, 512, 1024],
00136
00137):
00138
          cwd = os.getcwd()
00139
          if directory:
00140
             newpath = os.path.join(cwd, directory)
00141
          else:
            newpath = cwd
00142
00143
          os.chdir(newpath)
          files = next(os.walk("./"))[2]
00144
          filelist = []

for frame in FRAMES:
00145
00146
00147
              filelist.extend([k for k in files if frame in k and "tif" in k])
00148
          for fname in filelist:
             removeFF(fname, directory, roi)
00149
00150
00151
00152 if __name__ == "__main__":
00153
          import argparse
00154
00155
          parser = argparse.ArgumentParser()
00156
         parser.add_argument(
00157
              "-d", action="store", dest="directory", default="", help="VRST scan directory"
00158
00159
         parser.add_argument(
00160
              nargs="+",
00161
00162
              action="store"
00163
              dest="frames",
00164
              default="",
00165
              help="Frame numbers to process, eg. -f 2 3",
00166
00167
          args = parser.parse_args()
00168
          framelist = ["Frame_" + str(frame) for frame in args.frames]
00169
          generateFF(framelist, directory=args.directory)
00170
00171 """
00172 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00173 LLNL-CODE-838080
00175 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00176 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00177 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00178 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00179 be made under this license.
00180 """
```

# 7.45 C:/Users/hill35/git/camera\_python/nsCamera/utils/misc.py File Reference

#### Classes

· class nsCamera.utils.misc.fakeCA

#### **Namespaces**

- namespace nsCamera
- namespace nsCamera.utils
- namespace nsCamera.utils.misc

#### **Functions**

- nsCamera.utils.misc.makeLogLabels (logtag, label)
- nsCamera.utils.misc.getEnter (text)
- nsCamera.utils.misc.checkCRC (rval)
- nsCamera.utils.misc.str2bytes (astring)
- nsCamera.utils.misc.bytes2str (bytesequence)
- nsCamera.utils.misc.str2nparray (valstring)
- nsCamera.utils.misc.flattenlist (x)
- nsCamera.utils.misc.generateFrames (camassem, data, columns=1)
- nsCamera.utils.misc.loadDumpedData (filename="frames.txt", path=None, filetype="txt", sensor="daedalus", firstframe=None, lastframe=None, width=None, height=None, padToFull=None, firstrow=None, lastrow=None, maxwidth=None, maxheight=None, bytesperpixel=None, interlacing=None, columns=1)
- nsCamera.utils.misc.saveTiffs (self, frames, path=None, filename="Frame", prefix=None, index=None)
- nsCamera.utils.misc.plotFrames (self, frames, index=None)
- · nsCamera.utils.misc.partition (self, frames, columns)

## 7.46 misc.py

#### Go to the documentation of this file.

```
00001 # -*- coding: utf-8 00002 """
00003 Miscellaneous utilities, including batch processing of images acquired using the
00004 nsCamera. These are functions that don't require a cameraAssembler object to be
00005
       instantiated before use.
00007 Author: Jeremy Martin Hill (jerhill@llnl.gov)
00009 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00010 LLNL-CODE-838080
00012 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00013 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00014 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00015 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00016 be made under this license.
00017
00018 Version: 2.1.2 (February 2025)
00019 """
00020
00021 import binascii
```

```
00022 import collections
00023 import logging
00024 import os
00025 import sys
00026 from datetime import datetime
00027
00028 import numpy as np
00029 from matplotlib import pyplot as plt
00030 from past.builtins import raw_input
00031 from PIL import Image
00032
00033 # TODO: is the self-reference necessary?
00034 # import nsCamera.utils.misc
00035 from nsCamera.utils import crc16pure
00037
00038 # TODO: check error messages
00039 def makeLogLabels(logtag, label):
          if logtag is None:
00041
               logtag = ""
00042
00043
          logcritbase = "CRITICAL{logtag}: ".format(logtag=logtag)
          logerrbase = "ERROR{logtag}: ".format(logtag=logtag)
00044
00045
          logwarnbase = "WARNING{logtag}: ".format(logtag=logtag)
00046
          loginfobase = "INFO{logtag}: ".format(logtag=logtag)
          logdebugbase = "DEBUG{logtag}: ".format(logtag=logtag)
00047
00048
00049
          logcrit = "{base}{label}".format(base=logcritbase, label=label)
          logerr = "{base}{label}".format(base=logerrbase, label=label)
logwarn = "{base}{label}".format(base=logwarnbase, label=label)
00050
00051
          loginfo = "{base}{label}".format(base=loginfobase, label=label)
logdebug = "{base}{label}".format(base=logdebugbase, label=label)
00052
00053
00054
00055
          return logcrit, logerr, logwarn, loginfo, logdebug
00056
00057
00058 def getEnter(text):
00059
00060
          Wait for enter key to be pressed.
00061
00062
          text: message asking for keypress
00063
00064
          python, _, _, _, _ = sys.version_info
if python >= 3:
00065
00066
00067
               input(text)
00068
          else:
00069
               raw_input(text)
00070
00071
00072 def checkCRC(rval):
00073
00074
          Calculate CRC for rval[:-4] and compare with expected CRC in rval[-4:]
00075
00076
00077
              rval: hexadecimal string
00078
00079
          boolean, True if CRCs match, False if they don't match or the input is invalid
08000
00081
00082
          if not isinstance(rval, str) or len(rval) < 5:</pre>
00083
               logging.error("ERROR: checkCRC: Invalid input: {rval}".format(rval=rval))
00084
               return False
00085
          data_crc = int(rval[-4:], base=16)
00086
          CRC_calc = crc16pure.crc16xmodem(str2bytes(rval[:-4]))
00087
          return CRC_calc == data_crc
00088
00089
00090 def str2bytes(astring):
00091
00092
          Python-version-agnostic converter of hexadecimal strings to bytes
00093
00094
00095
              astring: hexadecimal string without '0x'
00096
00097
          Returns:
          byte string equivalent to input string
00098
00099
00100
00101
          python, _, _, _, = sys.version_info
if python >= 3:
00102
```

7.46 misc.pv 411

```
00103
              try:
00104
                  dbytes = binascii.a2b_hex(astring)
00105
                  logging.error(
00106
00107
                       "ERROR: str2bytes: invalid input: '{astring}'; returning zero"
00108
                      " byte".format(astring=astring)
00109
00110
                  dbytes = b" \x00"
00111
          else:
00112
00113
                  dbytes = astring.decode("hex")
00114
              except:
00115
                  logging.error(
                       "ERROR: str2bytes: invalid input: '{astring}'; returning zero "
00116
00117
                      " byte".format (astring=astring)
00118
00119
                  dbytes = b" \x00"
00120
          return dbytes
00121
00122
00123 def bytes2str(bytesequence):
00124
00125
          Python-version-agnostic converter of bytes to hexadecimal strings
00126
00127
00128
              bytesequence: sequence of bytes as string (Py2) or bytes (Py3)
00129
00130
          Returns:
          hexadecimal string representation of 'bytes' without '0x' """
00131
00132
00133
          try:
00134
              estring = binascii.b2a_hex(bytesequence)
00135
          except TypeError:
00136
              logging.error(
                  "ERROR: bytes2str: Invalid byte sequence: '{bytesequence}'; returning an"
00137
                  " empty string".format(bytesequence=bytesequence)
00138
00139
              )
              return ""
00140
          python, _, _, _, _ = sys.version_info
if python >= 3:
00141
00142
00143
              estring = str(estring)[2:-1]
00144
          return estring
00145
00146
00147 def str2nparray(valstring):
00148
00149
          Convert string into array of uint16s
00150
00151
00152
              valstring: string of hexadecimal characters
00153
00154
          Returns:
          numpy array of uint16
00155
00156
00157
          if not isinstance(valstring, str):
00158
              logging.error(
00159
                  "ERROR: str2nparray: Invalid input: {valstring} is not a string. Returning"
00160
                  " an empty array".format(valstring=valstring)
00161
00162
              return np.array([])
          stringlen = len(valstring)
00163
00164
          arraylen = int(stringlen / 4)
00165
          outarray = np.empty(int(arraylen), dtype="uint16")
00166
00167
          for i in range(0, arraylen):
00168
              try:
                 outarray[i] = int(valstring[4 * i : 4 * i + 4], 16)
00169
00170
              except ValueError:
00171
                  logging.error(
00172
                       "ERROR: str2nparray: input string does not represent a hexadecimal"
00173
                      " integer. Returning an empty array"
00174
00175
                  return np.array([])
00176
          return outarray
00177
00178
00179 def flattenlist(x):
00180
00181
          Flatten list of lists recursively into single list
00182
00183
          python, _, _, _ = sys.version_info
```

```
00184
          try:
00185
              if python >= 3:
00186
                   if isinstance(x, collections.abc.Iterable):
00187
                       return [a for i in x for a in flattenlist(i)]
00188
                   else:
00189
00190
00191
                  if isinstance(x, collections.Iterable):
00192
                       return [a for i in x for a in flattenlist(i)]
00193
                   else:
00194
00195
          except RecursionError:
00196
              logging.error(
                   "ERROR: flattenlist: input \{x\}' is pathological and cannot be flattened."
00198
                   " Attempting to return the input unchanged"
00199
00200
              return x
00201
00202
00203 class fakeCA:
00204
00205
          Fake 'cameraAssembler' object to use as a parameter object in offline functions.
00206
            Returned by generateFrames(), it contains the frame details required to save and
00207
            plot images.
00208
00209
00210
          def __init__(
00211
              self,
              sensorname="icarus2",
00212
00213
              firstframe=0.
00214
              lastframe=3.
00215
              width=512.
              height=1024,
00216
              padToFull=True,
00217
00218
               firstrow=0.
              lastrow=1023,
00219
00220
              maxwidth=512,
00221
              maxheight=1024,
00222
              bytesperpixel=2,
00223
              interlacing=None,
00224
              columns=1.
00225
              logtag=None,
00226
        ):
00227
              self.sensorname = sensorname
              self.boardname = None
self.padToFull = padToFull
00228
00229
              if logtag is None:
    self.logtag = ""
00230
00231
00232
00233
                   self.logtag = logtag
00234
              self.logcritbase = "CRITICAL" + self.logtag + ": "
00235
              self.logerrbase = "ERROR" + self.logtag + ": "
self.logwarnbase = "WARNING" + self.logtag + ": "
00236
00237
00238
               self.loginfobase = "INFO" + self.logtag + ": "
00239
               self.logdebugbase = "DEBUG" + self.logtag + ": "
00240
00241
              self.logcrit = self.logcritbase + "[FS] "
00242
              self.logerr = self.logerrbase + "[FS] "
00243
              self.logwarn = self.logwarnbase + "[FS] "
00244
              self.loginfo = self.loginfobase + "[FS] "
00245
              self.logdebug = self.logdebugbase + "[FS] "
00246
00247
              if self.sensorname == "icarus":
00248
                  import nsCamera.sensors.icarus as snsr
              elif self.sensorname == "icarus2":
00249
00250
                  import nsCamera.sensors.icarus2 as snsr
              elif self.sensorname == "daedalus":
00251
00252
                  import nsCamera.sensors.daedalus as snsr
00253
00254
              self.sensor = snsr(self)
00255
          def partition(self, parsed, columns):
00256
00257
               # return nsCamera.utils.misc.partition(self, parsed, columns)
00258
              return partition (self, parsed, columns)
00259
00260
00261 def generateFrames(camassem, data, columns=1):
00262
00263
          Processes data stream from board into frames and applies sensor-specific parsing.
00264
            Generates padded data for full-size option of setRows.
```

7.46 misc.py 413

```
If used for offline processing, replace the 'self' object with the parameter object
00265
00266
             returned by loadDumpedData().
           If the data stream is incomplete (e.g., from an interrupted download), the data is
00267
00268
             padded with zeros to the correct length.
00269
00270
00271
               self: during normal operation, is the CameraAssembler object. During offline
00272
                 operation, is a parameters object as returned by loadDumpedData()
00273
               data: text stream from board or loaded file, or numpy from loaded file
00274
               columns: 1 for single image, 2 for separate hemisphere images
00275
00276
           Returns: list of parsed frames
00277
00278
           logging.debug("DEBUG" + camassem.logtag + ": generateFrames")
00279
           if isinstance(data[0], str):
00280
               allframes = str2nparray(data)
00281
           else:
00282
               allframes = data
00283
           nframes = camassem.sensor.lastframe - camassem.sensor.firstframe + 1
00284
           frames = [0] * nframes
00285
           framesize = camassem.sensor.width * (
00286
               camassem.sensor.lastrow - camassem.sensor.firstrow + 1
00287
00288
           if hasattr(camassem, "ca"):
00289
               padIt = camassem.ca.padToFull
00290
           else:
               padIt = camassem.padToFull
00291
00292
           if padIt:
00293
               toprows = camassem.sensor.firstrow
00294
               botrows = (camassem.sensor.maxheight - 1) - camassem.sensor.lastrow
00295
               padtop = np.zeros(toprows * camassem.sensor.maxwidth, dtype=int)
00296
               padbot = np.zeros(botrows * camassem.sensor.maxwidth, dtype=int)
00297
               for n in range(nframes):
00298
                    thisframe = np.concatenate(
                        (padtop, allframes[n * framesize : ((n + 1) * framesize)], padbot)
00299
00300
00301
                    frames[n] = thisframe
00302
           else:
               for n in range(nframes):
00304
                    frames[n] = allframes[n * framesize : (n + 1) * framesize]
00305
           # self.clearStatus()
00306
           parsed = camassem.sensor.parseReadoff(frames, columns)
00307
           return parsed
00308
00309
00310 def loadDumpedData(
00311
           filename="frames.txt",
00312
           path=None,
           filetype="txt",
00313
00314
           sensor="daedalus",
00315
           firstframe=None,
00316
           lastframe=None,
00317
           width=None,
00318
           height=None,
00319
           padToFull=None,
00320
           firstrow=None,
00321
           lastrow=None,
00322
           maxwidth=None,
00323
           maxheight=None,
00324
           bytesperpixel=None,
00325
           interlacing=None,
00326
           columns=1,
00327):
00328
           """_summary_
00329
00330
               Output can be piped directly to saveTiffs:
                    pars,frames=loadDumpedData(filename="Dump.npy")
00331
00332
                    saveTiffs(pars, frames)
00333
           Args:
00334
               filename (str, optional): _description_. Defaults to "frames.txt".
00335
               path (_type_, optional): _description_. Defaults to None.
               filetype (str, optional): _description_. Defaults to "txt". sensor (str, optional): _description_. Defaults to "daedalus".
00336
00337
               firstframe (_type_, optional): _description_. Defaults to None. lastframe (_type_, optional): _description_. Defaults to None.
00338
00339
               width (_type_, optional): _description_. Defaults to None. height (_type_, optional): _description_. Defaults to None.
00340
00341
               padToFull (_type_, optional): _description_. Defaults to None.
00342
00343
               firstrow (_type_, optional): _description_. Defaults to None.
               lastrow (_type_, optional): _description_. Defaults to None. maxwidth (_type_, optional): _description_. Defaults to None.
00344
00345
```

```
00346
              maxheight (_type_, optional): _description_. Defaults to None.
00347
              interlacing (_type_, optional): _description_. Defaults to None.
00348
              columns (int, optional): _description_. Defaults to 1.
00349
00350
00351
          Tuple (parameter object, list of data frames)
00352
00353
00354
          logging.debug("DEBUG: loadDumpedData")
          if sensor.lower() == "daedalus":
00355
00356
              import nsCamera.sensors.daedalus as snsr
00357
          elif sensor.lower() == "icarus":
00358
              import nsCamera.sensors.icarus as snsr
          elif sensor.lower() == "icarus2":
00360
              import nsCamera.sensors.icarus2 as snsr
00361
00362
              logging.error(
00363
                   "ERROR loadDumpedData: invalid sensor type provided, defaulting to icarus2"
00364
00365
              import nsCamera.sensors.icarus2 as snsr
00366
00367
          def buildEmptvFrames():
00368
              cols = [0] * 512
00369
              frame = np.array([cols] * (lastrow - firstrow + 1))
00370
              return [frame] * (lastframe - firstframe + 1)
00371
00372
          # get defaults from class declarations if not specified as parameter
00373
          if firstframe is None:
00374
              firstframe = snsr.firstframe
00375
          if lastframe is None:
00376
              lastframe = snsr.lastframe
00377
          # TODO: first frame number
          if width is None:
    width = snsr.width
00378
00379
          if height is None:
00380
              height = snsr.height
00381
          if firstrow is None:
00382
              firstrow = snsr.firstrow
00383
          if lastrow is None:
00384
00385
              lastrow = snsr.lastrow
00386
          if maxwidth is None:
              maxwidth = snsr.maxwidth
00387
          if maxheight is None:
    maxheight = snsr.maxheight
00388
00389
00390
          if bytesperpixel is None:
00391
              bytesperpixel = snsr.bytesperpixel
00392
          if interlacing is None:
00393
              interlacing = snsr.interlacing
00394
00395
          parameters = fakeCA(
00396
              sensor,
00397
              firstframe,
00398
              lastframe,
00399
              width,
00400
              height,
00401
              padToFull,
00402
              firstrow,
00403
              lastrow,
00404
              maxwidth,
00405
              maxheight,
00406
              bytesperpixel,
00407
              interlacing,
00408
00409
00410
          if path is None:
              path = os.path.join(os.getcwd())
00411
00412
          datafile = os.path.join(path, filename)
00413
          if filename[-3:].lower() == "txt":
              filetype = "txt"
00414
          elif filename[-3:].lower() == "npy":
00415
00416
              filetype = "npy"
00417
          # TODO: return empty frames if error
00418
00419
          if filetype == "npy":
00420
              expectedlength = (lastframe - firstframe + 1) * (lastrow - firstrow + 1) * width
00421
                  f = np.load(datafile)
00422
00423
                  padding = expectedlength - len(f)
                  if padding:
00424
00425
                       logging.warning(
   "{logwarn}loadDumpedData: Payload is shorter than expected."
00426
```

7.46 misc.py 415

```
00427
                            " Padding with '0's".format(logwarn=parameters.logwarn)
00428
00429
                       f = np.pad(f, (0, padding), "constant", constant_values=(0))
00430
                   frames = generateFrames(parameters, f, columns)
00431
                   return parameters, frames
00432
               except OSError as err:
00433
00434
                   logging.error(
00435
                       "{logerr}loadDumpedData: OS error: {err}. Returning empty"
                       " frames.".format(logerr=parameters.logerr, err=err)
00436
00437
00438
                   return parameters, buildEmptyFrames()
00439
               except:
00440
                   logging.error(
00441
                        "{logerr}loadDumpedData: Unexpected error: {err}. Returning empty"
00442
                       " frames.".format(logerr=parameters.logerr, err=str(sys.exc_info()[0]))
00443
00444
                   return parameters, buildEmptyFrames()
00445
          # if filetype is not explicitly npy, try loading as text
00446
          else:
00447
               # Payload size as string implied by provided parameters
00448
               expectedlength = (
00449
                   4 * (lastframe - firstframe + 1) * (lastrow - firstrow + 1) * width
00450
00451
              try:
    f = open(datafile, "r")
00452
00453
                   s = f.read()
00454
00455
                   padding = expectedlength - len(s)
00456
00457
                   if padding:
00458
                       logging.warning(
                            "{logwarn}loadDumpedData: Payload is shorter than expected."
00459
                            " Padding with '0's".format(logwarn=parameters.logwarn)
00460
00461
00462
                       s = s.ljust(expectedlength, "0")
00463
00464
                   frames = generateFrames(parameters, s)
00465
                   return parameters, frames
00466
00467
               except OSError as err:
00468
                   logging.error(
                       "{logerr}loadDumpedData: OS error: {err}. Returning empty"
" frames.".format(logerr=parameters.logerr, err=err)
00469
00470
00471
00472
                   return parameters, buildEmptyFrames()
00473
               except ValueError:
00474
                   logging.error(
00475
                       "{logerr}loadDumpedData: Could not convert data to an integer."
00476
                       " Returning empty frames.".format(logerr=parameters.logerr)
00477
00478
                   return parameters, buildEmptyFrames()
00479
               except:
00480
                   logging.error(
00481
                       "{logerr}loadDumpedData: Unexpected error: {err}. Returning empty"
00482
                       " frames.".format(logerr=parameters.logerr, err=str(sys.exc_info()[0]))
00483
00484
                   return parameters, buildEmptyFrames()
00485
00486
00487 def saveTiffs(
00488
00489
          frames,
00490
          path=None,
          filename="Frame",
00491
00492
          prefix=None,
00493
          index=None,
00494):
00495
00496
          Save numpy array or list of numpy arrays or single array to disk as individual
00497
            tiffs, with frame number appended to filename. If used for standalone, use the
00498
            parameter object returned by loadDumpedData() as the first argument
00499
00500
          Args:
00501
               self: during normal operation, is cameraAssembler object. During offline
00502
                operation, is the parameter object returned by loadDumpedData()
               frames: numpy array or list of numpy arrays
00503
               path: save path, defaults to './output'
filename: defaults to 'Frame' followed by frame number
00504
00505
              prefix: prepended to 'filename', defaults to time/date (e.g. '160830-124704_')
00506
00507
```

```
00508
              index: number to start frame numbering
00509
00510
          Returns:
          Error string
00511
00512
00513
          # logging.info("INFO" + self.logtag + ": saveTiffs")
00514
          logging.info("{loginfo}: saveTiffs".format(loginfo=self.loginfo))
00515
00516
          if path is None:
00517
             path = os.path.join(os.getcwd(), "output")
00518
          if prefix is None:
00519
             prefix = datetime.now().strftime("%y%m%d-%H%M%S%f")[:-5] + "_"
00520
          if not os.path.exists(path):
00521
              os.makedirs(path)
00522
          if index is None:
00523
              firstnum = self.sensor.firstframe
00524
          else:
00525
              firstnum = index
00526
00527
          # if this is a text string from fast readoff, do the numpy conversion now
00528
          if isinstance(frames[0], str):
00529
              frames = generateFrames(frames)
00530
00531
          framestemp = np.copy(frames)
00532
          if np.issubdtype(type(framestemp[0]), np.number):
00533
              # if type(framestemp[0]) is np.uint16:
00534
              # single frame needs to be a list containing one frame
00535
              framestemp = [framestemp]
          for idx, frame in enumerate(framestemp):
00536
              if idx < len(framestemp) / 2:</pre>
00537
00538
                  interlacing = self.sensor.interlacing[0]
00539
              else:
00540
                  interlacing = self.sensor.interlacing[1]
00541
              try:
                  if self.padToFull:
00542
00543
                      frame = np.reshape(
                         frame, (self.sensor.maxheight // (interlacing + 1), -1)
00544
00545
                      )
00546
                  else:
                      frame = np.reshape(
00547
                          frame,
00548
00549
                               (self.sensor.lastrow - self.sensor.firstrow + 1)
00550
00551
                              // (interlacing + 1),
00552
                              -1,
00553
                          ),
00554
                  frameimg = Image.fromarray(frame, "I;16")
00555
00556
                  namenum = filename + "_%d" % firstnum
                  tifpath = os.path.join(path, prefix + namenum + ".tif")
00557
00558
                  frameimg.save(tifpath)
00559
                  firstnum += 1
00560
              except Exception:
00561
                  err = "saveTiffs: unable to save images"
                  # logging.error("ERROR" + self.logtag + ": " + err)
00562
00563
                  logging.error("{logerr}: {err}".format(logerr=self.logerr, err=err))
00564
          return err
00565
00566
00567 def plotFrames(self, frames, index=None):
00568
00569
          Plot frame or list of frames as individual graphs.
00570
00571
          Args:
00572
             self: during normal operation, is cameraAssembler object. During offline
00573
               operation, is the parameter object returned by loadDumpedData()
00574
              frames: numpy array or list of numpy arrays
00575
              index: number to start frame numbering
00576
00577
          Returns:
          Error string
00578
00579
00580
          # logging.info(self.loginfo + "plotFrames: index = " + str(index))
00581
          logging.info(
              "{loginfo}: plotFrames: index = {index}".format(
00582
00583
                  loginfo=self.loginfo, index=index
00584
00585
          )
          err = ""
00586
00587
          if index is None:
00588
              nframe = self.sensor.firstframe
```

7.46 misc.py 417

```
00589
00590
              nframe = index
00591
00592
          if not isinstance(frames, list):
00593
               frames = [frames]
00594
00595
          # if this is a text string from fast readoff, do the numpy conversion now
00596
          if isinstance(frames[0], str):
00597
               frames = generateFrames(frames)
00598
00599
          framestemp = np.copy(frames)
00600
          for idx, frame in enumerate(framestemp):
00601
               if idx < len(framestemp) / 2:</pre>
00602
                   interlacing = self.sensor.interlacing[0]
00603
               else:
00604
                  interlacing = self.sensor.interlacing[1]
00605
               try:
00606
                   if self.padToFull:
00607
                       frame = np.reshape(
00608
                           frame, (self.sensor.maxheight // (interlacing + 1), -1)
00609
                       )
00610
                   else:
00611
                       frame = np.reshape(
00612
                           frame,
00613
                            (
00614
                                (self.sensor.lastrow - self.sensor.firstrow + 1)
00615
                                // (interlacing + 1),
00616
                                -1,
00617
                           ),
00618
                       )
00619
               except:
00620
                   err = "{logerr}plotFrames: unable to plot frame".format(logerr=self.logerr)
00621
                   logging.error(err)
00622
               plt.imshow(frame, cmap="gray")
name = "Frame %d" % nframe
00623
00624
00625
               plt.title(name)
00626
               plt.show()
00627
               nframe += 1
00628
          return err
00629
00630
00631 \# TODO: separate images for hemispheres with different timing
00632
00633
00634 def partition(self, frames, columns):
00635
00636
          Extracts interlaced frames and divides images by hemispheres. If interlacing does
00637
            not evenly divide the height, remainder lines will be dropped
00638
00639
00640
              self: during normal operation, is sensor object. During offline
00641
                operation, is the parameter.sensor object returned by loadDumpedData()
00642
               frames: list of full-sized frames
00643
               columns: 1 for single image, 2 for separate hemisphere images
00644
00645
          Returns: list of deinterlaced frames
00646
00647
          logging.debug(
00648
               "{logdebug}partition: columns = {columns}, interlacing = {interlacing}".format(
00649
                   logdebug=self.logdebug, columns=columns, interlacing=self.sensor.interlacing
00650
00651
00652
00653
          def unshuffle(frames, ifactor):
00654
               warntrimmed = False
00655
               if self.padToFull:
00656
                   newheight = self.sensor.maxheight // (ifactor + 1)
00657
                   if newheight != (self.sensor.maxheight / (ifactor + 1)):
00658
                       warntrimmed = True
00659
               else:
00660
                   newheight = self.sensor.height // (ifactor + 1)
00661
                   if newheight != (self.sensor.height / (ifactor + 1)):
                       warntrimmed = True
00662
00663
00664
               if warntrimmed:
00665
                   logging.warning(
                       "{logwarn} partition: interlacing setting requires dropping of lines to"
" maintain consistent frame sizes ".format(logwarn=self.logwarn)
00666
00667
00668
00669
               delaced = []
```

```
for frame in frames:
                  for sub in range(ifactor + 1):
00672
                      current = np.zeros((newheight, self.sensor.width // columns), dtype=int)
00673
                      for line in range(newheight):
00674
                          current[line] = frame[(ifactor + 1) * line + sub]
00675
                     delaced.append(current)
00676
             nframes = self.sensor.lastframe - self.sensor.firstframe + 1
00677
             resorted = [None] * len(delaced)
00678
              for sub in range(ifactor + 1):
00679
                  for idx, frame in enumerate(frames):
                      resorted[sub * nframes + idx] = delaced[idx * (ifactor + 1) + sub]
00682
        if self.sensor.interlacing[0] != self.sensor.interlacing[1]:
             columns = 2 # true even if not explicitly requested by readoff
            columns == 1:
00686
              if self.sensor.interlacing == [0, 0]: # don't do anything
00687
                 return frames
00688
              else:
00689
                 return unshuffle(frames, self.sensor.interlacing[0])
00690
         else:
00691
              # reshape frame into the proper shape, then split horizontally
00692
              if self.padToFull:
00693
                  framesab = [
00694
                      np.hsplit(frame.reshape(self.sensor.maxheight, -1), 2)
00695
                      for frame in frames
00696
                 1
00697
              else:
00698
                  framesab = [
00699
                      np.hsplit(
00700
                          frame.reshape((self.sensor.lastrow - self.sensor.firstrow + 1), -1),
00701
00702
00703
                      for frame in frames
00704
                 1
              framesa = [hemis[0] for hemis in framesab]
00705
              framesb = [hemis[1] for hemis in framesab]
00706
00707
          if self.sensor.interlacing == [0, 0]:
00708
             return framesa + framesb
00709
          else:
00710
             return unshuffle(framesa, self.sensor.interlacing[0]) + unshuffle(
00711
                  framesb, self.sensor.interlacing[1]
00712
00713
00714
00715 """
00716 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00717 LLNL-CODE-838080
00718
00719 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00720 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00721 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00722 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00723 be made under this license.
00724 """
```

## 7.47 C:/Users/hill35/git/camera\_python/nsCamera/utils/Packet.py File Reference

### Classes

· class nsCamera.utils.Packet.Packet

#### **Namespaces**

- namespace nsCamera
- namespace nsCamera.utils
- namespace nsCamera.utils.Packet

7.48 Packet.py 419

### 7.48 Packet.py

```
Go to the documentation of this file.
```

```
00001 # -*- coding: utf-8 -*-
00003 Packet object for communication with boards
00005 Author: Brad Funsten (funsten1@llnl.gov)
00006 Author: Jeremy Hill (hill35@llnl.gov)
00007
00008 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00009 LLNL-CODE-838080
00011 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00012 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00013 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00014 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00015 be made under this license.
00016
00017 Version: 2.1.2 (February 2025)
00018 """
00019
00020 from __future__ import absolute_import
00021
00022 # mport binascii
00023 import sys
00024
00025 from nsCamera.utils import crc16pure
00026 from nsCamera.utils.misc import str2bvtes
00027
00028
00029 class Packet:
00030
00031
          Packet object for communication with boards. See ICD for details.
00032
00033
          Single Command/Response packet:
00034
          | 16 bits | 4 bits | 12 bits | 32 bits | 16 bits | Preamble | Command | Address | Data | CRC16 |
00035
00036
00037
00038
00039
          Read Burst Response packet:
00040
00041
          | 16 bits | 4 bits | 4 bits | 16 bits
00042
00043
          | Preamble | Command | Sub-command | Sequence ID %
00044
00045
                                                         16 bits | Variable | 16 bits
00046
                                                     용
00047
                                                     % Payload Length | Payload | CRC16
00048
00049
00050
00051
00052
00053
               # NOTE: 'numerical' components are handled as hex strings
00054
00055
               preamble="aaaa",
00056
               cmd="0",
00057
               addr=""
00058
               data="00000000",
               seqID="",
00059
00060
               payload_length="",
00061
              payload="",
               crc="",
00062
00063
00064
               self.PY3 = sys.version_info > (3,)
00065
               self.preamble = preamble # 16 bit packet preamble
               self.cmd = str(cmd) # 4 bit command packet
00066
               self.addr = addr.zfill(3) # 12 bit address packet
self.data = data.zfill(8) # 32 bit data packet
00067
00068
00069
               # 16 bit sequence ID packet (only Read Burst)
00070
               self.seqID = seqID
               # 16 bit payload packet (only Read Burst)
self.payload_length = payload_length
00071
00072
00073
               \mbox{\#} variable payload packet (only Read Burst) for now it's 16 bits
              self.payload = payload
# 16 bit CRC-CCIT (XModem) packet
00074
00075
00076
               self.crc = crc
```

```
00077
              self.type = ""
               if self.crc == "": # check if packet to be sent needs crc appended
00078
00079
                   self.crc = self.calculateCRC()
00080
00081
          def pktStr(self):
00082
00083
              Generate hexadecimal string form of packet
00084
00085
              packet as hexadecimal string without '0x' """
00086
00087
00088
              if self.seqID != "":
00089
                   # Read burst response
00090
                   packetparts = [
00091
                       self.preamble,
00092
                       self.cmd,
00093
                       self.seqID,
00094
                       self.payload_length,
00095
                       self.payload,
00096
                       self.crc,
00097
                  1
00098
              else:
00099
                  # Single Command/Response response
00100
                  packetparts = [self.preamble, self.cmd, self.addr, self.data, self.crc]
00101
              stringparts = [
                   part.decode("ascii") if isinstance(part, bytes) else part
00102
00103
                   \quad \quad \text{for part in packetparts} \quad
00104
              out = "".join(stringparts)
00105
00106
              return out
00107
00108
          def calculateCRC(self):
00109
00110
              Calculate CRC-CCIT (XModem) (2 bytes) from 8 byte packet for send and rcv
00111
00112
              Returns:
              CRC as hexadecimal string without '0x' """
00113
00114
              preamble = self.preamble
00115
              crc = self.crc
self.crc = ""
00116
00117
              self.preamble = ""
00118
00119
00120
              CRC_dec = crc16pure.crc16xmodem(str2bytes(self.pktStr()))
              \# input = int type decimal, output = hex string with 0x at the beginning CRC_hex_0x = "0x%0.4X" % CRC_dec
00121
00122
00123
               # make all hex letters lower case for comparison
00124
              CRC_hex = CRC_hex_0x.lower()
00125
               \# input = hex string with 0x at the beginning, output = hex str with 0x removed
00126
              CRC_hex = CRC_hex[2:]
00127
              self.preamble = preamble
00128
              self.crc = crc
00129
              return CRC_hex
00130
00131
          def checkCRC(self):
00132
00133
              Returns: boolean, True if CRC check passes
00134
00135
               return self.calculateCRC() == self.crc
00136
00137
          def checkReadPacket(self, resppkt):
00138
00139
               Confirm that Read Single occurred without error
00140
              Args:
00141
                  resppkt: response packet
00142
00143
              tuple (error string, response packet as string) _{\rm m\,m\,m}
00144
00145
              err = ""
00146
00147
              if int(resppkt.cmd.upper(), 16) - int(self.cmd.upper(), 16) != 0x8:
00148
                  err = "invalid command; "
00149
               if resppkt.addr.upper() != self.addr.upper():
00150
                  err += "invalid address; "
00151
              if resppkt.crc.upper() != resppkt.calculateCRC().upper():
                  err += "invalid CRC; "
00152
              return err, resppkt.pktStr()
00153
00154
00155
          def checkResponsePacket(self, resppkt):
00156
00157
              Confirm that Write Single occurred without error
```

```
00158
              Args:
00159
                 resppkt: response packet
00161
             Returns:
             tuple (error string, response packet as string)
00163
             err = ""
00165
             if int(resppkt.data, 16) & 1:
00166
                 err += "Checksum error; "
00167
              if int(resppkt.data, 16) & 2:
                 err += "Invalid command / command not executed; "
              err1, rval = self.checkReadPacket(resppkt)
00169
00170
             err += err1
00171
             return err, rval
00172
         def checkResponseString(self, respstr):
00174
00175
              Checks response string for error indicators
00176
00177
                 respstr: packet as hexadecimal string
00178
00179
             Returns:
             tuple (error string, response packet string)
00180
00181
00182
             respstring = respstr.decode(encoding="UTF-8")
00183
             resppkt = Packet(
00184
                 preamble=respstring[0:4],
00185
                 cmd=respstring[4],
00186
                 addr=respstring[5:8]
00187
                 data=respstring[8:16],
            )
00188
00189
00190
            if resppkt.cmd == "8":
00191
                 # verify response to write command
00192
             err, rval = self.checkResponsePacket(resppkt)
elif resppkt.cmd == "9":
00193
00194
                 err, rval = self.checkReadPacket(resppkt) # verify response to read command
00195
             else:
                 err = "Packet command invalid; "
00196
                 rval = ""
00197
00198
             return err, rval
00199
00200
00201 """
00202 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00203 LLNL-CODE-838080
00204
00205 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00206 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00207 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00208 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00209 be made under this license.
00210 """
```

# 7.49 C:/Users/hill35/git/camera\_python/nsCamera/utils/Subregister.py File Reference

#### Classes

· class nsCamera.utils.Subregister.SubRegister

#### **Namespaces**

- namespace nsCamera
- namespace nsCamera.utils
- · namespace nsCamera.utils.Subregister

### 7.50 Subregister.py

#### Go to the documentation of this file.

```
00001 # -*- coding: utf-8 -*-
00002 """
00003 Subregister object represents a subset of a full register
00004
00005 Author: Matthew Dayton (dayton5@llnl.gov)
00006 Author: Jeremy Martin Hill (jerhill@llnl.gov)
00007
00008 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00009 LLNL-CODE-838080
00010
00011 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00012 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00013 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00014 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00015 be made under this license.
00016
00017 Version: 2.1.2 (February 2025)
00018 """
00019
00020
00021 class SubRegister:
00022
00023
          Represents a subset of a 32-bit register [31..0] starting at 'start_bit' consisting
00024
            of 'width' bits. Consistent with the ICD usage, start_bit is MSB e.g., for [7..0],
00025
            the start_bit is '7'.
00026
00027
00028
          def __init__(
00029
              self,
00030
              board,
00031
              name,
00032
              register,
00033
              start_bit=31,
00034
              width=8,
00035
              writable=False,
00036
              value=255,
              minV=0,
00037
00038
              maxV=5,
00039
00040
              self.name = name
00041
              self.register = register
00042
              self.addr = board.registers[register]
00043
              self.start_bit = start_bit
00044
              self.width = width
              self.value = value
00045
00046
              self.max\_value = 2**width - 1 # used to normalize the input values to 1
00047
              self.min = 0
00048
              self.max = self.max_value
00049
              self.writable = writable
00050
              self.minV = minV
00051
              self.maxV = maxV
              \# resolution should be reset after init if actual min and max are different
00052
              self.resolution = (1.0 * maxV - minV) / self.max value
00053
00054
00055
00056 """
00057 Copyright (c) 2025, Lawrence Livermore National Security, LLC. All rights reserved.
00058 LLNL-CODE-838080
00060 This work was produced at the Lawrence Livermore National Laboratory (LLNL) under
00061 contract no. DE-AC52-07NA27344 (Contract 44) between the U.S. Department of Energy (DOE)
00062 and Lawrence Livermore National Security, LLC (LLNS) for the operation of LLNL.
00063 'nsCamera' is distributed under the terms of the MIT license. All new contributions must
00064 be made under this license. 00065 """
```

## Index

_CRT_SECURE_NO_WARNINGS	addr
Data.c, 304	nsCamera.utils.Packet.Packet, 196
UPnP.c, 344	nsCamera.utils.Subregister.SubRegister, 229
all	arch
nsCamera, 10	nsCamera.CameraAssembler.CameraAssembler, 89
nsCamera.boards, 10	args
nsCamera.comms, 13	nsCamera.utils.FlatField, 22
nsCamera.sensors, 14	arm
nsCamera.utils, 17	nsCamera.CameraAssembler.CameraAssembler, 47
init	nsCamera.comms.GigE.GigE, 125
nsCamera.boards.LLNL_v1.llnl_v1, 153	nsCamera.comms.RS422.RS422, 201
nsCamera.boards.LLNL_v4.llnl_v4, 173	armed
nsCamera.CameraAssembler.CameraAssembler, 43	nsCamera.CameraAssembler.CameraAssembler, 89
nsCamera.comms.GigE.GigE, 124	
nsCamera.comms.RS422.RS422, 199	batchAcquire
nsCamera.sensors.daedalus.daedalus, 99	nsCamera.CameraAssembler.CameraAssembler, 47
nsCamera.sensors.icarus.icarus, 139	baud
nsCamera.sensors.icarus2.icarus2, 146	nsCamera.comms.RS422.RS422, 208
nsCamera.sensors.sensorBase.sensorBase, 212	board
nsCamera.utils.misc.fakeCA, 119	nsCamera.CameraAssembler.CameraAssembler, 89
nsCamera.utils.Packet.Packet, 193	boardname
nsCamera.utils.Subregister.SubRegister, 229	nsCamera.CameraAssembler.CameraAssembler, 89
_crc16	nsCamera.utils.misc.fakeCA, 120
nsCamera.utils.crc16pure, 18	bytes2str
_fields_	nsCamera.CameraAssembler.CameraAssembler, 48
nsCamera.comms.GigE.GigE.ZESTETM1_CARD_INF	rO, nsCamera.utils.misc, 24
232	bytesperpixel
_ser	nsCamera.sensors.daedalus.daedalus, 113
nsCamera.comms.RS422.RS422, 208	nsCamera.sensors.icarus.icarus, 141
_strnicmp	nsCamera.sensors.icarus2.icarus2, 148
UPnP.c, 344	0.41. 4.410.
_zest	C:/Users/hill35/git/camera_python/nsCamera/initpy,
nsCamera.comms.GigE.GigE, 133	233
	C:/Users/hill35/git/camera_python/nsCamera/boards/initpy,
abort	234
nsCamera.CameraAssembler.CameraAssembler, 89	C:/Users/hill35/git/camera_python/nsCamera/boards/LLNL_v1.py,
abortReadoff	237, 238
nsCamera.CameraAssembler.CameraAssembler, 46	C:/Users/hill35/git/camera_python/nsCamera/boards/LLNL_v4.py,
action	249
nsCamera.utils.FlatField, 22	C:/Users/hill35/git/camera_python/nsCamera/CameraAssembler.py
ADC5_bipolar	260, 261
nsCamera.boards.LLNL_v1.llnl_v1, 169	C:/Users/hill35/git/camera_python/nsCamera/comms/initpy,
nsCamera.boards.LLNL_v4.llnl_v4, 189	235
ADC5_mult	C:/Users/hill35/git/camera_python/nsCamera/comms/GigE.py,
nsCamera.boards.LLNL_v1.llnl_v1, 169	287
net amora hoarde I I NII I V/I IInI V/I 180	

C:/Users/hill35/git/camera_python/nsCamera/comms/RS422	2.py, nsCamera.utils.Packet.Packet, 194 checkReadPacket
C:/Users/hill35/git/camera_python/nsCamera/comms/ZestE- 301, 315	TM1/Dæt6æmera.utils.Packet.Packet, 194 checkRegSet
C:/Users/hill35/git/camera_python/nsCamera/comms/ZestE <sup>-</sup> 326, 327	TM1/Ensocamera.CameraAssembler.CameraAssembler, 49
C:/Users/hill35/git/camera_python/nsCamera/comms/ZestE	
329, 330 C:/Users/hill35/git/camera_python/nsCamera/comms/ZestE	checkResponseString TM1/Pris@tentera.utils.Packet.Packet, 195
331, 340	checkSensorVoltStat
C:/Users/hill35/git/camera_python/nsCamera/comms/ZestE <sup>-</sup> 342, 354	
C:/Users/hill35/git/camera_python/nsCamera/comms/ZestE	nsCamera.sensors.icarus.icarus, 140
362, 375	checkStatus
C:/Users/hill35/git/camera_python/nsCamera/sensors/init_	
235, 236	nsCamera.boards.LLNL_v4.llnl_v4, 175
C:/Users/hill35/git/camera_python/nsCamera/sensors/daeda 377	checkStatus2
C:/Users/hill35/git/camera_python/nsCamera/sensors/icarus 387	s.py, nsCamera.boards.LLNL_v1.llnl_v1, 155 nsCamera.boards.LLNL_v4.llnl_v4, 175
C:/Users/hill35/git/camera_python/nsCamera/sensors/icarus	s2.py, nsCamera.CameraAssembler.CameraAssembler, $50$
390	clearStatus
C:/Users/hill35/git/camera_python/nsCamera/sensors/senso	
392	nsCamera.boards.LLNL_v4.llnl_v4, 175
C:/Users/hill35/git/camera_python/nsCamera/utils/initp 236, 237	by, nsCamera.CameraAssembler.CameraAssembler, 50 closecard
C:/Users/hill35/git/camera_python/nsCamera/utils/crc16pure	
401, 402	closeDevice
C:/Users/hill35/git/camera_python/nsCamera/utils/FlatField.  405, 406	py, nsCamera.CameraAssembler.CameraAssembler, 50 nsCamera.comms.GigE.GigE, 126
C:/Users/hill35/git/camera_python/nsCamera/utils/misc.py, 409	nsCamera.comms.RS422.RS422, 201 closesocket
C:/Users/hill35/git/camera_python/nsCamera/utils/Packet.py	
C:/Users/hill35/git/camera_python/nsCamera/utils/Subregist	
421, 422	columns
ca	nsCamera.sensors.daedalus.daedalus, 114
nsCamera.boards.LLNL_v1.llnl_v1, 169	nsCamera.sensors.icarus.icarus, 141
nsCamera.boards.LLNL_v4.llnl_v4, 189	nsCamera.sensors.icarus2.icarus2, 148
nsCamera.comms.GigE.GigE, 133	commname
nsCamera.comms.RS422.RS422, 208	nsCamera.CameraAssembler.CameraAssembler, 89
nsCamera.sensors.daedalus.daedalus, 113	comms
nsCamera.sensors.icarus.icarus, 141 nsCamera.sensors.icarus2.icarus2, 148	nsCamera.CameraAssembler.CameraAssembler, 90 configADCs
nsCamera.sensors.sensorBase.sensorBase, 227	nsCamera.boards.LLNL v1.llnl v1, 156
calculateCRC	nsCamera.boards.LLNL v4.llnl v4, 176
nsCamera.utils.Packet.Packet, 193	nsCamera.CameraAssembler.CameraAssembler, 51
CardInfo	Connection
nsCamera.comms.GigE.GigE, 133	nsCamera.comms.GigE.GigE, 134
CardInfoP	crc
nsCamera.comms.GigE.GigE, 133 checkCRC	nsCamera.utils.Packet.Packet, 197 CRC16_XMODEM_TABLE
nsCamera.CameraAssembler.CameraAssembler, 49	nsCamera.utils.crc16pure, 18
nsCamera.utils.misc, 25	crc16xmodem

nsCamera.utils.crc16pure, 18	dport
currtime	nsCamera.comms.GigE.GigE, 134
nsCamera.CameraAssembler.CameraAssembler, 90	dumpNumpy
cycle	nsCamera.CameraAssembler.CameraAssembler, 51
nsCamera.CameraAssembler.CameraAssembler, 90	dumpRegisters
	nsCamera.CameraAssembler.CameraAssembler, 52
daedalus_monitor_controls	dumpStatus
nsCamera.boards.LLNL_v1.llnl_v1, 169	nsCamera.boards.LLNL_v1.llnl_v1, 157
nsCamera.boards.LLNL_v4.llnl_v4, 189	nsCamera.boards.LLNL_v4.llnl_v4, 177
daedalus_subreg_aliases	nsCamera.CameraAssembler.CameraAssembler, 53
nsCamera.boards.LLNL_v1.llnl_v1, 169	dumpSubregisters
nsCamera.boards.LLNL_v4.llnl_v4, 189	nsCamera.CameraAssembler.CameraAssembler, 53
data	11.150
nsCamera.utils.Packet.Packet, 197	enableLED
Data.c	nsCamera.boards.LLNL_v1.llnl_v1, 158
_CRT_SECURE_NO_WARNINGS, 304	nsCamera.boards.LLNL_v4.llnl_v4, 178
ZestETM1_CloseConnection, 305	nsCamera.CameraAssembler.CameraAssembler, 53
ZESTETM1_COMMAND_MAILBOX_INT, 304	Error.c
ZESTETM1_COMMAND_READ_REG, 304	ZestETM1_ErrorHandler, 327
ZESTETM1_COMMAND_SPI, 304	ZestETM1_ErrorStrings, 327
ZESTETM1_COMMAND_WRITE_REG, 304	ZestETM1GetErrorMessage, 326
ZestETM1_OpenConnection, 305	ZestETM1RegisterErrorHandler, 326
ZestETM1_ReadData, 307	firstframe
ZestETM1_SendCommand, 308	nsCamera.sensors.daedalus.daedalus, 114
ZestETM1_SPIReadWrite, 308	nsCamera.sensors.icarus.icarus, 141
ZestETM1_WriteData, 309	nsCamera.sensors.icarus2.icarus2, 148
ZestETM1CloseConnection, 311	firstrow
ZestETM1OpenConnection, 311	nsCamera.sensors.daedalus.daedalus, 114
ZestETM1ReadData, 311	nsCamera.sensors.icarus.icarus, 141
ZestETM1ReadRegister, 312	nsCamera.sensors.icarus2.icarus2, 148
ZestETM1SetInterrupt, 313	flattenlist
ZestETM1SPIReadWrite, 313	nsCamera.CameraAssembler.CameraAssembler, 54
ZestETM1WriteData, 314	nsCamera.utils.misc, 25
ZestETM1WriteRegister, 314 datatimeout	FPGAboardtype
nsCamera.comms.RS422.RS422, 208	nsCamera.CameraAssembler.CameraAssembler, 90
default	FPGAinterfaces
nsCamera.utils.FlatField, 22	nsCamera.CameraAssembler.CameraAssembler, 90
defoff	FPGAinvalid
nsCamera.boards.LLNL_v4.llnl_v4, 189	nsCamera.CameraAssembler.CameraAssembler, 90
defsens	FPGANum
nsCamera.boards.LLNL v4.llnl v4, 189	nsCamera.CameraAssembler.CameraAssembler, 90
dest	fpganumID
nsCamera.utils.FlatField, 22	nsCamera.sensors.daedalus.daedalus, 114
detect	nsCamera.sensors.icarus.icarus, 142
nsCamera.sensors.daedalus.daedalus, 114	nsCamera.sensors.icarus2.icarus2, 149
nsCamera.sensors.icarus.icarus, 141	nsCamera.sensors.sensorBase.sensorBase, 227
nsCamera.sensors.icarus2.icarus2, 148	FPGArad
directory	nsCamera.CameraAssembler.CameraAssembler, 91
nsCamera.utils.FlatField, 23	FPGAsensor
disarm	nsCamera.CameraAssembler.CameraAssembler, 91
nsCamera.boards.LLNL_v1.llnl_v1, 156	FPGAVersion
nsCamera.boards.LLNL_v4.llnl_v4, 176	nsCamera.CameraAssembler.CameraAssembler, 91
nsCamera.CameraAssembler.CameraAssembler, 51	framelist
oamora.oamora losombior.oamera/issembior, 31	nsCamera.utils.FlatField, 23

generateFF	nsCamera.sensors.daedalus.daedalus, 114
nsCamera.utils.FlatField, 20	nsCamera.sensors.icarus.icarus, 142
generateFrames	nsCamera.sensors.icarus2.icarus2, 149
nsCamera.utils.misc, 26	help
getBoardInfo	nsCamera.utils.FlatField, 23
nsCamera.CameraAssembler.CameraAssembler, 54	HFW
getCardInfo	nsCamera.sensors.daedalus.daedalus, 115
nsCamera.comms.GigE.GigE, 126	nocamora.concoro.acoualac.acoualac, 110
getCardIP	icarus_monitor_controls
nsCamera.comms.GigE.GigE, 127	nsCamera.boards.LLNL_v1.llnl_v1, 169
	nsCamera.boards.LLNL_v4.llnl_v4, 189
getEnter	icarus_subreg_aliases
nsCamera.CameraAssembler.CameraAssembler, 55	nsCamera.boards.LLNL_v1.llnl_v1, 169
nsCamera.utils.misc, 27	nsCamera.boards.LLNL_v4.llnl_v4, 190
getFilenames	icarustype
nsCamera.utils.FlatField, 20	nsCamera.sensors.icarus.icarus, 142
getManualTiming	nsCamera.sensors.icarus2.icarus2, 149
nsCamera.CameraAssembler.CameraAssembler, 55	init_board_specific
nsCamera.sensors.daedalus.daedalus, 100	nsCamera.sensors.sensorBase.sensorBase, 216
nsCamera.sensors.sensorBase.sensorBase, 213	initBoard
getMonV	
nsCamera.CameraAssembler.CameraAssembler, 55	nsCamera.boards.LLNL_v1.llnl_v1, 160
getPot	nsCamera.boards.LLNL_v4.llnl_v4, 181
nsCamera.CameraAssembler.CameraAssembler, 56	nsCamera.CameraAssembler.CameraAssembler, 61
getPotV	initialize
nsCamera.CameraAssembler.CameraAssembler, 57	nsCamera.CameraAssembler.CameraAssembler, 61
getPressure	initPots
nsCamera.boards.LLNL_v1.llnl_v1, 159	nsCamera.boards.LLNL_v1.llnl_v1, 161
nsCamera.boards.LLNL_v4.llnl_v4, 179	nsCamera.boards.LLNL_v4.llnl_v4, 181
nsCamera.CameraAssembler.CameraAssembler, 58	nsCamera.CameraAssembler.CameraAssembler, 63
getRegister	initPowerCheck
nsCamera.CameraAssembler.CameraAssembler, 59	nsCamera.CameraAssembler.CameraAssembler, 63
getROIvector	initSensor
nsCamera.utils.FlatField, 20	nsCamera.boards.LLNL_v1.llnl_v1, 162
getSensorStatus	nsCamera.boards.LLNL_v4.llnl_v4, 181
nsCamera.sensors.sensorBase.sensorBase, 214	nsCamera.CameraAssembler.CameraAssembler, 64
getSensTemp	inittime
nsCamera.CameraAssembler.CameraAssembler, 59	nsCamera.CameraAssembler.CameraAssembler, 91
nsCamera.sensors.daedalus.daedalus, 101	interlacing
nsCamera.sensors.sensorBase.sensorBase, 214	nsCamera.sensors.daedalus.daedalus, 115
getSubregister	nsCamera.sensors.icarus.icarus, 142
nsCamera.CameraAssembler.CameraAssembler, 60	nsCamera.sensors.icarus2.icarus2, 149
getTemp	iplist
nsCamera.boards.LLNL_v1.llnl_v1, 159	nsCamera.CameraAssembler.CameraAssembler, 91
nsCamera.boards.LLNL_v4.llnl_v4, 180	lastframe
nsCamera.CameraAssembler.CameraAssembler, 61	nsCamera.sensors.daedalus.daedalus, 115
getTimer	nsCamera.sensors.icarus.icarus, 142
nsCamera.boards.LLNL_v1.llnl_v1, 160	nsCamera.sensors.icarus2.icarus2, 149
nsCamera.boards.LLNL_v4.llnl_v4, 180	lastrow
nsCamera.CameraAssembler.CameraAssembler, 61	nsCamera.sensors.daedalus.daedalus, 115
getTiming	nsCamera.sensors.icarus.icarus, 142
nsCamera.CameraAssembler.CameraAssembler, 61	nsCamera.sensors.icarus.icarus2, 149
nsCamera.sensors.sensorBase.sensorBase, 214	latchPots
	nsCamera.boards.LLNL_v1.llnl_v1, 162
height	nsCamera.boards.LLNL_v4.llnl_v4, 182
	1130amera.boards.cerve_v4.11111_v4, 102

nsCamera.CameraAssembler.CameraAssembler, 64	logwarn
loadDumpedData	nsCamera.boards.LLNL_v1.llnl_v1, 170
nsCamera.utils.misc, 27	nsCamera.boards.LLNL_v4.llnl_v4, 190
logcrit	nsCamera.CameraAssembler.CameraAssembler, 93
nsCamera.boards.LLNL_v1.llnl_v1, 169	nsCamera.comms.GigE.GigE, 135
nsCamera.boards.LLNL_v4.llnl_v4, 190	nsCamera.comms.RS422.RS422, 209
nsCamera.CameraAssembler.CameraAssembler, 91	nsCamera.sensors.sensorBase.sensorBase, 228
nsCamera.comms.GigE.GigE, 134	nsCamera.utils.misc.fakeCA, 121
nsCamera.comms.RS422.RS422, 209	logwarnbase
nsCamera.sensors.sensorBase.sensorBase, 227	nsCamera.CameraAssembler.CameraAssembler, 93
nsCamera.utils.misc.fakeCA, 120	nsCamera.utils.misc.fakeCA, 121
logcritbase	
nsCamera.CameraAssembler.CameraAssembler, 91	Main.c
nsCamera.utils.misc.fakeCA, 120	ZestETM1Close, 329
logdebug	ZestETM1Init, 329
nsCamera.boards.LLNL_v1.llnl_v1, 170	makeLogLabels
nsCamera.boards.LLNL_v4.llnl_v4, 190	nsCamera.utils.misc, 30
nsCamera.CameraAssembler.CameraAssembler, 92	max
nsCamera.comms.GigE.GigE, 134	nsCamera.utils.Subregister.SubRegister, 229
nsCamera.comms.RS422.RS422, 209	max_value
nsCamera.sensors.sensorBase.sensorBase, 227	nsCamera.utils.Subregister.SubRegister, 230
nsCamera.utils.misc.fakeCA, 120	maxframe
logdebugbase	nsCamera.sensors.daedalus.daedalus, 116
nsCamera.CameraAssembler.CameraAssembler, 92	nsCamera.sensors.icarus.icarus, 143
nsCamera.utils.misc.fakeCA, 120	nsCamera.sensors.icarus2.icarus2, 150
logerr	maxheight
nsCamera.boards.LLNL_v1.llnl_v1, 170	nsCamera.sensors.daedalus.daedalus, 116
nsCamera.boards.LLNL_v4.llnl_v4, 190	nsCamera.sensors.icarus.icarus, 143
nsCamera.CameraAssembler.CameraAssembler, 92	nsCamera.sensors.icarus2.icarus2, 150
nsCamera.comms.GigE.GigE, 134	maxV
nsCamera.comms.RS422.RS422, 209	nsCamera.utils.Subregister.SubRegister, 230
nsCamera.sensors.sensorBase.sensorBase, 227	maxwidth
nsCamera.utils.misc.fakeCA, 120	nsCamera.sensors.daedalus.daedalus, 116 nsCamera.sensors.icarus.icarus, 143
logerrbase	,
nsCamera.CameraAssembler.CameraAssembler, 92	nsCamera.sensors.icarus2.icarus2, 150 min
nsCamera.utils.misc.fakeCA, 121	nsCamera.utils.Subregister.SubRegister, 230
loginfo	minframe
nsCamera.boards.LLNL_v1.llnl_v1, 170	nsCamera.sensors.daedalus.daedalus, 116
nsCamera.boards.LLNL_v4.llnl_v4, 190	nsCamera.sensors.icarus.icarus, 143
nsCamera.CameraAssembler.CameraAssembler, 92	nsCamera.sensors.icarus2.icarus2, 150
nsCamera.comms.GigE.GigE, 134	minV
nsCamera.comms.RS422.RS422, 209	nsCamera.utils.Subregister.SubRegister, 230
nsCamera.sensors.sensorBase.sensorBase, 228	mmReadoff
nsCamera.utils.misc.fakeCA, 121	nsCamera.CameraAssembler.CameraAssembler, 64
loginfobase	mode
nsCamera.CameraAssembler.CameraAssembler, 92 nsCamera.utils.misc.fakeCA, 121	nsCamera.comms.GigE.GigE, 135
	nsCamera.comms.RS422.RS422, 209
loglabel	
nsCamera.sensors.daedalus.daedalus, 115 nsCamera.sensors.icarus.icarus, 142	name
nsCamera.sensors.icarus2.icarus2, 149	nsCamera.utils.Subregister.SubRegister, 230
logtag	nargs
nsCamera.CameraAssembler.CameraAssembler, 92	nsCamera.utils.FlatField, 23
nsCamera.utils.misc.fakeCA, 121	nframes
	nsCamera.sensors.daedalus.daedalus, 116

nsCamera.sensors.icarus.icarus, 143	ADC5_bipolar, 189
nsCamera.sensors.icarus2.icarus2, 150	ADC5_mult, 189
nsCamera, 9	ca, 189
all, 10	checkStatus, 175
nsCamera.boards, 10	checkStatus2, 175
all, 10	clearStatus, 175
nsCamera.boards.LLNL_v1, 11	configADCs, 176
nsCamera.boards.LLNL_v1.llnl_v1, 152	daedalus_monitor_controls, 189
init, 153	daedalus subreg aliases, 189
ADC5_bipolar, 169	defoff, 189
ADC5_mult, 169	defsens, 189
ca, 169	disarm, 176
checkStatus, 155	dumpStatus, 177
checkStatus2, 155	enableLED, 178
clearStatus, 155	getPressure, 179
configADCs, 156	getTemp, 180
daedalus_monitor_controls, 169	getTimer, 180
daedalus_subreg_aliases, 169	icarus_monitor_controls, 189
disarm, 156	icarus_subreg_aliases, 190
dumpStatus, 157	initBoard, 181
enableLED, 158	initPots, 181
getPressure, 159	initSensor, 181
getTemp, 159	latchPots, 182
getTimer, 160	logcrit, 190
icarus_monitor_controls, 169	logdebug, 190
icarus_subreg_aliases, 169	logerr, 190
initBoard, 160	loginfo, 190
initPots, 161	logwarn, 190
init dis, 101 initSensor, 162	readSRAM, 183
latchPots, 162	registers, 190
	reportEdgeDetects, 183
logdobus 170	reportStatus, 184
logdebug, 170 logerr, 170	resetTimer, 185
•	•
loginfo, 170	rs422_baud, 191
logwarn, 170	rs422_cmd_wait, 191
readSRAM, 163	setLED, 185
registers, 170	setPowerSave, 185
reportEdgeDetects, 163	setPPER, 186
reportStatus, 164	softReboot, 186
resetTimer, 164	startCapture, 187
rs422_baud, 170	subregisters, 191
rs422_cmd_wait, 170	subreglist, 191
setLED, 165	VREF, 191
setPowerSave, 165	waitForSRAM, 188
setPPER, 166	nsCamera.CameraAssembler, 11
softReboot, 166	nsCamera.CameraAssembler.CameraAssembler, 39
startCapture, 167	init, 43
subregisters, 171	abort, 89
subreglist, 171	abortReadoff, 46
VREF, 171	arch, 89
waitForSRAM, 168	arm, 47
nsCamera.boards.LLNL_v4, 11	armed, 89
nsCamera.boards.LLNL_v4.llnl_v4, 171	batchAcquire, 47
init, 173	board, 89

boardname, 89	logerrbase, 92
bytes2str, 48	loginfo, <mark>92</mark>
checkCRC, 49	loginfobase, 92
checkRegSet, 49	logtag, <mark>92</mark>
checkSensorVoltStat, 50	logwarn, <mark>93</mark>
checkStatus, 50	logwarnbase, 93
checkStatus2, 50	mmReadoff, 64
clearStatus, 50	oldtime, 93
closeDevice, 50	packageroot, 93
commname, 89	padToFull, 93
comms, 90	parsedtime, 93
configADCs, 51	parseReadoff, 65
currtime, 90	partition, 65
cycle, 90	payloaderror, 93
disarm, 51	platform, 94
dumpNumpy, 51	plotFrames, 65
dumpRegisters, 52	port, 94
dumpStatus, 53	powerCheck, 65
dumpSubregisters, 53	printBoardInfo, 66
enableLED, 53	PY3, 94
flattenlist, 54	pyth1, 94
FPGAboardtype, 90	pyth2, 94
FPGAinterfaces, 90	python, 94
FPGAinvalid, 90	read, 94
FPGANum, 90	readFrames, 67
FPGArad, 91	readImgs, 67
FPGAsensor, 91	readoff, 67
FPGAVersion, 91	readSerial, 68
getBoardInfo, 54	readSRAM, 68
getEnter, 55	reboot, 68
getManualTiming, 55	reinitialize, 68
getMonV, 55	reportEdgeDetects, 69
getPot, 56	reportStatus, 69
getPotV, 57	resetTimer, 69
getPressure, 58	resolveSubreg, 69
getRegister, 59	saveFrames, 70
getSensTemp, 59	saveHDF, 71
getSubregister, 60	saveNumpys, 72
getTemp, 61	saveTiffs, 74
getTimer, 61	savetime, 95
getTiming, 61	selectOscillator, 74
initBoard, 61	sendCMD, 74
initialize, 61	sensmanual, 95
initPots, 63	sensor, 95
initPowerCheck, 63	sensorname, 95
initSensor, 64	sensorSpecific, 74
inittime, 91	senstiming, 95
iplist, 91	setArbTiming, 74
latchPots, 64	setExtClk, 75
logcrit, 91	setFrames, 75
logcritbase, 91	setHighFullWell, 76
logdebug, 92	setInterlacing, 76
logdebugbase, 92	setLED, 77
logerr, 92	setManualShutters, 77
	,

setManualTiming, 77	sendCMD, 131
setPhiDelay, 77	skipError, 135
setPot, 77	writeSerial, 132
setPotV, 79	writeTimeout, 135
setPowerSave, 82	ZCountCards, 135
setPPER, 82	ZErrorDict, 136
setRegister, 82	ZESTETM1_CARD_INFO, 136
setRows, 83	ZOpenConnection, 136
setSubregister, 85	ZReadData, 136
setTiming, 86	ZWriteData, 136
setTriggerDelay, 86	nsCamera.comms.GigE.GigE.ZESTETM1_CARD_INFO,
setZeroDeadTime, 87	231
startCapture, 87	_fields_, 232
str2bytes, 87	ubyte4, 232
str2nparray, 87	ubyte6, 232
submitMessages, 87	nsCamera.comms.RS422, 13
timeout, 95	nsCamera.comms.RS422.RS422, 198
trigtime, 95	init, 199
unstringed, 96	_ser, 208
verblevel, 96	arm, 201
verbmap, 96	baud, 208
verbose, 96	ca, 208
version, 96	closeDevice, 201
waited, 96	datatimeout, 208
waitForSRAM, 88	logcrit, 209
writeSerial, 88	logdebug, 209
nsCamera.comms, 12	logerr, 209
all, 13	loginfo, 209
nsCamera.comms.GigE, 13	logwarn, 209
nsCamera.comms.GigE.GigE, 122	mode, 209
init, 124	par, 209
zest, 133	payloadsize, 210
arm, 125	port, 210
ca, 133	PY3, 210
CardInfo, 133	read_timeout, 210
CardInfoP, 133	readFrames, 202
closecard, 134	readoff, 203
closeDevice, 126	readSerial, 204
Connection, 134	sendCMD, 205
dport, 134	serialClose, 207
getCardInfo, 126	skipError, 210
getCardIP, 127	stop, 210
-	• •
logcrit, 134 logdebug, 134	write_timeout, 210
	writeSerial, 207
logerr, 134	nsCamera.sensors, 14
loginfo, 134	_all, 14
logwarn, 135	nsCamera.sensors.daedalus, 14
mode, 135	nsCamera.sensors.daedalus.daedalus, 97
openDevice, 127	init, 99
payloadsize, 135	bytesperpixel, 113
readFrames, 128	ca, 113
readoff, 129	columns, 114
readSerial, 130	detect, 114
readTimeout, 135	firstframe, 114

firstrow, 114	maxheight, 143
fpganumID, 114	maxwidth, 143
getManualTiming, 100	minframe, 143
getSensTemp, 101	nframes, 143
height, 114	padToFull, 143
HFW, 115	sens_registers, 143
interlacing, 115	sens_subregisters, 144
lastframe, 115	sensfam, 144
lastrow, 115	sensorSpecific, 140
loglabel, 115	specwarn, 144
maxframe, 116	width, 144
maxheight, 116	nsCamera.sensors.icarus2, 15
maxwidth, 116	nsCamera.sensors.icarus2.icarus2, 144
minframe, 116	init, 146
nframes, 116	bytesperpixel, 148
padToFull, 116	ca, 148
parseReadoff, 102	columns, 148
reportStatusSensor, 103	detect, 148
selectOscillator, 104	firstframe, 148
sens_registers, 116	firstrow, 148
sens_subregisters, 117	fpganumID, 149
sensfam, 117	height, 149
sensorSpecific, 104	icarustype, 149
setExtClk, 106	interlacing, 149
setHighFullWell, 107	lastframe, 149
setInterlacing, 107	lastrow, 149
setManualShutters, 109	loglabel, 149
setManualTiming, 109	maxframe, 150
setPhiDelay, 110	maxheight, 150
setTriggerDelay, 111	maxwidth, 150
setZeroDeadTime, 112	minframe, 150
	nframes, 150
specwarn, 117	
toffset, 117	padToFull, 150
tslope, 117	sens_registers, 150
width, 117	sens_subregisters, 151
ZDT, 117, 118	sensfam, 151
nsCamera.sensors.icarus, 15	sensorSpecific, 147
nsCamera.sensors.icarus.icarus, 137	specwarn, 151
init, 139	width, 151
bytesperpixel, 141	nsCamera.sensors.sensorBase, 16
ca, 141	nsCamera.sensors.sensorBase.sensorBase, 211
checkSensorVoltStat, 140	init, 212
columns, 141	ca, 227
detect, 141	checkSensorVoltStat, 213
firstframe, 141	fpganumID, 227
firstrow, 141	getManualTiming, 213
fpganumID, 142	getSensorStatus, 214
height, 142	getSensTemp, 214
icarustype, 142	getTiming, 214
interlacing, 142	init_board_specific, 216
lastframe, 142	logcrit, 227
lastrow, 142	logdebug, 227
loglabel, 142	logerr, 227
maxframe, 143	loginfo, 228

logwarn, 228	logcritbase, 120
parseReadoff, 216	logdebug, 120
reportStatusSensor, 217	logdebugbase, 120
selectOscillator, 217	logerr, 120
sensfam, 228	logerrbase, 121
setArbTiming, 218	loginfo, 121
setExtClk, 220	loginfobase, 121
setHighFullWell, 221	logtag, 121
setInterlacing, 221	logwarn, 121
setManualShutters, 221	logwarnbase, 121
setManualTiming, 222	padToFull, 121
setPhiDelay, 223	partition, 120
setTiming, 224	sensor, 122
setTriggerDelay, 226	sensorname, 122
setZeroDeadTime, 226	nsCamera.utils.Packet, 36
nsCamera.utils, 16	nsCamera.utils.Packet.Packet, 192
all, 17	init, 193
nsCamera.utils.crc16pure, 17	addr, 196
_crc16, 18	calculateCRC, 193
CRC16_XMODEM_TABLE, 18	checkCRC, 194
crc16xmodem, 18	checkReadPacket, 194
nsCamera.utils.FlatField, 19	checkResponsePacket, 194
action, 22	checkResponseString, 195
args, 22	cmd, 196
default, 22	crc, 197
dest, 22	data, 197
directory, 23	payload, 197
framelist, 23	payload_length, 197
generateFF, 20	pktStr, 196
getFilenames, 20	preamble, 197
getROlvector, 20	PY3, 197
help, 23	seqID, 197
nargs, 23	type, 198
parser, 23	nsCamera.utils.Subregister, 37
removeFF, 21	nsCamera.utils.Subregister.SubRegister, 228
removeFFall, 21	init, 229
tslopes, 22	addr, 229
nsCamera.utils.misc, 23	max, 229
bytes2str, 24	max value, 230
checkCRC, 25	maxV, 230
flattenlist, 25	min, 230
generateFrames, 26	minV, 230
getEnter, 27	name, 230
loadDumpedData, 27	register, 230
makeLogLabels, 30	resolution, 230
partition, 31	start_bit, 231
plotFrames, 32	value, 231
saveTiffs, 33	width, 231
str2bytes, 35	writable, 231
str2nparray, 35	aldtima
nsCamera.utils.misc.fakeCA, 118	oldtime
init, 119	nsCamera.CameraAssembler.CameraAssembler, 93
boardname, 120	openDevice
logcrit, 120	nsCamera.comms.GigE.GigE, 127

packageroot	ZestETM1_EraseFlashSector, 336
nsCamera.CameraAssembler.CameraAssembler, 93	ZESTETM1_ERROR, 333
padToFull	ZESTETM1_ERROR_CONN, 333
nsCamera.CameraAssembler.CameraAssembler, 93	ZESTETM1_ERROR_GENERAL, 334
nsCamera.sensors.daedalus.daedalus, 116	ZESTETM1_ERROR_STRING, 334
nsCamera.sensors.icarus.icarus, 143	ZestETM1_ErrorHandler, 340
nsCamera.sensors.icarus2.icarus2, 150	ZestETM1_ErrorStrings, 340
nsCamera.utils.misc.fakeCA, 121	ZestETM1_OpenConnection, 336
par	ZESTETM1_RATE_10MHz, 334
nsCamera.comms.RS422.RS422, 209	ZESTETM1_RATE_20MHz, 335
parsedtime	ZESTETM1_RATE_40MHz, 335
nsCamera.CameraAssembler.CameraAssembler, 93	ZESTETM1_REVERSE, 335
parser	ZestETM1_SendCommand, 338
nsCamera.utils.FlatField, 23	ZestETM1_SPIReadWrite, 338
parseReadoff	ZESTETM1_USER_DEVICE_ID, 335
nsCamera.CameraAssembler.CameraAssembler, 65	ZestETM1_WriteFlash, 339
nsCamera.sensors.daedalus.daedalus, 102	PY3
nsCamera.sensors.sensorBase.sensorBase, 216	nsCamera.CameraAssembler.CameraAssembler, 94
partition	nsCamera.comms.RS422.RS422, 210
nsCamera.CameraAssembler.CameraAssembler, 65	nsCamera.utils.Packet.Packet, 197
nsCamera.utils.misc, 31	pyth1
nsCamera.utils.misc.fakeCA, 120	nsCamera.CameraAssembler.CameraAssembler, 94
payload	pyth2
nsCamera.utils.Packet.Packet, 197	nsCamera.CameraAssembler.CameraAssembler, 94
payload_length	python
nsCamera.utils.Packet.Packet, 197	nsCamera.CameraAssembler.CameraAssembler, 94
payloaderror	
	wa a al
nsCamera.CameraAssembler.CameraAssembler, 93	read
	nsCamera.CameraAssembler.CameraAssembler, 94
payloadsize	
payloadsize nsCamera.comms.GigE.GigE, 135	nsCamera.CameraAssembler.CameraAssembler, 94
payloadsize	$ns Camera. Camera Assembler. Camera Assembler, \ 94\\ read\_time out$
payloadsize nsCamera.comms.GigE.GigE, 135 nsCamera.comms.RS422.RS422, 210	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67
payloadsize nsCamera.comms.GigE.GigE, 135 nsCamera.comms.RS422.RS422, 210 pktStr	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210 pktStr     nsCamera.utils.Packet.Packet, 196	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210 pktStr     nsCamera.utils.Packet.Packet, 196 platform	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94  plotFrames     nsCamera.CameraAssembler.CameraAssembler, 65	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202 readImgs
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94  plotFrames	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202 readImgs     nsCamera.CameraAssembler.CameraAssembler, 67
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94  plotFrames     nsCamera.CameraAssembler.CameraAssembler, 65     nsCamera.utils.misc, 32	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202 readImgs     nsCamera.CameraAssembler.CameraAssembler, 67 readoff
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94  plotFrames     nsCamera.CameraAssembler.CameraAssembler, 65     nsCamera.utils.misc, 32  port     nsCamera.CameraAssembler.CameraAssembler, 94	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202 readImgs     nsCamera.CameraAssembler.CameraAssembler, 67 readoff     nsCamera.CameraAssembler.CameraAssembler, 67
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94  plotFrames     nsCamera.CameraAssembler.CameraAssembler, 65     nsCamera.utils.misc, 32  port     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.comms.RS422.RS422, 210	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202 readImgs     nsCamera.CameraAssembler.CameraAssembler, 67 readoff     nsCamera.CameraAssembler.CameraAssembler, 67 nsCamera.CameraAssembler.CameraAssembler, 67 nsCamera.comms.GigE.GigE, 129
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94  plotFrames     nsCamera.CameraAssembler.CameraAssembler, 65     nsCamera.utils.misc, 32  port     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.comms.RS422.RS422, 210  powerCheck	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202 readImgs     nsCamera.CameraAssembler.CameraAssembler, 67 readoff     nsCamera.CameraAssembler.CameraAssembler, 67 nsCamera.CameraAssembler.CameraAssembler, 67 nsCamera.comms.GigE.GigE, 129 nsCamera.comms.RS422.RS422, 203
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94  plotFrames     nsCamera.CameraAssembler.CameraAssembler, 65     nsCamera.utils.misc, 32  port     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.comms.RS422.RS422, 210  powerCheck     nsCamera.CameraAssembler.CameraAssembler, 65	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202 readImgs     nsCamera.CameraAssembler.CameraAssembler, 67 readoff     nsCamera.CameraAssembler.CameraAssembler, 67 nsCamera.CameraAssembler.CameraAssembler, 67 nsCamera.comms.GigE.GigE, 129 nsCamera.comms.RS422.RS422, 203 readSerial
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94  plotFrames     nsCamera.CameraAssembler.CameraAssembler, 65     nsCamera.utils.misc, 32  port     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.comms.RS422.RS422, 210  powerCheck     nsCamera.CameraAssembler.CameraAssembler, 65  preamble	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202 readImgs     nsCamera.CameraAssembler.CameraAssembler, 67 readoff     nsCamera.CameraAssembler.CameraAssembler, 67 nsCamera.comms.GigE.GigE, 129     nsCamera.comms.RS422.RS422, 203 readSerial     nsCamera.CameraAssembler.CameraAssembler, 68
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94  plotFrames     nsCamera.CameraAssembler.CameraAssembler, 65     nsCamera.utils.misc, 32  port     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.comms.RS422.RS422, 210  powerCheck     nsCamera.CameraAssembler.CameraAssembler, 65  preamble     nsCamera.utils.Packet.Packet, 197	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202 readImgs     nsCamera.CameraAssembler.CameraAssembler, 67 readoff     nsCamera.CameraAssembler.CameraAssembler, 67 nsCamera.comms.GigE.GigE, 129     nsCamera.comms.RS422.RS422, 203 readSerial     nsCamera.CameraAssembler.CameraAssembler, 68 nsCamera.CameraAssembler.CameraAssembler, 68 nsCamera.comms.GigE.GigE, 130
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94  plotFrames     nsCamera.CameraAssembler.CameraAssembler, 65     nsCamera.utils.misc, 32  port     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.comms.RS422.RS422, 210  powerCheck     nsCamera.CameraAssembler.CameraAssembler, 65  preamble     nsCamera.utils.Packet.Packet, 197  printBoardInfo	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202 readImgs     nsCamera.CameraAssembler.CameraAssembler, 67 readoff     nsCamera.CameraAssembler.CameraAssembler, 67 nsCamera.comms.GigE.GigE, 129 nsCamera.comms.RS422.RS422, 203 readSerial     nsCamera.CameraAssembler.CameraAssembler, 68 nsCamera.comms.GigE.GigE, 130 nsCamera.comms.RS422.RS422, 204
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94  plotFrames     nsCamera.CameraAssembler.CameraAssembler, 65     nsCamera.utils.misc, 32  port     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.comms.RS422.RS422, 210  powerCheck     nsCamera.CameraAssembler.CameraAssembler, 65  preamble     nsCamera.utils.Packet.Packet, 197	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202 readImgs     nsCamera.CameraAssembler.CameraAssembler, 67 readoff     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 129     nsCamera.comms.RS422.RS422, 203 readSerial     nsCamera.CameraAssembler.CameraAssembler, 68     nsCamera.CameraAssembler.CameraAssembler, 68     nsCamera.Comms.GigE.GigE, 130     nsCamera.comms.RS422.RS422, 204 readSRAM
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94  plotFrames     nsCamera.CameraAssembler.CameraAssembler, 65     nsCamera.utils.misc, 32  port     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.comms.RS422.RS422, 210  powerCheck     nsCamera.CameraAssembler.CameraAssembler, 65  preamble     nsCamera.utils.Packet.Packet, 197  printBoardInfo     nsCamera.CameraAssembler.CameraAssembler, 66  Private.h	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202 readImgs     nsCamera.CameraAssembler.CameraAssembler, 67 readoff     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 129     nsCamera.comms.RS422.RS422, 203 readSerial     nsCamera.CameraAssembler.CameraAssembler, 68     nsCamera.CameraAssembler, 68     nsCamera.CameraAsse
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94  plotFrames     nsCamera.CameraAssembler.CameraAssembler, 65     nsCamera.utils.misc, 32  port     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.comms.RS422.RS422, 210  powerCheck     nsCamera.CameraAssembler.CameraAssembler, 65  preamble     nsCamera.utils.Packet.Packet, 197  printBoardInfo     nsCamera.CameraAssembler.CameraAssembler, 66  Private.h     closesocket, 333	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202 readImgs     nsCamera.CameraAssembler.CameraAssembler, 67 readoff     nsCamera.CameraAssembler.CameraAssembler, 67 nsCamera.comms.GigE.GigE, 129     nsCamera.comms.RS422.RS422, 203 readSerial     nsCamera.CameraAssembler.CameraAssembler, 68     nsCamera.comms.GigE.GigE, 130     nsCamera.comms.RS422.RS422, 204 readSRAM     nsCamera.boards.LLNL_v1.llnl_v1, 163     nsCamera.boards.LLNL_v4.llnl_v4, 183
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94  plotFrames     nsCamera.CameraAssembler.CameraAssembler, 65     nsCamera.utils.misc, 32  port     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.comms.RS422.RS422, 210  powerCheck     nsCamera.CameraAssembler.CameraAssembler, 65  preamble     nsCamera.utils.Packet.Packet, 197  printBoardInfo     nsCamera.CameraAssembler.CameraAssembler, 66  Private.h     closesocket, 333     SD_BOTH, 333	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202 readImgs     nsCamera.CameraAssembler.CameraAssembler, 67 readoff     nsCamera.CameraAssembler.CameraAssembler, 67 nsCamera.comms.GigE.GigE, 129     nsCamera.comms.RS422.RS422, 203 readSerial     nsCamera.CameraAssembler.CameraAssembler, 68     nsCamera.comms.GigE.GigE, 130     nsCamera.comms.RS422.RS422, 204 readSRAM     nsCamera.boards.LLNL_v1.llnl_v1, 163     nsCamera.boards.LLNL_v4.llnl_v4, 183     nsCamera.CameraAssembler.CameraAssembler, 68
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94  plotFrames     nsCamera.CameraAssembler.CameraAssembler, 65     nsCamera.utils.misc, 32  port     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.comms.RS422.RS422, 210  powerCheck     nsCamera.CameraAssembler.CameraAssembler, 65  preamble     nsCamera.utils.Packet.Packet, 197  printBoardInfo     nsCamera.CameraAssembler.CameraAssembler, 66  Private.h     closesocket, 333     SD_BOTH, 333     SOCKET, 335	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202 readImgs     nsCamera.CameraAssembler.CameraAssembler, 67 readoff     nsCamera.CameraAssembler.CameraAssembler, 67 nsCamera.comms.GigE.GigE, 129     nsCamera.comms.RS422.RS422, 203 readSerial     nsCamera.CameraAssembler.CameraAssembler, 68     nsCamera.comms.GigE.GigE, 130     nsCamera.comms.RS422.RS422, 204 readSRAM     nsCamera.boards.LLNL_v1.llnl_v1, 163     nsCamera.boards.LLNL_v4.llnl_v4, 183     nsCamera.CameraAssembler.CameraAssembler, 68 readTimeout
payloadsize     nsCamera.comms.GigE.GigE, 135     nsCamera.comms.RS422.RS422, 210  pktStr     nsCamera.utils.Packet.Packet, 196  platform     nsCamera.CameraAssembler.CameraAssembler, 94  plotFrames     nsCamera.CameraAssembler.CameraAssembler, 65     nsCamera.utils.misc, 32  port     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.CameraAssembler.CameraAssembler, 94     nsCamera.comms.RS422.RS422, 210  powerCheck     nsCamera.CameraAssembler.CameraAssembler, 65  preamble     nsCamera.utils.Packet.Packet, 197  printBoardInfo     nsCamera.CameraAssembler.CameraAssembler, 66  Private.h     closesocket, 333     SD_BOTH, 333	nsCamera.CameraAssembler.CameraAssembler, 94 read_timeout     nsCamera.comms.RS422.RS422, 210 readFrames     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 128     nsCamera.comms.RS422.RS422, 202 readImgs     nsCamera.CameraAssembler.CameraAssembler, 67 readoff     nsCamera.CameraAssembler.CameraAssembler, 67     nsCamera.comms.GigE.GigE, 129     nsCamera.comms.RS422.RS422, 203 readSerial     nsCamera.CameraAssembler.CameraAssembler, 68     nsCamera.comms.GigE.GigE, 130     nsCamera.comms.RS422.RS422, 204 readSRAM     nsCamera.boards.LLNL_v1.llnl_v1, 163     nsCamera.boards.LLNL_v4.llnl_v4, 183     nsCamera.CameraAssembler.CameraAssembler, 68 readTimeout     nsCamera.CameraAssembler.CameraAssembler, 68

nsCamera.utils.Subregister.SubRegister, 230 registers	nsCamera.CameraAssembler.CameraAssembler, 74 nsCamera.comms.GigE.GigE, 131
nsCamera.boards.LLNL_v1.llnl_v1, 170	nsCamera.comms.RS422.RS422, 205
nsCamera.boards.LLNL_v4.llnl_v4, 190	sens_registers
reinitialize	nsCamera.sensors.daedalus.daedalus, 116
nsCamera.CameraAssembler.CameraAssembler, 68	nsCamera.sensors.icarus.icarus, 143
removeFF	
	nsCamera.sensors.icarus2.icarus2, 150
nsCamera.utils.FlatField, 21 removeFFall	sens_subregisters
	nsCamera.sensors.daedalus.daedalus, 117 nsCamera.sensors.icarus.icarus, 144
nsCamera.utils.FlatField, 21	
reportEdgeDetects	nsCamera.sensors.icarus2.icarus2, 151
nsCamera.boards.LLNL_v1.llnl_v1, 163	sensfam
nsCamera.boards.LLNL_v4.llnl_v4, 183	nsCamera.sensors.daedalus.daedalus, 117
nsCamera.CameraAssembler.CameraAssembler, 69	nsCamera.sensors.icarus.icarus, 144
reportStatus	nsCamera.sensors.icarus2.icarus2, 151
nsCamera.boards.LLNL_v1.llnl_v1, 164	nsCamera.sensors.sensorBase.sensorBase, 228
nsCamera.boards.LLNL_v4.llnl_v4, 184	sensmanual
nsCamera.CameraAssembler.CameraAssembler, 69	nsCamera.CameraAssembler.CameraAssembler, 95
reportStatusSensor	sensor
nsCamera.sensors.daedalus.daedalus, 103	nsCamera.CameraAssembler.CameraAssembler, 95
nsCamera.sensors.sensorBase.sensorBase, 217	nsCamera.utils.misc.fakeCA, 122
resetTimer	sensorname
nsCamera.boards.LLNL_v1.llnl_v1, 164	nsCamera.CameraAssembler.CameraAssembler, 95
nsCamera.boards.LLNL_v4.llnl_v4, 185	nsCamera.utils.misc.fakeCA, 122
nsCamera.CameraAssembler.CameraAssembler, 69	sensorSpecific
resolution	nsCamera.CameraAssembler.CameraAssembler, 74
nsCamera.utils.Subregister.SubRegister, 230	nsCamera.sensors.daedalus.daedalus, 104
resolveSubreg	nsCamera.sensors.icarus.icarus, 140
nsCamera.CameraAssembler.CameraAssembler, 69	nsCamera.sensors.icarus2.icarus2, 147
rs422_baud	senstiming
nsCamera.boards.LLNL_v1.llnl_v1, 170	nsCamera.CameraAssembler.CameraAssembler, 95
nsCamera.boards.LLNL_v4.llnl_v4, 191	seqID
rs422_cmd_wait	nsCamera.utils.Packet.Packet, 197
nsCamera.boards.LLNL_v1.llnl_v1, 170	serialClose
nsCamera.boards.LLNL_v4.llnl_v4, 191	nsCamera.comms.RS422.RS422, 207
	setArbTiming
saveFrames	nsCamera.CameraAssembler.CameraAssembler, 74
nsCamera.CameraAssembler.CameraAssembler, 70	nsCamera.sensors.sensorBase.sensorBase, 218
saveHDF	setExtClk
nsCamera.CameraAssembler.CameraAssembler, 71	nsCamera.CameraAssembler.CameraAssembler, 75
saveNumpys	nsCamera.sensors.daedalus.daedalus, 106
nsCamera.CameraAssembler.CameraAssembler, 72	nsCamera.sensors.sensorBase.sensorBase, 220
saveTiffs	setFrames
nsCamera.CameraAssembler.CameraAssembler, 74	nsCamera.CameraAssembler.CameraAssembler, 75
nsCamera.utils.misc, 33	setHighFullWell
savetime	nsCamera.CameraAssembler.CameraAssembler, 76
nsCamera.CameraAssembler.CameraAssembler, 95	nsCamera.sensors.daedalus.daedalus, 107
SD_BOTH	nsCamera.sensors.sensorBase.sensorBase, 221
Private.h, 333	setInterlacing
selectOscillator	nsCamera.CameraAssembler.CameraAssembler, 76
nsCamera.CameraAssembler.CameraAssembler, 74	nsCamera.sensors.daedalus.daedalus, 107
nsCamera.sensors.daedalus.daedalus, 104	nsCamera.sensors.sensorBase.sensorBase, 221
nsCamera.sensors.sensorBase.sensorBase, 217	setLED
sendCMD	nsCamera.boards.LLNL_v1.llnl_v1, 165
	:

nsCamera.boards.LLNL_v4.llnl_v4, 185	nsCamera.sensors.icarus2.icarus2, 151
nsCamera.CameraAssembler.CameraAssembler, 77	start_bit
setManualShutters	nsCamera.utils.Subregister.SubRegister, 231
nsCamera.CameraAssembler.CameraAssembler, 77	startCapture
nsCamera.sensors.daedalus.daedalus, 109	nsCamera.boards.LLNL_v1.llnl_v1, 167
nsCamera.sensors.sensorBase.sensorBase, 221	nsCamera.boards.LLNL_v4.llnl_v4, 187 nsCamera.CameraAssembler.CameraAssembler, 87
setManualTiming nsCamera.CameraAssembler.CameraAssembler, 77	
nsCamera.sensors.daedalus.daedalus, 109	nsCamera.comms.RS422.RS422, 210
nsCamera.sensors.sensorBase.sensorBase, 222	str2bytes
setPhiDelay	nsCamera.CameraAssembler.CameraAssembler, 87
nsCamera.CameraAssembler.CameraAssembler, 77	nsCamera.utils.misc, 35
nsCamera.sensors.daedalus.daedalus, 110	str2nparray
nsCamera.sensors.sensorBase.sensorBase, 223	nsCamera.CameraAssembler.CameraAssembler, 87
setPot	nsCamera.utils.misc, 35
nsCamera.CameraAssembler.CameraAssembler, 77	submitMessages
setPotV	nsCamera.CameraAssembler.CameraAssembler, 87
nsCamera.CameraAssembler.CameraAssembler, 79	subregisters
setPowerSave	nsCamera.boards.LLNL_v1.llnl_v1, 171
nsCamera.boards.LLNL_v1.llnl_v1, 165	nsCamera.boards.LLNL_v4.llnl_v4, 191
nsCamera.boards.LLNL_v4.llnl_v4, 185	subreglist
nsCamera.CameraAssembler.CameraAssembler, 82	nsCamera.boards.LLNL_v1.llnl_v1, 171
setPPER	nsCamera.boards.LLNL_v4.llnl_v4, 191
nsCamera.boards.LLNL_v1.llnl_v1, 166	Atomic and
nsCamera.boards.LLNL_v4.llnl_v4, 186	timeout
nsCamera.CameraAssembler.CameraAssembler, 82	nsCamera.CameraAssembler.CameraAssembler, 95 toffset
setRegister	nsCamera.sensors.daedalus.daedalus, 117
nsCamera.CameraAssembler.CameraAssembler, 82	trigtime
setRows	nsCamera.CameraAssembler.CameraAssembler, 95
nsCamera.CameraAssembler.CameraAssembler, 83 setSubregister	tslope
nsCamera.CameraAssembler.CameraAssembler, 85	nsCamera.sensors.daedalus.daedalus, 117
setTiming	tslopes
nsCamera.CameraAssembler.CameraAssembler, 86	nsCamera.utils.FlatField, 22
nsCamera.sensors.sensorBase.sensorBase, 224	type
setTriggerDelay	nsCamera.utils.Packet.Packet, 198
nsCamera.CameraAssembler.CameraAssembler, 86	
nsCamera.sensors.daedalus.daedalus, 111	ubyte4
nsCamera.sensors.sensorBase.sensorBase, 226	nsCamera.comms.GigE.GigE.ZESTETM1_CARD_INFO
setZeroDeadTime	232
nsCamera.CameraAssembler.CameraAssembler, 87	ubyte6
nsCamera.sensors.daedalus.daedalus, 112	nsCamera.comms.GigE.GigE.ZESTETM1_CARD_INFO.
nsCamera.sensors.sensorBase.sensorBase, 226	unstringed
skipError	nsCamera.CameraAssembler.CameraAssembler, 96
nsCamera.comms.GigE.GigE, 135	UPnP.c
nsCamera.comms.RS422.RS422, 210	_CRT_SECURE_NO_WARNINGS, 344
SOCKET	strnicmp, 344
Private.h, 335	ZESTETM1_COMMAND_GET_SETTINGS, 344
softReboot	ZestETM1_GetAllAdapters, 345
nsCamera.boards.LLNL_v1.llnl_v1, 166	ZestETM1_GetCardInfo, 346
nsCamera.boards.LLNL_v4.llnl_v4, 186	ZestETM1_HTTPGet, 347
specwarn nsCamera.sensors.daedalus.daedalus, 117	ZestETM1_JoinGroup, 349
nsCamera.sensors.idaedalus.daedalus, 117	ZestETM1_LeaveGroup, 349
1130011110113.10011103.1001103, 144	ZestETM1_ReadSettings, 350

ZestETM1_SearchReq, 353	ZESTETM1_ILLEGAL_CONNECTION, 366
ZESTETM1_UPNP_ADDR, 344	ZESTETM1_ILLEGAL_PARAMETER, 366
ZESTETM1_UPNP_PORT, 344	ZESTETM1_ILLEGAL_STATUS_CODE, 366
ZestETM1CountCards, 350	ZESTETM1_INFO_BASE, 364
ZestETM1FreeCards, 352	ZESTETM1_INTERNAL_ERROR, 366
ZestETM1GetCardInfo, 353	ZESTETM1_INVALID_CONNECTION_TYPE, 366
· · · · · ·	ZESTETM1_MAX_ERROR, 366
value	ZESTETM1_MAX_INFO, 366
nsCamera.utils.Subregister.SubRegister, 231	ZESTETM1_MAX_WARNING, 366
verblevel	ZESTETM1_NULL_PARAMETER, 366
nsCamera.CameraAssembler.CameraAssembler, 96	ZESTETM1_OUT_OF_MEMORY, 366
verbmap	ZESTETM1_SOCKET_CLOSED, 366
nsCamera.CameraAssembler.CameraAssembler, 96	ZESTETM1_SOCKET_ERROR, 366
verbose	ZESTETM1_SPI_RATE, 365
nsCamera.CameraAssembler.CameraAssembler, 96	ZESTETM1_SPI_RATE_17_5, 366
version	ZESTETM1_SPI_RATE_35, 365
nsCamera.CameraAssembler.CameraAssembler, 96	ZESTETM1_SPI_RATE_8_75, 366
VREF	ZESTETM1_STATUS, 366
nsCamera.boards.LLNL_v1.llnl_v1, 171	ZESTETM1_SUCCESS, 366
nsCamera.boards.LLNL_v4.llnl_v4, 191	ZESTETM1_TIMEOUT, 366
waited	ZESTETM1_TYPE_TCP, 365
nsCamera.CameraAssembler.CameraAssembler, 96	ZESTETM1_TYPE_UDP, 365
waitForSRAM	ZESTETM1_VERSION_FALLBACK, 364
nsCamera.boards.LLNL_v1.llnl_v1, 168	ZESTETM1_WARNING_BASE, 364
nsCamera.boards.LLNL_v4.llnl_v4, 188	ZestETM1Close, 367
nsCamera.CameraAssembler.CameraAssembler, 88	ZestETM1CloseConnection, 367
width	ZestETM1CountCards, 367
	ZestETM1FreeCards, 369
nsCamera.sensors.daedalus.daedalus, 117 nsCamera.sensors.icarus.icarus, 144	ZestETM1GetCardInfo, 369
nsCamera.sensors.icarus2.icarus2, 151	ZestETM1GetErrorMessage, 370
nsCamera.utils.Subregister.SubRegister, 231	ZestETM1Init, 370
writable	ZestETM1OpenConnection, 370
nsCamera.utils.Subregister.SubRegister, 231	ZestETM1ReadData, 371
write_timeout	ZestETM1ReadRegister, 371
nsCamera.comms.RS422.RS422, 210	ZestETM1RegisterErrorHandler, 372
writeSerial	ZestETM1SetInterrupt, 372
nsCamera.CameraAssembler.CameraAssembler, 88	ZestETM1SPIReadWrite, 373
nsCamera.comms.GigE.GigE, 132	ZestETM1WriteData, 373
nsCamera.comms.RS422.RS422, 207	ZestETM1WriteRegister, 374
writeTimeout	ZESTETM1_CARD_INFO, 364
nsCamera.comms.GigE.GigE, 135	nsCamera.comms.GigE.GigE, 136
nodinicia.commo.dige.dige, 100	ZestETM1_CloseConnection
ZCountCards	Data.c, 305
nsCamera.comms.GigE.GigE, 135	Private.h, 336
ZDT	ZESTETM1_COMMAND_GET_SETTINGS
nsCamera.sensors.daedalus.daedalus, 117, 118	UPnP.c, 344
ZErrorDict	ZESTETM1_COMMAND_MAILBOX_INT
nsCamera.comms.GigE.GigE, 136	Data.c, 304
ZestETM1.h	ZESTETM1_COMMAND_READ_REG
ZESTETM1_CONNECTION, 365	Data.c, 304
ZESTETM1_CONNECTION_TYPE, 365	ZESTETM1_COMMAND_SPI
ZESTETM1_ERROR_BASE, 364	Data.c, 304
ZESTETM1_ERROR_FUNC, 365	ZESTETM1_COMMAND_WRITE_REG
ZESTETM1_HANDLE, 365	Data.c, 304

ZESTETM1_CONNECTION	UPnP.c, 349
ZestETM1.h, 365	ZESTETM1_MAILBOX_INT_CMD, 303
ZESTETM1_CONNECTION_HANDLE_MAGIC	ZESTETM1_MAILBOX_INT_RESPONSE, 303
Private.h, 333	ZESTETM1_MAX_ERROR
ZESTETM1_CONNECTION_STRUCT, 332	ZestETM1.h, 366
ZESTETM1_CONNECTION_TYPE	ZESTETM1_MAX_INFO
ZestETM1.h, 365	ZestETM1.h, 366
ZESTETM1_DEFAULT_TIMEOUT	ZESTETM1 MAX WARNING
Private.h, 333	ZestETM1.h, 366
ZestETM1 EraseFlashSector	ZESTETM1 NULL PARAMETER
Private.h, 336	ZestETM1.h, 366
ZESTETM1_ERROR	ZestETM1_OpenConnection
Private.h, 333	Data.c, 305
ZESTETM1 ERROR BASE	Private.h, 336
ZestETM1.h, 364	ZESTETM1_OUT_OF_MEMORY
ZESTETM1_ERROR_CONN	ZestETM1.h, 366
Private.h, 333	ZESTETM1 RATE 10MHz
ZESTETM1 ERROR FUNC	Private.h, 334
ZestETM1.h, 365	ZESTETM1_RATE_20MHz
ZESTETM1 ERROR GENERAL	
	Private.h, 335
Private.h, 334	ZESTETM1_RATE_40MHz
ZESTETM1_ERROR_STRING	Private.h, 335
Private.h, 334	ZESTETM1_READ_REG_CMD, 303
ZestETM1_ErrorHandler	ZESTETM1_READ_REG_RESPONSE, 303
Error.c, 327	ZestETM1_ReadData
Private.h, 340	Data.c, 307
ZestETM1_ErrorStrings	ZestETM1_ReadSettings
Error.c, 327	UPnP.c, 350
Private.h, 340	ZESTETM1_REVERSE
ZESTETM1_GET_SETTINGS_CMD, 343	Private.h, 335
ZESTETM1_GET_SETTINGS_RESPONSE, 343	ZestETM1_SearchReq
ZestETM1_GetAllAdapters	UPnP.c, 353
UPnP.c, 345	ZestETM1_SendCommand
ZestETM1_GetCardInfo	Data.c, 308
UPnP.c, 346	Private.h, 338
ZESTETM1_HANDLE	ZESTETM1_SOCKET_CLOSED
ZestETM1.h, 365	ZestETM1.h, 366
ZestETM1_HTTPGet	ZESTETM1_SOCKET_ERROR
UPnP.c, 347	ZestETM1.h, 366
ZESTETM1_ILLEGAL_CONNECTION	ZESTETM1_SPI_RATE
ZestETM1.h, 366	ZestETM1.h, 365
ZESTETM1_ILLEGAL_PARAMETER	ZESTETM1_SPI_RATE_17_5
ZestETM1.h, 366	ZestETM1.h, 366
ZESTETM1_ILLEGAL_STATUS_CODE	ZESTETM1_SPI_RATE_35
ZestETM1.h, 366	ZestETM1.h, 365
ZESTETM1_INFO_BASE	ZESTETM1_SPI_RATE_8_75
ZestETM1.h, 364	ZestETM1.h, 366
ZESTETM1_INTERNAL_ERROR	ZestETM1_SPIReadWrite
ZestETM1.h, 366	
ZESTETM1_INVALID_CONNECTION_TYPE	Private.h, 338
ZestETM1.h, 366	ZESTETM1_STATUS
ZestETM1_JoinGroup	ZestETM1.h, 366
UPnP.c, 349	ZESTETM1_SUCCESS
ZestETM1_LeaveGroup	ZestETM1.h, 366
_ ,	,

ZESTETM1_TIMEOUT	ZestETM1.h, 372
ZestETM1.h, 366	ZestETM1SetInterrupt
ZESTETM1_TYPE_TCP	Data.c, 313
ZestETM1.h, 365	ZestETM1.h, 372
ZESTETM1 TYPE UDP	ZestETM1SPIReadWrite
ZestETM1.h, 365	Data.c, 313
	•
ZESTETM1_UPNP_ADDR	ZestETM1.h, 373
UPnP.c, 344	ZestETM1WriteData
ZESTETM1_UPNP_PORT	Data.c, 314
UPnP.c, 344	ZestETM1.h, 373
ZESTETM1_USER_DEVICE_ID	ZestETM1WriteRegister
Private.h, 335	Data.c, 314
ZESTETM1_VERSION_FALLBACK	ZestETM1.h, 374
ZestETM1.h, 364	ZOpenConnection
ZESTETM1_WARNING_BASE	nsCamera.comms.GigE.GigE
ZestETM1.h, 364	ZReadData
ZESTETM1_WRITE_REG_CMD, 302	nsCamera.comms.GigE.GigE
ZESTETM1_WRITE_REG_RESPONSE, 302	ZWriteData
ZestETM1_WriteData	nsCamera.comms.GigE.GigE
Data.c, 309	
ZestETM1_WriteFlash	
Private.h, 339	
ZestETM1Close	
Main.c, 329	
ZestETM1.h, 367	
ZestETM1CloseConnection	
Data.c, 311	
ZestETM1.h, 367	
ZestETM1CountCards	
UPnP.c, 350	
ZestETM1.h, 367	
ZestETM1FreeCards	
UPnP.c, 352	
ZestETM1.h, 369	
ZestETM1GetCardInfo	
UPnP.c, 353	
ZestETM1.h, 369	
ZestETM1GetErrorMessage	
Error.c, 326	
ZestETM1.h, 370	
ZestETM1Init	
Main.c, 329	
ZestETM1.h, 370	
ZestETM1OpenConnection	
Data.c, 311	
ZestETM1.h, 370	
ZestETM1ReadData	
Data.c, 311	
ZestETM1.h, 371	
ZestETM1ReadRegister	
Data.c, 312	
ZestETM1.h, 371	
ZestETM1RegisterErrorHandler	
Error.c, 326	

136

136

136